



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com



Product Description

Item # 101382

You are bidding on one tooth polishing air prophylaxis unit. This mini, lightweight, air-powered tooth polisher can be used to remove plaque from a variety of surfaces. Using a combination of air pressure and polishing powder, our air prophylaxis system gives teeth back their natural whiteness. Whether polishing natural tooth enamel, porcelain, metal, or resin, our air prophylaxis unit effectively gets rid of stains from coffee, tea, wine, food, smoking, and etc., leaving you with a bright and beautiful smile!

Specifications & Features

- Compact ergonomic design
- Our 360 degree rotation effectively cleans even hidden tartar
- Especially designed to minimize splash-back
- Highly efficient jet of sand powder and pressurized air insures even polishing while reducing kinetic energy loss
- Easy operation, clean up, and maintenance
- Withstands high pressure and heat
- Working head of the polisher can be easily sterilized by submerging it in boiling water
- Water pressure: 0.05MPa - 0.5MPa
- Air pressure: 0.2MPa-0.5MPa
- Max. Water-flow: 30ml/min
- Max. Powder-flow: 2g-3g/min
- 4-Hole standard connection
- Weight: 1lbs

Packing List

- One Air Polisher
- Two Spare cover
- Two clean needles
- One operational manual in English