



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com



Apex Loccator Joypex 5

Item # 101954

Apex Locator Joypex 5 adopts the fifth generation circuit-T Surge Frequency & Low Frequency in the area of apex locator.

Distinctive function:

* **Battery operated – 5 x 1.5v AAA**

*Frequency operating system-0.1mm length determination

*Works in dry or wet canals

* **Adjustable LCD angle of view**

*Self calibrating

* **LCD panel could show the moving trace of file in the canal**

*Autoclavable file holder to avoid cross infection

*Audible warning system

*Adjustable volume of the sound indicator

*Excellence in working length determination comes with 5 different setting of apical adjustment around the minor diameter

*CE 0197



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Areentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

Standard accessory:

- *main unit 1 pc, testing cord 1 pc
- *file holder (autoclaveable) 2 pcs
- *contrary electrode 4 pcs, user's manual 1pc
- *root canal ruler (optional)

Standard



Optional

