



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com



Safendo Endo Motor

Item #102223

Safendo Endo Motor is a newly made root canal treatment device. It adopts brushless electric motor to more accurately control the speed and torque. It is more convenient for a dentist to operate.

Main Features:

- * Brushless electric motor
1:1 contra angle
Speed Range 1500-2000 rmp
Adjustable torque range 0.6-5.6 N.cm
- * **Apex locator included**
- * Adjustable LCD to make operation easier
- * In the process of root canal treatment, when the file reaches apex, the unit will stop working automatically
- * Saving Function: Then groups of program for pseed and torque can be saved
- * Auto reverse function: When torque exceeds the pre-set value, the file will reverse automatically to avoid the perforated and break-off