



# US DENTAL DEPOT

## US DENTAL DEPOT

6555 NW 9<sup>th</sup> Avenue suite 212  
Fort Lauderdale 33309  
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

[www.USDENTALDEPOT.com](http://www.USDENTALDEPOT.com)



### ECCO Digital Intra Oral Camera 3

Item # 102234

Simplest intra-oral camera ever! The camera weighs under a pound and fits into the palm of almost any hand. The camera is capable of accomplishing the full range of dental photography with ease from the 1:1 extreme close up of a posterior quadrant to anterior, lateral and all mirror photos including occlusals and linguals. This camera works great for portraits and personal photography as well.

#### Full Package Includes:

- Camera 8mp
- Direct connect downloading cord
- Imaging software
- Hard Carrying Case
- Cheek Refractors



# US DENTAL DEPOT

## US DENTAL DEPOT

6555 NW 9<sup>th</sup> Avenue suite 212  
Fort Lauderdale 33309  
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

[www.USDENTALDEPOT.com](http://www.USDENTALDEPOT.com)

- Intra oral mirror
- Lithium rechargeable Battery and Charger
- 2 gig Memory card [SD]
- Memory card reader
- DVD tutorial and written instructions

ECCO Digital Intraoral camera software program designed for image management, cosmetic imaging and printing. The program is capable of creating patient files with images. Users can record and cross reference data, run slide shows, and lay the images for printing.

Printing templates included: for before and afters, the typical orthodontic layout, letter writing, plus much more.

In addition the program includes: a full module for cosmetic imaging, allowing for such common imaging tasks as diastema closure, whitening, lengthening, smile gallery simulation and composite filling simulation.

The capability of the imaging software is virtually endless and is designed to rival the high-end dental imaging software found on the market.