



US DENTAL DEPOT

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com



DIGITAL RADIOGRAPHY SENSOR SYSTEM

THE DIGIREX Digital Dental Radiography System Features CCD Technology, providing high definition real-time Imaging for in-depth examination and precise diagnosis. Digital x-ray imaging is much safer compared to conventional x-ray film due to the dramatically reduced radiation exposure. The digital images can be easily stored, retrieved and transmitted electronically.

Technical Specifications:
Sensor Technology: CCD and Scintillator

Sensor Size 1

Image Area: 30 x 20mm
1,200 lines x 800 pixels

Sensor Size 2

Image Area: 36.5 x 25.9mm
1,872 lines x 1,330 pixels
Storage: -20 to 70 degrees Celsius

Sensor Cable

Length: 2.9 + 0.1m
Diameter: 3.45m + 0.15m

PC Minimum Requirements

Operating System: Windows 98 or better
CPU: 300MHz
RAM: 128MB
USB: 1.1 or 2.2
Video Card: SVGA
Monitor Resolution: 800 x 600