



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com



High Temperature High Pressure Autoclave Sterilizer

There are 3 sizes:

2-01-01 23L

2-01-01 18L

2-01-01 16L

All parts of apparatus are imported

- Italy UKLA water pump
- MOTOROLA chip and pressure sensor
- Imported vacuum pump
- PARKER valve housing

PRODUCT CHARACTERISTIC:

1. It adopts European B standard, with 3-times



**US
DENTAL
DEPOT**

US DENTAL DEPOT

6555 NW 9th Avenue suite 212
Fort Lauderdale 33309
Florida USA

USA - 1 1866 625 5558 – 1 954 874 6325

Colombia 57 1 489 9051

Peru 51 1 707 3444

Venezuela 58 212 729 1639

Argentina 54 11 5258 7148 Cba 54 351 569 1409

www.USDENTALDEPOT.com

pre-vacuum preceding vacuum drying. Residual

humidity is 0.1;vacuum measurement reaches 0.8bar.

It sterilizes miscellaneous instruments,

successfully solved the problem that autoclave

can not sterilize plastic instrument.

2. Computer control, pressure and temperature are displayed in digital.

**3. Quick spraying steam generator makes the inside temperature
of autoclave well-proportioned.**

**4. Waste water tank water-full alarm system avoiding
the possibility of circulated sterilizing with waste water.**

5. Automatic malfunction detection alarm.

6. The boiler is extended by 2.5mm stainless steel one shot.

**7. The scuppers at the bottom of the boiler make drying
more convenient.**