

# DualTrim™ Dry Model Trimmer



## High Performance Rotary Dry Model Trimming Machine That Is Convertible For Wet Trimming Operation



- **Fast and Accurate Trimming**

Dry trimming eliminates the potential ecological and accuracy problems associated with wet trimming

- **Convertible for Dry or Wet Operation**

Optional package allows the conversion to wet trimming operation for maximum versatility

- **New Exhaust Port Design**

1½" OD exhaust port promotes better capture of dry and wet grinding debris resulting in less cleaning

- **New Quick-Stick™ Abrasive Disk System**

Faster trimming and fast, "wrench-less" disk changes

- **Improved Water Spray Bar Design for Wet Trimming**

Moveable spray bar promotes improved cleaning of wheel and rotates out of the way during wheel or disk changes

- **New Noise Dampening Feet**

Quieter operation and easier installation

- **Ordering Information**

<u>Item Number</u>	<u>Description</u>
61850	<b>DualTrim Rotary Dry Model Trimmer, 1/3 HP, 10" Diameter, 120 V AC, 60 Hz</b> (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61855	<b>DualTrim Rotary Dry Model Trimmer, 1/3 HP, 10" Diameter, 220 V AC, 60 Hz</b> (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61860	<b>DualTrim Rotary Dry Model Trimmer, 1/2 HP, 10" Diameter, 120 V AC, 60 Hz</b> (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61865	<b>DualTrim Rotary Dry Model Trimmer, 1/2 HP, 10" Diameter, 220 V AC, 60 Hz</b> (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61875	<b>Water Conversion Kit For DualTrim</b> (Includes Water Spray Bar Assembly, Water Supply and Exhaust Hoses, Exhaust Hose Adapter, and Arbor Wrench Set)
78900	<b>Custom Safety Shield Kit</b> (Includes Clear Plastic Shield, Metal Goose-Neck, and Base that Clamps to Motor Base)
61798	<b>Rubber Feet For Model Trimmers, Set of 4</b>
61965	<b>Quick-Stick Abrasive Disks, Coarse Grit, 10" Diameter (w/ 2" Center Hole), Pkg. of 6</b>

Note: Dry model trimming requires the use of an appropriate dust collection system (not included) for removal of the dry grinding waste.

**Buffalo Dental Manufacturing Co., Inc.,**  
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA  
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751  
[www.buffalodental.com](http://www.buffalodental.com)

S61850 0607

**Please order today from  
your favorite Dealer!**