

# Shade of VITA 3D-Master vs. Vintage Halo Color Conversion Chart

with C,D color

VITA 3D-Master	Vintage Halo (Opaque)	Vintage Halo (Body)	Vintage Halo (Enamel)
0	0M1	W1O:1	W1B:1
	0M2		
	0M3	W2O:1	W2B:1
1	1M1	W3O:1	W3B:1
	1M2	A1O:3	A1B:3
		R1O:3	R1B:3
2	2L1.5	VOP:2	VDL:2
	2L2.5	B1O:1	VB1B:1
	2M1	A2O:1	A2B:1
	2M2	B2O:2	B2B:2
	2M3	D2O:1	D2B:1
	2R1.5	VB2B:1	VB3B:1
	2R2.5	A1O:3	A1B:3
		D2O:3	D2B:3
		VOP:2	VDL:2
		A2O:1	A2B:1
3	3L1.5	B3O:1	VB3B:1
	3L2.5	B4O:1	B4B:1
	3M1	C1O:1	C1B:1
	3M2	OM-P:5%	RED-D:5%
	3M3	A3O:1	A3B:1
	3R1.5	C2O:1	C2B:1
	3R2.5	OM-P:5%	RED-D:5%
		B4O:2	VB4B:2
		rootBO:1	rootBB:1
		A3.5O:1	A3.5B:1
4	4L1.5	C3O:1	C3B:1
	4L2.5	C4O:1	C4B:1
	4M1	rootBO:1	rootBB:1
	4M2	C4O:1	C4B:1
	4M3	OM-G:5%	RED-D:10%
	4R1.5	OM-P:5%	RED-D:5%
	4R2.5	A4O:1	VA4B:1
		rootAO:1	rootAB:1
		C4O:1	C4B:1
		OM-P:10%	RED-D:10%
5	5M1	R3.5O:1	R3.5B:1
	5M2	rootCO:1	rootCB:1
	5M3	OM-P:10%	RED-D:10%
		rootCO:1	rootCB:1

with Vm color

VITA 3D-Master	Vintage Halo (Opaque)	Vintage Halo (Body)	Vintage Halo (Enamel)
0	0M1	W1O:1	W1B:1
	0M2		
	0M3	W2O:1	W2B:1
1	1M1	W3O:1	W3B:1
	1M2	A1O:3	A1B:3
		R1O:3	R1B:3
2	2L1.5	VOP:2	VDL:2
	2L2.5	B1O:1	VB1B:1
	2M1	A2O:1	A2B:1
	2M2	B2O:2	B2B:2
	2M3	D2O:1	D2B:1
	2R1.5	VB2B:1	VB3B:1
	2R2.5	A1O:3	A1B:3
		D2O:3	D2B:3
		VOP:2	VDL:2
		A2O:1	A2B:1
3	3L1.5	B3O:1	VB3B:1
	3L2.5	B4O:1	B4B:1
	3M1	C1O:1	C1B:1
	3M2	OM-P:5%	RED-D:5%
	3M3	A3O:1	A3B:1
	3R1.5	C2O:1	C2B:1
	3R2.5	OM-P:5%	RED-D:5%
		B4O:2	VB4B:2
		rootBO:1	rootBB:1
		A3.5O:1	A3.5B:1
4	4L1.5	C3O:1	C3B:1
	4L2.5	C4O:1	C4B:1
	4M1	rootBO:1	rootBB:1
	4M2	C4O:1	C4B:1
	4M3	OM-G:5%	RED-D:10%
	4R1.5	OM-P:5%	RED-D:5%
	4R2.5	A4O:1	VA4B:1
		rootAO:1	rootAB:1
		C4O:1	C4B:1
		OM-P:10%	RED-D:10%
5	5M1	R3.5O:1	R3.5B:1
	5M2	rootCO:1	rootCB:1
	5M3	OM-P:10%	RED-D:10%
		rootCO:1	rootCB:1

**Color Conversion Chart**

VITA 3D-Master	Vintage AL & ZR
1M1	A1+1/4 W3
1M2	B2
2L1.5	B2+1/4 B3
2L2.5	B3+1/3 B2
2M1	B1+1/2 OA-Y <sup>1</sup>
2M2	A2
2M3	B3+1/4 B2
2R1.5	A1+1/3 OA-Y <sup>1</sup>
2R2.5	A3
3L1.5	D3
3L2.5	B3
3M1	C1
3M2	B3+A2
3M3	B4
3R1.5	B2+1/3 OA-Y <sup>1</sup>
3R2.5	A3+1/4 A3.5
4L1.5	C3
4L2.5	A4
4M1	D3+1/3 OA-R <sup>2</sup>
4M2	D3+1/3 A3.5
4M3	A4+1/3 OA-R <sup>2</sup>
4R1.5	C2+OA-Y <sup>1</sup>
4R2.5	A3.5+1/3 OA-V <sup>3</sup>
5M1	1/4 A3+ 3/4 C4+ 1/3 OA-V <sup>3</sup>
5M2	A4+OA-V <sup>3</sup>
5M3	3/4 A4+ 1/4B4+ OA-V <sup>3</sup> + OA-R <sup>2</sup>

OA-Y<sup>1</sup> - Opal Amber Yellow

OA-R<sup>2</sup> - Opal Amber Red

OA-V<sup>3</sup> - Opal Amber Violet

