

Lab. Accesories

Air Blow Guns



Air Blow Guns - KEYSTONE

Air Blow Guns - KEYSTONE

Code	Name	Price	+
446905	Keystone Air Gun Nozzle Only	US\$ 42.45	+Info

RUBBER GUN - HANDLER

RUBBER GUN - HANDLER

Metal, nickel plated blow gun activated by depressing lever located on top of gun.



Code	Name	Price	+
103829	#218 HANDLER Lever Blow Gun - 4.19in x 2in x 1.12in - (blow-off nozzle) Metal ...	US\$ 62.11	+Info
103830	#217 HANDLER Push Button Blow Gun -3.5in x 1.62in x 0.81in - (blow-off nozzle)...	US\$ 49.43	+Info
103831	#216 HANDLER Rubber Blow Gun- 3.062in x 0.62in - (blow-off nozzle) A rubber ai...	US\$ 18.00	+Info

Lab. Accesories

Bottles & Containers



BEAKER - ADS

[BAKER - ADS](#)

Features:

- * Can be plastic, disposable or glass
- * Easy to read graduations
- * Durable for repeated use
- * Ideal for caustics, corrosives and acids



US DENTAL DEPOT is an Authorized American Dental Supply dealer in US

DISPOSABLE

Code	Name	Price	+
117662	ADS Beaker disposable 1 oz. pk (25) - #B895-1	US\$ 21.43	+Info

GLASS BEAKER

Code	Name	Price	+
117671	ADS Beaker & Lid glass - 600 ml. ea - #B894-7	US\$ 30.36	+Info
117670	ADS Beaker & Lid glass - 400 ml. ea - #B894-6	US\$ 25.69	+Info
117669	ADS Beaker & Lid glass - 250 ml. ea - #B894-5	US\$ 18.82	+Info
117668	ADS Beakers glass - 150 ml. ea - #B894-10	US\$ 17.88	+Info
117667	ADS Beakers glass - 50 ml. ea. - #B894-8	US\$ 13.15	+Info

PLASTIC

Code	Name	Price	+
117666	ADS Beaker plastic - 600 ml. ea - #B894-4	US\$ 21.43	+Info
117665	ADS Beaker plastic - 400 ml. ea - #B894-3	US\$ 13.82	+Info
117664	ADS Beaker plastic - 250 ml. ea - #B894-2	US\$ 14.72	+Info
117663	ADS Beaker plastic - 100 ml. ea - #B894-1	US\$ 13.34	+Info

BOTTLE DROPPERS - ADS

BOTTLE DROPPERS - ADS

Features:

- * Brown and clear glass bottles and glass droppers
- * Plastic with cone top, spray head, pump, wash narrow mouth, wash wide mouth or applicator tubing



BROWN GLASS

Code	Name	Price	+
117679	ADS Bottle - BROWN Glass screw cap 2oz. - #E123	US\$ 3.70	+Info
117678	ADS Bottle - BROWN Glass dropper top 4oz. - #E115-BN	US\$ 4.96	+Info
117942	ADS Dropper / bottle glass BROWN 2oz. - #D958-2	US\$ 7.92	+Info
117943	ADS Dropper / bottle glass BROWN 4oz - #D958-7	US\$ 8.40	+Info

CLEAR GLASS

Code	Name	Price	+
117681	ADS Bottle - CLEAR Glass with brush cap round 1/2oz. - #E113-B	US\$ 3.70	+Info
117683	ADS Bottle - CLEAR Glass dropper top 4oz. - #E115	US\$ 4.96	+Info
117682	ADS Bottle - CLEAR Glass with brush cap square 1 oz. pk of 4 - #E113-2	US\$ 10.92	+Info
117680	ADS Bottle - CLEAR Glass with brush cap 1 oz. pk of 12 - #E113-A	US\$ 17.27	+Info
117944	ADS Droppers glass gross - #D958-51	US\$ 141.79	+Info

CONE TOP

Code	Name	Price	+
117684	ADS Bottle - Cone Top with RED tip 8oz. pk of 4 - #E119-A	US\$ 12.66	+Info
117685	ADS Bottle - Cone Top with RED tip 2oz. pk of 4 - #E119	US\$ 13.31	+Info

POLYSTOP

Code	Name	Price	+
117688	ADS Bottle - Dropper POLYSTOP 2oz. ea. - #B116-2	US\$ 21.58	+Info
117689	ADS Bottle - Dropper POLYSTOP 4oz. ea. - #B116-1	US\$ 30.33	+Info

SPRAY HEAD

Code	Name	Price	+
117691	ADS Bottle - Spray with spray head 4oz. (4) - #E112	US\$ 13.31	+Info
117692	ADS Bottle - Spray with spray head 8oz. - #E116	US\$ 4.96	+Info

W/ APPLICATOR TUBING

Code	Name	Price	+
117686	ADS Bottle - Dropper with applicator tubing 0.5oz. - #E111-A	US\$ 2.61	+Info
117687	ADS Bottle - Dropper with applicator tubing 1 oz. - #E118	US\$ 3.70	+Info

W/ PUMP

Code	Name	Price	+
117690	ADS Bottle - Pump with pump 16oz. pt (2) - #E111	US\$ 9.74	+Info

WASH NARROW MOUTH

Code	Name	Price	+
117694	ADS Bottles - Wash narrow mouth economy 16oz.- #E114	US\$ 7.88	+Info

WASH WIDE MOUTH

Code	Name	Price	+
117693	ADS Bottles - Wash wide mouth 8oz. (2) - #B116-4	US\$ 21.90	+Info

GLASS DROPPERS

Code	Name	Price	+
117947	ADS Droppers - 6in plastic graduated pk (10) - #D958-4	US\$ 12.34	+Info
117946	ADS Droppers - 3in glass pk (10) - #D958-5	US\$ 14.17	+Info
117945	ADS Droppers - 4 1/2in glass pk (10) - #D958-6	US\$ 21.47	+Info

Bottles & Containers - ADS

Bottles & Containers - ADS

Code	Name	Price	+
118311	Pmp-It Chemical Dispenser 8 02	US\$ 14.94	+Info

Bottles & Containers - JSP

Bottles & Containers - JSP

Code	Name	Price	+
461476	Jsp - Leakproof Glass Vial 1 3/3x9/	US\$ 13.87	+Info
461488	Jsp - Dial Container 11 Compartment	US\$ 1.89	+Info
461487	Jsp - Dial Container 22 Compartment	US\$ 4.55	+Info
461486	Jsp - Storage Box Crystal Clear 10	US\$ 0.98	+Info
461481	Jsp - Multi - Colored Aluminum Cont	US\$ 7.79	+Info
461480	Jsp - Leakproof Glass Vial 2 3/18x5	US\$ 14.34	+Info
461479	Jsp - Leakproof Glass Vial 144 Piec	US\$ 20.07	+Info
461478	Jsp - Leakproof Glass Vial 72 Piece	US\$ 17.59	+Info
461477	Jsp - Leakproof Glass Vial 144 Piec	US\$ 23.17	+Info
461463	Jsp - Stackable Jars 40mm Set of 7	US\$ 6.42	+Info
461475	Jsp - Leakproof Glass Vial 1 3/3x9/	US\$ 18.76	+Info
461474	Jsp - Leakproof Glass Vial 1 3/8x9/	US\$ 10.50	+Info
461473	Jsp - Leakproof Glass Vial 1 3/8x9/	US\$ 18.00	+Info
461472	Jsp - Leakproof Glass Vial 1x9.16in	US\$ 13.33	+Info
461471	Jsp - Leakproof Glass Vial 1x9.16in	US\$ 17.75	+Info
461468	Jsp - Stackable Jars 70mm Set of 4	US\$ 6.42	+Info
461467	Jsp - Stackable Jars 52mm 6 Sets of	US\$ 14.04	+Info
461464	Jsp - Stackable Jars 52mm Set of 5	US\$ 6.42	+Info

CYLINDER GRADUATER - ADS



DOUBLE SKALER

Code	Name	Price	+
117735	ADS Cylinder Graduated double scale 100 ml. - #C893-10	US\$ 26.04	+Info

GLASS

Code	Name	Price	+
117742	ADS Cylinder Graduated glass - 500 ml. ea - #C893-12	US\$ 112.90	+Info
117741	ADS Cylinder Graduated glass - 250 ml. ea - #C893-7	US\$ 48.40	+Info
117740	ADS Cylinder Graduated glass - 100 ml. ea - #C893-6	US\$ 41.88	+Info
117739	ADS Cylinder Graduated glass - 50 ml ea - #C893-5	US\$ 30.81	+Info
117738	ADS Cylinder Graduated glass - 25 ml ea - #C893-4	US\$ 30.81	+Info
117737	ADS Cylinder Graduated glass - 10 ml. ea - #C893-3	US\$ 26.36	+Info

PLASTIC

Code	Name	Price	+
117748	ADS Cylinder Graduated plastic - 1000 ml. ea - #C893-15	US\$ 63.70	+Info
117747	ADS Cylinder Graduated plastic - 250 ml. ea - #C893-11	US\$ 39.08	+Info

117746	ADS Cylinder Graduated plastic - 100 ml. ea - #C893-10	US\$ 26.04	+Info
117745	ADS Cylinder Graduated plastic - 50 ml. ea - #C893-9	US\$ 19.73	+Info
117744	ADS Cylinder Graduated plastic - 25 ml. ea - #C893-8	US\$ 19.53	+Info
117743	ADS Cylinder Graduated plastic - 10 ml. ea - #C893-2	US\$ 19.53	+Info

LAB LIQUID MEDIA - ADS

LAB LIQUID MEDIA - ADS



This product is available in 0.5 and 4 onzas sizes

Code	Name	Price	+
118122	ADS Lab-Liquid Media 4 oz. - # 410230	US\$ 69.46	+Info
118121	ADS Lab-Liquid Media 0.5oz. - #410229	US\$ 15.30	+Info

Lab. Accesories

Chemical Powders



CALCIUM HYDROXIDE POWDER - ADS

CALCIUM HYDROXIDE POWDER - ADS



Features:

* Chemically pure

CALCIUM CARBONATE

Code	Name	Price	+
117700	ADS Calcium Carbonate 4oz. - #S883-31	US\$ 15.14	+Info

CALCIUM HYDROXIDE

Code	Name	Price	+
117701	ADS Calcium Hydroxide USP 4oz. - #S883-10	US\$ 10.18	+Info
117702	ADS Calcium Hydroxide powder 1lb - #S883-10-5	US\$ 24.26	+Info

Lab. Accesories

General Accesories all Brands



General Accesories all Brands - BUFFALO

General Accesories all Brands - BUFFALO

Code	Name	Price	+
424457	Chip Syringe Versati	US\$ 18.18	+Info
424700	Chip Syringe Red Bul	US\$ 9.82	+Info

General Accesories all Brands - JSP

General Accesories all Brands - JSP

Code	Name	Price	+
437687	JSP Zip Lock PLASTIC BAG 3" x 5" 2	US\$ 13.12	+Info
461297	Jsp - Zip Lock Plastic Bag 10in x 8	US\$ 99.99	+Info
461630	Jsp - Plastic Bag - Self Sealing 2x	US\$ 5.32	+Info
461631	Jsp - Plastic Bag - Self Sealing 3x	US\$ 10.12	+Info
461632	Jsp - Plastic Bag - Self Sealing 3	US\$ 12.00	+Info
461633	Jsp - Plastic Bag - Self Sealing 8x	US\$ 30.00	+Info
461634	Jsp - Plastic Bag - Self Sealing 7x	US\$ 23.45	+Info
461635	Jsp - Plastic Bag - Self Sealing 1	US\$ 4.60	+Info
461636	Jsp - Plastic Bag - Self Sealing 4x	US\$ 16.91	+Info
461637	Jsp - Zip Seal Plastic Bags 3inx4in	US\$ 19.99	+Info

Lab. Accesories

Liquids & Solutions



ACETONE - ADS



Code	Name	Price	+
117625	ADS Acetone 1quart - #S883-1-6 (Acetone)	US\$ 23.04	+Info
117624	ADS Acetone 16oz. - #S883-1-5 (Acetone - (ship declaration as HAZ (like monom...	US\$ 16.42	+Info
117623	ADS Acetone 8oz. - #S883-1-4 (Acetone)	US\$ 10.95	+Info

ANTI GREENING LIQUID - ADS

ANTI GREENING LIQUID - ADS

Features:

- * Simple to use as a paint or dip
- * No extra firing steps required
- * Prevents greening and discoloration of porcelain margins when semi-precious alloy framework is used



- * Retards migration of oxides into the porcelain resulting in more life-like and vital appearing ceramic-metal restorations
- * This is the perfect partner for Furnace Decontaminator

US DENTAL DEPOT is an Authorized American Dental Supply dealer in US

Code	Name	Price	+
117657	ADS Anti GREENing Liquid refill pt 16oz. - #P672-5	US\$ 45.60	+Info
117656	ADS Anti GREENing Liquid refill 8oz. - #P672-4	US\$ 30.41	+Info
117655	ADS Anti GREENing Liquid with spray head 4oz. - #P672-3	US\$ 21.04	+Info
117654	ADS Anti GREENing Liquid dropper btl 1 oz. - #P672-1	US\$ 7.94	+Info

ACID NEUTRALIZER - ADS

[ACID NEUTRALIZER - ADS](#)

Features:

- * Osha compliance, eliminate acids including hydrofluoric
- * Powder instantly destroys all dangerous acid wastes from etching, cleaning, pickling, polishing, plating or stripping procedures
- * Specifically formulated for emergency clean up of laboratory acid spills and treatment of skin contact
- * Built-in, indicating color change shows you when its safe to rinse neutralized acids away
- * No waste, is a buffered blend of neutralizing agents effective against all common lab acids particularly Hydrofluoric and Stripping acids
- * Renders spilled acids safe for disposal
- * Store near your acids and keep it handy by your bench.



US DENTAL DEPOT is an Authorized American Dental Supply dealer in US

Code	Name	Price	+
117627	ADS Acid Neutralizer Container 2 lb. - #G864-9	US\$ 26.81	+Info
117626	ADS Acid Neutralizer Container 1lb. - #G864-91	US\$ 15.65	+Info
117934	ADS Dip-It NEUTRALIZER 10 per pk. - #A690-01	US\$ 22.14	+Info

ARMOR SHINE LIQUIDS - WHIPMIX

ARMOR SHINE LIQUIDS - WHIPMIX

Armor Shine improves case presentations by providing a protective, high gloss coating to gypsum models. A single coat prevents dirt and oils from staining the model and protects against shipping. With less than one micron in thickness, Armor Shine does not create occlusal interference. Ideal for producing diagnostic wax-ups, study, C&B and ortho models of the highest quality.



Code	Name	Price	+
120406	WHIP MIX Armor Shine Lq 3.5 Oz - # 09678	US\$ 10.25	+Info
120405	WHIP MIX Armor Shine Lq Liter - # 09679	US\$ 32.07	+Info

ETHYL ALCOHOL DENATURED - ADS

ETHYL ALCOHOL DENATURED 200 PROOF - ADS

Features:

- * Cleans
- * Removes oils & lubricants
- * Burning alcohol
- * Use in pressurized alcohol torches
- * General purpose solvent



Code	Name	Price	+
117760	ADS DenatuRED Alcohol (ETHYL ALCOHOL) gal. - #A966-8 (ships like monomer)	US\$ 41.01	+Info
117759	ADS DenatuRED Alcohol (ETHYL Alcohol) 1qt - # A966-6	US\$ 12.45	+Info

FOILCOTE LIQUID FOIL SUBSTITUTE - WHIPMIX

FOILCOTE LIQUID FOIL SUBSTITUTE - WHIPMIX

When applied to plaster or stone, it seals the pores and forms a smooth surface that gives excellent separation of the denture from the mold. With this product eliminates the time-consuming tin foiling and resulting creases.



Code	Name	Price	+
120409	WHIP-MIX Foilcote Lq 3.5 Oz - #09687	US\$ 10.96	+Info
120408	WHIP-MIX Foilcote Lq RED Gallon/3.75l - #09709	US\$ 78.62	+Info
120407	WHIP-MIX Foilcote 4gal Case - #27782	US\$ 269.50	+Info

Liquids & Solutions - ADS

Liquids & Solutions - ADS

Code	Name	Price	+
117677	Boric Acid Pwdr 4oz	US\$ 12.80	+Info

Liquids & Solutions - JSP

Liquids & Solutions - JSP

Code	Name	Price	+
461607	Jsp - Debubblizer - 32 Oz (946ml)	US\$ 16.89	+Info

Liquids & Solutions - WHIP MIX

Liquids & Solutions - WHIP MIX

Code	Name	Price	+
120433	Green Envy Tray Cleaner 1#Jar	US\$ 34.13	+Info

LUBRITEX N 12 - DIE LUBRICANT - WHIPMIX

LUBRITEX N 12 - DIE LUBRICANT - WHIPMIX

A water-soluble lubricant that is equally effective on all types of die materials. LUBRITEX will not react with die or pattern materials and will not affect investments.



Code	Name	Price	+
120419	WHIP-MIX Lubritex #12 3.5 Oz - #09784	US\$ 20.36	+Info
120418	WHIP MIX Lubritex #12 Pint/450ml - # 09792	US\$ 49.66	+Info

MODEL GLOW SOAP - WHIPMIX

MODEL GLOW SOAP - WHIPMIX

MODEL GLOW gives an incredibly high luster to gypsum models. Its special formulation does not yellow or flake and seals the pores to provide a protective coating. Its thin consistency does not affect occlusal accuracy. Long lasting and can be used full strength or diluted. Models can be buffed wet or dry.



Code	Name	Price	+
120421	WHIP-MIX Model Glow Gallon/3.75l - #29009	US\$ 57.61	+Info
120420	WHIP-MIX Model Glow 4gal Case - #29017	US\$ 190.40	+Info

PLASTIC WATER GRADUATES - W HIPMIX

PLASTIC WATER GRADUATES - WHIPMIX

The water/powder ratio of investments or die stones must be maintained precisely to secure the most accurate results. Whip Mix Water Graduates provide the means for measuring the water (or liquid) accurately. Water Graduates are plastic and available in two sizes, 26mL and 100mL, and are marked off in 1mL and 2mL intervals, respectively. The flanged base of each graduate provides stability.



Code	Name	Price	+
120457	WHIP-MIX Plastic Water Graduate 100ml - #05630	US\$ 6.99	+Info
120458	WHIP-MIX Plastic Water Graduate 26ml - #05622	US\$ 8.40	+Info

RING LINER - NON-ABESTOS CEL LULOSE - WHIPMIX

RING LINER - NON-ABESTOS CELLULOSE - WHIPMIX

A liner that has no asbestos and none of its inherent hazards. This material has the necessary cushioning, but more importantly the correct water absorption quality which feeds to the investment to provide and control setting expansion. During burnout it chars to a black ash which crumbles easily when the mold is pushed out of the ring.

Liners should be 3 (1/8") short of each end of the casting ring.



Code	Name	Price	+
120454	WHIP-MIX Ring Liner Roll Cellulose 3in - #14656	US\$ 25.43	+Info
120450	WHIP-MIX Ring Liner Roll Cellulose 50mm - #16195	US\$ 22.17	+Info
120453	WHIP-MIX Ring Liner Roll Cellulose 1-5/8in - #14648	US\$ 18.60	+Info
120452	WHIP-MIX Ring Liner Roll Cellulose 1-3/8 X60 - #14621	US\$ 16.24	+Info
120451	WHIP-MIX Ring Liner Roll Cellulose 1-1/8 X60 - #14605	US\$ 14.73	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US
DENTAL
DEPOT

+1 954 874 6325 
www.usdentaldepot.com 

SMOOTHX DEBUBBLING SOLUTION - WHIPMIX
SMOOTHX DEBUBBLING SOLUTION -WHIPMIX

A surface tension reducing liquid for treating patterns to permit better adherence and flow of investment.
Also excellent for treating synthetic impressions (Polysulfide, Silicone, or Polyether) before pouring casts.
Refillable, fine-spray bottles are easy to use and help keep debubbler clean by eliminating the need for dipping brushes or patterns.



Code	Name	Price	+
120431	WHIP-MIX Smoothex Debubblizer Gallon/3.75l - #12912	US\$ 302.27	+Info
120429	WHIP-MIX Smoothex Debubblizer Liter - #29394	US\$ 88.38	+Info
120428	WHIP-MIX Smoothex Pint/450ml - #09849	US\$ 47.60	+Info
120432	WHIP-MIX Smoothex Fine Spray Btl 8oz/225ml - #02313	US\$ 39.41	+Info
120430	WHIP-MIX Smoothex Fine Spray Bottle 4oz - #29386	US\$ 21.97	+Info

Catalog Index / Indice

Air Blow Guns	1
Bottles & Containers	1
Chemical Powders	5
General Accesories all Brands	5
Liquids & Solutions	6