

Acrylic Accessories

Acrylic Bonding



ACRYLIC BONDING - ADS

ACRYLIC BONDING - ADS

Features:

- * Enhances bond of acrylic to teeth
- * Prevents pop outs
- * Use as a bonding agent in reline procedures



US DENTAL DEPOT is an Authorized American Dental Supply dealer in US

Code	Name	Price	+
117630	ADS Acrylic Bonder 8oz. - #A130-4	US\$ 39.11	+Info
117629	ADS Acrylic Bonder 4oz. - #A130.3	US\$ 33.36	+Info

Acrylic Accessories

Stains for Acrylic



JET ADJUSTERS KIT- LIQUID RESIN SYSTEM - LANG

JET ADJUSTERS KIT- LIQUID RESIN SYSTEM - LANG

Features:

Liquid resin system for custom characterizing and staining temporary crowns and bridges.

This liquid resin system allows you to customize your newly created temporary crown or bridge for a more natural appearance to match the patient's dentition. We offer eight long-lasting, translucent and stable colors that bond directly to the fabricated tooth.

Jet Adjusters cure quickly (less than 2 minutes) and can be used chairside or in a lab prior to final cementation. Sealing properties eliminate porosity.

Liquid: Blue, Brown, Clear, Grey, Orange, Violet, White, Yellow, Thinner, and Brush Cleaner

Contain: Kit - bottle of each color 1x10 ml



REFILLS

Brush Cleaner

Brush Cleaner

Code	Name	Price	+
------	------	-------	---

114110	Jet Adjuster Liquid 10ml Brush	US\$ 24.57	+Info
--------	--------------------------------	------------	-----------------------

Colors
Colors

Code	Name	Price	+
114101	Jet Adjusters Liquid 10ml Blue	US\$ 24.57	+Info
114102	Jet Adjuster Liquid 10ml Brown	US\$ 24.57	+Info
114103	Jet Adjuster Liquid 10ml Clear	US\$ 24.57	+Info
114104	Jet Adjuster Liquid 10ml Grey	US\$ 24.57	+Info
114105	Jet Adjuster Liq 10ml Orange	US\$ 24.57	+Info
114106	Jet Adjuster Liq 10ml Violet	US\$ 24.57	+Info
114107	Jet Adjuster Liq 10ml White	US\$ 24.57	+Info
114108	Jet Adjuster Liquid 10ml Yello	US\$ 24.57	+Info

Thinner
Thinner

Code	Name	Price	+
114109	Jet Adjuster Liquid 10ml Thinn	US\$ 25.55	+Info

KIT

Code	Name	Price	+
114100	LANG Jet Adjusters Liquid Resin Kit 1x10 ml bottle of each color - Acrylic Stai...	US\$ 110.42	+Info

MINUTE STAIN - ACRYLIC RESIN STAINS - T AUB

MINUTE STAIN - ACRYLIC RESIN STAINS - TAUB

Minute Stains are quick setting colored acrylic liquids that are applied by brush to modify the shades of provisional and permanent acrylic and composite restorations. The stains cure and bond to all dental resins, including ethyl and methyl methacrylate, bis-acryl resins, polycarbonates, vinyl methacrylate copolymers, resin crowns and laminates, denture base, acrylic denture teeth and composite resins. Minute Stains can be applied in or out of the mouth.



Code	Name	Price	+
122415	TAUB - Minute Stain Kit - 3 Acrylic Stains Standard Color - #01-1000 YELLOW, C...	US\$ 120.38	+Info
122416	TAUB - Minute Stain Kit - 3 Color - Acrylic Stains Special - #01-1001 - # 01-1...	US\$ 120.38	+Info
475055	TAUB - Minute Stain Thinner 1/2 oz #01-4000	US\$ 7.22	+Info
475056	TAUB - Minute Stain Thinner 1/2 oz #01-4001	US\$ 9.81	+Info
122417	TAUB - Denture Base Stain Kit - 3 Color Acrylic Stains - #01-1010	US\$ 120.38	+Info
100436	Taub- Minute-Stains Kit 7 Colors Acrylic Stains (BLUE- BROWN - Cervical - GRAY ...	US\$ 218.36	+Info

Stains for Acrylic - OPTI-CRYL

Stains for Acrylic - OPTI-CRYL

Code	Name	Price	+
476928	Opti-Cryl Thermopolymerizable	US\$ 45.00	+Info

Catalog Index / Indice

Acrylic Bonding	1
Stains for Acrylic	1