

Acrylic Heat cured

Cross Linked Acrylics



OPTI-CRYL HEAT CURE ACRYLIC | NEW STETIC

OPTICRYL HEAT CURE ACRYLIC - NEW STETIC

Color stable Denture Base resin with excellent working conditions
 This Denture Acrylic resin will be compatible with all regular Denture Resin Acrylics in the market and can replace the great majority of them
 Its instructions will not surprise you, and make you move out of your usual way in order to process this Denture Resin.
 You will be able to process this material in the same way you are processing your Denture Acrylic resin today



Denture Resin Acrylic Heat curel, is available in 1 lb - 5 lb and 22 lb
 And as a Combo including the monomer

The proportions of use for the mixture are:
 By weight: Two parts of curing polymer Opti-Cryl and one part of curing monomer Opti-Cryl.
 By volume: Three parts of curing polymer Opti-Cryl and one part of curing monomer Opti-Cryl.

Opticryl Denture Resin Acrylic is FDA approved

EXCLUSIVE DISTRIBUTOR
 FDA Approved
 Made in Colombia

1lb & liquid

Clear
 Clear

Code	Name	Price	+
100404	Heat Cure Clear 1lb+8oz Mon	US\$ 31.21	+Info

Dark Pink Veined
 Dark Pink Veined

Code	Name	Price	+
100417	Heat Cure Dk Pink 1lb+8oz Mon	US\$ 31.21	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
100403	Heat Cure Lt Pink 1lb+8oz Mon	US\$ 31.21	+Info

Meharry Veined

Meharry Veined

Code	Name	Price	+
100416	Heat Cure Meharry 1lb+ 8oz Mon	US\$ 31.21	+Info

Original

Original

Code	Name	Price	+
100402	Heat Cure Original 1lb+8oz Mon	US\$ 31.21	+Info

1lb Powder Only

Clear

Clear

Code	Name	Price	+
105107	Heat Cure Clear 1lb Pwd.Only	US\$ 15.60	+Info

Dark Pink Veined

Dark Pink Veined

Code	Name	Price	+
105102	Heat Cure Dk Pink 1lb Pwd.Only	US\$ 15.60	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
105108	Heat Cure Lt Pink 1lb Pwd.Only	US\$ 15.60	+Info

Meharry Veined

Meharry Veined

Code	Name	Price	+
105109	Heat Cure Meharry 1lb Pwd.Only	US\$ 15.60	+Info

Original

Original

Code	Name	Price	+
103099	Heat Cure Orig. 1lb Pwd.Only	US\$ 15.60	+Info

5lb & liquid

Clear

Clear

Code	Name	Price	+
100414	Heat Cure Clear 5lb+1Qt Mon	US\$ 83.03	+Info

Dark Pink Veined

Dark Pink Veined

Code	Name	Price	+
100407	Heat Cure Dk Pink 5lb+1Qt Mon	US\$ 92.53	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
100413	Heat Cure Lt Pink 5lb+ 1Qt Mon	US\$ 92.53	+Info

Meharry Veined

Meharry Veined

Code	Name	Price	+
100405	Heat Cure Meharry 5lb+1Qt Mon	US\$ 83.98	+Info

Original

Original

Code	Name	Price	+
100412	Heat Cure Original 5lb+1Qt Mon	US\$ 92.53	+Info

5lb Powder Only

Clear

Clear

Code	Name	Price	+
105110	Heat Cure Clear 5lb Pwd.Only	US\$ 53.20	+Info

Dark Pink Veined

Dark Pink Veined

Code	Name	Price	+
105111	Heat Cure Dk Pink 5lb Pwd.Only	US\$ 62.70	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
105112	Heat Cure Lt Pink 5lb Pwd.Only	US\$ 62.70	+Info

Meharry Veined

Meharry Veined

Code	Name	Price	+
105113	Heat Cure Meharry 5lb Pwd.Only	US\$ 54.15	+Info

Original

Original

Code	Name	Price	+
105103	Heat Cure Orig. 5lb Pwd.Only	US\$ 62.70	+Info

22lb & liquid

Clear

Clear

Code	Name	Price	+
100421	Heat Cure Clear 22lb+ 4Qt Mon	US\$ 319.90	+Info

Dark Pink Veined

Dark Pink Veined

Code	Name	Price	+
100422	Heat Cure Dk Pink 22lb+4Qt Mon	US\$ 290.15	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
100419	Heat Cure Lt Pink 22lb+4Qt Mon	US\$ 290.15	+Info

Meharry Veined

Meharry Veined

Code	Name	Price	+
100418	Heat Cure Meharry 22lb+4Qt Mon	US\$ 290.15	+Info

Original

Original

Code	Name	Price	+
100420	Heat Cure Original 22lb+4Qt Mo	US\$ 290.15	+Info

22lb Powder Only

Clear

Clear

Code	Name	Price	+
105104	Heat Cure Clear 22lb Pwd.Only	US\$ 210.00	+Info

Dark Pink Veined

Dark Pink Veined

Code	Name	Price	+
105100	Heat Cure Dk Pink 22lb Pwd.Onl	US\$ 180.25	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
105105	Heat Cure Lt Pink 22lb Pwd.Onl	US\$ 180.25	+Info

Meharry Veined

Meharry Veined

Code	Name	Price	+
105106	Heat Cure Meharry 22lb Pwd.Onl	US\$ 180.25	+Info

Original

Original

Code	Name	Price	+
105101	Heat Cure Orig. 22lbo Pwd.Only	US\$ 180.25	+Info

Liquid Only

HEAT CURE MONOMER - OPTICRYL

Features:

optimal dimensional stability and structure free of porosities is achieved when the right liquid-to-powder ratio is used as indicated.

Note:

The Liquid can only be shipped by ground transportation, it can not be shipped by air.

Internationally maritime transportation is obligated. Minimum order is 700 us\$ + Hazardous Material Export license : 140 us\$



Code	Name	Price	+
100249	OPTI-CRYL Heat Curing Monomer Liquid Only 8oz -	US\$ 16.99	+Info
100250	OPTI-CRYL Heat Curing Monomer Liquid Only 32oz/ 1Quarter -	US\$ 33.35	+Info

Cross Linked Acrylics - COLTENE

Cross Linked Acrylics - COLTENE

Code	Name	Price	+
436292	Hygenic Denture Resin Material, Spe	US\$ 0.00	+Info
436302	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436305	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436311	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436320	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436322	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436323	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436331	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436333	Hygenic Denture Resin Material, Lab	US\$ 0.00	+Info
436336	Hygenic Denture Resin Material, Cle	US\$ 0.00	+Info

MILLENIUM Denture Acrylic | KEYSTONE

Millennium Heat cure Acrylic - KEYSTONE



Liquid only

Heat Cure Monomer - KEYSTONE

Code	Name	Price	+
446846	Millennium Acrylic - Liquid 8 oz, H	US\$ 14.03	+Info

Millenium

Millenium

Code	Name	Price	+
446235	Millennium Acrylic - Light Veined P	US\$ 38.55	+Info
446410	Millennium Acrylic Lt Veined P&L5lb	US\$ 125.57	+Info

SLEDGEHAMMER HEAT DENTURE ACRYLIC | KEYSTONE

SLEDGEHAMMER HEAT CURE ACRYLIC - KEYSTONE

- * FDA Compliant Quality System
- * Color Stable
- * No Bleaching, blanching or streaking
- * Manufactured under strict ISO and CGMP Guidelines for Guaranteed Quality and Product Reliability from batch to batch
- * Overnight or Shortcure Flexibility
- * High Impact Resistance ensures against breakage
- * 12-14 minute working times enables more units to be mixed simultaneously



20 MINUTES Liquid Only

20 MINUTES HEAT CURE SLEDGEHAMMER - KEYSTONE

An excellent product for the rapid curing of Sledgehammer Heat Cure Denture acrylics.

Kits with Sledgehammer Heat Cure powder and '20 minute' available upon request.

Note:

The Liquid can only be shipped by ground transportation, it can not be shipped by air.

Internationally maritime transportation is obligated. Minimum order is 700 us\$ + Hazardous Material Export license : 140 us\$



Code	Name	Price	+
114837	SLEDGEHAMMER 20-Minute Monomer, 8 Ounce Bottle. For rapid curing heat - #1000590	US\$ 19.50	+Info
034-1000591	SLEDGEHAMMER 20-Minute Monomer, 1 Quart Bottle. For rapid curing heat - #1000591	US\$ 37.20	+Info
34-1000517	SLEDGEHAMMER 20-Minute Monomer, 1 Quart Bottle. For rapid curing heat - #1000517	US\$ 41.10	+Info

Powder & Liquid

1lb

Light Fiber Plus

Light Fiber Plus

Code	Name	Price	+
34-1001965	Heat Cure Lig Fib 1lb+8oz P&L	US\$ 40.80	+Info

Light Fibered

Light Fibered

Code	Name	Price	+
034-1000535	Heat Cure Li Fiber 1lb+8oz P&L	US\$ 40.47	+Info

Orig.Light Redish

Orig.Light Redish

Code	Name	Price	+
34-1000470	Heat Cure LRP 1lb+8oz P&L	US\$ 40.80	+Info
34-1001945	Heat Cure LRP 1lb+8oz P&L	US\$ 40.80	+Info

Original

Original

Code	Name	Price	+
34-1000471	Heat Cure Original 1lb+8oz P&L	US\$ 40.80	+Info

Others

Others

Code	Name	Price	+
------	------	-------	---

445960	Sledgehammer Self Cure Denture Acry	US\$ 44.03	+Info
--------	-------------------------------------	------------	-----------------------

5lb

Clear

Clear

Code	Name	Price	+
34-1000475	Heat Cure Clear 5lbs+32oz P&L	US\$ 131.95	+Info

Light Fiber Plus

Light Fiber Plus

Code	Name	Price	+
34-1001966	Heat Cure Lig Fib 5lbs+32oz P&L	US\$ 130.41	+Info

Light Fibered

Light Fibered

Code	Name	Price	+
34-1000536	Heat Cure Lig Fib 5lbs/32oz P&L	US\$ 131.95	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
34-1000473	Heat Cure LRP 5lbs+32oz P&L	US\$ 120.07	+Info

Orig.Light Redish

Orig.Light Redish

Code	Name	Price	+
34-1001946	Heat Cure LRP 5lbs+32oz P&L	US\$ 131.95	+Info

Original

Original

Code	Name	Price	+
34-1000474	Heat Cure Ori 5lbs+32oz P&L	US\$ 131.95	+Info

25lb & 100lb

Clear

Clear

Code	Name	Price	+
34-1000478	Heat Cure Clear 25lbs+1gal P&L	US\$ 472.75	+Info
446336	Sledgehammer Acrylic Heat-Cure, Cle	US\$ 1,580.54	+Info

Light Fiber Plus

Light Fiber Plus

Code	Name	Price	+
34-1001967	Heat C. Lt.Fib 25lbs+1gal P&L	US\$ 464.17	+Info

Light Fibered

Light Fibered

Code	Name	Price	+
34-1000537	Heat Cure Li Fib 25lbs+1gal P&L	US\$ 472.75	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
34-1000476	Heat Cure LRP 25lbs+1gal P&L	US\$ 472.75	+Info
446834	Sledgehammer Acrylic Heat-Cure, Ori	US\$ 1,580.54	+Info

Orig.Light Redish

Orig.Light Redish

Code	Name	Price	+
34-1001947	Heat Cure LRP 25lbs+1gal P&L	US\$ 378.00	+Info

Original

Original

Code	Name	Price	+
34-1000477	Heat Cure Ori 25lbs+1gal P&L	US\$ 472.75	+Info
446844	Sledgehammer Acrylic Heat-Cure, Ori	US\$ 1,580.54	+Info

Powder only

1lb

Clear

Clear

Code	Name	Price	+
34-1000472	Heat Cure Clear 1lb+8oz P&L	US\$ 40.80	+Info
34-1000502	Heat Cure Clear 1lb Pwd only	US\$ 23.52	+Info

Light Fiber Plus

Light Fiber Plus

Code	Name	Price	+
34-1001961	Heat Cure Lig Fib 1lb Pwd only	US\$ 29.40	+Info

Light Fibered

Light Fibered

Code	Name	Price	+
34-1000543	Heat Cure Lig Fib 1lb Pwd only	US\$ 29.40	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
34-1000494	Heat Cure LRP 1lb Powder only	US\$ 29.40	+Info

34-1001941	Heat Cure LRP 1lb powder only	US\$ 29.40	+Info
------------	-------------------------------	------------	-----------------------

Original
Original

Code	Name	Price	+
34-1000498	Heat Cure Orig 1lb Powder only	US\$ 29.40	+Info

Others
Others

Code	Name	Price	+
446106	Sledgehammer Self-Cure Denture Acry	US\$ 36.68	+Info

5lb

Clear
Clear

Code	Name	Price	+
34-1000503	Heat Cure Clear 5lbs Pwd only	US\$ 90.16	+Info

Light Fiber Plus
Light Fiber Plus

Code	Name	Price	+
34-1001962	Heat Cure Lig Fib 5lbs Pwd only	US\$ 90.16	+Info

Light Fibered
Light Fibered

Code	Name	Price	+
34-1000544	Heat Cure Lig Fib 5lbs Pwd only	US\$ 86.55	+Info

Light Pink Veined
Light Pink Veined

Code	Name	Price	+
34-1000495	Heat Cure LRP 5lbs Powder only	US\$ 90.16	+Info

Orig.Light Redish
Orig.Light Redish

Code	Name	Price	+
34-1001942	Heat Cure LRP 5lbs powder only	US\$ 90.16	+Info

Original
Original

Code	Name	Price	+
34-1000499	Heat Cure Orig 5lbs Powder only	US\$ 90.16	+Info

25lb

Clear

Clear

Code	Name	Price	+
34-1000504	Heat Cure Clear 25lbs Pwd only	US\$ 383.31	+Info

Light Fiber Plus

Light Fiber Plus

Code	Name	Price	+
34-1001963	Heat Cure Lig Fib 25lbs Pwd onl	US\$ 383.31	+Info

Light Fibered

Light Fibered

Code	Name	Price	+
34-1000545	Heat Cure Lig Fib 25lbs Pwd onl	US\$ 383.31	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
34-1000496	Heat Cure LRP 25lbs Powder only	US\$ 383.31	+Info

Orig.Light Redish

Orig.Light Redish

Code	Name	Price	+
34-1001943	Heat Cure LRP 25lbs powder only	US\$ 295.52	+Info

Original

Original

Code	Name	Price	+
34-1000500	Heat Cure Ori 25lbs Powder only	US\$ 383.31	+Info

Acrylic Heat cured

HIGH-Impact Acrylic



DIAMOND D HEAT CURE ACRYLIC - KEYSTONE

DIAMOND D HEAT CURE ACRYLIC - KEYSTONE

Advanced natural-looking dentures require artistry as well as technical expertise. Diamond D was specially formulated for exceptional working and handling characteristics for anatomic sculpting.

This revolutionary new resin has just the right translucency and color variations found in natural tissue as well as ultra impact and flexural strength. The demand for function and aesthetics in dental fabrication is growing. Diamond D can help take a business to a new level.



Liquid only

DIAMOND D HEAT CURE MONOMER - KEYSTONE

To fabricate natural looking dentures it requires artistry, and technical expertise, as well as the materials to support those abilities. Diamond D Ultra Impact Denture Acrylic is now the premier denture acrylic in the industry.

The Notched Izod Impact Strength tests performed at a major polymer testing facility show that Diamond D Ultra Impact not only has a greater strength, but higher impact/fracture resistance than other Hi-Impact denture base acrylics tested.

With its exceptional stability Diamond D is offered as the base polymer for premium, impact resistant dentures offered by dentists and dental clinics who specialize in quality, good fitting dentures and dentures for family members of the dental community.

Note:

The Liquid can only be shipped by ground transportation, it can not be shipped by air.

Internationally maritime transportation is obligated. Minimum order is 700 us\$ + Hazardous Material Export license : 140 us\$



Code	Name	Price	+
34-1013014	KEYSTONE Diamond D Ultra Impact Acrylic, 8 ounce Heat Cure Monomer - #1013014	US\$ 21.38	+Info
34-1013015	KEYSTONE Diamond D Ultra Impact Acrylic, 1 quart Heat Cure Monomer - #1013015	US\$ 48.45	+Info
34-1013094	KEYSTONE Diamond D 20-Minute Monomer Liquid 32 oz/Bt. Heat Cure Ultra Impact Acr...	US\$ 58.58	+Info

Powder & Liquid

1lb

Code	Name	Price	+
34-1013061	Diamond-D Ultra Impact Acrylic, Heat Cure, Original Shade (PINK), 1lb & 8oz Liq...	US\$ 68.40	+Info
34-1013062	Diamond-D Ultra Impact Acrylic, Heat Cure, Light Shade, 1lb & 8oz Liquid monome...	US\$ 68.40	+Info
34-1013063	Diamond-D Ultra Impact Acrylic, Heat Cure, Light REDISH-PINK Shade, 1lb & 8oz ...	US\$ 68.48	+Info
34-1013064	Diamond-D Denture Acrylic, Heat Cure, Chroma-Essence, 1lb & 8oz Liquid monomer...	US\$ 68.40	+Info
34-1013112	Diamond-D Heat Cure Dark Veined 1lb & 8oz Liquid monomer P&L #1013112 -	US\$ 68.40	+Info

5lb

Code	Name	Price	+
1013065	Diamond-D Ultra Impact Acrylic, Heat Cure, Original Shade, 5lb powder & 1qt Liqu...	US\$ 228.48	+Info
1013066	Diamond-D Denture Acrylic, Heat Cure, Light Shade, 5lb powder & 1qt Liquid P&L...	US\$ 221.63	+Info
1013067	Diamond-D Denture Acrylic, Heat Cure, Light REDISH-PINK Shade, 5lb powder & 1qt ...	US\$ 226.46	+Info
1013068	Diamond-D Denture Acrylic, Heat Cure, Chroma-Essence Shade, 5lb powder & 1qt Li...	US\$ 228.48	+Info
1013113	Diamond-D Heat Cure Dark Veined 5lb powder & 1qt Liquid P&L -- P&L #1013113 -	US\$ 228.48	+Info

25lb

Code	Name	Price	+
1013069	Diamond-D Heat Cure Original 25lbs Powder & 4qts liquid Monomer - P&L #1013069 -	US\$ 756.41	+Info
1013070	Diamond-D Heat Cure Light 25lbs Powder & 4qts liquid Monomer - P&L #1013070 -	US\$ 705.65	+Info
1013071	Diamond-D Denture Acrylic, Heat Cure, Light REDISH-PINK Shade, 25lbs Powder & 4...	US\$ 755.69	+Info
1013072	Diamond-D Denture Acrylic, Heat Cure, Chroma-Essence Shade, 25lbs Powder & 4qts ...	US\$ 755.69	+Info
1013114	Diamond-D Heat Cure Powder and Liquid, Dark Vein, 25lbs Powder & 4qts liquid Mo...	US\$ 755.69	+Info

Powder only

1lb

Code	Name	Price	+
34-1013020	Diamond-D Denture Acrylic, Heat Cure, Original Shade, 1lbs Powder only - #1013020	US\$ 57.38	+Info
34-1013024	Diamond-D Denture Acrylic, Heat Cure, Light, 1lbs Powder only - #1013024	US\$ 57.38	+Info
34-1013028	Diamond-D Denture Acrylic, Heat Cure, Light REDISH-PINK, 1lbs Powder only - #10...	US\$ 57.38	+Info
34-1013032	Diamond-D Denture Acrylic, Heat Cure, Chroma-Essence, 1lbs Powder only - #1013032	US\$ 57.38	+Info
34-1013105	Diamond-D Heat Cure Dark Veined 1Lb Powder only #1013105 -	US\$ 57.38	+Info

5lb

Code	Name	Price	+
34-1013106	Diamond-D Heat Cure Dark Veined 5Lb Powder only #1013106 -	US\$ 210.60	+Info
34-1013033	Diamond-D Denture Acrylic, Heat Cure, Chroma-Essence, 5lbs Powder only - #1013033	US\$ 210.60	+Info
34-1013029	Diamond-D Denture Acrylic, Heat Cure, Light REDISH-PINK, 5lbs Powder only - #10...	US\$ 210.60	+Info
34-1013025	Diamond-D Denture Acrylic, Heat Cure, Light, 5lbs Powder only - #1013025	US\$ 210.60	+Info
34-1013021	Diamond-D Denture Acrylic, Heat Cure, Original Shade, 5lbs Powder only - #1013021	US\$ 210.60	+Info

25lb

Code	Name	Price	+
34-1013107	Diamond-D Heat-Cure Powder Only, Dark Vein, Package of 25lbs powder only - #101...	US\$ 696.54	+Info
34-1013034	Diamond-D Denture Acrylic, Heat Cure, Chroma-Essence, 25lbs Powder only - #1013034	US\$ 696.54	+Info
34-1013030	Diamond-D Denture Acrylic, Heat Cure, Light REDISH-PINK, 25lbs Powder only - #1...	US\$ 696.54	+Info
34-1013022	Diamond-D Denture Acrylic, Heat Cure, Original Shade, 25lbs Powder only - #1013022	US\$ 696.54	+Info

GC Nature-Cryl Hi-20ET Ultra High Impact - GC AMERICA

NATURE-CRYL HI-20ET

High Impact, Heat-Cured Denture Base Resin NATURE-CRYL HI-20ET is a high impact denture resins innovative polymer chemistry that allows extended working time and a 20-minute boiling water cure. This product permits computer-assisted shade matching, which provides aesthetically accurate upper and lower dentures, partial dentures, relines and repairs with all NATURE-CRYL system shades.



- Convenient 20-Minute cure - denture sections (up to 8.5mm) can be processed without exothermic porosity
 - Can also be processed using the conventional 2-hour or overnight curing cycles
 - Extended working time - 1 to 1.5 hours at room temperature, can be kept overnight if refrigerated - suitable for injection techniques

Liquid only

Liquid only

Code	Name	Price	+
421999	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 66.10	+Info
422389	Nature-Cryl Hi-20ET Ultra Liq 8oz	US\$ 26.88	+Info

Powder & Liquid

Powder & Liquid

Code	Name	Price	+
421742	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 59.72	+Info
421788	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 59.72	+Info
422043	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 59.72	+Info
422081	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 59.72	+Info
422381	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 59.72	+Info
422382	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 59.72	+Info
422384	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 59.72	+Info

Powder only

Powder only

Code	Name	Price	+
421706	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 178.79	+Info
422010	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 178.79	+Info
422373	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 178.79	+Info
422376	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 178.79	+Info
422383	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 178.79	+Info
422388	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 178.79	+Info
422394	Nature-Cryl Hi-20ET Ultra High Impa	US\$ 178.79	+Info

Shade Guide

Shade Guide

Code	Name	Price	+
422017	Nature-Cryl Shade Guide	US\$ 66.10	+Info

Lucitone 199 - DENTSPLY

Diamond-D Acrylic - DENTSPLY

Code	Name	Price	+
421077	LUCITONE 199 ORIGINAL	US\$ 1,158.02	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



**US
DENTAL
DEPOT**

+1 954 874 6325 
www.usdentaldepot.com 

Opticryl Denture High Impact Acrylic -
OPTI-CRYL

Opticryl High Impact - OPTI-CRYL



Code	Name	Price	+
112962	High Impact Orig. 1lb Pwd.Only	US\$ 39.00	+Info

Acrylic Heat cured
Microwave Cured Acrylic

Nature-CRYL MC MICROWAVE CURED Denture Acrylic | G
C AMERICA

NATURE-CRYL MC Microwave-Cured Denture Base Resin NATURE-CRYL MC is a denture acrylic resin specially created for microwave curing. This product assures that all high quality, precision dentures can be cured in just three minutes in a household microwave oven. Its unique plastic flask allows for more economical denture production. The product's computer-assisted shade matching assures aesthetically accurate upper and lower dentures, partial dentures, relines and repairs with all NATURE-CRYL system shades.



- [REDACTED]
- [REDACTED]
- [REDACTED]

Liquid only

Liquid only

Code	Name	Price	+
421664	Nature-Cryl MC Microwave Cured Dent	US\$ 74.75	+Info
422008	Nature-Cryl MC Microwave Cured Dent	US\$ 41.29	+Info

Powder & Liquid

Powder & Liquid

Clear

Clear

Code	Name	Price	+
422393	Nature-Cryl MC Microwave Cured Dent	US\$ 103.30	+Info

Heavy Shade

Heavy Shade

Code	Name	Price	+
422385	Nature-Cryl MC Microwave Cured Dent	US\$ 103.30	+Info

Light Fibered

Light Fibered

Code	Name	Price	+
421755	Nature-Cryl MC Microwave Cured Dent	US\$ 103.30	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
422371	Nature-Cryl MC Microwave Cured Dent	US\$ 103.30	+Info

Light Shade

Light Shade

Code	Name	Price	+
422375	Nature-Cryl MC Microwave Cured Dent	US\$ 103.30	+Info

Mild Shade

Mild Shade

Code	Name	Price	+
422042	Nature-Cryl MC Microwave Cured Dent	US\$ 103.30	+Info

Moderate Shade

Moderate Shade

Code	Name	Price	+
421693	Nature-Cryl MC Microwave Cured Dent	US\$ 103.30	+Info

Original

Original

Code	Name	Price	+
421754	Nature-Cryl MC Dent Resin Orig Kit	US\$ 103.30	+Info

Powder only

Powder only

Clear

Clear

Code	Name	Price	+
422374	Nature-Cryl MC Microwave Cured Dent	US\$ 288.12	+Info

Light Fibered

Light Fibered

Code	Name	Price	+
421993	Nature-Cryl MC Microwave Cured Dent	US\$ 288.12	+Info

Light Pink Veined

Light Pink Veined

Code	Name	Price	+
422395	Nature-Cryl MC Microwave Cured Dent	US\$ 288.12	+Info

Original

Original

Code	Name	Price	+
422391	Nature-Cryl MC Resin Powd Orig 5lb	US\$ 288.12	+Info

Powder only

Powder only

Code	Name	Price	+
421121	MC NATURE-CRYL Light 5lb Powder	US\$ 288.12	+Info
421597	Nature-Cryl MC Microwave Cured Dent	US\$ 288.12	+Info
421992	Nature-Cryl MC Microwave Cured Dent	US\$ 288.12	+Info
422372	Nature-Cryl MC Microwave Cured Dent	US\$ 288.12	+Info
422380	Nature-Cryl MC Microwave Cured Dent	US\$ 288.12	+Info

Catalog Index / Indice

Cross Linked Acrylics	1
HIGH-Impact Acrylic	11
Microwave Cured Acrylic	15