

## Acrylic Self Cured Diamond -D



### DIAMOND D SELF CURING | KEYSTONE ONE

#### SELF CURING ORIGINAL - DIAMOND D

Advanced natural-looking dentures require artistry as well as technical expertise.

Diamond D was specially formulated for exceptional working and handling characteristics for anatomic sculpting.

This revolutionary new resin has just the right translucency and color variations found in natural tissue as well as ultra impact and flexural strength. The demand for function and aesthetics in dental fabrication is growing. Diamond D can help you take your business to a new level.



#### Liquid only

#### DIAMOND D SELF CURE MONOMER - KEYSTONE

Advanced natural-looking dentures require artistry as well as technical expertise.

Diamond D was specially formulated for exceptional working and handling characteristics for anatomic sculpting.

This revolutionary new resin has just the right translucency and color variations found in natural tissue as well as ultra impact and flexural strength. The demand for function and aesthetics in dental fabrication is growing. Diamond D can help you take your business to a new level.

**Note:**

The Liquid can only be shipped by ground transportation, it can not be shipped by air.  
Internationally maritime transportation is obligated. Minimum order is 700 us\$ + Hazardous Material Export license : 140 us\$



Code	Name	Price	+
34-1013017	KEYSTONE Diamond D Ultra Impact Acrylic, 8 ounce Self Cure Monomer - #1013017	US\$ 27.15	+Info
34-1013018	KEYSTONE Diamond D Ultra Impact Acrylic, 1 quart Self Cure Monomer - #1013018	US\$ 50.10	+Info

## POWDER ONLY

### 1LB

#### Dark Pink

Dark Pink

Code	Name	Price	+
34-1013109	Diamond-d SC Dark Ve 1lb Pwd	US\$ 67.88	<a href="#">+Info</a>

#### Light Redish Pnk

Light Redish Pnk

Code	Name	Price	+
34-1013058	Diamond-d SC red Pi 1lb Powder	US\$ 65.16	<a href="#">+Info</a>

#### Light Shade

Light Shade

Code	Name	Price	+
34-1013054	Diamond-d SC Light 1lb Powder	US\$ 67.88	<a href="#">+Info</a>

#### Original

Original

Code	Name	Price	+
34-1013050	Diamond-d SC Orig 1lb Powder	US\$ 65.16	<a href="#">+Info</a>

### 5LB

#### Dark Pink

Dark Pink

Code	Name	Price	+
34-1013110	Diamond-d SC Dark Ve 5lb Pwd	US\$ 254.32	<a href="#">+Info</a>

#### Light Pink

Light Pink

Code	Name	Price	+
34-1013055	Diamond-d SC Light 5lb Powder	US\$ 254.32	<a href="#">+Info</a>

#### Light Redish Pnk

Light Redish Pnk

Code	Name	Price	+
34-1013059	Diamond-d SC red Pi 5lb Powder	US\$ 254.32	<a href="#">+Info</a>

#### Original

Original

Code	Name	Price	+
34-1013051	Diamond-d SC Orig 5lb Powder	US\$ 254.32	<a href="#">+Info</a>

## 25LB

### Dark Pink

Dark Pink

Code	Name	Price	+
34-1013111	DiamonD D SC Dark Ve 25lb pwd	US\$ 776.69	<a href="#">+Info</a>

### Light Pink

Light Pink

Code	Name	Price	+
34-1013056	DiamonD D SC Light 25lb Pwd	US\$ 776.69	<a href="#">+Info</a>

### Light Redish Pnk

Light Redish Pnk

Code	Name	Price	+
34-1013060	DiamonD D SC red Pi 25lb Powde	US\$ 776.69	<a href="#">+Info</a>

### Original

Original

Code	Name	Price	+
34-1013052	DiamonD D SC Orig 25lb Powder	US\$ 660.19	<a href="#">+Info</a>

## WITH MONOMER

### 1lb with Monomer

#### Light Redish Pnk

Light Redish Pnk

Code	Name	Price	+
34-1013075	Diamond-d SC Lig red 1lb P&L	US\$ 73.08	<a href="#">+Info</a>

#### Light Shade

Light Shade

Code	Name	Price	+
34-1013074	Diamond-d SC Light 1lb P&L	US\$ 73.08	<a href="#">+Info</a>

#### Original

Original

Code	Name	Price	+
34-1013073	Diamond-d SC Orig 1lb P&L	US\$ 73.08	<a href="#">+Info</a>

### 5lb with Monomer

#### Light Pink

Light Pink

Code	Name	Price	+
------	------	-------	---

34-1013078	Diamond-d SC Light 5lb P&L	US\$ 257.54	<a href="#">+Info</a>
<b>Light Redish Pnk</b> Light Redish Pnk			
Code	Name	Price	+
34-1013079	Diamond-d SC L Rd Pk 5lb P&L	US\$ 257.54	<a href="#">+Info</a>
<b>Original</b> Original			
Code	Name	Price	+
34-1013077	Diamond-d SC Orig 5lb P&L	US\$ 257.54	<a href="#">+Info</a>
<b>25lb with Monomer</b>			
<b>Light Pink</b> Light Pink			
Code	Name	Price	+
34-1013082	DiamonD D SC Lig 25lb P&L	US\$ 841.04	<a href="#">+Info</a>
<b>Light Redish Pnk</b> Light Redish Pnk			
Code	Name	Price	+
34-1013083	DiamonD D SC L Rd Pk 25lb P&L	US\$ 841.04	<a href="#">+Info</a>
<b>Original</b> Original			
Code	Name	Price	+
34-1013081	DiamonD D SC Ori 25lb P&L	US\$ 841.04	<a href="#">+Info</a>

**Acrylic Self Cured**  
**Self Cross Linked Acrylics**



**DUZ-ALL SELF CURE REPAIR ACRYLIC | KEYSTONE**



**Powder only**  
Powder only

**Clear**  
Clear

Code	Name	Price	+
445891	Duz-All Clear Standard Kit. All-Pur	US\$ 66.34	<a href="#">+Info</a>
446110	Duz-All 8 oz. Clear Powder. All-Pur	US\$ 49.80	<a href="#">+Info</a>
446367	Duz-All 1 lb. Clear Powder. All-Pur	US\$ 88.06	<a href="#">+Info</a>

## Pink

Pink

Code	Name	Price	+
446514	Duz-All 8 oz. Pink Powder. All-Purp	US\$ 49.80	<a href="#">+Info</a>
446850	Duz-All 1 lb. Pink Powder. All-Purp	US\$ 88.06	<a href="#">+Info</a>

## Pink Veined

Pink Veined

Code	Name	Price	+
446030	Duz-All 1 lb. White Powder. All-Pur	US\$ 75.73	<a href="#">+Info</a>
446446	Duz-All 8 oz. Vein Powder. All-Purp	US\$ 49.80	<a href="#">+Info</a>

## Liquid only

Liquid only

Code	Name	Price	+
446074	Duz-All 4 oz. Liquid. All Purpose S	US\$ 41.69	<a href="#">+Info</a>
446404	Duz-All 8 oz. Liquid. All Purpose S	US\$ 59.43	<a href="#">+Info</a>
446436	Duz-All 32 oz. Liquid. All Purpose	US\$ 102.27	<a href="#">+Info</a>

## Kits / Set

Kits / Set

## Pink

Pink

Code	Name	Price	+
446139	Duz-All Pink Standard Kit. All-Purp	US\$ 71.33	<a href="#">+Info</a>

## Pink Veined

Pink Veined

Code	Name	Price	+
446363	Duz-All Vein Standard Kit. All-Purp	US\$ 71.33	<a href="#">+Info</a>

## OPTI-CRYL SELF CURE ACRYLIC C | NEW STETIC

OPTICRYL DENTURE ACRYLIC RESIN - NEW STETIC

**Features:**

- \* Color stable Denture Base resin with excellent working conditions
- \* This Denture Acrylic resin will be compatible with all regular Denture Resin Acrylics in the market and can replace the great majority of them

The proportions of use for the mixture are:

By weight: Two parts of self curing polymer Opti-Cryl and one part of self curing monomer Opti-Cryl.

By volume: Three parts of self curing polymer Opti-Cryl and one part of self curing monomer Opti-Cryl.



### 1lb Powder Only

#### Clear

Clear

Code	Name	Price	+
105124	Self Cure Clear 1lb Pwd.Only	US\$ 18.04	<a href="#">+Info</a>

#### Dark Pink

Dark Pink

Code	Name	Price	+
105119	Self Cure Dk Pink 1lb Pwd.Only	US\$ 18.04	<a href="#">+Info</a>

#### Light Pink

Light Pink

Code	Name	Price	+
105120	Self Cure Lt Pink 1lb Pwd.Only	US\$ 18.04	<a href="#">+Info</a>

#### Meharry

Meharry

Code	Name	Price	+
105125	Self Cure Meharry 1lb Pwd.Only	US\$ 18.04	<a href="#">+Info</a>

#### Original

Original

Code	Name	Price	+
105121	Self Cure Orig. 1lb Pwd.Only	US\$ 18.04	<a href="#">+Info</a>

## 1lb with Monomer

### Clear

Clear

Code	Name	Price	+
100400	Self Cure Clear 1lb+8oz Mon	US\$ 28.76	<a href="#">+Info</a>

### Dark Pink

Dark Pink

Code	Name	Price	+
100406	Self Cure Dk Pink 1lb+8oz Mon	US\$ 28.76	<a href="#">+Info</a>

### Light Pink

Light Pink

Code	Name	Price	+
100428	Self Cure Lt Pink 1lb+8oz Mon	US\$ 28.76	<a href="#">+Info</a>

### Meharry

Meharry

Code	Name	Price	+
100399	Self Cure Meharry 1lb+ 8oz Mon	US\$ 28.76	<a href="#">+Info</a>

### Original

Original

Code	Name	Price	+
100398	Self Cure Original 1lb+8oz Mon	US\$ 28.76	<a href="#">+Info</a>

## 5lb Powder Only

### Clear

Clear

Code	Name	Price	+
105126	Self Cure Clear 5lb Pwd.Only	US\$ 60.80	<a href="#">+Info</a>

### Dark Pink

Dark Pink

Code	Name	Price	+
105122	Self Cure Dk Pink 5lb Pwd.Only	US\$ 68.40	<a href="#">+Info</a>

### Light Pink

Light Pink

Code	Name	Price	+
105127	Self Cure Lt Pink 5lb Pwd.Only	US\$ 68.40	<a href="#">+Info</a>

### Meharry

Meharry

Code	Name	Price	+
105128	Self Cure Meharry 5lb Pwd.Only	US\$ 68.40	<a href="#">+Info</a>

**Original**  
Original

Code	Name	Price	+
105123	Self Cure Orig. 5lb Pwd.Only	US\$ 68.40	<a href="#">+Info</a>

**5lb with Monomer**

**Clear**  
Clear

Code	Name	Price	+
100410	Self Cure Clear 5lb+1Qt Mon	US\$ 95.00	<a href="#">+Info</a>

**Dark Pink**  
Dark Pink

Code	Name	Price	+
100415	Self Cure Dk Pink 5lb+1Qt Mon	US\$ 102.60	<a href="#">+Info</a>

**Light Pink**  
Light Pink

Code	Name	Price	+
100411	Self Cure Lt Pink 5lb+1Qt Mon	US\$ 102.60	<a href="#">+Info</a>

**Meharry**  
Meharry

Code	Name	Price	+
100409	Self Cure Meharry 5lb+1Qt Mon	US\$ 102.60	<a href="#">+Info</a>

**Original**  
Original

Code	Name	Price	+
100408	Self Cure Original 5lb+1Qt Mon	US\$ 102.60	<a href="#">+Info</a>

**22lb Powder Only**

**Clear**  
Clear

Code	Name	Price	+
105114	Self Cure Clear 22lb Pwd.Only	US\$ 194.25	<a href="#">+Info</a>

**Dark Pink**  
Dark Pink

Code	Name	Price	+
105115	Self Cure Dk Pink 22lb Pwd.Onl	US\$ 194.25	<a href="#">+Info</a>



## Light Pink

Light Pink

Code	Name	Price	+
105116	Self Cure Lt Pink 22lb Pwd.Onl	US\$ 218.75	<a href="#">+Info</a>

## Meharry

Meharry

Code	Name	Price	+
105117	Self Cure Meharry 22lb Pwd.Onl	US\$ 194.25	<a href="#">+Info</a>

## Original

Original

Code	Name	Price	+
105118	Self Cure Orig. 22lb Pwd.Only	US\$ 194.25	<a href="#">+Info</a>

## 22lb with Monomer

### Clear

Clear

Code	Name	Price	+
100429	Self Cure Clear 22lb+4Qt Mon	US\$ 320.25	<a href="#">+Info</a>

### Dark Pink

Dark Pink

Code	Name	Price	+
100426	Self Cure Dk Pink 22lb+ 4Qt Mo	US\$ 320.25	<a href="#">+Info</a>

### Light Pink

Light Pink

Code	Name	Price	+
100401	Self Cure Lt Pink 22lb+4Qt Mon	US\$ 344.75	<a href="#">+Info</a>

### Meharry

Meharry

Code	Name	Price	+
100424	Self Cure Meharry 22lb+4Qt Mon	US\$ 320.25	<a href="#">+Info</a>

### Original

Original

Code	Name	Price	+
100423	Self Cure Original 22lb+4Qt M	US\$ 320.25	<a href="#">+Info</a>

## Only Liquid

### SELF CURE MONOMER - OPTICRYL

#### Features:

This self-cure acrylic resin -polymer and monomer- is used for repairing dental prosthesis and the making of orthodontic and orthopedic appliances.

The total amount of time required for repairing different acrylic structures or for making orthodontic appliances is 10 minutes approximately, with an average work time of 4 minutes.

This resin is highly resistant to fractures and does not require heat treatment for its polymerization.

#### Note:

The Liquid can only be shipped by ground transportation, it can not be shipped by air. Internationally maritime transportation is obligated. Minimum order is 700 us\$ + Hazardous Material Export license : 140 us\$



## Denture Opti-Cryl

### Denture Opti-Cryl

Code	Name	Price	+
100247	Monomer Self Cure 8oz	US\$ 11.68	<a href="#">+Info</a>
100248	Monomer Self Cure 32oz	US\$ 38.24	<a href="#">+Info</a>
100248-4	Monomer Self Cure 32oz x 4	US\$ 152.93	<a href="#">+Info</a>

## MILLENNIUM Self cure Acrylic | KEYSTONE

### Millennium Self cure Acrylic - KEYSTONE



Code	Name	Price	+
445985	Millennium Acrylic - Clear Powder 5	US\$ 137.36	<a href="#">+Info</a>
446831	Millennium Acrylic - Deep Veined Po	US\$ 137.36	<a href="#">+Info</a>
446835	Millennium Acrylic - Light Veined P	US\$ 123.62	<a href="#">+Info</a>

## Self Cross Linked Acrylics - COLTENE

Self Cross Linked Acrylics - COLTENE

Code	Name	Price	+
433488	Hygenic Cold Cure Resin, Cures with	US\$ 0.00	+Info
433501	Hygenic Cold Cure Resin, Cures with	US\$ 0.00	+Info
433605	Hygenic Cold Cure Resin, Cures with	US\$ 0.00	+Info
436304	Hygenic Cold Cure Resin Liquid, 1 q	US\$ 0.00	+Info
436318	Hygenic Cold Cure Resin, Cures with	US\$ 0.00	+Info
436321	Hygenic Cold Cure Resin, Cures with	US\$ 0.00	+Info

## SLEDGEHAMMER SELF DENTURE ACRYLIC | KEYSTONE

**SLEDGEHAMMER SELF CURE ACRYLIC - KEYSTONE**

Sledgehammer Multi-Pack Self Cure Acrylic is color stable and high impact resistance to prevent breakage. It has overnight and short cure flexibility with 12-14 minute working times enables more units to be mixed simultaneously. There's no bleaching, blanching or streaking.

FDA Compliant Quality System. Manufactured under strict ISO and CGMP Guidelines for Guaranteed Quality and Product Reliability from batch to batch.

**Note:**

The Liquid can only be shipped by ground transportation; it cannot be shipped by air.

Internationally maritime transportation is obligated. Minimum order is 700 us\$ + Hazardous Material Export license: 140 us\$



### Powder only

Powder only

### Clear

Clear

Code	Name	Price	+
34-1000513	Self Cure Clear 1lb Pwd only	US\$ 36.68	+Info
34-1000514	Self Cure Clear 5lbs Pwd only	US\$ 93.94	+Info
34-1000515	Self Cure Clear 25lbs Pwd only	US\$ 407.49	+Info

### Light Fiber Plus

Light Fiber Plus

Code	Name	Price	
34-1001970	Self Cure Lig 1lb Pwd only	US\$ 36.68	<a href="#">+Info</a>
34-1001971	Self Cure Lig Fib 5lbs Pwd only	US\$ 93.94	<a href="#">+Info</a>
34-1001972	Self Cure Li Fib 25lbs Pwd only	US\$ 346.37	<a href="#">+Info</a>

## Light Fibered

Light Fibered

Code	Name	Price	
34-1000547	Self Cure Lig Fib 1lb Pwd only	US\$ 36.68	<a href="#">+Info</a>
34-1000548	Self Cure Lig Fib 5lbs Pwd only	US\$ 93.94	<a href="#">+Info</a>
34-1000549	Self Cure Li Fib 25lbs Pwd only	US\$ 407.49	<a href="#">+Info</a>

## Light Red-Pnk

Light Red-Pnk

Code	Name	Price	
34-10001951	Self Cure Orig LRP 5lbs Pwd	US\$ 93.94	<a href="#">+Info</a>
34-1000507	Self Cure LRP 1lb Powder only	US\$ 36.68	<a href="#">+Info</a>
34-1000508	Self Cure LRP 5lbs Powder only	US\$ 93.94	<a href="#">+Info</a>
34-1000509	Self Cure LRP 25lbs Powder only	US\$ 407.49	<a href="#">+Info</a>
34-1001950	Self Cure Ori LRP 1lb Pwd only	US\$ 36.68	<a href="#">+Info</a>

## Original

Original

Code	Name	Price	
34-1000510	Self Cure Orig 1lb Powder only	US\$ 32.25	<a href="#">+Info</a>
34-1000511	Self Cure Orig 5lbs Powder only	US\$ 93.94	<a href="#">+Info</a>
34-1000512	Self Cure Orig 25lbs Pwd only	US\$ 407.49	<a href="#">+Info</a>
34-1001952	Self Cure Ori LRP 25lbs Pwd onl	US\$ 294.65	<a href="#">+Info</a>

## Liquid Only

Liquid Only

Code	Name	Price	
34-1000521	Monomer Self Cure 8oz Sledge	US\$ 15.43	<a href="#">+Info</a>
34-1000522	Monomer Self Cure 32oz Sledge	US\$ 47.48	<a href="#">+Info</a>

## Powder & Liquid

Powder & Liquid

### Clear

Clear

Code	Name	Price	
34-1000487	Self Cure Clear 5lbs+32oz P&L	US\$ 139.51	<a href="#">+Info</a>
34-1000490	Self Cure Clear 25lbs+1gal P&L	US\$ 452.27	<a href="#">+Info</a>

## Light Fiber Plus

Light Fiber Plus

Code	Name	Price	
34-1001973	Self Cure Lig Fib 1lb+8oz P&L	US\$ 44.33	<a href="#">+Info</a>
34-1001974	Self Cure Lig Fib 5lbs+32oz P&L	US\$ 139.51	<a href="#">+Info</a>

34-1001975 Self Cure Lig Fi 25lbs+1gal P&L

US\$ 384.43 [+Info](#)

## Light Fibered Light Fibered

Code	Name	Price	+
34-1000539	Self Cure Lig.Fib 1lb+8oz P&L	US\$ 44.03	<a href="#">+Info</a>
34-1000540	Self Cure Lig Fib 5lbs+32oz P&L	US\$ 139.51	<a href="#">+Info</a>
34-1000541	Self Cure Li Fib 25lbs+1gal P&L	US\$ 452.27	<a href="#">+Info</a>

## Light Red-Pnk Light Red-Pnk

Code	Name	Price	+
34-1000482	Self Cure LRP 1lb+8oz P&L	US\$ 44.03	<a href="#">+Info</a>
34-1000485	Self Cure LRP 5lbs+32oz P&L	US\$ 139.51	<a href="#">+Info</a>
34-1000488	Self Cure LRP 25lbs+1gal P&L	US\$ 452.27	<a href="#">+Info</a>

## Original Original

Code	Name	Price	+
34-1000483	Self Cure Original 1lb+8oz P&L	US\$ 44.03	<a href="#">+Info</a>
34-1000486	Self Cure Orig 5lbs+32oz mo P&L	US\$ 132.53	<a href="#">+Info</a>
34-1000489	Self Cure Orig 25lbs+1gal P&L	US\$ 452.27	<a href="#">+Info</a>
34-1001955	Self Cure Ori LRP 25lbs+1ga P&L	US\$ 330.07	<a href="#">+Info</a>

## TruRepair Self Cured - KEYSTONE | KEYSTONE

TruRepair Self Cured - KEYSTONE - KEYSTONE



## Kits / Set Kits / Set

Code	Name	Price	+
445906	TruRepair Standard Kit - Fast Bondi	US\$ 65.38	<a href="#">+Info</a>

## Liquid only Liquid only

Code	Name	Price	+
446188	TruRepair Liquid Refill - Fast Bond	US\$ 25.73	<a href="#">+Info</a>

## Powder only

Powder only

### Clear

Clear

Code	Name	Price	+
446469	TruRepair Powder Refill, Clear - Fa	US\$ 26.55	<a href="#">+Info</a>

### Pink

Pink

Code	Name	Price	+
446263	TruRepair Powder Refill, Pink - Fas	US\$ 26.55	<a href="#">+Info</a>

### Pink Veined

Pink Veined

Code	Name	Price	+
446418	TruRepair Powder Refill, Vein - Fas	US\$ 26.55	<a href="#">+Info</a>

## UNIFAST Self Cure Resin | GC AMERICA

GC AMERICA Unifast Trad - General Purpose Acrylic Resin Trad Liquid only 104ml  
#339291



- For pressure and brush techniques
- Quick setting
- High-dimensional stability
- Suitable as long-term temporary

Item Number Description  
 339291 UNIFAST Trad Liquid (100 g)  
 339292 UNIFAST Trad Liquid (250 g)  
 339102 UNIFAST Trad Powder, 3 Pink (Darker Pink, no fibers; 100 g)  
 339125 UNIFAST Trad Powder, 8 Live Pink (Pink, no fibers; 1 kg)  
 339105 UNIFAST Trad Powder, 8 Live Pink (Pink, no fibers; 100 g)  
 339115 UNIFAST Trad Powder, 8 Live Pink (Pink, no fibers; 250 g)  
 339124 UNIFAST Trad Powder, Ivory (1 kg)  
 339104 UNIFAST Trad Powder, Ivory (100 g)  
 339114 UNIFAST Trad Powder, Ivory (250 g)  
 339123 UNIFAST Trad Powder, Pink X (Orange Pink, no fibers; 1 kg)  
 339103 UNIFAST Trad Powder, Pink X (Orange Pink, no fibers; 100 g)  
 339113 UNIFAST Trad Powder, Pink X (Orange Pink, no fibers; 250 g)

### Liquid

Liquid

Code	Name	Price	+
421619	Unifast Trad Liquid 260ml	US\$ 67.99	<a href="#">+Info</a>

### Powder Only

Powder Only

### Dark Pink

Dark Pink

Code	Name	Price	+
421504	Unifast Trad Pwdr only #3 Pink 100g	US\$ 60.76	<a href="#">+Info</a>

## Live Pink

Live Pink

Code	Name	Price	+
421586	Unifast Trad Pwdr Live Pink 250g	US\$ 116.45	<a href="#">+Info</a>
421791	Unifast Trad Powder only, Live Pink	US\$ 60.76	<a href="#">+Info</a>
422392	Unifast Trad Powder only, Live Pink	US\$ 411.67	<a href="#">+Info</a>

## Pink X

Pink X

Code	Name	Price	+
421382	Unifast Trad Pwdr only Pink X 100gr	US\$ 60.76	<a href="#">+Info</a>
421389	Unifast Trad Pwdr only Pink X 250gr	US\$ 116.45	<a href="#">+Info</a>
422377	Unifast Trad Powder only Pink X 1kg	US\$ 403.44	<a href="#">+Info</a>

## Catalog Index / Indice

Diamond -D .....	1
Self Cross Linked Acrylics .....	4