

Adhesives & Glues

Cyanoacrylate Remover & Blocke



ACCUSOLV BOTTLE - ADS

[ACCLUSOLV BOTTLE - ADS](#)

Features:

- * Accusolv dissolves Cyanoacrylates cements
- * Bonded materials can be separated in seconds
- * Clean off instruments, counters and fingertips
- * Is a must for repositioning of bonds and removal of excess cement
- * Use to dissolve Acculite and Fill-It cements



US DENTAL DEPOT is an Authorized American Dental Supply dealer in US

Code	Name	Price	+
117622	ADS Accusolv qt. - #A743-6	US\$ 51.59	+Info
117621	ADS Accusolv pt 16oz. - #A743-5	US\$ 28.88	+Info
117620	ADS Accusolv Bottle 4oz. - #A743-2	US\$ 10.47	+Info

Cyanoacrylate Remover & Blocke - KEYSTONE

Cyanoacrylate Remover & Blocke - KEYSTONE

Code	Name	Price	+
034-1020100	Easy Removal Of Adhesive 4 oz	US\$ 19.62	+Info

GLUE BLOCKER - ADS

[GLUE BLOCKER - ADS](#)

Features:

- * Guards yourself against the effects of Cyanoacrylates, epoxy, acrylics, urethanes, rubber based materials, light cure resins, dental cements. Sealers, coatings, die materials and bonding agents
- * Protects skin from reactions to irritating chemicals
- * Reduces chapping, cracking and rough, dry skin
- * Invisible layer repels oily materials, cream contains penetrating emollients for smooth, soft, silky hands
- * Will not take the place of rubber gloves
- * Protects the skin from damage
- * Prevents Cyanoacrylates and composites from sticking to and building up on models, instruments and counter tops
- * Keeps your equipment looking like new and facilitates quick clean-ups



- * Pure, white, non-contaminating
- * Not a greasy silicone
- * Can be washed off with water

Code	Name	Price	+
118192	ADS Glue Blocker kit - #G910-1	US\$ 29.19	+Info
118193	ADS Glue Blocker HAND CREAM 8oz. - #G910-2	US\$ 24.33	+Info
118070	ADS Glue Blocker COMPOUND 2oz. - #G910-3	US\$ 13.59	+Info

Adhesives & Glues

Cyanoacrylates



ACCULITE - ADS

ACCULITE - ADS

Features:

- * Just a drop of Acculite will rapidly join dental materials for instant office repairs and constructions with a bond of un surpassed strength and durability
- * The pinpoint applicator lets you see and place Acculite exactly where you want it
- * Secure plastic articulators, mount models with just a few drops
- * Tack wire components together before soldering or welding
- * Repair broken models and impression trays
- * Accurate placement, less waste
- * Secure dowel pins in place
- * Join rubber, metals, woods, glass and most plastics in any combination
- * Surface hardener for all stone dies
- * Works very fast, lowest viscosity, initial setting time 5 seconds



US DENTAL DEPOT is an Authorized American Dental Supply dealer in US

FAST SET

Code	Name	Price	+
117609	ADS Acculite Fast Sct With Bottles 1lb. - #A745-4 - Cyanoacrylates Cement	US\$ 125.33	+Info
117607	ADS Acculite Fast Sct 1 oz. - #A745-2 - Cyanoacrylates Cement	US\$ 13.02	+Info
117606	ADS Acculite Fast Sct 0.5oz. - #A745-1 - Cyanoacrylates Cement	US\$ 10.10	+Info

GAP FILL

Code	Name	Price	+
117617	ADS Acculite Kit Gap Fill With Bottles 1lb. - #A747-5 - Cyanoacrylates Cement	US\$ 125.33	+Info
117616	ADS Acculite Gap Fill 2oz. (2 btls of 1oz) - #A747-3 - Cyanoacrylates Cement	US\$ 21.69	+Info
117615	ADS Acculite Gap Fill 1 oz. - #A747-2 - Cyanoacrylates Cement	US\$ 13.02	+Info
117614	ADS Acculite Gap Fill 0.5oz. - #A747-1 - Cyanoacrylates Cement	US\$ 10.10	+Info

KITS

Code	Name	Price	+
117618	ADS Acculite Kit 1 3 - 1/2oz. Btls. kit - #A744-1 - Cyanoacrylates Cement	US\$ 19.75	+Info
117619	ADS Acculite Kit 3 3 - 1/2oz. Btls. Plus Accusolv kit - #A744-3 - Cyanoacrylates...	US\$ 25.78	+Info

REGULAR

Code	Name	Price	+
117613	ADS Acculite Regular With Bottles 1lb. - #A746-4 - Cyanoacrylates Cement	US\$ 125.33	+Info
117612	ADS Acculite Regular 2oz. (2 - 1oz. btis) - #A746-3 - Cyanoacrylates Cement	US\$ 21.69	+Info
117610	ADS Acculite Regular 0.5oz. - #A746-1 - Cyanoacrylates Cement	US\$ 10.10	+Info
117611	ADS Acculite Regular 1 oz. - #A746-2 - Cyanoacrylates Cement	US\$ 13.02	+Info

Cyanoacrylates - BUFFALO

Cyanoacrylates - BUFFALO

Code	Name	Price	+
424655	Permabond 910 Cyanoacryl Adhesive	US\$ 8.93	+Info
425336	Permabond 910 Cyanoacrylate Adhesiv	US\$ 24.39	+Info

CYANOACRYLATE ADHESIVE - BESQUAL

SUPER CYANOACRYLATE ADHESIVE - BESQUAL

- * Formulated to quickly bond plastics, metal, and ceramics.
- * It offers excellent strength characteristics

Adhesive & Activators:



Code	Name	Price	+
102427	BESQUAL Adhesive (Cyanoacrylates) - Gel, industrial type for broken metals, 2 -...	US\$ 9.98	+Info
102424	BESQUAL Adhesive (Cyanoacrylates) - Thin for Plastic Articulators, 2 oz. Bottle...	US\$ 10.40	+Info
102425	BESQUAL Adhesive (Cyanoacrylates) - Medium for dowel pins, 2 oz. Bottle - #508-202	US\$ 10.40	+Info
102426	BESQUAL Cyanoacrylate Thick Adhesive for Broken Metals - 2oz - 1 Bottle #508-203	US\$ 10.40	+Info
102428	BESQUAL Activator Driers for Adhesive 1 Bottle 8oz - #509-008 - Dries Instantly ...	US\$ 15.20	+Info

FILL IT REFILL - ADS

FILL IT REFILL - ADS

Features:

- * Is a versatile gap filing cyanoacrylate available in high visibility colors
- * Repairs cracks, chips and bubbles in any material
- * Apply and then spray with Fill-It Activator for an instant cure allowing immediate handling
- * Increased gap filling ability means a new freedom to fill in hollow areas and surfaces that do not precisely match
- * Bond across gaps and joints
- * Tack wires and components into position before soldering or welding
- * Faster than acrylics or waxes
- * Bond porous and dissimilar materials together by spraying one surface with Fill-It Activator and applying Fill-It Base to the other surface, then joining
- * Protect fragile margins by painting with a thin cyanoacrylate base then spraying with Fill-It Activator creating a tough durable smooth finish



ACTIVATOR

Code	Name	Price	+
117981	ADS Fill It Activator refill qt btl. - #F680-6	US\$ 57.21	+Info
117980	ADS Fill It Activator refill 8oz. btl. - #F680-4	US\$ 18.04	+Info
117979	ADS Fill It Activator refill 4oz. btl. - #F680-3	US\$ 14.15	+Info
117978	ADS Fill It Activator refill 4oz. spray - #F680-3	US\$ 12.18	+Info
117977	ADS Fill It Activator refill 2oz. spray - #F680-0	US\$ 9.74	+Info

BASE & SOLVENT KITS

BLUE

Code	Name	Price	+
117958	ADS Fill - 1T BLUE 4 - 2oz. Act. - 4 - base & Solvent kit large lab - #F680-10-BE	US\$ 81.62	+Info
117956	ADS Fill - It BLUE 2oz. Act. - 1 Base & Solvent kit - #F680-1-BE	US\$ 23.37	+Info
117957	ADS Fill - It BLUE 1 oz. Act & 0.5oz. Base kit - #F680-11-BE	US\$ 13.56	+Info

CLEAR

Code	Name	Price	+
117972	ADS Fill - It CLEAR 4 - 2oz. Act. - 4 - base & Solvent kit large lab - #F680-10-CR	US\$ 81.62	+Info
117970	ADS Fill - It CLEAR 2oz. Act. - 1 Base & Sol. kit - #F680-1-CR	US\$ 23.37	+Info
117971	ADS Fill - It CLEAR 1 oz. Act. - 1/2oz. Base kit - #F680-11-CR	US\$ 13.56	+Info

ROSE

Code	Name	Price	+
117965	ADS Fill It Rose 4 - 2oz. Act. - 4 - base & Sol. kit large lab - #F680-10-RE	US\$ 81.62	+Info
117963	ADS Fill - It Rose 2oz. Act. - 1 Base & Sol. kit - #F680-1-RE	US\$ 23.37	+Info
117964	ADS Fill - 1T Rose 1 oz. Act. - 1/2oz. Base kit - #F680-11-RE	US\$ 13.56	+Info

BASE REFILL

BLUE

Code	Name	Price	+
117962	ADS Fill It BLUE Base bulk BLUE 1lb - #F680-8-BE	US\$ 155.41	+Info
117961	ADS Fill It BLUE Base refill BLUE 2oz. - #F680-9-BE	US\$ 26.69	+Info
117960	ADS Fill It BLUE Base refill BLUE 1 oz. - #F680-5-BE	US\$ 15.59	+Info
117959	ADS Fill It BLUE Base refill BLUE 1/2oz. - #F680-2-BE	US\$ 12.44	+Info

CLEAR

Code	Name	Price	+
117976	ADS Fill - It CLEAR Base refill CLEAR 1lb - #F680-8-CR	US\$ 155.41	+Info
117975	ADS Fill - It CLEAR Base refill CLEAR 2oz. - #F680-9-CR	US\$ 24.02	+Info
117974	ADS Fill - It CLEAR Base refill CLEAR 1 oz. - #F680-5-CR	US\$ 15.59	+Info
117973	ADS Fill - It CLEAR Base refill CLEAR 0.5oz. - #F680-2-CR	US\$ 12.44	+Info

ROSE

Code	Name	Price	+
117969	ADS Fill - It Rose Base refill rose 1lb - #F680-8-RE	US\$ 155.41	+Info
117968	ADS Fill - It Rose Base refill rose 2oz. - #F680-9-RE	US\$ 26.69	+Info
117967	ADS Fill - It Rose Base refill rose 1 oz. - #F680-5-RE	US\$ 15.59	+Info
117966	ADS Fill - It Rose Base refill rose 1/2oz. - #F680-2-RE	US\$ 12.44	+Info

SOLVENT REFILL

Code	Name	Price	+
117982	ADS Fill It Solvent refill 1 oz. - #F680-7	US\$ 5.99	+Info

TACKY LIQUID REFILL - ADS

ADS - TACKY LIQUID REFILL

Features:

- * For use as a pattern adhesive for luting wax and plastic patterns to models.
- * Can be used with impression trays to secure alginate impression materials.
- * May also be used on dies to secure retention beads before casting.
- * Non-flammable.



Code	Name	Price	+
118409	ADS Tacky Liquid refill pt 16oz. - #T 770-5	US\$ 50.98	+Info
118408	ADS Tacky Liquid refill 4oz. - #T 240-5	US\$ 23.08	+Info

118407	ADS Tacky Liquid 2 - 1 oz. btls - #T240-1	US\$ 28.33	+Info
--------	---	------------	-----------------------

Adhesives & Glues
Epoxy Glues



Glues - JSP

Glues - JSP

Code	Name	Price	+
461664	Epoxy 330 2 Tubes 1/2oz	US\$ 6.75	+Info

Catalog Index / Indice

Cyanoacrylate Remover & Blocker	1
Cyanoacrylates	2
Epoxy Glues	6