

Articulators and Occlusors

Accessories for articulators

Accessories for articulators - AD2 DENTAL

Accessories for articulators - AD2 DENTAL

Code	Name	Price	+
469671	Artex Adjustable Platform	US\$ 180.60	+Info
469658	EZ Bow Assembly	US\$ 236.60	+Info
469659	AD2/Panadent EZ Bow System	US\$ 390.60	+Info
469660	Stratos EZ Bow System	US\$ 390.60	+Info
469661	Artex EZ Bow System	US\$ 390.60	+Info
469662	SAM EZ Bow System	US\$ 390.60	+Info
469663	KaVo EZ Bow System	US\$ 390.60	+Info
469664	Denar Mark 300 Series EZ Bow System	US\$ 390.60	+Info
469665	Hanau Wide-Vue EZ Bow System	US\$ 390.60	+Info
469666	Whip Mix Quick Mount EZ Bow System	US\$ 390.60	+Info
469667	Denar Mark 2 EZ Bow System	US\$ 390.60	+Info
469668	AD2/Panadent Adjustable Platform	US\$ 180.60	+Info
469669	AD2/Panadent Fixed Platform	US\$ 138.60	+Info
469670	Stratos Adjustable Platform	US\$ 180.60	+Info
469657	Index Trays (bag of 50 trays)	US\$ 47.60	+Info
469672	SAM Adjustable Platform	US\$ 180.60	+Info
469673	KaVo Adjustable Platform	US\$ 180.60	+Info
469674	Denar Mark 300 Series Adjustable Pl	US\$ 180.60	+Info
469675	Hanau Wide-Vue Adjustable Platform	US\$ 180.60	+Info
469676	Whip Mix Quick Mount Adjustable Pla	US\$ 180.60	+Info
469677	Denar Mark 2	US\$ 180.60	+Info
469679	Mounting Table Assembly	US\$ 49.00	+Info
469680	Flat Bite Fork Assembly (includes b	US\$ 77.00	+Info
469681	Flat Bite Fork Only	US\$ 26.60	+Info
469682	Curved Bite Fork Only	US\$ 26.60	+Info
469683	Curved Bite Fork Assembly (includes	US\$ 77.00	+Info
469684	Bite Fork Stem Assembly	US\$ 63.00	+Info
474377	Artex Compatible Magnetic Mounting	US\$ 77.00	+Info
469644	AD2 Magnetic Mounting System (inclu	US\$ 138.60	+Info
469631	DynaLink	US\$ 0.49	+Info
469632	Straight Incisal Pin Module	US\$ 96.60	+Info
469633	Curved Incisal Pin Module	US\$ 152.60	+Info
469634	Incisal Table Assembly	US\$ 11.20	+Info
469635	Bite Fork Support	US\$ 96.60	+Info
469636	Articulator Support Leg	US\$ 12.60	+Info
469637	Cruzer Articulator with Mounting St	US\$ 773.50	+Info
469638	Cruzer Articulator (no mounting sta	US\$ 643.50	+Info
469639	Accu-Trays (box of 25)	US\$ 44.80	+Info
469640	Aesthetic Ruler	US\$ 5.60	+Info
469641	Magnet (neodymium)	US\$ 5.60	+Info

469642	AD2/Panadent/Whip Mix Compatible Sc	US\$ 70.00	+Info
469643	Denar/Hanau (Twin pin) Compatible S	US\$ 70.00	+Info
469630	Articulator Mounting Stand	US\$ 68.60	+Info
469645	Magnetic Mounting Base Assembly	US\$ 54.60	+Info
469646	Magnet Holder	US\$ 11.20	+Info
469647	AD2/Panadent Compatible Magnetic Mo	US\$ 77.00	+Info
469648	Panadent Compatible Magnetic Mounti	US\$ 88.20	+Info
469649	Whip Mix/Denar/Hanau Compatible Mag	US\$ 77.00	+Info
469650	Stratos Compatible Magnetic Mountin	US\$ 77.00	+Info
469651	SAM Compatible Magnetic Mounting PI	US\$ 77.00	+Info
469652	SAM Compatible Magnetic Mounting PI	US\$ 88.20	+Info
469653	Kavo Compatible Magnetic Mounting P	US\$ 77.00	+Info
469654	Denar Mark 300 Series Complatable P	US\$ 77.00	+Info
469655	Whip Mix/SAM QuickMount Plates (bag	US\$ 77.00	+Info
469656	Artex Compatible Magnetic Mounting	US\$ 77.00	+Info

Accessories for articulators - AMANN GIRRBACH

Accessories for articulators - AMANN GIRRBACH

Code	Name	Price	+
468012	Artex Arcon clip	US\$ 39.20	+Info
467932	Artex head support	US\$ 15.94	+Info
467933	Artex O-Ring for Set Support	US\$ 8.58	+Info
467936	Artex - semi axes NK/TK (pair)	US\$ 129.87	+Info
467968	Radial pin for Artex Carbon	US\$ 203.59	+Info
467970	Artex incisal table incl .218141+21	US\$ 58.94	+Info
467978	Micrometer pin Artex	US\$ 276.88	+Info
467979	Artex Top / Bottom Part	US\$ 207.05	+Info
467980	artex Cap for Head Support	US\$ 8.58	+Info
467981	Artex Allan key A	US\$ 56.35	+Info
467929	Screw for anterior guidance pin Art	US\$ 22.05	+Info
468014	Complete centric unit Artex type Ar	US\$ 274.43	+Info
468016	Artex CP/CR Condyle axes (pair)	US\$ 47.78	+Info
468017	Retrusion insert Artex CP	US\$ 79.03	+Info
468018	Axis and distractoin support Artex	US\$ 149.47	+Info
468019	Centric lever Artex CP/CR	US\$ 53.90	+Info
468022	Artex CR retraler Anschlag	US\$ 118.84	+Info
468024	Artex CR protrusion pointer (pair)	US\$ 302.61	+Info
468117	Artex separating spray	US\$ 49.55	+Info
467895	Screw for condyle setting Artex CN/	US\$ 26.95	+Info
467880	Artex-foot compl.	US\$ 25.74	+Info
467882	Artex Screw for Mode Plate	US\$ 29.40	+Info
467883	Artex O-Ring	US\$ 8.03	+Info
467887	Artex Anterior Guidance Screw	US\$ 25.74	+Info
467889	Artex Benett Setting Screw	US\$ 29.40	+Info
467890	Artex Condyle Setting Screw	US\$ 29.40	+Info
467892	Artex incisal needle	US\$ 22.05	+Info
467894	Artex CN condyles (pair)	US\$ 196.02	+Info
467874	Artex i-Tero model adapter	US\$ 586.84	+Info

467896	Artex CN/CT centric lock pair (set)	US\$ 29.40	+Info
467905	Artex C Anterior Guidance Screw Top	US\$ 25.74	+Info
467907	Upper part incl. plastic top Artex	US\$ 221.76	+Info
467909	Artex condylar ball MS	US\$ 22.05	+Info
467910	Artex Snap ring (10 pcs)	US\$ 22.05	+Info
467917	Artex click-axis	US\$ 8.58	+Info
467918	Artex click-ball	US\$ 8.58	+Info
467919	Artex click-spring	US\$ 8.58	+Info

Accessories for articulators - BIOART

Accessories for articulators - BIOART

Code	Name	Price	+
473868	Metallic Fork Replacement 1pk	US\$ 4.80	+Info
474948	Magnet disc for Mounting Plate	US\$ 4.32	+Info
475024	Screw X A7 Plus Articulator	US\$ 5.51	+Info

Accessories for articulators - KEYSTONE

Accessories for articulators - KEYSTONE

Code	Name	Price	+
446893	Magnet-X Base Former with Magnet	US\$ 14.60	+Info
466390	White Mounting Plate	US\$ 8.63	+Info

Accessories for articulators - OSUNG

Accessories for articulators - OSUNG

Code	Name	Price	+
449969	Osung Plane Plate, Occlusal, ARP1	US\$ 45.57	+Info

Accessories for articulators - PLASDENT

Accessories for articulators - PLASDENT

Code	Name	Price	+
423620	Plasdent Oblong Plastic Articulating	US\$ 67.79	+Info
424160	Plasdent Round Plastic Articulating	US\$ 79.75	+Info

Accessories for articulators - WHIP MIX

Accessories for articulators - WHIP MIX

Code	Name	Price	+
119954	Magnetic Converter Plates 2pk	US\$ 140.03	+Info
120059	Crossref Centerd Bite Record P	US\$ 240.24	+Info

ARTICULATORS ACCESORIES - BIOART

ARTICULATORS ACCESSORIES - BIOART

Features:

The new Elite Face-Bow brings an articulated and fixation system for the fork that allows a faster registration process for the professional and comfortable for the patient. It has a jig transfer assembly that eliminates the needs for face-bow in the assembling process of the casts in the articulator, turning the process faster, more practical and stable.

Its Naziun Relator with adjustable height, including an infra-orbital indicator, turns the registration process more precise. The face-bow is equipped with a sliding system which allows smooth lateral movement (opening/closing) and facilitates the overall use. The ear pieces are anatomical, removable and chemically sterilized.



Code	Name	Price	+
123175	BIOART Accessories Kit for Articulator Complete - Campers Plan - Fork support - ...	US\$ 59.25	+Info
123174	BIOART Adjustable Incisal Guiding Table (Metallic) each AMIM0069 (Accessories)	US\$ 32.59	+Info
124660	BIOART Campers Table Plane each APNC0446 (Accessories)	US\$ 11.70	+Info
122991	BIOART Rail Mounting Magnetic Plates units (10 pairs) - #APLP263	US\$ 61.05	+Info
416670	BIOART Elite Facebow toothless Fork Replacement 1pk - #AGCT0061(Accessories)	US\$ 12.16	+Info
417323	BIOART Edentulous Facebow Fork Replacement 1pk - #AGDE0495 (Accessories)	US\$ 4.75	+Info
416994	BIOART - FOX RULER PLANE (aluminum) each AJGS0544 (Accessories) - # ARFX0105	US\$ 27.79	+Info
416993	BIOART GARFO FORK- SUPPORT each ASUGO111 (Accessories)	US\$ 15.20	+Info
420787	BIOART A7 Plus Lock of upper frame - ORANGE piece #FTRA0765 (Accessories)	US\$ 2.00	+Info
461290	BIOART Transfer Set Support - SCTS0148	US\$ 25.28	+Info
466400	BIOART Broadrick Flag #ABAN0049	US\$ 48.64	+Info
466269	BIOART discs for the magnetic mounting plate(10 Pack) - #ADOM0066 (Accessories...	US\$ 7.50	+Info
466270	BIOART Magnetic Mounting PLATE ONLY EACH - #APMG0098 (Mounting Plates)	US\$ 30.32	+Info

ARTICULATORS ACCESSORIES ALL SERIES -WHIP MIX

ACCESSORIES ARTICULATORS - WHIP MIX



Features:

Whip Mix accessories enhance the capabilities of Whip Mix articulator systems and provide the clinician with a wide range of options for specialized needs.

ACCESSORIES 2000 SERIES

THE OTHERS ACCESSORIES FOR 2000 SERIES ARE IN: ACCESSORIES COMMON FOR ALL ARTICULATORS.

Code	Name	Price	+
119957	WHIP MIX Magnetic Disposable Mounting Plates (20 Pack) - #20002368 (Accessories...	US\$ 75.47	+Info
119956	WHIP MIX Magnetic Converter Plates - 2pk - #20002367 (Accessories Common for al...	US\$ 129.91	+Info
119955	WHIP MIX Magnetic Disposable Mounting Plates (100 Pack) - #20002369 (Accessorie...	US\$ 287.01	+Info
119872	WHIP MIX Quickmount Plates WHITE w/discs 100pk - #28710 (Accessories Common t...	US\$ 432.47	+Info

119855	WHIP MIX #8580b Plastic Mount Plates Box 500 - #28517 (Accessories for Articula...	US\$ 1,299.36	+Info
119854	WHIP MIX #8580b Plastic Mount Plates Box 100 - #28509 (Accessories common for a...	US\$ 272.54	+Info

ACCESSORIES 4000 SERIES

THE OTHERS ACCESSORIES FOR 2000 SERIES ARE IN: ACCESSORIES COMMON FOR ALL ARTICULATORS.

Code	Name	Price	+
119850	WHIP MIX #8780b Plastic Carrying Case - #18155 (Accessories common for all Arti...	US\$ 148.97	+Info
119853	WHIP MIX #8783 Universal Tote Bag - #18156 (Accessories common for all Articula...	US\$ 132.21	+Info

ACCESSORIES 8500 SERIES

Code	Name	Price	+
119862	WHIP MIX #8806b Spring Latch Assembly - #14524 (Accessories for 8500 Series Art...	US\$ 62.51	+Info

ACCESSORIES COMMON FOR ALL ARTICULATORS

THIS ACCESSORIES ARE COMMON FOR ALL TYPES ARTICULATORS

Code	Name	Price	+
119858	WHIP MIX #8526f Plastic Incisal Guide Flat - #04235 (Incisal Pins for Articulat...	US\$ 17.50	+Info
119859	WHIP MIX #8527 Screw For Plastic Incisal - #04243 (Incisal Pins for Articulators)	US\$ 7.97	+Info
119860	WHIP MIX #8528 Incisal Guide Pin Regular - #04251 (Incisal Pins for Articulators)	US\$ 42.80	+Info
119848	WHIP MIX #8533 Dovetail Incisal Block - #18870 (Accessories for Articulators)	US\$ 162.93	+Info
119929	WHIP MIX Hanau Adjustable Guide Table For Wm Articulator - #010631-000 (Accesso...	US\$ 252.88	+Info
120033	WHIP MIX Incisal Pin Sub Assy-p2- D5a - #110092 (Accessories for D5A Articulator)	US\$ 187.53	+Info
120027	WHIP MIX Lock Knob- Fossa- Labrelator (from Blank) - #20011061 (Accessories Com...	US\$ 22.90	+Info
119871	WHIP MIX Quickmount Dentoform Adapters 3 Pc - #28712 (Accessories common for al...	US\$ 363.23	+Info
120056	WHIP MIX Replacement Rubberbands-tmj - #650027 (Accessories for Articulators)	US\$ 10.48	+Info
120028	WHIP MIX Support Lockpin- Labrelator (from Blank) - #20011063 (Accessories Comm...	US\$ 22.90	+Info

DENAR - ARTICULATORS ACCESSORIES - WHIP MIX

DENAR - ARTICULATORS ACCESSORIES - WHIP MIX



Denar accessories enhance the capabilities of Denar articulator systems and provide the clinician with a wide range of options for specialized needs.
This accessories are common for all articulators.

ACCESSORIES COMMON FOR ALL ARTICULATORS

UNIVERSAL SPLIT CASE BASE ASSEMBLY

Code	Name	Price	+
119959	WHIP MIX Universal Split Cast Base Assembly 20 Mm - #20011745	US\$ 342.68	+Info
119960	WHIP MIX Universal Split Cast Base Assembly 10 Mm - #20011743	US\$ 302.82	+Info
119961	WHIP MIX Universal Split Cast Base Assembly 15 Mm - #20011744	US\$ 310.76	+Info

WRENCHS

Code	Name	Price	+
------	------	-------	---

120011	WHIP MIX Wrench .050 Hex Straight - #3261-2 (Accessories Common for Mark II Arti...	US\$ 26.22	+Info
120030	WHIP MIX Wrench- 1/16 Hex Straight - #3751	US\$ 24.79	+Info
120020	WHIP MIX Wrench 3/32 Hex Straight - #3261-1 (Accessories for all Articulators)	US\$ 33.35	+Info

OTHER ACCESSORIES

Code	Name	Price	+
119856	WHIP MIX Patient Record Cards pk 100 - #03891	US\$ 23.55	+Info
119897	WHIP MIX #8614 Hex Driver - #04545	US\$ 27.69	+Info
120022	WHIP MIX Fossa Insert-mkii - #300710-1 (Accessories for Mark II Articulators)	US\$ 37.39	+Info
119851	WHIP MIX #2225 Crvd Plate F/occlusal Std - #34213 (Accessories for Articulators)	US\$ 53.11	+Info
120057	WHIP MIX Replacement Inserts-tmj Tutor - #650026 (Accessories for Articulators)	US\$ 98.92	+Info

ADAPTABILITY

Code	Name	Price	+
120058	WHIP MIX Tmj Tutor Kit - #650000 (Accessories for Adaptability for Articulators)	US\$ 564.95	+Info
120055	WHIP MIX Best-bite Refill Kit - #20000012 (Accessories for Adaptability for Art...	US\$ 510.97	+Info
119843	WHIP MIX Fleximount pk 100 - #20011833 (Accessories for Adaptability for Artic...	US\$ 352.84	+Info
120053	WHIP MIX Best-bite Intro Kit - #20000011 (Accessories for Adaptability for Arti...	US\$ 273.43	+Info
120054	WHIP MIX Denar Deprogrammer (10 Pk) - #20000031 (Accessories for Adaptability f...	US\$ 95.18	+Info
119846	WHIP MIX Fleximount pk 20 - #20011832 (Accessories for Adaptability for Articul...	US\$ 70.71	+Info
120051	WHIP MIX Denar Bite Registration Wax - 25pk - #20000040 (Accessories for Adapta...	US\$ 67.85	+Info

D5A ARTICULATORS

Code	Name	Price	+
120012	WHIP MIX Medial Wall Insert-d5a-nyl-pair - #110065-4p (Accessories for D5A Arti...	US\$ 115.94	+Info
120023	WHIP MIX Super Wall Insert-d5a-nyl-strrt- Pair - #110066-1p (Accessories for D5A...	US\$ 79.03	+Info
120025	WHIP MIX Super Wall Insert-d5a-acr-3/8r- Pair - #110066-3p (Accessories for D5A...	US\$ 119.33	+Info

MARK 300 SERIES

Code	Name	Price	+
120026	WHIP MIX Labrelator Mark 300 Series - #20011380 (Accessories for Mark 300 Artic...	US\$ 1,036.50	+Info
120000	WHIP MIX Centri Check Mark 300 Series - #180036 (Accessories for Mark 300 Artic...	US\$ 805.93	+Info
119958	WHIP MIX Split Cast Plates 3 Piece Kit - #20011750 (Accessories for Mark 300 Ar...	US\$ 852.95	+Info
119953	WHIP MIX Mark 300 Magnetic Mounting Plate 100 Pack - #20011578 (Accessories for...	US\$ 258.31	+Info
120008	WHIP MIX Flex Tracks- Mark 310/320 - #20011395 (Accessories for Mark 300 Artic...	US\$ 291.27	+Info
120038	WHIP MIX Mark 300 Taper Incisal Pin w/support P6 - #20011397 - (Accessories for...	US\$ 185.78	+Info
120036	WHIP MIX Mark 300 Short Incisal Pin w/support-p4 - #20011386 - (Accessories for...	US\$ 185.78	+Info
119952	WHIP MIX Mark 300 Magnetic Mounting Plate 20 Pack - #20011577 (Accessories for ...	US\$ 75.47	+Info
120007	WHIP MIX Mark 300 Series Dust Cover - #20011392 (Accessories for Mark 300 Artic...	US\$ 49.69	+Info
120002	WHIP MIX Centri Check Stylus / Mark 300 Series - #180032 - (Accessories for Mar...	US\$ 126.49	+Info

VERIFICATION GAGES

Code	Name	Price	+
119998	WHIP MIX WHIP MIX Verification Gage - #20011670 (Accessories for Verification G...	US\$ 1,554.10	+Info
119999	WHIP MIX Mark 300 Series Verification Gage - #20011660 (Accessories common for ...	US\$ 1,554.05	+Info
120005	WHIP MIX Denar/hanau Verification Gage - #20011650 (Accessories for Verificatio...	US\$ 1,554.10	+Info

INCISAL PINS AND TABLES - WHIP MIX

INCISAL PINS AND TABLES - WHIP MIX

Features:

All incisal pins and tables are interchangeable between articulators. Incisal pins can be raised and lowered on an arc, simulating the natural arc of closure. The combinations shown must be used together on the articulator.



PINS

Code	Name	Price	+
120034	WHIP MIX Incisal Pin Assy-p2- D5a - #110093 (Accessories Common for D5A and Inc...	US\$ 271.06	+Info
120032	WHIP MIX Incisal Pin Assy-p6 - #110509 (Incisal Pins for Articulators)	US\$ 205.50	+Info
120037	WHIP MIX Denar Incisal Pin Assy-p4- BLACK - #112027 (Incisal Pins for Articulat...	US\$ 183.85	+Info
120031	WHIP MIX Incisal Pin-fillastre - #110448 (Incisal Pins for Articulators)	US\$ 120.63	+Info

TABLES

Code	Name	Price	+
120041	WHIP MIX T3 Mech Incisal Table D5a Assembly - #110109 (Accessories Common for D...	US\$ 549.78	+Info
120043	WHIP MIX Mechanical Table Assy- T6 - #110635 (Incisal Tables for Articulators)	US\$ 386.09	+Info
120042	WHIP MIX Incisal Table Flat T4 (12) - #110240-12 (Incisal Tables for Articulators)	US\$ 113.76	+Info
120040	WHIP MIX Incisal Table With Step-t2 - w/screw - #110241 (Incisal Tables for Art...	US\$ 28.32	+Info
120044	WHIP MIX Incisal Table Flat T4 - w/screw - #110240 (Incisal Tables for Articula...	US\$ 28.32	+Info

PARTS FOR ARTEX ARTICULATOR - AMANN GIRRBACH

Code	Name	Price	+
467830	Splitex key for Artex 116mm	US\$ 381.02	+Info
467831	Splitex key for Artex, 126mm, C-Ver	US\$ 381.03	+Info
467891	Artex Support Pin	US\$ 24.50	+Info
467898	Artex CT condyles (pair)	US\$ 476.58	+Info
467906	Artex plate 0 for C-version	US\$ 51.46	+Info
467911	NoPlast f. Artex Carbon	US\$ 763.50	+Info
467920	Artex click-locking bar (pair)	US\$ 22.05	+Info
467945	Artex Quickbite, 10 pcs.	US\$ 75.01	+Info
467963	Artex case	US\$ 195.56	+Info
467996	Artex Transfer Jig	US\$ 336.20	+Info

Articulators and Occlusors

Crown and Bridge Articulators

ARTICULATOR CROWN AND BRIDGE - BESQUAL

ARTICULATOR CROWN AND BRIDGE - BESQUAL

The articulator crown arch brass

The Besqual Crown and Bridge Articulator is specifically engineered for crown and bridgework.

Its easy to adjust and easy to separate parts.

These articulators are machined for quality and accuracy.



Code	Name	Price	+
100831	BESQUAL Articulator Crown and Bridge- Each -Brass- #603-100	US\$ 8.80	+Info
100832	BESQUAL Articulator Crown and Bridge- Each -Chrome #602-100	US\$ 5.95	+Info

Crown and Bridge Articulators - KEYSTONE

Crown and Bridge Articulators - KEYSTONE

Code	Name	Price	+
034-1050140	Art K-Chrome Crown&Bridge	US\$ 15.29	+Info

Crown and Bridge Articulators - RAY FOSTER

Crown and Bridge Articulators - RAY FOSTER

Code	Name	Price	+
470919	Cr & Bridge Articulators #57 Chrome	US\$ 25.20	+Info

Articulators and Occlusors

Denture Articulators



Apex (Denture) Articulators - KEYSTONE

Apex (Denture) Articulators - KEYSTONE



DELUXE ARTICULATOR

DELUXE ARTICULATOR

Code	Name	Price	+
034-1050095	Deluxe Articulator With Pin	US\$ 38.30	+Info

Other Models

Other Models

Code	Name	Price	+
446238	K-Chrome High Arch Articulator	US\$ 14.33	+Info
446113	Keystone Apex #2 Articulator, Chrom	US\$ 16.88	+Info
446154	Keystone Apex #2 Articulator, Bronz	US\$ 18.23	+Info

Apex (Denture) Articulators - RAY FOSTER

Apex (Denture) Articulators - RAY FOSTER

Code	Name	Price	+
470914	Apex #2 Articulator Chrome	US\$ 25.20	+Info
470915	Apex #2 Articulator Brass	US\$ 25.20	+Info
470916	Apex #8 Articulator Chrome	US\$ 24.19	+Info
470917	Apex #8 Articulator Brass	US\$ 25.20	+Info
470918	#3 Articulator w/ Anterior Pin	US\$ 78.40	+Info

Apex HIGH ARCH (Denture) Articulators - BESQUAL

High Arch articulators - BESQUAL



High Arch

High Arch

Code	Name	Price	+
470274	High Arch Denture Brass bronze	US\$ 10.64	+Info

Other Models

Other Models

Code	Name	Price	+
100827	High Arch Denture Chrome	US\$ 9.95	+Info

Apex LOW ARCH (Denture) Articulators - BESQUAL

ARTICULATOR APEX LOW ARCH DENTURE - BESQUAL

The articulator low arch bronze chrome

- * BESQUAL Articulator have been used by leading dental technicians and dentists for many years.
- * Constructed of bronze castings with a center support for durability and better plaster retention.
- * Machined to produce maximum accuracy for this type articulator, the Apex type articulator provides durable service.



Code	Name	Price	+
100830	BESQUAL Articulator (Apex) Low Arch Denture Each -Brass- #603-200 (Articulador ...	US\$ 10.64	+Info
100829	BESQUAL Articulator (Apex) Low Arch Denture Each -Chrome- (Articulador arco baj...	US\$ 10.32	+Info

Articulators and Occlusors

Disposable articulators



Disposable articulators - BESQUAL

Disposable articulators - BESQUAL

Code	Name	Price	+
100731	Disposable Art Black 100pk	US\$ 15.75	+Info
417391	Disposable Art Black Caps Only100pk	US\$ 14.32	+Info

Disposable articulators - DIGITECH

Disposable articulators - DIGITECH

Code	Name	Price	+
402380	3D ARTICULATOR plastic 100pk	US\$ 25.02	+Info

Disposable articulators - JSP

Disposable articulators - JSP

Code	Name	Price	+
437630	JSP Disposable Plastic Articulators	US\$ 15.66	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



**US
DENTAL
DEPOT**

+1 954 874 6325
www.usdentaldepot.com



MAINSTAY ARTICULATOR SYSTEM - WHIP MIX

MAINSTAY ARTICULATOR SYSTEM - WHIP MIX

The MainStay Disposable Articulator System was designed to save you time, money and waste. This system of disposable articulators does not require model formers or a gypsum base, and adjusts to accommodate any occlusal relationship. Its base was designed without external walls or raised edges to prevent interference with the seating of dies and model segments. This unique feature also means that it can easily be used with gypsums of all expansion ratios.

Choose between Pin Models in Full Arch, Anterior and Quadrant sizes. Pinless Models offered in a Quadrant size.

Features:

- * No External Walls- Flat table for easy removal and replacement of segments
- * Support Rods- Built-in to provide a vertical stop, ideal for edentulous models
- * Hinges- Durable and easy to separate. Resists breakage. Can be heated for custom occlusal relationship.

Pin System:

- * Stability of the MainStay? Pins
- * Numbered bases for quick, easy reference of die placement
- * No drilling- saves time
- * Raised anti-rotational spine prevents movement of die

Pinless System:

- * Excellent for scanning
- * No pinning, no drilling, no gluing- saves time
- * Built-In anti-rotational spine prevents movement of die
- * Numbered bases for quick, easy reference of die placement



PIN SYSTEM

Code	Name	Price	+
119899	WHIP MIX Mainstay Disposable Articulator Pin Quadrant 50 Sets - #08252 (Accesso...	US\$ 97.42	+Info
119902	WHIP MIX Mainstay Disposable Articulator Pin Anterior 25 Sets - #08253 (Accesso...	US\$ 59.54	+Info
119904	WHIP MIX Mainstay Disposable Articulator Pin Full Arch 25 Sets - #08251 (Access...	US\$ 59.54	+Info

PINLESS SYSTEM

Code	Name	Price	+
119901	WHIP MIX Mainstay Disposable Articulator Pinless Quadrant 50 Sets - #08254 (Acc...	US\$ 88.56	+Info

ACCESSORIES

Code	Name	Price	+
119903	WHIP MIX Mainstay Dowel Pins Only pk 1000 - #08256 (Incisal Pins for Articulat...	US\$ 98.25	+Info
119900	WHIP MIX Mainstay Stop Support Rods pk 20 - #08255 (Accessories for MainStay A...	US\$ 7.59	+Info

Articulators and Occlusors

Face bows

DENAR - SLIDEMATIC FACEBOWS

-WHIP MIX

DENAR - SLIDEMATIC FACEBOWS - WHIP MIX

Features:

Possibly the fastest facebow system you will ever use. A fast, easy and accurate means of transferring proper maxillary relationships.

Description:

- * Precision speed-glide gear mechanisms permits the facebow to open and close on the axis.
- * Convenient detachable transfer jig mounts to the articulator, freeing the facebow for the next patient.
- * Adaptable to most Denar, Hanau and Whip Mix Articulators.
- * Not compatible with Whip Mix 8500 or Hanau-Mate.
- * Comes complete with: one measuring bow, Denar transfer jig assembly, other needed accessories, instruction manual and carrying case.



SLIDEMATIC

Code	Name	Price	+
119989	WHIP MIX Slidematic With Qlt- 1 Denar Jig - #200012-1 (Accessories for Slidemat...	US\$ 1,166.43	+Info
119990	WHIP MIX Slidematic- 3 Jigs- Rev D Hole - #200011-6 (Accessories for Slidematic...	US\$ 1,123.69	+Info
119974	WHIP MIX Slidematic- Jig- Hanau Index & Support - #20009833 (Accessories for Sl...	US\$ 790.64	+Info
119985	WHIP MIX Slidematic- 1 Jig - HANAU #200045-1 (Accessories for Slidematic Facebo...	US\$ 783.40	+Info
119992	WHIP MIX Slidematic-rev D Hole - #200011-2 (Accessories for Slidematic Facebows...	US\$ 717.27	+Info

BITEFORK

Code	Name	Price	+
------	------	-------	---

119994	WHIP MIX Bitefork/shaft Assbly Slidematic w/qlt - #200008 (Accessories for Slid...	US\$ 709.49	+Info
119996	WHIP MIX Bitefork/shaft Assy- Dentulous - #200007 (Accessories for Slidematic F...	US\$ 150.22	+Info
119986	WHIP MIX Bitefork/slidematic - Qlt - #200038 (Accessories for Slidematic Facebo...	US\$ 49.15	+Info

INDEX JIG

Code	Name	Price	+
119993	WHIP MIX Trans Jig/index Assy-dentulous-rev D Hole - #200009-1 (Accessories for...	US\$ 323.54	+Info
119982	WHIP MIX Index/jig Assy-denar-molded - #200080 (Accessories for Slidematic Face...	US\$ 98.14	+Info

SHAFT

Code	Name	Price	+
119984	WHIP MIX Shaft Assbly/slidematic w/qlt - #200018-1 (Accessories for Slidematic ...	US\$ 619.62	+Info
119988	WHIP MIX Shaft Clamp Assembly - #200016-1 (Accessories for Slidematic Facebows ...	US\$ 182.47	+Info

ACCESSORIES

Code	Name	Price	+
119991	WHIP MIX Unit Assy - Slidematic - #20002950 (Accessories for Slidematic Facebow...	US\$ 395.80	+Info
119981	WHIP MIX Hanau Index/jig & Cast Support - #20009819 (Accessories for Slidematic...	US\$ 211.52	+Info
119979	WHIP MIX Index/jig Assy- WHIP MIX #200088 (Accessories for Slidematic Facebows ...	US\$ 143.24	+Info
119980	WHIP MIX WHIP MIX Index For Accumount Magnetic Mounting - #200088q (Accessories...	US\$ 145.95	+Info
119976	WHIP MIX Index/jig Assy- HANAU #200082 (Accessories for Slidematic Facebows for...	US\$ 126.56	+Info
119995	WHIP MIX Bitefork Assy - Dentulous - #200029 (Accessories for Slidematic Facebo...	US\$ 46.57	+Info
119987	WHIP MIX Bitefork Assy - Edentulous - #200030 (Accessories for Slidematic Faceb...	US\$ 46.57	+Info
119975	WHIP MIX Ref Plane Locator - #3266 (Accessories for Slidematic Facebows for Art...	US\$ 12.21	+Info

Face bows - AD2 DENTAL

Face bows - AD2 DENTAL

Code	Name	Price	+
469678	Anatomic Facebow (includes two bite	US\$ 539.50	+Info

Face bows - BIOART

Face bows - BIOART

Code	Name	Price	+
122989	BIOART ELITE Facebow each	US\$ 174.19	+Info
416745	STANDARD Facebow Only	US\$ 66.40	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



**US
DENTAL
DEPOT**

+1 954 874 6325 
www.usdentaldepot.com 

GENERAL FACEBOWS - WHIP MIX

GENERAL FACEBOWS - WHIP MIX

Whip Mix facebows provide a convenient, quick and accurate method of securing an axis location without the use of expensive equipment and time-consuming techniques.

**THESE FACEBOWS ARE COMPATIBLE WITH
2000, 3000, 4000, 8500 SERIES, HANAU
WIDE-VUE AND MODULAR ARTICULATORS**

Description:

- * Indirect mounting for easier access, increased stability, greater ease of use and optimum instrument efficiency
- * Detachable QuickLock Toggle (QLT) allows you to use the facebow on the next patient
- * The QLT helps prevent rotation of the bitefork and captures the patient's bite record by tightening one knob



GENERAL FACEBOWS

Code	Name	Price	+
119877	WHIP MIX #9155q Face Bow w/quick Lock/qms - #34344 (Accessories for Facebows)	US\$ 880.09	+Info
119874	WHIP MIX #9185 Quick Mount Face Bow Indirect - #34126 (Accessories for Facebows)	US\$ 583.69	+Info
119878	WHIP MIX #9185q Face Bow w/qms Transfer Base - #34358 (Accessories for Facebows)	US\$ 577.88	+Info

FACEBOWS ACCESSORIES

Code	Name	Price	+
119884	WHIP MIX #9650 QI Toggle Assy F/indirect Mtg - #34346 (Accessories for Facebows)	US\$ 722.46	+Info
119879	WHIP MIX #9650 QI Toggle Assy Direct/convert - #34341 (Accessories for Facebows)	US\$ 687.23	+Info
119983	WHIP MIX WHIP MIX Index Assembly - #200088qa (Accessories for Facebows)	US\$ 237.86	+Info
119867	WHIP MIX #9176qa Quickmount Transfer Base Assembly - #28705 (Accessories for Qu...	US\$ 134.74	+Info
119844	WHIP MIX #8579q Quickmount Remounting Jig - #03737 (Accessories for Quick Mount...	US\$ 108.87	+Info
119894	WHIP MIX #8605a Nasion Relator Assembly - #04405 (Accessories for Facebows)	US\$ 85.93	+Info
119888	WHIP MIX #8639 Toggle Assembly w/#8640 - #15954 (Accessories for Slidematic Fac...	US\$ 84.56	+Info
119852	WHIP MIX #8575 Remounting Jig - #03735 (Accessories for Articulators)	US\$ 79.30	+Info
119868	WHIP MIX #9174 Quickmount Facebow Fork Support - #28706 (Accessories for Quick ...	US\$ 73.51	+Info
119857	WHIP MIX #8587 Fork Support/indirect Bow - #33359 (Accessories for Facebows)	US\$ 70.27	+Info
119845	WHIP MIX #8585 Fork Support - #03743 (Accessories for Facebows)	US\$ 63.29	+Info
119887	WHIP MIX #8638 Horiz Slide Bar w/#8643 - #15946 (Accessories for Facebows)	US\$ 58.17	+Info
119889	WHIP MIX #9189 Support Block - #33370 (Accessories for Facebows)	US\$ 40.87	+Info
119892	WHIP MIX #8588 Fork Support Assbly For 200088q Index - #33371 (Accessories for ...	US\$ 39.60	+Info
119896	WHIP MIX #8631 Face Bow Fork Offset - #04480 (Accessories for Facebows)	US\$ 34.56	+Info

119895	WHIP MIX #8609 Face Bow Fork - #04464 (Accessories for Facebows)	US\$ 34.56	+Info
119890	WHIP MIX #9659 Face Bow Fork F/ql Toggle - #34340 (Accessories for Facebows)	US\$ 34.56	+Info
119891	WHIP MIX #9661 Fb Fork Offset/ql Toggle - #34333 (Accessories for Facebows)	US\$ 34.56	+Info
119880	WHIP MIX #9186 Support Bar F/HANAU #33375 (Accessories for Facebows)	US\$ 29.34	+Info
120018	WHIP MIX Face Bow Level 2 Pack - #20002951 (Accessories for Facebows)	US\$ 15.88	+Info
119893	WHIP MIX #8603 Plastic Ear Piece Set 2 - #04383 (Accessories for Facebows)	US\$ 38.42	+Info
119978	WHIP MIX Ref Plane Marker - #3265 (Accessories for Slidematic Facebows for Arti...	US\$ 12.21	+Info

HANAU FACEBOWS - WHIP MIX

HANAU FACEBOWS - WHIP MIX

Features:

- * Suspension from the external auditory meatus makes it a one operator bow
- * Meets the precise requirements for fixed restorative techniques
- * Available in several styles to comply with various methods of facebow registration and transfer and is adaptable to most makes of articulators
- * Scales are tilted forward for ease of reading and adjustment
- * The Bitefork/Biteplane clamps has an overhead integral lock - no separate wrenches are required. Finger pressure achieves locking without introducing torque into the box.



HANAU FACEBOWS

Code	Name	Price	+
119919	HANAU 153-16 Earpiece Facebow - Facebows Accesories - #008824-000 (Accessories ...	US\$ 560.10	+Info
119917	HANAU Spring Bow-drt/indrt Mt - Facebows Accesories - #010328-000 (Accessories ...	US\$ 476.67	+Info
119916	HANAU Spring Bow w/o Mtg Pltf - Facebows Accesories - #010351-000 (Accessories ...	US\$ 300.89	+Info

FACEBOWS ACCESSORIES

Code	Name	Price	+
119924	HANAU 182-412 Mounting Platform Assy - Facebows Accesories - #010342-000 (Acces...	US\$ 154.50	+Info
119923	HANAU 182-402 Transfer Clamp Asbly - Facebows Accesories - #010356-000 (Accesso...	US\$ 131.75	+Info
119938	HANAU 97-0 Maxillary Cast Support Assembly - Facebows Accesories - #003401-000 (...)	US\$ 78.98	+Info
119920	HANAU 132-104 Bite Fork - Facebows Accesories - #005726-000 (Accessories for Fa...	US\$ 46.01	+Info
119921	HANAU 132-110 Bite Plane - Facebows Accesories - #005727-000 (Accessories for F...	US\$ 45.55	+Info
119922	WHIP MIX 132-113 Biteplane - #005753-000 (Accessories for Facebows)	US\$ 41.42	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US
DENTAL
DEPOT

+1 954 874 6325

www.usdentaldepot.com



QUICK MOUNT FACEBOWS - WHIP MIX

QUICK MOUNT FACEBOWS - WHIP MIX

Features:

- * Direct mounting: facebow transfers directly to the articulator for mounting
- * Compatible to all Whip Mix articulators
- * Comes complete with: Whip Mix direct facebow, 8609 facebow fork



Code	Name	Price	+
119875	WHIP MIX #8645 Quick Mount Face Bow w/screws - #15911 (Accessories for Facebows)	US\$ 437.36	+Info
119881	WHIP MIX #9188 Transfer Assy - #34371 (Accessories for Facebows)	US\$ 222.02	+Info
119882	WHIP MIX #9197 Transfer Assy - #34436 (Accessories for Facebows)	US\$ 222.02	+Info
119885	WHIP MIX #9176a Transfer Base Assy - #34444 (Accessories for Facebows)	US\$ 122.00	+Info
119886	WHIP MIX #8679 Cross Bar - #06890 (Accessories for Facebows)	US\$ 37.76	+Info

Articulators and Occlusors

Fully-Adjustable articulators



Fully-Adjustable articulators - WHIP MIX

Fully-Adjustable articulators - WHIP MIX

Code	Name	Price	+
119962	D5a Kit - P2t3 Articulator	US\$ 4,888.69	+Info

Articulators and Occlusors

Magnetic Articulators



HATHO

POLISHING KIT

Thermoplastics
like *never* before



Packing Diamond D®
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory



FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



**US
DENTAL
DEPOT**

+1 954 874 6325

www.usdentaldepot.com



MAGNETIC ARTICULATORS | ECCO

Magnetic Articulators - ECCO



Code	Name	Price	+
102436	Art 3 Dimensional Magnetic	US\$ 34.20	+Info

MAGNETIC ARTICULATORS | KEYSTONE

Magnetic Articulators - KEYSTONE



ACCESSORIES

DELUXE ARTICULATOR

DELUXE ARTICULATOR

Code	Name	Price	+
034-1050076	Magnet for Articulator 4pk	US\$ 15.35	+Info

Other Models

Other Models

Code	Name	Price	+
419970	Flat Plate for Luxury 1pk	US\$ 9.86	+Info

DELUXE DENTURE MAGNETIC ARTICULATORS

DELUXE MAGNETIC ARTICULATORS - KEYSTONE

The Deluxe Magnetic Articulators comes with magnetic split casting to ensure easy removal and exact settings of models. These articulators are excellent for checking cases. It's completely interchangeable, lightweight and Teflon coated.



Code	Name	Price	+
034-1050074	KEYSTONE Deluxe Magnetic Articulator w/Pin - Light in weight and Teflon coated #...	US\$ 71.33	+Info
034-1050075	KEYSTONE Deluxe Magnetic Articulator with Magnetic Split Casting Ensures - #105...	US\$ 73.11	+Info

DELUXE MAGNETIC ARTICULATORS

DELUXE ARTICULATOR

Code	Name	Price	+
034-1050085	Luxury Art 3 Dimensional Artic	US\$ 114.10	+Info

PRECISION MAGNETIC ARTICULATOR

Other Models

Code	Name	Price	+
34-1050100	Precision Articulator	US\$ 79.17	+Info

Articulators and Occlusors

Non Adjustable articulators



HANAU - NON ADJUSTABLE ARTICULATORS - ARCON - WHIP MIX

HANAU - NON ADJUSTABLE ARTICULATORS - ARCON -WHIP MIX

These articulators are affordable, lightweight and easy to use. The non-adjustable articulators offer dental professionals, laboratories, and teaching facilities a choice of lightweight, easy to use articulators and accessories with the features they're looking for.



HANAU ARTICULATORS

Code	Name	Price	+
119972	WHIP MIX Automark Kit - P4t4 Articulator - #110407-4	US\$ 969.92	+Info
119915	HANAU Model Mate Plasterless Articulator - Articulator Non Adj - #20006899	US\$ 438.57	+Info
119913	HANAU Hanau Mate pk Assembly Articulator - Articulator Non Adj - #009856-000	US\$ 619.99	+Info
119914	HANAU Ulti-mate Articulator (1 Pk) - Articulator Non Adj -#20011209	US\$ 339.52	+Info

ARTICULATORS ACCESSORIES

Code	Name	Price	+
------	------	-------	---

119936	HANAU Ulti-mate Mounting Platform - Articulator Non Adj - #20011273 (Accessorie...	US\$ 101.21	+Info
119942	HANAU 154-110 Disposable Mounting Inserts (12 Pk) - Articulator Non Adj - #00911...	US\$ 144.94	+Info

HANDY ARTICULATORS - SHOFU

HANDY ARTICULATORS - SHOFU

Constructed on the basis of mean value of the jaw movement.

* When the centric screws are loosened, a spring located in the centre of the upper jaw permits movement and return to the original position

* The upper jaw maintains a horizontal position even when the Incisor Rod is removed

- This articulator comes with 2 easily removable, heat resistant plastic mounting plates to permit working with more than one model at a time



ACCESSORIES

Code	Name	Price	+
112052	SHOFU - 5072 - Metal Mounting Plate, Includes: (4) Universal Mounting Plates for...	US\$ 18.27	+Info
112055	SHOFU Incisor Rod lim - #5073	US\$ 10.62	+Info
112056	SHOFU Incisor Rod Screw - #5074	US\$ 5.03	+Info
112057	SHOFU Incisal Guide Pin - #5075	US\$ 5.65	+Info
112058	SHOFU Occlusal Plate lim - #5076	US\$ 24.07	+Info
112060	SHOFU Spring - #5078	US\$ 1.93	+Info
112061	SHOFU Centric Screw - #5079	US\$ 19.31	+Info

ARTICULATORS

Code	Name	Price	+
112051	SHOFU Handy IIM Artirculator With out Plate - #5070	US\$ 154.37	+Info
112053	SHOFU Handy 2M Articulator Complete - #5071	US\$ 169.67	+Info

Non Adjustable articulators - SHOFU

Non Adjustable articulators - SHOFU

Code	Name	Price	+
112024	Pro Arch I W/O Plate #4070	US\$ 223.70	+Info
112025	Pro Arch I Complete #4071	US\$ 246.90	+Info
112027	Incisal Rod Screw	US\$ 20.70	+Info
112042	Incisor Rod	US\$ 7.92	+Info
112046	Handy Iia Incisor Rd	US\$ 7.92	+Info

Articulators and Occlusors

Plasterless Articulators

ENZO UNIVERSAL PLASTERLESS ARTICULATOR - VISION

ENZO UNIVERSAL PLASTERLESS ARTICULATOR - VISION

Features:

- * The Enzo Articulator is proven to successfully integrate with many of the leading Dental CAD/CAM digital model systems.
- * The articulator continues to be a very popular choice with the use of standard models as it eliminates plaster and the resetting of articulations achieved in seconds
- * It replicates true lateral and protrusive excursions
- * Lightweight and perfectly balanced, it fits comfortably in the dental technicians hand



Code	Name	Price	+
103466	Enzo Universal Plasterless Articulator - Semiadjustable - 1001B-One Articulator ...	US\$ 191.83	+Info
113555	Enzo Disposable Plate 100 pk - Plastic universal disposable bases which slide in...	US\$ 57.58	+Info

GALLETTI TYPE ARTICULATORS - ECCO

GALLETTI TYPE ARTICULATOR - ECCO

Features:

- * Aluminum made and Sateinless steel working parts
- * Stone Free articulator
- * Also ocluser



Code	Name	Price	+
102434	Galletti type Articulator - Silver - Aluminum made and Sateinless steel working ...	US\$ 87.50	+Info
102435	ECCO - Galletti type Articulator - WHITE - Aluminum made and Sateinless steel wo...	US\$ 87.50	+Info

Plasterless Articulators - BESQUAL

Plasterless Articulators - BESQUAL

Code	Name	Price	+
103075	N-2 Plasterless Articulator	US\$ 15.92	+Info

Plasterless Articulators - KEYSTONE

Plasterless Articulators - KEYSTONE

Code	Name	Price	+
034-1050200	Galletti Art Stone Free	US\$ 544.00	+Info

Articulators and Occlusors
Plastic Articulators

Plastic Articulators - BUFFALO

Plastic Articulators - BUFFALO

Code	Name	Price	+
424605	Disposable plastic C&B 50pk	US\$ 62.62	+Info

Articulators and Occlusors
Semi-Adjustable articulators

SEMI-ADJUSTABLE ARTICULATORS | WHIP MIX



2000 SERIES ARTICULATORS - WHIP MIX

2000 SERIES ARTICULATORS - WHIP MIX

- * Interframe spacing of over 4" allowing for easy mounting of even the bulkiest casts
- * Excellent options for orthodontists and oral surgeons
- * Rounded edges and frame corners for enhanced visibility and posterior access
- * Open fossa design
- * Adjustable condylar inclination: 0-70
- * Fixed intercondylar distance: 110mm
- * Positive centric latch
- * "Accumount" system of interchangeability which allows accurate transfer of casts between articulators
- * Comes complete with: flat incisal table, dual-end incisal pin, 4 disposable mounting plates, instruction manual.



Code	Name	Price	+
119830	WHIP MIX #2220q Articulators w/stand - #34210	US\$ 1,528.51	+Info
119838	WHIP MIX #2340q Articulators w/qms - #34220	US\$ 1,333.27	+Info
119836	WHIP MIX #2240 Articulators - #02267	US\$ 1,225.25	+Info
119841	WHIP MIX #2240q Articulators w/qms MAGNETIC MOUNTING - #34219	US\$ 1,333.27	+Info

4000 SERIES ARTICULATORS - WHIP MIX

ARTICULATORS 4000 SERIES - WHIP MIX

Features:

- * 4000 Series Articulators:
- * Elimination of crossbar allows for enhanced visibility and lingual access during denture set-up
- * Easy-to-operate locking mechanisms allow:
- * Positive centric locking
- * Ability to move through excursive movements
- * Removal of upper member
- * Adjustable condylar inclination: 0-70
- * Adjustable progressive side shift: 0-25
- * Fixed intercondylar distance: 110mm
- * "Accumount" system of interchangeability which allows accurate transfer of casts between articulators
- * Support rod allows upper member stability when articulator is open
- * Comes complete with: flat incisal table, dual-end incisal pin, 4 disposable mounting plates, instruction manual



Code	Name	Price	+
119834	WHIP MIX #4640q Articulator w/qms - #02276	US\$ 1,393.08	+Info
119840	WHIP MIX #4641q Half Tracked Articulators w/ Qms - #02279	US\$ 1,393.08	+Info
119833	WHIP MIX #4641 Half Tracked Articulator - #02278	US\$ 1,284.66	+Info

8500 SERIES ARTICULATORS - WHIP MIX

8500 SERIES ARTICULATOR - WHIP MIX

Features:

Model 8500 Standard Articulators are designed to mechanically imitate the human jaw to assist in the fabrication of dentures, crowns, bridges, inlays, onlays and orthodontic appliances.

Description:

- * Adjustable condylar inclination: 0-70
- * Adjustable progressive side shift: 0-40
- * Adjustable intercondylar distance: S (96 mm); M (110 mm); L (124 mm)
- * Spring latch assembly holds articulator in centric position (when locked)
- * Spring latch assembly allows for easy separation of upper and lower members (when released)
- * Comes complete with: flat incisal table, dual-end incisal pin, 4 disposable mounting plates, instruction manual.



Code	Name	Price	+
119831	WHIP MIX #8500 Articulator - #03689 (Accessories for 8500 Serie Articulators)	US\$ 1,060.15	+Info
119869	WHIP MIX Qms Kit F/8500 Panadent Hanau Mod - #28707 (Accessories for Quick Moun...	US\$ 227.66	+Info

DENAR ARTICULATORS

DENAR - ANAMARK PLUS ARTICULATORS - WHIP MIX

DENAR - ANAMARK PLUS - WHIP MIX

Features:

The Anamark Plus Semi-Adjustable Articulator provides a fast, uncomplicated method of setting up a case for occlusal diagnosis and treatment. Side shift adjustments are simply a matter of selecting the desired fossae (.5mm or 1.0mm) and inserting them into the articulator.

Description:

- * 1 Pair of 1.0mm pre-programmed immediate side shift fossae
- * Removable upper member
- * Retaining springs hold upper and lower members together without being locked in centric
- * Positive centric latch
- * 110 mm intercondylar distance
- * Adjustable protrusive angle: 0-60
- * 7 progressive side shift
- * 3/4" (19 mm) curvature of superior wall
- * Black finish



Code	Name	Price	+
119964	WHIP MIX Anamark Kit Plus Articulator - #113000-201 (Accessories for Anamark Pl...	US\$ 1,610.38	+Info

DENAR - CADIAX COMPACT 2 - SEMI ADJUSTABLE ARTICULATORS - WHIP MIX

DENAR - CADIAX COMPACT 2 - SEMI ADJUSTABLE ARTICULATORS - WHIP MIX

System is designed for registration and display of mandibular hinge axis movements with data used to program the articulator.

- * Patient set-up and recordings are accomplished in less than 10 minutes
- * Easily fits into any teaching or practice philosophy
- * Visual aid and screening tool for TMJ anomalies

Features:

- * Real-time recording to visually locate the reference position
- * Increased resolution
- * Compatible with Windows
- * Uses new magnetic technology for optimized recording of biomechanical movements of the TMJ
- * Settings can be calculated for over 15 different models of articulators
- * Multiple recordings for comparison of tracings and accuracy verification
- * New Streamlined Cadiax unit
- * Space saving dimensions: 1.5"w x 6.5"h x 6.5"d
- * Portable
- * Plugs directly into computer USB port



Code	Name	Price	+
120046	WHIP MIX Cadiax Compact 2 w/o Gamma Software - #20002375 (Accessories for Cadi...	US\$ 16,892.72	+Info
120045	WHIP MIX Cadiax 2 Gamma Software (03-softe000c) - #20002916 (Accessories for Ca...	US\$ 5,384.13	+Info
120047	WHIP MIX Cadiax Compact 2 Recording Software Upgrade - #20002920 (Accessories f...	US\$ 1,226.98	+Info

DENAR - MARK 300 SERIES ARTICULATORS - WHIP MIX

DENAR - MARK 300 ARTICULATORS - WHIP MIX

Like all Denar products, the Mark 300 Series articulators are engineered for superior performance and tested thoroughly based on proven methodology. They bring a new standard of interchangeability with factory set accuracy within 20 microns. With their sleek, ergonomic design and powder-coated finish, these articulators have a contemporary look that's as durable as it is aesthetically pleasing.

Specifications:

- * Adjustable condylar inclination: 0-60
- * Adjustable immediate side shift: 0-2mm
- * Adjustable progressive side shift: 0-15
- * Compatible with the Denar Slidematic facebow
- * Comes complete with: Dust cover, 10 magnetic mounting plates, instructions and standard black carrying case



Code	Name	Price	+
119967	DENAR - Mark 330 Articulator - #20011389	US\$ 2,310.29	+Info
119966	DENAR - Mark 320 Articulator - #20011384	US\$ 1,915.99	+Info
119973	DENAR - Mark 310 Articulator - #20011382 (Accessories for Mark 300 Articulators)	US\$ 1,323.49	+Info

DENAR - MARK II SERIES ARTICULATORS - WHIP MIX

DENAR - MARK II SERIES - WHIP MIX

Mark II Series: Designed for simplicity and accuracy

- * Adjustments: protrusive angle: 0-60; immediate side shift: 0-4mm; progressive side shift: 5-15
- * Removable upper member
- * Retaining springs hold members together without being locked in centric
- * Clear finish
- * Positive centric latch
- * 110mm intercondylar distance
- * Rear wall posterior inclination of 25
- * Straight top wall



* Comes complete with: 10 disposable mounting plates, carrying case and instruction manual

Mark II Plus Articulator Features:

Adjustable incisal table with tapered incisal pin

Mark II Articulator Features:

Custom (flat) incisal table with short round incisal pin

Code	Name	Price	+
119969	WHIP MIX Mkii Kit P4t4 - #110500-2 (Accessories for Mark II Articulators)	US\$ 2,059.35	+Info

DENAR ACCESSORIES for ARTICULATORS

DENAR - MOUNTING PLATES ARTICULATOR S

DENAR - MOUNTING PLATES - WHIP MIX

Features:

- * Designed in conjunction with Dr. Tim Tremont and are an integral part of his Four Faces of Orthognatic Surgery system for accurate surgical outcomes
- * Upper and lower casts can be mounted to a "true horizontal" tansverse plane
- * Casts can easily be separated from the mounting plates for repositioning to the transverse plane
- * Enhances accuracy and efficiency of the cast mounting and model surgery
- * White surgical mounting plate attached to current screw type disposable mounting plate that fits Denar and Hanau? articulators



METAL MOUNTING PLATES

Code	Name	Price	+
119849	WHIP MIX #8580 Metal Mounting Plates Set 2 - #03867 (Accessories for Mounting	US\$ 21.06	+Info
119951	WHIP MIX 130-2012 Mtg Plates 1 Pair - #005057-000 (Accessories for Mounting Pla...	US\$ 22.16	+Info

MOUNTING PLATES WITH INSERT

Code	Name	Price	+
119947	WHIP MIX Mtg Plate w/insert-pk Of 10 - #1100026-10 (Accessories for Mounting P...	US\$ 42.01	+Info
119948	WHIP MIX Mtg Plate w/insert-pk Of 50 - #1100026-50 (Accessories for Mounting P...	US\$ 159.83	+Info
119946	WHIP MIX Montage Plates w/insert-pk Of 100 - #1100026-100 (Accessories for Mou...	US\$ 269.03	+Info

PLASTIC MOUNTING PLATES

Code	Name	Price	+
045-28495	WHIP MIX #8580B Plastic Mounting Plates - Box of 20 #28495 (Accessories for Mou...	US\$ 70.59	+Info
119861	WHIP MIX #8580b Plastic Mount Plates Box 20 - #28495 (Accessories for Articulat...	US\$ 70.59	+Info

DENAR - UNIVERSAL OCCLUSAL STAND ARTICULATORS

DENAR - UNIVERSAL OCCLUSAL STAND - WHIP MIX

The Denar Universal Occlusal Stand allows you to mount the working casts at Bonwills Equilateral Triangle on Denar, Hanau, and Whip Mix ("QM") articulators. Per lab owner, Al Fillastre, the new occlusal stand "allows the confident use of the articulator to capture all of its functional benefits even when a facebow is not available."

Advantages:

- * Convenient: Includes two interchangeable tables - flat and curved. The flat table accepts the Kois Index Tray.
- * Adjustable: Both the flat and curved tables may be set at 0 or 10 requiring only a small height adjustment.
- * Adaptable: Fits most Denar, Hanau and Whip Mix articulators.



ARTICULATOR

Code	Name	Price	+
120017	WHIP MIX Universal Occlusal Stand - #20011720 (Universal Occlusal Stand)	US\$ 411.89	+Info

ACCESSORIES

Code	Name	Price	+
120001	WHIP MIX Centri Check - #180035 (Universal Occlusal Stand)	US\$ 731.60	+Info
119997	WHIP MIX Centri Chk Patient Record Assy (10 Pair) - #180045 (Universal Occlusal...)	US\$ 172.07	+Info
120013	WHIP MIX Flat Table for Universal Occlusal Stand - #20011713 (Universal Occlusa...)	US\$ 64.57	+Info
120016	WHIP MIX Fence Assembly Occlusal Stand - #20011717 (Universal Occlusal Stand)	US\$ 46.08	+Info
120021	WHIP MIX D/h Thread Adaptor Occlusal Stand - #20011711 (Universal Occlusal Stand)	US\$ 48.32	+Info
120019	WHIP MIX Taper Pin - Occlusal Stand-mk300 Gage-split Cast - #20011665 (Universa...)	US\$ 11.43	+Info

HANAU ARTICULATORS

HANAU - 96H2 ARTICULATORS - WHIP MIX

HANAU - 96H2 ARTICULATORS - WHIP MIX

96H2 Semi-Adjustable, Non-Arcon Articulator:

- * Condylar guidances are associated with the lower member, which simulates the multiple functions of the Glenoid Fossa in centric, protrusive, right and left lateral positions
- * Enclosed track design prevents vertical displacement of the condyle element when waxing or making excursive movements with the instrument

96H2 Articulator Adjustments:

- * Progressive Side Shift: 0-30
- * Protrusive Angle: 0-75



ARTICULATOR

Code	Name	Price	+
119905	HANAU 96-h2 Articulator with adjustable Incisal Guide Articulator - Articulator ...	US\$ 1,423.27	+Info

ACCESSORIES

Code	Name	Price	+
119932	HANAU 158-115 Adj Guide Assembly - Accesories Serie 96H2 - #009551-000 (	US\$ 223.89	+Info
119937	HANAU 96-0 Articulator Plastering Stand - Accesories Serie 96H2 - #003368-000 (...	US\$ 66.67	+Info
119926	WHIP MIX - HANAU 96-83 Incisal Pin 130-3008/158-3001- Accesories Serie 96H2 - #0...	US\$ 36.40	+Info

HANAU - MODULAR SYSTEM ARTICULATORS - WHIP MIX

HANAU - MODULAR SYSTEM ARTICULATORS - WHIP MIX

Features:

Each Modular Semi-Adjustable Articulator features:

- * Interchangeable fossae provide maximum flexibility and precision, enabling proper occlusion for each patient
- * Programmed fossae with curved 19mm (3/4") radius simulate average anatomical movements
- * Positive centric latch
- * Interchangeability of accessories and casts
- * Wide lingual visibility
- * Removable upper member facilitate wax-ups
- * Sturdy construction and ease of use

Description:

- * Adjustable table
- * Dual-end incisal pin
- * 0-90 adjustable condylar inclination



- * 19 mm (3/4") superior radius
- * 0-30 progressive side shift
- * 0-4 mm adjustable immediate side shift.

ARTICULATORS

Code	Name	Price	+
119910	HANAU 194-611100 Radial Modular System Articulators - #014503-000 (Accessories ...	US\$ 1,404.31	+Info
119909	HANAU 190-291111 Modular System Articulators /adj. Bennett - #014102-000 (Acce...	US\$ 1,381.69	+Info
119908	HANAU 190-454110 Modular System Articulators -pr - #014116-000 (Accessories fo...	US\$ 1,344.95	+Info
119912	HANAU 192-051100 Modular System Articulators - #014300-000 (Accessories for Mod...	US\$ 1,156.66	+Info
119911	HANAU 192-254106 Modular System Articulators -side Shift - #014309-000 (Accesso...	US\$ 1,226.00	+Info

ACCESSORIES

Code	Name	Price	+
119927	HANAU 166-3009 Incisal Pin- Dual End - Modular System Accesories #010081-000 (A...	US\$ 57.63	+Info
119944	HANAU 166-107 Orbital Indicator Assembly (modular) - Accessories for Articulator...	US\$ 15.93	+Info

HANAU - WIDE-VUE ARTICULATORS - WHIP MIX

HANAU - WIDE VUE ARTICULATORS - WHIP MIX

Wide-Vue Semi-Adjustable Articulators:

- * Lingual access is unobstructed. You'll appreciate its enhanced view when positioning tooth-to-tooth contacts on models and final prostheses.
- * Closed tracking fossa and fixed upper member are ideal for fabricating prosthetic appliances.

Wide-Vue Articulator Features:

- * Progressive side shift: 0-30
- * Protrusive angle: -20+60
- * Comes complete with: orbitale indicator, dual-end incisal pin, 10 disposable mounting plates, instruction manual and carrying case

Maxillo Facial Articulator with Additional 25mm (1") Vertical Height Features:

- * Progressive side shift: 0-30
- * Protrusive angle: -20+60
- * Comes complete with: orbitale indicator, dual-end incisal pin, 2 pairs of metal mounting plates and instruction manual



ARTICULATORS

Code	Name	Price	+
119907	HANAU Package Assy- 183-2 Wide-Vue Articulator - #010885-000 (Accessories for W...	US\$ 1,376.99	+Info

ACCESSORIES

Code	Name	Price	+
119925	WHIP MIX 184-3043 Incisal Pin- Dual End 120mm - (Accessories for Wide-Vue Artic...	US\$ 56.52	+Info

119928	HANAU 183-402 Adj Incisal Pin Assy - Wide-Vue Accessories - #010670-000 (Access...	US\$ 508.51	+Info
119935	HANAU 184-404 Adjust. Incisal Guide - (Accessories for Wide-Vue Articulators)	US\$ 248.20	+Info
119941	HANAU 85-2 Remounting Record Jig - Wide-Vue Accessories #003102-000 - duplicate ...	US\$ 144.88	+Info
119943	HANAU 158-51 Remounting Record Jig - Wide-Vue Accessories - #009544-000 (Access...	US\$ 101.63	+Info

HANAU ACCESSORIES for ARTICULATORS

HANAU - ARTICULATORS ACCESSORIES - WHIP MIX



Features:

Hanau accessories enhance the capabilities of articulators systems and provide the clinician with a wide range of options for specialized needs.

ACCESSORIES COMMON FOR ALL ARTICULATORS

Code	Name	Price	+
119939	HANAU 142-1 Broadrick O P Analyzer - Accessories for Articulators - #006901-000...	US\$ 414.39	+Info
119931	HANAU 158-126 Anterior Tables-pk /12 - Accessories for Articulators - #009489-00...	US\$ 63.28	+Info
119940	HANAU 142-3401 Plastic Rec/pk 12 - Accessories for Articulators - #006916-000 -...	US\$ 50.74	+Info
119930	HANAU 158-134 Anterior Tables-pk /3 - Accessories for Articulators - #009487-000...	US\$ 17.49	+Info
119933	HANAU 184-3036 Thumbscrew - Accessories for Articulators - #010717-000 - duplic...	US\$ 7.55	+Info

ARTEX ARTICULATOR - AMANN GIRRIBACH



AMANN GIRRIBACH

ARTICULATORS

Code	Name	Price	+
467379	Ceramill Mind Artex CR	US\$ 1,800.94	+Info
467897	Artex Typ CT	US\$ 863.72	+Info
467908	Artex Typ BN	US\$ 443.50	+Info
468015	Artex Typ CP/CPR	US\$ 1,346.01	+Info
468020	Artex CP for university incl. Split	US\$ 1,593.21	+Info
468021	Artex Typ CR	US\$ 1,716.26	+Info
467899	Artex CT for university incl. Split	US\$ 1,204.31	+Info
467893	Artex Typ CN	US\$ 618.70	+Info
467982	Artex Facebow	US\$ 1,085.45	+Info
468023	Artex CR for university incl. Split	US\$ 2,012.47	+Info

QUICK MOUNT MAGNETIC SYSTEM - WHIP MIX

QUICK MOUNT MAGNETIC SYSTEM - WHIP MIX

Features:

The "QM" Magnetic System allows for quick, precise, and easy placement and removal of mounted casts on many of the most popular semi-adjustable and fully adjustable articulators. The end result is the saving of time and money. The "QM" Magnetic System is designed for use with Whip Mix, Denar or Hanau articulators.



Code	Name	Price	+
120050	WHIP MIX Quick Release Plate Assy-accuceph - #100890 (Instrument System for Rad...	US\$ 357.75	+Info
119870	WHIP MIX Qms Kit F/2240 3040 Hanau Wide Vue - #28708 (Accessories for Quick Mou...	US\$ 229.93	+Info
119864	WHIP MIX Quickmount Plates WHITE w/discs 20 - #28701 (Accessories for Quick M...	US\$ 105.60	+Info
119865	WHIP MIX Quickmount Plate Holder 13mm - #28703 (Accessories for Quick Mount Art...	US\$ 60.62	+Info
119866	WHIP MIX Quickmount Plate Holder 25mm - #28704 (Accessories for Quick Mount Art...	US\$ 60.62	+Info

Semi-Adjustable articulators - AD2 DENTAL

Semi-Adjustable articulators - AD2 DENTAL

Code	Name	Price	+
469626	Articulator with Straight Incisal P	US\$ 585.00	+Info
469627	Articulator with Curved Incisal Pin	US\$ 641.25	+Info
469628	Articulator with Straight Incisal P	US\$ 525.00	+Info
469629	Articulator w/ Curved Incisal Pin	US\$ 520.00	+Info
474818	Professional System, Magnetic	US\$ 1,950.00	+Info
474936	Professional System, Magnetic	US\$ 1,917.50	+Info

Semi-Adjustable articulators - ARTUS

Semi-Adjustable articulators - ARTUS

Code	Name	Price	+
469802	ARTUS 1000 Articulator w/Facebow	US\$ 115.05	+Info

HATHO

POLISHING KIT

Thermoplastics
like *never* before



Packing Diamond D®
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory



Semi-Adjustable articulators - BIOART

Semi-Adjustable articulators - BIOART



A7-PLUS Articulator

A7-PLUS Articulator

Code	Name	Price	+
122545	A7-PLUS Articulator w/Std Facebow	US\$ 213.82	+Info
474909	A7-PLUS E Articulator w/oFacebow	US\$ 219.73	+Info
474910	A7-PLUS E Articulat w/Std Facebow	US\$ 305.82	+Info
474911	A7-PLUS E Articulat w/Elite Facebow	US\$ 396.10	+Info

Other Models

Other Models

Code	Name	Price	+
122878	A7-PLUS Articulator w/o Facebook	US\$ 147.63	+Info
124661	4000-S Articulator	US\$ 173.69	+Info
124662	SURGICAL MT Articulator no/Fb	US\$ 619.20	+Info
416523	A7-PLUS Articulator w/ ELITE	US\$ 388.25	+Info
474410	A7-PLUS Articulator	US\$ 173.69	+Info

Catalog Index / Indice

Accessories for articulators	1
Crown and Bridge Articulators	7
Denture Articulators	8
Disposable articulators	10
Face bows	12
Fully-Adjustable articulators	16
Magnetic Articulators	16
Non Adjustable articulators	18
Plasterless Articulators	19
Plastic Articulators	21
Semi-Adjustable articulators	21