

Brushes

POLISHING BRUSHES

Bristle Brushes HP



BRISTLE BRUSHES - ECCO

BRISTLE BRUSHES - ECCO

Features:

- * Bristle Brushes are classic, high quality natural brushes ideal for cleaning and polishing of denture acrylic
- * For general laboratory polishing applications.
- * These brushes also excel at applying polishing materials and at polishing lab-processed composites, metal and ceramic crowns.
- * Durable design features permanently attached bristles for improved performance and extended life without flaring
- * Solid carbon-steel-plated HP mandrels for excellent fit and durability.

Abbott Robinson Brushes:

- * Abbott Robinson Brush Size #9-(14mm)
- * Abbott Robinson Brush Size #11-(17mm)
- * Abbott Robinson Brush Size #12-(25mm)



BRISTLE BRUSHES 17 mm (Ecco)

Code	Name	Price	+
101894	Bristle Brush Soft - 17mm - 11/16in Diameter #11 - Grey -12pk -Mounted 15000rpm ...	US\$ 5.90	+Info
101895	Bristle Brush Standard Stiff - 17mm - 11/16in Diameter #11SS - WHITE -12pk -Moun...	US\$ 5.90	+Info
101896	Bristle Brush Extra Stiff - 17mm - 11/16in Diameter #11XS - BLACK -12pk - Mounte...	US\$ 5.90	+Info

BRISTLE BRUSHES 14 mm (Ecco)

Code	Name	Price	+
117324	Bristle Brush Soft - 14mm - 9/16in Diameter #09 - Grey -12pk - Mounted 15000rpm ...	US\$ 5.84	+Info
117323	Bristle Brush Standard Stiff - 14mm - 9/16in Diameter #09SS - WHITE -12pk - Moun...	US\$ 5.90	+Info
117325	Bristle Brush Extra Stiff - 14mm - 9/16in Diameter #09 - BLACK -12pk - Mounted 1...	US\$ 5.84	+Info

BRISTLE BRUSHES 19 mm (Ecco)

Code	Name	Price	+
101897	Bristle Brush Soft - 19mm - 3/4in #12 Grey -12pk -Mounted 15000rpm #209-1201	US\$ 5.90	+Info
101898	Bristle Brush Standard Stiff - 19mm - 3/4in #12 WHITE -#12SS - 12pk -Mounted 150...	US\$ 5.90	+Info
101899	Bristle Brush Extra Stiff - 19mm - 3/4in #12 BLACK -12pk -Mounted 15000rpm #209-...	US\$ 5.90	+Info

Bristle Abbot Robinson Brushes - BUFFALO

Bristle Abbot Robinson Brushes - BUFFALO

Code	Name	Price	+
425300	Abbott-Robinson Stra	US\$ 149.13	+Info
424760	A-R HP #9 Soft Brushes 144pk	US\$ 138.81	+Info
424775	Abbott-Robinson Straight HP Mounted	US\$ 16.18	+Info
424806	Abbott-Robinson Straight HP Mounted	US\$ 16.18	+Info
424889	A-R Brushes #9 Soft Mounted 12pk	US\$ 16.62	+Info
425046	Abbott-Robinson Straight HP Mounted	US\$ 334.22	+Info
425120	Abbott-Robinson Small Inverted Cone	US\$ 13.57	+Info
425253	Nylon Washout Brush	US\$ 7.33	+Info
425256	Scouring Brush No 3P	US\$ 11.16	+Info
424752	A-R Straight HP Mounted #11 SS	US\$ 134.87	+Info
425301	Abbott-Robinson Straight HP Mounted	US\$ 144.78	+Info
425302	Abbott-Robinson Straight HP Mounted	US\$ 144.78	+Info
425304	A-R Straight HP Brush # 9 SS	US\$ 132.17	+Info
425305	A-R Brushes #11 XStiff Mounted 12pk	US\$ 148.99	+Info
425307	Abbott-Robinson Straight HP Mounted	US\$ 119.30	+Info
425308	Abbott-Robinson Straight HP Mounted	US\$ 149.13	+Info
425309	A-R HP Mounted Brushes #11 12x12pk	US\$ 141.92	+Info
424661	Scouring Brush, No.	US\$ 11.16	+Info
424526	Abbott-Robinson Straight HP Mounted	US\$ 149.13	+Info
424535	Abbott-Robinson Straight HP Mounted	US\$ 137.72	+Info
424557	A-R HP Mounted Brushes 12pk	US\$ 16.18	+Info
424573	Abbott-Robinson St HT #11 144pk	US\$ 138.81	+Info
424579	Abbott-Robinson St HT #12	US\$ 134.87	+Info
424583	Abbott-Robinson Str HT #12 stnd	US\$ 134.87	+Info
424598	Abbott-Robinson Str HP Large	US\$ 38.28	+Info
424608	Abbott-Robinson Str #11 Soft 12pk	US\$ 16.62	+Info
009-06194	#12 Abbott Rob X Stiff 144pk	US\$ 134.87	+Info
424695	A-R Mounted Brush #11 XS 12pk	US\$ 16.18	+Info
424703	Abbott-Robinson Straight HP Mounted	US\$ 15.79	+Info
424704	Abbott-Robinson Bristle Porte Polis	US\$ 18.06	+Info
424706	A-R Straight HP Mount SS#12-12pk	US\$ 16.35	+Info
424745	A-R HP Mounted Brushes #11 12pk	US\$ 16.18	+Info
424747	A-R Straight HP Mounted#11 XS 144pk	US\$ 134.87	+Info
424748	A-R Straight HP Mount #12 S 144pk	US\$ 140.21	+Info

Bristle Brushes - KEYSTONE

Bristle Brushes - KEYSTONE



Extra Hard

Extra Hard

Code	Name	Price	+
446937	Hatho Toro Miniature Brushes 3/4" E	US\$ 139.30	+Info

Extra Soft

Extra Soft

Code	Name	Price	+
446936	Hatho Toro Miniature Brushes 3/4",	US\$ 139.30	+Info

Extra Stiff

Extra Stiff

Code	Name	Price	+
034-1170070	#9 Bristle X Stiff 144pk	US\$ 142.26	+Info
034-1170071	#9 Bristle X Stiff 12pk	US\$ 17.63	+Info
034-1170080	#11 Bristle X Stiff 144pk	US\$ 142.26	+Info
034-1170081	#11 Bristle X Stiff 12pk	US\$ 17.63	+Info
034-1170090	#12 Bristle X Stiff 144pk	US\$ 142.26	+Info
034-1170091	#12 Bristle X Stiff 12pk	US\$ 17.63	+Info

Hard

Hard

Code	Name	Price	+
446932	Hatho Toro Mini Brush HARD 100pk	US\$ 139.30	+Info
446933	Hatho Toro Mini Br 11/16 Hard 100pk	US\$ 139.30	+Info

Soft

Soft

Code	Name	Price	+
034-1170010	#9 Bristle Soft 144pk	US\$ 138.25	+Info
034-1170011	#9 Bristle Soft 12pk	US\$ 19.88	+Info
034-1170020	#11 Bristle Soft 144pk	US\$ 138.25	+Info
034-1170021	#11 Bristle Soft 12pk	US\$ 16.65	+Info
034-1170030	#12 Bristle Soft 144pk	US\$ 155.47	+Info

034-1170031	#12 Bristle Soft 12pk	US\$ 16.65	+Info
420872	Hatho Goat hair Brush 12pk	US\$ 25.43	+Info
446934	Hatho Toro Miniature Brushes 3/4" S	US\$ 139.30	+Info
446935	Hatho Toro Miniature Brushes 11/16"	US\$ 139.30	+Info

Standard Stiff

Standard Stiff

Code	Name	Price	+
034-1170040	#9 Bristle Stnd Stiff 144pk	US\$ 140.56	+Info
034-1170041	#9 Bristle Stnd Stiff 12pk	US\$ 17.33	+Info
034-1170051	#11 Bristle Stnd Stiff 12pk	US\$ 17.33	+Info
034-1170060	#12 Bristle Stnd Stiff 144pk	US\$ 140.56	+Info
034-1170061	#12 Bristle Stnd Stiff 12pk	US\$ 17.33	+Info

Bristle Brushes - RENFERT

Bristle Brushes - RENFERT

Bison

Bison

Code	Name	Price	+
417518	Bison Polish brush 14mm 25pk	US\$ 28.94	+Info
417519	Bison Polish brushes 14mm 100pk	US\$ 102.05	+Info
417520	Bison Polish brush 18mm 25pk	US\$ 30.93	+Info
417521	Bison Polish brush 18mm 100pk	US\$ 112.36	+Info

Goathair

Goathair

Code	Name	Price	+
417441	Goathair Brush metal & acryl 12pk	US\$ 30.86	+Info

Slender

Slender

Code	Name	Price	+
417523	Slender brushes 36mm 12pk	US\$ 72.35	+Info

HATHO GEZA - KEYSTONE

HATHO GEZA - BRISTLE BRUSHES - KEYSTONE

The exclusive geared design and premium bristle materials consist of Goat Hair and Chung King. They offer the perfect miniature brush for deep, interproximal and occlusal polishing.

Benefits:

- * Chung King hear .
- * 3/32" shank .
- * Ideal speed is set at 15,000 RPM.
- * Available packages of 12 Only.



* Available in Grey, White and Black

US DENTAL DEPOT Inc is an Authorized KEYSTONE Dealer in the US.

Code	Name	Price	+
054-1670088	KEYSTONE Geza brush for deep Interproximal and Occlusal polishing, 3/32in shank...	US\$ 21.00	+Info
054-1670089	KEYSTONE Geza brush for deep Interproximal and Occlusal polishing, 3/32in shank...	US\$ 21.00	+Info
054-1670090	KEYSTONE Geza brush for deep Interproximal and Occlusal polishing, 3/32in shank...	US\$ 21.00	+Info

Condensing Brushes



Condensing Brushes - JSP

Condensing Brushes - JSP

Code	Name	Price	+
112200	Brush & Muslin Conv 3 in each	US\$ 10.73	+Info

Condensing Brushes - KEYSTONE

Condensing Brushes - KEYSTONE

Code	Name	Price	+
420873	Lathe Brush w/ Plast core 2.5 dia	US\$ 44.40	+Info
446938	Keystone Chungking Lathe Brush 1400	US\$ 28.28	+Info

Condensing Brushes - RENFERT

Condensing Brushes - RENFERT

Code	Name	Price	+
417522	Chungking brushes 80mm 12pk	US\$ 162.62	+Info

Lathe Brushes



1 row

1 row

Code	Name	Price	+
425338	Abbott-Robinson #B4 wood hub brush	US\$ 6.36	+Info

2 rows

2 rows

44mm

44mm

Plastic Hub

Plastic Hub

Code	Name	Price	+
023-788-1000	Slim Polishing Brushes- 44mm	US\$ 42.90	+Info
417524	SLIM Polish brush 44mm 100pk	US\$ 284.70	+Info

Wood Hub

Wood Hub

BESQUAL

BESQUAL

Code	Name	Price	+
101904	Wood B12/6C 47mmx 2row Conv 12pk	US\$ 20.72	+Info

BUFFALO

BUFFALO

Code	Name	Price	+
425342	Abbott-Robinson #B12 wood hub brush	US\$ 6.74	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1170100	Wood B12 47mmx 2row Conv 12pk	US\$ 34.43	+Info

50mm

50mm

Plastic Hub

Plastic Hub

Code	Name	Price	+
034-1170180	Plastic B12 47mmx 2r Conv 12pk	US\$ 27.08	+Info
034-1170190	Plastic B14 51mmx 2row Conv 12	US\$ 27.30	+Info

Wood Hub

Wood Hub

BUFFALO

BUFFALO

Code	Name	Price	+
425303	Abbott-Robinson #B14 wood hub brush	US\$ 5.98	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1170110	Wood B14 51mmx 2row Conv 12pk	US\$ 34.04	+Info

58mm

58mm

BESQUAL

BESQUAL

Code	Name	Price	+
101901	Wood B20/6B 63mmx 2row Conv 12pk	US\$ 22.32	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1170120	Wood B20 58mmx 2row Conv 12pk	US\$ 43.95	+Info

64mm

64mm

Plastic Hub

Plastic Hub

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1170200	Plastic B20 2.5in - 1/Pk	US\$ 28.65	+Info

RENFERT

RENFERT

Code	Name	Price	+
------	------	-------	---

417440

Special Brush Polishing 65mm 100pk

US\$ 367.63

[+Info](#)

Wood Hub

Wood Hub

Code	Name	Price	+
424464	Abbott-Robinson B20 Lathe Wood 12pk	US\$ 56.19	+Info
424688	A-R Lathe B20 Wood Brush 1pk	US\$ 8.80	+Info

3 rows

3 rows

44mm

44mm

Code	Name	Price	+
424728	Abbott-Robinson #B23 wood hub brush	US\$ 8.47	+Info

64mm

64mm

Plastic Hub

Plastic Hub

Code	Name	Price	+
034-1170210	Plastic B24 63mmx 3row Conv 12	US\$ 45.75	+Info

Wood Hub

Wood Hub

BESQUAL

BESQUAL

Code	Name	Price	+
101903	Wood B24/6A 63mmx 3row Conv 12pk	US\$ 23.20	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1170170	Wood B24 63mmx 3row Conv 12pk	US\$ 46.65	+Info

70mm

70mm

BUFFALO

BUFFALO

Code	Name	Price	+
425343	Abbott-Robinson #B25 wood hub brush	US\$ 9.17	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1170130	Wood B23 70mmx 3row Conv 12pk	US\$ 48.75	+Info
034-1170140	Wood B25 70mmx 3row Conv 12pk	US\$ 54.90	+Info

4 rows

4 rows

Plastic Hub

Plastic Hub



ATTRITOR

ATTRITOR

Code	Name	Price	+
475636	Polishing Brushes Plast Hub 4Rows	US\$ 4.40	+Info

BESQUAL

BESQUAL

Code	Name	Price	+
100778	Plastic B27 76mm 4row Conv 12p	US\$ 32.43	+Info
100778-1	Plastic B27 76mm 4 row Each	US\$ 4.60	+Info

Wood Hub

Wood Hub

BUFFALO

BUFFALO

Code	Name	Price	+
424721	Abbott-Robinson #B29 wood hub brush	US\$ 12.56	+Info
424755	Abbott-Robinson #B27 wood hub brush	US\$ 12.47	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1170150	Wood B27 76mmx 4row Conv 12pk	US\$ 77.14	+Info
034-1170160	Wood B29 76mmx 4row Conv 12pk	US\$ 57.75	+Info

Scotch Brite Brushes



ATTRITOR

ATTRITOR

Code	Name	Price	+
475627	Polishing Brushes Coarse 1pk	US\$ 1.62	+Info
475628	Polishing Brushes	US\$ 1.62	+Info
475630	Polishing Brushes Fine 1pk	US\$ 1.62	+Info
475631	Polishing Brushes Med Fine 1pk	US\$ 1.62	+Info

LATHE POLISHING DISCS SCOTCH BRITE - HATHO

POLISHING DISCS WITH SCOTCH BRITE - HATHO

Hatho Polishing Discs with Scotch Brite Provide superior performance to arbor bands or burs. The Scotch Brite line is ideal for working on softer materials such as dentures and removal of soft denture relined material. The core is wood encased in plastic, allowing easy mounting/dismounting from the lathe and ensures the disk remains true. The removable plastic side rings are a patented feature which exposes a fresh polishing area when needed.

The ideal speed for the 4" discs is around 2,800 RPM, while the 3" discs is around 1,400 rpm. Sold individually.



HATHO

HATHO

Code	Name	Price	+
054-1670008	Scotch Brite 100x25mm Med 1pk	US\$ 17.70	+Info
054-1670009	Scotch Brite 75x25mm Med 1pk	US\$ 16.01	+Info
054-1670010	Scotch Brite 100x13mm Med 1pk	US\$ 14.77	+Info
054-1670011	Scotch Brite 75x13mm Med 1pk	US\$ 12.16	+Info
054-1670012	Scotch Brite 100x25mm Fine 1pk	US\$ 17.70	+Info

054-1670014	Scotch Brite 100x13mm Fine 1pk	US\$ 13.94	+Info
-------------	--------------------------------	------------	-----------------------

KEYSTONE

KEYSTONE

Code	Name	Price	+
054-1670015	Scotch Brite 75x13mm Fine 1pk	US\$ 12.78	+Info
446566	Keystone Chungking Lathe Brush 2800	US\$ 7.79	+Info

Scotch Brite Brushes - BUFFALO

Scotch Brite Brushes - BUFFALO

Code	Name	Price	+
424537	Pre-Polishers 1-1/8"	US\$ 61.20	+Info

Scotch Brite Brushes - YETI

Scotch Brite Brushes - YETI

Code	Name	Price	+
417027	SATI BRUSH Polish wheel 50pk	US\$ 57.03	+Info

SCOTCH BRITE BRUSHES HP - HATH

O

SCOTCH BRITE BRUSHES - HATHO

Features:

- * 1in diameter
- * Available in Coarse, Medium, Fine
- * Available in 12/PK

Excellent for MOLLOPLAST and soft reliners grinding and polishing
 You can use the Hatho Miniature Scotch Brite wheel to polish services without ruining delicate anatomical surfaces on dental appliances. By using low speeds around 5,000 rpms and light pressure the acrylic should not overheat due to friction. These brushes come in coarse, medium and fine to let you control the type of polishing finish you need.



Code	Name	Price	+
054-1670016	KEYSTONE Scotch Brite Brushes - Coarse Polishing discs , Mounted on mandrel, 1...	US\$ 24.32	+Info
054-1670091	KEYSTONE Hatho Miniature HP Scotch Brite Brush for Final Buffing and Smoothing -...	US\$ 24.57	+Info
054-1670092	KEYSTONE Hatho Miniature HP Scotch Brite Brush for Final Buffing and Smoothing -...	US\$ 24.57	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US
DENTAL
DEPOT

+1 954 874 6325
www.usdentaldepot.com



Wire Brushes



MOUNTED WIRE BRUSHES - KEYSTONE

Features:

- * Abrasion - Free pre-polishing of metal occlusal surfaces
- * Ideal for maintaining the occlusal or approximal stops
- * Silky matt shine after polishing
- * Preparation for high-lustre polish
- * Steel/Brass Mounted Wire Brushes



US Dental Depot Is An Authorized Keystone Dealer Is
US

Code	Name	Price	+
034-1170245	KEYSTONE Rotary Wire Stiff Brass Brushes, Mounted on HP-sized Mandrels Designed...	US\$ 25.88	+Info
034-1170240	KEYSTONE HP Mounted Silver Wire Brushes 844B - 1in - 25mm - Stainless Steel - 12...	US\$ 25.88	+Info
034-1170230	KEYSTONE Rotary Wire Stiff Steel Brushes, Mounted on HP-sized Mandrels Designed...	US\$ 25.88	+Info

MOUNTED WIRE BRUSHES - RENFER

T

MOUNTED WIRE BRUSHES - RENFERT

Function & Performance

Abrasion-free pre-polishing of metal occlusal surfaces.

Advantages

- * Occlusal and approximal contact points are maintained.
- * Produces a silk-matt finish.
- * Ideal preparation for the high-lustre polish.
- * Max RPM 5000

**US DENTAL DEPOT is an Authorized
RENFERT dealer in US.**



Code	Name	Price	+
023-166-0000	RENFERT - HP Mounted Silver Wire Brushes 19mm 12pk-Stainless Steel- #1660000 #1...	US\$ 58.53	+Info
023-196-0000	RENFERT - Lathe Silver Wire Brushes 51mm 2 pieces-Stainless Steel- #1960000 #19...	US\$ 64.91	+Info

Wire Brushes - BUFFALO

Wire Brushes - BUFFALO

Code	Name	Price	+
425265	1" Steel Wire HP rot	US\$ 56.03	+Info
425266	3/4" Brass Wire HP r	US\$ 53.59	+Info

WIRE BRUSHES - SHARK

WIRE BRUSHES - SHARK

Brass is a softer wire than stainless steel or gold, and brass brushes are ideal for applications where you do not want to scratch other surfaces. Our brass brushes are made with 360 alloy, and are resistant to rust and corrosion. Brass brushes can be made with wire diameters that range from extremely soft to stiff, depending on your application. Abrasion-free pre-polishing of metal occlusal surfaces - bronze gold color

Specifications:

- * Lenth HP
- * Mx.Speed: 15.000
- * Diameter: 7/8 in = 22 mm - Shank 1/8



STAINLESS STEEL

+ Price	Code	Name		
	115801-1	2pk Mounted Stainless Steel Wire Brushes - Diam: 7/8 in=22mm - Shank 1/8 - Lenth...	US\$ 5.39	+Info
	115801	Mounted Stainless Steel Wire Brushes - Diam: 7/8 in=22mm - Shank 1/8 - Lenth HP ...	US\$ 12.29	+Info

GOLD BRASS

	Code	Name	Price	+
	115471-1	2pk Mounted Gold Brass Wire Brushes - Diam: 7/8 in=22mm - Shank 1/8 - Lenth HP -...	US\$ 3.35	+Info
	115471	Mounted Gold Brass Wire Brushes - Diam: 7/8 in=22mm - Shank 1/8 - Lenth HP - 12 ...	US\$ 11.61	+Info

Brushes

Porcelain Brushes

PORCELAIN BRUSHES | MPF



MPF

BRUSH CO

REPLACEMENT TIPS | MPF Brushes

HYBRID TIPS

	Code	Name	Price	+
	474234	MPF Hybrid Brush Tip # 2	US\$ 17.41	+Info
	474124	MPF Hybrid Brush Tip # 4	US\$ 23.88	+Info
	474125	MPF Hybrid Brush Tip # 6	US\$ 35.03	+Info
	474126	MPF Hybrid Brush Tip # 8	US\$ 43.44	+Info

TIPS | OPTIMUM

OPTIMUM SERIES

REPLACEMENT TIPS - MPF Brush Co

Features:

Replaceable premium quality tips for all your. Evolution tools.

#UsddMpf



Code	Name	Price	+
116494	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- size #2 (2/pk) Mfg.#...	US\$ 31.94	+Info
116495	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- size #4 (2/pk) Mfg.#...	US\$ 59.72	+Info
116496-1	MPF BRUSH Optimum (Nondas Design) ONE TIP Replacement Brush Tips- size #6 (1/p...	US\$ 45.60	+Info
116496	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- size #6 (2/pk TIPS)...	US\$ 91.30	+Info
116497	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- size #8 (2/pk) Mfg.#...	US\$ 142.80	+Info
116497-1	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- size #8 (1/pk) Mfg.#...	US\$ 80.38	+Info
116498	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- Stain (2/pk) Mfg.#10...	US\$ 32.04	+Info
116499	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- Opaque (2/pk) Mfg.#1...	US\$ 34.00	+Info
116500	MPF BRUSH Optimum (Nondas Design) Replacement Brush Tips- Contour (2/pk) Mfg.#...	US\$ 37.78	+Info

TIPS | EVOLUTION

Code	Name	Price	+
116544	MPF BRUSH Evolution Replaceable Brush Tips- size #8 (1/pk) Mfg.#100-1008 1161008...	US\$ 68.79	+Info
116545	MPF BRUSH Evolution Replaceable Brush Tips- Stain (1/pk) Mfg.#116-1009 1161009 -...	US\$ 11.48	+Info
116546	MPF BRUSH Evolution Replaceable Brush Tips- Glaze (1/pk) Mfg.#116-1010 1161010 o...	US\$ 20.05	+Info
116547	MPF BRUSH Evolution Replaceable Brush Tips- Contour (1/pk) Mfg.#116-1011 1161011...	US\$ 12.64	+Info

TIPS | REVOLUTION

Code	Name	Price	+
420876	MPF BRUSH Revolution Replaceable Brush Tips- size #6 (1/pk) Mfg.#12001006 200-10...	US\$ 37.89	+Info
420877	MPF BRUSH REvolution Replaceable Brush Tips- Glaze (1/pk) Mfg.#200-1010 2001010 ...	US\$ 24.89	+Info
421038	MPF BRUSH rEvolution Replaceable Brush Tips- size #2 (1/pk) Mfg.#200-1002 200100...	US\$ 23.99	+Info
421039	MPF BRUSH REvolution Replaceable Brush Tips- size #4 (1/pk) Mfg.#200-1004 200100...	US\$ 25.85	+Info
421041	MPF BRUSH REvolution Replaceable Brush Tips- size #8 (1/pk) Mfg.#200-1008 200100...	US\$ 48.24	+Info
421042	MPF BRUSH REvolution Replaceable Brush Tips- Stain (1/pk) Mfg.#200-1009 2001009 ...	US\$ 11.48	+Info
421044	MPF BRUSH REvolution Replaceable Brush Tips- Contour (1/pk) Mfg.#200-1011 200101...	US\$ 12.39	+Info

TIPS | STANDARD

Replaceable premium quality tips for all your Original Spring Activated Brushes and Original Non-Spring Activated Brushes.



Code	Name	Price	+
123902	MPF BRUSH Replacement Spring Brush Tips- size #2 - 1/pk 100-1002 - #MPF BRUSH sa...	US\$ 19.40	+Info
123903	MPF BRUSH Replacement Spring Brush Tips- size #4 - 1/pk 100-1004 - #MPF BRUSH sa...	US\$ 30.58	+Info
123904	MPF BRUSH Replacement Spring Brush Tips- size #6 - 1/pk 100-1006 - #MPF BRUSH sa...	US\$ 47.84	+Info
123905	MPF BRUSH Replacement Spring Brush Tips- size #8 - 1/pk 100-1008 - #MPF BRUSH sa...	US\$ 73.90	+Info

HYBRID PLUS BRUSHES

Code	Name	Price	+
475101	MPF Brush - Hybrid Plus Brush Kit 4 in 1 - # 700-0000	US\$ 76.97	+Info
475102	MPF Brush - Hybrid Plus Brush #2 - # 700-0002	US\$ 10.20	+Info
475103	MPF Brush - Hybrid Plus Brush #4 - # 700-0004	US\$ 14.61	+Info
475104	MPF Brush - Hybrid Plus Brush #6 - # 700-0006	US\$ 16.51	+Info
475105	MPF Brush - Hybrid Plus Brush #8 - # 700-0008	US\$ 19.00	+Info
475106	MPF Brush - Hybrid Plus Brush #9 (STAIN) - # 700-0009	US\$ 7.27	+Info

EVOLUTION BRUSHES

EVOLUTION BRUSHES - MPF Brush Co

Brushed aluminum and carbon fiber. Slick.

Inspired by Peter Pizzi, MDT, Evolution was carefully crafted and developed with the highest quality materials.

The instrumentation is machined out of premium, lightweight brushed aluminum. The grip is thoughtfully made of carbon fibre for elegance sake.

- * Finest kolinsky hair in the world.
- * Exclusive hair treatment & ferrule bonding process.
- * Designed by master ceramist Peter Pizzi.



#UsddMpf

BRUSHES

Code	Name	Price	+
116536	MPF BRUSH Evolution Brush- Size #4 1pk Mfg.#116-0004 1160004 - #BRUSH co MPF BRU...	US\$ 64.98	+Info
116537	MPF BRUSH Evolution Brush- Size #6 1pk Mfg.#116-0006 1160006 - #MPF BRUSH co	US\$ 47.63	+Info
116538	MPF BRUSH Evolution Brush Size #8 1pk Mfg.#116-0008 1160008 - #MPF BRUSH co	US\$ 100.06	+Info
116539	MPF BRUSH Evolution Stain Brush 1pk Mfg.#116-0009 1160009 (Replaced by Revoluti...	US\$ 24.00	+Info
116540	MPF BRUSH Evolution Glaze Brush 1pk Mfg.#116-0010 1160010 - #MPF BRUSH co	US\$ 23.99	+Info
116541	MPF BRUSH Evolution Contour Brush 1pk DISCONTINUED SEE REPLACEMENT Mfg.#116-0011...	US\$ 43.22	+Info

EVOLUTION MASTER KITS

EVOLUTION MASTER KITS BRUSHES - MPF Brush Co

Premium line of application specific kolinsky

Natural kolinsky hair:

Evolution boasts only the rarest and most precious male Siberian kolinsky hair in the world. From the coldest regions on earth - for the ceramist, this means an exclusive feeling when building ceramic - a superior capacity for moisture retention, paired with the perfect balance of resiliency, tension and strength.

#UsddMpf



Code	Name	Price	+
------	------	-------	---

116548	MPF BRUSH Evolution Master Brush Kit- size #4 Mfg.#116-2004 1162004 - #MPF BRUSH	US\$ 117.02	+Info
116549	MPF BRUSH Evolution Master Brush Kit- size #6 Mfg.#116-2006 1162006 - #MPF BRUSH	US\$ 188.21	+Info
116550	MPF BRUSH Evolution Master Brush Kit- size #8 Mfg.#116-2008 1162008 - #MPF BRUSH	US\$ 193.37	+Info

OPTIMUM SPRING BRUSHES LINE

MPF OPTIMUM Brushes - MPF BRUSH

OPTIMUM SPRING BRUSHES

SPRING ACTIVATED BRUSHES - MPF Brush Co

Innovative spring activation system. The perfect tip with a simple flick'. Replaceable premium Kolinsky sable brushtips.

The Dental M.P.F. Spring Activated Ceramic Brush features the unique and world-renowned 'flick' action, which reshapes and contours the brush hair with one simple movement. High-end replaceable hair tips are made of the finest 100% male Siberian Kolinsky hair. The high-quality, light aluminum handle is also available in five brilliant colors.

#UsddMpf



BRUSHES

Code	Name	Price	+
123000	MPF BRUSH Optimum (Nondas Design) Spring Brush- Size #4 1pk - Silver (line) ...	US\$ 33.31	+Info
123001	MPF BRUSH Optimum (Nondas Design) Spring Brush- Size #6 1pk - LIGHT BLUE (lin...	US\$ 78.55	+Info
123002	MPF BRUSH Optimum (Nondas Design) Spring Brush- Size #8 1pk - BLACK (line) M...	US\$ 80.28	+Info

KITS

Code	Name	Price	+
123003	MPF BRUSH Optimum (Nondas Design) Spring 3-in-1 Brush Kit- Silver (line) Mfg...	US\$ 139.45	+Info
123004	MPF BRUSH Optimum (Nondas Design) Spring 3-in-1 Brush Kit- LIGHT BLUE (line)...	US\$ 145.79	+Info
123005	MPF BRUSH Optimum (Nondas Design) Spring 3-in-1 Brush Kit- BLACK (line) Mfg....	US\$ 136.28	+Info

OPTIMUM LADIES SPRING BRUSHES

OPTIMUM LADIES SPRING BRUSHES - MPF Brush Co

The latest line of veneering kolinsky tools from European brush maker M.P.F.

Stunning new handle desing.

Optimum sets new design standards in balance, structure and ergonomics. The handle is precision machined out of lightweight aluminum, and features a sleek curvature profile for total comfort and control.

The comfort of flick innovation.

With one single movement, the brush hair is reshaped and contoured to a sharp, working point. With just a flick, excess water is expelled from the brush without wiping or handling the hair in any way



#UsddMpf

BRUSHES

Code	Name	Price	+
116505	MPF BRUSH Optimum (Nondas Design) Lady Spring Brush- Size #4 1pk - PINK (	US\$ 44.68	+Info
116506	MPF BRUSH Optimum (Nondas Design) Lady Spring Brush- Size #6 1pk - PURPLE (	US\$ 63.74	+Info
116507	MPF BRUSH Optimum (Nondas Design) Lady Spring Brush- Size #8 1pk - Midnight BL...	US\$ 92.45	+Info

KITS

Code	Name	Price	+
116508	MPF BRUSH Optimum (Nondas Design) Spring Lady- 3-in-1 Brush Kit - PINK Mfg.#10...	US\$ 161.93	+Info
122967	MPF BRUSH Optimum (Nondas Design) Spring Lady- 3-in-1 Brush Kit - PURPLE	US\$ 153.83	+Info
122968	MPF BRUSH Optimum (Nondas Design) Spring Lady- 3-in-1 Brush Kit - Midnight BLU...	US\$ 153.83	+Info

OPTIMUM MEN SPRING BRUSHES

OPTIMUM MEN SPRING BRUSHES - MPF Brush Co

Spring activated tip. Perfectly balanced. Replaceable premium Kolinsky sable brushtips.

Integrated thorn allows to release pinned die-segments from model base.

Colors: Champagne, Midnight and Citrus

Size Availables: 4-6-8



#UsddMpf

KITS

Code	Name	Price	+
116501	MPF BRUSH Optimum (Nondas Design) Spring 3-in-1 Brush Kit- Midnight BLUE (men...	US\$ 161.93	+Info
122965	MPF BRUSH Optimum (Nondas Design) Spring 3-in-1 Brush Kit- ORANGE (mens line	US\$ 150.59	+Info
122966	MPF BRUSH Optimum (Nondas Design) Spring 3-in-1 Brush Kit- Champagne (mens li...	US\$ 142.50	+Info

BRUSHES

Code	Name	Price	+
116488	MPF BRUSH Optimum (Nondas Design) Spring Brush Mens Line Size #4 1pk - Citrus ...	US\$ 81.82	+Info
116489	MPF BRUSH Optimum (Nondas Design) Spring Brush Mens Line Size #6 1pk - Champag...	US\$ 63.74	+Info
116490	MPF BRUSH Optimum (Nondas Design) Spring Brush Mens Line Size #8 1pk - Midnigh...	US\$ 89.59	+Info

OPTIMUM SPRING MASTER KITS

OPTIMUM MASTER KITS - MPF Brush Co

Designed by Nondas

M.P.F Optimum ceramic brush designed by Nondas, this truly unique design will revolutionize the dental ceramic brush industry ,with its balanced lightweight aluminum structure and ergonomic design, this brush handle is by far the greatest brush design ever to come onto the dental ceramic market, hence the name "optimum" .



With its streamlined curves, the M.P.F Optimum ceramic brush positions into the hand so well, this allows for total control of your build up technique.

Attention to detail such as the metal point on the end of the handle so that you can punch your dowel pins from the base of your models without having to change tools. The posterior shaft with its perforated design acts as a mini hammer for condensing your ceramic.

The M.P.F Brush Company Philosophy is to insist on high quality brush hair. Our high quality Kolinsky hair being a natural product is far superior to any synthetic hair the market has to offer.

Kit Contains:

- 1 Ceramic Brush, Stain
- 1 Ceramic Brush, Opaque
- 1 Ceramic Brush, Contour
- 1 Spring Ceramic Brush in your choice

#UsddMpf

Code	Name	Price	+
116502	MPF BRUSH Optimum (Nondas Design) Master Kit- size #4 (mens line) Mfg.#102-5...	US\$ 148.14	+Info
116503	MPF BRUSH Optimum (Nondas Design) Master Kit- size #6 (mens line) Mfg.#102-5...	US\$ 151.64	+Info
116504	MPF BRUSH Optimum (Nondas Design) Master Kit- size #8 (mens line) Mfg.#102-5...	US\$ 179.03	+Info

THOMAS SINGH BRUSHES

THOMAS SINGH BRUSHES - MPF Brush Co

The M.P.F. Stay Close to Nature Thomas Singh Signature features a smaller hair bundle which allows for the highly accurate and controlled placement of white or pink ceramic providing the esthetic excellence demands. Comfortable handle with flick innovation.

Machined out of premium, lightweight aluminum, the handle features spring activation for a super sharp working brush tip. The red handle color was chosen to exemplify confidence and boldness while in the hand of the Smile Maker. The replaceable brush tip system offers a nice finishing touch.



#UsddMpf

Code	Name	Price	+
116563	MPF BRUSH Thomas Singh Signature Brush 1P.K. Mfg.#118-2000 1182000 - #MPF BRUSH	US\$ 40.19	+Info
116564	MPF BRUSH Thomas Singh Signature Replacement Brush Tips (2/P.K.) Mfg.#118-2001 ...	US\$ 42.65	+Info

NONDAS DESIGN | OPTIMUM MASTER BRUSHES Non-Spring

OPTIMUM MASTER BRUSHES - MPF Brush Co

Designed by Nondas.

M.P.F Optimum ceramic brush non-Spring designed by Nondas, this truly unique design will revolutionize the dental ceramic brush industry ,with its balanced lightweight aluminum structure and ergonomic design, this brush handle is by far the greatest brush design ever to come onto the dental ceramic market, hence the name "optimum" .

With its streamlined curves, the M.P.F Optimum ceramic brush positions into the hand so well, this allows for total control of your build up technique. Attention to detail such as the metal point on the end of the handle so that you can punch your dowel pins from the base of your models without having to change tools. The posterior shaft with its perforated design acts as a mini hammer for condensing your ceramic.

The M.P.F Brush Company Philosophy is to insist on high quality brush hair. Our high quality Kolinsky hair being a natural product is far superior to any synthetic hair the market has to offer.



#UsddMpf

Code	Name	Price	+
116491	MPF BRUSH Optimum (Nondas Design) Stain Brush 1pk Mfg.#102-0009 1020009 - Brus...	US\$ 26.41	+Info
116492	MPF BRUSH Optimum (Nondas Design) Opaque Brush 1pk Mfg.#102-0010 1020010 - #MP...	US\$ 36.76	+Info
116493	MPF BRUSH Optimum (Nondas Design) Contour Brush 1pk Mfg.#102-0011 1020011 - #M...	US\$ 25.00	+Info

JUNGO ENDO BRUSHES

JUNGO ENDO BRUSHES - MPF Brush Co

The M.P.F. Stay Close to Nature Jungo Endo Signature features a smaller hair bundle which allows for the highly accurate and controlled placement of white or pink ceramic providing the esthetic excellence Jungo demands. Comfortable handle with flick innovation

Machined out of premium, lightweight aluminum, the handle features spring activation for a super sharp working brush tip. The red handle color was chosen to exemplify confidence and boldness while in the hand of the Smile Maker. The replaceable brush tip system offers a nice finishing touch.



#UsddMpf

Code	Name	Price	+
116561	MPF BRUSH Jungo Endo Signature Brush- RED 1pk Mfg.#118-1000 1181000 - #MPF BRUSH co	US\$ 35.48	+Info
122970	MPF BRUSH Jungo Endo Signature Brush- BLACK 1pk Mfg.#118-0000 1180000 - #MPF BRU...	US\$ 35.48	+Info
116562	MPF BRUSH Jungo Endo Singature Replacement Brush Tips (2/pk) Mfg.#118-1001 11810...	US\$ 36.63	+Info
124541	MPF BRUSH Jungo Endo Signature Brush- Grey1pk Mfg.#118-1500 1181500 - #MPF BRUSH co	US\$ 35.48	+Info

NON-SPRING BRUSHES

NON-SPRING BRUSHES - MPF Brush Co

The latest line of veneering kolinsky tools from European brush maker M.P.F. Kolinsky hair - the real difference 5 Handle colors , you pick

Highlights

- * Latest brush line from brush maker M.P.F.
- * Fully weighted handle with ergonomic curvature
- * Replaceable kolinsky tip system.
- * Finest kolinsky hair in the world.
- * Exclusive hair treatment & ferrule bonding process.
- * Designed by master ceramist Nondas Vlachapoulos.



Available Colors: Blue, orange, green, black and red.

#UsddMpf

Code	Name	Price	+
122964	MPF BRUSH Original Non-Spring 3-in-1 Brush Kit - BLUE Mfg.#200-5001 2005001 - #M...	US\$ 131.35	+Info
122963	MPF BRUSH Original Non-Spring 3-in-1 Brush Kit - ORANGE Mfg.#200-4001 2004001 - ...	US\$ 131.35	+Info
122962	MPF BRUSH Original Non-Spring 3-in-1 Brush Kit - GREEN Mfg.#200-3001 2003001 - #...	US\$ 141.56	+Info
122961	MPF BRUSH Original Non-Spring 3-in-1 Brush Kit - BLACK Mfg.#200-2001 2002001 - #...	US\$ 143.02	+Info
122960	MPF BRUSH Original Non-Spring 3-in-1 Brush Kit - RED Mfg.#200-1001 2001001 - #MP...	US\$ 128.43	+Info

REVOLUTION BRUSHES

REVOLUTION BRUSHES - MPF BRUSH Co

Light-weight, ergonomic aluminium handle.
Replaceable premium Kolinsky sable brushtips.



Brush + Tips KITS

Code	Name	Price	+
400088	MPF BRUSH MPF REVOLUTION BRUSH KIT SAME 3 TIPS size #2 MIDNIGHT BLACK 200-3002	US\$ 69.63	+Info
400090	MPF BRUSH REVOLUTION Brush Tips #4 3pk #200-3004	US\$ 90.61	+Info
400091	MPF BRUSH REVOLUTION Brush Tips #6 3pk #200-3006	US\$ 125.06	+Info
400092	MPF BRUSH REVOLUTION Brush Tips #8 3pk #200-3008	US\$ 166.99	+Info

Brushes

Code	Name	Price	+
400077	MPF BRUSH REVOLUTION Brush #2 Midnight BLACK #200-0002	US\$ 34.34	+Info
400078	MPF BRUSH REVOLUTION Brush #4 Silver Moon #200-0004	US\$ 64.98	+Info
400079	MPF BRUSH REVOLUTION Brush #6 Lead Grey #200-0006	US\$ 48.11	+Info
400080	MPF BRUSH REVOLUTION Brush #8 Sea BLUE #200-0008	US\$ 100.06	+Info
400081	MPF BRUSH REVOLUTION Brush Stain #00 #200-0009	US\$ 23.99	+Info
400082	MPF BRUSH REVOLUTION Brush Glaze #200-0010	US\$ 22.55	+Info
400083	MPF BRUSH REVOLUTION Brush Contour #200-0011	US\$ 23.99	+Info

REVOLUTION MASTER BRUSH

REVOLUTION MASTER BRUSH - MPF BRUSH Co

Light-weight, ergonomic aluminium handle.
Replaceable premium Kolinsky sable brushtips.



Code	Name	Price	+
400084	REVOLUTION Master Brush KIT #2 #200-2002	US\$ 110.83	+Info
400085	REVOLUTION Master Brush KIT #4 #200-2004	US\$ 117.02	+Info
400086	REVOLUTION Master Brush KIT #6 #200-2006	US\$ 126.90	+Info
400087	REVOLUTION Master Brush KIT #8 #200-2008	US\$ 146.99	+Info

CLASSIC KOLINSKY BRUSHES

CLASSIC KOLINSKY BRUSHES - MPF Brush Co

Classic Kolinsky Brushes M.P.F. are produced using the highest quality kolinsky hair.

Inside these brush heads are several metallic beads, which act as a condenser when shaken.

This action sends the moisture trapped within the brush head, directly to the brush tip, giving a very sharp point to work with.

The handles are made from wood, varnished in a matt blue colour.



#UsddMpf

Kit

Code	Name	Price	+
124255	MPF BRUSH Classical Kolinsky Kit 4 pcs - #105-0001	US\$ 158.74	+Info

BRUSHES

Code	Name	Price	+
400003	MPF BRUSH Classic Brush 1pc Size #2 Premium kolinsky sable brushtips. New plasti...	US\$ 30.06	+Info
122969	MPF BRUSH Classic Stain Brush Mfg #105-0009 1050009 - # (handle BLUE or BLACK m...	US\$ 12.64	+Info
122974	MPF BRUSH Classic Kolinsky Ceramic Opaque Brush Mfg #105-0010 1050010 - # (han...	US\$ 15.50	+Info
122975	MPF BRUSH Classic Kolinsky Ceramic Contour Brush Mfg #105-0011 1050011 - # (han...	US\$ 15.50	+Info
122971	MPF BRUSH Classic Kolinsky Ceramic Brush #4 Mfg #105-0004 1050004 - # (handle B...	US\$ 34.34	+Info
122972	MPF BRUSH Classic Kolinsky Ceramic Brush #6 Mfg #105-0006 1050006 - # (handle B...	US\$ 51.57	+Info
122973	MPF BRUSH Classic Kolinsky Ceramic Brush #8 Mfg #105-0008 1050008 - # (handle B...	US\$ 68.51	+Info
122976	MPF BRUSH Classic Kolinsky Ceramic Occlusal Brush Mfg #105-0012 1050012 - # (ha...	US\$ 15.34	+Info

SYNTHESIS BRUSHES

SYNTHESIS BRUSHES - MPF Brush Co

The latest brush in our range, a practical, affordable addition to your ceramic build up department, the MPF Synthesis brush fuses both the humid properties of Kolinsky and the maneuverability of synthetic hair, giving you the best of both elements.

#UsddMpf



Code	Name	Price	+
123906	MPF BRUSH -Synthesis Brush- 4 in 1 Kit - #105-0000 - #MPF BRUSH co	US\$ 46.96	+Info
122994	MPF BRUSH Synthesis Brush #2 Mfg.#105-2000 1052000 - #MPF BRUSH co	US\$ 22.86	+Info

122995	MPF BRUSH Synthesis Brush #4 Mfg.#105-4000 1054000 - #MPF BRUSH co	US\$ 24.12	+Info
122996	MPF BRUSH Synthesis Brush #6 Mfg.#105-6000 1056000 - #MPF BRUSH co	US\$ 33.59	+Info
122997	MPF BRUSH Synthesis Brush #8 Mfg.#105-8000 1058000 - #MPF BRUSH co	US\$ 40.77	+Info

ZIRCONIA STAINING BRUSHES

ZIRCONIA STAINING BRUSH - MPF Brush Co

Zirconia brush kit - Brush for coloring of zirconia before sintering.



Code	Name	Price	+
123438	MPF BRUSH Zirconia Staining Brush Kit 3 Brushes Kit - Mfg.#105-1000 1051000 - #...	US\$ 45.94	+Info
421148	MPF BRUSH Zirconia Staining Brush size #1 Brushes Kit - Mfg.#105-1001 - #MPF BRU...	US\$ 15.50	+Info
421149	MPF BRUSH Zirconia Staining Brush size #2 - Mfg.#105-1002 - #MPF BRUSH co	US\$ 17.23	+Info
421150	MPF BRUSH Zirconia Staining Brush size #3 Brushes Kit - Mfg.#105-1003 - #MPF BRU...	US\$ 18.38	+Info

PORCELAIN BRUSHES | SMILE LINE



REPLACEMENTS TIPS | SMILE LINE



Other Models

Other Models

Code	Name	Price	+
123623	RSPCT Aluminum handle only	US\$ 25.00	+Info
123774	Compo Brush #000 Style Italian	US\$ 70.00	+Info
123776	#000 Brush tips Posterior 50pk	US\$ 83.13	+Info
123842	FlexMount Connector 2pk	US\$ 40.00	+Info
123846	N.era #8 w/ Flexible Connector	US\$ 94.28	+Info

REPLACEMENT TIPS | CLEAR HANDLE COLLECTION

Code	Name	Price	+
123607	SMILE LINE * Spare tips #3/0 / 5pcs 1951	US\$ 70.00	+Info
123605	SMILE LINE * Spare tips #1 / 5pcs 1941	US\$ 65.00	+Info
123603	SMILE LINE * Spare tips #4 / 3pcs 1931	US\$ 90.00	+Info
123601	SMILE LINE * Spare tips #6 / 2pcs 1921	US\$ 130.00	+Info
123598	SMILE LINE * Spare tips #8-LT / 2pcs 1911-A	US\$ 81.69	+Info
123596	SMILE LINE * Spare tips #8-R / 2pcs 1911	US\$ 147.26	+Info
123609	SMILE LINE * Spare tip for Dust brush / 1pc 1961	US\$ 38.75	+Info

REPLACEMENT TIPS | KOLINSKY BRUSHES

Code	Name	Price	+
123722	SMILE LINE Brush tip #8-R module - Kolinsky hair - stainless steel connection 41...	US\$ 89.10	+Info
123721	SMILE LINE Brush tip #8-LT module - Kolinsky hair - stainless steel connection 4...	US\$ 47.99	+Info

REPLACEMENT TIPS | N ERA BRUSH

Code	Name	Price	+
123757	SMILE LINE N.era brush tips #1 / 2pcs 4201-1	US\$ 86.62	+Info
123758	SMILE LINE N.era brush tips #4 / 2pcs 4201-4	US\$ 102.24	+Info
123759	SMILE LINE N.era brush tips #6 / 2pcs 4201-6	US\$ 119.28	+Info
123760	SMILE LINE N.era brush tips #8 / 2pcs 4201-8	US\$ 136.32	+Info
123761	SMILE LINE N.era build-up brush #8- Short Ferrule - Spare tips (2pcs) 4201SF-8	US\$ 126.45	+Info

REPLACEMENT TIPS | PASTE OPAQUE BRUSHES

Each combination requires more than fifty operations before being carefully controlled and packed.
The mechanical parts are manufactured in Switzerland by specialists in watch and machinery making and all tips are of the highest quality and of certified European production.



Code	Name	Price	+
123644	SMILE LINE Paste opaque brush tip module- large 4000-LO	US\$ 45.00	+Info
123645	SMILE LINE Paste opaque brush tip module- small 4000-M	US\$ 38.75	+Info
123678	SMILE LINE Paste opaque brush tips refill (large) / 5pcs 4001-LO	US\$ 45.94	+Info
123679	SMILE LINE Paste opaque brush tips refill / 5pcs 4001-M	US\$ 43.13	+Info

REPLACEMENT TIPS | SYNTEC

Code	Name	Price	+
------	------	-------	---

123689	SMILE LINE SynTec brush tip #0 refill / 5pcs 4001-Sy0	US\$ 43.75	+Info
123690	SMILE LINE SynTec brush tip #3/0 refill / 5pcs 4001-Sy3/0	US\$ 43.75	+Info
123667	SMILE LINE SynTec brush tip #0 module 4000-Sy0	US\$ 28.13	+Info
123668	SMILE LINE SynTec brush tip #3/0 module 4000-Sy3/0	US\$ 28.13	+Info

REPLACEMENT TIPS | MODULE

Code	Name	Price	+
123673	SMILE LINE Brush tip #4 module 4000-W	US\$ 51.88	+Info
123672	SMILE LINE Dust brush module 4000-VD	US\$ 58.13	+Info
123671	SMILE LINE Brush tip #6 module 4000-V6	US\$ 76.25	+Info
123643	SMILE LINE Brush tip #5 module 4000-LC	US\$ 45.94	+Info
123648	SMILE LINE Stain brush tip #000 module 4000-O	US\$ 35.94	+Info

REPLACEMENT TIPS | STANDARD

REPLACEMENT TIPS | STANDARD

Code	Name	Price	+
123692	SMILE LINE Brush tips #4 refill / 3pcs 4001-W	US\$ 90.00	+Info
123677	SMILE LINE Brush tips #5 refill / 3pcs 4001-LC	US\$ 65.00	+Info
123680	SMILE LINE Brush tips #2 refill / 5pcs - 4001-N	US\$ 65.00	+Info
123681	SMILE LINE Brush tips #3/0 refill / 5pcs 4001-O	US\$ 55.94	+Info

CLEAR HANDLE COLLECTION - SMILE LINE

SMILE LINE - CLEAR HANDLE COLLECTION KOLINSKY BRUSHES

Thousands of ceramists throughout the world have been faithful to this iconic product line. The collection of brushes ?Clear Handle? is made from first choice and severely controlled Kolinsky sable hair. This offers you the confidence of the best possible quality: a **very sharp and stable point** as well as a **superior hair resilience**. With a concern for durability and rationalization, Smile Line has developed all sizes with **exchangeable brush tips**.

ADAPTED TO YOUR WORKING STYLE

The shape of your brush tip matches your build-up technique or style and this is the reason why we propose not one but two different shapes of #8 brush tips: the ?Regular? one, which is the best seller, and the ?Long Tip? for who prefers a tip even longer and sharper.

The Cut-Back Instrument of the ?Clear Handle? collection allows you all options of reduction, sculpture and separation of unfired porcelain. It automatically comes along with its assortment of x-thin tools.



CLEAR HANDLE COLLECTION - SMILE LINE


Code	Name	Price	+
123594	SMILE LINE CLEAR Handle brushes- assortment of 6 1900	US\$ 430.00	+Info
123606	SMILE LINE * CLEAR Handle- brush #000 1950	US\$ 43.00	+Info
123604	SMILE LINE * CLEAR Handle- brush #1 1940	US\$ 45.00	+Info
123602	SMILE LINE * CLEAR Handle- brush #4 1930	US\$ 55.00	+Info
123600	SMILE LINE * CLEAR Handle- brush #6 1920	US\$ 106.25	+Info
123597	SMILE LINE CLEAR Handle- brush #8-LT 1910-A	US\$ 71.12	+Info
123595	SMILE LINE * CLEAR Handle- brush #8-R 1910	US\$ 130.00	+Info
123608	SMILE LINE * CLEAR Handle- Dust brush 1960	US\$ 70.00	+Info

COMPO BRUSH - SMILE LINE

COMPO BRUSH - SMILE LINE

- * Smile Line's black synthetic fiber allows superior accuracy with the best contrast
- * # 000 Rear: you can easily fold them for ergonomic access to the back areas
- * Extend its life: clean with acetone
- * Brush tips for indirect use

Available: Flat - Set



Code	Name	Price	+
123773	SMILE LINE *Compo Brush Flat Brush Style Italiano for anteriors complete with t...	US\$ 70.00	+Info

COMPO BRUSH - SMILE LINE

Smile Line Replacement Tips - SMILE LINE

Code	Name	Price	+
466408	Compo Brush Set	US\$ 103.12	+Info

KOLINSKY BRUSHES - SMILE LINE

KOLINSKY BRUSHES - SMILE LINE

Features:
Smile Line introduces ?FlexMount?, a flexible rubber connector for combining the handle and the tool. This new concept is adapted to Instrument by Smile Line, more specifically to our ZR mixing spatula and our Porcelain Brushes.



BUILD-UP

Code	Name	Price	+
123660	SMILE LINE Kolinsky build-up brush #6- Short Ferrule - Spare tip 4000-SF-6	US\$ 70.18	+Info
123662	SMILE LINE Kolinsky build-up brush #8LT- Short Ferrule - Spare tip 4000-SF-8LT	US\$ 85.47	+Info

FLEXMOUNT

Code	Name	Price	+
123843	SMILE LINE Kolinsky Brush #6 w/FlexMount Connector - Module only FL-4020SF-6	US\$ 80.00	+Info
123845	SMILE LINE -Kolinsky Brush #8LT w/FlexMount Connector - module tip only FL-4020S...	US\$ 98.10	+Info
123844	SMILE LINE *Kolinsky Brush #8R w/FlexMount Connector - Module tip only FL-4020SF-8R	US\$ 98.10	+Info

Other Models

Other Models

Code	Name	Price	+
473479	RSPCT, spare brush tips #6 Regular	US\$ 129.36	+Info
473462	Njoy Nera Brush #6 Sex on/the Beach	US\$ 69.09	+Info
473463	Njoy N.era Brush #4 Sex on the Beac	US\$ 58.80	+Info
473464	Njoy N.era Brush #1 Sex on the Beac	US\$ 50.35	+Info
473466	N.era Build-up brush #8, Short Ferr	US\$ 147.00	+Info
473470	RSPCT, assortment of 1 complete bru	US\$ 254.31	+Info
473471	RSPCT, Complete Brush #8R Regular	US\$ 109.37	+Info
473472	RSPCT, Complete Brush #8LT Long Tip	US\$ 109.37	+Info
473473	RSPCT, Complete Brush #6R Regular	US\$ 94.08	+Info
473474	RSPCT, Complete Brush #6LT Long Tip	US\$ 94.08	+Info
473475	RSPCT - Complete Brush #3	US\$ 82.32	+Info
473476	RSPCT, Complete Brush #1	US\$ 76.44	+Info
473477	RSPCT, spare brush tips #8 Regular	US\$ 141.12	+Info
473478	RSPCT, spare brush tips #8 Long tip	US\$ 141.12	+Info
473461	Njoy N.era Brush #8 Sex on the Beac	US\$ 77.54	+Info
473480	RSPCT, spare brush tips #8 Long Tip	US\$ 129.36	+Info
473481	RSPCT, spare brush tips #3 (2pcs)	US\$ 105.84	+Info
473482	RSPCT, spare brush tips #1 (2pcs)	US\$ 88.20	+Info
473524	Kolinsky Build-up brush #8R, Short	US\$ 117.60	+Info
473525	Kolinsky Build-up brush #8LT, Short	US\$ 124.58	+Info
473526	N.era Build-up brush #8R, Short Fer	US\$ 105.84	+Info
473527	Instrument Set of 2 (Kolinsky brush	US\$ 282.24	+Info
473528	Instrument Set of 2 (Kolinsky brush	US\$ 294.00	+Info
473529	Instrument Set of 2 (N.era brush #8	US\$ 294.00	+Info
473530	Instrument Set of 2 (Kolinsky brush	US\$ 311.64	+Info
473531	Instrument Set of 2 (Kolinsky brush	US\$ 329.28	+Info
473532	Instrument Set of 2 (N.era brush #8	US\$ 317.52	+Info
473448	N.era bRRRush Assortment of 4	US\$ 352.80	+Info
473367	Gingiva surface texture brush set (	US\$ 170.52	+Info
473382	Clear Handle Brush #8R Complete + t	US\$ 164.20	+Info
473387	N.era bRRRush handle only	US\$ 44.10	+Info
473391	Gingiva surface texture hard Brush	US\$ 38.96	+Info
473392	Gingiva surface texture soft Brush	US\$ 38.96	+Info
473397	Kolinsky Build-up brush #6, Short F	US\$ 78.65	+Info
473398	Kolinsky Build-up brush #8R, Short	US\$ 105.84	+Info
473399	Kolinsky Build-up brush #8LT, Short	US\$ 105.84	+Info
473444	N.era bRRRush #8	US\$ 117.60	+Info
473445	N.era bRRRush #6	US\$ 105.84	+Info
473446	N.era bRRRush #4	US\$ 94.08	+Info
473447	N.era bRRRush #1	US\$ 82.32	+Info
123610	Clear Handle Cut Back w/Blades	US\$ 86.25	+Info
473449	Njoy N.era Brush #8 Mojito	US\$ 77.54	+Info
473450	Njoy N.era Brush #6 Mojito	US\$ 69.09	+Info
473451	Njoy N.era Brush #4 Mojito	US\$ 58.80	+Info
473452	Njoy N.era Brush #1 Mojito	US\$ 50.35	+Info

473453	Njoy N.era Brush #8 Yellow Cab	US\$ 77.54	+Info
473454	Njoy N.era Brush #6 Yellow Cab	US\$ 69.09	+Info
473455	Njoy N.era Brush #4 Yellow Cab	US\$ 58.80	+Info
473456	Njoy N.era Brush #1 Yellow Cab	US\$ 50.35	+Info
473457	Njoy N.era Brush #8 Cosmopolitan	US\$ 77.54	+Info
473458	Njoy N.era Brush #6 Cosmopolitan	US\$ 69.09	+Info
473459	Njoy N.era Brush #4 Cosmopolitan	US\$ 58.80	+Info
473460	Njoy N.era Brush #1 Cosmopolitan	US\$ 50.35	+Info

N.ERA BRRRUSH - SMILE LINE


SMILE LINE - N.ERA BRRRUSH THE VIBRATION THAT MAKES A DIFFERENCE

N.era bRRRush is an idea of Roberto Devolio, MDT in Curitiba, Brasil. While bRRRush is using the revolutionary N.era tips, a unique vibe zone has been machined on its handle.

When letting the nail of your forefinger stroke back and forth on the surface, you create a **micro-vibration which propagates up to the very end of your brush tip**. The small portion of porcelain on your brush will almost detach while transforming itself into a small drop very easy to apply to your build-up. The short video clip below shows you how impressive it is.

The advantage of this technique lies in the fact that **much less air is integrated in the porcelain during build-up**, synonymous with less porosities in the porcelain after firing - spectacular especially in the incisal area.

N.era bRRRushes are available as a set of four sizes, n 8, 6, 4 and 1 and can of course be ordered also as single instruments. The replacement tips are the same as N.era tips. It is also possible to order the handle of bRRRush alone, in case you wish to mount a traditional Kolinsky hair brush tip for example.



Code	Name	Price	+
123732	SMILE LINE N.era bRRRush- set of 4 4200-BRRR-Asst (Vibrating Porcelain Brush)	US\$ 356.59	+Info
123728	SMILE LINE N.era Brush #1 4200-BRRR-1 (Vibrating Porcelain Brush)	US\$ 76.82	+Info
123729	SMILE LINE bRRRush #4 4200-BRRR-4 (Vibrating Porcelain Brush)	US\$ 89.10	+Info
123730	SMILE LINE bRRRush #6 4200-BRRR-6 (Vibrating Porcelain Brush)	US\$ 102.35	+Info
123731	SMILE LINE bRRRush #8 4200-BRRR-8 (Vibrating Porcelain Brush)	US\$ 117.78	+Info

N.ERA BRUSH - SMILE LINE

SMILE LINE - N.ERA BRUSH

N.era brushes are manufactured out of a space-age filament, never seen yet in the dental industry. It is a unique polyimide fiber that undergoes several special chemical treatments.

The final properties (resilience, water absorption, stability at the tip) are better than the highest quality of Kolinsky hair. N.era is far superior because it last months longer.

Not only is the N.era the highest quality, but it is also constant. Because of the material's nature and production process, the N.era brush batch #1 is exactly the same as batch #999.



BRUSHES

Code	Name	Price	+
123727	SMILE LINE N.era brush- assortment of 4 4200	US\$ 270.00	+Info
123726	SMILE LINE N.era brush #8- Macchiato 4200-8	US\$ 102.24	+Info
123725	SMILE LINE N.era brush #6- Cappuccino 4200-6	US\$ 76.34	+Info
123724	SMILE LINE N.era brush #4- Moka 4200-4	US\$ 79.52	+Info
123723	SMILE LINE N.era brush #1- Ristretto 4200-1	US\$ 68.16	+Info

N.ERA NJOY BRUSHES - SMILE LINE

SMILE LINE - N.ERA - NJOY - BRUSHES THE FEATHERLIGHT VERSION OF THE N.ERA BRUSH

On request of some customers, Smile Line has decided to complete the N.era assortment with **Njoy: colorful, full of joy and good mood? but also extra-light N.era brushes !**

[Njoy] handle is hollow, injected out of ABS-thermoplastic and overmoulded with a soft-touch silicone for a very nice grip.

YOUR CUSTOMIZED SET :

The concept [Njoy] allows you **to compose your own set of brushes according to your own tastes**. Enough to determine the color you wish to receive for each of the four sizes in the set. Your Smile Line partner will be glad to compose your personal and customized set ! Some possibilities are shown on the images above.

IMPORTANT :

N.era [Njoy] brushes have no interchangeable tips. Smile Line chose that option mainly in order



to diminish the weight of the instrument. Of course, this has a positive impact on its retail price while keeping all the advantages of N.era!

Bloody Mery

Code	Name	Price	+
123733	SMILE LINE N.era brush 4200-Njoy-BLO-1	US\$ 50.01	+Info
123734	SMILE LINE N.era brush 4200-Njoy-BLO-4	US\$ 50.01	+Info
123735	SMILE LINE N.era brush 4200-Njoy-BLO-6	US\$ 50.01	+Info
123736	SMILE LINE N.era brush 4200-Njoy-BLO-8	US\$ 50.01	+Info

Blue Lagoon

Code	Name	Price	+
123741	SMILE LINE N.era brush 4200-LAG-1	US\$ 50.01	+Info
123742	SMILE LINE N.era brush 4200-Njoy-LAG-4	US\$ 50.01	+Info
123743	SMILE LINE N.era brush 4200-Njoy-LAG-6	US\$ 50.01	+Info
123744	SMILE LINE N.era brush 4200-Njoy-LAG-8	US\$ 66.22	+Info

Cosmopolitan

Code	Name	Price	+
123737	SMILE LINE N.era brush 4200-Njoy-COS-1	US\$ 50.01	+Info
123738	SMILE LINE N.era brush 4200-Njoy-COS-4	US\$ 50.01	+Info
123739	SMILE LINE N.era brush 4200-COS-6	US\$ 50.01	+Info
123740	SMILE LINE N.era brush 4200-Njoy-COS-8	US\$ 50.01	+Info

Mojito

Code	Name	Price	+
123745	SMILE LINE N.era brush 4200-Njoy-MOJ-1	US\$ 50.01	+Info
123746	SMILE LINE N.era brush 4200-Njoy-MOJ-4	US\$ 50.01	+Info
123747	SMILE LINE N.era brush 4200-MOJ-6	US\$ 50.01	+Info
123748	SMILE LINE N.era brush 4200-MOJ-8	US\$ 50.01	+Info

S.O The Beach

Code	Name	Price	+
123749	SMILE LINE N.era brush 4200-SOB-1	US\$ 50.01	+Info
123750	SMILE LINE N.era brush 4200-SOB-4	US\$ 50.01	+Info
123751	SMILE LINE N.era brush 4200-SOB-6	US\$ 50.01	+Info
123752	SMILE LINE N.era brush 4200--Njoy.SOB-8	US\$ 50.01	+Info

Yellow Cab

Code	Name	Price	+
123753	SMILE LINE N.era brush 4200-YEC-1	US\$ 50.01	+Info
123754	SMILE LINE N.era brush 4200-YEC-4	US\$ 50.01	+Info
123755	SMILE LINE N.era brush 4200-YEC-6	US\$ 50.01	+Info
123756	SMILE LINE N.era brush 4200-YEC-8	US\$ 50.01	+Info

SYNTEC COLLECTION BRUSH - SMILE LINE SMILE LINE - SYNTEC COLLECTION BRUSH

The SYNTEC synthetic Collection offers clients the beauty and ancient elegance of Smile Line within a more affordable price range. Each piece in the Syntec Collection is set with only high quality synthetic gold fibers with a unique interlock technique. For those contemplating owning their first piece of art jewelry, the Classical Collection provides a gentle, yet grand step into the world of hand-crafted splendor. The line is most likely the best offering of synthetic brushes available to the Laboratory industry.

The fiber used is conical at its end, and this combined with the interlock construction technique allows for a consistent and stable brush tip. A huge advantage for the ceramist. The handle, made out of transparent acrylic rods is smooth, light and offers an excellent grip to its user.



Code	Name	Price	+
123587	SMILE LINE SynTec brushes- assortment of 6 1800	US\$ 60.91	+Info
123589	SMILE LINE SynTec brush #000 / 2pcs - 1800-3/0	US\$ 29.42	+Info
123588	SMILE LINE SynTec brush #0 / 2pcs 1800-0	US\$ 16.51	+Info
123590	SMILE LINE SynTec brush #2 / 2pcs 1800-2	US\$ 16.51	+Info
123591	SMILE LINE *SynTec brush #4 / 2pcs 1800-4	US\$ 17.29	+Info
123592	SMILE LINE *SynTec brush #6 / 2pcs 1800-6	US\$ 23.81	+Info
123593	SMILE LINE *SynTec brush #8 / 2pcs 1800-8	US\$ 30.99	+Info

PORCELAIN BRUSHES | RENFE RT



BASIC LINE BRUSHES - RENFERT

BASIC LINE BRUSHES - RENFERT

The Main Difference Between Basic Line And Other Low Cost Dental Brushes Is The Choice Of Material And Manufacturing Method.

Manufacturing process

- * Optimal length of bristles
- * Sorting of unsuitable bristles
- * Manual, careful depilation
- * Special adhesive technique
- * High-quality ferrule

Bristle quality & Standard Red Sable

Up to approx. 80% of low cost bristles used, which are usually meant for the cosmetic or painting market.

Manufacturing process

Bristles too short, no sorting out. Chemical acidification
Poor quality retention of bristles which are often too short. Collar with sharp edges.

The top bristle quality and the complex manufacturing process provide the Basic line brush with all the essential requirements for working with porcelain. The Basic line offers very high quality at a fair and attractive price.



Code	Name	Price	+
023-1717-0000	RENFERT - BASIC-LINE Brushes Set of 6-#1717-0000 #17170000	US\$ 79.86	+Info
023-1717-0010	RENFERT - BASIC-LINE Brushes size 0/1-2pk-#1717-0010 #17170010	US\$ 13.31	+Info
023-1717-0001	RENFERT - BASIC-LINE Brushes size #1-2pk-#1717-0001 #17170001	US\$ 16.10	+Info
023-1717-0002	RENFERT - BASIC-LINE Brushes size #2-2pk-#1717-0002 #17170002	US\$ 16.64	+Info
023-1717-0004	RENFERT - BASIC-LINE Brushes size #4-2pk-#1717-0004 #17170004 DISCONTINUED BY...	US\$ 20.94	+Info
023-1717-0006	RENFERT - BASIC-LINE Brushes size 6-2pk-#1717-0006 #17170006	US\$ 39.61	+Info
120685	RENFERT - Basic-line Brushes- Size 07 - #17170007 - 2pk	US\$ 29.79	+Info
120686	RENFERT - - Basic-line Brushes- Size 08 - # 17170008	US\$ 30.40	+Info

CERAMICUS BRUSHES - RENFERT

CERAMICUS BRUSHES - RENFERT

With properties superior to both types of brush known so far. This special material and the manufacturing process to create the Ceramicus finally provide.

- * What specialists in porcelain aesthetics really need:
- * A brush tip which remains perfectly formed.
- * Sufficient tension force for picking up larger amounts of porcelain.
- * Moistening the porcelain powder optimally.
- * Build up interdental regions without the use of a metal instrument.
- * Direct working procedures; no need to rinse continually, saves time and allows concentrated working.

The custom grain of the wooden handle makes every Ceramicus unique.

- * Ceramicus Size 01
- * Ceramicus Size 02
- * Ceramicus Size 04
- * Ceramicus Size 06
- * Ceramicus Size 08
- * Ceramicus Big Brush
- * Ceramicus Size 1/0
- * Ceramicus Size 2/0
- * Ceramicus Opaquing brush



Code	Name	Price	+
120684	RENFERT - Ceramicus Staining Brush Set 00/000 - #17030000	US\$ 47.89	+Info
023-1716-0000	RENFERT - Ceramicus Brush Set of 6-#171-60000 #17160000 (#0 #1 #2 #4 #6 #8)	US\$ 112.36	+Info
023-1716-0120	RENFERT - Ceramicus Brush Size #00 2pcs-#1716-0120 #17160120	US\$ 32.93	+Info
023-1716-0110	RENFERT - Ceramicus Brush Size #0-2pcs-#1716-0110 #17160110	US\$ 32.57	+Info
023-1716-0001	RENFERT - Ceramicus Brush Size #01 2pcs-#1716-0001 #17160001	US\$ 27.93	+Info
023-1716-0002	RENFERT - Ceramicus Brush Size #02 2pcs-#1716-0002 #17160002	US\$ 30.86	+Info
023-1716-0004	RENFERT - Ceramicus Brush Size #04 2pcs-#1716-0004 #17160004	US\$ 31.92	+Info
023-1716-0006	RENFERT - Ceramicus Brush Size #06 Each-#1716-0006 #17160006	US\$ 34.05	+Info
023-1716-0008	RENFERT - Ceramicus Brush Size #08 Each-#1716-0008 #17160008	US\$ 46.82	+Info
023-1716-0088	RENFERT - Ceramicus Big Brush Each-(#1716-0088) #17160088	US\$ 74.48	+Info
023-1716-0200	RENFERT - Ceramicus Brush Opaque 2pcs-(#1716-0200) #17160200	US\$ 43.63	+Info

GENIUS PORCELAIN BRUSHES - RENFE

RT

GENIUS PORCELAIN BRUSHES - RENFERT

Function & Performance

High-quality porcelain brush with a multifunctional handle. The benefits, function and design set new standards. The natural hair consists of 1A Kolinsky quality, which means an ideal combination of stability and elasticity for porcelain processing. A fine brush tip is formed due to thick hair.

Advantages

- * The multifunction handle protects the brush tip.
- * Optimum storage position.
- * Non-slip texture.
- * Comfortable handling

The protection mechanism

- * Hood for the brush tip during work breaks.
- Protective sleeve in the handle is slipped easily over the tip.
- * Ideal storage position - standing on its tip. This prevents porcelain deposits and splaying of tip bristles.



Code	Name	Price	+
023-1715-0000	RENFERT - Genius Brushes (4 Pcs)-SeT - #1715-0000 #17150000 High quality porcel...	US\$ 186.60	+Info
023-1715-0002	RENFERT - Genius Brush-size #2-Each- #1715-0002 #17150002 High quality porcelai...	US\$ 39.90	+Info
023-1715-0004	RENFERT - Genius Brush-size #4-Each- #1715-0004 #17150004 High quality porcelai...	US\$ 58.53	+Info
023-1715-0006	RENFERT - Genius Brush-Size 6-Each- #1715-0006 #17150006 High quality porcelain...	US\$ 73.42	+Info
023-1715-0008	RENFERT - Genius Brush-Size 8-Each- #1715-0008 #17150008 High quality porcelain...	US\$ 104.11	+Info

KOLINSKY BRUSHES - RENFERT

KOLINSKY - NATURAL BRISTLE BRUSHES FOR PORCELAIN - RENFERT

Natural bristle brush made from Kolinsky red sable bristles - quality and aesthetics combined.

Advantages

- * Fine, flexible brush tip.
- * Optimum liquid reservoir.
- * 6 different sizes for individual uses.

Brush forms Size 1 and 0/1 - stain brushes

Long, fine bristle tips for deeper spots in occlusal or interdental areas.

Size 2 - correction brush

Ideal for applying small amounts of porcelain portions after corrective firing.



Sizes 4 and 6

For porcelain modelling of large surfaces.

Size 8

Large liquid reservoir.

US DENTAL DEPOT is an Authorized RENFERT dealer in US.

#UsddRenfert

Code	Name	Price	+
023-1713-1000	RENFERT - Kolinsky Brushes (6 Pcs)-Set kiT - Natural bristle brush made from Ko...	US\$ 183.46	+Info
023-1713-1110	RENFERT - Kolinsky Brushes-Size 0/1-2 Pcs-Natural bristle brush made from Kolin...	US\$ 34.92	+Info
417576	RENFERT - LTD Edition V Musella build-up EACH - #17251500	US\$ 149.15	+Info
023-1713-1001	RENFERT - Kolinsky Brushes-size #1-2 Pcs-Natural bristle brush made from Kolins...	US\$ 39.90	+Info
023-1713-1002	RENFERT - Kolinsky Brushes-size #2-2 Pcs-Natural bristle brush made from Kolins...	US\$ 53.20	+Info
023-1713-1004	RENFERT - Kolinsky Brushes-size #4-2 Pcs-Natural bristle brush made from Kolins...	US\$ 71.30	+Info
023-1713-1006	RENFERT - Kolinsky Brushes-Size 6-2 Pcs-Natural bristle brush made from Kolinsk...	US\$ 96.89	+Info
023-1713-1008	RENFERT - Kolinsky Brushes-Size 8-2 Pcs-Natural bristle brush made from Kolinsk...	US\$ 144.54	+Info

LAY-ART BRUSHES STYLE - RENFERT

LAY-ART BRUSHES STYLE - RENFERT

Lay:Art style natural hair brushes provide premium quality for the highest demands. An innovative manufacturing technique achieves maximum values in terms of tension, dimensional stability and durability of the brush tip. The ergonomic handle design with Soft Touch is optimally balanced and fits very well in the hand.

Advantages

- * Top quality with natural hair for consistent working comfort at the highest level.
- * 10 sizes with different shape concepts, perfectly coordinated to individual requirements.
- * The Big Brush technique with sizes 8, 6 and 4 increases the moisture reserve and tension.

The optimal addition

Size 2 for small amounts, opaque brush for uniform application and staining brush for subtle shade accents - available as a practical set or individually.

2:

Size 2 is ideal for the exact application of small amounts.

Opaque:



The perfect combination of thick and long hairs (7.5 mm /0.29"). It nestles on the framework surface and ensures an effective, uniform opaque and liner application.

Color:

The staining brush: High strength right into the tip for targeted application of perfect shade accents and specific details.

#UsddRenfert

SET

Code	Name	Price	+
1725-1000	LAYART Natural Brush Set 3pieces-#1725-1000 #17251000 - Includes - x1 #2 Brush -...	US\$ 82.57	+Info

BOLD

Code	Name	Price	+
1725-0018	LAYART Natural Brush-Style size #8 Bold-Sold Individually- #1725-0018 #17250018	US\$ 160.72	+Info
1725-0016	LAYART Natural Brush-Style size #6 Bold-Sold Individually- #1725-0016 #17250016	US\$ 109.84	+Info
1725-0014	LAYART Natural Brush-Style size #4 Bold-Sold Individually- #1725-0014 #17250014	US\$ 55.06	+Info

SLIM

Code	Name	Price	+
1725-0008	LAYART Natural Brush-Style size #8 Slim-Sold Individually- #1725-0008 #17250008	US\$ 160.72	+Info
1725-0006	LAYART Natural Brush-Style size #6 Slim-Sold Individually- #1725-0006 #17250006	US\$ 114.42	+Info
1725-0004	LAYART Natural Brush-Style size #4 Slim-Sold Individually- #1725-0004 #17250004	US\$ 55.06	+Info

CONE

Code	Name	Price	+
1725-0040	LAYART Natural Brush-Style size #4 Cone-Sold Individually- #1725-0040 #17250040	US\$ 72.35	+Info

HAired

Code	Name	Price	+
1725-0002	LAYART Natural Brush-Style size #2-Sold In Package of 2- #1725-0002 #17250002	US\$ 73.42	+Info

OPAKER

Code	Name	Price	+
1725-0001	LAYART Natural Brush Opaker-Sold In Package of 2- #1725-0001 #17250001	US\$ 56.87	+Info

STAIN

Code	Name	Price	+
1725-0000	LAYART Natural Stain/ Color Brush #00-Sold In Package of 2- #1725-0000 #17250000	US\$ 58.53	+Info

PROFI LINE BRUSHES - RENFERT

PROFI BRUSHES - RENFERT

High-End Natural Bristle Brushes

Function & Performance

High quality porcelain brush made from 1A quality Kolinsky natural bristles.

Advantages

- * Weight incorporated in shaft gives optimum balance.
- * Comfortable to hold due to ergonomic handle.

Staining brush

Very high degree of elasticity for ideal take-up of large and small amounts of porcelain material.

Glazing brush

Ideal for "rubbing in" glaze.

Cleaning brush

For condensing the model.

US DENTAL DEPOT IS AN AUTHORIZED RENFERT DEALER IN US.



Code	Name	Price	+
023-1711-0100	RENFERT - Profi Brushes (6 Pcs)-Set - High quality porcelain brush made from A1...	US\$ 93.88	+Info
023-1711-0004	RENFERT - Profi Brushes-size #4-2 Pcs-High quality porcelain brush made from A1...	US\$ 32.32	+Info
023-1711-0006	RENFERT - Profi Brush-Size 6-Each-High quality porcelain brush made from A1 qua...	US\$ 36.34	+Info
023-1711-0008	RENFERT - Profi Brush-Size 8-Each-High quality porcelain brush made from A1 qua...	US\$ 73.63	+Info
023-1711-0120	RENFERT - Profi Staining Brushes-#2 X 2 Pcs High quality porcelain brush made f...	US\$ 18.24	+Info
023-1711-0106	RENFERT - Profi Brushes-BlunT - 2 Pcs-High quality porcelain brush made from A1...	US\$ 18.10	+Info
23-1711-0200	RENFERT - Profi Opaque Brush SET/2pcs-High quality porcelain brush made from A1...	US\$ 17.92	+Info

Renfert Genius Brushes - RENFERT

Renfert Genius Brushes - RENFERT

Code	Name	Price	+
417571	Genius Refill #2 2pk	US\$ 39.90	+Info
417572	Genius Refill #4 2pk	US\$ 63.84	+Info
417573	Genius Refill #6 2pk	US\$ 92.77	+Info
417574	Genius Refill #8 2pk	US\$ 114.11	+Info

TAKANISHI BRUSHES - RENFERT

TAKANISHI BRUSHES - RENFERT

Function & Performance

Very high quality synthetic bristle brush which has been well proven for many years for the optimal porcelain build-up.

Advantages

- * The brush tip remains stable.
- * Highly flexible brush bristles.
- * Ideal for reducing porcelain.



Code	Name	Price	+
023-1714-0000	RENFERT - Takamishi Brushes-(6 Pcs)-SeT - (#1714-0000 #17140000)-	US\$ 71.13	+Info
023-1714-0120	RENFERT - Takamishi Brushes-Size 00-2 Pcs-(#1714-0120 #17140120)-	US\$ 26.61	+Info
023-1714-0110	RENFERT - Takamishi Brushes-Size 0-2 Pcs-(#1714-0110 #17140110)	US\$ 26.61	+Info
023-1714-0001	RENFERT - Takamishi Brushes-size #1-2 Pcs-(#1714-0001 #17140001)-	US\$ 23.41	+Info
023-1714-0002	RENFERT - Takamishi Brushes-size #2-2 Pcs-(#1714-0002 #17140002)-	US\$ 23.41	+Info
023-1714-0003	RENFERT - Takamishi Brushes-size #3-2 Pcs-(#1714-0003 #17140003)-	US\$ 23.41	+Info
023-1714-0004	RENFERT - Takamishi Brushes-Size #4-2 Pcs-(#1714-0004 #17140004)-	US\$ 28.74	+Info
023-1714-0005	RENFERT - Takamishi Brushes-Size #5-2 Pcs-(#1714-0005 #17140005)-	US\$ 30.86	+Info
023-1714-0006	RENFERT - Takamishi Brushes-Size 6-2 Pcs-(#1714-0006 #17140006)-	US\$ 37.25	+Info
023-1714-0008	RENFERT - Takamishi Brushes-Size #8-2 Pcs-(#1714-0008 #17140008)-	US\$ 46.82	+Info
023-1714-0012	RENFERT - Takamishi Brush-size #12-Each-(#1714-0012 #17140012)-	US\$ 63.52	+Info
023-1701-0000	RENFERT - Takamishi Brushes-2Ea-00/000-4 Pcs-(#1701-0000 #1701-0000)-	US\$ 47.89	+Info
023-1714-0200	RENFERT - Takamishi Opaque Brush-(#1714-0200 #17140200)-	US\$ 31.92	+Info

PORCELAIN BRUSHES | KEYSTONE



CERAMIST BRUSHES - CLEAR HANDLE

CERAMIST BRUSHES - CLEAR HANDLE - KEYSTONE

Features:

- * Keystone offers an ALL-NATURAL high quality Kolinsky sable bristles.
- * The elasticity of hair recovers its original shape after stroke.
- * These brushes give the ability to control a porcelain flow, while holding and dispensing material predictably.
- * Sable Bristles
- * Available Sizes: Big Brush (8, 6, 4, 2, 1, 0)



Code	Name	Price	+
100787	KEYSTONE Kolinsky Ceramist Brush, #0, Made with all natural high quality - #117...	US\$ 11.40	+Info
100789	KEYSTONE Kolinsky Ceramist Brush, #2, Made with all natural high quality - #117...	US\$ 11.66	+Info
100790	KEYSTONE Kolinsky Ceramist Brush, #4, Made with all natural high quality - #117...	US\$ 12.16	+Info
100791	KEYSTONE Kolinsky Ceramist Brush, #6, Made with all natural high quality - #117...	US\$ 16.01	+Info
100792	KEYSTONE Kolinsky Ceramist Brush, #8, Made with all natural high quality - #117...	US\$ 17.10	+Info
100793	KEYSTONE Kolinsky Ceramist Brush, Big Brush, Made with all natural high quality...	US\$ 37.86	+Info

KOLINSKY BRUSHES

KEYSTONE - KOLINSKY BRUSHES

These are legendary natural Kolinsky Brushes. These brushes have a small bead inside that helps manage the humidity on the brush. You'll see the difference in your work immediately after using Kolinsky Brushes.

Available Sizes: #000, #00, #0, #1, #2, #3, #4, #5, #6, #7, #8, #2 Flat, #4 Flat, Big Brush



Code	Name	Price	+
124701	KEYSTONE Kolinsky Brush, #000, Made from natural Kolinsky sable, single brush - ...	US\$ 5.42	+Info
124702	KEYSTONE Kolinsky Brush, #00, Made from natural Kolinsky sable, single brush - ...	US\$ 5.42	+Info
124703	KEYSTONE Kolinsky Brush, #0, Made from natural Kolinsky sable, single brush - #...	US\$ 5.42	+Info
124704	KEYSTONE Kolinsky Brush, #1, Made from natural Kolinsky sable, single brush - #...	US\$ 7.51	+Info
124705	KEYSTONE Kolinsky Brush, #2, Made from natural Kolinsky sable, single brush - #...	US\$ 8.55	+Info
124706	KEYSTONE Kolinsky Brush, #3, Made from natural Kolinsky sable, single brush - #...	US\$ 9.03	+Info
124707	KEYSTONE Kolinsky Brush, #4, Made from natural Kolinsky sable, single brush - #...	US\$ 7.13	+Info
124708	KEYSTONE Kolinsky Brush, #5, Made from natural Kolinsky sable, single brush - #...	US\$ 7.44	+Info
124709	KEYSTONE Kolinsky Brush, #6, Made from natural Kolinsky sable, single brush - #...	US\$ 14.36	+Info
124710	KEYSTONE Kolinsky Brush, #7, Made from natural Kolinsky sable, single brush - #...	US\$ 11.21	+Info
124711	KEYSTONE Kolinsky Brush, #8, Made from natural Kolinsky sable, single brush - #...	US\$ 12.69	+Info
124712	KEYSTONE Kolinsky Brush, #2 Flat, Made with all natural Kolinsky sable, Single ...	US\$ 5.61	+Info
124713	KEYSTONE Kolinsky Brush, #4 Flat, Made with all natural Kolinsky sable, Single ...	US\$ 6.18	+Info

RED SABLE BRUSHES

RED SABLE BRUSHES - KEYSTONE

- * Keystone offers an ALL-NATURAL high quality Red Sable bristles.
- * This hair is red in color and is sometimes mixed or passed off as Kolinsky.
- * There are many different grades of Red Sable, depending on the region of the world from which the animal comes.
- * Available Sizes: Right Angle, 00 Right Angle, 000 Right Angle, 2 Flat, 4 Flat, 6 Flat, 000, 00, 0, 1, 2, 3, 4, 5, 6, 7, 14 Whipper



Code	Name	Price	+
034-1170280	KEYSTONE RED Sable Brush #0, single brush - #1170280	US\$ 4.70	+Info
034-1170370	KEYSTONE RED Sable Brush #00 Right Angle, single brush - #1170370	US\$ 3.99	+Info
034-1170270	KEYSTONE RED Sable Brush #00, single brush - #1170270	US\$ 4.70	+Info
034-1170360	KEYSTONE RED Sable Brush #000 Right Angle, single brush - #1170360	US\$ 3.91	+Info
034-1170260	KEYSTONE RED Sable Brush #000, single brush - #1170260	US\$ 4.70	+Info
034-1170290	KEYSTONE RED Sable Brush #1, single brush - #1170290	US\$ 4.85	+Info
034-1170390	KEYSTONE RED Sable Brush #1, single brush - #1170290	US\$ 3.99	+Info
034-1170355	KEYSTONE RED Sable Brush #14 Whipper, single brush - #1170355	US\$ 15.08	+Info
034-1170410	KEYSTONE RED Sable Brush #2 Flat, single brush - #1170410	US\$ 4.84	+Info
034-1170300	KEYSTONE RED Sable Brush #2, single brush - #1170300	US\$ 5.23	+Info
034-1170310	KEYSTONE RED Sable Brush #3, single brush - #1170310	US\$ 5.89	+Info
034-1170420	KEYSTONE RED Sable Brush #4 Flat, single brush - #1170420	US\$ 7.13	+Info
034-1170320	KEYSTONE RED Sable Brush #4, single brush - #1170320	US\$ 6.56	+Info
034-1170330	KEYSTONE RED Sable Brush #5, single brush - #1170330	US\$ 10.07	+Info
034-1170430	KEYSTONE RED Sable Brush #6 Flat, single brush - #1170430	US\$ 5.89	+Info
034-1170340	KEYSTONE RED Sable Brush #6, single brush - #1170340	US\$ 10.26	+Info
034-1170350	KEYSTONE RED Sable Brush #7, single brush - #1170350	US\$ 10.74	+Info
034-1170770	KEYSTONE SYNTHETIC RED Sable Big Brush Straight Round Big Brush #8 - (with smal...	US\$ 20.81	+Info

SYNTHETIC CERAMIC BRUSHES

SYNTHETIC CERAMIST BRUSHES - KEYSTONE

Features:

- * These Synthetic Ceramist Brushes combine excellent quality with an affordable price.
- * All of these brushes are available for private label.
- * Minimum quantities are required.



SYNTHETIC BIG BRUSH

Code	Name	Price	+
------	------	-------	---

034-1170700

KEYSTONE #9 Big Brush - SYNTHETIC Brush - Small bead inside for vibration #1170700

US\$ 41.32

[+Info](#)

FLAT

Code	Name	Price	+
034-1171220	KEYSTONE Synthetic Ceramist Brush, #2 Flat, single brush - #1171220	US\$ 4.09	+Info
034-1171230	KEYSTONE Synthetic Ceramist Brush, #4 Flat, single brush - #1171230	US\$ 9.12	+Info
034-1171240	KEYSTONE Synthetic Ceramist Brush, #6 Flat, single brush - #1171240	US\$ 9.12	+Info

ROUND

Code	Name	Price	+
034-1171000	KEYSTONE Synthetic Ceramist Brush, #000, single brush - #1171000	US\$ 5.39	+Info
034-1171010	KEYSTONE Synthetic Ceramist Brush, #00, single brush - #1171010	US\$ 3.99	+Info
034-1171020	KEYSTONE Synthetic Ceramist Brush, #0, single brush - #1171020	US\$ 3.99	+Info
034-1171030	KEYSTONE Synthetic Ceramist Brush, #1, single brush - #1171030	US\$ 4.18	+Info
034-1171040	KEYSTONE Synthetic Ceramist Brush, #2, single brush - #1171040	US\$ 6.56	+Info
034-1171050	KEYSTONE Synthetic Ceramist Brush, #3, single brush - #1171050	US\$ 8.84	+Info
034-1171060	KEYSTONE Synthetic Ceramist Brush, #4, single brush - #1171060	US\$ 9.03	+Info
034-1171070	KEYSTONE Synthetic Ceramist Brush, #5, single brush - #1171070	US\$ 9.02	+Info
034-1171080	KEYSTONE Synthetic Ceramist Brush, #6, single brush - #1171080	US\$ 9.69	+Info
034-1171090	KEYSTONE Synthetic Ceramist Brush, #7, single brush - #1171090	US\$ 10.07	+Info
034-1171200	KEYSTONE Synthetic Ceramist Brush, #8, single brush - #1171200	US\$ 10.17	+Info
034-1171205	KEYSTONE SYNTHETIC Brush #10 - Straight Round #1171205 -	US\$ 14.19	+Info

PORCELAIN BRUSHES | HARVEST DENTAL



FUSION BRUSH KOLINSKY

Brushes

Brushes

Code	Name	Price	+
475508	Fusion Brush #2	US\$ 72.08	+Info
475509	Fusion Brush #4	US\$ 81.49	+Info
475510	Fusion Brush #6	US\$ 94.67	+Info
475511	Fusion Brush #8	US\$ 97.36	+Info
475512	Fusion Brush #Stain	US\$ 69.92	+Info
475513	Fusion Brush #Opaque/Glaze	US\$ 68.41	+Info
475520	Fusion Brush - 4in1 Kit	US\$ 230.34	+Info

BrushTips

BrushTips

Code	Name	Price	+
475514	Fusion Brushtip #2 - 2 pc.	US\$ 80.68	+Info
475515	Fusion Brushtip #4 - 2 pc.	US\$ 95.35	+Info
475516	Fusion Brushtip #6 - 2 pc.	US\$ 100.09	+Info
475517	Fusion Brushtip #8 - 2 pc.	US\$ 134.32	+Info
475519	Fusion Brush Tip Opaque/Glaze	US\$ 77.46	+Info

Master Kits / Set

Master Kits / Set

Code	Name	Price	+
475521	Fusion Brush - Master Kit #2	US\$ 168.58	+Info
475522	Fusion Brush - Master Kit #4	US\$ 173.09	+Info
475523	Fusion Brush - Master Kit #6	US\$ 183.92	+Info
475524	Fusion Brush - Master Kit #8	US\$ 203.73	+Info

NYLO BRUSH SYNTHETIC

Brushes

Brushes

Code	Name	Price	+
475538	NyLo Brush Synthetic #2 Build	US\$ 45.36	+Info
475539	NyLo Brush Synthetic #4 Build	US\$ 47.78	+Info
475540	NyLo Brush Synthetic #6 Build	US\$ 56.32	+Info
475541	NyLo Brush Synthetic #8 Build	US\$ 59.73	+Info
475542	NyLo Brush Synthetic Stain	US\$ 38.51	+Info
475543	NyLo Brush Synthetic Glaze	US\$ 44.65	+Info

BrushTips

BrushTips

Code	Name	Price	+
475544	NyLo Brush Synthetic #2 Tips 2/	US\$ 59.16	+Info
475545	NyLo Brush Synthetic #4 Tips 2/	US\$ 66.70	+Info

475546	NyLo Brush Synthetic #6 Tips 2/	US\$ 76.37	+Info
475547	NyLo Brush Synthetic #8 Tips 2/	US\$ 87.13	+Info
475548	NyLo Brush Synthetic Stain Tips	US\$ 41.95	+Info
475549	NyLo Brush Synthetic Glaze Tips	US\$ 58.09	+Info

Kits / Set

Kits / Set

Code	Name	Price	+
475550	NyLo Brush Synthetic 4 in 1 Kit	US\$ 141.54	+Info
475551	NyLo Brush Synthetic #2 Master K	US\$ 109.98	+Info
475552	NyLo Brush Synthetic #4 Master K	US\$ 123.51	+Info
475553	NyLo Brush Synthetic #6 Master K	US\$ 127.10	+Info
475554	NyLo Brush Synthetic #8 Master K	US\$ 133.42	+Info

PURA BRUSH KOLINSKY

Brushes

Brushes

Code	Name	Price	+
475491	Pura Brush #2	US\$ 74.23	+Info
475492	Pura Brush #4	US\$ 77.46	+Info
475493	Pura Brush #6	US\$ 115.12	+Info
475494	Pura Brush #8	US\$ 130.49	+Info
475495	Pura Brush #Stain	US\$ 58.09	+Info
475496	Pura Brush #Opaque/Glaze	US\$ 68.85	+Info

Brushtips

Brushtips

Code	Name	Price	+
475497	Pura Brushtip #2 - 2 pc.	US\$ 92.51	+Info
475498	Pura Brushtip #4 - 2 pc.	US\$ 104.36	+Info
475499	Pura Brushtip #6 - 2 pc.	US\$ 138.66	+Info
475500	Pura Brushtip #8 - 2 pc.	US\$ 187.78	+Info
475501	Pura Brushtip #Stain - 2 pc.	US\$ 57.01	+Info
475502	PuraBrushtip #Opaque/Glaze - 2 pc.	US\$ 77.46	+Info

Kits / Set

Kits / Set

4 in 1 Set

4 in 1 Set

Code	Name	Price	+
475503	Pura Brush 4in1 Kit	US\$ 218.40	+Info

Master Kits

Master Kits

Code	Name	Price	+
475504	Pura Brush Master Kit #2	US\$ 168.58	+Info

475505	Pura Brush Master Kit #4	US\$ 172.18	+Info
475506	Pura Brush Master Kit #6	US\$ 185.70	+Info
475507	Pura Brush Master Kit #8	US\$ 205.54	+Info

PORCELAIN BRUSHES | TAKAGUTTI



Code	Name	Price	+
100783	TAKAGUTTI - Kolinsky Brushes Straight Round #2 RED k-Sable Pincel Porcelanas # BR752	US\$ 5.00	+Info
100786	TAKAGUTTI - Kolinsky Brushes Straight Round #8 RED k-Sable Pincel Porcelanas # BR758	US\$ 6.11	+Info

PORCELAIN BRUSHES | YETI

CARAMISTA CERAMIC BRUSH - YETI

YETI - CARAMISTA Brush for Ceramic

Description

* Storage of the water reservoir - the ceramic is kept evenly moist during layering

- * The natural spring force of the brush hair allows a razor-sharp separation, as well as targeted and accurate application
- * Noticeable extension of the service life of the "Caramista" brush hairs
- * The hair always returns to a fine point, even after prolonged use, so disheveled and splayed brush heads in the dry state are a thing of the past
- * Exclusive brush hair surpasses previously known properties
- * Pleasant handling thanks to ergonomically shaped brush handles
- * Handmade in Germany



Code	Name	Price	+
416917	YETI CARAMISTA-Brush-SET for ceramic (5 pc) No.1, no.O/B, no.6,no.8, no.10 - #31...	US\$ 97.33	+Info
416914	YETI CARAMISTA-Brush No.0 for ceramic (2 pc) - #312-3000	US\$ 18.02	+Info
416915	YETI CARAMISTA-Brush No.1 for ceramic (2 pc) - #312-3001	US\$ 24.57	+Info
416910	YETI CARAMISTA-Brush No.4 for ceramic - #312-2004	US\$ 15.14	+Info
416911	YETI CARAMISTA-Brush No.6 for ceramic - #312-2006	US\$ 24.57	+Info
416912	YETI CARAMISTA-Brush No.8 for ceramic - #312-2008	US\$ 24.57	+Info
416916	YETI CARAMISTA-Opaquer/Bonder Brush for ceramic (2 pc) - #312-4000	US\$ 28.65	+Info

Staining brush - YETI

Staining brush - YETI

Code	Name	Price	+
416923	CARAMISTA-PURE Set	US\$ 60.47	+Info

PORCELAIN BRUSHES | SHOFU

PORCELAIN BRUSHES - SHOFU

Shape-retentive, well-balanced, ergonomic

Natural-hair and synthetic porcelain brushes with convenient sculpting properties and an extremely long lifetime.

The use of high-quality porcelain brushes is basic to a precise and time-saving buildup of natural tooth shapes. SHOFU natural-hair brushes feature fine, shape-retentive tips and a continuous release of moisture.

The natural-hair brush A7 (A = animal hair) is therefore ideal for the application of porcelain in large quantities and the creation of delicate details.

The synthetic porcelain brush no. 5, made of very fine nylon bristles, has proved itself as an "all-round brush" for any application in porcelain layering. It maintains these properties even when using mixing liquids that contain glycerin.

Brushes types:

- * OP-O Flat natural-hair brush for an efficient application of opaquer pastes and opaque liners
- * S-O Fine liner for staining, designed to create very fine natural color accents
- * S-1 Ideal for the characterization of pressed or milled porcelain objects
- * P-3 All-round brush for daily use
- * P-5 Classic porcelain brush for all requirements
- * P-7 Big Brush design with a fine tip



Code	Name	Price	+
118815	SHOFU Porcelain Brush Set (6 Pcs) #6985	US\$ 69.60	+Info
118814	Porcelain Brush Op-0 #6986	US\$ 11.40	+Info
118813	Porcelain Brush S-0 #6988	US\$ 11.40	+Info
118812	Porcelain Brush S-1 #6989	US\$ 11.40	+Info
118811	Porcelain Brush P-3 #6991	US\$ 14.88	+Info
118810	Porcelain Brush P-5 #6992	US\$ 19.34	+Info
118809	Porcelain Brush P-7 #6993	US\$ 28.37	+Info

Brushes

Smile Line Replacement Tips

Smile Line Replacement Tips - SMILE LINE

Smile Line Replacement Tips - SMILE LINE

Code	Name	Price	+
123661	Kolinsky Buil Up #8R Spare Tip	US\$ 85.47	+Info

Catalog Index / Indice

POLISHING BRUSHES	1
Porcelain Brushes	14
Smile Line Replacement Tips	46