

Burs Diamond

FG Diamond Burs



FG ROBOT POINTS DIAMOND - SHOFU

FG ROBOT POINTS DIAMOND - SHOFU

The computerized production process guarantees an uniform distribution of a precise amount of natural diamond particles.

The diamond abrasives are available in the five different grit sizes superfine, fine, regular, coarse and supercoarse.

Six different shapes are available for minimally invasive therapy. Compared with conventional diamonds the working parts are small; the narrow and slender instrument neck allows an unimpeded view of the surface that is to be prepared.



Indication

In the dental practice

For cavity and crown preparations and for adjusting direct restorative materials

In the dental laboratory

E. G. Fine adjustments to full anatomically designed occlusal surfaces on porcelain crowns, inlays, onlays and approximal areas

Range of use

Suitable for the following materials: tooth enamel, porcelain, composites, filling materials such as GIC, cements

Benefits:

- * Exact, uniform distribution of the diamond particles
- * Shank and working end are produced from a one-piece blank
- * Long lasting
- * Precision running
- * Vibration-free action

Material properties

Working end

- * Uniformly coated with natural diamond particles
- * 5 grit sizes: supercoarse, coarse, standard, fine, superfine

Shapes and sizes

77 shapes for a wide range of preparations and procedures as well as 6 shapes for minimally invasive therapy

- * Supercoarse grit (X, black banded): 10
- * Coarse grit (C, green-banded shank): 6
- * Regular (plain shank): 77
- * Fine (F, red-banded shank): 11
- * Superfine (V, yellow-banded shank): 19

Shank

Stainless steel shank

Working end and shank as one piece and exactly centered

Shank shapes

FG = Turbine

Recommended speed range

120,000 - 450,000 min-1

The recommended speed range varies according to the type of shank and the shape and size of the working end.

The maximum permitted speed should not be exceeded!

Note

Robot Points FG can be cleaned in an ultrasonic bath and sterilized in an autoclave.

TAPERED CYLINDER (FLAT END)

Code	Name	Price	+
107-824	SHOFU FG Robot Points Diamond-1pk #0824-1 Bur	US\$ 7.49	+Info
107-820F	SHOFU FG Robot Points Diamond-1pk #0820F-1 Bur	US\$ 7.49	+Info
107-0830	SHOFU FG Robot Points Diamond-1pk #0830-1 Bur	US\$ 7.49	+Info
107-0822	SHOFU FG Robot Points Diamond-1pk #0822-1 Bur	US\$ 5.99	+Info

TAPERED CYLINDER (NEEDLE)

Code	Name	Price	+
107-850F	SHOFU FG Robot Points Diamond-1pk #0850F-1 Bur	US\$ 7.49	+Info
107-850	SHOFU FG Robot Points Diamond-1pk #0850-1 Bur	US\$ 7.49	+Info
107-848	SHOFU FG Robot Points Diamond-1pk #0848-1 Bur	US\$ 7.49	+Info

INVERTED CONE

Code	Name	Price	+
107-0866	SHOFU FG Robot Points Diamond-1pk #0866-1 Bur	US\$ 7.49	+Info
107-0864	SHOFU FG Robot Points Diamond-1pk #0864-1 Bur	US\$ 7.49	+Info
107-0860	SHOFU FG Robot Points Diamond-1pk #0860-1 Bur	US\$ 7.49	+Info

STRAIGHT CYLINDER

Code	Name	Price	+
107-0800	SHOFU FG Robot Diamond Points-1pk -1 Bur #0800	US\$ 7.49	+Info
107-803-1	SHOFU FG Robot Points Diamond-1pk #0803-1 Bur	US\$ 7.49	+Info

ROUND

Code	Name	Price	+
107-0874	SHOFU FG Robot Points Diamond-1pk #0874-1 Bur	US\$ 7.42	+Info
107-0875	SHOFU FG Robot Points Diamond-1pk #0875-1 Bur	US\$ 7.49	+Info

BUD

Code	Name	Price	+
107-883	SHOFU FG Robot Points Diamond-1pk #0883-1 Bur	US\$ 7.42	+Info
470534	SHOFU FG Robot Points Diamond-6pk #0860-1 Bur	US\$ 44.89	+Info
470535	SHOFU FG Robot Points Diamond-6pk #0848-1 Bur	US\$ 44.89	+Info
470537	SHOFU FG Robot Points Diamond-6pk # 0850-1 Piedra Diamante	US\$ 44.89	+Info
470538	SHOFU FG Robot Points Diamond-6pk # 0850F-1 Bur	US\$ 44.89	+Info

FLAME

Code	Name	Price	+
107-840F	SHOFU FG Robot Points Diamond-1pk #0840F-1 Bur	US\$ 5.99	+Info

PEAR

Code	Name	Price	+
107-0878	SHOFU FG Robot Points Diamond-1pk #0878-1 Pear Bur	US\$ 7.49	+Info

Burs Diamond

HP Diamond Burs



HP Diamond Burs - BESQUAL

HP Diamond Burs - BESQUAL

Code	Name	Price	+
470746	5X LOT of BESQUAL HP Diamond Coated	US\$ 5.80	+Info
470747	5X LOT of BESQUAL HP Diamond Coated	US\$ 5.80	+Info

HP Diamond Burs - JSP

HP Diamond Burs - JSP

Code	Name	Price	+
461400	Jsp - Hp Diamond Burs 6 Pack Round	US\$ 6.00	+Info
461401	Jsp - Diamond Burs 6 Pack Tapered	US\$ 6.00	+Info
461402	Jsp - Diamond Burs 6 Pack Cylinder	US\$ 6.00	+Info

HP DIAMOND COATED BUR - BESQUAL

HP DIAMOND COATED BUR - BESQUAL

We are proud to present to you our range of custom made diamond drills products, using the latest manufacturing techniques, resulting in a high quality precision product with excellent cutting qualities and attractive pricings.



Code	Name	Price	+
102163	BESQUAL HP Diamond Coated Bur - DC1 - Small Round #311-001	US\$ 3.95	+Info
102164	BESQUAL HP Diamond Coated Bur - DC2 - Large Round #311-002	US\$ 3.95	+Info
102165	BESQUAL HP Diamond Coated Bur - DC3 - Pointed Taper #311-003	US\$ 3.95	+Info
102168	BESQUAL HP Diamond Coated Bur - DC4 - Small Taper #311-004	US\$ 3.95	+Info
102169	BESQUAL HP Diamond Coated Bur - DC5 - Large Taper #311-005	US\$ 3.95	+Info
102170	BESQUAL HP Diamond Coated Bur - DC6 - Flat Taper Imported diamond powder - #311-006	US\$ 3.95	+Info
102171	BESQUAL HP Diamond Coated Bur - DC7 - Round Taper #311-007	US\$ 3.95	+Info
102172	BESQUAL HP Diamond Coated Bur - DC8 - Round Cylinder #311-008	US\$ 3.95	+Info
102173	BESQUAL HP Diamond Coated Bur - DC9 - Small Inverted - #311-009	US\$ 3.95	+Info
102174	BESQUAL HP Diamond Coated Bur - DC10 - Large Inverted #311-010	US\$ 3.95	+Info
102175	BESQUAL HP Diamond Coated Bur - DC11 - Small Pointed #311-011	US\$ 3.95	+Info
102176	BESQUAL HP Diamond Coated Bur - DC12 - Large Pointed Imported diamond powder - #...	US\$ 3.95	+Info
102177	BESQUAL HP Diamond Coated Bur - DC13 - Large Cylinder Imported diamond powder - ...	US\$ 3.95	+Info
102178	BESQUAL HP Diamond Coated Bur - DC14 - Small Cylinder - #311-014	US\$ 3.95	+Info
102179	BESQUAL - HP Diamond Coated Bur - DC15 - Small Wheel -	US\$ 3.95	+Info
102180	BESQUAL HP Diamond Coated Bur - DC16 - Large Wheel - #311-016	US\$ 3.95	+Info

Catalog Index / Indice

FG Diamond Burs	1
HP Diamond Burs	3