

Casting Accessories

Alloy rings



ALLOY CASTING RINGS - KEYSTONE

ALLOY CASTING RINGS - KEYSTONE

Features:

- * Stainless Steel.
- * High-Temperature resistance.
- * Shapes Available: Round and Oval.
- * Diameters Available: 35mm, 40mm, 45mm, 50mm, 65mm, 80mm, and 90mm.

* Even after repeated heatings they remain precision fitting and stable in size

* They do not flake, corrode, or peel away, thus they have an unusually long life.



US Dental Depot Is An Authorized Keystone Dealer In US.

Code	Name	Price	+
1070050	KEYSTONE Metal Casting Rings - 1.25in x 1.25in - 35mm x35mm Stainless Steel - Hi...	US\$ 7.35	+Info
034-1070070	KEYSTONE Metal Casting Rings - 1.25in x 1.5in - 35mm x 40mm Stainless Steel - Hi...	US\$ 7.14	+Info
1070060	KEYSTONE Metal Casting Rings - 1.25in x 1.75in - 35mm x45mm - Stainless Steel - ...	US\$ 7.35	+Info
034-1070080	KEYSTONE Metal Casting Rings - 1.5in x 1.5in - 40mm x 40mm Stainless Steel - HI ...	US\$ 9.30	+Info
1070085	KEYSTONE Metal Casting Rings - 1.5in x 1.75in - 40mm x 45mm Stainless Steel - Hi...	US\$ 10.27	+Info
034-1070090	KEYSTONE Metal Casting Rings - 1.75in x 1.5in - 45mm x 40mm Stainless Steel - HI...	US\$ 10.42	+Info
034-1700100	KEYSTONE Metal Casting Rings - 1.75in x 1.75in - 45mm x 45mm Stainless Steel - H...	US\$ 9.60	+Info
1070040	KEYSTONE Metal Casting Rings - 1.75in x 2in - 45mm x 50mm Stainless Steel - Hi T...	US\$ 11.62	+Info
034-1700110	KEYSTONE Casting Ring 2in x 2in , Stainless Steel - Hi Temp Resistant - #1700110	US\$ 14.17	+Info
034-1700120	KEYSTONE Casting Ring 2.5in x 2in , Stainless Steel - Hi Temp Resistant - #170...	US\$ 15.30	+Info
034-1700130	KEYSTONE Casting Ring 3in x 2.5in , Stainless Steel - Hi Temp Resistant - #170...	US\$ 21.53	+Info
034-1700140	KEYSTONE Casting Ring 3.5in x 2.5in , Stainless Steel - Hi Temp Resistant - #1...	US\$ 31.20	+Info

Alloy rings - JSP

Alloy rings - JSP

Code	Name	Price	+
437749	JSP Casting Ring & Base, Clear PVC	US\$ 8.98	+Info


Casting Rings - BUFFALO

Code	Name	Price	+
424846	AccuCast No. 17 Ring Flask 1-3/4" x	US\$ 10.08	+Info
424871	AccuCast No. 12 Ring Flask 1-1/4" x	US\$ 9.52	+Info
424877	AccuCast No. 35 Ring Flask 3-1/2" x	US\$ 21.63	+Info
425424	AccuCast No. 45 Ring Flask 3-1/2" x	US\$ 41.47	+Info
425426	AccuCast No. 25 Ring Flask 2-1/2" x	US\$ 11.20	+Info

CASTING RINGS - WHIPMIX

CASTING RINGS -WHIPMIX

All Whip Mix Casting Rings are made of a special Pyronel alloy metal or Stainless Steel. Special length rings in the diameters listed are available upon request. When ordering special rings always list the diameter first and then the height.



STAINLESS STEEL

Code	Name	Price	+
120287	WHIP-MIX #4391 Casting Ring 2.5 x 2.25 (63.5x57 mm) 1 pk - #11754	US\$ 46.01	+Info
120279	WHIP-MIX #4214 Casting Ring 2 x 1.87 (51 x 48 mm) 1 pk - #29319	US\$ 42.22	+Info
120285	WHIP-MIX #4089 Casting Ring 1.50 x 2.18 (38x55 mm) 1 pk - #11770	US\$ 33.43	+Info
120284	WHIP-MIX #4090 Metallic Casting Ring 1.75 x 2.2 (45x55 mm) 1pk - #11789	US\$ 38.22	+Info
120286	WHIP MIX #4177 Casting Ring 1.75 (44.45mm) 1 pk - #11967	US\$ 29.91	+Info
120283	WHIP-MIX #4181 Casting Ring 1.50 x 1.56 (38x40 mm) - 1 pk #11738	US\$ 27.06	+Info
120288	WHIP-MIX #4088 Casting Ring 1.25 x 2.2 (32x55 mm) 1 pk - #11762	US\$ 26.46	+Info
120282	WHIP-MIX #4087 Casting Ring 1.25 x 1.87 (32x48 mm) - 1 pk #11711	US\$ 24.75	+Info
120281	WHIP-MIX #4086 Casting Ring 1.25 x 1.62 (32x41 mm) - 1 pk #11703	US\$ 21.79	+Info
120280	WHIP-MIX #4085 Casting Ring 1.25 x 1.37 (32x35 mm) - 1 pk #11681	US\$ 20.17	+Info

PLASTIC CASTING RING

Code	Name	Price	+
120289	WHIP MIX #4093 Plastic Casting Ring 90g 1 pk - #32999	US\$ 39.35	+Info
120290	WHIP MIX #4092 Plastic Casting Ring 60g 1 pk - #35734	US\$ 39.37	+Info
120291	WHIP MIX #4095 Plastic Casting Ring 270g 1 pk - #33014	US\$ 39.35	+Info
120292	WHIP MIX #4094 Plastic Casting Ring 180g 1 pk - #33006	US\$ 29.42	+Info

Casting Accessories

Casting Accessories



Casting Accessories - BUFFALO

Casting Accessories - BUFFALO

Code	Name	Price	+
------	------	-------	---

425269	Replacement Shield w	US\$ 23.10	+Info
425272	Buffalo No. 3 Solder	US\$ 10.48	+Info

Casting Accesories - JSP

Casting Accesories - JSP

Code	Name	Price	+
461496	Jsp - Melting Dish 200Gm	US\$ 6.00	+Info
461498	Jsp - Melting Dish Handle	US\$ 8.13	+Info
461499	Jsp - Melting Dish + Melting Dish H	US\$ 16.56	+Info
461500	Jsp - Melting Dish	US\$ 8.44	+Info
461501	Jsp - Mold Plate(Pair) 5 x 7in	US\$ 7.80	+Info
461502	Jsp - Kevlar Mold Cuttng Glove Smal	US\$ 9.75	+Info

Casting Accesories - SHARK

Casting Accesories - SHARK

Code	Name	Price	+
251-215109	Pencil Refractory Mater ea	US\$ 6.50	+Info

EASY RINGS - VANIMAN

EASY RINGS - VANIMAN

Tapered for easy removal of the investment. A wide range of sizes is available to meet any casting requirement.



Code	Name	Price	+
122664	VANIMAN Easy Rings - 1in/4in D X 15/8in H - Rings And Formers - Each #Vmc-f125	US\$ 6.08	+Info
122665	VANIMAN Easy Rings - 1in/2in D X 13/4in H - Rings And Formers - Each #Vmc-f15	US\$ 6.08	+Info
122666	VANIMAN Easy Rings - 13/4in D X 19/16in H - Rings And Formers - Each #Vmc-f175	US\$ 5.81	+Info
122667	VANIMAN Easy Rings - 2in D X 2in H - Rings And Formers - Each #Vmc-f20	US\$ 6.60	+Info
122668	VANIMAN Easy Rings - 2in/2in D X 2in/4in H - Rings And Formers - #Vmc-f25	US\$ 8.45	+Info
122669	VANIMAN Easy Rings - D X 2in/4in H - Rings And Formers - Each #Vmc-f30	US\$ 10.56	+Info

Casting Accessories

Casting tongs



Casting tongs - BEGO

Casting tongs - BEGO

Code	Name	Price	+
10-11599	Mould Tongs 64cm Each	US\$ 180.76	+Info
10-39754	Mould Tongs 55cm 1pk	US\$ 168.94	+Info
470582	Mould Tongs 64cm 1pk	US\$ 133.22	+Info

Casting tongs - KEYSTONE

Casting tongs - KEYSTONE

Code	Name	Price	+
034-1620037	Tongs- Casting 15 Scissor Type	US\$ 34.72	+Info

Casting Accessories

Crucible formers



CRUCIBLE FORMER - KEYSTONE

CRUCIBLE FORMERS - KEYSTONE

Features:

- * The Crucible Formers have a tapered ring design for easy removal of mold.
- * These are made of rubber for long lasting and reliable performance.
- * These eliminate the need for ring liner.

Available Diameter Sizes:

- 1-1/4"
- 1-1/2"
- 1-3/4"
- 2"
- 2-1/2"
- 3"
- 3-1/2"

* Sold individually.

US Dental Depot Is An Authorized Keystone Dealer In US.



Code	Name	Price	+
034-1700141	KEYSTONE Rubber Crucible Former - BLACK - 1.25inch in Diameter - 35mm #1700141 -	US\$ 9.38	+Info
034-1700142	KEYSTONE Rubber Crucible Former - BLACK 1.5in Diameter - 40mm #1700142 -	US\$ 12.08	+Info
034-1700143	KEYSTONE Rubber Crucible Former, BLACK, 1-3/4in Diameter. Single Former - #170...	US\$ 13.08	+Info
034-1700144	KEYSTONE Rubber Crucible Former, BLACK, 2in Diameter. Single Former - #1700144	US\$ 13.72	+Info
034-1700145	KEYSTONE Rubber Crucible Former, BLACK, 2-1/2in Diameter. Single Former - #170...	US\$ 14.17	+Info
034-1700146	KEYSTONE Rubber Crucible Former, BLACK, 3in Diameter. Single Former - #1700146	US\$ 10.95	+Info
034-1700147	KEYSTONE Rubber Crucible Former - BLACK - 3.5inch in Diameter - 90mm #1700147 -	US\$ 15.60	+Info

Crucible formers - BUFFALO

Crucible formers - BUFFALO

Code	Name	Price	+
424720	No. 35R Rubber Sprue	US\$ 8.48	+Info
424898	No. 17R Rubber Sprue	US\$ 7.36	+Info
425398	No. 25R Rubber Sprue	US\$ 8.48	+Info
425399	No. 12R Rubber Sprue	US\$ 7.36	+Info

Crucible formers - NOBILIUM

Crucible formers - NOBILIUM

Code	Name	Price	+
N-1007	Plastic Sprue Former [blue] ea	US\$ 4.76	+Info
N-62167	Plastic Sprue Former white ea	US\$ 8.90	+Info

CRUCIBLE FORMERS - WHIPMIX

CRUCIBLE FORMERS - WHIPMIX

Whip Mix CRUCIBLE FORMERS are available in several diameters and large and small oval shapes made of durable plastic.



Code	Name	Price	+
120293	WHIP-MIX #4081 Crucible Former 1-1/4in - #11908	US\$ 26.03	+Info
120294	WHIP-MIX #4179 Crucible Former 1-1/2in - #11940	US\$ 26.03	+Info
120295	WHIP-MIX #4177 Crucible Former 1-3/4in - #11967	US\$ 26.03	+Info
120296	WHIP-MIX #4390 Crucible Former 2-1/2in - #11975	US\$ 26.03	+Info
120297	WHIP-MIX #4215 Crucible Former Conical 2in - #35149	US\$ 26.03	+Info

FUNNEL FORMERS - BEGO

FUNNEL FORMERS - BEGO

For combination crucible

Funnel former with reservoir for combination crucible



Code	Name	Price	+
10-52075	BEGO Funnel Former with Reservoir 10pk #52075	US\$ 94.05	+Info
10-52060	BEGO Funnel Former Standard 10pk #52060	US\$ 78.04	+Info
419867	BEGO UNIVERSAL FUNNEL FORMERS - 100 PIECES - #52068	US\$ 52.29	+Info
419866	BEGO FUNNEL FORMERS FOR NAUTILUS - 10 PIECES - #52066	US\$ 90.58	+Info

SPRUE FORMERS - VANIMAN

[SPRUE FORMER - VANIMAN](#)

A cylindrical silicon base which holds the sprues upright in the casting ring during investing. It forms a funnel in the investment to guide the molten material into the mold during casting.



Code	Name	Price	+
111391	VANIMAN Spectrum Sprue Former 1.25in - Purple - Each (#VMC-SB125)	US\$ 8.65	+Info
111392	VANIMAN Spectrum Sprue Former 1.50in - Dark BLUE - Each (#VMC-SB150)	US\$ 6.09	+Info
111390	VANIMAN Spectrum Sprue Former 1.75in - Medium BLUE - Each (#VMC-SB175)	US\$ 6.60	+Info
111393	VANIMAN - Spectrum Sprue Former 2 in - LIGHT BLUE - Each # VMC-SB200	US\$ 7.61	+Info

111394	VANIMAN Spectrum Sprue Former 2.5in - YELLOW - Each (#VMC-SB250)	US\$ 7.39	+Info
111395	VANIMAN Spectrum Sprue Former 3in - ORANGE - Each (#VMC-SB300)	US\$ 7.66	+Info

Casting Accessories

Crucible Ringless formers



Crucible Ringless formers - BEGO

Crucible Ringless formers - BEGO

Code	Name	Price	+
419875	Rubber Base w/ Funnel 4pk	US\$ 51.23	+Info

CRUCIBLE RINGLESS FORMERS - ECCO

CRUCIBLE RINGLESS FORMERS - ECCO



CASTING RINGLESS

Code	Name	Price	+
045-32999	ECCO - Casting Ringless and Crucible Former (base & ring) (#4093) - 1.5in Di...	US\$ 14.40	+Info
045-33006	ECCO - Castings Ringless and Crucible Formers (base & ring) #4094 2inch Diamet...	US\$ 14.40	+Info
101093	ECCO - Castings Ringless and Crucible Formers (base & ring) - #409 2.5in Diame...	US\$ 14.40	+Info
114469	ECCO - Castings Ringless and Crucible Formers (base & ring) Set of 3 pcs with ...	US\$ 54.68	+Info

PRESS RINGLESS

Code	Name	Price	+
100467	ECCO - Press Ringless And Crucible Formers (base & ring) - 90gr - Ringless 1.5...	US\$ 21.60	+Info
100468	ECCO - Press Ringless And Crucible Formers (base & ring) - 180gr - Ringless 2i...	US\$ 20.21	+Info

Crucible Ringless formers - JSP

Crucible Ringless formers - JSP

Code	Name	Price	+
437750	JSP Casting Ring & Base, Clear PVC	US\$ 8.00	+Info
461509	Jsp - Clear Castng Ring & Base 1 1/	US\$ 8.00	+Info
461510	Jsp - Clear Castng Ring&Base 2in	US\$ 8.00	+Info
461511	Silicone Model Formers - Set of 5	US\$ 11.69	+Info

MOULD FORMERS - BEGO

MOULD FORMERS - BEGO

- * Eliminates fixing and grinding of the investment models when the BEGO combination duplicating flask is used.
- * Both mould formers can also be used with all other duplicating systems.

US DENTAL DEPOT Inc is an Authorized BEGO Dealer in the US.



Code	Name	Price	+
10-52400	BEGO Mould Formers Large BLUE 4pk #52400	US\$ 88.09	+Info
10-52390	BEGO Mould Formers Small RED 4pk #52390	US\$ 62.98	+Info

MOULD RINGS - BEGO

MOULD RINGS - BEGO

- * For crown and bridge work
- * Metal mould ring size 1 for 60 g of investment material.
- * Long service life thanks to special steel design.
- * Suitable for all BEGO crown and bridge investment materials.

US DENTAL DEPOT Inc is an Authorized BEGO Dealer in the US.



Code	Name	Price	+
10-52424	BEGO Mould Rings Size 9 4pk #52424	US\$ 158.12	+Info
10-52423	BEGO Mould Rings Size 6 4pk #52423	US\$ 136.89	+Info
10-52422	BEGO Mould Rings size #3 4pk #52422	US\$ 158.12	+Info
10-52419	BEGO - Mould Rings Size 1 4pk #52419	US\$ 62.56	+Info

PERFECT RINGS - ADS

PERFECT RINGS - ADS

Introducing all new Perfect Ring specifically engineered for use with lithium disilicate restorations. Perfect Ring is made from the same clear material our traditional ringless system for casting alloys. The height and volume dimensions have been specifically designed to fit all current pressing furnaces.

- * Eliminates investment cracking
- * Increases green strength of investment
- * No leaking investment
- * No excess investment flash
- * More convenient and saves time
- * Works great and gives a perfect fitting ring



Code	Name	Price	+
118269	ADS Perfect Ring Hot PINK for Lithium Disilicate 100 gms 13mm diameter plunger ...	US\$ 41.62	+Info
118270	ADS Perfect Ring Hot PINK for Lithium Disilicate 200 gms 13mm diameter plunger ...	US\$ 41.62	+Info

RINGLESS RING AND FORMERS - BEGO

RINGLESS RING AND FORMERS - BEGO

Universally applicable for many casting systems.

Description:

- * Time savings in relation to mould systems with foil sleeve, iron ring, etc..
- * Compatible with Rapid Wax System.
- * Minimal wear, thus lower costs than with comparable systems.
- * Universally applicable for many casting systems.
- * Easy separation of mould and mould ring.
- * For all BEGO crown and bridge investment materials.



US DENTAL DEPOT Inc is an Authorized BEGO Dealer in the US.

Code	Name	Price	+
10-52665	BEGO Ringless Ring & Former 1.25in 1pk (#52665) Universally applicable for ma...	US\$ 59.58	+Info
10-52665A	BEGO - Ringless Ring & Former 1.5in 1pk (#52665A) Universally applicable for m...	US\$ 94.58	+Info
10-52666	BEGO Ringless Ring & Former 2in 1pk (#52666) Universally applicable for many c...	US\$ 59.58	+Info
10-52667	BEGO Ringless Ring & Former 2.5in 1pk (#52667) Universally applicable for many...	US\$ 59.58	+Info

Casting Accessories

Crucibles all types



Crucibles all types - BEGO

Crucibles all types - BEGO

Code	Name	Price	+
10-52477	Crucible Spray/Lubricant 100ml	US\$ 55.50	+Info

Crucibles all types - BUFFALO

Crucibles all types - BUFFALO

Code	Name	Price	+
------	------	-------	---

009-28012	Crucible Jumbo Clay Round Hole	US\$ 21.56	+Info
425396	Type 'A' Quartz cruc	US\$ 25.48	+Info
425397	Type 'A' clay crucib	US\$ 21.56	+Info

Crucibles all types - ECCO

Crucibles all types - ECCO

Code	Name	Price	+
421028	Hooded Crucible 7oz	US\$ 27.83	+Info

Crucibles all types - EDG

Crucibles all types - EDG

Code	Name	Price	+
102295	Ceramic Crucible for Powercast	US\$ 265.00	+Info

Crucibles all types - JSP

Crucibles all types - JSP

Code	Name	Price	+
437640	JSP HOODED QUARTZ CRUCIBLE, 2/Pk. A	US\$ 8.65	+Info

Crucibles all types - KERR

Crucibles all types - KERR

Code	Name	Price	+
101453	Zircon Alumina Crucible Yellow	US\$ 30.87	+Info

Crucibles all types - NOBILIUM

Crucibles all types - NOBILIUM

Code	Name	Price	+
112257	Fused Silica Crucible 3 pk	US\$ 37.35	+Info

CRUCIBLES FOR FORNAX - BEGO

CRUCIBLES FOR FORNAX - BEGO

The new Bego ceramic crucibles for Fornax set the highest standards. An innovative manufacturing procedure for crucibles that are resistant high temperature, developed in cooperation with scientific institutes produces an extremely homogenous material structure contributing to consistently reproducible accuracy of castings.

US DENTAL DEPOT is an Authorized BEGO dealer in US.



Code	Name	Price	+
10-52455	BEGO Crucible Ceramic Inserts for Fornax 6pk #52455 #52455	US\$ 131.67	+Info
10-52843	BEGO NEW Crucible for Fornax 6pk #52483	US\$ 207.58	+Info
10-52482	BEGO - Melting Ceramic Crucible for Fornax T & G models 6pk #52482 #52483	US\$ 201.23	+Info
10-52454	BEGO Crucible Graphite Inserts for Fornax 6pk #52454 #52454	US\$ 170.88	+Info

CRUCIBLES FOR FUNDOR-T - BEGO


CRUCIBLES FOR FUNDOR-T - BEGO

Ceramic Crucibles for cobalt-chrome partial denture alloys and non-precious alloys.

The melting crucibles are made of special ceramic material and have a long service life through a high heat change resistance.

For Melting crucible, deep trough, shallow trough, for precious metal, for cobalt-chrome partial denture alloys and non-precious alloys that can be used with Fundor T.

US DENTAL DEPOT is an Authorized BEGO dealer in US.



Code	Name	Price	+
10-52426	BEGO Crucible for Co-Cr + Non Precious Alloys for Fundor T 6pk #52426 #52426	US\$ 198.91	+Info
10-52425	BEGO Crucible for Precious Alloys for Fundor T 6pk #52425 #52425	US\$ 198.91	+Info

CRUCIBLES TYPE - KEYSTONE

CRUCIBLES TYPES - KEYSTONE

Features:

* Most popular Crucibles Types that fit standard and different types of casting machines.

* Types :

Centrifugal

Standard Hooded

Standard Slotted

Two Stage Slotted

* Generally Sold Individually.

US Dental Depot Is An Authorized Keystone Dealer In US.



CENTRIFUGAL

KEYSTONE - CENTRIFUGAL

Features:

* This Centrifugal Type Crucible fits standard Dentsply type casting machines. Sold individually.

Code	Name	Price	+
034-5920250	KEYSTONE Crucible Centrifical Type - #5920250	US\$ 23.55	+Info

STANDARD HOODED

KEYSTONE - STANDARD HOODED

Features:

- * This high quality Standard Hooded Type Crucible offers optimum performance in induction furnaces.
- * This crucible fits standard type Kerr, Lucas and Emesco casting machines.

Code	Name	Price	+
034-5920235	KEYSTONE Crucible Hooded Standard - #5920235	US\$ 22.95	+Info

STANDARD SLOTTED

KEYSTONE - STANDARD SLOTTED

Features:

- * These high quality Standard Slotted Type Crucibles offer optimum performance.
- * This crucible fits standard type Kerr and Emesco casting machines.

Code	Name	Price	+
034-5920200	KEYSTONE Crucible Slotted Standard - #5920200	US\$ 22.95	+Info
034-5920201	KEYSTONE #32 Standard Slotted Type Crucible - for European Machines Same As #30 ...	US\$ 24.11	+Info

TWO STAGE SLOTTED

KEYSTONE - TWO STAGE SLOTTED CRUCIBLE

Features:

- * This Two State Slotted Crucible fits standard type Kerr and Emesco casting machines.

Code	Name	Price	+
034-5920230	KEYSTONE Crucible 2 Stage Slotted - #5920230	US\$ 21.75	+Info

HI-HEAT QUARTZ CRUCIBLES - B ESQUAL

HI-HEAT QUARTZ CRUCIBLES - BESQUAL

The large, flat melting area allows more alloy of various shapes to be heated evenly, an important factor in the proper melting of nickel based alloys. The design of this crucible reflects heat for more uniform melting. For broken arm casting machines.



Code	Name	Price	+
------	------	-------	---

100126	BESQUAL Hi Heat Quartz Standard Crucible Slotted (Broken Arm Casting) . - Kerr /...	US\$ 9.41	+Info
100044	BESQUAL Hi Heat Quartz Standard Hooded Crucible Round Orifice (QHN-2/719-000)- F...	US\$ 10.12	+Info
116093	BESQUAL Hi Heat Quartz Mini Crucible Hooded - Round Orifice - For Broken Arm Cas...	US\$ 9.34	+Info

INDUCTION CRUCIBLE - KEYSTONE

INDUCTION CRUCIBLE - KEYSTONE



Features:

* This high quality Induction Crucible offers optimum performance in induction furnaces.

US Dental Depot Is An Authorized Keystone Dealer In US.

Code	Name	Price	+
113656	KEYSTONE Induction Crucible II each Fits standard type Howmedica and electronic ...	US\$ 24.60	+Info
113655	KEYSTONE Induction Crucible Hooded #49A - Induction Crucible Hooded #49A #5920245	US\$ 25.28	+Info
113654	KEYSTONE Induction Crucible - Fits standard type Tycomatic and Unitek casting ma...	US\$ 22.35	+Info

NAUTILUS CRUCIBLES - BEGO

NAUTILUS CRUCIBLES - BEGO

BEGO has pushed the envelope once again with the newly designed Nautilus Crucibles. Automatic casting for consistent results now has new options such as adjustable melt hold times and optional longer pressure settings after casting - these features give the technician greater control in casting larger implant bars and full arch C&B frames. The new design has a large color touch screen for easy operation, energy saving "Eco-Mode" and integrated cooling system making the new Nautilus Crucibles the intelligent environmentally friendly choice. The induction heating guarantees a shorter melting cycle and minimizes oxidation which facilitates finishing time. Let BEGO help you reach your C&B goal of consistent castings by learning more about the proven performance of vacuum/pressure, split crucible technology of the Nautilus Crucible induction casting machine.



Code	Name	Price	+
10-52473	BEGO Nautilus Glass Inserts for crucibles 4pk #52473 #52473	US\$ 317.30	+Info

10-52488	BEGO Crucible Ceramic for all Nautilus casting models 4pk #52488 #-52488	US\$ 184.90	+Info
10-52468	BEGO Nautilus Graphite Inserts for crucibles 6pk #52468 #52468	US\$ 130.75	+Info

Casting Accessories

Debubblizers & Wetting Agents



CERPRESS SL TENSION RELIEF SPRAY - ADS

CERPRESS SL TENSION RELIEF SPRAY - ADS

- * Cerpress SL Pressable Ceramic System is used by thousands of dental laboratories throughout the country.
- * It has been on the market over 6 years and been clinically tested by several organizations.



Code	Name	Price	+
117712	ADS Cerpress SI Tension Relief Spray 4oz. - #201045	US\$ 18.17	+Info

CLEARLY DEBUBBLIZER - ADS

CLEARLY DEBUBBLIZER - ADS

Features:

- * Advanced formula Debubblizer for the new resin-toughened stone as recommended by their manufacturers
- * Completely water soluble, crystal clear and alcohol-free
- * Produces perfect castings
- * Investments flow smoothly across wax patterns
- * Essential on the new alginate/silicone hybrid and chromatic impression materials
- * Clearly Debubblizer can be applied by painting, spraying or dipping
- * Thinnest coating gives the best results



Code	Name	Price	+
117753	ADS Debubblizer CLEARLY gal. - #D926-8	US\$ 41.26	+Info
117752	ADS Debubblizer CLEARLY qt - #D926-6	US\$ 20.68	+Info
117751	ADS Debubblizer CLEARLY pt 16oz. - #D926-5	US\$ 14.77	+Info
117750	ADS Debubblizer CLEARLY / SPRAY 8oz. - #D926-4	US\$ 11.35	+Info
117749	ADS Debubblizer CLEARLY - SPRAY 4oz. - #D926-3	US\$ 8.76	+Info

DEBUBBLIZER AND WAX PATTERN CLEAR - ADS

[DEBUBBLIZER AND WAX PATTERN CLEAR - ADS](#)



Features:

- * Formulated for investment to flow evenly over the wax pattern
- * Eliminates bubbles and reduces the surface tension of waxes
- * Produces castings with a satin smooth finish
- * Spray on all impression materials before pouring in stone
- * Cleans epoxy dies of oils and particles

Code	Name	Price	+
117758	ADS Debubblizer & Wax PATTERN CLEANER gal. - #D916-8	US\$ 52.23	+Info
117757	ADS Debubblizer & Wax PATTERN CLEANER qt - #0916-6	US\$ 20.66	+Info
117756	ADS Debubblizer & Wax PATTERN CLEANER pt 16oz. - #D916-5	US\$ 14.18	+Info
117755	ADS Debubblizer & Wax PATTERN CLEANER 8oz. - spray - #D916-4	US\$ 10.56	+Info
117754	ADS Debubblizer & Wax PATTERN CLEANER 4oz. - spray - #0916-3	US\$ 8.78	+Info

Debubblizers & Wetting Agents - ALMORE

Debubblizers & Wetting Agents - ALMORE

Code	Name	Price	+
423080	Almore Debubblizer Debubblizer/Surf	US\$ 35.70	+Info
423081	Almore Debubblizer Debubblizer/Surf	US\$ 17.00	+Info

Debubblizers & Wetting Agents - BEGO

Debubblizers & Wetting Agents - BEGO

Code	Name	Price	+
105064	Aurofilm Wetting Agent 3.5oz	US\$ 25.62	+Info

Debubblizers & Wetting Agents - JSP

Debubblizers & Wetting Agents - JSP

Code	Name	Price	+
437629	JSP Debubblizer Surfactant, Blue -	US\$ 10.00	+Info
437654	JSP Debubblizer Surfactant - 8 oz.	US\$ 8.95	+Info
437663	JSP Debubblizer Surfactant, Blue -	US\$ 29.54	+Info
437734	JSP Debubblizer Surfactant - 16 oz.	US\$ 13.66	+Info
437735	JSP Debubblizer Surfactant - 128 oz	US\$ 29.54	+Info
437736	JSP Debubblizer Surfactant - 32 oz.	US\$ 18.79	+Info
437737	JSP Debubblizer Surfactant, Blue -	US\$ 13.66	+Info

Debubblizers & Wetting Agents - KEYSTONE

Debubblizers & Wetting Agents - KEYSTONE

Code	Name	Price	+
034-6121900	Blue Debubblizer 1 Gal	US\$ 31.88	+Info
445889	Mizzy Wax and Impression Debubblize	US\$ 12.84	+Info

Debubblizers & Wetting Agents - ZHERMACK

Debubblizers & Wetting Agents - ZHERMACK

Code	Name	Price	+
115863	Tensilab 120ml	US\$ 17.65	+Info

PICOSILK SPRAY BOTTLE - RENFER

T

PICOSILK BOTTLE - RENFERT

Picosilk is a wetting agent for wax models within crown, bridge and model casting technique. For tension-free wax patterns

Advantages

- * Spray on finely using the atomizer before investing.
- * Wax does not take on water-repellent characteristics after use.
- * Improves the flow of investment material.



Code	Name	Price	+
023-1550-0100	RENFERT Picosilk Refill Bottle-500ml	US\$ 45.42	+Info
105083	RENFERT Picosilk Agent w/spray 75ml- #1550-0000 #15500000	US\$ 28.26	+Info

SIMPLY SURFACTANT - ADS

SIMPLY SURFACTANT - ADS

The superior wetting agent. A modern surface tension eliminator designed for the new resin toughened stones and silicone impression materials. Formulated for investment to flow smoothly over wax patterns and hydrophobic impressions. Easy to use, simply paint, spray or dip wax patterns and impressions. Drain or blow off excess liquid. Produces defect-free pours with a fine, smooth finish and improved surface detail.



Code	Name	Price	+
------	------	-------	---

118377	ADS Simply Surfactant LIQUID Wetting AGENT 8oz. - #D 936-5	US\$ 14.77	+Info
118378	ADS Simply Surfactant LIQUID Wetting AGENT pt 16oz. - #D 936-6	US\$ 20.66	+Info
118379	ADS Simply Surfactant LIQUID Wetting AGENT qt - #S 650-2	US\$ 33.50	+Info

Casting Accessories

Liners 4 Casting Rings



Liners 4 Casting Rings - GC AMERICA

Liners 4 Casting Rings - GC AMERICA

Code	Name	Price	+
423004	GC America New Casting Liner - Line	US\$ 86.75	+Info

LINING TAPE - ADS

LINING TAPE - ADS

Type 380 ceramic fiber tape is asbestos-free and made from high purity refractory fibers (not cellulose) with a melting point of 3400F. Lining Tape is used as a direct replacement for asbestos ring liners. Available in convenient sizes to fit all casting rings. Avoids asbestos hazard. Use with all investments including high temperature types.



Code	Name	Price	+
118134	ADS Lining Tape Ceramic 1 3/4in x 1/16in x 50in roll - #T380-6 casting ring line...	US\$ 34.97	+Info
118135	ADS Lining Tape Ceramic 2in x 1/32in x 50in roll - #T380-7 casting ring liners]	US\$ 27.09	+Info
118132	ADS Lining Tape Ceramic 1 1/2in x 1/16in x 50in roll - #T380-4 casting ring line...	US\$ 25.74	+Info
118133	ADS Lining Tape Ceramic 1 3/4in x 1/32in x 50in roll - #T380-5 casting ring line...	US\$ 25.58	+Info
118130	ADS Lining Tape Ceramic 1 1/4in x 1/16in x 50in roll - #T380-2 casting ring line...	US\$ 23.93	+Info
118131	ADS Lining Tape Ceramic 1 1/2in x 1/32in x 50in roll - #T380-3 casting ring line...	US\$ 22.42	+Info
118129	ADS Lining Tape Ceramic 1 1/4in x 1/32in x 50in roll - #T380-1 casting ring line...	US\$ 18.07	+Info

RING LINERS - KEYSTONE

RING LINERS - KEYSTONE

Our Casting Ring Liner is made of refractory fibers that have a melting point in excess of 3,200F. The Liner is asbestos free and sold in 50-foot rolls. Various available sizes.

Features:

- * Asbestos free
- * 50 foot rolls

Available Sizes: 1-1/4" H x 1/32" Thick, 1-1/4" H x 1/16" Thick, 1-1/2" H x 1/32" Thick, 1-1/2" H x 1/16" Thick, 1-3/4" H x 1/32" Thick, 1-3/4" H x 1/16" Thick, 2" H x 1/32" Thick, 2" H x 1/16" Thick, 3" H x 1/16" Thick



Code	Name	Price	+
101147	KEYSTONE Casting Ring Liner 1 1/4in Ht. x 1/32in Thick. Asbestos-Free, Single...	US\$ 17.55	+Info
101148	KEYSTONE Casting Ring Liner 1 1/4in Ht. x 1/16in Thick. Asbestos-Free, Single...	US\$ 20.25	+Info
101149	KEYSTONE Casting Ring Liner 1 1/2in Ht. x 1/32in Thick. Asbestos-Free, Single...	US\$ 19.80	+Info
101150	KEYSTONE Casting Ring Liner 1 1/2in Ht. x 1/16in Thick. Asbestos-Free, Single...	US\$ 27.97	+Info
101151	KEYSTONE Casting Ring Liner 1 3/4in Ht. x 1/32in Thick. Asbestos-Free, Single...	US\$ 28.65	+Info
101152	KEYSTONE Casting Ring Liner 1 3/4in Ht. x 1/16in Thick. Asbestos-Free, Single...	US\$ 37.12	+Info
101153	KEYSTONE Casting Ring Liner 2in Ht. x 1/32in Thick. Asbestos-Free, Single 50 ...	US\$ 33.60	+Info
101154	KEYSTONE Casting Ring Liner 2in Ht. x 1/16in Thick. Asbestos-Free, Single 50 ...	US\$ 39.68	+Info
101155	KEYSTONE Casting Ring Liner 3in Ht. x 1/16in Thick. Asbestos-Free, Single 50 ...	US\$ 52.27	+Info

TAPE EEZ - ADS

TAPE EEZ - ADS

Super concentrated liquid for wetting ADS casting ring lining tapes. One pint makes 5 quarts. Designed to create the workability and handling texture of asbestos. Expansion properties of the investment will not be altered. Burns out clean in the furnace without charring.



Code	Name	Price	+
118411	ADS Tape Eez qt - #N 928-6 (wetting agent for ring lanners)	US\$ 23.48	+Info
118410	ADS Tape Eez 16oz. pt - #T 770.6 (wetting agent for ring lanners)	US\$ 18.88	+Info

Catalog Index / Indice

Alloy rings	1
Casting Accesories	2
Casting tongs	3
Crucible formers	4
Crucible Ringless formers	6
Crucibles all types	8
Debubblizers & Wetting Agents	13
Liners 4 Casting Rings	16