

## Casting Accessories

### Alloy rings



### ALLOY CASTING RINGS - KEYSTONE

#### ALLOY CASTING RINGS - KEYSTONE

**Features:**

- \* Stainless Steel.
- \* High-Temperature resistance.
- \* Shapes Available: Round and Oval.
- \* Diameters Available: 35mm, 40mm, 45mm, 50mm, 65mm, 80mm, and 90mm.

- \* Even after repeated heatings they remain precision fitting and stable in size
- \* They do not flake, corrode, or peel away, thus they have an unusually long life.



**US Dental Depot Is An Authorized Keystone Dealer In US.**

| Code        | Name  | Price      | +                     |
|-------------|---|------------|-----------------------|
| 1070050     | KEYSTONE Metal Casting Rings - 1.25in x 1.25in - 35mm x35mm Stainless Steel - Hi... | US\$ 9.31  | <a href="#">+Info</a> |
| 034-1070070 | KEYSTONE Metal Casting Rings - 1.25in x 1.5in - 35mm x 40mm Stainless Steel - Hi... | US\$ 11.38 | <a href="#">+Info</a> |
| 1070060     | KEYSTONE Metal Casting Rings - 1.25in x 1.75in - 35mm x45mm - Stainless Steel - ... | US\$ 8.94  | <a href="#">+Info</a> |
| 034-1070080 | KEYSTONE Metal Casting Rings - 1.5in x 1.5in - 40mm x 40mm Stainless Steel - HI ... | US\$ 11.66 | <a href="#">+Info</a> |
| 1070085     | KEYSTONE Metal Casting Rings - 1.5in x 1.75in - 40mm x 45mm Stainless Steel - Hi... | US\$ 11.99 | <a href="#">+Info</a> |
| 034-1070090 | KEYSTONE Metal Casting Rings - 1.75in x 1.5in - 45mm x 40mm Stainless Steel - HI... | US\$ 12.78 | <a href="#">+Info</a> |
| 034-1700100 | KEYSTONE Metal Casting Rings - 1.75in x 1.75in - 45mm x 45mm Stainless Steel - H... | US\$ 10.86 | <a href="#">+Info</a> |
| 1070040     | KEYSTONE Metal Casting Rings - 1.75in x 2in - 45mm x 50mm Stainless Steel - Hi T... | US\$ 13.56 | <a href="#">+Info</a> |
| 034-1700110 | KEYSTONE Casting Ring 2in x 2in , Stainless Steel - Hi Temp Resistant - #1700110    | US\$ 15.60 | <a href="#">+Info</a> |
| 034-1700120 | KEYSTONE Casting Ring 2.5in x 2in , Stainless Steel - Hi Temp Resistant - #170...   | US\$ 16.73 | <a href="#">+Info</a> |
| 034-1700130 | KEYSTONE Casting Ring 3in x 2.5in , Stainless Steel - Hi Temp Resistant - #170...   | US\$ 21.53 | <a href="#">+Info</a> |
| 034-1700140 | KEYSTONE Casting Ring 3.5in x 2.5in , Stainless Steel - Hi Temp Resistant - #1...   | US\$ 31.20 | <a href="#">+Info</a> |

### Casting Rings - BUFFALO

| Code   | Name                                | Price      | +                     |
|--------|-------------------------------------|------------|-----------------------|
| 424846 | AccuCast No. 17 Ring Flask 1-3/4" x | US\$ 11.34 | <a href="#">+Info</a> |
| 424871 | AccuCast No. 12 Ring Flask 1-1/4" x | US\$ 11.78 | <a href="#">+Info</a> |
| 424877 | AccuCast No. 35 Ring Flask 3-1/2" x | US\$ 21.63 | <a href="#">+Info</a> |
| 425424 | AccuCast No. 45 Ring Flask 3-1/2" x | US\$ 41.47 | <a href="#">+Info</a> |
| 425426 | AccuCast No. 25 Ring Flask 2-1/2" x | US\$ 14.40 | <a href="#">+Info</a> |

### CASTING RINGS - WHIP MIX

#### CASTING RINGS -WHIP MIX

All Whip Mix Casting Rings are made of a special Pyronel alloy metal or Stainless Steel. Special length rings in the diameters listed are available upon request. When ordering special rings always list the diameter first and then the height.



#### STAINLESS STEEL

| Code   | Name  | Price      | +     |
|--------|---|------------|-------|
| 120287 | WHIP MIX #4391 Casting Ring 2.5 x 2.25 ( 63.5x57 mm) 1 pk - #11754      | US\$ 53.92 | +Info |
| 120279 | WHIP MIX #4214 Casting Ring 2 x 1.87 (51 x 48 mm) 1 pk - #29319         | US\$ 49.50 | +Info |
| 120285 | WHIP MIX #4089 Casting Ring 1.50 x 2.18 (38x55 mm) 1 pk - #11770        | US\$ 39.19 | +Info |
| 120284 | WHIP MIX #4090 Metallic Casting Ring 1.75 x 2.2 (45x55 mm) 1pk - #11789 | US\$ 46.97 | +Info |
| 120283 | WHIP MIX #4181 Casting Ring 1.50 x 1.56 (38x40 mm) - 1 pk #11738        | US\$ 31.72 | +Info |
| 120288 | WHIP MIX #4088 Casting Ring 1.25 x 2.2 (32x55 mm) 1 pk - #11762         | US\$ 34.64 | +Info |
| 120282 | WHIP MIX #4087 Casting Ring 1.25 x 1.87 (32x48 mm) - 1 pk #11711        | US\$ 29.00 | +Info |
| 120281 | WHIP MIX #4086 Casting Ring 1.25 x 1.62 (32x41 mm) - 1 pk #11703        | US\$ 25.55 | +Info |
| 120280 | WHIP MIX #4085 Casting Ring 1.25 x 1.37 (32x35 mm) - 1 pk #11681        | US\$ 25.26 | +Info |

### EASY RINGS - VANIMAN

#### EASY RINGS - VANIMAN

Tapered for easy removal of the investment. A wide range of sizes is available to meet any casting requirement.



| Code   | Name   | Price      | +     |
|--------|--|------------|-------|
| 122664 | VANIMAN Easy Rings - 1in/4in D X 15/8in H - Rings And Formers - Each #Vmc-f125 | US\$ 6.26  | +Info |
| 122665 | VANIMAN Easy Rings - 1in/2in D X 13/4in H - Rings And Formers - Each #Vmc-f15  | US\$ 6.26  | +Info |
| 122666 | VANIMAN Easy Rings - 13/4in D X 19/16in H - Rings And Formers - Each #Vmc-f175 | US\$ 6.16  | +Info |
| 122667 | VANIMAN Easy Rings - 2in D X 2in H - Rings And Formers - Each #Vmc-f20         | US\$ 7.00  | +Info |
| 122668 | VANIMAN Easy Rings - 2in/2in D X 2in/4in H - Rings And Formers - #Vmc-f25      | US\$ 8.96  | +Info |
| 122669 | VANIMAN Easy Rings - D X 2in/4in H - Rings And Formers - Each #Vmc-f30         | US\$ 11.21 | +Info |

#### Casting Accessories

### Casting Accessories



### Casting Accesories - BUFFALO

#### Casting Accesories - BUFFALO

| Code   | Name                 | Price      | +     |
|--------|----------------------|------------|-------|
| 425269 | Replacement Shield w | US\$ 23.10 | +Info |
| 425272 | Buffalo No. 3 Solder | US\$ 11.79 | +Info |

### Casting Accesories - SHARK

Casting Accesories - SHARK

| Code       | Name                       | Price     | +                     |
|------------|----------------------------|-----------|-----------------------|
| 251-215109 | Pencil Refractory Mater ea | US\$ 6.50 | <a href="#">+Info</a> |

### Casting Accessories

### Casting tongs



### Casting tongs - BEGO

Casting tongs - BEGO

| Code     | Name                 | Price       | +                     |
|----------|----------------------|-------------|-----------------------|
| 10-39754 | Mould Tongs 55cm 1pk | US\$ 168.94 | <a href="#">+Info</a> |
| 470582   | Mould Tongs 64cm 1pk | US\$ 133.22 | <a href="#">+Info</a> |

### Casting tongs - KEYSTONE

Casting tongs - KEYSTONE

| Code        | Name                           | Price      | +                     |
|-------------|--------------------------------|------------|-----------------------|
| 034-1620037 | Tongs- Casting 15 Scissor Type | US\$ 44.85 | <a href="#">+Info</a> |

### Casting Accessories

### Crucible formers



### CRUCIBLE FORMER - KEYSTONE

#### CRUCIBLE FORMERS - KEYSTONE

#### Features:

- \* The Crucible Formers have a tapered ring design for easy removal of mold.
- \* These are made of rubber for long lasting and reliable performance.
- \* These eliminate the need for ring liner.

#### Available Diameter Sizes:

- 1-1/4"
- 1-1/2"
- 1-3/4"
- 2"
- 2-1/2"
- 3"
- 3-1/2"

\* Sold individually.

**US Dental Depot Is An Authorized Keystone Dealer In US.**



| Code | Name | Price | + |
|------|------|-------|---|
|------|------|-------|---|

|             |   |            |                       |
|-------------|---|------------|-----------------------|
| 034-1700141 | KEYSTONE Rubber Crucible Former - BLACK - 1.25inch in Diameter - 35mm #1700141 -  | US\$ 11.81 | <a href="#">+Info</a> |
| 034-1700142 | KEYSTONE Rubber Crucible Former - BLACK 1.5in Diameter - 40mm #1700142 -          | US\$ 13.28 | <a href="#">+Info</a> |
| 034-1700143 | KEYSTONE Rubber Crucible Former, BLACK, 1-3/4in Diameter. Single Former - #170... | US\$ 14.69 | <a href="#">+Info</a> |
| 034-1700144 | KEYSTONE Rubber Crucible Former, BLACK, 2in Diameter. Single Former - #1700144    | US\$ 15.10 | <a href="#">+Info</a> |
| 034-1700145 | KEYSTONE Rubber Crucible Former, BLACK, 2-1/2in Diameter. Single Former - #170... | US\$ 15.59 | <a href="#">+Info</a> |
| 034-1700146 | KEYSTONE Rubber Crucible Former, BLACK, 3in Diameter. Single Former - #1700146    | US\$ 16.17 | <a href="#">+Info</a> |
| 034-1700147 | KEYSTONE Rubber Crucible Former - BLACK - 3.5inch in Diameter - 90mm #1700147 -   | US\$ 15.44 | <a href="#">+Info</a> |

### Crucible formers - BEGO

Crucible formers - BEGO

| Code   | Name                  | Price      | +                     |
|--------|-----------------------|------------|-----------------------|
| 419893 | Silicone mould former | US\$ 61.90 | <a href="#">+Info</a> |

### Crucible formers - BUFFALO

Crucible formers - BUFFALO

| Code   | Name                 | Price      | +                     |
|--------|----------------------|------------|-----------------------|
| 424720 | No. 35R Rubber Sprue | US\$ 10.49 | <a href="#">+Info</a> |
| 424898 | No. 17R Rubber Sprue | US\$ 9.11  | <a href="#">+Info</a> |
| 425398 | No. 25R Rubber Sprue | US\$ 10.49 | <a href="#">+Info</a> |
| 425399 | No. 12R Rubber Sprue | US\$ 9.11  | <a href="#">+Info</a> |

### CRUCIBLE FORMERS - WHIP MIX

CRUCIBLE FORMERS - WHIP MIX

Whip Mix CRUCIBLE FORMERS are available in several diameters and large and small oval shapes made of durable plastic.



| Code   | Name  | Price      | +                     |
|--------|---|------------|-----------------------|
| 120293 | WHIP MIX #4081 Crucible Former 1-1/4in - #11908     | US\$ 30.50 | <a href="#">+Info</a> |
| 120294 | WHIP MIX #4179 Crucible Former 1-1/2in - #11940     | US\$ 30.50 | <a href="#">+Info</a> |
| 120295 | WHIP MIX #4177 Crucible Former 1-3/4in - #11967     | US\$ 30.50 | <a href="#">+Info</a> |
| 120296 | WHIP MIX #4390 Crucible Former 2-1/2in - #11975     | US\$ 30.50 | <a href="#">+Info</a> |
| 120297 | WHIP MIX #4215 Crucible Former Conical 2in - #35149 | US\$ 33.55 | <a href="#">+Info</a> |

### FUNNEL FORMERS - BEGO

FUNNEL FORMERS - BEGO

For combination crucible

Funnel former with reservoir for combination crucible



| Code     | Name  | Price      | +                     |
|----------|---|------------|-----------------------|
| 10-52075 | BEGO Funnel Former with Reservoir 10pk #52075 | US\$ 94.05 | <a href="#">+Info</a> |

|          |   |            |                       |
|----------|---|------------|-----------------------|
| 10-52060 | BEGO Funnel Former Standard 10pk #52060               | US\$ 78.04 | <a href="#">+Info</a> |
| 419867   | BEGO UNIVERSAL FUNNEL FORMERS - 100 PIECES - #52068   | US\$ 52.29 | <a href="#">+Info</a> |
| 419866   | BEGO FUNNEL FORMERS FOR NAUTILUS - 10 PIECES - #52066 | US\$ 90.58 | <a href="#">+Info</a> |

## SPRUE FORMERS - VANIMAN

### [SPRUE FORMER - VANIMAN](#)



A cylindrical silicon base which holds the sprues upright in the casting ring during investing. It forms a funnel in the investment to guide the molten material into the mold during casting.

| Code   | Name   | Price     | +                     |
|--------|--|-----------|-----------------------|
| 111391 | VANIMAN Spectrum Sprue Former 1.25in - PURPLE - Each ( #VMC-SB125 )      | US\$ 9.18 | <a href="#">+Info</a> |
| 111392 | VANIMAN Spectrum Sprue Former 1.50in - Dark BLUE - Each ( #VMC-SB150 )   | US\$ 6.39 | <a href="#">+Info</a> |
| 111390 | VANIMAN Spectrum Sprue Former 1.75in - Medium BLUE - Each ( #VMC-SB175 ) | US\$ 6.60 | <a href="#">+Info</a> |
| 111394 | VANIMAN Spectrum Sprue Former 2.5in - YELLOW - Each ( #VMC-SB250 )       | US\$ 7.84 | <a href="#">+Info</a> |
| 111395 | VANIMAN Spectrum Sprue Former 3in - ORANGE - Each ( #VMC-SB300 )         | US\$ 8.12 | <a href="#">+Info</a> |

## Casting Accessories

### Crucible Ringless formers



## Crucible Ringless formers - BEGO

Crucible Ringless formers - BEGO

| Code   | Name                      | Price      | +                     |
|--------|---------------------------|------------|-----------------------|
| 419875 | Rubber Base w/ Funnel 4pk | US\$ 51.23 | <a href="#">+Info</a> |

## CRUCIBLE RINGLESS FORMERS - ECCO

### CRUCIBLE RINGLESS FORMERS - ECCO



## CASTING RINGLESS

| Code      | Name  | Price      | +                     |
|-----------|---|------------|-----------------------|
| 045-32999 | ECCO - Casting Ringless and Crucible Former ( base & ring ) ( #4093 ) - 1.5in Di... | US\$ 21.00 | <a href="#">+Info</a> |
| 045-33006 | ECCO - Castings Ringless and Crucible Formers ( base & ring ) #4094 2inch Diamet... | US\$ 21.00 | <a href="#">+Info</a> |
| 101093    | ECCO - Castings Ringless and Crucible Formers ( base & ring ) - #409 2.5in Diame... | US\$ 21.00 | <a href="#">+Info</a> |
| 114469    | ECCO - Castings Ringless and Crucible Formers ( base & ring ) Set of 3 pcs with ... | US\$ 72.90 | <a href="#">+Info</a> |

## PRESS RINGLESS

| Code   | Name  | Price      | +                     |
|--------|---|------------|-----------------------|
| 100467 | ECCO - Press Ringless And Crucible Formers ( base & ring ) - 90gr - Ringless 1.5... | US\$ 21.00 | <a href="#">+Info</a> |
| 100468 | ECCO - Press Ringless And Crucible Formers ( base & ring ) - 180gr - Ringless 2i... | US\$ 21.00 | <a href="#">+Info</a> |

### MOULD FORMERS - BEGO

#### MOULD FORMERS - BEGO

- \* Eliminates fixing and grinding of the investment models when the BEGO combination duplicating flask is used.
- \* Both mould formers can also be used with all other duplicating systems.

**US DENTAL DEPOT Inc is an Authorized BEGO Dealer in the US.**



| Code     | Name                                     | Price      | +     |
|----------|--|------------|-------|
| 10-52400 | BEGO Mould Formers Large BLUE 4pk #52400 | US\$ 88.09 | +Info |
| 10-52390 | BEGO Mould Formers Small RED 4pk #52390  | US\$ 62.98 | +Info |

### MOULD RINGS - BEGO

#### MOULD RINGS - BEGO

- \* For crown and bridge work
- \* Metal mould ring size 1 for 60 g of investment material.
- \* Long service life thanks to special steel design.
- \* Suitable for all BEGO crown and bridge investment materials.

**US DENTAL DEPOT Inc is an Authorized BEGO Dealer in the US.**



| Code     | Name                                | Price       | +     |
|----------|-------------------------------------|-------------|-------|
| 10-52424 | BEGO Mould Rings size #9 4pk #52424 | US\$ 158.12 | +Info |
| 10-52423 | BEGO Mould Rings size #6 4pk #52423 | US\$ 136.89 | +Info |
| 10-52422 | BEGO Mould Rings size #3 4pk #52422 | US\$ 158.12 | +Info |

### PERFECT RINGS - ADS

#### PERFECT RINGS - ADS

Introducing all new Perfect Ring specifically engineered for use with lithium disilicate restorations. Perfect Ring is made from the same clear material our traditional ringless system for casting alloys. The height and volume dimensions have been specifically designed to fit all current pressing furnaces.

- \* Eliminates investment cracking
- \* Increases green strength of investment
- \* No leaking investment
- \* No excess investment flash
- \* More convenient and saves time
- \* Works great and gives a perfect fitting ring





| Code   | Name   | Price      | +                     |
|--------|--|------------|-----------------------|
| 118269 | ADS Perfect Ring Hot PINK for Lithium Disilicate 100grs 13mm diameter plunger n... | US\$ 38.53 | <a href="#">+Info</a> |
| 118270 | ADS Perfect Ring Hot PINK for Lithium Disilicate 200grs 13mm diameter plunger n... | US\$ 38.53 | <a href="#">+Info</a> |

## RINGLESS RING AND FORMERS - BEGO

### RINGLESS RING AND FORMERS - BEGO

Universally applicable for many casting systems.

#### Description:

- \* Time savings in relation to mould systems with foil sleeve, iron ring, etc..
- \* Compatible with Rapid Wax System.
- \* Minimal wear, thus lower costs than with comparable systems.
- \* Universally applicable for many casting systems.
- \* Easy separation of mould and mould ring.
- \* For all BEGO crown and bridge investment materials.



**US DENTAL DEPOT Inc is an Authorized BEGO Dealer in the US.**

| Code     | Name  | Price      | +                     |
|----------|---|------------|-----------------------|
| 10-52665 | BEGO Ringless Ring & Former 1.25in 1pk ( #52665 ) Universally applicable for ma...  | US\$ 59.58 | <a href="#">+Info</a> |
| 10-52666 | BEGO Ringless Ring & Former 2in 1pk ( #52666 ) Universally applicable for many c... | US\$ 59.58 | <a href="#">+Info</a> |
| 10-52667 | BEGO Ringless Ring & Former 2.5in 1pk ( #52667 ) Universally applicable for many... | US\$ 59.58 | <a href="#">+Info</a> |

## Casting Accessories

### Crucibles all types

## Crucibles all types - BEGO

Crucibles all types - BEGO

| Code     | Name                           | Price      | +                     |
|----------|--------------------------------|------------|-----------------------|
| 10-52477 | Crucible Spray/Lubricant 100ml | US\$ 55.50 | <a href="#">+Info</a> |

## Crucibles all types - ECCO

Crucibles all types - ECCO

| Code   | Name                | Price      | +                     |
|--------|---------------------|------------|-----------------------|
| 421028 | Hooded Crucible 7oz | US\$ 32.16 | <a href="#">+Info</a> |

### Crucibles all types - BUFFALO

Crucibles all types - BUFFALO



| Code      | Name                           | Price      | +                     |
|-----------|--------------------------------|------------|-----------------------|
| 009-28012 | Crucible Jumbo Clay Round Hole | US\$ 22.12 | <a href="#">+Info</a> |
| 425396    | Type 'A' Quartz cruc           | US\$ 25.48 | <a href="#">+Info</a> |
| 425397    | Type 'A' clay crucib           | US\$ 22.12 | <a href="#">+Info</a> |

### Crucibles all types - KERR

Crucibles all types - KERR

| Code   | Name                           | Price      | +                     |
|--------|--------------------------------|------------|-----------------------|
| 101453 | Zircon Alumina Crucible Yellow | US\$ 30.87 | <a href="#">+Info</a> |

### CRUCIBLES FOR FORNAX - BEGO

#### CRUCIBLES FOR FORNAX - BEGO

The new Bego ceramic crucibles for Fornax set the highest standards. An innovative manufacturing procedure for crucibles that are resistant high temperature, developed in cooperation with scientific institutes produces an extremely homogenous material structure contributing to consistently reproducible accuracy of castings.



**US DENTAL DEPOT is an Authorized BEGO dealer in US.**

| Code     | Name  | Price       | +                     |
|----------|---|-------------|-----------------------|
| 10-52455 | BEGO Crucible Ceramic Inserts for Fornax 6pk #52455 #52455  | US\$ 131.67 | <a href="#">+Info</a> |
| 10-52843 | BEGO NEW Crucible for Fornax 6pk #52483                     | US\$ 148.28 | <a href="#">+Info</a> |
| 10-52454 | BEGO Crucible Graphite Inserts for Fornax 6pk #52454 #52454 | US\$ 170.88 | <a href="#">+Info</a> |




## CRUCIBLES FOR FUNDOR-T - BEGO

### CRUCIBLES FOR FUNDOR-T - BEGO

Ceramic Crucibles for cobalt-chrome partial denture alloys and non-precious alloys.

The melting crucibles are made of special ceramic material and have a long service life through a high heat change resistance.

For Melting crucible, deep trough, shallow trough, for precious metal, for cobalt-chrome partial denture alloys and non-precious alloys that can be used with Fundor T.



**US DENTAL DEPOT is an Authorized BEGO dealer in US.**

| Code     | Name   | Price       | +                     |
|----------|--|-------------|-----------------------|
| 10-52426 | BEGO Crucible for Co-Cr + Non Precious Alloys for Fundor T 6pk #52426 #52426 | US\$ 198.91 | <a href="#">+Info</a> |
| 10-52425 | BEGO Crucible for Precious Alloys for Fundor T 6pk #52425 #52425             | US\$ 198.91 | <a href="#">+Info</a> |

## CRUCIBLES TYPE - KEYSTONE

### CRUCIBLES TYPES - KEYSTONE

**Features:**

\* Most popular Crucibles Types that fit standard and different types of casting machines.

\* Types :

Centrifugal

Standard Hooded

Standard Slotted

Two Stage Slotted

\* Generally Sold Individually.



**US Dental Depot Is An Authorized Keystone Dealer In US.**

### CENTRIFUGAL

#### KEYSTONE - CENTRIFUGAL

**Features:**

\* This Centrifugal Type Crucible fits standard Dentsply type casting machines. Sold individually.

| Code        | Name  | Price      | +                     |
|-------------|---|------------|-----------------------|
| 034-5920250 | KEYSTONE Crucible Centrifugal Type - #5920250 | US\$ 23.55 | <a href="#">+Info</a> |

### STANDARD HOODED

#### KEYSTONE - STANDARD HOODED

**Features:**

- \* This high quality Standard Hooded Type Crucible offers optimum performance in induction furnaces.
- \* This crucible fits standard type Kerr, Lucas and Emesco casting machines.

| Code        | Name   | Price      | +     |
|-------------|--|------------|-------|
| 034-5920235 | KEYSTONE Crucible Hooded Standard - #5920235 | US\$ 22.95 | +Info |

### STANDARD SLOTTED

#### KEYSTONE - STANDARD SLOTTED

**Features:**

- \* These high quality Standard Slotted Type Crucibles offer optimum performance.
- \* This crucible fits standard type Kerr and Emesco casting machines.

| Code        | Name  | Price      | +     |
|-------------|---|------------|-------|
| 034-5920200 | KEYSTONE Crucible Slotted Standard - #5920200                                       | US\$ 22.95 | +Info |
| 034-5920201 | KEYSTONE #32 Standard Slotted Type Crucible - for European Machines Same As #30 ... | US\$ 24.60 | +Info |

### TWO STAGE SLOTTED

#### KEYSTONE - TWO STAGE SLOTTED CRUCIBLE

**Features:**

- \* This Two State Slotted Crucible fits standard type Kerr and Emesco casting machines.

| Code        | Name   | Price      | +     |
|-------------|--|------------|-------|
| 034-5920230 | KEYSTONE Crucible 2 Stage Slotted - #5920230 | US\$ 21.75 | +Info |

### HI-HEAT QUARTZ CRUCIBLES - B

### ESQUAL

#### HI-HEAT QUARTZ CRUCIBLES - BESQUAL

The large, flat melting area allows more alloy of various shapes to be heated evenly, an important factor in the proper melting of nickel based alloys. The design of this crucible reflects heat for more uniform melting. For broken arm casting machines.



| Code   | Name  | Price     | +     |
|--------|---|-----------|-------|
| 100126 | BESQUAL Hi Heat Quartz Standard Crucible Slotted (Broken Arm Casting) . - Kerr /... | US\$ 7.20 | +Info |
| 100044 | BESQUAL Hi Heat Quartz Standard Hooded Crucible Round Orifice (QHN-2/719-000)- F... | US\$ 7.20 | +Info |
| 116093 | BESQUAL Hi Heat Quartz Mini Crucible Hooded - Round Orifice - For Broken Arm Cas... | US\$ 7.20 | +Info |

**INDUCTION CRUCIBLE - KEYSTONE**  
**INDUCTION CRUCIBLE - KEYSTONE**



**Features:**

\* This high quality Induction Crucible offers optimum performance in induction furnaces.

**US Dental Depot Is An Authorized Keystone Dealer In US.**

| Code   | Name  | Price      | +     |
|--------|---|------------|-------|
| 113656 | KEYSTONE Induction Crucible II each Fits standard type Howmedica and electronic ... | US\$ 24.60 | +Info |
| 113655 | KEYSTONE Induction Crucible Hooded #49A - Induction Crucible Hooded #49A #5920245   | US\$ 24.77 | +Info |
| 113654 | KEYSTONE Induction Crucible - Fits standard type Tycomatic and Unitek casting ma... | US\$ 22.35 | +Info |

**NAUTILUS CRUCIBLES - BEGO**  
**NAUTILUS CRUCIBLES - BEGO**

BEGO has pushed the envelope once again with the newly designed Nautilus Crucibles. Automatic casting for consistent results now has new options such as adjustable melt hold times and optional longer pressure settings after casting - these features give the technician greater control in casting larger implant bars and full arch C&B frames. The new design has a large color touch screen for easy operation, energy saving "Eco-Mode" and integrated cooling system making the new Nautilus Crucibles the intelligent environmentally friendly choice. The induction heating guarantees a shorter melting cycle and minimizes oxidation which facilitates finishing time. Let BEGO help you reach your C&B goal of consistent castings by learning more about the proven performance of vacuum/pressure, split crucible technology of the Nautilus Crucible induction casting machine.



| Code     | Name   | Price       | +     |
|----------|--|-------------|-------|
| 10-52473 | BEGO Nautilus Glass Inserts for crucibles 4pk #52473 #52473              | US\$ 317.30 | +Info |
| 10-52488 | BEGO Crucible Ceramic for all Nautilus casting models 4pk #52488 #-52488 | US\$ 184.90 | +Info |
| 10-52468 | BEGO Nautilus Graphite Inserts for crucibles 6pk #52468 #52468           | US\$ 130.75 | +Info |

### Casting Accessories

## Debubblizers & Wetting Agents



### CERPRESS SL TENSION RELIEF SPRAY - ADS

#### CERPRESS SL TENSION RELIEF SPRAY - ADS

- \* Cerpress SL Pressable Ceramic System is used by thousands of dental laboratories throughout the country.
- \* It has been on the market over 6 years and been clinically tested by several organizations.



| Code   | Name  | Price      | +     |
|--------|---|------------|-------|
| 117712 | ADS Cerpress SI Tension Relief Spray 4oz. - #201045 | US\$ 18.17 | +Info |

### CLEARLY DEBUBBLIZER - ADS

#### CLEARLY DEBUBBLIZER - ADS

#### Features:

- \* Advanced formula Debubbler for the new resin-toughened stone as recommended by their manufacturers
- \* Completely water soluble, crystal clear and alcohol-free
- \* Produces perfect castings
- \* Investments flow smoothly across wax patterns
- \* Essential on the new alginate/silicone hybrid and chromatic impression materials
- \* Clearly Debubbler can be applied by painting, spraying or dipping
- \* Thinnest coating gives the best results



| Code   | Name   | Price      | +     |
|--------|--|------------|-------|
| 117753 | ADS Debubbler CLEARLY 1 gallon . - #D926-8   | US\$ 57.63 | +Info |
| 117752 | ADS Debubbler CLEARLY 1qt ( 32oz ) - #D926-6 | US\$ 20.68 | +Info |
| 117751 | ADS Debubbler CLEARLY pt 16oz. - #D926-5     | US\$ 13.67 | +Info |
| 117750 | ADS Debubbler CLEARLY / SPRAY 8oz. - #D926-4 | US\$ 11.32 | +Info |
| 117749 | ADS Debubbler CLEARLY - SPRAY 4oz. - #D926-3 | US\$ 8.76  | +Info |

### DEBUBBLIZER AND WAX PATTERN CLEAR - ADS

#### DEBUBBLIZER AND WAX PATTERN CLEAR - ADS

#### Features:

- \* Formulated for investment to flow evenly over the wax pattern
- \* Eliminates bubbles and reduces the surface tension of waxes
- \* Produces castings with a satin smooth finish
- \* Spray on all impression materials before pouring in stone
- \* Cleans epoxy dies of oils and particles



| Code   | Name   | Price      | +     |
|--------|--|------------|-------|
| 117758 | ADS Debubblizer & Wax PATTERN CLEANER gal. - #D916-8         | US\$ 52.23 | +Info |
| 117757 | ADS Debubblizer & Wax PATTERN CLEANER qt - #0916-6           | US\$ 20.66 | +Info |
| 117756 | ADS Debubblizer & Wax PATTERN CLEANER pt 16oz. - #D916-5     | US\$ 13.67 | +Info |
| 117755 | ADS Debubblizer & Wax PATTERN CLEANER 8oz. - spray - #D916-4 | US\$ 11.32 | +Info |
| 117754 | ADS Debubblizer & Wax PATTERN CLEANER 4oz. - spray - #0916-3 | US\$ 8.78  | +Info |

### Debubblizers & Wetting Agents - BEGO

Debubblizers & Wetting Agents - BEGO

| Code   | Name                         | Price      | +     |
|--------|------------------------------|------------|-------|
| 105064 | Aurofilm Wetting Agent 3.5oz | US\$ 25.36 | +Info |

### Debubblizers & Wetting Agents - KEYSTONE

Debubblizers & Wetting Agents - KEYSTONE

| Code        | Name                             | Price      | +     |
|-------------|----------------------------------|------------|-------|
| 034-6121900 | Blue Debubblizer 1 Gal           | US\$ 31.88 | +Info |
| 445889      | Mizzy Wax&Impress Debubbliz 32oz | US\$ 16.44 | +Info |

### PICOSILK SPRAY BOTTLE - RENFER

#### PICOSILK BOTTLE - RENFERT

Picosilk is a wetting agent for wax models within crown, bridge and model casting technique. For tension-free wax patterns

#### Advantages

- \* Spray on finely using the atomizer before investing.
- \* Wax does not take on water-repellent characteristics after use.



\* Improves the flow of investment material.

| Code          | Name  | Price      | +     |
|---------------|---|------------|-------|
| 023-1550-0100 | RENFERT - Picosilk Refill Bottle-500ml                      | US\$ 48.94 | +Info |
| 105083        | RENFERT - Picosilk Agent w/spray 75ml- #1550-0000 #15500000 | US\$ 23.41 | +Info |

## SIMPLY SURFACTANT - ADS

### SIMPLY SURFACTANT - ADS

The superior wetting agent. A modern surface tension eliminator designed for the new resin toughened stones and silicone impression materials. Formulated for investment to flow smoothly over wax patterns and hydrophobic impressions. Easy to use, simply paint, spray or dip wax patterns and impressions. Drain or blow off excess liquid. Produces defect-free pours with a fine, smooth finish and improved surface detail.



| Code   | Name   | Price      | +     |
|--------|--|------------|-------|
| 118377 | ADS Simply Surfactant LIQUID Wetting AGENT 8oz. - #D 936-5         | US\$ 14.77 | +Info |
| 118378 | ADS Simply Surfactant LIQUID Wetting AGENT pt 16oz. - #D 936-6     | US\$ 20.66 | +Info |
| 118379 | ADS Simply Surfactant LIQUID Wetting AGENT 1qt ( 32oz ) - #S 650-2 | US\$ 35.50 | +Info |

## SMOOTHX DEBUBBLING SOLUTION - WHIP MIX

### SMOOTHX DEBUBBLING SOLUTION -WHIP MIX

A surface tension reducing liquid for treating patterns to permit better adherence and flow of investment.

Also excellent for treating synthetic impressions (Polysulfide, Silicone, or Polyether) before pouring casts.

Refillable, fine-spray bottles are easy to use and help keep debubblizer clean by eliminating the need for dipping brushes or patterns.



| Code   | Name  | Price       | +     |
|--------|---|-------------|-------|
| 120431 | WHIP MIX Smoothex Debubblizer Gallon/3.75l - #12912 | US\$ 330.25 | +Info |
| 120429 | WHIP MIX Smoothex Debubblizer Liter - #29394        | US\$ 103.57 | +Info |
| 120428 | WHIP MIX Smoothex Pint/450ml - #09849               | US\$ 55.80  | +Info |
| 120432 | WHIP MIX Smoothex Fine Spray Btl 8oz/225ml - #02313 | US\$ 46.21  | +Info |
| 120430 | WHIP MIX Smoothex Fine Spray Bottle 4oz - #29386    | US\$ 25.75  | +Info |

## Casting Accessories

### Liners 4 Casting Rings



### LINING TAPE - ADS

#### LINING TAPE - ADS

Type 380 ceramic fiber tape is asbestos-free and made from high purity refractory fibers (not cellulose) with a melting point of 3400F. Lining Tape is used as a direct replacement for asbestos ring liners. Available in convenient sizes to fit all casting rings. Avoids asbestos hazard. Use with all investments including high temperature types.



| Code   | Name  | Price      | +     |
|--------|---|------------|-------|
| 118134 | ADS Lining Tape Ceramic 1 3/4in x 1/16in x 50in roll - #T380-6 casting ring line... | US\$ 29.36 | +Info |
| 118135 | ADS Lining Tape Ceramic 2in x 1/32in x 50in roll - #T380-7 casting ring liners ]    | US\$ 27.09 | +Info |
| 118132 | ADS Lining Tape Ceramic 1 1/2in x 1/16in x 50in roll - #T380-4 casting ring line... | US\$ 25.74 | +Info |
| 118133 | ADS Lining Tape Ceramic 1 3/4in x 1/32in x 50in roll - #T380-5 casting ring line... | US\$ 25.58 | +Info |
| 118130 | ADS Lining Tape Ceramic 1 1/4in x 1/16in x 50in roll - #T380-2 casting ring line... | US\$ 23.93 | +Info |
| 118131 | ADS Lining Tape Ceramic 1 1/2in x 1/32in x 50in roll - #T380-3 casting ring line... | US\$ 22.42 | +Info |
| 118129 | ADS Lining Tape Ceramic 1 1/4in x 1/32in x 50in roll - #T380-1 casting ring line... | US\$ 18.07 | +Info |

### RING LINER - NON-ABESTOS CELLULOSE - WHIP MIX

#### RING LINER - NON-ABESTOS CELLULOSE - WHIP MIX

A liner that has no asbestos and none of its inherent hazards. This material has the necessary cushioning, but more importantly the correct water absorption quality which feeds to the investment to provide and control setting expansion. During burnout it chars to a black ash which crumbles easily when the mold is pushed out of the ring.

Liners should be 3 (1/8") short of each end of the casting ring.



| Code   | Name  | Price      | +     |
|--------|---|------------|-------|
| 120454 | WHIP MIX Ring Liner Roll Cellulose 3in - #14656       | US\$ 28.94 | +Info |
| 120450 | WHIP MIX Ring Liner Roll Cellulose 50mm - #16195      | US\$ 25.23 | +Info |
| 120453 | WHIP MIX Ring Liner Roll Cellulose 1-5/8in - #14648   | US\$ 21.18 | +Info |
| 120452 | WHIP MIX Ring Liner Roll Cellulose 1-3/8 X60 - #14621 | US\$ 18.48 | +Info |
| 120451 | WHIP MIX Ring Liner Roll Cellulose 1-1/8 X60 - #14605 | US\$ 16.77 | +Info |

**FREE SHIPPING**  
FOR ALL DOMESTIC ACCOUNTS!



**US  
DENTAL  
DEPOT**

+1 954 874 6325   
www.usdentaldepot.com 

### RING LINERS - KEYSTONE

#### RING LINERS - KEYSTONE

Our Casting Ring Liner is made of refractory fibers that have a melting point in excess of 3,200F. The Liner is asbestos free and sold in 50-foot rolls. Various available sizes.

**Features:**

- \* Asbestos free
- \* 50 foot rolls

**Available Sizes:** 1-1/4" H x 1/32" Thick, 1-1/4" H x 1/16" Thick, 1-1/2" H x 1/32" Thick, 1-1/2" H x 1/16" Thick, 1-3/4" H x 1/32" Thick, 1-3/4" H x 1/16" Thick, 2" H x 1/32" Thick, 2" H x 1/16" Thick, 3" H x 1/16" Thick



| Code   | Name   | Price      | +     |
|--------|--|------------|-------|
| 101147 | KEYSTONE Casting Ring Liner 1 1/4in Ht. x 1/32in Thick. Asbestos-Free, Single... | US\$ 17.55 | +Info |
| 101148 | KEYSTONE Casting Ring Liner 1 1/4in Ht. x 1/16in Thick. Asbestos-Free, Single... | US\$ 23.03 | +Info |
| 101149 | KEYSTONE Casting Ring Liner 1 1/2in Ht. x 1/32in Thick. Asbestos-Free, Single... | US\$ 20.55 | +Info |
| 101150 | KEYSTONE Casting Ring Liner 1 1/2in Ht. x 1/16in Thick. Asbestos-Free, Single... | US\$ 27.98 | +Info |
| 101151 | KEYSTONE Casting Ring Liner 1 3/4in Ht. x 1/32in Thick. Asbestos-Free, Single... | US\$ 28.65 | +Info |
| 101152 | KEYSTONE Casting Ring Liner 1 3/4in Ht. x 1/16in Thick. Asbestos-Free, Single... | US\$ 37.13 | +Info |
| 101153 | KEYSTONE Casting Ring Liner 2in Ht. x 1/32in Thick. Asbestos-Free, Single 50 ... | US\$ 33.60 | +Info |
| 101154 | KEYSTONE Casting Ring Liner 2in Ht. x 1/16in Thick. Asbestos-Free, Single 50 ... | US\$ 39.68 | +Info |
| 101155 | KEYSTONE Casting Ring Liner 3in Ht. x 1/16in Thick. Asbestos-Free, Single 50 ... | US\$ 52.28 | +Info |

### TAPE EEZ - ADS

#### TAPE EEZ - ADS

Super concentrated liquid for wetting ADS casting ring lining tapes. One pint makes 5 quarts. Designed to create the workability and handling texture of asbestos. Expansion properties of the investment will not be altered. Burns out clean in the furnace without charring.



| Code   | Name   | Price      | +     |
|--------|--|------------|-------|
| 118411 | ADS Tape Eez qt - #N 928-6 ( wetting agent for ring liners )       | US\$ 23.48 | +Info |
| 118410 | ADS Tape Eez 16oz. pt - #T 770.6 ( wetting agent for ring liners ) | US\$ 18.88 | +Info |

## Catalog Index / Indice

|                                     |    |
|-------------------------------------|----|
| Alloy rings .....                   | 1  |
| Casting Accessories .....           | 2  |
| Casting tongs .....                 | 3  |
| Crucible formers .....              | 3  |
| Crucible Ringless formers .....     | 5  |
| Crucibles all types .....           | 7  |
| Debubblizers & Wetting Agents ..... | 12 |
| Liners 4 Casting Rings .....        | 14 |