

## Chissels electric and neumatic

### Laboratory Chissels & Acc



#### Laboratory Chissels & Acc - BUFFALO

Laboratory Chissels & Acc - BUFFALO

Code	Name	Price	+
425262	Pearson No. 18 Chise	US\$ 7.76	+Info
425263	Pearson No. 16 Chise	US\$ 7.76	+Info
425264	Pearson No. 17 Chise	US\$ 7.76	+Info

#### Laboratory Chissels & Acc - ESSENCE

Laboratory Chissels & Acc - ESSENCE

Code	Name	Price	+
100776	Neumatic Chisel + 2 tips	US\$ 241.50	+Info

#### Laboratory Chissels & Acc - RENFERT

Laboratory Chissels & Acc - RENFERT

Code	Name	Price	+
417644	Fastener Power Pillo	US\$ 26.51	+Info
417663	Pillo Hose clamp	US\$ 5.65	+Info
417673	Valve for new pillos	US\$ 40.06	+Info
417674	US air fitting for Pillos	US\$ 21.12	+Info
417675	Rubber Grip for Pillo 4000	US\$ 65.55	+Info
417677	Valve for power pillo SET	US\$ 142.12	+Info
421065	Quick Connect -For Power-Pillo	US\$ 46.65	+Info

#### PILLO CHISEL - RENFERT

##### PILLO CHISEL - RENFERT

##### Function & Performance

Stresses in bridges, cracks in acrylic dentures and splintered teeth are things of the past.

Laboratory-fabricated restorations are gently, easily and quickly deflasked using the Power pillo and Pillo pneumatic chissels.

##### Advantages

- \* Save time and money with powerful and gentle devesting of sensitive restorations (castings, dentures).
- \* Absolutely maintenance-free due to special oil-free mechanism.
- \* Maximum health protection provided by minimal transmission of vibration to the joints



### Features

The Power pillo delivers a higher performance, smoother running and a quick-release coupling connection

### Details

- \* Infinitely variable power setting.
- \* Four different tungsten carbide chisel inserts are optionally available.
- \* The ergonomic form provides an optimum, reliable grip and a high degree of working comfort.

### #UsddRenfert

Code	Name	Price	+
114041	RENFERT Power Pillo deflasking Chisel Pneumatic chisel- #5022-5000 #50225000	US\$ 678.38	+Info
120734	RENFERT Pillo- Pneumatic Chisel - #50224000	US\$ 491.04	+Info

## POINTED CHISEL - RENFERT

### POINTED CHISEL - RENFERT

#### Feature

Very long service life and highly durable



Code	Name	Price	+
023-5022-0200	RENFERT Flat Chisel-Wide-For Power-Pillo-Each #5022-0200 #50220200	US\$ 80.37	+Info
023-5022-0100	RENFERT Fiat Chisel-Narrow-For Power-Pillo-Each- #5022-0100 #50220100	US\$ 81.21	+Info
023-5022-0400	RENFERT Pointed Chisel For Power-Pillo-Each- #5022-0400 #50220400	US\$ 82.03	+Info
023-5022-0300	RENFERT Channeled Chisel For Power-Pillo-Each-R #5022-0300 #50220300	US\$ 82.03	+Info
023-5021-0400	RENFERT Fiat Chisel-Extra Wide-For Old Style Pillo-Each-#5021-0400 #50210400	US\$ 60.25	+Info

## Catalog Index / Indice

Laboratory Chissels & Acc ..... 1