

Die Materials

Die Hardenders



CLEAR COAT - ADS

CLEAR COAT HARDENER & SILVER ? ADS

Features:

- * Provides instant protection for stone models and die materials
- * Penetrates dies and creates a clear, hard, glossy, tough surface
- * Prevents abrasion and shipping of dies
- * Protects margins and critical areas
- * Dries instantly and is waterproof
- * One bottle protects hundreds of dies
- * Paint on or dip, will not change die dimensions
- * Compatible with ADS Die Spacer
- * Removes Clear Coat from delicate dies without distortion or damage



HARDENER

Code	Name	Price	+
117717	ADS CLEAR Coat Hardener refill qt - #C760-6	US\$ 90.77	+Info
117715	ADS CLEAR Coat Hardener refill 4oz. - #C760-2	US\$ 25.09	+Info

SOLVENT

Code	Name	Price	+
117721	ADS CLEAR Coat Solvent refill 1 oz. - #C760-12	US\$ 4.23	+Info

Die Hardenders - BEGO

Die Hardenders - BEGO

Code	Name	Price	+
10-52008	Durofluid Model Spray 100ml	US\$ 29.46	+Info

Die Hardenders - KEYSTONE

Die Hardenders - KEYSTONE

Code	Name	Price	+
446232	Keystone Die Hardener Yeti, Complet	US\$ 17.02	+Info
446444	Keystone Die Hardener Yeti, Complet	US\$ 133.29	+Info

Die Hardenders - RENFERT

Die Hardenders - RENFERT

Die:Master

Die:Master

Code	Name	Price	+
1956-0600	DieMaster Duo Hardener 15ml	US\$ 15.96	+Info

Other Models

Other Models

Code	Name	Price	+
113841	Hardening Bath Sealing 125ml	US\$ 35.12	+Info

Die Hardenders - TAUB

Die Hardenders - TAUB

Code	Name	Price	+
464573	Taub Stone and Die Hardener Kit	US\$ 36.84	+Info

Die Hardenders - YETI

Die Hardenders - YETI

Code	Name	Price	+
34-1860046	Yeti Die Hardener 18ml	US\$ 20.01	+Info

P.D.Q. KIT - WHIP MIX

P.D.Q. KIT - WHIP MIX

A Complete System for Preparing Dies

The **P.D.Q. Kit** is designed to facilitate the production of conventional crown and bridges and refractory dies for all-porcelain restorations. Each kit contains half-ounce bottles of stone hardener, blue and grey spacer liquids and lubricant, each with a convenient brush applicator and packaged in a functional plastic organizer for convenient access. The kit also includes a one-ounce bottle of thinner formulated for use with both the space liquids and the hardener. Each item in the kit is available separately.

P.D.Q. Hardener a liquid polymer that is painted on the stone die. It soaks into the stone pores to seal, harden and waterproof. It protects margins, yet does not change the dimensions. This hardener also prevents water absorption from porcelain during build-up.



P.D.Q. Spacer liquids each provide a 12 micron space in a single coating. Alternating blue and grey colors ensure thickness -- two applications provide the optimum 24 microns of space to relieve dies. The liquids are extremely fast setting to maximize efficiency and are highly visible.

P.D.Q. Lubricant is a water soluble lubricant with low surface tension that will form an extremely thin, even film. It facilitates the removal of wax patterns.

P.D.Q. Thinner may be used with both the Spacer liquids and the Hardener as required to restore optimum handling properties. Thinner replenishes the liquid that evaporates if the bottle is left uncovered.

Code	Name	Price	+
120422	WHIP MIX Pdq Kit - #10839	US\$ 127.16	+Info
120425	WHIP MIX Pdq Spacer GRAY .5 Oz/14ml - #10804	US\$ 33.03	+Info
120426	WHIP MIX Pdq Spacer BLUE .5 Oz/14ml - #10790	US\$ 33.03	+Info
120427	WHIP MIX Pdq Hardener .5 Oz/14ml - #10626	US\$ 33.03	+Info
120423	WHIP MIX Pdq Thinner 28ml - #10820	US\$ 17.96	+Info
120424	WHIP MIX Pdq Lubricant .5 Oz/14ml - #10812	US\$ 13.22	+Info

Die Materials

Epoxi & Plastic die materials



Epoxi & Plastic die materials - KEYSTONE

Epoxi & Plastic die materials - KEYSTONE



Code	Name	Price	+
446280	Tri-Dynamics Tri-Epoxy Bulk 1lb	US\$ 13.61	+Info
446588	Keystone Tri-Dynamics Tri-Epoxy Bul	US\$ 13.61	+Info
446928	Keystone Tri-Dynamics Tri-Epoxy Bul	US\$ 13.61	+Info
446929	Keystone Tri-Dynamics Tri-Epoxy Bul	US\$ 11.41	+Info
446930	Keystone Tri-Dynamics Tri-Epoxy Bul	US\$ 13.61	+Info
474236	Tri-Dynamics Hardener -	US\$ 15.13	+Info

EPOXY RESIN - ADS

EPOXY RESIN - ADS

- * Provides excellent dimensional stability
- * Can be trimmed and sawed with ease, yet is extremely hard
- * Holds delicate details for consistent accuracy
- * Bubble reducing formula is available in all popular colors
- * Color are carefully selected to facilitate reading of margin details and prevent eye strain and fatigue
- * Dies can be demolded and used in 3 to 3 1/2 hours
- * No need to centrifuge



10 GRAMS

Code	Name	Price	+
117779	ADS Die Epoxy 10 gram kit - BLUE 12 units (120grams in total) - #8000-10-BE	US\$ 24.72	+Info
117784	ADS Die Epoxy 10 gram kit - WHITE 12 units - #8000-10-WE	US\$ 22.25	+Info

BULK EPOXY RESIN ONLY

1 LB

Code	Name	Price	+
118445	ADS Die Epoxy Bulk 1lb. - #1500-110	US\$ 14.75	+Info

EPOXI ACCESORIES

Code	Name	Price	+
117837	ADS Die Epoxy Hardener JAR - refill 4oz. - #8000-220	US\$ 5.21	+Info

PLASTICAL MODEL & DIE MATERIAL KIT - ADS

PLASTICAL MODEL & DIE MATERIAL KIT - ADS

A low viscosity liquid plastic, that cures from liquid to an extremely hard plastic model in three minutes. Excellent for small pours. Can be demolded and worked on in 15 minutes. No detectable shrinkage, smell or gumming of burs. No messy syringes needed with Plactical.



COLOR GREEN

Code	Name	Price	+
118295	ADS Plactical - GREEN 2 - 4oz. Btls. 8oz. kit - #PL-80-GN	US\$ 32.77	+Info

COLOR TAN

Code	Name	Price	+
118297	ADS Plactical - Tan 2 - 4oz. Btls. 8oz. kit - #PL-80-TN	US\$ 32.77	+Info

COLOR WHITE

Code	Name	Price	+
118300	ADS Plactical - WHITE 2 - 8oz. Btls. 16oz. kit - #PL-16-WE	US\$ 66.26	+Info
118299	ADS Plactical - WHITE 2 - 4oz. Btls. 8oz. kit - #PL-80-WE	US\$ 40.12	+Info

Die Materials

Gloss Dip for Models



HI-GLOSS - MODEL GLOSSING LIQUID - TAUB

HI-GLOSS - MODEL GLOSSING LIQUID - TAUB

A liquid-wax glaze you paint on any wax surface to create a smooth, high gloss. Apply to wax up veneers, cast crowns, etc. These products Eliminates scratches and burns out clean. Makes castings smoother. Use on Wax Denture Try-Ins and Bite Blocks, Stone or Plaster.



8 OZ

Code	Name	Price	+
122470	TAUB - Hi-gloss - Liquid Wax Glaze - 8 Oz - #16-2000	US\$ 22.97	+Info

16 OZ

Code	Name	Price	+
122471	TAUB - Hi-gloss - Liquid Wax Glaze - 16 Oz - #16-2001	US\$ 41.06	+Info

32 OZ

Code	Name	Price	+
122472	TAUB - Hi-gloss - Liquid Wax Glaze - 32 Oz - #16-2002	US\$ 63.02	+Info

1 GALLON

Code	Name	Price	+
122473	TAUB - Hi-gloss - Liquid Wax Glaze - 1gal - #16-2003	US\$ 198.91	+Info

MODEL GLOW SOAP - WHIP MIX

MODEL GLOW SOAP - WHIP MIX

MODEL GLOW gives an incredibly high luster to gypsum models. Its special formulation does not yellow or flake and seals the pores to provide a protective coating. Its thin consistency does not affect occlusal accuracy. Long lasting and can be used full strength or diluted. Models can be buffed wet or dry.



Code	Name	Price	+
120421	WHIP MIX Model Glow Gallon/3.75l - #29009	US\$ 74.29	+Info
120420	WHIP MIX Model Glow 4gal Case - #29017	US\$ 215.73	+Info

HATHO
POLISHING KIT



Thermoplastics
like *never* before

Packing Diamond D®
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US
DENTAL
DEPOT

+1 954 874 6325
www.usdentaldepot.com



SATIN DIP - ADS

SATIN DIP - ADS

Use as a dip or paint-on. Apply Satin Dip to study models and buff them with a soft cloth to an impressive, high satiny finish for presentations.



Code	Name	Price	+
118371	ADS Satin Dip MODEL GLOSS 16oz. pt - #S401-5	US\$ 11.83	+Info
118372	ADS Satin Dip MODEL GLOSS qt - #S401-6	US\$ 13.64	+Info
118373	ADS Satin Dip MODEL GLOSS gal - #S401-8	US\$ 31.52	+Info

Die Materials

Wax Glaze Polishing



WAX SHINE - ADS

WAX SHINE - ADS

Spreads easily on all dental waxes. Dries in minutes to a high shine without buffing. Masks fine scratches in the wax surface. The practical alternative to manual polishing of wax-ups, bite blocks, and wax try-ins. Forms a crystal clear, durable finish. Strengthens and hardens thin wax copings. Safe, ready to use, cleans up with water.



Code	Name	Price	+
118439	ADS Wax Shine Finishing 16oz. pt - #W 220-5	US\$ 23.10	+Info
118440	ADS Wax Shine Finishing qt - #W 220-6	US\$ 32.98	+Info

Catalog Index / Indice

Die Hardeners	1
Epoxi & Plastic die materials	3
Gloss Dip for Models	5
Wax Glaze Polishing	7