

Die Spacer

Light Cure Die Spacer



LC DIE SPACERS - YETI

YETI - LIGHT CURING - LIGHT SPACER



- Advantages of the light curing die spacer LIGHT SPACER:
- * homogeneous curing
 - * creates a very hard, scratch-resistant and durable surface
 - * extremely smooth surface - no problem when lifting the coping
 - * excellent adhesion on plaster base
 - * unlimited application time cures only when put into the light curing unit

Curing in any available standard unit

Curing time: approx. 2 - 5 min
Delivery: in black UVA protection bottles of 10 ml

Die LIGHT SPACER can only be used with each other! Not to be mixed with other Spacers!

Code	Name	Price	+
34-1860213	YETI Light Curing liquid Die Spacer - 3 micron CLEAR - 10ml bottle - K#1860213 #...	US\$ 27.99	+Info
34-1860214	YETI Light Curing liquid - RED Die Spacer - 10 micron RED - 10ml bottle - #520-2...	US\$ 27.99	+Info
34-1860215	YETI Light Curing liquid - BLUE Die Spacer - 10 micron BLUE - 10ml bottle - #ART...	US\$ 26.67	+Info

Die Spacer

Paint on Die Spacer



DIE:MASTER - SYSTEMATIC DIE PREPARATION - RENFERT

DIE:MASTER - SYSTEMATIC DIE PREPARATION - RENFERT

Function & Performance

With the die:master system, Renfert offers air-drying die varnish in four layer thicknesses (10, 13, 15 and 20 um) and six opaque colors.

Die:Master Set

- * Faster and clearer access to all components of the system thanks to the stand.
- * Practical, firm base for drying all types of die.
- * Quick and clean storage with supplied high quality separating brush with magnetic fixation.

Die:Master Die Spacers



- * Extremely steam-proof and scratch-resistant thanks to high density, glass-like die spacer layer.
- * Generates maximum bond strength between hardener and die spacer due to primer technology.
- * Die spacer will not puddle on the brush or die, yet it sticks well to line angles due to its extremely thixotropic flowability.

Die:Master Duo

- * Cross links the die stone structure for in-depth hardness and high edge-strength.
- * Maximum bond strength thanks to chemical bond to die spacer layer.
- * 0 um thickness for a 100% dimensionally-stable die.

Die:Master Aqua Thinner

- * Long-term, sustainable quality,- even after frequent use.
- * Controlled regulation of the viscosity because it can be dispensed accurately drop-by-drop.

Picosep

Silicone based separator completes the die:master system components in an optimal way.

#UsddRenfert

STAND KITS

Die:Master

Die:Master

Code	Name	Price	+
120705	Die:master Set (ivory)	US\$ 147.41	+Info
1956-0800	DieMaster Stand Kit - red	US\$ 57.86	+Info
1956-0900	DieMaster Stand Kit - Blue	US\$ 57.86	+Info
1956-1000	DieMaster Stand Kit Gold/Silve	US\$ 81.44	+Info

Other Models

Other Models

Code	Name	Price	+
417582	Die:master set 2xAqua SET	US\$ 160.72	+Info

Refills

Die:Master

Die:Master

Code	Name	Price	+
1956-0700	DieMaster Thinner 15ml	US\$ 15.96	+Info
1956-0400	DieMaster Silver 15ml	US\$ 19.95	+Info
1956-0200	DieMaster red 15ml	US\$ 19.95	+Info
1956-0100	DieMaster Gray 15ml	US\$ 13.84	+Info
1956-0500	DieMaster Gold 15ml	US\$ 19.95	+Info
1956-0300	DieMaster Blue 15ml	US\$ 19.95	+Info

417581	Die:master Aqua Thinner 30ml	US\$ 15.96	+Info
417580	Die:master Aqua Dentine 15ml	US\$ 23.87	+Info

PICO-FIT | RENFERT

PICO-FIT - RENFERT

Function & Performance

High quality die spacer varnish with good contrast on the most popular plasters.

Advantages

- * Short drying period.
- * Flow property individually controlled using the thinner.
- * Builds an even layer thickness and a smooth surface.



US DENTAL DEPOT is an Authorized RENFERT dealer in US.

Code	Name	Price	+
023-1954-0600	RENFERT - Pico-Fit Silver 15ml-#1954-0600 #19540600 -(Espaciador de Modelos-Dow...	US\$ 23.04	+Info

TRU-FIT DIE SPACER - TAUB

TRU-FIT DIE SPACER - TAUB

A controlled technique for relieving dies before wax-up, by alternating coatings of gold and silver colored liquids over the die for uniform coverage of 25 microns space. Each coat is approximately 6 microns. Liquids are fast setting. Use on stone, amalgam, plated or epoxy dies. It is a high spot indicator. Tru-Fit's gold and silver colors can be used as a visual aid on models for patient education. Gold solution can be used to simulate crowns, inlays and bridgework. Silver solution can be used to simulate amalgam restorations, partial denture framework.

According to studies at a leading university, precision castings will not seat, be as much as 200 microns, unless relieved. To provide adequate space for cement and common internal discrepancies for ALL crowns and inlays, four coats of fast-drying Tru-Fit will provide much better marginal fit by allowing castings to seat more completely. Tru-Fit also greatly reduces the need for occlusal adjusting.



Refill

Code	Name	Price	+
469817	Tru-Fit Spacer Gold 4oz	US\$ 65.01	+Info
469818	Tru-Fit Spacer Silver 4oz	US\$ 65.01	+Info

469819	Tru-Fit Thinner 8oz	US\$ 20.37	+Info
469820	Empty Bottle W/Brush & Cap Ass.	US\$ 4.85	+Info

Tru-Fit

Tru-Fit

Code	Name	Price	+
138-138	Tru-Fit Kit Cement Spacer 1/2	US\$ 43.05	+Info
138-139	Tru-Fit Kit Cement Sp. 1oz G/S	US\$ 61.79	+Info
138-139-2OZ	Tru-Fit Kit Cement Spacer 2Oz	US\$ 59.09	+Info
138-139-4OZ	Tru-Fit Kit Cement Spacer 4Oz	US\$ 156.08	+Info

DIE SPACER - KEYSTONE

DIE SPACER - KEYSTONE

Features:

- * Fast drying
- * Heat resistant
- * 1oz bottles
- * Provides excellent uniform seating of the restoration



Code	Name	Price	+
034-1028525	KEYSTONE Die Spacer, Silver, 1oz Bottle - #1028525	US\$ 20.33	+Info
034-1028523	KEYSTONE BLUE Die Spacer - 1oz bottle #1028523 -	US\$ 20.33	+Info
034-1028524	KEYSTONE YELLOW Die Spacer - 1oz bottles #1028524 -	US\$ 20.33	+Info
034-1028527	KEYSTONE Thinner Die Spacer - 1oz bottle #1028527 -	US\$ 20.33	+Info

DIE SPACER - YETI

DIE SPACERS - YETI KEYSTONE

Description:

Yeti Die Spacers seals plaster and serves to stimulate the space needed for fixing cement between the crown and the die. This item is able to fulfill even the highest demands.

Features:

- * Available in Clear and Silver.
- * Comes in 18ml bottles.



Code	Name	Price	+
100718	YETI CLEAR Die Spacer 18ml bottle - #1860006 - #501-0018 (New Packaging 18ml)	US\$ 26.67	+Info
417079	YETI COLOR SPACER RED 7my 1x18ml #500-0001 replaced by 500-0018 (New Packaging ...	US\$ 26.67	+Info
417082	YETI DIE SPACER gold 13my 1x18ml #504-0018 (New Packaging 18ml) #504-0018	US\$ 26.67	+Info
186022	YETI Silver Die Spacer 18ml Bottle #ART 503-0018 #1860225 (New Packaging 18ml)	US\$ 26.67	+Info

DIE SPACER REFILL - ADS
DIE SPACER REFILL - ADS

Features:

- * Creating the exact spacing before wax-up is essential.
- * Produce controlled relief of dies & consistent film thickness.
- * Die Spacer provides excellent uniform seating of the restorations.
- * Heat resistant - lies flat.
- * Easy on your eyes.
- * One even coat provides approximately 12 microns relief.
- * Dries in minutes



Specify color

- * Blue, Clear, Cream, Gray, Green, Orange, Red, White, Yellow, Silver and Thinner

DIE SPACER AND THINNER KIT

ECONOMY 4 BOTTLE KIT

Code	Name	Price	+
117878	ADS Die Spacer economy - YELLOW kit (4 btl.) - #D750-6-YW	US\$ 34.99	+Info

DIE SPACER KITS STANDARD

Code	Name	Price	+
102404	ADS Die Spacer Kit - Standard - 1oz ADS Die Spacer + 1oz Thinner - ORANGE - #D75...	US\$ 14.80	+Info

FLASH DIE SPACER - ADS



KIT

Code	Name	Price	+
102491	ADS Flash Die Spacer Kit - Sliver 1oz + thinner #D751-1	US\$ 15.41	+Info

SILVER

Code	Name	Price	+
117925	ADS Die Spacer Flash Dry SILVER (spacer 1oz & thinner 1oz) kit - #D751-1-SR	US\$ 15.41	+Info

THINNER

Code	Name	Price	+
117933	ADS Die Spacer Flash Dry THINNER qt - #D751-12	US\$ 34.40	+Info

**TrueFIT Die Spacer - HARVEST DENT
AL**



Code	Name	Price	+
475331	HARVEST DENTAL - Harvest Surefit Die Spacer, Blue, 1/2oz - # 7151	US\$ 26.46	+Info
475332	HARVEST DENTAL - Harvest Surefit Die Spacer, Grey, 1/2oz - # 7152	US\$ 26.46	+Info
475333	HARVEST DENTAL - Harvest Surefit Die Spacer, Thinner, 1/2oz - # 7153	US\$ 26.46	+Info
475334	HARVEST DENTAL - Harvest Surefit Die Spacer, Dentine, 1/2oz - # 7156	US\$ 26.46	+Info
475335	HARVEST DENTAL - Harvest Surefit Die Spacer, Red, 1/2oz - # 7158	US\$ 26.46	+Info
475336	HARVEST DENTAL - Harvest Surefit Die Spacer, Gold, 1/2oz - # 7159	US\$ 26.46	+Info
475337	HARVEST DENTAL - Harvest Surefit Die Spacer, Silver, 1/2oz - # 7160	US\$ 26.46	+Info
475339	HARVEST DENTAL - Harvest Surefit Die Spacer Kit, Gold/Silver/Thinner, 1/2oz ea...	US\$ 54.00	+Info
475340	HARVEST DENTAL - Harvest Surefit Die Spacer Kit, Red/Blue/Thinner, 1/2oz each -...	US\$ 54.00	+Info

Die Spacer
Tooth Color Die Spacer

SENSATION SL TOOTH COLOR REFILLS - ADS
AMERICAN DENTAL SUPPLY - SENSATION SL TOOTH



- * Sensation SL Perfect Match Divex Liquid reduces the surface reaction during the pressing procedure on Lithium Disilicate (e.max)
- * Sensation SL Perfect Match Glaze Paste is available in 2 formulas: "MPC" for porcelain-to-metal and leucite-reinforced glass pressable ceramics and "EZ" for zirconia and lithium disilicate (e.max)
- * Perfect Match MPC Glaze Paste is a ready-to-use premixed glaze specifically designed for porcelain-to-metal and leucite-reinforced glass pressable ceramics (Cerpress SL)
- * Is a smooth glaze that allows the restoration to provide a beautiful glaze in one firing
- * The pre-mix provides the exact consistency for every application
- * Perfect Match EZ Glaze Paste is a ready-to-use premixed glaze specifically designed for zirconia and lithium disilicate all-ceramic restorations
- * Is a smooth glaze that allows the restoration to provide a beautiful glaze in one firing
- * The pre-mix provides the exact consistency for every application, is wear resistant, creating a protective coat and may be applied to polished zirconia
- * Perfect Match MPC and EZ Glaze Paste will not clump, get bumpy or leave streaks and provides a smooth glaze every time, allowing the ceramist to achieve a perfect glaze without losing the surface anatomy on the restoration

REFILL

Code	Name	Price	+
115022	ADS Die Aid Refill 0.5oz Natural shades- Sensation SL Perfect Match #PM101-12 PM...	US\$ 13.62	+Info

HATHO
POLISHING KIT

Thermoplastics
like *never* before



HATHO
POLISHING KIT

Thermoplastics
like *never* before



Packing Diamond D®
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory



Packing Diamond D®
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory



TruePrep Tooth Color Die Spacer - HA
RVEST DENTAL



Code	Name	Price	+
475345	HARVEST DENTAL - Harvest TruePrep Spacer Kit, (9 shades-7ml each) - # 7500	US\$ 266.18	+Info
475346	HARVEST DENTAL - Harvest TruePrep Spacer, HD1, 7ml - # 7501	US\$ 40.22	+Info
475347	HARVEST DENTAL - Harvest TruePrep Spacer, HD2, 7ml - # 7502	US\$ 40.22	+Info
475348	HARVEST DENTAL - Harvest TruePrep Spacer, HD3, 7ml - # 7503	US\$ 40.22	+Info
475349	HARVEST DENTAL - Harvest TruePrep Spacer, HD4, 7ml - # 7504	US\$ 40.22	+Info
475350	HARVEST DENTAL - Harvest TruePrep Spacer, HD5, 7ml - # 7505	US\$ 40.22	+Info
475351	HARVEST DENTAL - Harvest TruePrep Spacer, HD6, 7ml - # 7506	US\$ 40.22	+Info
475352	HARVEST DENTAL - Harvest TruePrep Spacer, HD7, 7ml - # 7507	US\$ 40.22	+Info
475353	HARVEST DENTAL - Harvest TruePrep Spacer, HD8, 7ml - # 7508	US\$ 40.22	+Info
475354	HARVEST DENTAL - Harvest TruePrep Spacer, HD9, 7ml - # 7509	US\$ 40.22	+Info

Catalog Index / Indice

Light Cure Die Spacer	1
Paint on Die Spacer	1
Tooth Color Die Spacer	6