

Wax pots & Dipping

Dipping pots



Dipping pots - BESQUAL

Dipping pots - BESQUAL

Code	Name	Price	+
101900	Digital Dipping Pot 110v	US\$ 111.25	+Info
469840	Single Well Wax Pot	US\$ 111.25	+Info
469841	Digital 4 Compartment Wax Pot	US\$ 101.24	+Info

Dipping pots - ECCO

Dipping pots - ECCO

Code	Name	Price	+
102470	Digital Dipping Pot 110v/220v	US\$ 89.25	+Info

Dipping pots - KEYSTONE

Dipping pots - KEYSTONE

Code	Name	Price	+
034-7000335	Wax Dipp Pot/Blue - 110Volt	US\$ 63.00	+Info

Dipping pots - YETI

Dipping pots - YETI

Code	Name	Price	+
417041	Dippy PRO digital 115v	US\$ 177.61	+Info

HOTTY - WAX COPINGS - RENFERT

HOTTY - WAX COPINGS - RENFERT

Function & Performance

The wax dipping pots hotty and hotty LED are the optimal units for producing dimensionally stable wax copings with a uniform thickness in only a few seconds.

Advantages

- * Quick, continuous heating phase.
- * Accurate results due to permanent compensation of minimal temperature fluctuations.
- * Uniform temperature distribution ensures consistent layer thickness of the wax copings.

Features

The membrane keypad of the hotty LED enables the temperature to be set accurately to within a degree. The



temperature status is permanently indicated by the LED display.
The hotty allows for temperature input and display via a potentiometer.

Details

- * Actual value indication in the dipping area of the dies.
- * Infinitely variable temperature setting from 60 - 110C (140 - 230F).
- * Ergonomically contoured housing provides optimal support for the fingers.
- * No contamination of the unit by wax residue thanks to the specially designed wax runback pot.
- * The integrated lid provides protection against dirt and dust.

HOTTY LED

Code	Name	Price	+
023-1461-1000	RENFERT Hotty Led 110 Volts Each- #14611000 #14611000 Calentador Cera Electric W...	US\$ 299.49	+Info
023-1461-0000	RENFERT Hotty Led 220 Volts Each- #1461-0000 #14610000 Calentador Cera Electric ...	US\$ 347.68	+Info

HOTTY

Code	Name	Price	+
023-1459-1000	RENFERT Hotty 110 Volts-Each- #1459-1000 #14601000 Calentador Cera Electric Waxer	US\$ 240.75	+Info
023-1459-0000	RENFERT Hotty 220 Volts Each- #1460-0000 #14600000	US\$ 277.28	+Info

Wax pots & Dipping

Wax pot



MULTI DIGITAL WAXER UNITS - BESQUAL

MULTI DIGITAL WAXER UNITS - BESQUAL

This Multi Wax Pot is designed to melt a variety of wax with a liquid condition.

This unit can heat some amount of 3 different type of waxes quickly and keep its liquid state. So it make the waxing works easy and convenient



Code	Name	Price	+
118734	BesQual S900 Multi Digital Waxer unit 3 in 1 - (No-Flame + Wax Dipping Pot + El...	US\$ 306.25	+Info
115634	BesQual S800 Multi Digital Wax Pot with LCD Display #115634 #723-800	US\$ 148.50	+Info

Wax pot - ECCO


Wax pot - ECCO

Code	Name	Price	+
116115	Waxing Unit Digital 3 box 220v	US\$ 106.73	+Info
205-47715	Waxing Unit Digital 3 box 110v	US\$ 106.73	+Info

THE HEAT ZONE - WHIPMIX

THE HEAT ZONE - WHIPMIX

The Heat Zone heats wax modeling tools without flame using an electromagnetic induction process and is sensor activated when a tool is inserted into the "zone". It heats the tool safely and extremely fast, replacing the traditional Bunsen Burner. The induction form of heating ensures complete safety for non-metallic items such as fingers and clothing eliminating the potential for burns or fire.



Description:

- * Heats wax tools fast without flame.
- * No carbon deposits on instruments or on wax ups.
- * Energy efficient as power is used only at the moment of heating.
- * No gas lines required, making unit mobile and operational anywhere an electric outlet is available.
- * Indicator light glows green for ready mode and red for when in use.
- * Safety light blinks when instrument is left in too long
- * Replaceable inserts keeps unit clean of wax residue (5 included)
- * Two Year warranty
- * The Heat Zone provides a safe, clean, convenient way to heat wax tools without compromising productivity

Code	Name	Price	+
120299	WHIP-MIX Heat Zone 110v - #10147	US\$ 455.63	+Info
120300	WHIP-MIX Heat Zone 220vol - #10247	US\$ 33.97	+Info
120298	WHIP-MIX Replacement Liners Heat Zone Silicone Pk/2 - #10149	US\$ 13.35	+Info

VARIO E-WAXER - RENFERT

VARIO E-WAXER - RENFERT

Function & Performance

The Vario E allows you to keep your sculpting wax at the required working temperature. Protect your wax against too high temperatures, your sculpting tips against mechanical loading and your wax pattern against high stresses.

Advantages

- * 30% saving in time due to the preheating function.
- * Multi-functionality provided by three wax pots with separate control.
- * Individual setting of the wax consistency (from creamy to fluid) using infinitely variable temperature setting of 45 - 110C (113 - 230F).

Details

- * The permanent temperature control of the actual value input is indicated by the LED display.



- * Also suitable for processing light-curing waxes.
- * Ergonomic working height ensures relaxed working.
- * The deep wax pot is ideal for dipping copings.

Code	Name	Price	+
112256	RENFERT Vario E-Waxer 110 Volts-#1452-1000 #14521000 -The Vario E heats modellin...	US\$ 460.99	+Info
112215	RENFERT Vario E-Waxer 220 Volts- #1452-0000 #14520000 -The Vario E heats modelli...	US\$ 473.35	+Info

Wax pot - RAY FOSTER

Wax pot - RAY FOSTER

Code	Name	Price	+
115235	Electric Wax Heatr w/Thrm 110v	US\$ 197.10	+Info
470921	WH40 Standard Wax Heater	US\$ 106.65	+Info

WAXPROFI - RENFERT

WAXPROFI - RENFERT

Function & Performance

Your individual skill is required to give a denture freshness and vitality. Our Waxprofi is required to implement this efficiently. Enjoy the pleasure of how easy setting-up has become.

Advantages

- * 50% time saving due to the preheating function.
- * Individual setting of the wax consistency (from creamy to fluid) by infinitely variable temperature setting from 40 - 110C (104 - 230F).
- * Work without continually refilling due to a large volume pot (120 ml / 4.08 fl.oz.).

Details

- * Ergonomic height ensures easy uptake of wax portions.
- * Clean handling due to functional pot design.
- * Integrated lid protects the wax.
- * Suitable for the use of light-curing waxes.



Code	Name	Price	+
23-14401000	RENFERT Waxprofi 110 Volts-178 x 54 x 120mm (7.01 x 2.13 x 4.73 inch)-110 V-50 ...	US\$ 365.84	+Info
23-14400000	RENFERT Waxprofi 220 Volts-178 x 54 x 120mm (7.01 x 2.13 x 4.73 inch)-220 V/50 ...	US\$ 399.59	+Info

Catalog Index / Indice

Dipping pots	1
Wax pot	2