

Duplicating Material

All purpose Duplicating Mat.



All purpose Duplicating Mat. - GC AMERICA

All purpose Duplicating Mat. - GC AMERICA

Code	Name	Price	+
421361	Dupli-Coe-Loid 16Lb/2gal	US\$ 117.18	+Info
421720	Dupli-Coe-Loid 8 Lb. Jar (1 Gallon)	US\$ 76.26	+Info

All purpose Duplicating Mat. - NOBILIUM

All purpose Duplicating Mat. - NOBILIUM

Code	Name	Price	+
N-1009	Nobilflex Dup Coll. 3.5 gal	US\$ 125.00	+Info

ALL PURPOSE DUPLICATING MATERIAL - CFS

DUPLICATING MATERIAL - ALL PURPOSE GEL - CFS

CFS Duplicating Gel is one of our more highly regarded secrets.
We take pride in our product



Code	Name	Price	+
235106	CFS Duplicating Material - All purpose - Gel - 1 gal - a highly regarded secret...	US\$ 53.00	+Info
294106	CFS Duplicating Material - All purpose - Gel - 5 gal - a highly regarded secret...	US\$ 154.85	+Info

ALL PURPOSE DUPLICATION MATERIAL - BEGO

ALL-PURPOSE DUPLICATING MATERIAL - BEGO

This All-Purpose Duplicating Material can be used with all stone and investment models also for pour acrylic dentures.

Features:

- * A quality, reversible duplicating material for use with all investment types
- * Ideal for use with all duplicating conditioners
- * All-Purpose
- * No wash required



CastoGel
CastoGel

Code	Name	Price	+
106940	Castogel Rvrsibl Dupl Gel 6kg	US\$ 111.55	+Info
Wirogel			
Wirogel			
Code	Name	Price	+
419892	WiroGel C- 10 kg	US\$ 148.09	+Info
WIROGEL M			
WIROGEL M			
Code	Name	Price	+
10-54352	WiroGel M Dupl Gel 10kg- #54354	US\$ 112.07	+Info

ALL PURPOSE DUPLICATION MATERIAL - KEYSTONE

ALL-PURPOSE DUPLICATING MATERIAL - KEYTONE

This All-Purpose Duplicating Material can be used with all stone and investment models also for pour acrylic dentures.

Features:

- * A quality, reversible duplicating material for use with all investment types
- * Ideal for use with all duplicating conditioners
- * All-Purpose
- * No wash required



Code	Name	Price	+
034-1310300	KEYSTONE Duplicating Material - All Purpose - 1gal - Dual - Gypsum and Investmen...	US\$ 72.17	+Info
034-1310330	KEYSTONE All-Purpose Duplicating Material, PINK. 5 Gallon. Can be used with all...	US\$ 267.05	+Info

ALL PURPOSE DUPLICATION MATERIAL - NOBILIU M



Instaloid

Code	Name	Price	+
N-7220	NOBILIUM Instaloid Concentrated Duplicating Material - All Purpose Can mix smal...	US\$ 48.00	+Info

GI Mask gingival - COLTENE



Code	Name	Price	+
426608	GI Mask Automix 2x50ml	US\$ 88.33	+Info
426708	GI Mask Automix New Formula - Silic	US\$ 136.76	+Info
426849	GI Mask Silicone-based Gingival Mas	US\$ 53.14	+Info
426928	GI Mask Silicone-based Gingival Mas	US\$ 132.12	+Info
427052	GI Mask Silicone-based Gingival Mas	US\$ 25.94	+Info
436272	GI Mask Silicone-based Gingival Mas	US\$ 234.36	+Info

GINGIFAST Gingival Mask - ZHERMACK

GINGIFAST - ZHERMACK

Laboratory A-silicone for gingival reproduction on models.Rigid type

Use:

* Gingival reproduction for implant prosthetic models. Indicated for direct technique

* Uses Pistol 1:1 or 4:1

Characteristics:

* Addition silicone formulation in a Cartridge

Advantages:

* Easy to Trim: permits accurate margin definition

* Indicated for implant prosthetics

* Precision and dimensional stability



GINGIFAST CAD

ELASTIC

Code	Name	Price	+
118463	ZHERMACK Gingifast CAD Elastic - Standard Pack - Z#C203227 (2 x 50 cartridges + ...	US\$ 66.67	+Info

RIGID

Code	Name	Price	+
118465	ZHERMACK Gingifast CAD rigid - Standard pack - Z#C203232 (2 x 50 cartridges + 12...	US\$ 66.67	+Info

ELASTIC

Code	Name	Price	+
111389	ZHERMACK Gingifast Elastic - Silicone Removable Gingival Mask: 2 - 50 mL Cartrid...	US\$ 58.95	+Info

RIGID

Code	Name	Price	+
111388	ZHERMACK Gingifast Rigid - Silicone Removable Gingival Mask: 2 - 50 mL Cartridge...	US\$ 58.95	+Info

SEPARATOR

Code	Name	Price	+
106929	zhermack Gingifast Separator - 10ml bottle (Z#C400887) - Z#C400888 For gingiv...	US\$ 10.96	+Info

Gingival Mask - COLTENE

Gingival Mask - COLTENE

Code	Name	Price	+
436275	GI Mask Silicone Trimmer	US\$ 32.44	+Info
436280	GI Mask Silicone-based Gingival Mas	US\$ 0.00	+Info

Mask Silicone - HENRY SCHEIN

Mask Silicone - HENRY SCHEIN



Code	Name	Price	+
420278	Gingival Mask Silicone Std Kit	US\$ 38.67	+Info

Duplicating Material

Lab Putty



Lab Putty - GC AMERICA

Lab Putty - GC AMERICA



Code	Name	Price	+
421352	GCLT Laboratory Putty B&C	US\$ 286.60	+Info

ZETALABOR - ZHERMACK

ZETALABOR - ZHERMACK

C-Silicone Laboratory Putty

Use:

- * Duplication of denture moulds
- * Control matrix for the design of metal structures
- * Impressions for copy denture technique
- * Matrix for use in making temporary crowns/bridges
- * Matrix for the reproduction of gingival tissue (e.g. Zhermack's gingifast)
- * Insulation of teeth/gypsum in muffle

Chacarcteristics:

- * To be mixed with gel or liquid catalyst (not included)
- * Final hardness (85 Shore A)
- * Non-pigmented

Adventages:

- * High initial flowability
- * Good final hardness
- * High detail reproduction
- * No deformation in flask. Good breaking strength
- * Heat resistance 140°C
- * High elastic resistance



CATALYST



Code	Name	Price	+
470855	ZHERMACK Indurent LAB, NEW catalyst for Zetalabor and Titanium condensation sili...	US\$ 10.80	+Info

ECONOMY

Code	Name	Price	+
114170	ZHERMACK Zetalabor Silicone Lab Putty, Regular Type, GRAY, 5 Kg - #C400811	US\$ 86.87	+Info
114172	ZHERMACK Zetalabor Economy Pack (1x 10 kg canister (22 lbs) / 1x spoon) High Qua...	US\$ 149.20	+Info
118525	ZHERMACK Zetalabor Eco Pack (25 kg) - Z#C400812	US\$ 307.48	+Info

STANDARD

Code	Name	Price	+
114171	ZHERMACK Zetalabor Standard Pack (1x 2.6 kg tub (5.7 lbs) High Quality Putty 85...	US\$ 55.55	+Info

Lab Putty - 3M ESPE

Lab Putty - 3M ESPE

Code	Name	Price	+
428221	Express Putty Spoon, Package of 2.	US\$ 26.32	+Info

Lab Putty - COLTENE WHALADENT

Lab Putty - COLTENE WHALADENT

Code	Name	Price	+
426467	Lab Putty Polysil Modelling 1.35kg	US\$ 77.22	+Info
426794	Coltene Lab-Putty Hard Std Set	US\$ 109.13	+Info
426801	Lab Putty Polysilox Hard 40ml	US\$ 40.24	+Info
430793	Coltene/Whaledent Lab Putty Polysil	US\$ 366.54	+Info
431665	Coltene/Whaledent Lab-Putty Hard Ba	US\$ 196.23	+Info

Lab Putty - PAC-DENT

Lab Putty - PAC-DENT

Code	Name	Price	+
449027	Pac-Dent Vinyl polysiloxane Putty,	US\$ 61.97	+Info

Duplicating Material

Silicone 4 duplicating (Material)



DUPLICATION SILICONE ELITE - ZHERMACK

ZHERMACK - DUPLICATION SILICONE ELITE

A-Silicone for laboratory model duplication. Ideal for duplicating stumps using the all ceramic technique and models with accentuated undercuts and investments. Use of a plastic duplication flask is recommended.

Advantages

- * Rapidity: a duplicate can be obtained in just 10 minutes (23)
- * Mixing facilitated by the 1:1 base to catalyst ratio
- * High fluidity
- * Absolute precision for faithful reproduction of detail
- * Dimensional stability over time and non-deformability, allowing a number of duplications to be made
- * Compatible with all plasters, polyurethane resins, phosphate investments and acrylic resins



Standard Pack

Fast

Code	Name	Price	+
115864	Elite Double 22 Fast 1+1 Kg #C400834 - Z#C400834 Duplicating Silicone	US\$ 119.49	+Info

Regular

Code	Name	Price	+
118526	ZHERMACK Elite Double 8 Standard Pack (1 kg) - Z#C400830 Duplicating Silicone	US\$ 119.49	+Info
118527	ZHERMACK Elite Double 22 Standard Pack (1 kg) - Z#C400832 Duplicating Silicone	US\$ 119.49	+Info
118528	ZHERMACK Elite Double 32 Standard Pack (1 kg) - Z#C400833 Duplicating Silicone	US\$ 119.49	+Info

Economy Pack

This is the Economy Value Pack from Zhermack.

Is a 5 kg package

Code	Name	Price	+
118534	ZHERMACK Elite Double 16 FAST Eco Pack (5 kg) - Z#C400847 Duplicating Silicone	US\$ 516.36	+Info
118530	ZHERMACK Elite Double 22 Eco Pack (5 kg) - Z#C400840 Duplicating Silicone	US\$ 516.36	+Info
118532	ZHERMACK Elite Double 22 FAST Eco Pack (5 kg) - Z#C400842 Duplicating Silicone	US\$ 516.36	+Info

DUPLICATION SILICONE REDU-IT - ASD

DUPLICATION SILICONE REDUT-IT - ADS

Premium model duplicating silicone. Great dimensional stability & accuracy. Won't shrink on standing. Unlimited duplication can be made from master model at any time. Easy release from any model material.



Code	Name	Price	+
118317	ADS redu It GREEN Duplicating Silicone Material 2 pint kit - #R774-5	US\$ 119.89	+Info
118318	ADS redu It GREEN Duplicating Silicone Material - 2qt kit - Part A & Part B - ...	US\$ 171.77	+Info
118319	ADS redu It GREEN Duplicating Silicone Material 2 gallon kit - #R774-8	US\$ 483.78	+Info

PLATINIUM 95 & 85 - ZHERMACK

PLATINIUM - ZHERMACK

High precision laboratory A-Silicone putty material

Use:

- * Duplication of denture moulds
- * Ceramics technique printed on verticulator
- * Die in full dentures
- * Control matrix for the design of metal structures
- * Denture repair
- * Bite recording key for positioning in articulator



- * Matrix for use in making temporary crowns/bridges
- * Matrix support for the working of metal free aesthetic materials
- * Insulation of teeth/plaster in flask

Code	Name	Price	+
118521	ZHERMACK Platinum 85 Eco Pack (5 kg) - Z#C400723 (Silicone material for duplica...	US\$ 324.54	+Info
118519	ZHERMACK Platinum 95 Eco Pack (5 kg) - Z#C400710 (Silicone material for duplica...	US\$ 381.23	+Info
118522	ZHERMACK Platinum 85 Standard Pack (1 kg) - Z#C400725 (Silicone material for du...	US\$ 114.88	+Info
118518	ZHERMACK Platinum 95 Standard Pack (1 kg) - Z#C400700 (Silicone material for du...	US\$ 114.88	+Info

Silicone 4 duplicating (Material) - BEGO

Silicone 4 duplicating (Material) - BEGO

Code	Name	Price	+
471456	Wirosil Duplicating Silicone, 2	US\$ 381.54	+Info

Silicone 4 duplicating (Material) - BUFFALO

Silicone 4 duplicating (Material) - BUFFALO

Code	Name	Price	+
424524	Flexistone Plus Rigid Insul Kit	US\$ 209.22	+Info
424759	Flexistone for Insulating Acrylic o	US\$ 91.95	+Info

Silicone 4 duplicating (Material) - ZHERMACK

Silicone 4 duplicating (Material) - ZHERMACK

Platinum

Platinum

Code	Name	Price	+
420817	Platinum 85 Eco Pack 5kg	US\$ 381.23	+Info

Titanium

Titanium

Code	Name	Price	+
118513	Titanium Standard Pack (2.6 k	US\$ 54.07	+Info
118514	Titanium Eco Pack (5 kg)	US\$ 86.27	+Info

Catalog Index / Indice

All purpose Duplicating Mat.	1
Gingival Mask	2
Lab Putty	4
Silicone 4 duplicating (Material)	6