

## Flasks & Compresses

### Compresses




### BIG BOY FLASK PRESS - HANDLER

#### BIG BOY FLASK PRESS - HANDLER

A Heavy Duty Flask Press manufactured of malleable/ductile cast iron. With alloy steel shaft, grey iron handle and malleable/ductile iron pressure plate. Holes in frame allow for easy bench mounting, or 38A Clamps may be utilized for easy bench mounting and/or removal.

Both 38 and 38B Big Boy Flask Presses are designed for the busy laboratory.



Code	Name	Price	+
103359	#38B HANDLER Big Boy Press #38B Heavy Duty Two Flask Press manufactured of malle...	US\$ 288.33	<a href="#">+Info</a>
103402	#38 HANDLER Big Boy Flask Press - 10in x 4.5in x 8.5in - For single Flask - Cons...	US\$ 265.18	<a href="#">+Info</a>

### Compresses - BESQUAL

Compresses - BESQUAL

Code	Name	Price	+
100695	Compress #2 Spring Bronze ea	US\$ 81.00	<a href="#">+Info</a>

### Compresses - HANDLER

Compresses - HANDLER

Code	Name	Price	+
417399	Clamp for Press for 2 flask	US\$ 75.60	<a href="#">+Info</a>

### FLASK COMPRESS - HANDLER

FLASK COMPRESS - HANDLER

Our most popular denture compress is built to hold one, two or three conventional flasks for simultaneous curing. Sturdily built of bronze alloy, its rustproof cast body will provide years of service. The stainless steel springs give equalized compensatory pressure during curing operation. A suitable Allen wrench, wood capped holding handle and brass lifting handle are supplied with each compress. Flasks not included.



Code	Name	Price	+
103420	#61B HANDLER Flask Compress - 4.12in x 2.5in x 7.5in - Two Flask	US\$ 198.88	<a href="#">+Info</a>

### HANAU - FLASK COMPRESS - WHIP MIX

#### HANAU - FLASK COMPRESS - WHIP MIX

**Features:**

Durable, one-piece casting of special high tensile manganese bronze, notable for its strength. Special wrench and holding handle are supplied. Convenient lifting loop is attached, permitting easy compress placement or removal from hot curing baths.



#### ACCESSORIES

Code	Name	Price	+
120088	WHIP MIX 44-1017 Wrench 9/16 In - #001015-000 ( Accessories for Flask Compress)	US\$ 37.63	<a href="#">+Info</a>

#### Flasks & Compresses

### Crown & Bridge Flasks



### Crown & Bridge Flasks - DENTAL FARM

Crown & Bridge Flasks - DENTAL FARM

Code	Name	Price	+
474542	PHOTO FLASK	US\$ 420.47	<a href="#">+Info</a>

#### Flasks & Compresses

### Denture Flasks



### DENTURE FLASK BRONZE - BESQUAL

DENTURE FLASK BRONZE - BESQUAL



- \* Excellent quality
- \* Very strong and durable
- \* Polished insides and easy cleaning

Code	Name	Price	+
100691	BESQUAL Bronze Denture Flask Upper #611-100	US\$ 81.00	<a href="#">+Info</a>
100692	BESQUAL Lower Denture Flask, Ejector Type. Made of high quality bronze alloy - ...	US\$ 81.00	<a href="#">+Info</a>

### Denture Flasks - HANDLER

Denture Flasks - HANDLER

Code	Name	Price	+
103366	Lower Flask 4x3.75x2.25in	US\$ 122.65	<a href="#">+Info</a>
103367	Upper Flask 4x3.75x2.25in	US\$ 122.65	<a href="#">+Info</a>

103369	De Flasking Disc 3x2.75x0.37in	US\$ 32.83	<a href="#">+Info</a>
--------	--------------------------------	------------	-----------------------

### Denture Flasks - KEYSTONE

Denture Flasks - KEYSTONE

Code	Name	Price	+
034-1502400	Keystone Reline Jig	US\$ 29.64	<a href="#">+Info</a>

### Denture Flasks - LANG

Denture Flasks - LANG

Code	Name	Price	+
473742	ReFlex Reline Jig 3 Post Jig	US\$ 76.98	<a href="#">+Info</a>

### HANAU - VARSITY FLASK EJETOR TYPE - W HIP MIX

HANAU - VARSITY FLASK EJETOR TYPE - WHIP MIX

- \* Made of forged brass, interior and exterior are both finished
- \* Tapered, stainless steel guide pins securely and accurately align flask sections

**Dimensions:** 35/8"L x 4"W x 21/3"H



### FLASKS

Code	Name	Price	+
120084	HANAU 58-0 Lower Flask - Laboratory Accesories - #001702-000 ( Accessories for L...	US\$ 239.22	<a href="#">+Info</a>
120086	WHIP MIX - HANAU 57-0 Upper Flask - Laboratory Accesories - #001602-000 ( Access...	US\$ 229.65	<a href="#">+Info</a>

### ACCESSORIES

Code	Name	Price	+
120082	WHIP MIX 59-1648 Pry Bar - #001814-000 ( Accessories for Varsity Flask)	US\$ 136.82	<a href="#">+Info</a>
120087	WHIP MIX 57-1587 Disc- Flask - #001613-000 ( Accessories for Varsity Flask)	US\$ 30.79	<a href="#">+Info</a>

### Flasks & Compresses

## Duplicating Flasks

### Duplicating Flasks - BEGO

Duplicating Flasks - BEGO

Code	Name	Price	+
------	------	-------	---

10-52090	Kombi Dupl Flask 2 size	US\$ 76.93	<a href="#">+Info</a>
----------	-------------------------	------------	-----------------------

### DUPLICATING FLASK - KEYSTONE

#### DUPLICATING FLASK - KEYSTONE

These Duplicating Flasks are used in producing full or partial dentures. High quality and long lasting. Available in raised bottom and flat bottom



Code	Name	Price	+
034-1350080	KEYSTONE 568 Flask - Duplicating - Flat Bottom #1350080 #DF78	US\$ 40.13	<a href="#">+Info</a>
034-1350090	KEYSTONE 568A Flask - Duplicating Raised Bottom #1350090 - D78A	US\$ 40.13	<a href="#">+Info</a>

### DUPLICATING FLASK - RAY FOSTER

#### RAY FOSTER Duplicating Flasks

##### Features:

- \* Ray Foster Duplicating Flasks are made of durable die cast aluminum and finished with a smooth baked-on powder coat.
- \* The finish facilitates the easy removal of the duplicating material.
- \* They are available with either a flat or raised base.



Code	Name	Price	+
101587	RAY FOSTER - Duplicating Flask - Raised Base - Base: 4in ( 101mm ) x 3.12in ( 96...	US\$ 35.88	<a href="#">+Info</a>
101586	RAY FOSTER - Duplicating Flask - Flat Base - Base: 4in ( 101mm ) x 3.12in ( 96mm...	US\$ 36.40	<a href="#">+Info</a>

### Duplicating Flasks - MESTRA

#### Duplicating Flasks - MESTRA

Code	Name	Price	+
417395	Plastic Duplicating Flasks	US\$ 37.80	<a href="#">+Info</a>
417396	Duplicating Tray for Silicones	US\$ 30.94	<a href="#">+Info</a>

## Flasks & Compresses

### Flasks Accessories

### Flasks Accessories - HANDLER

#### Flasks Accessories - HANDLER



Code	Name	Price	+
103358	#38A Bench Clamps	US\$ 46.26	<a href="#">+Info</a>
103364	#44D Bronze K.O. 2.43x0.18in	US\$ 31.20	<a href="#">+Info</a>
103418	Spring Filler 3-1/2x2-1/2in	US\$ 82.73	<a href="#">+Info</a>

### Flasks Accessories - WHIP MIX

Flasks Accessories - WHIP MIX

Code	Name	Price	+
120080	Flask Band 3 Pk	US\$ 40.90	<a href="#">+Info</a>

### Flasks & Compresses

## Flasks for Flexibles thermopla



### Flasks for Flexibles thermopla - CFS

Flasks for Flexibles thermopla - CFS

Code	Name	Price	+
202106	Alum Flask for Nylon Inject ea	US\$ 176.25	<a href="#">+Info</a>

### Flasks for Flexibles thermopla - FLEXITE

Flasks for Flexibles thermopla - FLEXITE

Code	Name	Price	+
101824	Alum Flask for Nylon Inject ea	US\$ 185.00	<a href="#">+Info</a>

### Flasks for Flexibles thermopla - GFM

Flasks for Flexibles thermopla - GFM

Code	Name	Price	+
122292	Alum Flask for Nylon Inject ea	US\$ 180.00	<a href="#">+Info</a>

### Flasks for Flexibles thermopla - SABLEX

Flasks for Flexibles thermopla - SABLEX

Code	Name	Price	+
421047	NEW Alum Flask for Nylon Inject ea	US\$ 80.00	<a href="#">+Info</a>

### Flasks for Flexibles thermopla - SNOW ROCK

Flasks for Flexibles thermopla - SNOW ROCK

Code	Name	Price	+
474223	Jet Flask ONLY Extra Standard Flask	US\$ 280.01	<a href="#">+Info</a>

## Flasks for Flexibles thermopla - VALPLAST

Flasks for Flexibles thermopla - VALPLAST

Code	Name	Price	+
829-20582	Bolts for flasks	US\$ 5.63	<a href="#">+Info</a>
829-205FL03	Vpst Aluminum Flask	US\$ 276.27	<a href="#">+Info</a>

### Flasks & Compresses

## Flasks for Microwave



## Flasks for Microwave - BUFFALO

Flasks for Microwave - BUFFALO

Code	Name	Price	+
424494	Microwavable Flask,	US\$ 154.44	<a href="#">+Info</a>

## Flasks for Microwave - GC AMERICA

Flasks for Microwave - GC AMERICA

Code	Name	Price	+
421622	Nature-Cryl MC FRP Microwave Flask	US\$ 425.17	<a href="#">+Info</a>
422013	Nature-Cryl MC FRP Flask Guard	US\$ 93.10	<a href="#">+Info</a>
422889	Nature-Cryl MC FRP Microwave Replac	US\$ 130.68	<a href="#">+Info</a>

## TECNOFLASK KIT MUFFLE - KEYSTON E

### TECNOFLASK KIT MUFFLE - KEYSTONE

Keystone Industries could not be more excited for a product that allows Diamond D and other acrylics to be processed in much quicker time and easier fashion inside of a microwave.

The TF-100 muffle will give your laboratory a cleaner environment and make the process more enjoyable for the production of your work. Included with the kit and available separately is the PL-100 protection plate, which extends the muffles service life with extra needed protection.



Code	Name	Price	+
116597	KEYSTONE Tecnoflask it Muffle + Board Microweave Flasks - microweavable mfg #100...	US\$ 261.32	<a href="#">+Info</a>
116598	KEYSTONE Tecnoflask Muffle only mfg #1009200 -	US\$ 197.83	<a href="#">+Info</a>

### Flasks & Compresses

## Pour Flasks



### Pour Flasks - GC AMERICA

Pour Flasks - GC AMERICA

Code	Name	Price	+
123015	POUR-N-CURE FLASK	US\$ 114.95	<a href="#">+Info</a>

### Pour Flasks - MESTRA

Pour Flasks - MESTRA

Code	Name	Price	+
417397	Pouring Flask KIT	US\$ 110.27	<a href="#">+Info</a>

## Catalog Index / Indice

Compresses .....	1
Crown & Bridge Flasks .....	2
Denture Flasks .....	2
Duplicating Flasks .....	3
Flasks Accessories .....	4
Flasks for Flexibles thermopla .....	5
Flasks for Microwave .....	6
Pour Flasks .....	6