

Gypsums & Die Stone Acc

Lab Plasters TYPE II



GYPSUMS LABORATORY PLASTER TYPE 2 - WHIP MIX

GYPSUMS LABORATORY PLASTER TYPE 2 - WHIP MIX

An economical white plaster formulated for general purpose laboratory use. Improved to provide a creamier mix consistency. Available in both regular set and fast set.



ISO Type 2

Code	Name	Price	+
119734	WHIP MIX Laboratory Plaster Regular Set White 33lb/15kg - #330086 TYPE-II (TYPE...	US\$ 34.54	+Info
119736	WHIP MIX Laboratory Plaster Fast set 33lb/15kg - #330085 TYPE-II (TYPE-2)	US\$ 35.09	+Info

LAB PLASTER - ECCO

LAB PLASTER - ECCO

Features:

- * Available in 25lb and 50lb
- * White
- * A fine low expansion plaster
- * Recommended for models- casts and flasking

Physical Properties:

- * 47cc
- * 100 grams
- * 6 - 8 minutes initial set
- * With a compression strength of up to 2,300 lbs.
- * PSI/15.86 MPA after 1 hour - up to 3,500 lbs.
- * PSI/24.13 MPA dry



FAST LAB PLASTER

Code	Name	Price	+
112421	ECCO - PLASTER FAST SET 50lb - #900-4035 low expansion TYPE-II (TYPE-2)	US\$ 42.43	+Info

REGULAR LAB PLASTER

Code	Name	Price	+
112420	ECCO - PLASTER REGULAR SET 50lb - #900-4037 TYPE-II (TYPE-2) A fine low expan...	US\$ 41.06	+Info
112419	ECCO - PLASTER REGULAR SET 25lb -TYPE-II - A fine low expansion plaster recommen...	US\$ 25.57	+Info

Gypsums & Die Stone Acc

Lab Stones TYPE III



LAB STONE TYPE-3 - ECCO

LAB STONE- ECCO

Features:

- * Available in 25lb and 50lb
- * Strong Casts
- * Comparable with Quickstone

Physical Properties:

- * 30 cc
- * 100 grams
- * 8 - 10 minutes initial set
- * With a compression strength of up to 4-000 lbs.
- * PSI/27.58 MPA after 1 hour - up to 8-500 lbs.
- * PSI/58.61 MPA dry.



Code	Name	Price	+
112416	ECCO - LABORATORY STONE 25lb Buff - LABSTONE TYPE-III (TYPE-3) FAST SET	US\$ 34.09	+Info
112415	ECCO - LABORATORY STONE 47.5lb Buff - LABSTONE TYPE-III (TYPE-3) #900-4029 FA...	US\$ 51.52	+Info
474372	ECCO - LABORATORY STONE 50lb BLUE- TYPE-III - FAST SET .	US\$ 52.36	+Info

ELITE MODEL - ZHERMACK

ELITE MODEL FAST - ZHERMACK

Type 3 fast dental stone for antagonists and models in general

Use:

- * Fast production of models in general, antagonists and flasking of removable prostheses

Characteristics:

- * Thixotropic
- * Low expansi
- * High mechanical characteristics
- * Universal use
- * Available in three colours: Sky Blue, Light Cream, White

Advantages:

- * Fast casting of the model in only 15:00 min maintaining an adequate working time
- * Thixotropic for simplified model production
- * Extreme versatility due to its superior characteristics compared to other stones of the same class
- * Optimal compressive strength



55 LIBRAS CARTON ECO

Code	Name	Price	+
118566	ZHERMACK Elite Model 55lb Carton Eco Pack Steel BLUE - Z#C410221 TYPE-III (TYPE...	US\$ 90.82	+Info
118567	ZHERMACK Elite Model 55lb Carton Eco Pack WHITE - Z#C410222 TYPE-III (TYPE-3)	US\$ 90.82	+Info

6.6 LIBRAS BAG

Code	Name	Price	+
118550	ZHERMACK Elite Model 6.6lb Bag Steel BLUE - Z#C410070 TYPE-III (TYPE-3)	US\$ 15.98	+Info

MODEL FAST

55 LIBRAS CARTON ECO

Code	Name	Price	+
118569	ZHERMACK Elite Model FAST 55lb Carton Eco Pack Light Cream - Z#C410226 TYPE-III ...	US\$ 79.23	+Info
118570	ZHERMACK Elite Model FAST 55lb Carton Eco Pack WHITE - Z#C410228 TYPE-III (TYPE...	US\$ 79.23	+Info

6.6 LIBRAS BAG

Code	Name	Price	+
118547	ZHERMACK Elite Model FAST 6.6lb Bag WHITE - Z#C410063 TYPE-III (TYPE-3)	US\$ 17.00	+Info
118548	ZHERMACK Elite Model FAST 6.6lb Bag Sky BLUE - Z#C410065 TYPE-III (TYPE-3)	US\$ 16.97	+Info

GYPSUMS MICROSTONE TYPE III - WHIP MIX

GYPSUMS MICROSTONE TYPE III - WHIP MIX

MICROSTONE's fine particle size and dense surface make it our finest general purpose laboratory stone. It is especially suited for processing acrylic dentures as it accepts tinfoil substitutes readily.

ISO Type 3



GOLDEN

Code	Name	Price	+
119704	WHIP MIX Microstone Golden 33lb/15kg - #330054 TYPE-III (TYPE-3)	US\$ 57.67	+Info

Lab Stones TYPE III - ECCO

Lab Stones TYPE III - ECCO

Code	Name	Price	+
474371	LAB STONE 25lb White	US\$ 48.33	+Info

GYPSUMS QUICKSTONE TYPE III - WHIP MIX

GYPSUMS QUICKSTONE TYPE III - WHIP MIX

QUICKSTONE is an economical dental stone with consistent handling properties. Its smooth creamy mix makes it ideal for a variety of lab procedures.

ISO Type 3



BUFF

Code	Name	Price	+
119706	WHIP MIX Quickstone Buff 33lb - #330062 TYPE-III (TYPE-3)	US\$ 42.30	+Info

Lab Stones TYPE III - GC AMERICA

Lab Stones TYPE III - GC AMERICA

Code	Name	Price	+
421607	Coe-Cal Stone Buff Hard 50lb	US\$ 60.38	+Info
421648	Coe-Cal Stone Buff Hard 25lb	US\$ 45.18	+Info
421713	Coecal Stone White TYPE-III 50lb	US\$ 60.38	+Info

Lab Stones TYPE III - HOUSE BRAND

Lab Stones TYPE III - HOUSE BRAND

Code	Name	Price	+
469264	Orthodontic Stone White 25lb	US\$ 55.59	+Info

Lab Stones TYPE III - WHIP MIX

Lab Stones TYPE III - WHIP MIX

Code	Name	Price	+
119700	Rapidflask Buff 33lb/15kg	US\$ 45.73	+Info
469830	Buffstone Buff 33lb/15kg	US\$ 44.54	+Info

Gypsums & Die Stone Acc

Die Stones TYPE IV

NATURAL Die Stones TYPE-4



DIE STONE TYPE IV | ECCO
DIE STONE TYPE IV - ECCO

Available in 25lb and 50lb
A fine hard- accurate stone recommended for models- crown and bridge- casts and dies

25Lb

Code	Name	Price	+
112408	ECCO - DIE STONE 25lb GREEN- TYPE-4 Reg Set - Ultra fine die stone	US\$ 41.72	+Info
112410	ECCO - DIE STONE 25lb Violet- TYPE-4 Reg Set -	US\$ 41.73	+Info
112409	ECCO - DIE STONE 25lb BLUE - TYPE-4 Reg Set - 900-4025	US\$ 41.73	+Info

50Lb

Code	Name	Price	+
112412	ECCO - DIE STONE 50lb BLUE - TYPE-4 Reg Set - 900-4020	US\$ 65.67	+Info
112411	ECCO - DIE STONE 50lb GREEN- TYPE-4 Reg Set- #900-4021	US\$ 65.67	+Info
112413	ECCO - DIE STONE 50lb Violet- TYPE-4 Reg Set - 900-4023	US\$ 65.67	+Info
112414	ECCO - XL DIE STONE 50lb BLUE- TYPE-4 Reg Set - 900-4046 -	US\$ 69.07	+Info

ELITE ROCK Die Stone | ZHERMACK

ELITE ROCK - ZHERMACK

Type 4 extra hard dental die stone for fixed prostheses

Use:

- * Production of models and removable dies in fixed and removable prostheses

Characteristics:

- * High mechanical characteristics
- * Thixotropic
- * Accurate reproduction of detail
- * Available in five different colours: Sandy Brown, Light Amber, Silver Grey, Cream and White

Advantages:

- * Simplified mixing and balanced thixotropy for easier model production
- * High compressive strength
- * Reduced abrasion and higher precision of model
- * The silver grey colour is optimised for detection by modern laser and optical scanners
- * Low expansion, even after 48 hours



55 Lbs Die Stone

ELITE MASTER Stone 55lbs

Code	Name	Price	+
118576	ZHERMACK Elite Master - 55lb Carton - Economy Pack - Desert Sand - Z#C410403	US\$ 123.69	+Info
118580	ZHERMACK Elite Master-55lb Carton - Economy Pack - Sandy BROWN - Z#C410411	US\$ 123.69	+Info

ELITE Rock 55Lbs

Elite Rock

Elite Rock

Code	Name	Price	+
118555	Elite Rock 55 lb Carton Eco	US\$ 98.89	+Info
118556	Elite Rock 55 lb Carton Eco	US\$ 98.89	+Info
118557	Elite Rock 55 lb Carton Eco	US\$ 99.00	+Info
118558	Elite Rock 55 lb Carton Eco	US\$ 98.89	+Info
118561	Elite Rock FAST 55 lb Carton	US\$ 108.68	+Info

Other Models

Other Models

Code	Name	Price	+
118562	Elite Stone 55 lb Carton Eco	US\$ 100.09	+Info
118563	Elite Stone 55 lb Carton Eco	US\$ 107.73	+Info

6.6lbs Stone

Elite Rock

Elite Rock

Code	Name	Price	+
118542	Elite Rock 6.6 lb Bag Sandy	US\$ 22.01	+Info
118543	Elite Rock 6.6 lb Bag Brown	US\$ 22.01	+Info

Other Models

Other Models

Code	Name	Price	+
118545	Elite Stone 6.6 lb Bag Navy	US\$ 22.00	+Info

PRIMA ROCK Gypsum Type-4 | WHIP MIX
GYPSUMS PRIMA ROCK TYPE 4 - WHIP MIX

PRIMA-ROCK has the highest compressive strength of any Whip Mix die material. Vacuum spatulation provides a smooth mix which flows readily. It's compatible with all types of impression materials and upon setting, PRIMA-ROCK develops a dense surface with superior hardness.



ISO Type 4

YELLOW

Code	Name	Price	+
119640	WHIP MIX Prima Rock YELLOW 120/70g - #16632 - TYPE-4 -	US\$ 150.81	+Info

VIOLET

Code	Name	Price	+
119638	WHIP MIX Prima Rock Violet 55lb/25kg - #31712 - TYPE-4 -	US\$ 294.93	+Info
119641	WHIP MIX Prima Rock Violet 120/70g - #16624 - TYPE-4 -	US\$ 150.81	+Info

SILKY ROCK Gypsum Type-4 | WHIP MIX
GYPSUMS SILKY ROCK TYPE 4 - WHIP MIX

SILKY-ROCK is characterized by exceptional hardness and density throughout, and an unusual "silky-smooth" pouring quality that makes it flow freely into any type of impression. It works well with all types of impression materials. These properties, together with its low setting expansion, make SILKY-ROCK one of the world's most popular die materials.



ISO Type 4

VIOLET

Code	Name	Price	+
------	------	-------	---

119673	WHIP MIX Silky Rock Violet 33lb/15kg #330026	US\$ 58.26	+Info
119677	WHIP MIX Silky Rock Violet 120/70g - #01945-TYPE-4	US\$ 112.50	+Info
119671	WHIP MIX Silky Rock Violet 80/140g - #01953-TYPE-4	US\$ 112.50	+Info

YELLOW

Code	Name	Price	+
119680	WHIP MIX Silky Rock YELLOW 120/70g - #02038-TYPE-4	US\$ 112.50	+Info
119681	WHIP MIX Silky Rock YELLOW 80/140g - #02046-TYPE-4	US\$ 112.50	+Info
119679	WHIP MIX Silky Rock YELLOW 25#/11kg - #01996-TYPE-4	US\$ 40.37	+Info

WHITE

Code	Name	Price	+
119684	WHIP MIX Silky Rock WHITE 33lb/15kg - #330027	US\$ 58.26	+Info
119683	WHIP MIX Silky Rock WHITE 80/140g - #02119-TYPE-4	US\$ 112.50	+Info
119682	WHIP MIX Silky Rock WHITE 120/70g - #02100-TYPE-4	US\$ 112.50	+Info

IVORY

Code	Name	Price	+
119678	WHIP MIX Silky Rock Ivory 33lb/15kg - #330023	US\$ 66.10	+Info

LIGHT BROWN

Code	Name	Price	+
119676	WHIP MIX Silky Rock Light BROWN 33lb/15kg - #330025	US\$ 58.26	+Info

ULTI ROCK Gypsum Type-4 | WHIP MIX

GYPSUMS ULTI ROCK TYPE 4 - WHIP MIX

UltiRock is an exceptionally creamy, hard premium stone for every type of gypsum application. Its smooth non-reflective surface and high durability makes trimming full models, ditching dies and scanning easy and reliable.

In addition, UltiRock has a silky smooth consistency while mixing and the stones flowability and precision enables it to reproduce minute detail.

ISO Type 4



GOLDEN BROWN

Code	Name	Price	+
119696	WHIP MIX Ultirock Golden BROWN 11lb/5kg - #35025-TYPE-4	US\$ 58.33	+Info
119692	WHIP MIX Ultirock Golden BROWN 33lb - #330040 TYPE-4	US\$ 100.77	+Info

RESIN Die Stones TYPE-4



RESIN DIE STONE Type 4 | ECCO
RESIN DIE STONE - ECCO

ISO Type 4

The unique blend of synthetic resin and alpha gypsum results in improved surface smoothness and increased resistance to abrasion. Its low setting expansion and dimensional stability make it ideal for implant and complex restorative cases where precision and accuracy are required.



25 Lb

Code	Name	Price	+
112422	ECCO - RESIN DIE STONE 25lb BLUE- TYPE-IV (TYPE-4) -#900-4039 inComparable wi...	US\$ 66.86	+Info
112423	ECCO - RESIN DIE STONE 25lb GREEN - TYPE-IV (TYPE-4) -inComparable with Royal...	US\$ 52.00	+Info

50 Lb

Code	Name	Price	+
112424	ECCO - RESIN DIE STONE 50lb BLUE- TYPE-IV (TYPE-4) #900-4041	US\$ 88.50	+Info
112425	ECCO - RESIN DIE STONE 50lb GREEN - TYPE-IV (TYPE-4) - #900-4042	US\$ 119.67	+Info

Resin Die Stones TYPE 4 - GC AMERICA

Resin Die Stones TYPE 4 - GC AMERICA



12kg

12kg

Code	Name	Price	+
421750	GC Fujirock EP G.Brown 12kg	US\$ 135.56	+Info

22.7kg

22.7kg

Code	Name	Price	+
470492	New Fujirock Alabaster Ivory 50lb	US\$ 149.00	+Info
474517	Fujirock G. Brown Fast Set 22.7Kg	US\$ 149.00	+Info
476499	New Fujirock IMP IMP Golden Brown	US\$ 142.20	+Info

3kg

3kg

Code	Name	Price	+
421394	Fujirock Alabast IVORY Reg Set 3kg	US\$ 23.93	+Info
421689	Fujirock Golden Brown Fast Set 3kg	US\$ 23.13	+Info
421790	New Fujirock Alabaster Ivory Regula	US\$ 23.13	+Info
421817	New Fujirock IMP Golden Brown Reg	US\$ 23.13	+Info
421818	Fujirock Golden BROWN Reg Set 3kg	US\$ 23.13	+Info

5kg

5kg

Code	Name	Price	+
421376	GC Fujirock EP Golden Brown 5kg	US\$ 53.43	+Info
421811	GC Fujirock EP Type IV 5kg	US\$ 53.43	+Info

RESINROCK Resin Die Stone Type 4 | WHIP MIX

GYPSUMS RESINROCK TYPE 4 - WHIP MIX

ResinRock, a new, resin fortified die stone produces models of the highest quality. The unique blend of synthetic resin and alpha gypsum results in improved surface smoothness and increased resistance to abrasion. Its low setting expansion and dimensional stability make it ideal for implant and complex restorative cases where precision and accuracy are required. ResinRock is easy to work with - its creamy, thixotropic mix pours easily under vibration and stacks readily.



Available in seven colors: peach, blue, grey, ivory, white, green and golden brown. Peach, blue, green and golden brown provide excellent contrast to casting waxes. The neutral tones of grey, ivory and white complement all ceramic materials.

Please note the use of Gypsum Hardener is not recommended with this product.

ISO Type 4

IVORY

Code	Name	Price	
119656	WHIP MIX Resinrock Ivory 25 bags / 70gr Box - #34973	US\$ 53.18	+Info
119660	WHIP MIX Resinrock Ivory 5kg 2/2.5kg - #34970	US\$ 61.75	+Info
119658	WHIP MIX Resinrock Ivory 25/70g 4 Box Case - #34972	US\$ 186.22	+Info

PEACH

Code	Name	Price	
119663	WHIP MIX Resinrock Peach 25/70g Box - #34983	US\$ 53.18	+Info
119667	WHIP MIX Resinrock PEACH 5kg 2/2.5kg - #34916	US\$ 56.14	+Info
119661	WHIP MIX Resinrock Peach 25/70g 4 Box Case - #34991	US\$ 186.22	+Info

BLUE

Code	Name	Price	
119665	WHIP MIX Resinrock BLUE 5kg (2x2.5kg) - #34924	US\$ 61.75	+Info

SNAP STONE Resin Die Stone Type 4

| WHIP MIX

GYPSUMS SNAP STONE - WHIP MIX

SNAP-STONE is an exceptionally fast setting stone with an early high compressive strength, which allows separation of the cast from the impression in 5 minutes. Because SNAP-STONE is fast setting, hand mixing is recommended. SNAP-STONE is compatible with all impression materials and its versatile, light pink color makes it suitable for many uses.



Procedures which can be performed faster with SNAP-STONE include:

- * Bleaching trays
- * Mouth guards
- * Vacuum formed splints
- * Custom trays
- * Orthodontic appliances
- * Temporaries
- * Denture repairs
- * Study casts
- * Articulator cast mounting
- * CAD-CAM restorations

ISO Type 4

Physical Properties:

- * Water Powder Ratio: 23 ml / 100 g
- * Working Time: 60 - 90 seconds
- * Setting Time: 2 minutes
- * Setting Expansion: 0.15 %
- * Compressive Strength Wet (1 hr): 6,000 psi
- * Compressive Strength Dry (48 hrs): 14,000 psi

Code	Name	Price	+
119730	WHIP MIX Snap-Stone PINK 33lb/15kg (Gypsums Snap Stone) #330080	US\$ 82.27	+Info
470331	WHIP MIX Snap-Stone PINK 80-140 g Package (Gypsums Snap Stone) #26514	US\$ 108.23	+Info

Gypsums & Die Stone Acc
Die Stones TYPE V



GYPsums HARD ROCK TYPE 5 - WHIP MIX

GYPsums HARD ROCK TYPE 5 - WHIP MIX

The material is compatible with all types of impression materials. Its high expansion makes it especially suited for polyvinyl or polyether impression materials. This material has a creamy and thixotropic consistency for good flow under vibration, yet is very stackable. Its extended working time facilitates multiple pours. The stone has a smooth non-brittle surface so dies trim easily without shipping.



ISO Type 5

Code	Name	Price	+
119645	WHIP MIX Hard Rock GREEN 120/70g - #32077-TYPE-4	US\$ 150.57	+Info

Gypsums & Die Stone Acc
MOUNTING Stone Plasters



GYPsums - MOUNTING PLASTER - WHIP MIX

GYPsums - MOUNTING PLASTER - WHIP MIX

These products are specially formulated plaster which offers a low setting expansion and fast setting time. Its low compressive strength allows mounting plates to be easily cleaned. * Incorporate 15 seconds, hand spatulate 30 seconds.



ISO Type 1

Code	Name	Price	+
119740	WHIP MIX Mounting Plaster 33lb/15kgs - #330090 (Gypsums Mounting Plaster)	US\$ 54.26	+Info
119739	WHIP MIX Mounting Plaster 100/100g - #15550 (Gypsums Mounting Plaster)	US\$ 106.33	+Info

GYPSUMS - MOUNTING STONE - WHIP MIX

GYPSUMS - MOUNTING STONE - WHIP MIX

Low-Expansion, Fast-Setting Stone

A very hard, fast setting gypsum material formulated expressly for mounting casts to articulators. Its very low setting expansion prevents inaccuracies in occlusal cast relationships that result when conventional model stone or plaster is used. Although MOUNTING STONE has excellent stacking qualities, it moves very freely under the slightest pressure during its working time. It can be used as a model stone when a rapid set is required.



ISO Type 3

WHITE

Code	Name	Price	+
119728	WHIP MIX Mounting Stone WHITE 33lb/15kg - #330059(Gypsums Mounting Stone)	US\$ 75.13	+Info

BLUE

Code	Name	Price	+
119727	WHIP MIX Mounting Stone BLUE 33lb/15kg - #330058 (Gypsums Mounting Stone)	US\$ 64.28	+Info

Gypsums & Die Stone Acc

CAD/CAM Scannable Die Stone



CAD-CAM Scannable Die Stone - PEMACO

CAD-CAM Scannable Die Stone - PEMACO

Code	Name	Price	+
111386	CAD-CAM Stone-Plaster 5kg	US\$ 45.00	+Info

HATHO
POLISHING KIT

Thermoplastics
like *never* before



Packing Diamond D®
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory



GYPSUMS LEAN ROCK TYPE 5 - WHIP MIX
GYPSUMS LEAN ROCK TYPE 5 - WHIP MIX

Lean Rock, ISO Type 5 die stone, offers increased expansion to compensate for shrinkage of vinyl polysiloxane and other rubber-based impression materials. LeanRock has the same great qualities as Lean Rock which is specifically formulated for laboratories who want to pour small batches (about 6) and have working models ready in 15 minutes.

Small batch processing is based on lean principals to improve work flow and reduce wait time because models can be trimmed and pinned after only 15 minutes, decreasing downtime waiting for stone to set. This enables smaller batches to be sent down the line for trimming, eliminating the "bottleneck" associated with large batches hitting all at once, giving you more controlled release of product.

Description:

- * Saves time - Models can be separated and trimmed after only 15 minutes, producing better process flow and less downtime.
- * Scannable - Models may be used with Cad Cam systems eliminating the need for a specialty scan stone or spray.
- * Strength - Lean Rock is hard with superior abrasion resistance.
- * Clear Margins- Rich Ivory color provides clean, clear margins.

The decreased processing time also speeds up remakes and those emergency repairs needed in a hurry.

Directions for use - Hand incorporate for 10-15 seconds then mechanically mix under vacuum for 30 seconds.

ISO Type 5



IVORY

Code	Name	Price	+
119688	WHIP MIX Lean Rock Ivory 33lb/15kg - #330030 (Cad Cam Scannable Die Stone gypsum)	US\$ 87.35	+Info

BLUE

Code	Name	Price	+
119687	WHIP MIX Lean Rock BLUE 33lb/15kg - #330029 (Cad Cam Scannable Die Stone gypsum)	US\$ 79.40	+Info

Gypsums & Die Stone Acc

POURABLE Base Stone



POURABLE Base Stone - HOUSE BRAND

POURABLE Base Stone - HOUSE BRAND

Code	Name	Price	+
112417	POUR BASE STONE 50lb Buff	US\$ 52.70	+Info
112418	POUR BASE STONE 50lb White	US\$ 52.70	+Info

POURABLE Base Stone - WHIP MIX

POURABLE Base Stone - WHIP MIX

Code	Name	Price	+
119713	Flowstone Buff 33lb/15kg	US\$ 71.49	+Info
119714	Flowstone Blue 50lb /22kg	US\$ 103.81	+Info
119718	Flowstone Fs White 33lbs/15kg	US\$ 80.06	+Info

Gypsums & Die Stone Acc

Gypsum Accessories



Accessories for Plasters - ZHERMACK

Accessories for Plasters - ZHERMACK



Code	Name	Price	+
118594	Gypsum Measuring Set	US\$ 4.62	+Info

Gypsums & Die Stone Acc

Gypsum Hardners



FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US DENTAL DEPOT

+1 954 874 6325
www.usdentaldepot.com



Gypsum Hardners - BUFFALO

Gypsum Hardners - BUFFALO



Code	Name	Price	+
424570	Stalite Solut Gypsum Hardener 1ga	US\$ 41.33	+Info

GYPSUMS HARDENER - WHIP MIX

GYPSUMS HARDENER - WHIP MIX

Replaces water as the gauging liquid for mixing with gypsum materials where increased hardness and reduced brittleness are desired; will also improve the "running" or flowing quality of the mix. It does not clinically affect the other properties of the die or model materials or plasters. This formula is not damaged by freezing.

This product is not recommended for gypsums that contain resin materials, such as Resin Rock



Code	Name	Price	+
120415	WHIP MIX Gypsum Hardener 4gal Case - #27774	US\$ 301.71	+Info
120416	WHIP MIX Gypsum Hardener Gallon/3.75l - #09725	US\$ 84.61	+Info
120417	WHIP MIX Gypsum Hardener Pint/450ml - #09717	US\$ 22.56	+Info

MODEL EEZ - ADS

MODEL EEZ - ADS

Use Model-Eez, instead of water, to produce models that are smoother, harder and stronger. Avoid broken, chipped or cracked models. Model-Eez represents a major breakthrough in creating perfect models with all commercially available gypsums. Conditions the gypsum and produces a low-viscosity, high-density, smooth, bubble-free mix. Model-Eez plaster mix pours smoothly and easily into impressions to capture the fine oral detail.



Code	Name	Price	+
118215	ADS Model Eez qt - #M756-6 (mixing liquid for harder stone results)	US\$ 14.12	+Info
118216	ADS Model Eez gal - (mixing liquid for harder stone results) #M756-8	US\$ 32.23	+Info

ARTICULATOR & FLASK CLEANER - ADS ARTICULATOR & FLASK CLEANER - ADS

Features:

- * Removes plaster, stone or investment without ultrasonic
- * Is a heavy duty laboratory tested concentrate
- * Formulated to penetrate, separate and disintegrate gypsum based materials
- * Surfactant action formula gets under thick deposits which then lift off, leaving a clean bright and tarnish free flask and articulator
- * Non-hazardous, safe for all hard surfaces
- * Caustic free solution is mild to hands
- * Biodegradable, economical overnight soak
- * Removes lubricants and stone from articulators and flasks without scrubbing or scraping
- * Solubilizing action penetrates joints, crevices and delicate areas
- * Rinse clean with water



US DENTAL DEPOT is an Authorized American Dental Supply dealer in US

Code	Name	Price	+
117660	ADS Articulator & Flask CLEANER qt - #A946-6	US\$ 18.29	+Info
117661	ADS Articulator & Flask CLEANER gal - #A946-8	US\$ 45.02	+Info

GYP SUMS REMOVERS - GYP STRIPS - WHIP MIX GYP SUMS REMOVERS - GYP STRIPS - WHIP MIX

A highly effective gypsum remover, not only keeps equipment sparkling clean, it saves time and energy, too. GYP-STRIP'S penetrating action removes plaster, stone and gypsum investments from mixing bowls, acrylic dentures, denture flasks, casting rings and non-aluminum articulators. GYP-STRIP is versatile, you can use it in an ultrasonic, or as an overnight soaking agent. Available in one or four gallon jugs, GYP-STRIP requires no mixing or diluting.



Code	Name	Price	+
120413	WHIP MIX Gyp Strip Gypsum Remover Gal/3.75l - #27014	US\$ 96.73	+Info
120414	WHIP MIX Gyp Strip Gypsum Remover 4gal Case - #27022	US\$ 309.39	+Info

PLASTER & STONE REMOVER - ADS

PLASTER & STONE REMOVER - ADS

Designed to give you the safest most effective product. Solubilizing action rapidly penetrates joints, crevices and delicate areas better than any other cleaner. Leaves denture base clean, and looking like new. Safe - no hazard shipping restrictions. Non-caustic, non-corrosive. Safe for use on aluminum. Use directly in the ultrasonic tank, or in glass or plastic beaker. Safe for hands, no rubber gloves needed. Can be shipped in any weather. Biodegradable product can be flushed down lab drain for disposal. Removes tarnish and oxidation and leaves instruments shining. High capacity for calcium and magnesium compounds. Dissolves large quantities of gypsum, lasts longer, and saves time. Plaster deposits disintegrate and dissolve in minutes.



Code	Name	Price	+
118290	ADS Plaster & Stone Remover Solution qt - #P926-6	US\$ 11.99	+Info
118291	ADS Plaster & Stone Remover Solution gal - #P926-8	US\$ 30.33	+Info
118287	ADS Plaster & Stone REMOVER powder 1Lb (to make 1gal) - #P926-9	US\$ 30.33	+Info
118288	ADS Plaster & Stone REMOVER powder 5Lb (to make 5 gal) - #P926-91	US\$ 78.44	+Info
118289	ADS Plaster & Stone REMOVER powder 10Lb (to make 10 gal) - #P926-92	US\$ 120.27	+Info

Plasters & Gypsum Removers - COLTENE

Plasters & Gypsum Removers - COLTENE

Code	Name	Price	+
427074	BioSonic Plaster and Stone Remover,	US\$ 41.25	+Info

Plasters & Gypsum Removers - JSP

Plasters & Gypsum Removers - JSP

Code	Name	Price	+
437700	Plaster & Stone Remover 32 oz	US\$ 14.11	+Info

Plasters & Gypsum Removers - KEYSTONE

Plasters & Gypsum Removers - KEYSTONE

Code	Name	Price	+
034-1850065	DP2400 Plaster & Stone Remover	US\$ 30.60	+Info

Plasters & Gypsum Removers - L&R

Plasters & Gypsum Removers - L&R

Code	Name	Price	+
402004	Plaster & Stone Remover 1gal 4cs	US\$ 113.27	+Info

Plasters & Gypsum Removers - RENFERT

Plasters & Gypsum Removers - RENFERT

Code	Name	Price	+
417583	GO 2011 Speed Remover 2lt	US\$ 63.84	+Info

Catalog Index / Indice

Lab Plasters TYPE II	1
Lab Stones TYPE III	2
Die Stones TYPE IV	4
Die Stones TYPE V	12
MOUNTING Stone Plasters	12
CAD/CAM Scannable Die Stone	13
POURABLE Base Stone	14
Gypsum Accessories	15
Gypsum Hardners	15
Plasters & Gypsum Removers	16