

Inlay Resins

Light Cure Lab Composite

CERAMAGE Light Cure Lab Composite | SHOFU



CERAMAGE KITS CERAMAGE KITS - SHOFU

Ceramage contains 73% by weight of micro-fine ceramic particle, which are embedded in an organic polymer matrix and combines the advantages of ceramic and composite to give a light-curing "esthetic" composite for use in any indication. This composition increases the homogenous structure resulting in optimised physical properties and enhanced esthetics.

Indications:

- * Fully veneered crowns and bridges
- * Fully veneered telescopic and implant-supported restorations
- * Precision attachments
- * Occlusal facings
- * Inlays / onlays
- * Veneers
- * Long-term temporary restorations
- * Modifying conventional acrylic teeth

Benefits:

- * Excellent flexural and compressive strength
- * Excellent abrasion resistance of occlusal facings
- * Highest aesthetics with a light refractive index similar to natural teeth
- * Extensive range of translucent effect colors
- * High color stability
- * Perfect viscosity and modelling properties



Code	Name	Price	+
117314	CERAMAGE Intro Set A2 # 1807	US\$ 318.35	+Info
116836	CERAMAGE Intro Set A3 # 1808	US\$ 318.35	+Info
112457	CERAMAGE Standard Set Kit #1806 (As ddiscontinued from factory, we will provid...	US\$ 1,625.80	+Info

CERAMAGE BODY

CERAMAGE BODY - SHOFU

Ceramage is a zirconium silicate micro ceramic which creates indirect restorations that exhibit virtually the same light transmission as natural teeth, maintaining a remarkable translucency. Ceramage bonds to a variety of substructures, including non-precious and high noble alloys.

With ceramage, an extensive shade selection for natural tooth and gum color reproduction, in addition to new bleach white shades, ceramage can be used to create superior anterior and posterior crowns, veneers, implant-supported restorations and inlays and onlays. Ceramage full set of gum colors are also ideal for replicating gingival anatomy.

Ceramage Advantages:

- * Microfine structure gives high compressive strength
- * Does not vitrify, maintaining excellent color stability
- * Outstanding tensile strength and elasticity
- * Superior wear characteristics

Authorized SHOFU Dealer



Code	Name	Price	+
111378	CERAMAGE Body A1 #1901 - 2.6ml / 4.6g Syringe	US\$ 31.35	+Info
110547	Ceramage Body A2 #1902 - 2.6ml / 4.6g Syringe #1902	US\$ 45.60	+Info
110548	CERAMAGE Body A3 #1903 - 2.6ml / 4.6g Syringe #1903	US\$ 45.60	+Info
110549	CERAMAGE Body A3.5 - 2.6ml / 4.6g Syringe #1904	US\$ 45.60	+Info
110550	CERAMAGE Body A4 - 2.6ml / 4.6g Syringe #1905	US\$ 45.60	+Info
110551	CERAMAGE Body B1 - 2.6ml / 4.6g Syringe #1906	US\$ 45.60	+Info
110552	CERAMAGE Body B2 - 2.6ml / 4.6g Syringe #1907	US\$ 45.60	+Info
110553	CERAMAGE Body B3 - 2.6ml / 4.6g Syringe #1908	US\$ 44.69	+Info
110554	CERAMAGE Body B4 - 2.6ml / 4.6g Syringe #1909	US\$ 44.69	+Info
110555	CERAMAGE Body C1 - 2.6ml / 4.6g Syringe #1910	US\$ 45.60	+Info
110556	CERAMAGE Body C2 - 2.6ml / 4.6g Syringe #1911	US\$ 45.14	+Info
110557	CERAMAGE Body C3 - 2.6ml / 4.6g Syringe #1912	US\$ 45.60	+Info
110558	CERAMAGE Body C4 - 2.6ml / 4.6g Syringe #1913	US\$ 45.60	+Info
110559	CERAMAGE Body D2 - 2.6ml / 4.6g Syringe #1914	US\$ 44.69	+Info
110560	CERAMAGE Body D3 - 2.6ml / 4.6g Syringe #1915	US\$ 44.69	+Info
110561	CERAMAGE Body D4 - 2.6ml / 4.6g Syringe #1916	US\$ 45.60	+Info
110562	CERAMAGE RTA Body #1917 - 2.6ml / 4.6g Syringe #1917	US\$ 45.60	+Info
110563	CERAMAGE R2 Body #1918 - 2.6ml / 4.6g Syringe #1918	US\$ 45.60	+Info
110564	CERAMAGE R3 Body #1919 - 2.6ml / 4.6g Syringe #1919	US\$ 45.60	+Info
110565	CERAMAGE R3.5 Body #1920 - 2.6ml / 4.6g Syringe #1920	US\$ 45.60	+Info

CARAMAGE UP



Code	Name	Price	+
470779	Ceramage UP Body A1B	US\$ 23.58	+Info
470793	Ceramage UP Incisal 59 5gr	US\$ 24.08	+Info

CERAMAGE ACCESSORIES

CERAMAGE ACCESSORIES - SHOFU

Using the experience and knowledge gained over many years, SHOFU has developed Ceramage.

Ceramage contains 73% by weight of micro-fine ceramic particle, which are embedded in an organic polymer matrix and combines the advantages of ceramic and composite to give a light-curing "esthetic" composite for use in any indication. This composition increases the homogenous structure resulting in optimised physical properties and enhanced esthetics.

Ceramage combines natural color reproduction with extraordinary strength and elasticity and is recommended for conventional cases as well as implant restorations.

Indications:

- * Fully veneered crowns and bridges
- * Fully veneered telescopic and implant-supported restorations
- * Precision attachments
- * Occlusal facings
- * Inlays / onlays
- * Veneers
- * Long-term temporary restorations
- * Modifying conventional acrylic teeth

Benefits:

- * Excellent flexural and compressive strength
- * Excellent abrasion resistance of occlusal facings
- * Highest aesthetics with a light refractive index similar to natural teeth
- * Extensive range of translucent effect colors
- * High color stability



- * Perfect viscosity and modelling properties
- * Superb polishability
- * Outstanding plaque resistance
- * Excellent biocompatibility

Code	Name	Price	+
111380	CERAMAGE Metal Photo Primer 5ml #1740	US\$ 54.06	+Info
111381	CERAMAGE Modeling Liquid 6ml #1991	US\$ 33.50	+Info
111382	CERAMAGE Separator Liquid 7ml #1993	US\$ 33.50	+Info
111383	CERAMAGE Spacer 7ml #1994	US\$ 22.57	+Info
116873	CERAMAGE Cleaner #1995	US\$ 33.50	+Info
113540	CERAMAGE Oxybar #1992 - 10 ml Syringe Oxibarrier	US\$ 32.83	+Info

CERAMAGE CERVICAL

CERAMAGE CERVICAL - SHOFU

Ceramage is a zirconium silicate micro ceramic which creates indirect restorations that exhibit virtually the same light transmission as natural teeth, maintaining a remarkable translucency. Ceramage bonds to a variety of substructures, including non-precious and high noble alloys.

With ceramage, an extensive shade selection for natural tooth and gum color reproduction, in addition to new bleach white shades, ceramage can be used to create superior anterior and posterior crowns, veneers, implant-supported restorations and inlays and onlays. Ceramage full set of gum colors are also ideal for replicating gingival anatomy.

Ceramage Advantages:

- * Microfine structure gives high compressive strength
- * Does not vitrify, maintaining excellent color stability
- * Outstanding tensile strength and elasticity
- * Superior wear characteristics

Raves And Recommendations:

- * Ceramage exhibits perfect bond strength to any alloy, even to titanium.
- * Ceramage can be used on existing porcelain restorations.
- * For the best glaze on ceramage, use dura polish and dura polish-dia

"As a dental laboratory owner, it's important to me to have materials that meet the demands of the ever-changing dental industry and ceramage, far exceeds other composite materials on the market in versatility, strength, and esthetics.

Us Dental Depot is an Authorized SHOFU Dealer in US.



Code	Name	Price	+
110566	CERAMAGE Cervical AC1 #1881 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110567	CERAMAGE Cervical AC2 #1882 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info

110568	CERAMAGE Cervical BC1 #1883 - 2.6ml / 4.6g Syringe	US\$ 31.92	+Info
110569	CERAMAGE Cervical BC2 #1884 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110570	CERAMAGE Cervical CC1 #1885 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110571	CERAMAGE Cervical CC2 #1886 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
116840	CERAMAGE DC1 2.6ml / 4.6g Syringe #1887	US\$ 45.60	+Info
110572	CERAMAGE Cervical DC2 #1888 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info

CERAMAGE DENTIN BLEACH SHADES

CERAMAGE DENTIN BLEACH SHADES - SHOFU



Using the experience and knowledge gained over many years, SHOFU has developed Ceramage.

Ceramage contains 73% by weight of micro-fine ceramic particle, which are embedded in an organic polymer matrix and combines the advantages of ceramic and composite to give a light-curing "esthetic" composite for use in any indication. This composition increases the homogenous structure resulting in optimised physical properties and enhanced esthetics.

Ceramage combines natural color reproduction with extraordinary strength and elasticity and is recommended for conventional cases as well as implant restorations.

Indications:

- * Fully veneered crowns and bridges
- * Fully veneered telescopic and implant-supported restorations
- * Precision attachments
- * Occlusal facings
- * Inlays / onlays
- * Veneers
- * Long-term temporary restorations
- * Modifying conventional acrylic teeth

Benefits:


- * Excellent flexural and compressive strength
- * Excellent abrasion resistance of occlusal facings
- * Highest aesthetics with a light refractive index similar to natural teeth
- * Extensive range of translucent effect colors
- * High color stability
- * Perfect viscosity and modelling properties
- * Superb polishability
- * Outstanding plaque resistance
- * Excellent biocompatibility

Code	Name	Price	+
112879	CERAMAGE Whitening Super WHITE 2.6ml / 4.6g Syringe #1924 W0B -	US\$ 45.60	+Info
112880	CERAMAGE Whitening WHITE Shade 1 2.6ml / 4.6g Syringe #1925 W1B -	US\$ 45.60	+Info
112881	CERAMAGE Whitening WHITE Shade 2 2.6ml / 4.6g Syringe #1926 W2B -	US\$ 45.60	+Info
112882	CERAMAGE Whitening WHITE Shade 3 2.6ml / 4.6g Syringe #1927 W3B -	US\$ 45.60	+Info

CERAMAGE FLOWABLE RESIN

CERAMAGE FLOWABLE RESIN - SHOFU

- * Ceramage is an easy-to-apply paste
- * Accurate reproduction of the natural tooth shade
- * Its abrasion resistance ensures molars are protected from the opposing dentition
- * Flowable Composite Resin suitable for adjustments to the build-up and repairing small areas
- * In combination with M.L. Primer produces a strong bond to any type of metal, including precious alloys



SALE OF INDIVIDUAL SYRINGE FOR REPLACEMENT

FLOWABLE

Code	Name	Price	+
116841	CERAMAGE Flowable A3B #1921 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
116839	CERAMAGE Flowable-Opacious Dentin A3 ODA3 #1871 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
116842	CERAMAGE Flowable F-Rootab #1922 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110600	CERAMAGE Flowable Enamel F-59 #1897 - 2.6ml / 4.6g Syringe same as incisal	US\$ 44.69	+Info

FLOW GUM

Code	Name	Price	+
110641	CERAMAGE Flowable F-T Glass #1941 - 2.6ml / 4.6g Syringe	US\$ 41.50	+Info
110643	CERAMAGE Flowable F-Gum BROWN #1946 - 2.6ml / 4.6g Syringe	US\$ 41.04	+Info
110644	CERAMAGE Flowable F-Gum RED #1948 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
116847	CERAMAGE Flowable F-Gum-V #1947 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110642	CERAMAGE Flowable F-Gum WHITE #1943 - 2.6ml / 4.6g Syringe -	US\$ 45.60	+Info

CERAMAGE GUM COLORS

CERAMAGE GUM COLORS - SHOFU

Ceramage is a zirconium silicate micro ceramic which creates indirect restorations that exhibit virtually the same light transmission as natural teeth, maintaining a remarkable translucency. Ceramage bonds to a variety of substructures, including non-precious and high noble alloys.

With ceramage, an extensive shade selection for natural tooth and gum color reproduction, in addition to new bleach white shades, ceramage can be used to create superior anterior and posterior crowns, veneers, implant-supported restorations and inlays and onlays. Ceramage full set of gum colors are also ideal for replicating gingival anatomy.

Ceramage Advantages:

- * Microfine structure gives high compressive strength
- * Does not vitrify, maintaining excellent color stability
- * Outstanding tensile strength and elasticity



*** Superior wear characteristics**

Gum Colors Ideal For Gingival Modification, Implant Cases And Previously Fixed Porcelain Restorations

Shofu Raves And Recommendations

- * Ceramage exhibits perfect bond strength to any alloy, even to titanium.
- * Ceramage can be used on existing porcelain restorations.
- * For the best glaze on ceramage, use dura polish and dura polish-dia

"As a dental laboratory owner, it's important to me to have materials that meet the demands of the ever-changing dental industry and ceramage, far exceeds other composite materials on the market in versatility, strength, and esthetics.

Us Dental Depot is an Authorized SHOFU Dealer in US.

CERAMAGE GUM

Code	Name	Price	+
110637	CERAMAGE Gum Light #1959 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110638	CERAMAGE Gum Dark #1960 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110639	CERAMAGE Gum ORANGE #1961 - 2.6ml / 4.6g Syringe	US\$ 44.69	+Info
112180	CERAMAGE Gum Translucent - 2.6ml / 4.6g Syringe #1940	US\$ 42.86	+Info

CERAMAGE INCISAL / ENAMEL

CERAMAGE INCISAL / ENAMEL - SHOFU



Ceramage is a zirconium silicate micro ceramic which creates indirect restorations that exhibit virtually the same light transmission as natural teeth, maintaining a remarkable translucency. Ceramage bonds to a variety of substructures, including non-precious and high noble alloys.

Shofu Raves and Recommendations:

- * Ceramage exhibits perfect bond strength to any alloy, even to titanium.
- * Ceramage can be used on existing porcelain restorations.
- * For the best glaze on Ceramage, use Dura Polish and Dura Polish-Dia

Authorized SHOFU Dealer

Code	Name	Price	+
110594	CERAMAGE Enamel 56 #1891 - 2.6ml / 4.6g Syringe same as incisal	US\$ 45.60	+Info
110595	CERAMAGE Enamel 57 #1892 - 2.6ml / 4.6g Syringe same as incisal	US\$ 44.69	+Info
110596	CERAMAGE Enamel 58 #1893 - 2.6ml / 4.6g Syringe same as incisal	US\$ 45.60	+Info
110597	CERAMAGE Enamel 59 #1894 - 2.6ml / 4.6g Syringe same as incisal	US\$ 45.60	+Info
110598	CERAMAGE Enamel 60 #1895 - 2.6ml / 4.6g Syringe same as incisal	US\$ 45.60	+Info
110599	CERAMAGE Enamel 61 #1896 - 2.6ml / 4.6g Syringe same as incisal	US\$ 45.60	+Info

CERAMAGE INCISAL EFFECT MODIFIERS

CERAMAGE INCISAL EFFECT MODIFIERS - SHOFU

- * Ceramage is an easy-to-apply paste
- * Accurate reproduction of the natural tooth shade
- * Its abrasion resistance ensures molars are protected from the opposing dentition
- * Flowable Composite Resin suitable for adjustments to the build-up and repairing small areas
- * In combination with M.L. Primer produces a strong bond to any type of metal, including precious alloys



SALE OF INDIVIDUAL SYRINGE FOR REPLACEMENT

Code	Name	Price	+
116848	CERAMAGE - ENAMEL EFFECTS Milky Syringe 4.6g #1951	US\$ 45.60	+Info
116849	CERAMAGE - ENAMEL EFFECTS WHITE Enamel Syringe 4.6gr #1952 -	US\$ 45.60	+Info
116850	CERAMAGE - ENAMEL EFFECTS Occlusal Syringe 4.6g #1953	US\$ 42.86	+Info
116851	CERAMAGE - ENAMEL EFFECTS Amber YELLOW Syringe 4.6g #1954	US\$ 45.60	+Info
116852	CERAMAGE - Concentrated ENAMEL EFFECTS Amber RED Syringe 4.6g #1955 AM-R	US\$ 28.73	+Info
116853	CERAMAGE - ENAMEL EFFECTS Amber Violet Syringe 4.6g #1956	US\$ 45.60	+Info
116854	CERAMAGE - ENAMEL EFFECTS YELLOW Syringe 4.6g #1957	US\$ 45.60	+Info
116855	CERAMAGE - ENAMEL EFFECTS Mamelon PINK Syringe 4.6g #1958	US\$ 45.60	+Info

CERAMAGE OPACIOUS DENTIN

CERAMAGE OPACIOUS DENTIN - SHOFU

Ceramage is a zirconium silicate micro ceramic which creates indirect restorations that exhibit virtually the same light transmission as natural teeth, maintaining a remarkable translucency. Ceramage bonds to a variety of substructures, including non-precious and high noble alloys. With ceramage, an extensive shade selection for natural tooth and gum color reproduction, in addition to new bleach white shades, ceramage can be used to create superior anterior and posterior crowns, veneers, implant-supported restorations and inlays and onlays. Ceramage full set of gum colors are also ideal for replicating gingival anatomy.



Ceramage Advantages:

- * Microfine structure gives high compressive strength
- * Does not vitrify, maintaining excellent color stability
- * Outstanding tensile strength and elasticity
- * Superior wear characteristics
- * Gum colors ideal for gingival modification, implant cases and previously fixed porcelain restorations

Shofu Raves and Recommendations:

- * Ceramage exhibits perfect bond strength to any alloy, even to titanium.
- * Ceramage can be used on existing porcelain restorations.
- * For the best glaze on Ceramage, use Dura Polish? and Dura Polish-Dia?

"As a dental laboratory owner, it's important to me to have materials that meet the demands of the ever-changing dental industry and Ceramage, far exceeds other composite materials on the market in versatility, strength, and esthetics."

Us Dental Depot is an Authorized SHOFU Dealer in US.

Code	Name	Price	+
110573	CERAMAGE Opacious Dentin A1 #1851 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110574	CERAMAGE Opacious Dentin A2 #1852 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110575	CERAMAGE Opacious Dentin A3 #1853 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110576	CERAMAGE Opacious Dentin A3.5 #1854 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110577	CERAMAGE Opacious Dentin A4 #1855 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110579	CERAMAGE Opacious Dentin B1 #1856 - 2.6ml / 4.6g Syringe	US\$ 44.69	+Info
110580	CERAMAGE Opacious Dentin B2 #1857 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110581	CERAMAGE Opacious Dentin B3 #1858 - 2.6ml / 4.6g Syringe	US\$ 40.13	+Info
110582	CERAMAGE Opacious Dentin B4 #1859 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
110583	CERAMAGE Opacious Dentin C1 #1860 - 2.6ml / 4.6g Syringe #1860	US\$ 45.60	+Info
110584	CERAMAGE Opacious Dentin C2 #1861 - 2.6ml / 4.6g Syringe #1861	US\$ 45.60	+Info
110585	CERAMAGE Opacious Dentin C3 #1862 - 2.6ml / 4.6g Syringe #1862	US\$ 45.60	+Info
110586	CERAMAGE Opacious Dentin C4 #1863 - 2.6ml / 4.6g Syringe #1863	US\$ 45.60	+Info
110587	CERAMAGE Opacious Dentin D2 #1864 - 2.6ml / 4.6g Syringe #1864	US\$ 45.60	+Info
110588	CERAMAGE Opacious Dentin D3 #1865 - 2.6ml / 4.6g Syringe #1865	US\$ 45.60	+Info
110589	CERAMAGE Opacious Dentin D4 #1866 - 2.6ml / 4.6g Syringe #1866	US\$ 45.60	+Info
110590	CERAMAGE Opacious Root A - 2.6ml / 4.6g Syringe #1867	US\$ 45.60	+Info
110591	CERAMAGE Opacious R2 #1868 - 2.6ml / 4.6g Syringe #1868	US\$ 45.60	+Info
110593	CERAMAGE Opacious Dentin R3.5 #1870 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info

CERAMAGE OPAQUE

CERAMAGE OPAQUE - SHOFU

Ceramage is a zirconium silicate micro ceramic which creates indirect restorations that exhibit virtually the same light transmission as natural teeth, maintaining a remarkable translucency. Ceramage bonds to a variety of substructures, including non-precious and high noble alloys. With Ceramage, an extensive shade selection for natural tooth and gum color reproduction, in addition to new bleach white shades, Ceramage can be used to create superior anterior and posterior crowns, veneers, implant-supported restorations and inlays and onlays. Ceramage full set of gum colors are also ideal for replicating gingival anatomy.

Advantages:

- * Microfine structure gives high compressive strength
- * Does not vitrify, maintaining excellent color stability
- * Outstanding tensile strength and elasticity
- * Superior wear characteristics
- * Gum colors ideal for gingival modification, implant cases and previously fixed porcelain restorations

Raves and Recommendations:



- * Ceramage exhibits perfect bond strength to any alloy, even to titanium.
- * Ceramage can be used on existing porcelain restorations.
- * For the best glaze on Ceramage, use Dura Polish? and Dura Polish-Dia?

As a dental laboratory owner, it's important to me to have materials that meet the demands of the ever-changing dental industry and Ceramage, far exceeds other composite materials on the market in versatility, strength, and esthetics.

Code	Name	Price	+
112885	CERAMAGE Pre-Opaque - 2.6ml / 4.6g Syringe #1820	US\$ 31.35	+Info
110607	CERAMAGE Opaque A4 #1825 - 2.6ml Syringe	US\$ 31.35	+Info
110610	CERAMAGE Opaque B3 #1828 - 2.6ml Syringe	US\$ 31.35	+Info
110611	CERAMAGE Opaque B4 #1829 - 2.6ml Syringe	US\$ 31.35	+Info
110612	CERAMAGE Opaque C1 #1830 - 2.6ml Syringe	US\$ 31.35	+Info
110613	CERAMAGE Opaque C20 #1831 - 2.6ml Syringe	US\$ 31.35	+Info
110614	CERAMAGE Opaque C3 #1832 - 2.6ml Syringe	US\$ 31.35	+Info
110615	CERAMAGE Opaque C4 #1833 - 2.6ml Syringe	US\$ 31.35	+Info
110616	CERAMAGE Opaque D2 #1834 - 2.6ml Syringe	US\$ 31.35	+Info
110618	CERAMAGE Opaque D4 #1836 - 2.6ml Syringe	US\$ 31.35	+Info
110619	CERAMAGE Root AO #1837 - 2.6ml Syringe	US\$ 31.35	+Info
110620	CERAMAGE R2 #1838 - 2.6ml Syringe	US\$ 31.35	+Info
110621	CERAMAGE R3 #1839 - 2.6ml Syringe	US\$ 31.35	+Info
110622	CERAMAGE R3.5 #1840 - 2.6ml Syringe	US\$ 31.35	+Info

CERAMAGE OPAQUE UNIVERSAL

CERAMAGE OPAQUE UNIVERSAL - SHOFU

The new light-cured opaquer pastes of the Universal Opaque system offer excellent masking properties and bond strengths on all commonly used dental alloys, thanks to specially optimized micro-fillers.

Improved thixotropic properties have now further increased softness and flowability.

The pastes are extremely easy and quick to apply. This ready-to-use opaquer system, consisting of one Pre-Opaque and 31 Opaque materials, can be used as a base for the reproduction of any shade of the Ceramage and Solidex composite systems, and also in combination with other crown and bridge composites.

Indication

Masking of frameworks made of precious and non-precious alloys

Benefits

- * Universal system



- * Easy and quick to apply
- * Excellent masking properties
- * Durable bond

Code	Name	Price	+
116874	CERAMAGE Universal Opaquers A1O #2080	US\$ 32.02	+Info
116875	CERAMAGE Universal Opaquers A2O # 2081	US\$ 30.10	+Info
116876	CERAMAGE Universal Opaquers A3O #2082	US\$ 31.35	+Info
116877	CERAMAGE Universal Opaquers A3.5O #2083	US\$ 32.02	+Info
116878	CERAMAGE Universal Opaquers A4O #2084	US\$ 32.02	+Info
116879	CERAMAGE Universal Opaquers B1O #2085	US\$ 32.02	+Info
116880	CERAMAGE Universal Opaquers B2O #2086	US\$ 32.02	+Info
116881	CERAMAGE Universal Opaquers B3O #2087	US\$ 32.02	+Info
116882	CERAMAGE Universal Opaquers B4O #2088	US\$ 32.02	+Info
116883	CERAMAGE Universal Opaquers C1O #2089	US\$ 32.02	+Info
116884	CERAMAGE Universal Opaquers C2O #2090	US\$ 32.02	+Info
116885	CERAMAGE Universal Opaquers C3O #2091	US\$ 32.02	+Info
116886	CERAMAGE Universal Opaquers C4O #2092	US\$ 32.02	+Info
116887	CERAMAGE Universal Opaquers D2O #2093	US\$ 32.02	+Info
116888	CERAMAGE Universal Opaquers D3O #2094	US\$ 32.02	+Info
116889	CERAMAGE Universal Opaquers D4O #2095	US\$ 32.02	+Info
116891	CERAMAGE Universal Opaquers R2O #2097	US\$ 32.02	+Info
116892	CERAMAGE Universal Opaquers R3O #2098	US\$ 32.02	+Info
116893	CERAMAGE Universal Opaquers R3.5O #2099	US\$ 32.02	+Info
116895	CERAMAGE Universal Opaquers W2O #2101	US\$ 32.02	+Info
116896	CERAMAGE Universal Opaquers W3O #2102	US\$ 32.02	+Info
116899	CERAMAGE Universal Opaquers Go #2105	US\$ 32.02	+Info
116898	CERAMAGE Universal Opaquers Bgo #2104	US\$ 32.02	+Info
116901	CERAMAGE Universal Opaquers Ino #2107	US\$ 32.02	+Info
116900	CERAMAGE Universal Opaquers Mo #2106	US\$ 32.02	+Info
116904	CERAMAGE Universal Opaquers Po #2110	US\$ 32.02	+Info
116905	CERAMAGE Universal Opaquers Pre-Opaquers #2111	US\$ 31.35	+Info
116890	CERAMAGE Universal Opaquers Root Ao #2096	US\$ 32.02	+Info
116902	CERAMAGE Universal Opaquers Vpo #2108	US\$ 32.02	+Info
116903	CERAMAGE Universal Opaquers Vro #2109	US\$ 32.02	+Info
116894	CERAMAGE Universal Opaquers Wo #2100	US\$ 31.35	+Info

CERAMAGE TRANSPARENT

CERAMAGE TRANSPARENT - SHOFU

Using the experience and knowledge gained over many years, SHOFU has developed Ceramage.

Ceramage contains 73% by weight of micro-fine ceramic particle, which are embedded in an organic polymer matrix and combines the advantages of ceramic and composite to give a light-curing "esthetic" composite for use in any indication. This composition increases the homogenous structure resulting in optimised physical properties and enhanced esthetics.



TRANSPARENT INCISAL MODIFIERS

Code	Name	Price	+
114931	CERAMAGE Translucent T #1931 - 2.6ml / 4.6g Syringe	US\$ 45.14	+Info
114932	CERAMAGE Translucent HVT #1932 - 2.6ml / 4.6g Syringe	US\$ 44.23	+Info
114933	CERAMAGE Translucent LVT #1933 - 2.6ml / 4.6g Syringe	US\$ 37.85	+Info
114934	CERAMAGE Translucent T-GLASS - TG #1934 - 2.6ml / 4.6g Syringe	US\$ 44.69	+Info
114935	CERAMAGE Translucent B-GLASS - BG #1935 - 2.6ml / 4.6g Syringe	US\$ 45.60	+Info
116843	CERAMAGE Grey Trans #1936	US\$ 45.60	+Info

CERVICAL TRANSPARENT

Code	Name	Price	+
116844	CERAMAGE Cervical Trans A #1937	US\$ 45.60	+Info
116845	CERAMAGE Cervical Trans B #1938	US\$ 45.60	+Info
116846	CERAMAGE Cervical Trans R #1939	US\$ 45.60	+Info

OPAQUE MODIFIER

OPAQUE MODIFIER - SHOFU

Using the experience and knowledge gained over many years, SHOFU has developed Ceramage.

Ceramage contains 73% by weight of micro-fine ceramic particle, which are embedded in an organic polymer matrix and combines the advantages of ceramic and composite to give a light-curing "esthetic" composite for use in any indication. This composition increases the homogenous structure resulting in optimised physical properties and enhanced esthetics.

Ceramage combines natural color reproduction with extraordinary strength and elasticity and is recommended for conventional cases as well as implant restorations.

Indications:

- * Fully veneered crowns and bridges
- * Fully veneered telescopic and implant-supported restorations
- * Precision attachments
- * Occlusal facings
- * Inlays / Onlays
- * Veneers
- * Long-Term temporary restorations
- * Modifying conventional acrylic teeth

Benefits:

- * Excellent flexural and compressive strength
- * Excellent abrasion resistance of occlusal facings
- * Highest aesthetics with a light refractive index similar to natural teeth
- * Extensive range of translucent effect colors
- * High color stability
- * Perfect viscosity and modelling properties



- * Superb polishability
- * Outstanding plaque resistance
- * Excellent biocompatibility

Code	Name	Price	+
110623	CERAMAGE MO #1841 - 2.6ml / 4.6g Syringe	US\$ 34.64	+Info
110624	CERAMAGE BLUE GRAY Opaque #1842 - 2.6ml / 4.6g Syringe	US\$ 34.64	+Info
110625	CERAMAGE GRAY Opaque #1843 - 2.6ml / 4.6g Syringe	US\$ 34.64	+Info
110626	CERAMAGE WHITE Opaque #1844 - 2.6ml / 4.6g Syringe Ceramage -	US\$ 34.64	+Info
110627	CERAMAGE Gum Opaque #1845 - PINK - 2.6ml / 4.6g Syringe Ceramage silicate -	US\$ 34.64	+Info
116837	CERAMAGE Opaque 2ml #1848	US\$ 34.64	+Info
116838	CERAMAGE Opaque 2ml #1849	US\$ 34.64	+Info

Light Cure LAB Composite - GC AMERI CA

Light Cure Lab Composite - GC



GRADIA

GRADIA(R) Indirect Composite Restoration System
Product Information
 GRADIA(R) is an indirect restoration system for crown and bridge, inlays, onlays and veneers. This product is an advanced 'MFR composite'. Its internal structure consists of reinforced bonds between an organic/inorganic filler and a lightly-filled matrix resin which together impart a very high mechanical strength while being 'gentle' on opposing teeth. This makes GRADIA(R) particularly suitable for posterior high-wear, high-pressure restorations that are prone to chipping or cracking when made with porcelain. The composite is also simple to handle, biocompatible, very technique tolerant and can be easily polished to a beautiful natural gloss. When you want to use a composite restorative and the quality of your work demands using only the best, there is no finer system than GRADIA(R).

Years of experience have shown that with GRADIA(R) you can create beautiful, life-like restorations without the 'paleness' of traditional composites as well as restorations with all the toughness of natural teeth. Another benefit with GRADIA(R) is that GRADIA(R) gum shades, our gingiva build-up composite, is an integral part of the GRADIA(R) system. Used separately, each product gives superb results. But if a patient requires teeth and gum restorations where their 'white' and 'red' esthetics need to be finely balanced, then the complete system is ideal. The shades and accessories of both products are fully harmonized so you can precisely color match your restorations to the patient and to each other, regardless of age or ethnicity.

Accessories for Gradia

Accessories for Gr

Code	Name	Price	+
422045	Gradia Separator, 5 ml Bottle	US\$ 38.56	+Info
422046	Gradia Air Barrier, 10 ml Bottle	US\$ 54.93	+Info
422125	Gradia Die Hardener, 5 ml Bottle	US\$ 38.56	+Info
475212	Metal Primer Z 5ml	US\$ 114.24	+Info

Cervical

Cervical

Code	Name	Price	
422929	Gradia Cervical Translucent Light C	US\$ 54.93	+Info
422950	Gradia Cervical Translucent Light C	US\$ 54.93	+Info

Dentine

Dentine

Code	Name	Price	+
422936	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
422937	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
422956	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
422965	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
422967	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
422975	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
422985	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
422996	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info
423002	Gradia Dentin Light Cured Micro Cer	US\$ 54.93	+Info

Enamel

Enamel

Code	Name	Price	+
421374	Gradia Enamel Light Cured Micro Cer	US\$ 54.93	+Info
422923	Gradia Enamel Light Cured Micro Cer	US\$ 54.93	+Info
422931	Gradia Enamel Intensive Light Cured	US\$ 54.93	+Info
422952	Gradia Halo Enamel Light Cured Mic	US\$ 54.93	+Info
422954	Gradia Enamel Light Cured Micro Cer	US\$ 54.93	+Info
422968	Gradia Pearl Enamel Light Cured Mic	US\$ 54.93	+Info
422988	Gradia Enamel Intensive Light Cured	US\$ 54.93	+Info
422997	Gradia Enamel Light Cured Micro Cer	US\$ 54.93	+Info

Intensive Color

Intensive Color

Code	Name	Price	+
422916	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422922	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422925	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422930	Gradia Direct Gradia Intensive Colo	US\$ 44.70	+Info
422935	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422944	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422945	Gradia Intensive Color Light Cured	US\$ 44.70	+Info
422964	Gradia Direct Gradia Intensive Colo	US\$ 46.35	+Info
422969	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422977	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422993	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
422995	Gradia Intensive Color Light Cured	US\$ 46.35	+Info
423000	Gradia Direct Gradia Intensive Colo	US\$ 46.35	+Info
423001	Gradia Intensive Color Light Cured	US\$ 46.35	+Info

Kits /Set

Kits /Set

Code	Name	Price	+
421926	Gradia Intensive Color Set Contains	US\$ 447.02	+Info

Mamelon

Mamelon

Code	Name	Price	+
422921	Gradia Mamelon Satin Light Cured Mi	US\$ 46.35	+Info
422938	Gradia Mamelon Satin Light Cured Mi	US\$ 46.35	+Info
422989	Gradia Mamelon Satin Light Cured Mi	US\$ 46.35	+Info

Opaque

Opaque

Gradia Indirect

Gradia Indirect

Code	Name	Price	+
421375	Gradia Opaque Lt Cured MC OB1	US\$ 50.99	+Info
421995	Gradia Opaque Lt Cured MC OA2	US\$ 46.35	+Info
422060	Gradia Opaque Lt Cured MC OA1	US\$ 50.99	+Info
422120	Gradia Opaque Lt Cured MC OB2	US\$ 50.99	+Info
422126	Gradia Foundation Opaque (FO) Light	US\$ 46.35	+Info
422915	Gradia Opaque Light Cured Micro Cer	US\$ 46.35	+Info
422918	Gradia Opaque Light Cured Micro Cer	US\$ 46.35	+Info
422940	Gradia Opaque Lt Cured OA3.5	US\$ 46.35	+Info
422947	Gradia Opaque Light Cured Micro Cer	US\$ 46.35	+Info
422963	Gradia Opaque Light Cured Micro Cer	US\$ 46.35	+Info
422981	Gradia Opaque Lt Cured OC3	US\$ 50.99	+Info
422982	Gradia Opaque Light Cured Micro Cer	US\$ 46.35	+Info
422990	Gradia Opaque Light Cured Micro Cer	US\$ 46.35	+Info
422998	Gradia Opaque Light Cured Micro Cer	US\$ 46.35	+Info

Opaque Dentin

Opaque Dentin

Code	Name	Price	+
422086	Gradia Opaque Dentin Light Cured Mi	US\$ 54.93	+Info
422914	Gradia Opaque Dentin Intensive Ligh	US\$ 54.93	+Info
422919	Gradia Opaque Dentin Intensive Ligh	US\$ 54.93	+Info
422924	Gradia Opaque Dentin Light Cured Mi	US\$ 54.93	+Info
422933	Gradia Opaque Dentin Light Cured Mi	US\$ 54.93	+Info
422934	Gradia Opaque Dentin Lt Cured Micro	US\$ 54.93	+Info
422951	Gradia Opaque Dentin Intensive Ligh	US\$ 54.93	+Info
422955	Gradia Opaque Dentin Light Cured Mi	US\$ 54.93	+Info
422957	Gradia Opaque Dentin Light Cured Mi	US\$ 54.93	+Info
422959	Gradia Opaque Dentin Intensive Ligh	US\$ 52.97	+Info
422979	Gradia Opaque Dentin Intensive Ligh	US\$ 54.93	+Info
422980	Gradia Opaque Dentin Light Cured Mi	US\$ 54.93	+Info
422994	Gradia Opaque Dentin Light Cured Mi	US\$ 54.93	+Info
422976	Gradia Direct Gradia Opaqus Dentin	US\$ 54.93	+Info

Shade Guide

Shade Guide

Code	Name	Price	+
422118	Gradia Air Shade Guide Kit	US\$ 26.10	+Info
422131	Gradia Shade Guide Plate	US\$ 26.10	+Info

Shoulder Dentin

Shoulder Dentin

Code	Name	Price	+
422917	Gradia Shoulder Dentin Light Cured	US\$ 54.93	+Info
422928	Gradia Shoulder Dentin Light Cured	US\$ 54.93	+Info
422958	Gradia Shoulder Dentin Light Cured	US\$ 54.93	+Info
422986	Gradia Shoulder Dentin Light Cured	US\$ 54.93	+Info

Translucent

Translucent

Code	Name	Price	+
422966	Gradia Translucent Light Cured Micr	US\$ 54.93	+Info
422971	Gradia Translucent Light Cured Micr	US\$ 54.93	+Info
422973	Gradia Translucent Light Cured Micr	US\$ 54.93	+Info

GRADIA GUM

Gradia Gum

Gum Body shades

Gum

Gradia Gum

Gradia Gum

Code	Name	Price	+
421784	Gradia Gum Light Cured Gingival Sha	US\$ 70.51	+Info
422001	Gradia Gum Shade Gum G21, 2.9 ml Sy	US\$ 55.70	+Info
422073	Gradia Gum Fiber Shade GF7, Light-C	US\$ 26.48	+Info

Other Models

Other Models

Code	Name	Price	+
475215	GC Gradia Gum Body Shade G23 2.9 ML	US\$ 45.76	+Info

Gum Modifier

Code	Name	Price	+
421732	Gradia Gum Modifier Shade GM33, 2.4	US\$ 47.14	+Info
421802	Gradia Gum Modifier Shade GM35, 2.4	US\$ 47.14	+Info
421803	Gradia Gum Modifier Shade GM36, 2.4	US\$ 47.14	+Info
422000	Gradia Direct Gradia Gum Modifier R	US\$ 47.14	+Info
422025	Gradia Gum Modifier Shade G24, 2.9	US\$ 55.70	+Info
422064	Gradia Gum Modifier Shade GM31, 2.4	US\$ 45.45	+Info

422065	Gradia Gum Modifier Shade GM32, 2.4	US\$ 47.14	+Info
422066	Gradia Gum Modifier Shade GM30, 2.4	US\$ 47.14	+Info

Gum Opaque

Gum Opaque

Code	Name	Price	+
421804	Gradia Gum Opaque GO12 Shade, 2.4 m	US\$ 47.14	+Info
422002	Gradia Gum Opaque Shade GO22, 2.9 m	US\$ 55.70	+Info
422057	Gradia Gum Opaque Shade GO20, 2.9 m	US\$ 55.70	+Info
422067	Gradia Gum Opaque GO11 2.4 ml	US\$ 46.67	+Info
422069	Gradia Gum Opaque Shade GO13, 2.4 m	US\$ 49.29	+Info
422124	Gradia Gum Light Cured Gingival Sha	US\$ 47.14	+Info

Kits /Set

Kits /Set

Code	Name	Price	+
421311	Gradia Gum Gingival Starter Kit	US\$ 456.31	+Info

GRADIA PLUS

Kits / Set

Kits / Set

Code	Name	Price	+
475225	GRADIA Plus Composite Kit - Complet	US\$ 311.00	+Info
475226	GRADIA PLUS Layer Set	US\$ 1,170.99	+Info
475227	GRADIA PLUS Paint Set	US\$ 797.02	+Info
475228	GRADIA PLUS Gum Shades Set	US\$ 881.16	+Info
475229	GC GRADIA PLUS One Body Set	US\$ 739.37	+Info
475230	GRADIA PLUS Accessory Set	US\$ 722.99	+Info
475300	GRADIA PLUS Layer Pro Set	US\$ 1,196.68	+Info

LB Light Body

LB Light Body

Code	Name	Price	+
475268	GRADIA PLUS Paste LB-Blue, 2ml 1pc	US\$ 62.32	+Info
475279	GRADIA PLUS Paste GLB-CL, 2ml 1pc	US\$ 68.55	+Info
475278	GRADIA PLUS Paste GLB-3, 2ml 1pc	US\$ 62.32	+Info
475277	GRADIA PLUS Paste GLB, 2ml 1pc	US\$ 62.32	+Info
475276	GRADIA PLUS ONE BODY LB-W, 2.0ml 1p	US\$ 62.32	+Info
475275	GRADIA PLUS ONE BODY LB-D, 2.0ml 1p	US\$ 62.32	+Info
475274	GRADIA PLUS ONE BODY LB-C, 2.0ml 1p	US\$ 62.32	+Info
475273	GRADIA PLUS ONE BODY LB-B, 2.0ml 1p	US\$ 62.32	+Info
475272	GRADIA PLUS ONE BODY LB-A, 2.0ml 1p	US\$ 62.32	+Info
475271	GRADIA PLUS Paste LB-Inlay TD, 2ml	US\$ 62.32	+Info
475270	GRADIA PLUS Paste LB-Inlay E, 2ml 1	US\$ 62.32	+Info
475269	GRADIA PLUS Paste LB-Milky, 2ml 1pc	US\$ 62.32	+Info
475216	GRADIA PLUS Paste GLB-1 2ml	US\$ 62.32	+Info
475267	GRADIA PLUS Paste LB-Grey, 2ml 1pc	US\$ 62.32	+Info

475266	GRADIA PLUS Paste LB-Red, 2ml 1pc	US\$ 62.32	+Info
475265	GRADIA PLUS Paste LB-Yellow, 2ml 1p	US\$ 62.32	+Info
475264	GRADIA PLUS Paste LB-Orange, 2ml 1p	US\$ 62.32	+Info
475263	GRADIA PLUS Paste LB-Base Opal, 2ml	US\$ 62.32	+Info
475262	GRADIA PLUS Paste LB-Base ODW, 2ml	US\$ 62.32	+Info
475261	GRADIA PLUS Paste LB-Base OD, 2ml 1	US\$ 62.32	+Info
475260	GRADIA PLUS Paste LB-Base DW, 2ml 1	US\$ 62.32	+Info
475259	GRADIA PLUS Paste LB-Base D, 2ml 1p	US\$ 62.32	+Info
475258	GRADIA PLUS Paste LB-Base CLF 2ml,	US\$ 62.32	+Info
475257	GRADIA PLUS Paste LB-Base E, 2ml 1p	US\$ 62.32	+Info

HB Heavy Body

HB Heavy Body

Code	Name	Price	+
475246	GRADIA PLUS Paste HB-ED, 3.3ml, 1pc	US\$ 60.38	+Info
475256	GRADIA PLUS Paste GHB-CL, 3.3ml, 1p	US\$ 60.38	+Info
475255	GRADIA PLUS Paste GHB-3, 3.3ml, 1pc	US\$ 60.38	+Info
475254	GRADIA PLUS Paste HB-DW, 3.3ml, 1pc	US\$ 60.38	+Info
475253	GRADIA PLUS Paste HB-ODW, 3.3ml, 1p	US\$ 60.38	+Info
475252	GRADIA PLUS Paste HB-ODD, 3.3ml, 1p	US\$ 60.38	+Info
475251	GRADIA PLUS Paste HB-ODC, 3.3ml, 1p	US\$ 60.38	+Info
475250	GRADIA PLUS Paste HB-ODB, 3.3ml, 1p	US\$ 60.38	+Info
475249	GRADIA PLUS Paste HB-ODA, 3.3ml, 1p	US\$ 60.38	+Info
475248	GRADIA PLUS Paste HB-PE, 3.3ml, 1pc	US\$ 60.38	+Info
475247	GRADIA PLUS Paste HB-CLF, 3.3ml, 1p	US\$ 60.38	+Info
475213	GRADIA PLUS Paste GHB-1, 3.3ml, 1pc	US\$ 60.38	+Info
475245	GRADIA PLUS Paste HB-EL, 3.3ml, 1pc	US\$ 60.38	+Info
475244	GRADIA PLUS Paste HB-DD2, 3.3ml, 1p	US\$ 60.38	+Info
475243	GRADIA PLUS Paste HB-DC3, 3.3ml, 1p	US\$ 60.38	+Info
475242	GRADIA PLUS Paste HB-DB3, 3.3ml, 1p	US\$ 60.38	+Info
475241	GRADIA PLUS Paste HB-DB1, 3.3ml, 1p	US\$ 60.38	+Info
475240	GRADIA PLUS Paste HB-DA3.5, 3.3ml,	US\$ 60.38	+Info
475239	GRADIA PLUS Paste HB-DA3, 3.3ml, 1p	US\$ 60.38	+Info
475238	GRADIA PLUS Paste HB-DA2, 3.3ml, 1p	US\$ 60.38	+Info
475237	GRADIA PLUS Paste HB-DA1, 3.3ml, 1p	US\$ 60.38	+Info
475214	GRADIA PLUS Paste GHB-2, 3.3ml	US\$ 60.38	+Info

LP Lustre Paste

LP Lustre Paste

Code	Name	Price	+
475280	GRADIA PLUS LP-A, 0.8ml, 1pc	US\$ 62.32	+Info
475281	GRADIA PLUS LP-B, 0.8ml, 1pc	US\$ 62.32	+Info
475282	GRADIA PLUS LP-C, 0.8ml, 1pc	US\$ 62.32	+Info
475283	GRADIA PLUS LP-D, 0.8ml, 1pc	US\$ 62.32	+Info
475284	GRADIA PLUS LP-CLF, 0.8ml, 1pc	US\$ 62.32	+Info
475285	GRADIA PLUS LP-Blue, 0.8ml, 1pc	US\$ 62.32	+Info
475286	GRADIA PLUS LP-Grey, 0.8ml, 1pc	US\$ 62.32	+Info
475287	GRADIA PLUS LP-Cream, 0.8ml, 1pc	US\$ 62.32	+Info
475288	GRADIA PLUS LP-Lavender, 0.8ml, 1pc	US\$ 62.32	+Info

475289	GC GRADIA PLUS LP-CL, 0.8ml, 1pc	US\$ 62.32	+Info
475290	GRADIA PLUS GLP-Bright Red, 0.8ml,	US\$ 62.32	+Info
475291	GRADIA PLUS GLP-Violet, 0.8ml, 1pc	US\$ 62.32	+Info
475292	GRADIA PLUS LP Diluting Liquid, 3ml	US\$ 77.52	+Info

Opaque

Opaque

Code	Name	Price	+
475218	GRADIA PLUS Opaque GO-1 2.0ml	US\$ 60.00	+Info
475219	GRADIA PLUS Opaque GO-2, 2.0ml, 1 p	US\$ 60.00	+Info
475231	GRADIA PLUS Opaque O-Base 2.0ml	US\$ 60.00	+Info
475232	GRADIA PLUS Opaque OA, 2.0ml, 1 pc.	US\$ 60.00	+Info
475233	GRADIA PLUS Opaque OB, 2.0ml, 1 pc.	US\$ 60.00	+Info
475234	GRADIA PLUS Opaque OC, 2.0ml, 1 pc.	US\$ 60.00	+Info
475235	GRADIA PLUS Opaque OD 2.0ml, 1 pc.	US\$ 60.00	+Info
475236	GRADIA PLUS Opaque OA1, 2.0ml, 1 pc	US\$ 60.00	+Info

Accessories

Accessories

Dispensing Tips

Dispensing Tips

Code	Name	Price	+
475296	GRADIA PLUS DTip PTWide LP Cover 5p	US\$ 9.74	+Info
475297	GRADIA PLUS DTip Ntype LP Cover 5pc	US\$ 8.96	+Info
475298	GRADIA PLUS DTip LP Cover 5pc	US\$ 21.42	+Info
475299	GRADIA PLUS Mixotip 10pcs	US\$ 56.08	+Info

Liquids

Liquids

Code	Name	Price	+
475211	GRADIA PLUS AIR BARRIER, 10ml, 1 pc	US\$ 35.83	+Info
475293	GRADIA PLUS Modelling Liquid 3ml, 1	US\$ 37.00	+Info
475294	GRADIA PLUS SEPARATOR. 5ml, 1 pc.	US\$ 64.28	+Info
475295	GRADIA PLUS DIE-HARDNER, 5ml, 1 pc.	US\$ 39.35	+Info

Rubber Tips

Rubber Tips

Code	Name	Price	+
475221	GRADIA Plus Comp. Rubber Tip #1 Chi	US\$ 95.04	+Info
475222	GRADIA Plus Comp. Rubber Tip #2 Con	US\$ 95.04	+Info

Spatulas

Spatulas

Code	Name	Price	+
475223	GRADIA Plus Composite Spatula #2	US\$ 95.04	+Info
475224	GRADIA Plus Composite Spatula #3	US\$ 95.04	+Info

LITE ART COLOR STAIN | CERAMAGE | SHOFU

LITE ART COLOR CERAMAGE - SHOFU

Lite Art light-curing coloring Stain pastes for Ceramage are designed to recreate the internal effects and characteristics of natural teeth.

Works in composites, and all light cure compounds. In order to be used in regular **Acrylic Dentures, you need to use CERARESIN as bonding agent**



Characteristics:

- * Excellent curing characteristics due to newly developed multifunctional monomer and photo initiator.
- * Low viscosity - easy brush technique.
- * Fluorescence of natural teeth.
- * Wide range of expressive colors for individual characteristics.
- * Supreme handling properties.

Authorized SHOFU Dealer.

Colors

Colors

Code	Name	Price	+
116857	Lite At Col Paste White	US\$ 28.65	+Info
116858	Lite Art Col Paste Black	US\$ 28.08	+Info
116859	Lite Art Col Paste red	US\$ 28.65	+Info
116860	Lite Art Col Paste A Shade	US\$ 28.65	+Info
116861	Lite Art Col Paste B Shade	US\$ 31.35	+Info
116862	Lite Art Col Paste Khaki	US\$ 28.65	+Info
116863	Lite Art Col Paste Org	US\$ 28.65	+Info
116864	Lite Art Col Paste Blue Gray	US\$ 25.50	+Info
116865	Lite Art Col Paste Violet	US\$ 28.65	+Info
116866	Lite Art Col Paste Org Brwn	US\$ 28.65	+Info
116867	Lite Art Col Paste Dk red Brwn	US\$ 28.65	+Info
116868	Lite Art Col Paste Blasck Brwn	US\$ 28.65	+Info
116869	Lite Art Col Paste Cyan	US\$ 28.65	+Info
116870	Lite Art Col Paste Magenta	US\$ 28.65	+Info
116871	Lite Art Col Paste Yellow	US\$ 28.65	+Info

Kits / Set

Kits / Set

Code	Name	Price	+
116856	Lite Art Full Set	US\$ 383.78	+Info

Liquid

Liquid

Code	Name	Price	+
116872	Lite Art Clear Liq	US\$ 28.65	+Info

PARAFIL LAB Light Cure Lab Composite | PRIM E DENTAL

PARAFIL LAB - PRIME DENTAL

Features:

Parafil LAB is a light-cured indirect restoration system for crown and bridge, inlays, onlays and veneers. It has unsurpassed durability, natural opalescence and excellent lifelike esthetics. The unique chemistry of Parafil LAB couples nano particle/ZIRCONIUM filler with a urethane dimethacrylate matrix to produce a superior composite with unsurpassed high strength, wear resistance and superior polishability. Parafil LAB material is excellently polishable and allows highly esthetic results to be achieved. The material structure of Parafil LAB contributes to the excellent grinding and polishing properties. Parafil LAB offers Parafil



Gingiva for gingiva build-up. Used separately, each product gives superb results. If a patient requires teeth and gum restorations where their ?white? and ?red? esthetics need to be finely balanced, then the complete system is ideal. The shades and accessories of both products are fully harmonized so you can precisely color match your restoration to the patient and to each other, regardless of age or ethnicity.

Ease Of Cure

- * Does not require an oven for a final cure
- * Cure composite with any LED light

Superior Physical

- * Properties
- * High wear resistance
- * Unsurpassed strength
- * Can be used for a wide range of applications

Optimized Filler Particle Size

- * Paste has a consistency that makes it easier to mold and shape when compared to porcelain and ceramic.
- * Holds shape without slumping
- * Tooth like light reflection
- * No problem polishing

COMPOSITE

Code	Name	Price	+
115722	PARAFIL Lab Composite 4.5g Syringe A1 #9999-70A1	US\$ 16.00	+Info
115723	Parafil LAB A2 - 4.5g Syr - Zirconium Composite #9999-70A2	US\$ 16.00	+Info
115724	Parafil LAB A3 - 4.5g Syr - Zirconium Composite #9999-70A3	US\$ 16.00	+Info
123893	Parafil LAB A3.5 - 4.5gr Syr - Zirconium Composite - #9999-70A3.5	US\$ 16.00	+Info
115725	Parafil LAB A4 - 4.5g Syr - Zirconium Composite #9999-70A4	US\$ 16.00	+Info
123465	Parafil LAB B3 - 4.5gr Syr - Zirconium Composite - #9999-70B3	US\$ 16.00	+Info
115727	Parafil LAB Natural Translucent - 4.5g Syr - Zirconium Composite - # 9999-70NT	US\$ 16.00	+Info
123894	Parafil LAB OA3 - 4.5gr Syr - Zirconium Composite - #9999-70OA3	US\$ 16.00	+Info

123897

Parafil LAB Opal UD - Opal Universal dentine - Zirconium Composite - #9999-71UD

US\$ 35.20

[+Info](#)

GINGIVAL

Body

Body

Code	Name	Price	+
123480	Parafil LAB Ging 4.5gr G21	US\$ 16.00	+Info
420269	Parafil LAB Ging 4.5gr G23	US\$ 16.00	+Info
474943	Parafil LAB Gingiva G22 4.5gr	US\$ 11.50	+Info

Modifiers

Modifiers

Code	Name	Price	+
420266	Parafil LAB Ging Opaq 4.5gr GM21	US\$ 12.80	+Info
420267	Parafil LAB Gingiva Modifier GM22	US\$ 16.00	+Info
420268	Parafil LAB Ging Modifier GM23 2gr	US\$ 16.00	+Info
420270	Parafil LAB Ging Opaq 2gr GOM2	US\$ 16.00	+Info

Opaque

Opaque

Code	Name	Price	+
123479	Parafil LAB Ging Opaq 2gr GO	US\$ 18.00	+Info

INCISAL

Code	Name	Price	+
123474	Parafil LAB IU Incisal 4.5gr Syr - Zirconium Composite - #9999-70IU	US\$ 16.00	+Info

LIQUIDS

Code	Name	Price	+
123895	Parafil LAB Modeler 7ml bottle - Zirconium Composite - #9999-7260	US\$ 16.00	+Info
123896	Parafil LAB Composite Primer 7ml bottle - Zirconium Composite - #9999-7230	US\$ 16.00	+Info
123898	Parafil LAB Separator 3gr Syringe - Zirconium Composite - #9999-7240	US\$ 14.40	+Info
123900	Parafil LAB Metal Primer 7ml bottle - Zirconium Composite - #9999-7220	US\$ 16.00	+Info

OPALESCENT

Code	Name	Price	+
123476	Parafil LAB OTB Opalescent 4.5gr Syr - Zirconium Composite - #9999-70OTB	US\$ 16.00	+Info

OPAQUER

Code	Name	Price	+
124490	PARAFIL LAB Foundation Universal Opaque - Syringe #9999-70F-UO	US\$ 16.00	+Info

Other Models

Other Models

Code	Name	Price	+
------	------	-------	---

124347	Parafil LAB compo KIT 10-Syrin	US\$ 400.00	+Info
124348	Parafil LAB compo KIT 4-Syrin	US\$ 126.40	+Info

STAIN

Code	Name	Price	+
124340	Parafil LAB Stain Syringe BEIGE1.6g - Zirconium Composite - #9999-71BE	US\$ 16.00	+Info
124341	Parafil LAB Stain Syringe BLUE- Zirconium Composite - #9999-71BL	US\$ 16.00	+Info
124342	Parafil LAB Stain Syringe BROWN1.6g - Zirconium Composite - #9999-71BR	US\$ 16.00	+Info
124343	Parafil LAB Stain Syringe YELLOW 1.6g - Zirconium Composite - #9999-71YE	US\$ 16.00	+Info
124344	Parafil LAB Stain Syringe ORANGE 1.6g - Zirconium Composite - #9999-71OR	US\$ 16.00	+Info
124345	Parafil LAB Stain Syringe GRAY 1.6g - Zirconium Composite - #9999-71GR	US\$ 16.00	+Info
124346	Parafil LAB Stain Syringe PINK 1.6g - Zirconium Composite - #9999-71PI	US\$ 16.00	+Info
123899	Parafil LAB Stain Syringe WHITE 1.6g - Zirconium Composite - #9999-71WH	US\$ 16.00	+Info

Inlay Resins

Resin pattern

DURALAY RESIN PATTERN- RED - RELIANCE

DURALAY RESIN PATTERN - RED - RELIANCE

Forms a dense, accurate, non-shrinking pattern that can be trimmed, handled, and invested without fear of breakage, shipping, distortion or temperature changes.

Burns out the same temperature as wax, leaving no residue.



Kits /Set

Kits /Set

Code	Name	Price	+
071-2244	Duralay red Standard Pack 2oz	US\$ 30.53	+Info
071-2249A	Duralay Laboratory KIT 8oz	US\$ 88.52	+Info
071-2249D	Duralay Student Kit With 1oz.	US\$ 24.83	+Info

Powder refills

Powder refills

Code	Name	Price	+
071-2245	Powder Only -2oz Shade red	US\$ 16.73	+Info
071-2245-1	Powder Only -1oz Shade red	US\$ 18.98	+Info
071-2249B	Duralay red Powder Only - 8oz	US\$ 56.34	+Info

Resin pattern - DURACRYL

Resin pattern - DURACRYL

Code	Name	Price	+
100427	Duracryl Flow Kit Pattern 55ml/60gr	US\$ 11.36	+Info

Resin pattern - GC AMERICA

Resin pattern - GC AMERICA



Pattern Resin

Pattern Resin



Kits /Set

Kits /Set

Code	Name	Price	+
233-396_A	Pattern Resin LS Kit	US\$ 107.80	+Info

Liquid refills

Liquid refills

Code	Name	Price	+
113792	Pattern Resin Liquid 100ml	US\$ 49.47	+Info
421509	Pattern Resin LS 250 g Liquid Acryl	US\$ 81.12	+Info
421683	Pattern Resin LS 100 g Liquid Acryl	US\$ 49.47	+Info

Powder refills

Powder refills

Code	Name	Price	+
113791	Pattern Resin Powder 100gr	US\$ 58.05	+Info
421350	Pattern Resin LS 100 g Powder Acryl	US\$ 58.05	+Info
421476	Pattern Resin LS 1 kg Powder	US\$ 343.70	+Info

Unifast Trad

Unifast Trad

Liquid only

Liquid only

Code	Name	Price	+
421341	UNIFAST Trad Liquid 104ml	US\$ 34.67	+Info

Powder refills

Powder refills

Code	Name	Price	+
421343	Unifast Trad Pwdr only Ivory 100gr	US\$ 62.15	+Info
421688	Unifast Trad Pwdr only Ivory 250gr	US\$ 116.45	+Info
421879	Unifast Trad Powder only, Ivory, Me	US\$ 411.67	+Info

Catalog Index / Indice

Light Cure Lab Composite	1
Resin pattern	23