

Model Trimmers

Arch Model Trimmers



Arch Model Trimmers - RENFERT

Arch Model Trimmers - RENFERT



Code	Name	Price	+
23-1805-0000	Millo Pro Arch Trim w/bur 220v	US\$ 1,151.19	+Info
23-1805-1000	Millo Pro Arch Trim w/bur 120v	US\$ 1,113.88	+Info
417579	Millo Arch Trimmer 120v	US\$ 861.18	+Info

Arch Model Trimmers - BESQUAL

Arch Model Trimmers - BESQUAL

Code	Name	Price	+
100469	Arch Model Trimmer 110V	US\$ 178.50	+Info
116071	Arch Model Trimmer 220V	US\$ 189.52	+Info

Arch Model Trimmers - INTEGRAL SYS

Arch Model Trimmers - INTEGRAL SYS

Code	Name	Price	+
115054	ROUTER ISI-898	US\$ 699.46	+Info

Arch Model Trimmers - MESTRA

Arch Model Trimmers - MESTRA

Code	Name	Price	+
420201	Interior trimmer Trim w/bur 220v	US\$ 1,062.45	+Info
420207	Interior trimmer Trim w/bur 110v	US\$ 1,062.45	+Info

Arch Model Trimmers - RAY FOSTER

Arch Model Trimmers - RAY FOSTER

Code	Name	Price	+
470895	Model Arch Trimmer MR20	US\$ 1,034.48	+Info
470896	Model Arch Trimmer w/ CDC1 Dust	US\$ 2,000.62	+Info

Model Trimmers

Disc Model Trimmers



Double Disc Model Trimmers

Double Disc



Ortho Wet Double Timmers

Ortho Wet Model Timmers

Code	Name	Price	+
112855	Dual Disc Ortho Trimmer 220v	US\$ 4,830.05	+Info
112857	Dual Disc Ortho Trimmer 110v	US\$ 4,303.98	+Info

Wet Double Model trimmers

Wet Double Model trimmers

Code	Name	Price	+
118758	Twin Mod Trim. 10in 110v	US\$ 675.00	+Info

Dry Model Trimmers

Dry Model Trimmers

BUFFALO

BUFFALO

Code	Name	Price	+
008-61650	DryTrim Premium Trimmer 1/3HP 110v	US\$ 1,219.40	+Info
008-61660	DryTrim Premium Trimmer 1/3Hp 220v	US\$ 1,256.92	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
116011	Dry Mod Trim 10in 110v w/DC	US\$ 2,203.50	+Info
116012	Dry Mod Trim 10 in 110v	US\$ 1,232.80	+Info

Ortho Wet Model Timmers

Ortho Wet Model Timmers

Carborundum disc

Carborundum disc

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
101585	Ortho Mod Trim 12in 110v Full	US\$ 2,233.40	+Info
112772	Ortho Mod Trim 12in 220v Std	US\$ 2,233.40	+Info
112846	Ortho Mod Trim 12in 110v Std	US\$ 1,989.90	+Info

WHIP MIX

WHIP MIX

Code	Name	Price	+
120216	Dlx Ortho Trim 12in 115v White	US\$ 1,907.70	+Info

Diamond disc

Diamond disc

Code	Name	Price	+
470888	12 Ortho Trimmer w Diamond + Access	US\$ 3,047.20	+Info
470925	12in Ortho Trimmer w/ Diamond Wheel	US\$ 2,737.80	+Info

Wet & Dry Model Trimmer

Wet & Dry Model Trimmer

Carborundum disc

Carborundum disc

Code	Name	Price	+
424490	DualTrim Dry Trimmer 10in 1/2H 120v	US\$ 890.40	+Info
424491	DualTrim Dry Trimmer 10in 1/3H 120v	US\$ 917.00	+Info
425204	DualTrim Dry Trimmer 10in 1/2H 220v	US\$ 919.80	+Info
425207	DualTrim Dry Trimmer 10in 1/3H 220v	US\$ 924.00	+Info

Diamond disc

Diamond disc

Code	Name	Price	+
420203	Model Trimmer Diamond Coated 220V	US\$ 1,324.35	+Info
420209	Model Trimmer Diamond Coated 110V	US\$ 1,525.50	+Info

Wet Model trimmers

Wet Model trimmers

Adhesive discs

Adhesive discs

Code	Name	Price	+
113519	MT3 Trimmer w/Klettfix 110v WET	US\$ 1,068.44	+Info

Carborundum disc

Carborundum disc

10in Discs Model Trimmers

10in Discs

BUFFALO

BUFFALO

Code	Name	Price	+
008-61790A	Model Trimmer 10in 1/3HP 110v	US\$ 821.80	+Info
008-61792	Model Trimmer 10in 1/2HP 110v	US\$ 968.80	+Info
008-61792A	Model Trimmer 10in 1/2HP 220v	US\$ 1,041.11	+Info
008-61810	Model Trimmer 10in 1/3HP 220v	US\$ 859.60	+Info
425208	10" diameter Carboru	US\$ 1,109.52	+Info

DENTAL FARM

DENTAL FARM

Code	Name	Price	+
474535	MT-WET	US\$ 1,185.72	+Info

HANDLER

HANDLER

Code	Name	Price	+
101767	Model Trimmer 10in 1/3HP	US\$ 871.09	+Info
103344	Model Trimmer 10in 1/4HP	US\$ 823.30	+Info
115151	Model Trimmer 10in W/Sol-Valve	US\$ 953.89	+Info
115153	Mod-Trimmer 10in 1/4HP W/S/V	US\$ 942.79	+Info

MESTRA

MESTRA

Code	Name	Price	+
419522	Model trimmer MULHACEN	US\$ 761.40	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
101583	MT10 Model Trimmer 10in 1/3HP	US\$ 1,121.58	+Info
470884	MT10A Trimmer w/ water & shield	US\$ 1,278.36	+Info

WHIP MIX

WHIP MIX

Code	Name	Price	+
120209	Model Trimmer 10in 115v White	US\$ 1,265.11	+Info

12in Discs Model Trimmers

12in Discs

BUFFALO

BUFFALO

Code	Name	Price	+
008-61820	Model Trimmer 12in 1/2HP 110v	US\$ 1,016.40	+Info
008-61825	Model Trimmer 12in 1/3HP 110v	US\$ 835.80	+Info
008-61830	Model Trimmer 12in 1/2HP 220v	US\$ 1,073.80	+Info
424782	12" Carborundum Trimmer Wheel	US\$ 39.51	+Info
425202	12" diameter Carboru	US\$ 1,130.96	+Info

HANDLER

HANDLER

Code	Name	Price	+
101765	Model Trimmer 12in 3/4HP	US\$ 988.40	+Info
101766	Model Trimmer 12in 1/2HP	US\$ 988.40	+Info
101769	Model Trimmer 12in 1/3HP	US\$ 979.88	+Info
103354	Model Trimmer 12in 3/4HP	US\$ 1,210.00	+Info
115164	Model Trimmer 12in 1/2HP W/S/V	US\$ 1,103.48	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
101584	MT12 Model Trimmer 12in 1/2HP	US\$ 1,522.24	+Info
470886	MT12 Trimmer w/ Water & Shield	US\$ 1,838.48	+Info

WHIP MIX

WHIP MIX

Code	Name	Price	+
120211	Model Trimmer 12in 115v White	US\$ 1,293.16	+Info
120214	Dlx Model Trim 12in 115v White	US\$ 1,567.03	+Info
120217	Model Trimmer 12in 115v 3/4 Hp	US\$ 2,102.84	+Info

Diamond disc

Diamond disc

10in Discs Model Trimmers

10in Discs

HANDLER

HANDLER

Code	Name	Price	+
115152	Model Trimmer 10in W/Diam-Wh	US\$ 1,246.34	+Info
115156	Mod Trimmer 10in 1/4HP W/D/Wh	US\$ 1,210.00	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
470885	MT10 Model Trimmer w/ Diamond wheel	US\$ 1,867.96	+Info
470887	MT12 Model Trimmer w/ Diamond wheel	US\$ 2,256.80	+Info
470923	MT10 Trimmer w/ Diamond + Accessory	US\$ 2,131.94	+Info
470924	MT12 Trimmer w/ Diamond + Accessory	US\$ 2,562.30	+Info

12in Discs Model Trimmers

12in Discs

HANDLER

HANDLER

Code	Name	Price	+
115158	Model Trimmer 12in 1/2HP W/D/W	US\$ 1,512.80	+Info
115163	Model Trimmer 12in 1/3HP W/D/W	US\$ 1,355.34	+Info

WHIP MIX

WHIP MIX

Code	Name	Price	+
120215	Diamond Dlx Trim 12in 115v 1/2	US\$ 2,496.29	+Info

Partially Diamond

Partially Diamond

Code	Name	Price	+
1808-1500	MT3 Trimmer w/Marathon 110v WET	US\$ 1,692.21	+Info

Model Trimmers

Discs for model trimmers



10in Discs

10in Discs



Adhesive Discs

Adhesive Discs

BUFFALO

BUFFALO

Code	Name	Price	+
008-61960	Mod Trimer Wheels Adhes #10	US\$ 36.54	+Info
424504	Quick-Stick Adhesive Disc 10in	US\$ 35.84	+Info
424668	Plastic Back Up Whee	US\$ 20.44	+Info

DEDECO

DEDECO

Code	Name	Price	+
031-7060	Stickies-Mod Trim Dis 10in 6pk	US\$ 89.50	+Info
031-7061	Stickies-Mod Trim Dis 10in 6pk	US\$ 51.55	+Info
031-7062	Stickies-Mod Trim Dis 10in Kt	US\$ 77.03	+Info
031-7063	Stickies-Mod Trim Dis 10in 6pk	US\$ 71.10	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
103832	Model Trim Disc Set 10in	US\$ 48.00	+Info
103833	Model Trim Replace discs 10in	US\$ 39.45	+Info
103834	Model Trim Back Wheel 10in	US\$ 14.70	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
101711	Aluminum Back Plate for 10in	US\$ 117.45	+Info
101712	Alumi Back Plate Removal Tool	US\$ 36.12	+Info

RENFERT

RENFERT

Code	Name	Price	+
120690	Klettfix For Mt2	US\$ 286.56	+Info
23-1803-1000	KLETTFIX MTplus Dia. 234m 5pk	US\$ 193.05	+Info
23-1803-1100	KLETTFIX MTplus 80gt 5pk	US\$ 93.81	+Info
23-1803-1200	KLETTFIX MTplus 120gt 5pk	US\$ 90.06	+Info

Carborundum disc

Carborundum disc

BUFFALO

BUFFALO

Code	Name	Price	+
424619	10" Carborundum Trim Wheel	US\$ 33.34	+Info

HANDLER

HANDLER

Code	Name	Price	+
103333	Carbo-dent Wheel 10in X Coarse	US\$ 53.97	+Info
103334	Carbo-dent Wheel 10in Coarse	US\$ 29.94	+Info
103335	Carbo-dent Wheel 10in Medium	US\$ 45.27	+Info

103337	Reinforced Wheel 10in Coarse	US\$ 29.94	+Info
103338	Reinforced Wheel 10in Fine	US\$ 53.97	+Info
103339	Reinforced Wheel 10in Medium	US\$ 45.27	+Info
115165	Carbo-dent Wheel 10in X Fine	US\$ 44.60	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1900620	Wheel Carb Coarse 10in 1 Hole	US\$ 30.06	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
101713	Replacement Wheel 10in Coarse	US\$ 80.52	+Info
101714	Replacement Wheel 10in XCoarse	US\$ 73.44	+Info
470892	Coated Abrasive 10in Thick plastic	US\$ 34.04	+Info
470894	Trimmer Disc 10in 3pk + mouniting	US\$ 108.23	+Info

WHIP MIX

WHIP MIX

Code	Name	Price	+
120247	#30293 Coarse Trim Wheel 10 In	US\$ 66.46	+Info
120249	#30102 Med Abrasive Wheel 10in	US\$ 66.46	+Info

Diamond disc

Diamond disc

BESQUAL

BESQUAL

Code	Name	Price	+
102213	Wheel Diam Coated 10in Medium	US\$ 262.50	+Info

BUFFALO

BUFFALO

Code	Name	Price	+
425281	10" Diamond model t	US\$ 441.54	+Info

HANDLER

HANDLER

Code	Name	Price	+
103342	Wheel Diam Coated 10in	US\$ 566.05	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1900642	Wheel Diam Coated 10in Medium	US\$ 371.52	+Info

RAY FOSTER
RAY FOSTER

Code	Name	Price	+
400666	10" Diamond Abrasive MT Wheel	US\$ 668.25	+Info

RENFERT
RENFERT

Code	Name	Price	+
120691	Infinity Trimmer Disc For Mt3	US\$ 1,497.95	+Info

Partially Diamond
Partially Diamond

RENFERT
RENFERT

Code	Name	Price	+
023-1803-2001	Marathon Disc For MTs	US\$ 799.74	+Info

WHIP MIX
WHIP MIX

Code	Name	Price	+
120245	#30430 Blue Diamond Wheel 10in	US\$ 899.77	+Info

12in Discs
12in Discs



Adhesive Discs
Adhesive Discs

BUFFALO
BUFFALO

Code	Name	Price	+
424663	Plastic Back Up Whee	US\$ 22.68	+Info
425279	12in Trimm. Disc Adhes Coarse 4pk	US\$ 35.84	+Info

DEDECO

DEDECO

Code	Name	Price	+
031-7070	Stickies-Mod Trim Dis 12in Kt	US\$ 79.70	+Info
031-7071	Stickies-Mod Trim Dis 12in 4pk	US\$ 51.22	+Info
031-7072	Stickies-Mod Trim Dis 12 Kt	US\$ 79.70	+Info
031-7073	Stickies-Mod Trim Dis 12 4pk	US\$ 59.16	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1900680	Model Trim Disc Kit 12in	US\$ 41.40	+Info
034-1900690	Model Trim Replace Disc 12in	US\$ 34.50	+Info
034-1900695	Model Trim Refill Disc 12in	US\$ 179.86	+Info
103835	Model Trim Washer 10/12in	US\$ 6.80	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
101710	Aluminum Back Plate for 12in	US\$ 125.55	+Info

Carborundum disc

Carborundum disc

HANDLER

HANDLER

Code	Name	Price	+
103336	Carbo-dent Wheel 12in Fine	US\$ 66.09	+Info
103347	Carbo-dent Wheel 12in Coarse	US\$ 34.29	+Info
103348	Carbo-dent Wheel 12in Medium	US\$ 58.15	+Info
103349	Reinforced Wheel 12in Coarse	US\$ 34.29	+Info
103350	Reinforced Wheel 12in Fine	US\$ 66.09	+Info
115169	Carbo-dent Wheel 12in X Fine	US\$ 54.65	+Info
115170	Carbo-dent Wheel 12in X Coarse	US\$ 64.77	+Info

KEYSTONE

KEYSTONE

Code	Name	Price	+
034-1900630	Wheel Carb Coarse 12in 1 hole	US\$ 34.56	+Info
034-1900635	Wheel. Carb Coarse 12in 2 Hole	US\$ 33.90	+Info
034-1900700	Model Trim Back Wheel 12in	US\$ 16.61	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
101715	Replacement Wheel 12in Coarse	US\$ 73.00	+Info
101716	Replacement Wheel 12in Fine	US\$ 77.38	+Info

101717	Replacement Wheel 12in XCoarse	US\$ 78.84	+Info
470889	Coated Abrasive 12in Coarse 3pk	US\$ 101.25	+Info
470890	Coated Abrasive 12in Coarse 5pk	US\$ 101.25	+Info
470891	Coated Abrasive 12in Thick plastic	US\$ 37.44	+Info
470893	Trimmer Disc 12in 3pk + mouniting	US\$ 112.88	+Info

WHIP MIX

WHIP MIX

Code	Name	Price	+
120248	#30097 Med Abrasive Wheel 12in	US\$ 81.76	+Info
120250	#30138 Extra Coarse Wheel 12in	US\$ 81.76	+Info
120251	#30292 Coarse Abrasive Wheel 1	US\$ 81.76	+Info
120252	#30433 Psa Disc 50grit 12in B	US\$ 634.14	+Info
120253	#30433 Psa Disc 50grit 12in B	US\$ 121.91	+Info
120254	Conversion Kit W/plst Disc 12i	US\$ 98.52	+Info

Diamond disc

Diamond disc

BESQUAL

BESQUAL

Code	Name	Price	+
114020	Diam wheel Coated 12in Medium	US\$ 271.95	+Info

BUFFALO

BUFFALO

Code	Name	Price	+
425280	12" Diamond model tr	US\$ 521.07	+Info

HANDLER

HANDLER

Code	Name	Price	+
103343	Diam wheel Coated 12in	US\$ 674.13	+Info

RAY FOSTER

RAY FOSTER

Code	Name	Price	+
400661	12" Diamond Abrasive MT Wheel	US\$ 733.60	+Info

Partially Diamond

Partially Diamond

Code	Name	Price	+
120246	#30432 Blue Diamond Wheel 12in	US\$ 1,053.64	+Info

Model Trimmers

Model Trimmer Accessories

ACCESSORIES FOR MODEL TRIMMERS - WHIP MIX

ACCESSORIES FOR MODEL TRIMMERS - WHIP MIX

Orthodontic Work Table

A precision-calibrated, easily adjustable work table offering near vertical cutting force. Comes with durable Vertisquare, Angle Guide and Template. Fits 12-inch size.

Electric Water Valve

Regulates water supply in conjunction with the trimmers on/off switch, saving water and eliminating constant flow adjustment.

Plastic Sediment Tray

Designed to fit under all Whip Mix Model Trimmers and most other major brands, this high-impact polystyrene sediment tray is extremely strong yet lightweight and easy to clean. A central platform keeps trimmer clean and dry. Left and right side drain holes make for convenient placement on either side of sink. Unit measures 18" long x 16-1/2" wide x 3" high (45.7 cm x 42 cm x 7.6 cm).

Water Spray Attachment

A lever-operated water spray head for easy cleaning of casts.

Splash Shield

Clear plastic shield provides protection from water spray/debris. Model Trimmer Wheels 12" abrasive wheel available in medium, coarse and extra coarse grits. 10" available in medium and coarse grit.

Blue Wonder Diamond Wheels

Diamond particles embedded in fiberglass substrate provide fast, aggressive cutting action. One year guarantee. Available in 10-inch and 12-inch sizes.

12-Inch PSA Disc

Pressure-sensitive discs available in 50 grit.

Conversion Kit

Allows conversion from abrasive wheel to Pressure Sensitive Adhesive Disc System (PSA). Standard kit contains three 50 grit discs, plus a 12" Plastic Backing Disc.



Code	Name	Price	+
120223	WHIP MIX #36082 Splash Shield Assy F/3/4 Mt - #09669	US\$ 151.43	+Info
120225	WHIP MIX #36124 Cleancut Wheel Cleaner - #09602	US\$ 27.61	+Info
120231	WHIP MIX #35173 Splash Shield F/trimmer - #22373	US\$ 136.68	+Info
120232	WHIP MIX Sediment Tray Plastic WHITE - #22590	US\$ 296.38	+Info
120236	WHIP MIX #30043 Orthodontic Angle Guide - #21792	US\$ 72.37	+Info
120482	WHIP MIX #10225 Rubber Feet Large pk -4 - #96011	US\$ 28.82	+Info

Burs & Cutters for Model Trimmers

Burs & Cutterd for Model Trimmers



RAY FOSTER

RAY FOSTER

Code	Name	Price	+
470902	A420 SaburTooth Carbide Cylinder	US\$ 44.70	+Info
470903	A421 SaburTooth Carb Taper Cylinder	US\$ 46.19	+Info
470904	SaburTooth Carb Cylinder Fine	US\$ 42.47	+Info
470905	SaburTooth Carb Taper Cylin Fine	US\$ 46.19	+Info

RENFERT

RENFERT

Code	Name	Price	+
23-1806-0001	Millo Bur - Tap Cross-cut Fine	US\$ 153.59	+Info
23-1806-0002	Millo Bur - Tap Cross-cut 6bl	US\$ 150.52	+Info
23-1806-0003	Millo Bur - Cyl Cross-cut 6bl	US\$ 153.59	+Info

Model Trimmer Accesories - BUFFALO

Model Trimmer Accesories - BUFFALO

Code	Name	Price	+
425205	Water Conversion Kit For DualTrim	US\$ 55.97	+Info

Catalog Index / Indice

Arch Model Trimmers	1
Disc Model Trimmers	2
Discs for model trimmers	6
Model Trimmer Accessories	11