

Models

BASE MODEL Systems



PIN-CAST BASE FORMERS SYSTEM - RENFERT

PIN-CAST BASE FORMERS - RENFERT

Function & Performance

Model base production kit for effective formation of model bases for plaster arches.

Advantages

- * Saves time and material.
- * Straight forward base system with re-usable parts.
- * Variable base height and size.

Details

- * Accurately fitting primary and secondary bases.
- * Easy removal and precise relocation of the models.
- * Corrosion resistant discs and magnet holders.
- * Strong magnets for secure hold.

Use

- * Separating dental arches: Plaster against plaster.
- * For the dental arch base with pins.
- * Removal of the base plate.
- * Isolation of the model base: Plaster against plaster.
- * Filling out the split cast form with base plaster.
- * Removal of the rubber ring after the plaster has set.
- * Articulating the models in the articulator or occludator parallel to the table level.



REFILLS

Code	Name	Price	+
113934	RENFERT - PIN-CAST Retention Rings Zinc Plated 1000pk-For secure fixation of no...	US\$ 32.99	+Info
417459	RENFERT - Pin-Cast magnets and magnet holders 30pk- #4130000	US\$ 54.27	+Info
417460	RENFERT - Pin-Cast retention discs P120 - #4140000	US\$ 146.37	+Info

RUBER Rings

Other Models

Other Models

Code	Name	Price	+
417455	Large rubber rings for Bi-V-Pin 5pk	US\$ 34.05	+Info
417456	Small rubber rings for Bi-Pin 5pk	US\$ 34.05	+Info

417458	Small rubber rings for Bi-V-Pi 5pk	US\$ 34.05	+Info
--------	------------------------------------	------------	-----------------------

PIN-Cast

PIN-Cast

Code	Name	Price	+
417453	Pin-Cast Rubber Rings L short 5pk	US\$ 29.79	+Info
417454	Pin-Cast rubber Rings L 27mm 5pk	US\$ 29.79	+Info
417457	Small rubber rings for Bi-Pin L 5pk	US\$ 34.05	+Info

Sets / Kit

Code	Name	Price	+
023-4100326	RENFERT - PIN-CAST Set for for SHORT Bi-Pin w/ sleeve - BASE FORMER 13.5m - ...	US\$ 198.06	+Info
023-4100346	RENFERT - PIN-CAST Set for for LONG Bi-Pin w/ sleeve - BASE FORMER former 17.5...	US\$ 231.88	+Info
023-4100366	RENFERT - PIN-CAST Set for BI-V-PIN w/ sleeve - BASE FORMER 12.5m #4100-336 #4...	US\$ 241.54	+Info

PIN BASE SYSTEMS - READYBASE

PIN BASE SYSTEMS - READYBASE



MAXILAR LOWER

NO PINS

BESQUAL

BESQUAL

Code	Name	Price	+
421081	Opposing / No Pins	US\$ 67.50	+Info
421087	Opposing / No Pins	US\$ 78.00	+Info

READYBASE

READYBASE

Code	Name	Price	+
421083	Opposing / No Pins	US\$ 67.50	+Info
421097	Opposing / No Pins	US\$ 36.92	+Info

WITH STONE BASE AND PINS

BESQUAL

BESQUAL

Code	Name	Price	+
421088	S-Base Pinned Stone Low FA 100pk	US\$ 213.00	+Info

READYBASE

READYBASE

Code	Name	Price	+
421091	with Stone Base and Pins	US\$ 84.00	+Info
421096	with Stone Base and Pins	US\$ 109.50	+Info

MAXILAR UPPER

NO PINS

Code	Name	Price	+
421089	BESQUAL S-Base No Pin Blanks (Opposing), Upper Full Arch, 100/Box, BROWN - #539...	US\$ 78.00	+Info

WITH STONE BASE AND PINS

BESQUAL

BESQUAL

Code	Name	Price	+
421080	with Stone Base and Pins	US\$ 168.00	+Info
421082	with Stone Base and Pins	US\$ 168.00	+Info

READYBASE

READYBASE

Code	Name	Price	+
421086	with Stone Base and Pins	US\$ 213.00	+Info
421090	with Stone Base and Pins	US\$ 84.00	+Info
421095	with Stone Base and Pins	US\$ 109.50	+Info

QUADRANT

Code	Name	Price	+
421093	READYBASE - BesQual Master Base Partial No Pins / Opposing for MP1 - Partial / Q...	US\$ 45.44	+Info
421085	READYBASE - BROWN S-Base Partial No Pins / Opposing for BP2 - Partial / Quadrant...	US\$ 50.76	+Info
421092	READYBASE - BesQual Master Base Partial with Stone Plate & Pins - Partial / Quad...	US\$ 108.00	+Info
421084	READYBASE - BROWN S-Base Partial with Stone Base & Pins - Partial / Quadrant wit...	US\$ 124.50	+Info

ADAPTOR

Code	Name	Price	+
421094	READYBASE - BesQual Master Base Partial Articulator Adaptor - Jig - Partial / Qu...	US\$ 21.60	+Info
421098	READYBASE - BesQual Master Base Full Arch or Anterior Adaptor Jig - Full or Ante...	US\$ 21.60	+Info

Bi-PIN - RENFERT



Other Models

Other Models

Code	Name	Price	+
417550	Pin Guide for Bi-Pins 1pk	US\$ 41.54	+Info

Models

Model Formers



Model Formers - BESQUAL

Model Formers - BESQUAL



	OUTER DIMENSIONS	INNER DIMENSIONS
Small:	2.55x3.00x0.75 in	(2.20x2.65x0.60 in)
Medium:	2.70x3.25x0.75 in	(2.35x2.90x0.60 in)
Large:	2.80x3.50x0.75 in	(2.45x3.15x0.60 in)
Quadrant Left:	2.70x1.80x0.75 in	(2.40x1.35x0.60 in)
Quadrant Right:	2.70x1.80x0.75 in	(2.40x1.35x0.60 in)

Code	Name	Price	+
475019	Base Mold FULL ARCH LARGE	US\$ 6.50	+Info
475020	Base Mold FULL ARCH MEDIUM	US\$ 6.50	+Info
475021	Base Mold FULL ARCH SMALL	US\$ 6.43	+Info
101674	Base Mold S-M-L-2Q 5pk	US\$ 18.40	+Info
421054	Base Mold Quad Right Each	US\$ 6.50	+Info
421055	Base Mold Quad Left Each	US\$ 6.50	+Info

Model Formers - BEGO

Model Formers - BEGO

Code	Name	Price	+
------	------	-------	---

419883	MODEL BASE FORMERS U1	US\$ 73.78	+Info
419884	MODEL BASE FORMERS U2 5 PIECES	US\$ 73.78	+Info
419885	MODEL BASE FORMER	US\$ 79.38	+Info
419886	MODEL BASE FORMER 02- Pk OF 5	US\$ 79.38	+Info
471443	Model Base Former, Assorted 1 pack	US\$ 65.84	+Info

Model Formers - PINDEX




Code	Name	Price	+
426555	COLTENE Pindex Flexible Rubber Base Mold - Full Arch, Single Mold - # PX147	US\$ 61.87	+Info
436467	COLTENE Pindex Flexible Rubber Base Molds - Right and Left Quadrants - # PX148	US\$ 88.38	+Info
436479	COLTENE Pindex Flexible Rubber Base Molds, Box of 3 Molds: Full Arch, Right and ...	US\$ 112.21	+Info

MODEL FORMERS - PLASDENT

MODEL FORMERS - PLASDENT

Model formers are the easiest way to create bases when pouring up impressions for models



DISPOSABLE MODEL FORMERS

Assorted

Assorted

Code	Name	Price	+
001-MFD-AST-120	Disp. Model Formers 120pk	US\$ 26.16	+Info

Disposable M. Formers

Disposable M. Formers

Code	Name	Price	+
001-MFD-AST-60	Disp. Model Formers 60pk Ea	US\$ 12.29	+Info
001-MFD-DL	Disp. Mdl Formers Lg/dp Cl 60p	US\$ 14.83	+Info
001-MFD-L	Disp. Mdl Formers Lg Cl 60pk	US\$ 14.84	+Info
001-MFD-M	Disp. Mdl Formers Med Cl 60pk	US\$ 12.99	+Info

SILICONE MODEL FORMERS

Code	Name	Price	+
001-MF-A	PLASDENT Model Formers Kit-#MF-A-6 Pcs ASSORTED Size/Box	US\$ 6.17	+Info
001-MF-S	PLASDENT Small Model Formers-#MF-S-Color: YELLOW-3 Pcs Upper & 3 Pcs Lower/Box	US\$ 4.36	+Info
001-MF-M	PLASDENT Medium Model Formers-#MF-M-Color: GREEN-3 pcs upper & 3 pcs lower/box	US\$ 6.39	+Info
001-MF-L	PLASDENT Large Model Formers-#MF-L-Color: BLUE-3 Pcs Upper & 3 Pcs Lower/Box	US\$ 8.01	+Info

Model Formers - RENFERT

Model Formers - RENFERT

Code	Name	Price	+
417452	Pin-Cast base former Small 3pk	US\$ 88.05	+Info

Catalog Index / Indice

BASE MODEL Systems	1
Model Formers	4