

## Dowel Pins

### Accu trays



### Accu-Trac - COLTENE

Accu-Trac - COLTENE



Code	Name	Price	+
063-ATR_12	Accu-Trac Refill Kit	US\$ 150.45	+Info
417363	Accu-Trac Refill Kit	US\$ 314.31	+Info
426671	Accu-Trac Precision Die System Intr	US\$ 129.64	+Info
436469	Accu-Trac Precision Die System - Ar	US\$ 212.50	+Info
436470	Accu-Trac Precision Die System - Ar	US\$ 103.61	+Info
436472	Accu-Trac Precision Die System - Mo	US\$ 84.08	+Info
436473	Accu-Trac Precision Die System - Ar	US\$ 25.51	+Info
436474	Accu-Trac Precision Die System - Lo	US\$ 104.21	+Info
436476	Accu-Trac Precision Die System - Re	US\$ 2,562.52	+Info
436478	Accu-Trac Precision Die System - In	US\$ 1,535.14	+Info

## Dowel Pins

### Aluminum Dowel pins



### DOWEL PIN ALUMINUM - ECCO

DOWEL PIN ALUMINUM - ECCO

**Features:**

- \* Aluminum, which is light, ductile, conductive, and very resistance to oxidation.
- \* Compatible with regular pindex drill
- \* Dimensions: 23mm length - 18mm taper - 5mm head
- \* Available in packages of 100pk or 1000pk

For more related products please visit [www.usdentaldepot.com](http://www.usdentaldepot.com)



**ECCO**  
ECCO

Code	Name	Price	+
100447-1	Dowel Pin Med Alum 100pk	US\$ 4.00	<a href="#">+Info</a>

**JSP**  
JSP

Code	Name	Price	+
100447	Dowel Pin Med Alum 1000pk	US\$ 12.78	<a href="#">+Info</a>

**Dowel Pins**

**Brass dowel pins ( bronze )**



**Brass dowel pins ( bronze ) - DIGITECH**

Brass dowel pins ( bronze ) - DIGITECH

Code	Name	Price	+
100026	Dowel Pin Med Brass #2 1000pk	US\$ 21.50	<a href="#">+Info</a>

**Brass dowel pins ( bronze ) - KEYSTONE**

Brass dowel pins ( bronze ) - KEYSTONE

Code	Name	Price	+
446294	Dowel Pins #2 Medium 100pk	US\$ 9.12	<a href="#">+Info</a>

**DOWEL PIN BRASS - KEYSTONE**

**DOWEL PIN BRASS - KEYSTONE**

These Dowel brass pins are made of top quality Brass  
Dowel pins are available in a variety of sizes and quality materials.

- \* Available packages: 100pk or 1000pk
- \* Available Size: Small or Medium



Code	Name	Price	+
034-1300015	KEYSTONE- Dowel Pins Medium Brass #2 - 1000pk - K # 1300015 Dowel Pin	US\$ 54.18	<a href="#">+Info</a>
034-1310011	KEYSTONE- Dowel Pins Small Brass Anterior #1 - 100pk -	US\$ 9.12	<a href="#">+Info</a>
034-1310010	KEYSTONE Swiss Dowel Pins #1, Small, 1000/Pk. Made of top quality brass - #1310010	US\$ 64.82	<a href="#">+Info</a>

Dowel Pins  
Bronze Dowel Pins



Bronze Dowel Pins - KEYSTONE

Bronze Dowel Pins - KEYSTONE

Code	Name	Price	+
034-1310210	805-Js J - Pins -Universal M	US\$ 74.83	<a href="#">+Info</a>

Dowel Pins  
Ceramic Dowel Pin



Ceramic Dowel Pin - GC AMERICA

Ceramic Dowel Pin - GC AMERICA

Code	Name	Price	+
422054	G-Cera Ceramic Zirconia Refractory	US\$ 86.02	<a href="#">+Info</a>

Dowel Pins  
Double Dowel Pins



DR. PIN DOUBLE DOWEL PIN - ECCO

DR. PIN DOUBLE DOWEL PIN - ECCO

Features:

Precise Fitting Metal Pins & Sleeves with Rubber Caps [F-210]

Quantity: 100pk and 1000pk



Code	Name	Price	+
100032	ECCO - DR PIN Double Dowel Pins with Metalic Sleeves - 1000pk- ( D.R. ) ( MM325 ...	US\$ 51.81	<a href="#">+Info</a>
100032-1	ECCO - DR PIN Double Dowel Pins with Metalic Sleeves - 100pk- Precise Fitting Me...	US\$ 8.03	<a href="#">+Info</a>

EURO PIN DOWEL w/ PLASTIC SLEEVE 2 LEGS - DIGITECH

EURO PIN DOWEL w/ PLASTIC SLEEVE 2 LEGS

Available in 100pk & 1000pk

Features

- \* More accurate Pin
- \* Smooth, when male gets fitted into female sleeve.
- \* After insertion, male & female are adhered together accurately.
- \* Two solid metal pins.
- \* Its female sleeve is separated into two parts for its better grip and accuracy.
- \* Precision joining in the portion of retention head.



- \* With Plastic Sleeve
- \* Nickel Plated Zinc
- \* Extended pin head allows for increased stability
- \* Able to compensate for the expansion of anhydrite while hardening

Code	Name	Price	+
114800	EURO-PIN Double dowel pin with Plastic Sleeves 1000pk Precise Fitting Metal Pins...	US\$ 56.28	<a href="#">+Info</a>
114801	EURO-PIN Double dowel pin with Plastic Sleeves 100pk Precise Fitting Metal Pins...	US\$ 9.00	<a href="#">+Info</a>

### MASTER PIN I DOWEL PIN with METAL SLEEVE EVE 1 1/2 LEGS - ECCO

#### MASTER PIN I DOWEL PIN WITH METAL SLEEVE 1 1/2 LEGS - ECCO

**Features:**

Master Double Twin Pin with Metal Sleeves and Rubber Caps.  
Die guards for added stability.  
Box of 100 or 1000 sets - 1 1/2 Leg



Code	Name	Price	+
113749	MASTER II Double Dowel Pin with Metalic Sleeves ( 1.5 legs ) and Rubber Caps Die...	US\$ 71.28	<a href="#">+Info</a>
113749-1	MASTER II Double Dowel Pin with Metalic Sleeves ( 1.5 legs ) and Rubber Caps Die...	US\$ 11.36	<a href="#">+Info</a>

### MASTER PIN II DOUBLE DOWEL PIN w/ PLASTIC SLEEVE EVES 1 1/2 LEGS - ECCO

#### MASTER PIN II DOUBLE DOWEL PIN w/ PLASTIC SLEEVES 1 1/2 LEGS - ECCO

**Features:**

Master Double Twin Pin with Plastic Sleeves and Rubber Caps.  
Perfect fit without movements.  
Die guards for added stability.  
-Box of 100 or 1000 sets - 1 1/2 Leg



Code	Name	Price	+
113750	MASTER II Double Dowel Pins with Plastic Sleeves 1000pk Precise Fitting Metal Pi...	US\$ 81.00	<a href="#">+Info</a>
113750-1	MASTER II Double Dowel Pins with Plastic Sleeves 100pk Twin Pin - Precise Fittin...	US\$ 10.50	<a href="#">+Info</a>

**BI-FIX MODEL PIN- RENFERT**



Code	Name	Price	+
417445	RENFERT - Bi-Fix Pin Model pin Dual guidance model pin with fixing wire for exa...	<b>US\$ 36.18</b>	<a href="#">+Info</a>
417446	RENFERT - Bi-Fix Pin Model pin Dual guidance model pin with fixing wire for exa...	<b>US\$ 293.27</b>	<a href="#">+Info</a>

**BI-PIN LONG WITH SLEEVE - RENFERT**

**BI-PIN LONG WITH SLEEVE - RENFERT**

**Model pin with dual guidance for precise, rotation-free positioning of the die.**

- \* Pin and sleeve are high-lustre nickel plated - no oxidisation.
- \* Only one drill-hole per die necessary.

**Contains:** 100pk or 1000pk

**US Dental Depot is an Authorized RENFERT Dealer in US.**

**#USDDRENFERT**



Code	Name	Price	+
113825	RENFERT - Bi-Pin Long with sleeve 1000pk-#346-2000 #3462000 -Double Dowel Pins-...	<b>US\$ 260.09</b>	<a href="#">+Info</a>
113824	RENFERT - Bi-Pin Long with sleeve 100pk-#346-1000 #3461000 -Double Dowel Pins- ...	<b>US\$ 31.92</b>	<a href="#">+Info</a>

## BI-PIN LONG WITH SLEEVE AND SPIKE - RENFERT

### BI-PIN LONG WITH SLEEVE AND SPIKE - RENFERT

**Features:**

High gloss nickel-plated modeling pin with double guide for exact and anti-rotational positioning of the stump.

**Advantage:**

- \* Shorten the time of elaboration of the model.
- \* Two guides as anti-torsion protection.
- \* Optimal friction thanks to the conical shape of the guide.

Available in 100 & 1000 per pack



Code	Name	Price	+
417443	RENFERT - Bi-Pin long with sleeve and spike P100 - High-lustre nickel-plated mo...	US\$ 37.91	<a href="#">+Info</a>
417444	RENFERT - Bi-Pin long with sleeve and spike #3432000	US\$ 503.75	<a href="#">+Info</a>

## BI-PIN w/o SLEEVES - RENFERT

### BI-PIN W/O SLEEVES - RENFERT

High-lustre nickel plated model pin with dual guidance.

**Features**

- \* Only one drill-hole per die necessary.
- \* Dual guidance prevents rotation.
- \* Precision fit with gentle friction between pin and sleeve.
- \* No oxidisation of pins even with time
- \* Narrow pin shape is suitable for cramped die situations.

**US Dental Depot is an Authorized RENFERT Dealer in US.**



Code	Name	Price	+
113838	RENFERT - Bi-Pin Short No Sleeves 1000pk-#336-2000 #3362000 -Double Dowel Pins-...	US\$ 168.44	<a href="#">+Info</a>
113837	RENFERT - Bi-Pin Long No Sleeves 1000pk-#353-3000 #3533000 -Double Dowel Pins-H...	US\$ 183.54	<a href="#">+Info</a>

## BI-PIN WITH SKIPE - RENFERT

### BI-PIN WITH SPIKE - RENFERT

**Features:**

High gloss nickel-plated modeling pin with double guide for exact and anti-rotational positioning of the stump.

**Advantage:**

- \* Shorten the time of elaboration of the model.
- \* Two guides as anti-torsion protection.
- \* Optimal friction thanks to the conical shape of the guide.

**Available in 100 & 1000 per pack**



Code	Name	Price	+
417447	RENFERT - Bi-Pin with spike NO sleeve P100 - High-lustre nickel-plated model pi...	US\$ 29.93	+Info
417448	RENFERT - Bi-Pin with spike NO sleeves 1000pk - #3523000	US\$ 207.99	+Info

## BI-PINS SHORT w/ SLEEVE & SPIKE - RENFERT

### BI-PIN SHORT WITH SLEEVE AND SPIKE - RENFERT

High-lustre nickel plated model pin with dual guidance. A special spike for fixation in the impression shortens the production time of the model.

**Feature**

- \* Long spike ensures secure hold in the impression.
- \* Plaster arch and model base can be made one after the other in the impression.
- \* Dual purpose model:
  - \* As a saw cut die.
  - \* As an anatomical model.
- \* Pin and sleeve are high lustre nickel plated ? no oxidization.

**Package:** 100 or 1000pk

**US Dental Depot is an Authorized RENFERT Dealer in US.**



Code	Name	Price	+
113819	RENFERT - Bi-Pin Short with sleeve and spike 1000pk-#323-2000 #3232000 -Double ...	US\$ 398.37	+Info
113822	RENFERT - Bi-Pin Short with sleeve and spike 100pk-#323-1000 #3231000 -Double D...	US\$ 18.96	+Info

### BI-PINS SHORT w/ SLEEVE - RENFERT

#### BI-PIN SHORT WITH SLEEVE - RENFERT

Model pin with dual guidance for precise, rotation-free positioning of the die.

#### Features

- \* Pin and sleeve are high-lustre nickel plated - no oxidation.
- \* Only one drill-hole per die necessary

**Package:** 100pk or 1000pk



Code	Name	Price	+
113835	RENFERT - Bi-Pin Short with sleeve 1000pk-#326-2000 #3262000 Double Dowel Pins-...	US\$ 221.98	<a href="#">+Info</a>
113834	RENFERT - Bi-Pin Short with sleeve 100pk-#326-1000 #3261000 -Double Dowel Pins-...	US\$ 31.92	<a href="#">+Info</a>

### BI-V-PINS DOUBLE DOWEL PIN w/ METAL SLEEVES - RENFERT

#### BI-V-PINS DOUBLE DOWEL PIN w/ METAL SLEEVES - RENFERT

Narrow model pin with two dowel pins, ideal for positioning in limited space.

- \* Metal sleeve
- \* Nickel-plated pin and sleeve - no oxidation.
- \* Precision fit with slight friction between pin and sleeve.

**Package:** 100 or 1000pk

**US Dental Depot is an Authorized RENFERT Dealer in US.**

#UsddRenfert



Code	Name	Price	+
113818	RENFERT - Dowel Bi-V-Pin with Metalic sleeve-1000pk-#328-2000 #3282000 -Double ...	US\$ 188.26	<a href="#">+Info</a>
113821	RENFERT - Dowel Bi-V-Pin with metal sleeve-100pk-#328-1000 #3281000 -Double Dow...	US\$ 31.92	<a href="#">+Info</a>

### Double Dowel Pins - DIGITECH

Double Dowel Pins - DIGITECH

Code	Name	Price	+
------	------	-------	---



114451	TWIN PIN Double 1000pk	US\$ 61.60	<a href="#">+Info</a>
--------	------------------------	------------	-----------------------

## Double Dowel Pins - KEYSTONE

Double Dowel Pins - KEYSTONE

Code	Name	Price	+
034-1310055	Double Dowel Pin w/Slv 1000pk	US\$ 119.77	<a href="#">+Info</a>

## MAINSTAY DOWEL PIN SYSTEM - WHIP MIX

### MAINSTAY DOWEL PIN SYSTEM - WHIP MIX

The MainStay Dowel Pin System is a uniquely designed system for fabricating crown and bridge working models with sectioned, removable dies. The protective sleeve and tapered pin's unique rectangular shape prevents the rotation of sectioned dies utilizing a single pin application.

For the other side, the Mainstay II dowel pin system is designed for labs that want the added stability of a dual pin while maintaining an easy to use pinning system.


The Mainstay II pins are designed for use in crown and bridge working models with sectioned removable dies. Like the original Mainstay, Mainstay II utilizes a protective sleeve and rubber end caps. Its protective sleeve and tapered, dual pin design prevents reverse placement of sectioned dies. In addition the protective sleeve grips the dowel pin creating positive indexing while protecting the base stone material from abrasion, once abraded the base stone can alter the fit of the die causing costly rework.

Finally, the rubber end cap stabilizes quadrants or full arch models and eliminates the use of a dowel pin bur.

The Mainstay Drill Bit can be used for drilling holes for either Mainstay II or the original MainStay single pins.

- \* Asymmetric, dual pin design eliminates reverse placement of dies
- \* Ribbed plastic sleeve holds tight while pin removes easily.
- \* Rubber end cap stabilizes case while pouring base stone
- \* Metal to stone contact at bottom of pin for added stability

And finally, MainStay glue is specifically formulated for the dental laboratory. It offers an advanced cyanoacrylate adhesive that exhibits versatility, reliability, precision and extended shelf life. It



Instantly bonds plaster, metal, acrylic or porcelain.

### DOWEL PINS

Code	Name	Price	+
120445	WHIP MIX Mainstay li Dual Dowel Pin System pk /1000 - #08240	US\$ 232.76	<a href="#">+Info</a>
120447	WHIP MIX Mainstay Dowel Pin System pk /1000 - #08236	US\$ 232.76	<a href="#">+Info</a>

### BI-V-PINS DOUBLE DOWEL PIN w/ PLASTIC SLEEVES - RENFERT

#### BI-V-PIN WITH PLASTIC SLEEVE - RENFERT

Narrow model pin with two dowel pins, ideal for positioning in limited space.

#### Description

- \* V-shaped profile suitable for use with close standing dies such as lower anteriors.
- \* Only one drilling needed for each die.
- \* Tapered shape glides easily into final position.
- \* Open-ended sleeve prevents clogging and simplifies cleaning.
- \* Reproducible positioning without vertical adjustment of bite.
- \* Small base thanks to reduced pin length (12.5 mm).



**Package:** 100 or 1000pk

**US Dental Depot is an Authorized RENFERT Dealer in US.**

Code	Name	Price	+
113817	RENFERT - Dowel Bi-V-Pin with plastic sleeve-1000p-#329-2000 #3292000 -Double D...	US\$ 142.87	<a href="#">+Info</a>
113820	RENFERT - Dowel Bi-V-Pin with plastic sleeve-100pk-#329-1000 #3291000 -Double D...	US\$ 28.94	<a href="#">+Info</a>

### PINDEX Double Dowel Pins | COLTENE



#### Self-Articulating Pins

Code	Name	Price	+
427075	COLTENE Pindex Self-Articulating Pins and Plastic Sleeves, Refill Pack of 100 Pi...	US\$ 80.59	<a href="#">+Info</a>
431666	COLTENE Pindex Self-Articulating Pins and Plastic Sleeves, Bulk Pack of 1000 Pin...	US\$ 332.85	<a href="#">+Info</a>


436477	COLTENE Pindex Plastic Sleeves for Self Articulating Pins, Box of 1000 Sleeves - ...	US\$ 209.17	<a href="#">+Info</a>
--------	--	-------------	-----------------------

## MR PINS DOUBLE DOWEL PIN w/ METAL SLEEVES - SELECT DENT AL

### MR PINS DOUBLE DOWEL PIN w/ METAL SLEEVES - SELECT

**Features:**

- \* Pins and sleeves
- \* Compatible with virtually all pin setting systems
- \* Very economical
- \* Precise fitting metal pin & metal sleeve with improved retention
- \* Die guards for clean professional models



**Package:** 100 or 1000pk

Code	Name	Price	+
4001998	SELECT M.R. PIN Double Dowel Pins with Metalic Sleeve - 1000pk	US\$ 129.56	<a href="#">+Info</a>
4001998-1	SELECT M.R. PIN Double Dowel Pins with Metalic Sleeve - 100pk ( from 1000 bag S#...	US\$ 13.28	<a href="#">+Info</a>

**Dowel Pins**

## Reverse Dowel Pins



## Reverse Dowel Pins - BESQUAL

Reverse Dowel Pins - BESQUAL

Code	Name	Price	+
100257	Reverse Dowel Pin 1000pk	US\$ 14.32	<a href="#">+Info</a>

## Reverse Dowel Pins - KEYSTONE

Reverse Dowel Pins - KEYSTONE

Code	Name	Price	+
034-1310070	Reverse Dowel Pin 1000pk	US\$ 22.50	<a href="#">+Info</a>

**Dowel Pins**

## Rubber Caps for pins



## Rubber Caps for pins - BESQUAL

Rubber Caps for pins - BESQUAL

Code	Name	Price	+
115865	Dowel pin locator 1000pk	US\$ 14.32	<a href="#">+Info</a>

## Rubber Caps for pins - RENFERT

Rubber Caps for pins - RENFERT

Code	Name	Price	+
113828	Rubber caps for BI-PINS 500pk	US\$ 36.18	<a href="#">+Info</a>

### Dowel Pins

## Single Dowel Pin



## PRO-FIX WITH PLASTIC SLEEVE - RENFERT

PRO-FIX WITH PLASTIC SLEEVE - RENFERT

### Function & Performance

Model pin with Easy-Glide-Function and a guidance pin.



### Advantages

- \* Stress-free removal through gentle die removal without jolting movements.
- \* Short friction path of 1 mm for safe removal of diverging dies.
- \* Guarantees even friction strength with all dies.

Code	Name	Price	+
113829	RENFERT - Pro-Fix Precision pin 1000pk-#367-2000 #3672000 -Easy-Glide Now also ...	US\$ 136.30	<a href="#">+Info</a>
113830	RENFERT - Pro-Fix Precision pin 100pk-#367-1000 #3671000-Easy-Glide Now also fo...	US\$ 23.94	<a href="#">+Info</a>

### Dowel Pins

## Single Dowels w/ Sleeves




## ROUND Single Pins w/ Sleeves



**Single Dowels w/ Sleeves - DIGITECH**

Single Dowels w/ Sleeves - DIGITECH



Code	Name	Price	+
101549	MAINSTAY Type Dowels 1000pk	US\$ 41.58	+Info

**Single Pins w/ Sleeves - COLTENE**

Single Dowels w/ Sleeves - COLTENE



Code	Name	Price	+
063-PX114	PX114 Long Pins w/Slve 1000pk	US\$ 133.87	+Info
436471	Pindex Short (10 mm) Pins and Plast	US\$ 33.75	+Info
436466	Pindex Medium (16 mm) Pins with Syn	US\$ 106.90	+Info
063-PX115	Short Pins&Plastic Sleeves1000pk	US\$ 138.56	+Info

### MASTER PIN III SINGLE DOWEL PIN

#### MASTER PIN III SINGLE DOWEL PIN - ECCO


Diamond shape provides added stability and prevents from bending, rotation and movements  
Slim pin good for anterior



Available Packages: 100Pcs or 1000Pcs

Code	Name	Price	+
113751-1	MASTER PIN II 100pk Single Dowel Pin Diamond shape with Plastic Sleeves and Rubb...	US\$ 9.90	+Info
113751	MASTER PIN II 1000pk Single Dowel Pin Diamond shape with Plastic Sleeves and Rub...	US\$ 64.50	+Info

### FLAT PINS - Flat dowel pins



### MR FLAT PIN AND SLEEVE - SELECT

#### MR FLAT PIN AND SLEEVE - SELECT

#### Advantages of the MR Flats:

- \* New concept in Dowel Pins
- \* Flat side prevents rotation associated with round pins
- \* Tapered side assures pins are always replaced properly
- \* Great for anterior and posterior dies
- \* Plastic sleeve has die guard built in (No need to assemble- just snap off end)



### DIGITECH

#### DIGITECH

Code	Name	Price	+
4001991	MR Flat Pin w/ Plast- 100pk	US\$ 9.39	+Info

### SELECT DENTAL

#### SELECT DENTAL

Code	Name	Price	+
4001990	MR Flat Pin w/ Plast- 1000pk	US\$ 98.00	+Info

### TELESCOP FLAT DOWEL PIN - ECCO

#### TELESCOP FLAT DOWEL PIN - ECCO

Female sleeve of Telescop Soft Pin is made of soft Plastics. This soft female part provides users feeling tight in the pin-process.

#### The Advantages

- \* Female sleeve of Telescop Soft Pin is made of soft Plastics.
- \* This soft female part provides users feeling tight in the pin-process.
- \* Female sleeve is united shape of sleeve-body and bottom-cap,
- \* so there is no need to put the cap on the pin botherly.
- \* Users can separate Female's cap easily.
- \* Females' keeping ball can catch the male.

**Total length:** 19.3mm

**Available Packages:** 100pk or 1000pk



### DIGITECH

#### DIGITECH

Code	Name	Price	+
114798	Telescop Flat Pin w/ P 1000pk	US\$ 44.76	+Info


### ECCO

#### ECCO

Code	Name	Price	+
114799	Telescop Flat Pin w/ P 100pk	US\$ 5.49	+Info

### Plastic Single Dowels w/ Sleeves - ECCO

Single Dowels w/ Sleeves - ECCO



Code	Name	Price	+
100397	Plastic Pin w/ Sleeve 1000pk	US\$ 24.26	+Info

Dowel Pins  
Stick Pins



Stick Pins - BESQUAL

Stick Pins - BESQUAL

Code	Name	Price	+
114452	Spike Dwl Pin Brass #2 100pk	US\$ 20.00	+Info

Stick Pins - KEYSTONE

Stick Pins - KEYSTONE

Code	Name	Price	+
034-1310040	Stick Dowel Pin Brass 100pk	US\$ 15.30	+Info

Stick Pins - RENFERT

Stick Pins - RENFERT

Code	Name	Price	+
417450	Smart-Pin short 1000pk	US\$ 80.86	+Info
417451	Sleeves for Smart-Pin 1000pk	US\$ 80.86	+Info



## Catalog Index / Indice

Accu trays .....	1
Aluminum Dowel pins .....	1
Brass dowel pins ( bronze ) .....	2
Bronze Dowel Pins .....	3
Ceramic Dowel Pin .....	3
Double Dowel Pins .....	3
Reverse Dowel Pins .....	11
Rubber Caps for pins .....	11
Single Dowel Pin .....	12
Single Dowels w/ Sleeves .....	12
Stick Pins .....	16