

Electro Polishers

Electro polisher Units

Electro polisher Units - AMERICAN

Electro polisher Units - AMERICAN

| Code | Name | Price | + |
|--------|---------------------------|-------------|-----------------------|
| 101375 | Electro-Polisher 1Fr 110v | US\$ 639.44 | +Info |

Electro polisher Units - BESQUAL

Electro polisher Units - BESQUAL



S730

S730

| Code | Name | Price | + |
|--------|----------------------------------|-------------|-----------------------|
| 123440 | Electro Polishing Unit S730-220V | US\$ 575.25 | +Info |

S730 Chrome Cobalt

S730 Chrome Cobalt

| Code | Name | Price | + |
|--------|-----------------------------|-------------|-----------------------|
| 114037 | Electro Polishing Unit 110v | US\$ 575.25 | +Info |

ELTROPOL - ELECTROLYTIC POLISHING UNITS - BEGO

ELTROPOL - ELECTROLYTIC POLISHING UNITS - BEGO

- * Automatically calculates polishing time for up to two objects.
- * Innovative new heating concept quickly brings the unit up to operating temperature.
- * Major time-saving due to simultaneous polishing of two Co-Cr partial dentures.
- * New, user-friendly operating panel with display and soft keys.
- * Indicator to show when the solution in the polishing bath is in need of change.
- * Easy drainage of tank via drain hose.
- * Excellent polishing results ensured by uniform movement of the polishing bath.
- * Includes a supplementary cathode for objects with a deep palate.
- * Uniform polishing effect achieved by automatic current stabilization.



| Code | Name | Price | + |
|--------|------------------------------------------------------------------------------------|---------------|-------|
| 101348 | BEGO Eltropol 300 Control Unit 110vol ts 50Hz - Each #26310 (electropolishing ... | US\$ 2,966.75 | +Info |

Electro Polishers

Electropolishers Solutions



Electropolishers Solutions - BEGO

Electropolishers Solutions - BEGO

| Code | Name | Price | + |
|--------|---------------------------|------------|-------|
| 416534 | Wirolyt Polish Liquid 1lt | US\$ 43.44 | +Info |

Electropolishers Solutions - CFS

Electropolishers Solutions - CFS

| Code | Name | Price | + |
|--------|--------------------------------|------------|-------|
| 114911 | Electro Polisher Solution 1 ga | US\$ 72.00 | +Info |

Catalog Index / Indice

| | |
|----------------------------------|---|
| Electro polisher Units | 1 |
| Electropolishers Solutions | 2 |