

Rubber Polishers

Acrylic & Flexlble Polishers



Acrylic & Flex. Polishers - ATTRITOR

Acrylic & Flex. Polishers - ATTRITOR



Code	Name	Price	+
416688	Acry-Caps Coars 5pk	US\$ 8.63	+Info
ACRYCAPS-KIT	Acry-Caps Coars 15pk	US\$ 18.73	+Info
416689	ATTRITOR Acry-Caps Med 5pk	US\$ 8.63	+Info
416690	ATTRITOR Acry-Caps Fine 5pk	US\$ 8.63	+Info
416692	ACRY-CAP POLISHING KIT	US\$ 34.59	+Info
416694	Acryl polisher Kit 8 pcs	US\$ 20.91	+Info

SILICON POINTS FOR FLEXIBLES - CFS

SILICONE MOUNTED POINTS FOR FLEXIBLE MATERIALS - CFS

HP mounted

Used for finishing & polishing & retouching and adjustments on all unbreakable material like Splint material, flexible materials like Valplast, Flexite, Deflex, GFM



Code	Name	Price	+
245106	CFS Rubber Points Small - 10pk - (mounted used for finishing & polishing & reto...	US\$ 9.00	+Info
246106	CFS Rubber Points Large 10pk - D#R103M - U11 U12 (mounted used for finishing & ...	US\$ 8.10	+Info
259106	CFS Diamont Caps - 5pk - #Rh101M (mounted used for finishing & polishing & reto...	US\$ 10.43	+Info
115868	CFS KENDA type - 5pk - #Rh106M] (mounted used for finishing & polishing & reto...	US\$ 8.64	+Info

Acrylic & Flexlble Polishers - SHOFU

Acrylic & Flexlble Polishers - SHOFU

Code	Name	Price	+
111659	Acrylic Contour/Finishing Kit	US\$ 187.43	+Info

Acrylic & Flex. Polishers - CFS

Acrylic & Flex. Polishers - CFS



Code	Name	Price	+
220106	Brown Rubber Wheels 100pk	US\$ 35.69	+Info

ACRY POINT - SHOFU

ACRY POINT POLISHERS - SHOFU

Advantajes:

Incorporating the latest polishing technology from SHOFU, AcryPoints are designed for adjusting and polishing to a high luster all acrylic appliances, either in the laboratory or Clinic.

- * Available in three silicon carbide grit
- * Use Dark Grey for adjustments, Brown for finishing and Light Grey for polishing
- * Use AcryPoint for: acrylic resin denture material

Us Dental Depot is an Authorized SHOFU Dealer in US.



ACRYPOINT BP1 - LARGER

Code	Name	Price	+
106766	SHOFU Acrypoint Polisher - Fine - Light Grey - BP1 - 6pk - #0428	US\$ 34.68	+Info
106764	SHOFU Acrypoint Polisher - Medium - BROWN - BP1 - 6pk - #0427	US\$ 35.04	+Info
106763	SHOFU Acrypoint Polisher - Coarse - Dark Grey - BP1 - 6pk - #0426	US\$ 34.68	+Info

ACRYPOINT PC2 - SMALLER

Code	Name	Price	+
106757	SHOFU Acrypoint Polisher - Coarse - Dark Grey - PC2-12pk - #0423	US\$ 35.39	+Info
106759	SHOFU Acrypoint Polisher - Medium - BROWN - PC2-12pk - #0424	US\$ 31.85	+Info
106761	SHOFU Acrypoint Polisher - Fine - Light Grey - PC2-12pk - #0425	US\$ 35.39	+Info

ACRYLIC POLISHING KITS - SHOFU

ACRYLIC POLISHING KIT HP - SHOFU

Features:

A complete chairside or lab system for adjusting and polishing acrylic materials.

- * Handpiece kit
- * Carbide burs for smooth, fast adjustment
- * AcryPoints are impregnated with high quality abrasives
- * Your choice of a sterilizable bur block or plastic kit
- * Use Acrylic Adjusting and Polishing kit for: acrylic resin denture material



Use Acrylic Adjusting and Polishing kit for: acrylic resin denture material

Us Dental Depot is an Authorized SHOFU Dealer in US.

Code	Name	Price	+
111653	SHOFU Acrylic Classic Plastic Hp Polishing Kit 0321 - #84T and 70 A - 2 each BP1...	US\$ 83.84	+Info

RUBBER POINT FOR FLEXIBLES - BESQUAL

RUBBER POINT FOR FLEXIBLES - BESQUAL

Used for finishing & polishing & retouching and adjustments on all unbreakable material like Splint material, flexible materials like Valplast, Flexite, Deflex, GFM



GREEN RUBBER POINT

Code	Name	Price	+
100779	BESQUAL Mounted HP Rubber Point Medium GREEN - for acrilycs and flexibles - 100p...	US\$ 52.50	+Info
101140	BESQUAL Mounted HP Rubber Point Large GREEN - for acrilycs and flexibles - 100pk...	US\$ 67.50	+Info
103931	BESQUAL Mounted HP Rubber Point Large GREEN - for acrilycs and flexibles - 10pk	US\$ 10.15	+Info
103929	BESQUAL Mounted HP Rubber Point Medium GREEN - for acrilycs and flexibles - 10pk	US\$ 8.22	+Info

RED RUBBER POINT

Code	Name	Price	+
101141	BESQUAL Mounted HP Rubber Point Large red- for acrilycs and flexibles - 100pk #2...	US\$ 67.50	+Info
100780	BESQUAL Mounted HP Rubber Point Medium RED - for acrilycs and flexibles - 100pk ...	US\$ 52.50	+Info
103932	BESQUAL Mounted HP Rubber Point Large red- for acrilycs and flexibles - 10pk	US\$ 10.57	+Info
103930	BESQUAL Mounted HP Rubber Point Medium RED - for acrilycs and flexibles - 10pk	US\$ 7.23	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US DENTAL DEPOT

+1 954 874 6325 
www.usdentaldepot.com 

SoftGrind Sleeves & Mandrel - BUFFALO



SoftGrind Grinding Sleeves

Code	Name	Price	+
425372	SoftGrind Grinding Sleeve Refills f	US\$ 35.93	+Info
425374	Soft Grinding Kit	US\$ 62.62	+Info
425101	Soft Grind Mandrel	US\$ 26.74	+Info

Rubber Polishers

Clasp polisher (Rubber)



RUBBER POINTS - KEYSTONE

RUBBER POINTS - KEYSTONE

The Superior Rubber Clasp Polishers are made of pure fine silicon carbide abrasives in synthetic rubber. This combination makes for excellent finishes on chrome with long wear.

Features:

- * Available Styles: Superior Clasp Polisher (Standard), Superior Rubber Points (Hard), Superior Blue, Flex Rubber Point, Superior Blue Clasp Polishers, Superior Green Points, Superior Clasp Polishers Black, Superior Red Clasp Polishers, Superior Special Hard Red Rubber Point
- * Available Size: 15/16"x1/4" (23.8 x 6.3mm)
- * Made with pure silicon carbide
- * Available in 100/Box




Code	Name	Price	+
034-1640050	KEYSTONE Superior Rubber Clasp Polishers BLACK, Medium, 100pk. Silicon - #1640050	US\$ 16.28	+Info
034-1640070	KEYSTONE Superior Rubber Clasp Polishers BLUE, Fine, 100pk. Silicon Carbide - #...	US\$ 16.28	+Info
034-1640075	KEYSTONE Superior Rubber Clasp Polishers BLUE, Flex, 100pk. Silicon Carbide - #...	US\$ 16.28	+Info
034-1640010	KEYSTONE Superior Rubber Clasp Polishers Dark Grey, Regular, 100pk. Silicon - #...	US\$ 16.28	+Info
034-1640060	KEYSTONE Superior Rubber Clasp Polishers GREEN, Medium Fine, 100pk. Silicon - #...	US\$ 16.28	+Info
034-1640080	KEYSTONE Superior Rubber Clasp Polishers red, Fine, 100pk. Silicon Carbide - #1...	US\$ 16.28	+Info
034-1640081	KEYSTONE Superior Rubber Clasp Polishers red, Superior Hard, 100pk. Silicon - #...	US\$ 16.28	+Info
034-1640020	KEYSTONE Wheels Rubber Hard Clasp Polish, Superior Hard Same Quality - #1640020	US\$ 16.28	+Info

RUBBER POINTS CLASP - BESQUAL

RUBBER POINTS CLASP - BESQUAL

Rubber wheel for super luster polishing; for Porcelains and Metals



Code	Name	Price	+
100490	BESQUAL RED HARD Rubber Clasp polisher 100pk - 240 Grit - fine hard #204-1035	US\$ 10.56	+Info
100745	BESQUAL BLACK Rubber Clasp polisher 100pk - 150 Grit Coarse oarse - #204-1024	US\$ 12.00	+Info
100744	BESQUAL Unmounted Rubber Wheels, Square Edge - GREEN, MEDIUM 100pk. Silicon - #...	US\$ 14.50	+Info
100743	BESQUAL BLUE FLEX Rubber Clasp polisher 100pk - 220 Grit fine - #204-1070	US\$ 11.88	+Info
100742	BESQUAL GRAY HARD Rubber Clasp polisher - 100pk - 240 Grit -X- fine stiff - #2...	US\$ 12.00	+Info

Rubber Polishers

Rubber wheels



INTERPROXIMAL RUBBER WHEELS - KEYSTONE

INTERPROXIMAL RUBBER WHEELS - KEYSTONE

We offer several different types of Interproximal Rubber Wheels. We have White Wheels that are fine and non-contamination, great for extra fine polishing. Green Wheels are medium fine grit, Red Wheels are fine grit, Blue Wheels are fine (non-precious metals), and Black Wheels are medium coarse.



Features:

- * Made of aluminum oxide
- * Available Wheels: White, Green, Red, Black
- * Size: 22 x 0.40 mm
- * Available in 25/Box

Code	Name	Price	+
034-1900850	KEYSTONE Interproximal Rubber Wheels, BLACK, Medium Coarse for Finishing 25pk ...	US\$ 15.38	+Info
034-1900855	KEYSTONE Interproximal Rubber Wheels, GREEN, Medium Fine All Purpose - #1900855	US\$ 15.38	+Info
034-1900860	KEYSTONE Interproximal Rubber Wheels, BLUE, Fine for Polishing Non-Precious - #...	US\$ 15.38	+Info
034-1900865	KEYSTONE Interproximal Rubber Wheels, red, Fine and Extremely Flexible, Silicon...	US\$ 15.38	+Info
034-1900870	KEYSTONE Interproximal Rubber Wheels, WHITE, Fine and Non-Contaminating - #1900870	US\$ 15.38	+Info

RUBBER WHEELS - BESQUAL



Code	Name	Price	+
100747	BESQUAL Rubber Square Edge Wheel Polishers, 100/Box. FINE HARD RED silicon car...	US\$ 14.50	+Info
100748	BESQUAL SQUARE EDGE RED SUPERIOR red- FINE (Alum Oxide) Rubber Wheel - For Fle...	US\$ 14.50	+Info
100749	BESQUAL SQUARE EDGE GRAY Rubber Wheels 100pk -240 gri FINE for metals #202-1005	US\$ 14.50	+Info
100750	BESQUAL SQUARE EDGE BLUE FLEX Rubber Wheels 100pk - 220 Grit FINE SOFT for met...	US\$ 14.50	+Info
100751	BESQUAL Unmounted Rubber Wheels, Square Edge - GREEN, MEDIUM Silicon carbide -...	US\$ 14.50	+Info
100752	BESQUAL Unmounted Rubber Wheels, Square Edge - BLACK, COARSE silicon carbide ...	US\$ 14.50	+Info

RUBBER WHEELS KNIFE EDGE - KEYSTONE

RUBBER WHEELS KNIFE EDGE - KEYSTONE

Knife Edge Rubber Wheels are made of aluminum carbide and come in several different styles for different purposes.



- * The Green wheel is medium fine which is an "All Purpose" polisher great for metal.
- * Blue Wheel is fine for posing non-precious, semi-precious and chrome cobalt.
- * The Black Wheel is medium coarse for finishing and polishing chrome cobalt metals.
- * The Red Wheel is fine and extremely flexible and made of silicon carbide.
- * The White Wheel is fine and non-contaminating for extra fine polishing. Its also extremely flexible.

Features:

- * Made of aluminum carbide
- * Red Wheels are made of silicon carbide
- * Available Sizes: 5/8" (15.8mm), 7/8" (22.2mm)
- * Available in 100/Box

Code	Name	Price	+
034-1900610	KEYSTONE KNIFE EDGE rubber Wheels Dk Grey 5/8 -15.8mm 100box #1900610	US\$ 26.80	+Info
034-1900830	KEYSTONE Knife Edge Rubber Wheels, GREEN, Medium Fine All Purpose For Polishing...	US\$ 26.80	+Info
034-1900835	KEYSTONE Knife Edge Rubber Wheels, BLUE, Fine For Polishing Non-Precious - #190...	US\$ 26.80	+Info
034-1900840	KEYSTONE Knife Edge Rubber Wheels, red, Fine and Extremely Flexible, Silicon - ...	US\$ 26.80	+Info
034-1900845	KEYSTONE Knife Edge Rubber Wheels, WHITE, Fine and Non-Contaminating For Extra ...	US\$ 26.80	+Info

034-1900875	KEYSTONE Knife Edge Rubber Wheels, BLACK, Medium Coarse 100pk - #1900875	US\$ 26.80	+Info
034-1900880	KEYSTONE Knife Edge Rubber Wheels, GREEN, Medium Fine All Purpose For Polishing...	US\$ 25.46	+Info
034-1900885	KEYSTONE Knife Edge Rubber Wheels, BLUE, Fine For Polishing Non-Precious - #190...	US\$ 26.80	+Info
034-1900990	KEYSTONE Knife Edge Rubber Wheels, red, Fine and Extremely Flexible, Silicon - ...	US\$ 23.58	+Info
034-1900995	KEYSTONE Knife Edge Rubber Wheels, WHITE, Fine and Non-Contaminating For Extra ...	US\$ 26.80	+Info


RUBBER WHEELS SUP HARD - KEYSTONE

SUPERIOR RUBBERS CLAPS POLISHERS - KEYSTONE

The Superior Rubber Clasp Polishers are made of pure fine silicon carbide abrasives in synthetic rubber. This combination makes for excellent finishes on chrome with long wear.

Features:

- * Available Styles: Superior Clasp Polisher (Standard), Superior Rubber Points (Hard), Superior Blue, Flex Rubber Point, Superior Blue Clasp Polishers, Superior Green Points, Superior Clasp Polishers Black, Superior Red Clasp Polishers, Superior Special Hard Red Rubber Point
- * Available Size: 15/16"x1/4" (23.8 x 6.3mm)
- * Made with pure silicon carbide
- * Available in 100/Box



Code	Name	Price	+
100489	KEYSTONE Rubber Polishing Wheels- 7/8in x 1/8in (22.2 x 3.1 mm), Superior Har...	US\$ 19.87	+Info
034-1900820	KEYSTONE Superior Rubber Polishing Wheels, 100pk, WHITE, 7/8 x 1/8in - #1900820	US\$ 22.05	+Info
034-1900815	KEYSTONE Rubber Polishing Wheels- 7/8in x 1/8in (22.2 x 3.1 mm), Superior RED...	US\$ 13.42	+Info
034-1900810	KEYSTONE Superior Rubber Polishing Wheels- 7/8in x 1/8in (22.2 x 3.1 mm) - #1...	US\$ 19.87	+Info
034-1900805	KEYSTONE Rubber Polishing Wheels- 7/8in x 1/8in (22.2 x 3.1 mm), Superior GRE...	US\$ 19.87	+Info
034-1900800	KEYSTONE Rubber Polishing Wheels 7/8in x 1/8in (22.2 x 3.1 mm) 100pk Superior...	US\$ 19.87	+Info
034-1900795	KEYSTONE Rubber Polishing Wheels- 7/8in x 1/8in (22.2 x 3.1 mm), Superior BLU...	US\$ 19.87	+Info
034-1900580	KEYSTONE Superior Rubber Polishing Wheels, 100pk, Coarse, GREEN Utility, 7/8 x ...	US\$ 17.63	+Info
034-1900560	KEYSTONE Super Rubber Polishing Wheels, 100pk, Superior Hard, Dark GRAY, 7/8 x ...	US\$ 15.10	+Info
034-1900550	KEYSTONE Rubber Polishing Wheels, Superior Regular, 7/8 x 1/8in (22.2 x 3.1 mm...	US\$ 19.86	+Info
034-1900540	KEYSTONE Superior Rubber Polishing Wheel, 100pk, Medium Soft, Light Grey, 7/8 -...	US\$ 15.10	+Info

Rubber Polishers

Rubber Wheels for Lathes

Lathe Rubber Wheel - CFS



Code	Name	Price	+
242106	CFS BROWN Rubber Wheels - Large - 3 inches used before giving polish to unbreaka...	US\$ 20.89	+Info

Rubber Wheels for Lathes - KEYSTONE

Rubber Wheels for Lathes - KEYSTONE

Code	Name	Price	+
034-1900050	Rubber Wheel Gray Regular 1pk	US\$ 11.88	+Info

Catalog Index / Indice

Acrylic & Flexible Polishers	1
Clasp polisher (Rubber)	4
Rubber wheels	5
Rubber Wheels for Lathes	7