

Separators & Diluent

Wax Separators

DIE LUBES



[DIE LUBE Separators | ADS](#)

[DIE EEZ ORIGINAL - ADS](#)

[DIE EEZ ORIGINAL DIE LUBES - ADS](#)

Features

- * Recommended for all types of dies
- * Flows smoothly, covering any type of surface
- * Contains no mineral oil
- * Forms a uniform thin film never losing detail or fit, and is specially important to deep critical areas
- * Does not dissolve or distort wax
- * Water soluble formula
- * When rinsed off with water Die Eez carries away surface contamination
- * May be painted on, sprayed or dipped



Code	Name	Price	+
117778	ADS DIE-Eez (Die Lubricant) regular qt - #D896-6 (Die-lubes or Die & wax sepa...	US\$ 45.60	+Info
117777	ADS DIE-Eez (Die Lubricant) regular 16oz. - #D896-5 (Die-lubes or Die & wax s...	US\$ 30.33	+Info
117776	ADS DIE-Eez (Die Lubricant) regular 8oz. - #D896-4 (Die-lubes or Die & wax se...	US\$ 20.05	+Info

117775	ADS DIE-Eez (Die Lubricant) regular 4oz. - #D896-3 (Die-lubes or Die & wax se...	US\$ 14.26	+Info
117774	ADS DIE-Eez (Die Lubricant) regular 2oz. - #D896-2 (Die-lubes or Die & wax se...	US\$ 9.24	+Info

DIE SEP LUBRICANT - ADS [DIE SEP LUBRICANT - ADS](#)

Features:

- * Lubricates stone, metal, plastic or porcelain dies
- * Water soluble - use on wet or dry stone
- * Separates wax patterns from all types of dental dies
- * Not recommended for lacquer die spacers
- * Discontinued by original manufacturer and now available from American Dental Supply, Inc



Code	Name	Price	+
111133	ADS DIE-Sep Lubricant - 1Quart #J896-6	US\$ 45.60	+Info
117859	ADS DIE-Sep 16oz. pt - #J896-5	US\$ 35.76	+Info
117858	ADS DIE-Sep 8oz. - #J896-4	US\$ 20.81	+Info
117857	ADS DIE-Sep 4oz. - #J896-3	US\$ 12.65	+Info
117856	ADS DIE-Sep 2oz. - #J896-2 (simil to Die Lube)	US\$ 8.40	+Info

MICRO LUBE - ADS [MICRO LUBE - ADS](#)

Easily releases completed wax patterns. Forms a fine, thin, non-oily film. Especially useful in the hydrocolloid technique. Completely water base and non-flammable. An excellent lubricant for crown and bridge work. Microlube enhances polishing of waxed cases and smooth wax ups. Cleans die surfaces and wax when rinsed off with water. Can be painted on or dipped.



Code	Name	Price	+
118205	ADS Microlube 4oz. - #M896-3	US\$ 14.26	+Info
118206	ADS Microlube 8oz. - #M896-4	US\$ 20.88	+Info
111134	ADS MicroLube - 1Pint #M896-5 easily releases completed wax patterns	US\$ 31.38	+Info
118207	ADS Microlube pt 16oz. - #M896-5	US\$ 31.38	+Info
118208	ADS Microlube qt - #M896-6	US\$ 50.97	+Info

DIE LUBE - NEY DIE LUBE - NEY

An excellent separating medium for wax patterns. It is non-oily & will not evaporate from the die, so it can be used sparingly. Since it is water soluble, it does not have to be washed off the pattern before investing. Acts as a surface tension reducer for the wax pattern.

Specifications:

- * For all dies: stone, epoxy and metal
- * Features a non-oily, non-evaporating formula



Code	Name	Price	+
9995105	Die Lube 8oz.	US\$ 22.00	+Info
9995106	THE CLASSIC NEY Die Lube 16oz.	US\$ 48.00	+Info

Die Lubes - SMILE LINE Die Lubes - SMILE LINE

Code	Name	Price	+
123782	Wax Separator JP52 20ml	US\$ 48.75	+Info

LUBRITEX N 12 - DIE LUBRICANT - WHIP MIX LUBRITEX N? 12 - DIE LUBRICANT - WHIP MIX

A water-soluble lubricant that is equally effective on all types of die materials. LUBRITEX will not react with die or pattern materials and will not affect investments.



Code	Name	Price	+
120419	WHIP MIX Lubritex #12 3.5 Oz - #09784	US\$ 23.87	+Info

Separators & Diluent Acrylic separators



Acrylic separators - BUFFALO Acrylic separators - BUFFALO

Code	Name	Price	+
424453	Foil Cote Liquid Separator Pink 1Qt	US\$ 14.13	+Info

LIQUID TIN-FOILS



Tin-Foil - JSP

Tin-Foil - JSP

Code	Name	Price	+
437628	Tin-Foil-Cote Separating	US\$ 12.00	+Info
437661	Tin-Foil-Cote Separating 32oz	US\$ 32.03	+Info
461603	Jsp - Tin-Foil Cote High	US\$ 12.00	+Info
461605	Tin-Foil Cote HV Separat. 32oz	US\$ 16.00	+Info
461606	Tin-Foil Cote Separat HV 1gl	US\$ 30.36	+Info
437632	Tin-Foil-Cote Separat Liq 1gallon	US\$ 35.94	+Info

Tin-Foils - ADS

LIQUI-FOIL PINK - ADS

LIQUI-FOIL PINK - ADS

For high production labs that demand superior performance and quality. Ideal performance. Ready to use. Adheres to stone, gels quickly, penetrates, and forms a fine tough film. Ideal viscosity results in smooth application whether painted on or poured. Requires only one coat. Film doesn't lift from edges. Smooth flow on stone and around teeth. Will not stick to teeth or metal. Rinses off cleanly with water. Two super stable formulas. Contains no alcohol to evaporate. Shipped year round in any temperature.



Code	Name	Price	+
118200	ADS Liqui-Foil PINK 1qt Quart - #T550-1pk - Acrylic Separators Liquid Tin-Foil	US\$ 13.41	+Info
118201	ADS Liqui-Foil PINK gallon - #T550-8-pk - Acrylic Separators Liquid Tin-Foil	US\$ 34.32	+Info
118199	ADS Liqui-Foil PINK 8oz. - #T550-4-pk - Acrylic Separators Liquid Tin-Foil	US\$ 8.00	+Info

LIQUI-FOIL Premiun - ADS

LIQUI-FOIL - ADS

For high production labs that demand superior performance and quality. Ideal performance. Ready to use. Adheres to stone, gels quickly, penetrates, and forms a fine tough film. Ideal viscosity results in smooth application whether painted on or poured. Requires only one coat. Film doesn't lift from edges. Smooth flow on stone and around teeth. Will not stick to teeth or metal. Rinses off cleanly with water. Two super stable formulas. Contains no alcohol to evaporate. Shipped year round in any temperature.



Code	Name	Price	+
116040	ADS LIQUI-FOIL CLEAR Tin Foil substitute 1lts / 1qts - #T550-6-CR - Acrylic Sep...	US\$ 13.41	+Info
118198	ADS Liqui-Foil CLEAR gal - #T550-8-cr - Acrylic Separators Liquid Tin-Foil	US\$ 34.32	+Info
118196	ADS Liqui-Foil CLEAR 8oz. - #T550-4-cr - Acrylic Separators Liquid Tin-Foil	US\$ 8.00	+Info

X-FOIL PINK - ADS

X - FOIL PINK - ADS

For high production labs that demand superior performance and quality. Ideal performance. Ready to use. Adheres to stone, gels quickly, penetrates, and forms a fine tough film. Ideal viscosity results in smooth application whether painted on or poured. Requires only one coat. Film doesn't lift from edges. Smooth flow on stone and around teeth. Will not stick to teeth or metal. Rinses off cleanly with water. Two super stable formulas. Contains no alcohol to evaporate. Shipped year round in any temperature.



Code	Name	Price	+
118450	ADS X-Foil Liqui-Foil Premium - PINK 1qt - #T 555-6-pk - Acrylic Separators Liqu...	US\$ 17.96	+Info
118451	ADS X-Foil Liqui-Foil Premium - PINK gallon - #S 823-1 - Acrylic Separators Liqu...	US\$ 52.11	+Info

X-FOIL PREMIUM - ADS

X-FOIL PREMIUM - ADS

For high production labs that demand superior performance and quality. Ideal performance. Ready to use. Adheres to stone, gels quickly, penetrates, and forms a fine tough film. Ideal viscosity results in smooth application whether painted on or poured. Requires only one coat. Film doesn't lift from edges. Smooth flow on stone and around teeth. Will not stick to teeth or metal. Rinses off cleanly with water. Two super stable formulas. Contains no alcohol to evaporate. Shipped year round in any temperature.



Code	Name	Price	+
118448	ADS X-Foil Liqui-Foil Premium - CLEAR qt - #T 555-8-CR - Acrylic Separators Liqu...	US\$ 19.03	+Info
118449	ADS X-Foil Liqui-Foil Premium - CLEAR gal - #T 555 - Acrylic Separators Liquid T...	US\$ 55.23	+Info

Tin-Foils - KEYSTONE

Code	Name	Price	+
034-1360025	KEYSTONE Kefoil Specially Formulated Multi-purpose Liquid Separator and Tin Foil...	US\$ 37.46	+Info

Tin-Foil - BUFFALO

Code	Name	Price	+
424587	BUFFALO Foil Cote Liquid Separator - Ready to use PINK Alginate Protective Coati...	US\$ 29.33	+Info

FOILCOTE FOIL SUBSTITUTE - WHIP MIX

FOILCOTE LIQUID FOIL SUBSTITUTE - WHIP MIX

When applied to plaster or stone, it seals the pores and forms a smooth surface that gives excellent separation of the denture from the mold. With this product eliminates the time-consuming tin foiling and resulting creases.



Code	Name	Price	+
120409	WHIP MIX Foilcote Lq 3.5 Oz - #09687	US\$ 13.72	+Info
120408	WHIP MIX Foilcote Lq RED Gallon/3.75l - #09709	US\$ 94.83	+Info
120407	WHIP MIX Foilcote 4gal Case - #27782	US\$ 360.21	+Info

Tin-Foil - GC AMERICA

Tin-Foil - GC AMERICA

Code	Name	Price	+
421809	Coe-Sep Tin-Foil Substitute.	US\$ 44.02	+Info

Acrylic separators - CFS

Acrylic separators - CFS

Code	Name	Price	+
237106	CFS Acrylic Separator 1Gallon	US\$ 39.90	+Info

Acrylic separators - ECCO

Acrylic separators - ECCO

Code	Name	Price	+
102360	ECCO Acrylic Separator 120ml	US\$ 6.36	+Info

Acrylic separators - GC AMERICA

Acrylic separators - GC AMERICA

Code	Name	Price	+
421651	Coe-Sep Tin-Foil Substitute.	US\$ 85.27	+Info
423003	Coe-Sep Tin-Foil Substitute.	US\$ 11.69	+Info
470708	Multi-Sep Universal Separating 8oz	US\$ 89.77	+Info

Acrylic separators - KEYSTONE

Acrylic separators - KEYSTONE



Diamond D

Diamond D

Code	Name	Price	+
470546	Diamond D Ultra Sep 1qt	US\$ 19.95	+Info

Acrylic separators - LANG

Acrylic separators - LANG

Code	Name	Price	+
114120	Release Agent Liquid Foil 4oz	US\$ 13.89	+Info
114121	Release Agent Liquid Foil 0.5G	US\$ 41.73	+Info

DENTURE SEP - ADS

DENTURE SEP - ADS

Features

- * Designed to function as a premium denture silicone insulating paste
- * Provides clean and fast devesting of all types of acrylic dentures
- * Prevents investment from sticking to plastic or porcelain teeth
- * Inert to denture base resin
- * Protects teeth and resin material from fracturing
- * Ideal viscosity spreads easily across the model and teeth
- * Thixotropic paste will cling to surfaces and not run or slump
- * Gingival and interproximal areas are formed smooth and clean with minimal finishing and polishing needed
- * Exceptional retention of fine details
- * Prevents damage and distortion during processing



Code	Name	Price	+
117764	ADS Denture Sep CLEANER 2oz. - #D400-3	US\$ 9.54	+Info
117763	ADS Denture Sep CATALYST 1/4oz. - #D400-2	US\$ 11.52	+Info
117765	ADS Denture Sep INSULATING PASTE 14oz. - #D400-6	US\$ 50.38	+Info
117762	ADS Denture Sep INSULATING PASTE kit 14oz + 1/4oz catalyst, mixing cups and stic...	US\$ 69.10	+Info
117761	ADS Denture Kit ASSORTED SAMPLE KIT - #D9999	US\$ 108.63	+Info

RELEASE-SEP PLUS - ADS

RELEASE-SEP PLUS - ADS

RELEASE-SEP PLUS is a high performance synthetic lubricant and separating agent designed for the needs of the dental industry.

A light coating cures in seconds to form a highly repellent mold releasing finish. Can be sprayed or painted on virtually any surface. The separating film is long lasting and durable. Will not wash off with water. Does not have to be reapplied for multiple releases.

RELEASE-SEP PLUS is an excellent stone to plaster separator. Apply to impression materials to release stone and create smooth perfect models. Apply to impression trays to release alginate. Use on rubber base molds for clean release of epoxy, polyurethanes and stone. Apply to articulators, casting rings and flasks for easy cleanup. No more messy Vaseline.

RELEASE-SEP PLUS releases waxes from dies and models. Releases thermoformed materials and splints from models. Will not alter dimensions of dies.

RELEASE-SEP PLUS eliminates scrubbing and scraping hardened impression and investment residues from sinks, trays and bowls. Saves time, labor and protects equipment. Finger pump spray bottles are available which produce a fine mist for even coating and no waste. Lasts longer than aerosol cans. Special dental formulation contains no chlorofluorocarbons or dimethyl silicone oils.



Code	Name	Price	+
118332	ADS Release Separator Plus spray 8oz. - #R811-4	US\$ 18.37	+Info
118333	ADS Release Sep Plus bottle 16oz. p1 - #R811-5	US\$ 35.31	+Info

Separators & Diluent

Die Sealers



DIE SEALER REFILL - ADS

[DIE SEALER REFILL - ADS](#)

Features:

- * Creates a smooth, releasing surface for intricate direct porcelain build-up on stone dies
- * Burns out clean and will not cause discoloration
- * Closes the open pores of stone dies
- * Protects dies from shipping or breaking
- * Provides a hard, smooth working surface without distorting margins



DIE SEALER REFILL

Code	Name	Price	+
117850	ADS Die Prep Dryer ea - #0790-1	US\$ 14.08	+Info
117851	ADS Die Sealer refill 1 oz. - #0235-2	US\$ 13.87	+Info
117852	ADS Die Sealer refill 4oz. - #0235-3	US\$ 31.22	+Info
117853	ADS Die Sealer refill 16oz. pt - #D235-61	US\$ 71.31	+Info

THINNER

Code	Name	Price	+
117854	ADS Die Sealer THINNER 1 oz. - #D235-7	US\$ 7.76	+Info

Die Sealers - ADS

Die Sealers - ADS

Code	Name	Price	+
117855	Die Sealer & Sep.1oz ea sealer	US\$ 30.83	+Info
118217	Model Giant Sealant Demonstrat	US\$ 60.75	+Info

DIE SHIELD SPRAY - ADS
[DIE SHIELD - ADS](#)

Features:

- * Hardens and toughens delicate refractory dies and models before duplication
- * Helps prevent abrasion and shipping damage
- * Ceramic liquid impregnates pores in refractory investment, dries quickly and hardens the surface
- * Supplied with spray head for easy application
- * Can also be applied by brushing or dipping toughens margins and critical areas
- * Effective to 3000F



Code	Name	Price	+
117861	ADS Die Shield SPRAY 8oz. - #D510-4	US\$ 39.45	+Info
117860	ADS Die Shield SPRAY 4oz. - #D510-3	US\$ 26.25	+Info

ISOFIX 2000 - RENFERT
ISOFIX 2000 - RENFERT

Function & Performance

Film-free isolation of plaster against plaster from a practical spray bottle. Ideal for producing saw cut models.

Advantages

- * Soaks into the plaster surface.
- * Seals the open pores of the plaster model.
- * No layer build-up.



Code	Name	Price	+
101468	RENFERT - Iso-Fix 2000 RefillLts Pack unit 2 x 1 Lts- #1720-2000 #17202000 resis...	US\$ 105.14	+Info

Separators & Diluent

Flexible Separators



Flexible Separators - VALPLAST

Flexible Separators - VALPLAST

Code	Name	Price	+
829-20256	V-Sep Separator Gallons - Each	US\$ 80.54	+Info

Separators & Diluent

Mold Release & Lubricant Silicones



EEZ OFF - ADS

EZ OFF - ADS

Features:

- * Non-oily, non-toxic universal anti-sticking agent
- * Lasts 3 times longer than can sprays
- * Forms a fine, super effective releasing and lubricating film
- * Pleasant aroma and unlimited applications
- * Contains no silicones or chlorofluorocarbons to harm environment
- * Spray alginate impressions and flasks to release stone
- * Releases die pins from investment
- * Apply to trays, bowls, articulators, flasks, sprue formers, sinks and counters for quick clean ups



Code	Name	Price	+
117948	ADS Eez Off spray bottle 8oz. - #E806-4 (Die-lubes or Die & wax separators sili...	US\$ 20.60	+Info
117949	ADS Eez Off refill qt - #E806-6 (Die-lubes or Die & wax separators silicones)	US\$ 59.26	+Info

GLYCERINE USP - ADS

GLYCERINE U.S.P - ADS

ADS Glycerin U.S.P. is 99.5%.



Code	Name	Price	+
118072	ADS Glycerine USP qt - #S883-20-6	US\$ 24.80	+Info
118071	ADS Glycerine USP pt 16oz. - #S883-20-5	US\$ 16.24	+Info

Lubricant Silicones - YETI

Lubricant Silicones - YETI

Code	Name	Price	+
417071	Cooling Lubricant for Wet	US\$ 47.30	+Info

Mold Release & Lubricant Silicones - BUFFALO

Mold Release & Lubricant Silicones - BUFFALO

Code	Name	Price	+
424455	Sta-Lube Silicone Inst Spray 10.5oz	US\$ 17.89	+Info

Mold Release Agent - CFS

Mold Release Agent - CFS



Code	Name	Price	+
211106	CFS Mold Release Spray ea	US\$ 15.92	+Info

Mold Release Agent - KEYSTONE

Mold Release Agent - KEYSTONE



Code	Name	Price	+
034-5940310	Tri-Epoxy Mold Release Spray	US\$ 13.72	+Info

Mold Release Agent - VALPLAST

Mold Release Agent - VALPLAST



Code	Name	Price	+
829-20213	Mold Release Spray	US\$ 20.16	+Info

PETROSIL SEPARATOR - ADS

PETROSIL SEPARATOR - ADS

Release agent for stone, plaster and epoxy. Save time by using two of the best known release agents, Liquid Petroleum Jelly and Pure Silicone. Petro-Sil sprays on evenly and leaves a barely detectable film. Non-toxic formula is fast drying. Save time with a guaranteed release on every model.



Code	Name	Price	+
118274	ADS Petro- Separator Agent 8oz. - #PS-08	US\$ 27.91	+Info
118275	ADS Petro- Sil Separator Agent 16oz (2x 8oz.) pk . - #PS-160	US\$ 46.27	+Info

PLASTICAL RELEASE AGENT- ADS

PLASTICAL RELEASE AGENT- ADS

Specifically designed to work with synthetic Die Materials. A terrific separator for resins, silicones, epoxies, and composites. Easy spray on or paint. Non-flammable formula leaves a fine, thin film. Spray on steel instruments to create a non-stick surface for composite material.



Code	Name	Price	+
102442	ADS Plactical Release Agent - Spray - 4oz #P510-3	US\$ 18.18	+Info
118292	ADS Plactical Release Agent 4oz. - #P510-3	US\$ 18.18	+Info
118293	ADS Plactical Release Agent 8oz. - #P510-4	US\$ 28.32	+Info
118294	ADS Plactical Release Agent 16oz. refill pt - #P510-5	US\$ 40.92	+Info

Separators & Diluent

Plasters Separators

N-SEP - ADS

N-SEP - ADS

Newly formulated for denture, partial and orthodontic laboratories. Provides clean separation between plaster and stone without film build-up or flaking. Non-flammable, non-hazardous. Applied as a dip or paint-on to stone or plaster. Brown tinted for visibility.



Code	Name	Price	+
118222	ADS N - Sep PLASTER SEPARATOR qt - #N210-6	US\$ 14.28	+Info
118223	ADS N - Sep PLASTER SEPARATOR gal - #N210-8	US\$ 31.47	+Info

Separators & Diluent

Porcelain Separators



DIAX-100 PORCELAIN SEPARATOR - ADS

DIAX-100 PORCELAIN SEPARATOR - ADS

Features:

- * Is specifically formulated to guarantee release of porcelain from stone model or die
- * Evaporates instantly and leaves a barely detectable, dry film
- * Build porcelain directly to die without any sticking or discoloration
- * Excellent for porcelain butt margins and bridge lapping
- * Non-contaminating, non-toxic and odorless



Code	Name	Price	+
117773	ADS DiAx - 100 PORCELAIN SEPARATOR qt - #DX-100-6	US\$ 101.81	+Info
117772	ADS DiAx - 100 PORCELAIN SEPARATOR 8oz. bottle - #DX-100-8	US\$ 43.65	+Info
117771	ADS DiAx - 100 PORCELAIN SEPARATOR 4oz. - spray - #DX-100-4	US\$ 27.03	+Info
117770	ADS DiAx - 100 PORCELAIN SEPARATOR 2oz. spray - #DX-100-2	US\$ 18.18	+Info

DIE SEALER/PORCELAIN SEPARATOR - ADS

DIE SEALER/PORCELAIN SEPARATOR - ADS

Features:

- * Creates a smooth, releasing surface for intricate direct porcelain build-up on stone dies
- * Burns out clean and will not cause discoloration
- * DIE SEALER closes the open pores of stone dies and protects from shipping or breaking providing a hard smooth working surface without distorting margins
- * PORCELAIN SEPARATOR releases porcelain from sealed dies and eliminates the need for messy tissue paper which causes the porcelain to dry out



Code	Name	Price	+
118302	ADS Porcelain Separator With Brush Cap 1 oz. - #D235-4	US\$ 15.11	+Info

118303	ADS Porcelain Separator Refill 4oz. - #D235-5	US\$ 44.15	+Info
--------	---	------------	-----------------------

Porcelain Separators - RENFERT

Porcelain Separators - RENFERT

Code	Name	Price	+
023-1552-0030	Picosep Isolation Agent-30ml	US\$ 19.15	+Info
105084	Iso-Stick each	US\$ 40.43	+Info
468958	IsoFix 2000 Spray Bottle Only	US\$ 9.88	+Info

Porcelain Separators - YETI

Porcelain Separators - YETI

Code	Name	Price	+
102101	Yeti Lube 18ml Bottle	US\$ 18.54	+Info

Separators & Diluent

Rubber Separators



FLEX-SEP - ADS

[FLEX-SEP - ADS](#)

Features:

- * Separates and protects, dries in seconds to a thin rubbery film
- * Paint on stone dies for clean separation from composites and acrylics
- * Paint on plastic teeth as a protective separation in denture processing
- * Paint on epoxy and metal dies as a die lube/separator
- * Can be peeled off most materials when you are done



KITS

Code	Name	Price	+
118004	ADS Flex Sep 1oz. BLUE/THINNER kit - #F896-12	US\$ 15.41	+Info
118005	ADS Flex Sep 1oz. red/1oz. THINNER kit - #F896-13	US\$ 15.41	+Info
118006	ADS Flex Sep 1oz.CLEAR/THINNER kit - #F896-11	US\$ 15.41	+Info

REFILL

BLUE

Code	Name	Price	+
118015	ADS Flex Sep refill BLUE 8oz. - #F896-42	US\$ 55.40	+Info

118012	ADS Flex Sep refill BLUE 4oz. - #F896-32	US\$ 36.96	+Info
118007	ADS Flex Sep BRUSH BTL - BLUE 1 oz. - #F896-22	US\$ 12.99	+Info

CLEAR

Code	Name	Price	+
118014	ADS Flex Sep refill CLEAR 8oz. - #F896-41	US\$ 55.40	+Info
118010	ADS Flex Sep refill CLEAR 4oz. - #F896-31	US\$ 36.96	+Info
118008	ADS Flex Sep BRUSH BTL - CLEAR 1 oz. - #F896-21	US\$ 12.99	+Info

RED

Code	Name	Price	+
118013	ADS Flex Sep refill RED 8oz. - #F896-43	US\$ 55.40	+Info
118011	ADS Flex Sep refill RED 4oz. - #F896-33	US\$ 36.96	+Info
118009	ADS Flex Sep BRUSH BTL - RED 1 oz. - #F896-23	US\$ 12.99	+Info

THINNER

Code	Name	Price	+
118016	ADS Flex Sep THINNER 1 oz. - #F896-5	US\$ 5.71	+Info

RUBBER SEP - SEPARATOR AND RELEASE AGENT - TAUB

RUBBER SEP - SEPARATOR AND RELEASE AGENT - TAUB

Rubber Sep is a latex liquid that creates a peelable film. It is a die spacer, release agent and separator for acrylic and composite resin restorations. Paint a thin coat onto a preparation and let it dry which takes a few minutes. The pink Rubber Sep dries translucent red and white Rubber Sep dries clear. The prep now has a 12 micron coating that allows the resin to cure onto the surface, making the cured resin easier to remove.



PINK

4 OZ

Code	Name	Price	+
122438	TAUB - Rubber Sep - PINK (red) - 4oz (96 Cc) - Wide Mouth Bottle - #10-2002	US\$ 61.76	+Info

8 OZ

Code	Name	Price	+
122440	TAUB - Rubber Sep - PINK - 8oz - #10-2003	US\$ 84.24	+Info

16 OZ

Code	Name	Price	+
122442	TAUB - Rubber Sep - PINK - 16oz - #10-2004	US\$ 121.25	+Info

32 OZ

Code	Name	Price	+
122444	TAUB - Rubber Sep - PINK - 32oz - #10-2006	US\$ 236.06	+Info

1 GALLON

Code	Name	Price	+
122443	TAUB - Rubber Sep - PINK - 1gal - #10-2005	US\$ 856.21	+Info

WHITE

6 CC

Code	Name	Price	+
122457	TAUB - Rubber Sep - WHITE - 6 Cc - w/ Brush And Cap - #10-2057	US\$ 0.00	+Info

4 OZ

Code	Name	Price	+
122452	TAUB - Rubber Separator - WHITE - 4oz (96 Cc) - Wide Mouth Bottle - #10-2052	US\$ 61.76	+Info

8 OZ

Code	Name	Price	+
122454	TAUB - Rubber Sep - WHITE - 8oz - #10-2053	US\$ 84.81	+Info

16 OZ

Code	Name	Price	+
122456	TAUB - Rubber Sep - WHITE - 16oz - #10-2054	US\$ 122.11	+Info

BLUE

4 OZ WITH BRUSH

Code	Name	Price	+
122462	TAUB - Rubber Sep - BLUE - 4oz (96 Cc) - Wide Mouth Bottle - #10-2072	US\$ 61.76	+Info

Rubber Separators - ADS

Rubber Separators - ADS

Code	Name	Price	+
118345	Rubber Dip Nin CoatClr 8oz	US\$ 25.89	+Info
118346	Rubber Dip Nin Coat red 8oz	US\$ 26.69	+Info

Separators & Diluent

Tin-Foil

Tin-Foil - KEYSTONE

Tin-Foil - KEYSTONE



Code	Name	Price	+
446884	Keystone Tin-Foil - .002" Thick x 6	US\$ 98.89	+Info

Tin-Foil - BUFFALO

Tin-Foil - BUFFALO



Code	Name	Price	+
424452	Tin-Foil .001" Roll, 1 pound	US\$ 91.17	+Info
424578	Tin Foil .002" x 6" Wide Roll	US\$ 94.86	+Info
424580	Tin-Foil .003" Roll, 1 pound	US\$ 91.17	+Info

Separators & Diluent
Zirconia Separators



Zirconia Separators - AMANN GIRRBACH

Zirconia Separators - AMANN GIRRBACH

Code	Name	Price	+
760561	Ceramill sep (lubricant)	US\$ 26.79	+Info

Catalog Index / Indice

Wax Separators	1
Acrylic separators	3
Die Sealers	9
Flexible Separators	10
Mold Release & Lubricant Silicones	10
Plasters Separators	13
Porcelain Separators	14
Rubber Separators	15
Tin-Foil	17
Zirconia Separators	18