

## Silicon wheels and points

### Silicon Amalgam Polishers



### Silicon Amalgam Polishers - SHOFU

Silicon Amalgam Polishers - SHOFU



#### Brownie

Brownie

Code	Name	Price	+
107-0413	Brownie Mini Points Fg	US\$ 24.97	+Info

#### Other Models

Other Models

Code	Name	Price	+
112041	Um Supergreenie Whl 12pk	US\$ 16.99	+Info

### Silicon wheels and points

## Silicon Hard Alloy Polishers



### Amalgam & Alloy Polishers - SHOFU

Amalgam & Alloy Polishers - SHOFU



Code	Name	Price	+
111632	HP Gold Polish Kit #0303	US\$ 32.98	+Info

### HP HARDIE POLISHERS - SHOFU

HP HARDIE POLISHERS - SHOFU

#### Advantage!!!

For hard alloys. These silicone polishers contain high-quality abrasive materials for polishing hard palladium, nickel-chrome, NonPrecious Alloys, Model Castings and Cobalt-chrome alloys. Use

SHAPE CATEGORY	FL2	PC2	WH6	Floppie	BP1	Unm. Wheel #8
SHAPE						
ISO #	030	050	120	220	100	230
<b>BLACK</b>						
HP (12 pk)	PN	-	0362	-	0366	0368 (6 pk)
CA (12 pk)	PN	0390	-	-	-	-
<b>BROWN</b>						
HP (12 pk)	PN	0361	0363	0365	0367	0369 (6 pk)
CA (12 pk)	PN	0391	-	-	-	-

Hardie Black for Contouring and Finishing. Use Hardie Brown For Polishing.

US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.

### Black

Code	Name	Price	+
111671	SHOFU HP Hardie BLACK Pc2 12pk - #0362	US\$ 27.20	+Info
111674	SHOFU HP Hardie BLACK Flame Floppie 12pk - #0366	US\$ 27.20	+Info
111676	SHOFU HP Hardie BLACK Bp1 6pk - #0368	US\$ 27.20	+Info
111679	SHOFU Unmounted Hardie BLACK Wheel #8 12pk - #0384	US\$ 16.22	+Info

### Brown

Code	Name	Price	+
111670	SHOFU Hp Hardie BROWN Flame Fl2 12pk - #0361 Mfg #0361	US\$ 27.20	+Info
111672	SHOFU HP Hardie BROWN Point PC2 12pk - #0363	US\$ 27.20	+Info
111673	SHOFU HP Hardie BROWN Wheel WH6 12pk - #0365	US\$ 27.20	+Info
111675	SHOFU HP Hardie BROWN Flame Floppie 12pk - #0367	US\$ 27.20	+Info
111677	SHOFU HP Hardie BROWN Bp1 6pk - #0369	US\$ 27.20	+Info

### Silicon wheels and points

## Silicone Points & Wheels



### POLISOFT - QUICK AND EASY PRE-POLISHING - RENFERT

#### POLISOFT - QUICK AND EASY PRE-POLISHING - RENFERT

#### Function & Performance

Effective polishing of precious and non-precious metal alloys.

#### Advantages

- \* Optimum alternative to rubber or silicone polishers.
- \* Quickly smoothes all metal surfaces.
- \* Virtually no dust or odour produced.



Code	Name	Price	+
023-91-0000	RENFERT - Polisoft Wheel 22X3mm-50 Pcs- #910-000 #910000 Pano Pulir Control Rueda	US\$ 54.27	+Info
023-90-0000	RENFERT - Polisoft A Wheel 22X3mm-50 Pcs- #900-000 #900000 Pano Pulir Control R...	US\$ 42.56	+Info

### SILICONE POINT POLISHERS SHOFU

#### POLISHERS - SHOFU



Use Brownies- Greenies - Super Greenies - Hardies to achieve a smooth surface, Greenies for a lustrous polish and Super greenies for a super-polish.  
Specially selected Brown and Green Silicone Points in HP shank for polishing and super polishing of precious alloys.

#### Brownie

Code	Name	Price	+
111603	SHOFU HP BROWNIE Wheel WH6-12pk - #0273	US\$ 26.28	+Info
111602	SHOFU HP BROWNIE Point PC2-12pk - #0272	US\$ 26.28	+Info
111601	SHOFU HP BROWNIE Knife Kn7-12pk - #0271	US\$ 26.28	+Info
111713	SHOFU HP BROWNIE Flame Floppie-12pk - #0458	US\$ 26.28	+Info

#### Greenie

Code	Name	Price	+
111714	SHOFU Hp GREENIE Flame Floppie 12pk - #0459	US\$ 24.70	+Info
111606	SHOFU HP GREENIE Wheel WH6 12pk - #0276	US\$ 25.75	+Info
111605	SHOFU HP GREENIE Point PC2 12pk - #0275	US\$ 26.28	+Info
111604	SHOFU HP GREENIE Knife Kn7 12pk - #0274	US\$ 26.28	+Info

#### Super Greenie

Code	Name	Price	+
111609	SHOFU HP SUPERGREENIE Wheel WH6 12pk - #0279	US\$ 26.28	+Info
111608	SHOFU HP SUPERGREENIE Point PC2 12pk - #0278	US\$ 26.28	+Info
111607	SHOFU HP SUPERGREENIE Knife Kn7 12pk - #0277	US\$ 26.28	+Info

### Silicone Points & Wheels - JSP

Silicone Points & Wheels - JSP

Code	Name	Price	+
100476	Blue Pnt Fine 22x 3mm 100p	US\$ 31.49	+Info

### Silicone Points & Wheels - RENFERT

Silicone Points & Wheels - RENFERT

Code	Name	Price	+
417438	Silicone Polishers 22x3.2mm 100pk	US\$ 81.44	+Info

**FREE SHIPPING**  
FOR ALL DOMESTIC ACCOUNTS!



**US  
DENTAL  
DEPOT**

+1 954 874 6325   
www.usdentaldepot.com 

### SILICONE POINTS AND WHEELS - KEYSTONE ONE

#### SILICONE POINTS AND WHEELS - KEYSTONE

Our Deluxe Silicone Points and Wheels are specially formatted with the finest grade of silicone. These items provide excellent results and are long-lasting, making them a cost-effective choice.

Available in Square Edge, Knife Edge and Pointed.

#### These points and wheels also come in various grits:

- \* Black is a medium grit for all metals
- \* White is a coarse grit made for all metals, amalgam, porcelain and composites.
- \* Blue is a fine grit for final polish of porcelain or metals.
- \* Pink is an X-fine grit works well for porcelain.

#### Features:

- \* Square Edge - 7/8" x 1/8" (22.2 x 3.1mm)
- \* Knife Edge - 7/8" (22.2mm)
- \* Made of fine-grade silicone
- \* Available in 100/box



### SILICONE WHEELS AND POINTS 100PK

#### Silicone Knife Edge Wheels

Code	Name	Price	+
034-1909005	KEYSTONE Deluxe Silicone Knife Edge Wheels, BLACK, Medium Grit , For All Metal...	US\$ 59.63	<a href="#">+Info</a>
034-1909020	KEYSTONE Deluxe Silicone Knife Edge Wheels, WHITE, Coarse Grit , For All Metal...	US\$ 59.63	<a href="#">+Info</a>
034-1909035	KEYSTONE Deluxe Silicone Knife Edge Wheels, BLUE, Fine Grit , For Final Polish...	US\$ 59.63	<a href="#">+Info</a>

#### Silicone Points

Code	Name	Price	+
034-1909010	KEYSTONE Deluxe Silicone Points, BLACK, For All Metals, Medium Grit , box of 1...	US\$ 59.63	<a href="#">+Info</a>
034-1909025	KEYSTONE Deluxe Silicone Points, WHITE, Coarse Grit , For All Metals, Amalgam ...	US\$ 52.47	<a href="#">+Info</a>
034-1909055	KEYSTONE Deluxe Silicone Points, PINK, X-Fine Grit , For Porcelain, box of 100...	US\$ 59.63	<a href="#">+Info</a>
034-1909040	KEYSTONE Deluxe Silicone Points, BLUE, Fine Grit , For Final Polish of Porcela...	US\$ 47.70	<a href="#">+Info</a>

#### Silicone Square Edge Wheels

Code	Name	Price	+
034-1909000	KEYSTONE Deluxe Silicone Square Edge Wheels, BLACK, Medium Grit , For All - #1...	US\$ 47.70	<a href="#">+Info</a>
034-1909015	KEYSTONE Deluxe Silicone Square Edge Wheels, WHITE, Coarse Grit , For All - #1...	US\$ 59.63	<a href="#">+Info</a>
034-1909030	KEYSTONE Deluxe Silicone Square Edge Wheels, BLUE, Fine Grit , For Final Polis...	US\$ 57.24	<a href="#">+Info</a>

034-1909045

KEYSTONE Deluxe Silicone Square Edge Wheels, PINK, X-Fine Grit , For Porcelain...

US\$ 59.63

[+Info](#)

### SILICONE WHEELS AND POINTS 10PK ( Keystone )

#### Silicone Knife Edge Wheels

Code	Name	Price	+
034-1909020-10	KEYSTONE WHITE Coarse Grit - 22mm x 3mm - Deluxe Silicone Wheels Knife edge - ...	US\$ 6.35	<a href="#">+Info</a>
034-1909035-10	KEYSTONE BLUE fine Grit - 22mm - Deluxe Silicon Wheel Knife edge - 10pk K #190...	US\$ 5.78	<a href="#">+Info</a>

#### Silicone Points 10pk ( Keystone )

Code	Name	Price	+
034-1909010-10	KEYSTONE - BLACK Medium Grit Deluxe Silicone Points - 10pk 1909010-10	US\$ 5.19	<a href="#">+Info</a>

#### Silicone Square Edge Wheels

Code	Name	Price	+
034-1909015-10	KEYSTONE WHITE Coarse Grit - 22mm x 3mm - Deluxe Silicone Square Edge - 10pk	US\$ 6.35	<a href="#">+Info</a>
034-1909030-10	KEYSTONE BLUE fine Grit - 22mm x 3mm - Deluxe Silicon Wheel Square Edge - 10p...	US\$ 6.41	<a href="#">+Info</a>

## Catalog Index / Indice

Silicon Amalgam Polishers .....	1
Silicon Hard Alloy Polishers .....	1
Silicone Points & Wheels .....	2