

Solders



Chrome Cobalt Welding Rods

Chrome Cobalt Welding Rods - E

Chrome Cobalt Welding Rods - E

Code	Name	Price	+
100965	E-Chrome-2 - Solder 1oz	US\$ 21.60	+Info

SOLDERING RODS - BEGO

SOLDERING RODS - BEGO

- * Melting interval 1100 - 1150 C
- * Soldering temperature 1180 C



US DENTAL DEPOT Inc is an Authorized BEGO Dealer in US.

Code	Name	Price	+
10-52622	BEGO Soldering Rods for Wirobond - 8pk #52622 Melting interval: 1150-1100 C- Sol...	US\$ 60.83	+Info
10-52520	BEGO Soldering Rods Co-Cr - 5pk #52520 Melting interval: 1150-1100 C- Soldering ...	US\$ 50.16	+Info

Solders

Fluxes Solders



CASTING FLUX - ADS

CASTING FLUX - ADS

Features:

- * High efficiency, reducing formula
- * No carbon particles to form inclusions
- * Fluoride free formula
- * For gold and silver alloy castings
- * Wets the molten alloy surface and increases the fluidity of metal, thus minimizing porosity
- * Flux cleans up old metal for recasting, reduces metal oxides and inhibits formation of new oxides
- * Use to glaze the interior of crucibles



ALL HIGH

Code	Name	Price	+
118067	ADS Flux - All High FUSING SOLDERING 1/2oz. - #499911	US\$ 15.30	+Info

BROWN FLUX

Code	Name	Price	+
118136	ADS Flux - BROWN Sd - 1 HIGH TEMP. SYRINGE 20 cc paste - #F110-2	US\$ 18.55	+Info
118069	ADS Flux - BROWN Sd - 1 HIGH TEMP. JAR 40 cc paste - #F110-3	US\$ 21.88	+Info

CASTING

Code	Name	Price	+
118140	ADS Flux - Casting WHITE powder 3oz.powder - #C400-2	US\$ 15.81	+Info
118141	ADS Flux - Casting w/ Fluoride WHITE powder 4oz. powder - #C400-3	US\$ 15.81	+Info

DEOXIDIZER

Code	Name	Price	+
118139	ADS Flux - With Deoxidizer BROWN powder 4oz. powder - #C400-1	US\$ 15.81	+Info

WHITE FLUX

Code	Name	Price	+
118138	ADS Flux - WHITE Sd - 2 NON PRECIOUS SYRINGE 20 cc paste - #F111-3	US\$ 15.81	+Info
118137	ADS Flux - WHITE Sd - 2 NON PRECIOUS JAR 40 cc paste - #F111-4	US\$ 19.46	+Info

Fluxes Solders - BEGO

Fluxes Solders - BEGO

Code	Name	Price	+
010-52530	Minoxid Flux - Bottle - 80gr	US\$ 55.50	+Info
10-52526	Wiromelt Melting Powder 80g	US\$ 86.84	+Info

Fluxes Solders - DENTALALLOY

Fluxes Solders - DENTALALLOY

Code	Name	Price	+
102341	SuperFlux Jar 4oz	US\$ 23.15	+Info

Solders

Heat Protectors & Blockers



HEAT BARRIER GEL - ADS

HEAT BARRIER GEL - ADS

Ready to use, asbestos-free Heat Barrier Gel protects acrylic and plastic by absorbing the heat during soldering and brazing procedures. Prevents warping and discoloration. Cover the surface next to heated area, wash off with water when work is done. Use to 2500F. Will not lift or curl away from the edges of the heated area. Squeeze bottle with dispensing spout



Code	Name	Price	+
118078	ADS Heat Barrier Gel 4oz. - #H703-3	US\$ 18.55	+Info
118079	ADS Heat Barrier Gel 8oz. - #H703-4	US\$ 26.75	+Info

Heat Protectors & Blockers - ADS

Heat Protectors & Blockers - ADS

Code	Name	Price	+
118080	Heat Transfer Oil #20 qt	US\$ 27.91	+Info

Solders

Non Precious Ceramic Solders



Non Precious Ceramic Solders - BEGO

Non Precious Ceramic Solders - BEGO

Code	Name	Price	+
010-52625	Wiolot Soldering Rods 4gr	US\$ 60.83	+Info

Non Precious Ceramic Solders - BESQUAL

Non Precious Ceramic Solders - BESQUAL

Code	Name	Price	+
102340	Solder Nickel Chrome 1oz	US\$ 14.40	+Info

Non Precious Ceramic Solders - DENTALALLOY

Non Precious Ceramic Solders - DENTALALLOY

Code	Name	Price	+
100736	SUPER-1 PRE-SOLDER 1oz	US\$ 24.00	+Info

Solders

Soldering Accessories



PLACE-IT - ADS

PLACE-IT - ADS

Safe and asbestos free. Just a dab holds your assembly during soldering, welding, and firing. Place-It will not rob heat from joints, insuring an excellent bond. Heat shield protects acrylic during soldering. Use to repair furnace parts and to create custom supports for laminate veneers. Use to 3200F. No shrinkage on setting. Sets hard in any desired form or shape. Create custom sagger trays, pegs, furnace fixtures, and custom refractory supports for porcelain. Instantly secure multiple unit bridgework for soldering and welding. Apply easily to soldering board or coping interior with handy syringe. Hardens quickly. Ready to use -no mixing, setting or high heat required. Easily change shades, characterize, or add on to your porcelain restorations. Just a dab of Place-It



supports laminate veneers, all ceramic crowns, inlays, and onlays in position for refiring. Use Place-It instead of investment. Apply directly to your sagger tray creating a refractory die in seconds. No degassing, no mixing, no waiting. Perfect for your busy lab, saves hours of work with investment.

TYPE 2

Code	Name	Price	+
118283	ADS Place-It Type 2 TYPE 2 small - #P702-6	US\$ 16.44	+Info
118282	ADS Place-It Type 2 TYPE 2 medium - #P702-5	US\$ 22.36	+Info
118281	ADS Place-It Type 2 TYPE 2 large - #P702-4	US\$ 30.25	+Info
118284	ADS Place-It Type 2 TYPE 2 PRO kit - #P702-11	US\$ 111.29	+Info

TYPE 3

Code	Name	Price	+
118286	ADS Place-It Type 2 TYPE 3 small - #P702-8	US\$ 19.66	+Info
118285	ADS Place-It Type 2 TYPE 3 medium - #P702-9	US\$ 26.25	+Info

Soldering Accessories - BUFFALO

Soldering Accessories - BUFFALO

Code	Name	Price	+
425275	Buffalo No. 2A Solde	US\$ 9.35	+Info
425276	Buffalo No. 2 Solder	US\$ 8.81	+Info

SOLDERING BOARD - ADS

SOLDERING BOARD - ADS

Economical protection for work tables or benches. Withstands quick heating and cools in seconds. Use in front of furnace to cool sagger trays. Use under bunsen burner to prevent scorching of formica or wood counters.



Code	Name	Price	+
118389	ADS Soldering Board 6in x 6in ea - #S 701-2	US\$ 22.29	+Info
118390	ADS Soldering Board 6in x 12in ea - #S 701-3	US\$ 37.95	+Info
118391	ADS Soldering Board 12in x 12in ea - #S 924-6	US\$ 52.55	+Info

TURN N' SOLDER TABLE - ADS

TURN N' SOLDER TABLE - ADS

Revolutionary Turn N' Solder is an economical soldering board, designed for light duty, mounted in a stable, heat-resistant turntable. Super convenience and soldering ease. Pin your workpiece in place or use Place-It for support. Use Turn N' Solder directly in front of



porcelain funace for cooling and manipulation of sagger trays. Resistant to molten metals and most chemicals. economical lightweight ceramic. Temperature tolerance of 3000F. Withstands rapid heating and cools to the touch in seconds.

Code	Name	Price	+
118435	ADS TurnNTable Solder Frame 6in x 6in x 0.5in boards kit - #S705-5	US\$ 53.59	+Info
118434	ADS TurnNTable Solder Frame 6in x 6in x 0.5in board table - #T703-1	US\$ 39.98	+Info
118436	ADS Turntable support frame ea - #W 946-3-25	US\$ 26.82	+Info

Catalog Index / Indice

Chrome Cobalt Welding Rods	1
Fluxes Solders	1
Heat Protectors & Blockers	2
Non Precious Ceramic Solders	3
Soldering Accessories	3