

## Vacuum Forming Thermoplastic Materials

### Base Strenghteners



### Base Strenghteners - KEYSTONE

Base Strenghteners - KEYSTONE

Code	Name	Price	+
801100	Base Strenghtener .100 25pk	US\$ 58.27	<a href="#">+Info</a>

## Vacuum Forming Thermoplastic Materials

### Bleaching / Whitening Dual laminates



### Bleaching / Whitening Dual laminates - KEYSTONE

Bleaching / Whitening Dual laminates - KEYSTONE

Code	Name	Price	+
474316	Dual BleachLamin. w/Foam round 50pk	US\$ 78.68	<a href="#">+Info</a>

### BLEACHING LAMINATES - KEYSTONE

#### BLEACHING LAMINATES - KEYSTONE

Its become an acceptable dental technique to bleach teeth using carbamide peroxide and a bleaching tray. The Bleaching Laminate is ideally suited for making the bleaching tray and withstanding bleaching products. It is thin, only .020 (.5mm) inches thick once its formed over the model. It has a foam liner to absorb the bleach, thus decreasing the bleaching time for the patient. The Bleaching Laminate is our best-selling item for bleaching processes.



#### Features:

- \* Clear material
- \* .020 (.5mm) thick
- \* Available in 5? x 5? (127 x 127mm) sheets
- \* Also available in round
- \* Available in 10 or 50 pieces per pack

Code	Name	Price	+
034-9505780-1	KEYSTONE Dual Bleaching Laminates with Foam - 10pk - ( With foam liner ) K#96057...	US\$ 79.80	<a href="#">+Info</a>
034-9505780	PRO-FORM - Dual Bleaching Laminates with Foam - 50pk - With foam liner #9605780 ...	US\$ 78.68	<a href="#">+Info</a>

## Vacuum Forming Thermoplastic Materials

### Coping Material



### Coping Material - BUFFALO

Coping Material - BUFFALO

Code	Name	Price	+
424454	.020" Coping and Tem	US\$ 13.30	+Info

### COPING MATERIAL - KEYSTONE

#### COPING MATERIAL - KEYSTONE

These cloudy, stiff polypropylene sheets are excellent for fabricating temporary crowns, bridges and copings. The material gives good definition and will not adhere to acrylics.



#### Specifications:

- \* Excellent for Temporary Bridges
- \* Get Good Definition
- \* Vacuum Forming Material
- \* Three Thicknesses : 0.20in .030in .040in
- \* Available in 5" x 5" ( 127x127mm ) sheets
- \* Round
- \* Package: 50 or 625 pieces per pack

**Us Dental Depot Is An Authorized Keystone Dealer In Us**

Code	Name	Price	+
419987	PRO-FORM - Coping Material Round - .030in -0.7mm) , Round Sheets 125mm, - 50pk ...	US\$ 14.73	+Info

#### 50pk

Code	Name	Price	+
157017	KEYSTONE Coping Material - .020in ( 0.5mm ) - 5x5 sheets - 50pk - ( K#9613470 )...	US\$ 11.86	+Info
101160	KEYSTONE Coping Material - .030in (0.7mm) - 50pk - 5x5 sheets - ( K#9613480 ) ( ...	US\$ 12.75	+Info
700015	KEYSTONE Coping Material .040in (1mm) - 5x5 sheets - 50pk - ( K#9613510 ) ( 5X5in )	US\$ 17.15	+Info

#### 625pk

Code	Name	Price	+
157018	PRO-FORM - Coping Material - .020in - 0.5mm - 5x5 sheets - 625p #9613540 #70001...	US\$ 122.36	+Info
700012	PRO-FORM - Coping Material - .030in - 0.7mm - 5x5 sheets - 625pk - 5X5in #96135...	US\$ 152.81	+Info
700017	PRO-FORM - Coping Material .040in - 1mm - 5x5 sheets - 625p #9613530 -	US\$ 209.44	+Info

### CROWN & BRIDGE - KEYSTONE

#### CROWN & BRIDGE - KEYSTONE

Our crown & bridge material is thin, clear rigid material that works great for temporary bridges. The material is easy to work with while it does not bond to acrylic.



#### Specifications:

- \* Acrylic Does Not Bond to Material
- \* Packages of 50 or 625 Sheets
- \* 5 x 5 inch
- \* For Temporary Bridges
- \* Vacuum Forming Material

Code	Name	Price	+
034-1570160	PRO-FORM - Crown & Bridge Material .020in 0.5mm 625pk 5x5 Sheet #9613700	US\$ 296.58	+Info
034-1570150	KEYSTONE Crown & Bridge Material - .020in ( 0.5mm ) - 50pk - 5x5 sheets - ( K#96...	US\$ 27.07	+Info
034-1570150-1	KEYSTONE Crown & Bridge Material - .020in ( 0.5mm ) - 25pk - 5x5in sheets - ( K#...	US\$ 14.37	+Info

### Vacuum Forming Thermoplastic Materials

## Cutters and Knives



### Cutters and Knives - BUFFALO

#### Cutters and Knives - BUFFALO

Code	Name	Price	+
008-80500	Thermaknife Intro 120v AC	US\$ 281.05	+Info
122329	ThermaPlier - 110v	US\$ 285.12	+Info
424445	ThermaKnife Repl Pen w/ Knife Edge	US\$ 49.84	+Info
424508	ThermaKnife Replac Cord Pen Only	US\$ 35.87	+Info
424522	ThermaPlier Hot tipped thermal plie	US\$ 285.12	+Info
424523	ThermaKnife Intro 220v AC	US\$ 264.00	+Info
424586	ThermaKnife Repl Pen w/ Concave Tip	US\$ 54.24	+Info
424636	ThermaKnife Repl Pen w/ Waxing Tip	US\$ 49.84	+Info
425201	ThermaKnife Console Only 120V AC.	US\$ 215.28	+Info

### Cutters and Knives - KEYSTONE

#### Cutters and Knives - KEYSTONE

Code	Name	Price	+
035-P360-001	Proform Accessories Utility Sc	US\$ 15.45	+Info

### Vacuum Forming Thermoplastic Materials

## Dual laminates 4 Bruxism



### Dual laminates 4 Bruxism - KEYSTONE

Dual laminates 4 Bruxism - KEYSTONE

Code	Name	Price	+
112003	NightGuards Soft/Soft 12pk 3mm	US\$ 57.45	<a href="#">+Info</a>

### Dual laminates 4 Bruxism - PRO

Dual laminates 4 Bruxism - PRO

Code	Name	Price	+
420169	NightGuards Rnd Soft/Soft 12pk 3mm	US\$ 72.66	<a href="#">+Info</a>

## NIGHTGUARDS & E-GASKET DUAL LAMINATES - KEYSTONE

### NIGHTGUARDS & E-GASKET DUAL LAMINATES - KEYSTONE

#### Features:

This laminate or layered soft/soft material was made specifically for the making of niteguards. The layering gives the thermoplastic a more abrasion-resistant characteristic that will withstand grinding better than traditional soft sheet materials.

- \* Available in 5"x 5" (127 x127mm) sheets
- \* .120 (3mm) thickness
- \* Also available in round
- \* Available in 12/Pkg
- \* Available in three thicknesses 0.080 (2mm) , 0.120 (3mm), and 0.200 (5mm)
- \* Available in 5"x 5" (127x127mm) sheets
- \* Also available in round
- \* Available in 12/Pkg

#### Dual Laminates (E-Gasket)

The Pro-form Dual Laminate (E-gasket) is a unique combination of two materials formed into a single laminated sheet. It was made with bruxing cases in mind. The soft liner assures patient comfort and the hard surface assures long-lasting wear. Also, acrylic will bond to the hard surface to allow the dental technician to modify the occlusal surface. Available in 3 different thicknesses: 0.080 (2mm), 0.120, (3mm), and .200 (5mm).



### HARD/SOFT

Code	Name	Price	+
113005	PRO-FLEX- PRO-FORM E-Gasket 5x5in .080 (2mm) thick, CLEAR, ( Dual Laminate 2...	US\$ 57.45	<a href="#">+Info</a>
113003	PRO-FLEX - PRO-FORM Dual Laminate 5x5in .120 (3mm) thick, ( E-Gasket 50 Hard...	US\$ 57.45	<a href="#">+Info</a>
113004	PRO-FLEX - E-Gasket - Dual Laminates 5mm 5x5 in (E-Gasket 50 Hard / 50 Soft ) ...	US\$ 80.36	<a href="#">+Info</a>

### ROUND

Code	Name	Price	+
419978	PRO-FLEX - E-Gasket - Dual Laminates .080 (2mm) E-Gasket CLEAR, #9604700R1	US\$ 72.66	+Info
419979	PRO-FLEX - E-Gasket - Dual Laminates .080 -2mm ( 50/50 ) E-Gasket CLEAR, #96047...	US\$ 72.66	+Info
419980	PRO-FLEX - E-Gasket - Dual Laminates .120 (3mm) Round 125mm ( E-Gasket 50 Hard ...	US\$ 72.66	+Info
419981	PRO-FLEX - E-Gasket - Dual Laminates .120 (3mm) ( E-Gasket 50 Hard / 50 Soft ) C...	US\$ 72.66	+Info
419982	PRO-FLEX - E-Gasket - Dual Laminates .200 5mm Round 125mm ( E-Gasket 50 Hard / ...	US\$ 80.36	+Info
419983	PRO-FLEX - E-Gasket - Dual Laminates .200 -5mm Round 120mm (E-Gasket 50 Hard / 5...	US\$ 80.36	+Info

### Vacuum Forming Thermoplastic Materials

## Mouthguards

## FUN MOUTHGUARDS - PRO-FORM

### FUN MOUTHGUARDS - PRO-FORM

#### Features:

PRO-FORM Tie Dye Mouthguards designs are a new line added to its extensive line of laminated mouthguards. The Pro-form mouthguard tensile strength of the laminate is excellent due to the laminate's two layers of materials making up the mouthguard. Density is maintained at the pressure laminating process so shrinkage of the mouthguard is uniformly controlled.

All Pro-form Mouthguards are made using the same laminating procedure to insure the same Pro-form quality.

- \* Numerous color options to meet the needs of your patients
- \* 5x5 - 6pk - 0.150- thick
- \* Available in 6 different styles



### 6 PK

Code	Name	Price	+
034-7954000	PRO-FORM Tie-Dye Multi-color Mouthguard Laminate - Patriot, 6/Pk. 5x5in .157in...	US\$ 29.78	+Info
113909	PRO-FORM Tie-Dye Multi-color Mouthguard Laminate - Wilderness, 6/Pk. 5x5in - #...	US\$ 29.78	+Info
113910	PRO-FORM Tie-Dye Multi-color Mouthguard Laminate - Midnight, 6/Pk. 5x5in .157i...	US\$ 29.78	+Info
113911	PRO-FORM Tie-Dye Multi-color Mouthguard Laminate - Shamrock, 6/Pk. 5x5in .157i...	US\$ 29.78	+Info
113912	PRO-FORM Tie-Dye Multi-color Mouthguard Laminate - Flamingo, 6/Pk. 5x5in .157i...	US\$ 29.78	+Info
113913	PRO-FORM Tie-Dye Multi-color Mouthguard Laminate - Rainbow, 6/Pk. 5x5in .157in...	US\$ 29.78	+Info
034-7974000	PRO-FORM Tie-Dye Multi-color Mouthguard Laminate, ASSORTED - Patriot 12pk - #7...	US\$ 29.78	+Info

## Mouthguards - BUFFALO

### Mouthguards - BUFFALO

Code	Name	Price	+
424959	.020" Coping and Tem	US\$ 80.54	+Info
424960	.020" Coping and Tem	US\$ 147.68	+Info

425197	BUFF-TUFF #49 pliers, Serrated, tw	US\$ 8.40	<a href="#">+Info</a>
--------	------------------------------------	-----------	-----------------------

### GLITTER MOUTHGUARDS - PRO-FORM

#### PRO-FORM GLITTER GUARDS

Pro-form athletic mouthguard line got a little more flash with Glitter Guards. The Glitter Guard material adds a bit of sparkle to the Pro-form lauded line, which delivers protection, affordability and style.

- \* .160 (4mm) Thick
- \* 5"x 5" (127x127mm ) squares
- \* Package of 6



By laminating two layers of soft EVA together, the material has excellent tensile strength. All of these features make the Pro-form mouthguard laminate the material of choice for custom-made mouthguards

Code	Name	Price	+
034-7965070	PRO-FORM Glitter Guard Mouthguard Laminate - ASSORTED: BLUE, YELLOW, Purple - #...	US\$ 29.78	<a href="#">+Info</a>
034-7965010	PRO-FORM Glitter Guard Mouthguard Laminate - BLUE, 6/Pk. 5x5in .157in (4mm) -...	US\$ 29.78	<a href="#">+Info</a>
034-7965020	PRO-FORM Glitter Guard Mouthguard Laminate - ORANGE, 6/Pk. 5x5in .157in (4mm)...	US\$ 29.78	<a href="#">+Info</a>
034-7965030	PRO-FORM Glitter Guard Mouthguard Laminate - YELLOW, 6/Pk. 5x5in .157in (4mm)...	US\$ 29.78	<a href="#">+Info</a>
034-7965040	PRO-FORM Glitter Guard Mouthguard Laminate - GREEN, 6/Pk. 5x5in .157in (4mm) ...	US\$ 29.78	<a href="#">+Info</a>
034-7965050	PRO-FORM Glitter Guard Mouthguard Laminate - Purple, 6/Pk. 5x5in .157in (4mm)...	US\$ 29.78	<a href="#">+Info</a>
034-7965060	PRO-FORM Glitter Guard Mouthguard Laminate - PINK, 6/Pk. 5x5in .157in (4mm) -...	US\$ 29.78	<a href="#">+Info</a>

### Mouthguards - KEYSTONE

Mouthguards - KEYSTONE

Code	Name	Price	+
470596	Mthgrd TriColor Gr/Cl/Yw 12 pk	US\$ 2.44	<a href="#">+Info</a>

### Mouthguards - PLASDENT

Mouthguards - PLASDENT

Code	Name	Price	+
424366	Plasdent Dental Mouth Guard, Adult,	US\$ 8.15	<a href="#">+Info</a>
424367	Plasdent Dental Mouth Guard, Child,	US\$ 8.10	<a href="#">+Info</a>

### Mouthguards - PRO-FORM

Mouthguards - PRO-FORM

Code	Name	Price	+
115434	PRO-FORM Complete Sports Pack	US\$ 144.43	<a href="#">+Info</a>
421144	Youth Mouthguard - White -300pk	US\$ 264.74	<a href="#">+Info</a>

## ONE COLOR MOUTHGUARD - LAMINATED - PRO-FORM

### ONE COLOR MOUTHGUARD - LAMINATED - PRO-FORM

#### Pro-form Mouthguard Laminates

Pro-form manufactures the only pre-made laminated Pro-form mouthguard material in the industry. This material has long been considered the leader in custom-made mouthguards. By laminating two layers of soft EVA together, the laminate has excellent tensile strength. All of these features make the Pro-form mouthguard laminate the material of choice for custom-made mouthguards.

#### Features:

- \* Available in .160 (4mm) 5"x5" (127 x 127mm) squares
- \* Also available in round 120mm and 125mm diameter
- \* Available in 12/Pkg and 1 Pk
- \* Available Colors: Clear, Orange, Yellow, Maroon, White, Green, Black, Red, Blue, Purple, Pink
- \* Available in assorted color packs too

These custom-made athletic Pro-form mouthguards are proven to reduce the incidence of tooth breakage by properly distributing the stress of a hard impact to the entire length of the tooth. A double layer of laminated sheet vinyl and a lingual plate imbedded behind the incisors gives the piece long-lasting durability. The Pro-form mouthguard maintains its form because of the heat and pressure laminating process. The tensile strength is excellent due to the two layers of laminate while the density is maintained during the pressure laminated process for controlled, uniformed shape. Overall, the risk of injuries is lessened by about 90% when wearing a custom-made mouthguard.

Used at all levels of competition, Pro-form has proved itself as a top performer on the field when athletes need it the most

#### The following is a list of direct benefits of wearing a custom Pro-form mouthguard

- \* Allows athletes to breathe during competition and enhances performance.
- \* Prevents tongue, lips and cheeks from lacerations from sharp ends of the upper teeth.
- \* Lessens risk of damage to the back of the teeth following a blow delivered to the underside of the lower jaw. Without a mouthguard, such impacts can



cause cusp fractures and complete root fractions.

\* Lessens the risk of concussions following an impact around the mouth. This is due to the prevention of full posterior translation of the condyles to the base of the skull.

\* Lessens the risk of jaw fractures by absorbing the energy of a traumatic blow.

\* Improves confidence of players so the athletes can concentrate more on the efforts and execution of their sport

\* Can be used for multiple sports.\* With proper care, the mouthguard can be used for more than a year.

\* Extremely beneficial for athletes who are undergoing orthodontic treatment.

\* Myriad color options to meet the needs and wants of patients.

### 1 PK SQUARE 5 X 5 in

Code	Name	Price	+
9596920-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate CLEAR.160in ( 4mm ) 1pk 5x5 Sheet (...)	US\$ 10.27	+Info
9598060-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate BLUE .160in ( 4mm ) 1pk 5x5 Sheet (...)	US\$ 11.29	+Info
9598730-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate BLACK .160in ( 4mm ) 1pk 5x5 Sheet ...	US\$ 8.78	+Info
9600800-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate GREEN .160in ( 4mm ) 1pk 5x5 Sheet ...	US\$ 8.78	+Info
9600820-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate Maroon .160in ( 4mm ) 1pk 5x5 Sheet...	US\$ 8.78	+Info
9600850-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate ORANGE.160in ( 4mm ) 1pk 5x5 Sheet ...	US\$ 8.78	+Info
9600870-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate Purple .160in ( 4mm ) 1pk 5x5 Sheet...	US\$ 8.78	+Info
9600910-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate RED .160in ( 4mm ) 1pk 5x5 Sheet ( ...)	US\$ 8.78	+Info
9601020-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate WHITE .160in ( 4mm ) 1pk 5x5 Sheet ...	US\$ 11.22	+Info
9601040-1	PRO-FORM - SINGLE-COLOR Mouthguards Laminate YELLOW .160in ( 4mm ) 1pk 5x5 Sheet...	US\$ 8.78	+Info

### 12 PK SQUARE 5 X 5 in

Code	Name	Price	+
9600180	PRO-FORM SINGLE-COLOR Mouthguards Laminate ASSORTED .160in - 4mm 12pk 5x5 Sheet...	US\$ 69.09	+Info
9596920	PRO-FORM SINGLE-COLOR Mouthguard Laminate - CLEAR, 12/Box. 5x5in .160in (4mm)...	US\$ 68.39	+Info
9598730	PRO-FORM SINGLE-COLOR Mouthguard Laminate - BLACK, 12pk. 5x5in .160in (4mm) - ...	US\$ 69.79	+Info
9600800	PRO-FORM SINGLE-COLOR Mouthguard Laminate - GREEN, 12/Box. 5x5in .160in (4mm)...	US\$ 69.79	+Info
9600820	PRO-FORM SINGLE-COLOR Mouthguard Laminate - Maroon, 12pk. 5x5in .160in (4mm) ...	US\$ 69.79	+Info
9600870	PRO-FORM SINGLE-COLOR Mouthguard Laminate - Purple, 12pk. 5x5in .160in (4mm) ...	US\$ 69.79	+Info
9600886	PRO-FORM SINGLE-COLOR Mouthguard Laminate - PINK, 12/Box. 5x5in .160in (4mm) ...	US\$ 69.79	+Info
9600910	PRO-FORM SINGLE-COLOR Mouthguard Laminate - RED, 12/Box. 5x5in .160in (4mm) T...	US\$ 69.79	+Info
9601020	PRO-FORM SINGLE-COLOR Mouthguard Laminate - WHITE, 12/Box. 5x5in .160in (4mm)...	US\$ 69.79	+Info
9601040	PRO-FORM SINGLE-COLOR Mouthguard Laminate - YELLOW, 12pk. 5x5in .160in (4mm) ...	US\$ 55.13	+Info
9598060	PRO-FORM SINGLE-COLOR Mouthguard Laminate - BLUE, 12/Box. 5x5in .160in (4mm) ...	US\$ 69.79	+Info
9600850	PRO-FORM SINGLE-COLOR Mouthguard Laminate - ORANGE, 12pk. 5x5in .160in (4mm) ...	US\$ 69.79	+Info

### 12 PK ROUND 125MM

Code	Name	Price	+
9600820R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate Maroon .160in - 4mm 12pk Round 120mm...	US\$ 69.79	+Info
9601040R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate YELLOW .160in - 4mm 12pk Round 120mm...	US\$ 69.79	+Info
9601040R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate YELLOW.160 - 4mm Round 125mm Sheet 1...	US\$ 69.79	+Info



9601020R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate WHITE .160in - 4mm 12pk Round 120mm ...	US\$ 57.23	+Info
9601020R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate WHITE .160 - 4mm Round 125mm Sheet 1...	US\$ 69.79	+Info
9600910R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate RED .160in - 4mm 12pk Round 120mm Sh...	US\$ 69.79	+Info
9600910R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate RED .160 - 4mm Round 125mm Sheet 12p...	US\$ 69.79	+Info
9600870R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate Purple .160in - 4mm 12pk Round 120mm...	US\$ 68.39	+Info
9600870R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate Purple .160 - 4mm Round 125mm Sheet ...	US\$ 69.79	+Info
9600850R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate ORANGE .160in - 4mm 12pk Round 120mm...	US\$ 69.79	+Info
9600850R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate ORANGE.160 - 4mm Round 125mm Sheet 1...	US\$ 69.79	+Info
9596920R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate CLEAR .160 - 4mm Round Clear 125mm ...	US\$ 69.79	+Info
9600820R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate Maroon .160 - 4mm Round 125mm Sheet ...	US\$ 69.79	+Info
9600800R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate GREEN .160in - 4mm 12pk Round 120mm ...	US\$ 69.79	+Info
9600800R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate GREEN .160 - 4mm Round 125mm Sheet 1...	US\$ 69.79	+Info
9600180R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate ASSORTED .160in - 4mm 12pk Round 120...	US\$ 69.79	+Info
9600180R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate ASSORTED 12pk .160 - 4mm Round 125m...	US\$ 69.79	+Info
9598730R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate BLACK .160in - 4mm 12pk Round 120mm ...	US\$ 67.00	+Info
9598730R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate BLACK .160 - 4mm Round 125mm Sheet 1...	US\$ 69.79	+Info
9598060R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate BLUE .160in - 4mm 12pk Round 120mm S...	US\$ 69.79	+Info
9598060R1	PRO-FORM SINGLE-COLOR Mouthguards Laminate BLUE .160 - 4mm Round 125mm Sheet 12...	US\$ 69.79	+Info
9596920R2	PRO-FORM SINGLE-COLOR Mouthguards Laminate CLEAR .160in - 4mm 12pk Round Clear ...	US\$ 67.70	+Info

## PREFORMED HOME MOUTHGUARDS

Best Athletic Mouthguards in 2023: Shock Doctor Gel Max Mouthguard.  
 Best Lightweight Design: SISU Aero NextGen Mouthguard.  
 Best Value: Shock Doctor Pro Mouthguard.  
 Best for Fighters: Venum Challenger Mouthguard.  
 Best for Braces: Shock Doctor Double Braces Mouthguard.



Code	Name	Price	+
------	------	-------	---

425405	Polyshok GAMEguard Small RED SSD pr	US\$ 10.15	+Info
822710	Instant Mouthguards Clear 2pk	US\$ 9.00	+Info
425416	PolyShoK GAMEguard PLUS Small Yello	US\$ 11.62	+Info
425415	Polyshok GAMEguard Large Neon Yello	US\$ 10.15	+Info
425414	PolyShoK GAMEguard PLUS Large BLACK	US\$ 11.62	+Info
425413	PolyShoK GAMEguard PLUS Small RED P	US\$ 11.62	+Info
425412	Polyshok GAMEguard Small BLUE SSD p	US\$ 10.15	+Info
425411	PolyShoK GAMEguard PLUS Large RED P	US\$ 11.62	+Info
425410	PolyShoK GAMEguard PLUS Large GREEN	US\$ 11.62	+Info
425409	Polyshok GAMEguard Small Neon Yello	US\$ 10.15	+Info
425408	PolyShoK GAMEguard PLUS Small BLUE	US\$ 11.62	+Info
425407	PolyShoK GAMEguard PLUS Small Gree	US\$ 11.62	+Info
425406	PolyShoK GAMEguard PLUS Large YELL	US\$ 11.62	+Info
116135	Instant Mouthguards Blue 2pk	US\$ 8.50	+Info
425404	Polyshok GAMEguard Small BLACK SSD	US\$ 10.15	+Info
425403	Polyshok GAMEguard Small Neon Green	US\$ 10.15	+Info
425401	Polyshok GAMEguard Large RED SSD pr	US\$ 10.15	+Info
425400	PolyShoK GAMEguard PLUS Small Black	US\$ 11.62	+Info
424884	Polyshok GAMEguard Large BLUE SSD p	US\$ 10.15	+Info
424835	Polyshok GAMEguard Large BLACK SSD	US\$ 10.15	+Info
34-9601301	Model White PF2 each	US\$ 41.55	+Info
34-9601300	Model Black PF2 each	US\$ 33.24	+Info
116139	Instant Mouthguards green 2pk	US\$ 8.50	+Info
116138	Instant Mouthguards Black 2pk	US\$ 8.33	+Info
116137	Instant Mouthguards Orange 2pk	US\$ 6.80	+Info
116136	Instant Mouthguards Yellow 2pk	US\$ 6.80	+Info

## ONE COLOR RESIN MOUTHGUARD - PRO-FORM

### ONE COLOR RESIN MOUTHGUARD - PRO-FORM

#### Features:

- \* 160 (4mm) thick
- \* 5"x 5" (127 x 127mm) sheets
- \* Also available in round
- \* Available Colors: Clear, Red, Blue, Green, Yellow, Black, Orange, Eva White, Purple

#### Available in Packages of:

- \* 01 Sheet
- \* 12 Sheets
- \* 25 Sheets
- \* 300 Sheets

The 160 (4mm) thick Pro-form Mouthguard Resin Sheets are soft for energy absorption, adding more protection. The resin is perfect for mouthguards and niteguards.



### 300pk

Code	Name	Price	+
7000625	PRO-FORM Mouthguard Resin Sheets, BLUE, 5x5in .160 thick, Soft for energy 300p...	US\$ 298.58	+Info
7000626	PRO-FORM Mouthguard Resin Sheets, GREEN, 5x5in .160 thick - 300pk - Soft for ...	US\$ 298.58	+Info
7000627	PRO-FORM Mouthguard Resin Sheets, RED, 5x5in .160 thick, - 300pk - Soft for e...	US\$ 298.58	+Info
7000628	PRO-FORM Mouthguard Resin Sheets, YELLOW, 5x5in .160 thick - 300 pk - Soft for...	US\$ 298.58	+Info
7000632	PRO-FORM Mouthguard Resin Sheets, WHITE, 5x5in .160 thick - 300pk - Soft for ...	US\$ 298.58	+Info
7000656	PRO-FORM Mouthguard Resin Sheets, CLEAR, 5x5in 0.160in thick - 300 pk - Soft ...	US\$ 292.61	+Info
700063	PRO-FORM Mouthguard Resin Sheets, BLACK, 5x5in .160 thick - 300 pk - Soft for...	US\$ 318.48	+Info
7000629	PRO-FORM Mouthguard Resin Sheets, ORANGE, 5x5in .160 thick, Soft for energy - ...	US\$ 298.58	+Info

### 25pk

Code	Name	Price	+
7000654	PRO-FORM Mouthguard Resin Sheets, CLEAR, 5x5in 0.160in ( 4mm ) thick 25pk ( K#...	US\$ 30.56	+Info
101161	PRO-FORM - RESIN Mouthguards RED.160in ( 4mm ) 25pk 5x5 Sheet ( K#9597910/25 ) (...	US\$ 0.28	+Info
101162	PRO-FORM - RESIN Mouthguards BLUE.160in ( 4mm ) 25pk 5x5 Sheet ( K#9597680/25 ) ...	US\$ 0.25	+Info
101163	PRO-FORM - RESIN Mouthguards GREEN.160in ( 4mm ) 25pk 5x5 Sheet ( K#9597700/25 )...	US\$ 0.28	+Info
101164	PRO-FORM - RESIN Mouthguards YELLOW.160in 25 pk 5x5 Sheet ( 5x5in ) # 9597760/25	US\$ 0.28	+Info
101165	PRO-FORM - RESIN Mouthguards ORANGE.160in ( 4mm ) 25pk 5x5 Sheet ( K#9597720/25 ...	US\$ 0.27	+Info
101166	PRO-FORM - RESIN Mouthguards WHITE.160in ( 4mm ) 25pk 5x5 Sheet ( K#9597750/25 )...	US\$ 0.27	+Info
959766	PRO-FORM - RESIN Mouthguards BLACK.160in ( 4mm ) 25pk 5x5 Sheet ( K#9597660/25 )...	US\$ 0.27	+Info

### 2pk

Code	Name	Price	+
959766-2	PRO-FORM - RESIN Mouthguards BLACK.160in ( 4mm ) 2pk ( K#9597660/2 ) 5x5 square ...	US\$ 3.94	+Info
101166-2	PRO-FORM - RESIN Mouthguards WHITE.160in ( 4mm ) 2pk 5x5 square Sheet ( K#959775...	US\$ 0.05	+Info
101165-2	PRO-FORM - RESIN Mouthguards ORANGE.160in ( 4mm ) 2pk 5x5 square Sheet( K#959772...	US\$ 3.94	+Info
101164-2	PRO-FORM - RESIN Mouthguards YELLOW.160in ( 4mm ) 2pk ( K#9597760/2 ) 5x5 square...	US\$ 4.06	+Info
101163-2	PRO-FORM - RESIN Mouthguards GREEN.160in ( 4mm ) 2pk 5x5 square Sheet( K#9597700/...	US\$ 2.95	+Info
101162-2	PRO-FORM - RESIN Mouthguards BLUE.160in ( 4mm ) 2pk ( K#9597680/2 ) 5x5 square S...	US\$ 3.94	+Info
101161-2	PRO-FORM - RESIN Mouthguards RED.160in ( 4mm ) 2pk ( K#9597910/2 ) 5x5 square Sh...	US\$ 0.05	+Info

## TRI-COLOR MOUTHGUARDS LAMINATES - PRO-FORM

### PRO-FORM TRI-COLOR MOUTHGUARDS

Pro-form Mouthguard Laminate is a strong two layer material made of high-energy absorbing EVA vinyls.

Pro-form MouthGuards are the Official Mouthguard of the U.S Women's Water Polo team for the 2012 Olympics.

ProForm has the advantage over boil-and-bite types of mouthguards because of their superior fit and retention. In addition, ProForm mouthguards are thin, comfortable, custom fitted and will retain their shape long after other mouthguards have worn out.

#### ADVANTAGES Pro-Form Tri-Color Laminate mouthguards

- \* Increased oxygen intake for maximum performance and endurance
- \* Improved speech for clear communication on the field



- \* Natural comfort with a secure fit
- \* Pro-Form Tri-Color Laminate mouthguards Provides maximum impact dispersion and protection against tooth damage
- \* Custom fit for any size mouth
- \* Numerous color options to meet the needs of your patients
- \* fv0.150" thick in tri-color sheets.
- \* We offer Tri-Color Mouthguards in packages of 1 and 12 units

ProForm Tri-Color Mouthguards has the advantage over boil-and-bite types of ProForm Tri-Color Mouthguards because of their superior fit and retention. In addition, ProForm Tri-Color Mouthguards are thin, comfortable, custom fitted and will retain their shape long after other mouthguards have worn out.

Tensile strength, softness, compression percentage and uniform density are all important characteristics of Pro-form mouthguard materials. ProForm Tri-Color Mouthguards maintains these characteristics using a heat/pressure laminating process. The Pro-form mouthguard tensile strength of the laminate is excellent due to the laminate's two layers of materials making up the mouthguard. Density is maintained at the pressure laminating process so shrinkage of the mouthguard is uniformly controlled

### 12PK

Code	Name	Price	+
113583	PRO-FORM TRI-COLOR Mouthguards ASSORTED 5-5 5x5 12pk 0.160in (4mm ) Thick #96001...	US\$ 57.93	+Info
113582	PRO-FORM TRI-COLOR Mouthguards ORANGE/ WHITE/ ORANGE 5x5in 12pk 0.160in - 4mm Th...	US\$ 55.13	+Info
113581	PRO-FORM TRI-COLOR Mouthguards Purple/ BLUE/ YELLOW 5x5in 12pk 0.150in thic #959...	US\$ 69.79	+Info
113580	PRO-FORM TRI-COLOR Mouthguards RED/ YELLOW/ BLUE 5x5in 12pk 0.150in thick	US\$ 69.79	+Info
113579	PRO-FORM TRI-COLOR Mouthguard Laminate - RED/BLACK/RED, 12pk. 5x5in .160in - #...	US\$ 69.79	+Info
113578	PRO-FORM TRI-COLOR Mouthguard Laminate - RED/WHITE/BLACK, 12/Box. 5x5in .160in ...	US\$ 69.79	+Info
113577	PRO-FORM TRI-COLOR Mouthguard Laminate - RED/WHITE/BLUE, 12/Box. 5x5in .160in ...	US\$ 69.79	+Info
034-9599707	PRO-FORM TRI-COLOR Mouthguards BLACK/YELLOW/GREEN 5x5in 12pk 0.160in ( 4mm ) thi...	US\$ 61.46	+Info
034-9599531	PRO-FORM TRI-COLOR Mouthguards BLACK/WHITE/YELLOW 5x5in 12pk .160 - 4mm Thic #95...	US\$ 69.79	+Info
034-9599280	PRO-FORM TRI-COLOR Mouthguard Laminate - Purple/WHITE/Purple, 12/Box. 5x5in - #...	US\$ 67.00	+Info
034-9598920	PRO-FORM TRI-COLOR Mouthguard Laminate - BLUE/YELLOW/BLUE, 12/Box. 5x5in .160in...	US\$ 69.79	+Info

### 1PK

#### Package of 1

Code	Name	Price	+
99060-1	Pro-Form TRI-COLOR Mouthguards GREEN/CLEAR/YELLOW5x5in 1pk 0.160in ( 4mm ) thick...	US\$ 2.44	+Info
98803-1	PRO-FORM TRI-COLOR Mouthguards BLACK/WHITE/GREEN 5x5in 1pk 0.160in ( 4mm ) thick...	US\$ 0.05	+Info
98920-1	PRO-FORM TRI-COLOR Mouthguards BLUE/YELLOW/BLUE 5x5in 1pk 0.160in ( 4mm ) thick ...	US\$ 5.63	+Info
113582-1	PRO-FORM TRI-COLOR Mouthguards ORANGE/WHITE/ORANGE 5x5in 1pk 0.160in ( 4mm ) Thi...	US\$ 5.18	+Info
113581-1	PRO-FORM TRI-COLOR Mouthguards Purple/BLUE/YELLOW 5x5in 1pk 0.160in thick #95991...	US\$ 5.67	+Info
113580-1	PRO-FORM TRI-COLOR Mouthguards RED/YELLOW/BLUE 5x5in 1pk 0.160in thick ##9599205-1	US\$ 5.18	+Info
113579-1	PRO-FORM TRI-COLOR Mouthguards RED/BLACK/RED 5x5in 1pk 0.160in thick ##9598940-1	US\$ 5.18	+Info
113578-1	PRO-FORM TRI-COLOR Mouthguards RED/WHITE/BLACK 5x5in 1pk 0.160in thick ##9598800-1	US\$ 5.63	+Info

113577-1	PRO-FORM TRI-COLOR Mouthguards RED/WHITE/BLUE 5x5in 1pk 0.160in thick ##9598840-1	US\$ 0.70	<a href="#">+Info</a>
----------	-----------------------------------------------------------------------------------	-----------	-----------------------

## TWO COLOR MOUTHGUARDS + LAMINATES - PRO-FORM

### PRO-FORM DUAL COLOR MOUTHGUARDS

#### Feature

Proform Mouthguard Laminate is a strong two layer material made of high-energy absorbing EVA vinyls. Dual-color Mouth guard Laminate in assorted two-color combos. The Pro-form Mouthguards Soft EVA Resin is a rubber-like clear material easily formed and trimmed for bleaching trays, soft bruxing appliances, indirect bonding and TMJ splints. 12 sheets per pack.

#### ProForm Advantages

- \* Increased oxygen intake for maximum performance and endurance
- \* Improved speech for clear communication on the field
- \* Natural comfort with a secure fit
- \* Provides maximum impact dispersion and protection against tooth damage
- \* Custom fit for any size mouth
- \* Numerous color options to meet the needs of your patients
- \* Dual-color Mouthguard Laminate in assorted two-color combos.
- \* Energy-absorbing mouthguard sheets
- \* 0.150" thick in dual-color sheets.

Laboratory processed professional grade ProForm mouthguards offer the highest level of protection, retention, comfort and fit, without interfering with speech or breathing. ProForm has the advantage over boil-and-bite types of mouthguards because of their superior fit and retention.

In addition, ProForm mouthguards are thin, comfortable, custom fitted and will retain their shape long after other mouthguards have worn out. Pro-form mouthguard Tensile strength, softness, compression percentage and uniform density are all important characteristics of mouthguard materials. ProForm maintains these characteristics using a heat/pressure laminating process. Tensile strength of the laminate is excellent due to the laminate's two layers of materials making up the mouthguard. Density is maintained at the pressure laminating process so shrinkage of the mouthguard is uniformly controlled.

**We offer Dual-Color Mouthguards in packages of 1 and 12 units**



### MOUTHGUARDS 12 PK

Code	Name	Price	+
------	------	-------	---

113506	PRO-FORM DUAL-COLOR Mouthguard Laminate - GREEN/YELLOW, 12/Box. 5x5in .160in ...	US\$ 69.79	+Info
113507	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLACK/RED, 12/Box. 5x5in .160in (4m...	US\$ 69.79	+Info
113508	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLACK/YELLOW, 12/Box. 5x5in .160in ...	US\$ 60.72	+Info
113509	PRO-FORM DUAL-COLOR Mouthguard Laminate - RED/YELLOW, 12pk. 5x5in .160in - #9...	US\$ 69.79	+Info
113510	PRO-FORM DUAL-COLOR Mouthguard Laminate - ASSORTED Colors, 12/Box. 5x5in .160i...	US\$ 66.30	+Info
113685	PRO-FORM DUAL-COLOR Mouthguard Laminate - GREEN/BLUE, 12/Box. 5x5in .160in - ...	US\$ 69.79	+Info
113686	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLACK/WHITE, 12/Box. 5x5in .160in -...	US\$ 57.23	+Info
113687	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLUE/RED, 12/Box. 5x5in .160in (4mm...	US\$ 69.79	+Info
113688	PRO-FORM DUAL-COLOR Mouthguard Laminate - Purple/YELLOW, 12/Box. 5x5in .160in ...	US\$ 69.79	+Info
113689	PRO-FORM DUAL-COLOR Mouthguard Laminate - Maroon/YELLOW, 12pk. 5x5in .160in -...	US\$ 69.79	+Info
113690	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLUE/YELLOW, 12/Box. 5x5in .160in -...	US\$ 69.79	+Info
113768	PRO-FORM DUAL-COLOR Mouthguard Laminate - RED/WHITE, 12pk. 5x5in .160in (4mm)...	US\$ 69.79	+Info
113769	PRO-FORM DUAL-COLOR Mouthguard Laminate - ORANGE/WHITE, 12/Box. 5x5in .160in ...	US\$ 69.79	+Info
113770	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLUE/WHITE, 12/Box. 5x5in .160in - ...	US\$ 69.79	+Info
113771	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLUE/ORANGE, 12pk. 5x5in .160in - #...	US\$ 69.79	+Info
113772	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLACK/ORANGE, 12/Box. 5x5in .160in ...	US\$ 69.79	+Info
113773	PRO-FORM DUAL-COLOR Mouthguard Laminate - WHITE/Purple, 12/Box. 5x5in .160in ...	US\$ 69.79	+Info
113774	PRO-FORM DUAL-COLOR Mouthguard Laminate - GREEN/WHITE, 12/Box. 5x5in .160in -...	US\$ 69.79	+Info
113775	PRO-FORM DUAL-COLOR Mouthguard Laminate - Maroon/WHITE, 12pk. 5x5in .160in - ...	US\$ 69.79	+Info
113776	Pro-Form Dual-Color Mouthguard Laminate - RED/CLEAR, 12/Box. 5x5in .160in (4m...	US\$ 57.48	+Info
113777	PRO-FORM DUAL-COLOR Mouthguard Laminate - BLACK/BLUE, 12/Box. 5x5in .160in - ...	US\$ 69.79	+Info
113778	PRO-FORM DUAL-COLOR Mouthguards Laminate ORANGE/CLEAR 5x5 12pk 0.160in thick #9...	US\$ 69.79	+Info

### MOUTHGUARDS 1 PK

Code	Name	Price	+
113506-1	PRO-FORM DUAL-COLOR Mouthguards Laminate GREEN/YELLOW 5x5 1pk 0.160in thick #R...	US\$ 5.24	+Info
113507-1	PRO-FORM DUAL-COLOR Mouthguards Laminate BLACK/RED 5x5 1pk 0.160in thick	US\$ 5.24	+Info
113508-1	<del>#RR-1</del> PRO-FORM DUAL-COLOR Mouthguards Laminate BLACK/ YELLOW 5x5 1pk 0.160in thick	US\$ 0.06	+Info
113509-1	<del>#</del> PRO-FORM DUAL-COLOR Mouthguards Laminate RED/YELLOW 5x5 1pk 0.160in thick #RR-...	US\$ 5.12	+Info
113685-1	PRO-FORM DUAL-COLOR Mouthguards Laminate GREEN/BLUE 5x5 1pk 0.160in thick	US\$ 5.24	+Info
113686-1	<del>#RR</del> PRO-FORM DUAL-COLOR Mouthguards Laminate BLACK/WHITE 5x5 1pk 0.160in thick	US\$ 5.63	+Info
113687-1	<del>#RR</del> PRO-FORM DUAL-COLOR Mouthguards Laminate BLUE/RED 5x5 1pk 0.160in thick #RR-11...	US\$ 5.63	+Info
113688-1	PRO-FORM DUAL-COLOR Mouthguards Laminate Purple/YELLOW 5x5 1pk 0.160in thick #...	US\$ 5.12	+Info
113689-1	PRO-FORM DUAL-COLOR Mouthguards Laminate Maroon/YELLOW 5x5 1pk 0.160in thick #...	US\$ 4.62	+Info
113690-1	PRO-FORM DUAL-COLOR Mouthguards Laminate BLUE/YELLOW 5x5 1pk 0.160in thick #RR...	US\$ 5.24	+Info
113768-1	PRO-FORM DUAL-COLOR Mouthguards Laminate RED/WHITE 5x5 1pk 0.160in thick	US\$ 5.24	+Info
113769-1	<del>#RR-1</del> PRO-FORM DUAL-COLOR Mouthguards Laminate ORANGE/WHITE 5x5 1pk 0.160in thick #R...	US\$ 5.24	+Info
113770-1	PRO-FORM DUAL-COLOR Mouthguards Laminate BLUE/WHITE 5x5 1pk 0.160in thick #RR-...	US\$ 0.06	+Info
113771-1	PRO-FORM DUAL-COLOR Mouthguards Laminate BLUE/ORANGE 5x5 1pk 0.160in thick #RR...	US\$ 5.12	+Info
113772-1	PRO-FORM DUAL-COLOR Mouthguards Laminate BLACK/ORANGE 5x5 1pk 0.160in thick #R...	US\$ 5.46	+Info
113773-1	PRO-FORM DUAL-COLOR Mouthguards Laminate WHITE/Purple 5x5 1pk 0.160in thick ...	US\$ 5.24	+Info
113774-1	PRO-FORM DUAL-COLOR Mouthguards Laminate GREEN/WHITE 5x5 1pk 0.160in thick #RR...	US\$ 5.24	+Info

113775-1	PRO-FORM DUAL-COLOR Mouthguards Laminate Maroon/WHITE 5x5 1pk 0.160in thick #R...	US\$ 5.57	<a href="#">+Info</a>
113776-1	Pro-Form DUAL-COLOR Mouthguards Laminate RED/CLEAR 5x5 1pk 0.160in thick #RR-1...	US\$ 7.28	<a href="#">+Info</a>
113777-1	PRO-FORM DUAL-COLOR Mouthguards Laminate BLACK/BLUE 5x5 1k 0.160in thick #	US\$ 5.63	<a href="#">+Info</a>
113778-1	PRO-FORM DUAL-COLOR Mouthguards Laminate ORANGE/CLEAR 5x5 1pk 0.160in thick #R...	US\$ 5.24	<a href="#">+Info</a>

## YOUTH RESIN MOUTHGUARD - PRO-FORM

### YOUTH RESIN MOUTHGUARD - PRO-FORM

The .120 thick Youth Mouthguard Laminates are great for little kids who are playing sports. Made of soft EVA, for energy absorption and is thin enough to fit the smaller mouth of a child. It remains comfortable for the child who is just getting into athletic competition.



#### Features:

- \* .120 thick
- \* 5" x 5" square sheets
- \* Also available in round
- \* Package of 25
- \* Available in 5 other colors
- \* Available in 6 Colors

### 25pk

Code	Name	Price	+
101169	PRO-FORM - YOUTH RESIN Mouthguards CLEAR .120in ( 3mm ) 5x5 square Sheet 25pk ( ...	US\$ 22.99	<a href="#">+Info</a>
101170	PRO-FORM - YOUTH RESIN Mouthguards RED .120in ( 3mm ) 5x5 square Sheet 25pk ( K#...	US\$ 31.80	<a href="#">+Info</a>
101171	PRO-FORM - YOUTH RESIN Mouthguards BLUE .120in ( 3mm ) 5x5 square Sheet 25pk ( K...	US\$ 21.64	<a href="#">+Info</a>
101172	PRO-FORM - YOUTH RESIN Mouthguards GREEN .120in ( 3mm ) 5x5 square Sheet 25pk (...	US\$ 18.83	<a href="#">+Info</a>
101173	PRO-FORM - YOUTH RESIN Mouthguards YELLOW .120in ( 3mm ) 5x5 square Sheet 25pk ...	US\$ 27.27	<a href="#">+Info</a>
420851	PRO-FORM YOUTH RESIN Mouthguards WHITE .120in - 3mm 5x5 square Sheet 25p #9611...	US\$ 22.58	<a href="#">+Info</a>
101174	PRO-FORM - YOUTH RESIN Mouthguards BLACK .120in ( 3mm ) 5x5 square Sheet 25pk ( ...	US\$ 22.99	<a href="#">+Info</a>

### 300pk

Code	Name	Price	+
157088	PRO-FORM - YOUTH RESIN Mouthguards CLEAR .120in ( 3mm ) 5x5 square Sheet 300pk (...	US\$ 236.80	<a href="#">+Info</a>
700527	PRO-FORM Youth Mouthguard Material, BLUE, 5x5in .120 thick, Soft EVA - #9611282	US\$ 255.01	<a href="#">+Info</a>
700528	PRO-FORM Youth Mouthguard Material, BLACK, 5x5in .120 thick, Soft EVA - #9611285	US\$ 255.01	<a href="#">+Info</a>
700531	PRO-FORM Youth Mouthguard Material, RED, 5x5in .120 thick, Soft EVA for energy...	US\$ 255.01	<a href="#">+Info</a>
700538	PRO-FORM Youth Mouthguard Material, YELLOW, 5x5in .120 thick, Soft EVA - #9611284	US\$ 255.01	<a href="#">+Info</a>
700539	PRO-FORM Youth Mouthguard Material, GREEN, 5x5in .120 thick, Soft EVA - #9611283	US\$ 255.01	<a href="#">+Info</a>

**RETAINER MATERIAL - KEYSTONE**  
**RETAINER MATERIAL - KEYSTONE**

His stiff but clear material is ideal for orthodontic retainer.  
The long-lasting material enables acrylics to stick, but also provides the long-lasting form that retainers need.

**Specifications:**

- \* Clear material
- \* 5" x 5" ( 127x127mm ) sheets
- \* Also available in round
- \* Available in 25, 50 and 625 pieces per pack
- \* .030 (75mm)



**Us Dental Depot Is An authorized Keystone dealer In Us**

**25pk**

Code	Name	Price	+
101181	PRO-FORM Retainer Material, 5x5in .030 Super CLEAR, pkg of 25. Stiff but CLEAR...	<b>US\$ 58.43</b>	<a href="#">+Info</a>

**50pk**

Code	Name	Price	+
101182	PRO-FORM Retainer Material, 5x5in .030 Super CLEAR, pkg of 50. Stiff but CLEAR...	<b>US\$ 88.07</b>	<a href="#">+Info</a>

**625pk**

Code	Name	Price	+
101183	PRO-FORM - Retainer Material Super CLEAR - .030in - 0.7mm - 625 pk #9614875 -	<b>US\$ 824.59</b>	<a href="#">+Info</a>

**Vacuum Forming Thermoplastic Materials**  
**Soft Eva**



**Soft Eva - BUFFALO**

Code	Name	Price	+
424441	.040" Clear soft EVA 25pk	<b>US\$ 20.51</b>	<a href="#">+Info</a>
424476	.040 Soft EVA Clear 100pk	<b>US\$ 59.72</b>	<a href="#">+Info</a>
424509	.040 Soft EVA Clear 300pk	<b>US\$ 146.65</b>	<a href="#">+Info</a>
424520	.040 Soft EVA Clear 12pk	<b>US\$ 11.23</b>	<a href="#">+Info</a>

**FREE SHIPPING FOR ALL DOMESTIC ACCOUNTS!**  **US DENTAL DEPOT** +1 954 874 6325   
www.usdentaldepot.com 



**SOFT EVA - KEYSTONE**  
**SOFT EVA - KEYSTONE**

Soft EVA is great for the fabrication of customer fluoride trays, bleaching trays and model duplication. The material is soft, clear and easily forms. Also the material can be trimmed effortlessly. The material is ideal for soft bruxing patients.

**Features:**

- \* Soft, Clear Material
- \* Easily Formed and Trimmed
- \* Custom Fluoride Trays, Model Duplication
- \* 5 x 5 inches
- \* Bleaching Trays, Soft Bruxing
- \* Packages of 10, 25, or 300 Sheets
- \* Vacuum Forming Material



**10pk**

Code	Name	Price	+
103921	KEYSTONE Soft Eva - .040in ( 1mm ) - 10pk - 5in x5in sheets - CLEAR ( K #9597040 ...	<b>US\$ 6.86</b>	<a href="#">+Info</a>
103922	KEYSTONE Soft Eva - .060in ( 1.5mm ) - 10pk - 5in x5in sheets - CLEAR ( K#959713...	<b>US\$ 7.57</b>	<a href="#">+Info</a>
103923	KEYSTONE Soft Eva - .080in ( 2mm ) - 10pk - 5x5in sheets - CLEAR ( K#9597240 ) (...	<b>US\$ 11.06</b>	<a href="#">+Info</a>

**25pk**

Code	Name	Price	+
100613	KEYSTONE Soft Eva - .040in ( 1mm ) - 25pk - 5x5 sheets - CLEAR ( K#9596980 ) ( 5...	<b>US\$ 13.37</b>	<a href="#">+Info</a>
101158	KEYSTONE Soft Eva - .060in ( 1.5mm ) - 25pk - 5x5 sheets - CLEAR ( K#9597090 ) (...	<b>US\$ 18.32</b>	<a href="#">+Info</a>
9597150	KEYSTONE Soft Eva - .080in ( 2mm ) - 25pk - 5x5 sheets - CLEAR ( K#9597150 ) ( 5...	<b>US\$ 20.15</b>	<a href="#">+Info</a>
9611240	KEYSTONE Soft Eva - .120in ( 3mm ) - 25pk - 5x5 sheets - CLEAR ( K#9611240 ) ( ...	<b>US\$ 21.20</b>	<a href="#">+Info</a>

**300pk**

Code	Name	Price	+
034-1570182	PRO-FORM .040 Soft EVA Tray Material 5x5in 300pk. Soft, CLEAR, easily formed -...	<b>US\$ 152.66</b>	<a href="#">+Info</a>
034-1570184	PRO-FORM .060 Soft EVA Tray Material 5x5in 300pk. Soft, CLEAR, easily formed -...	<b>US\$ 219.59</b>	<a href="#">+Info</a>
034-1570186	PRO-FORM .080 Soft EVA Tray Material 5x5in 300pk. Soft, CLEAR, easily formed -...	<b>US\$ 241.57</b>	<a href="#">+Info</a>
101159	PRO-FORM Soft EVA CLEAR 5x5in .120 thick 300pk Soft EVA for energy - #9611280	<b>US\$ 255.01</b>	<a href="#">+Info</a>

**ROUND**

Code	Name	Price	+
417353	PRO-FORM - Soft Eva - .060in (1.5mm) 125mm - 300pk - Round sheets - CLEAR #95971...	<b>US\$ 219.59</b>	<a href="#">+Info</a>
469770	KEYSTONE .080 (2mm) Soft EVA Round Clear 120mm - 25pk (#9597150R2)	<b>US\$ 28.05</b>	<a href="#">+Info</a>
469771	KEYSTONE .060 (1.5mm) Soft EVA Round Clear 120mm - 300/Pkg (#9597130R2)	<b>US\$ 219.59</b>	<a href="#">+Info</a>
469772	KEYSTONE .040 (1mm) Soft EVA Round Clear 120mm - 300/Pkg (#9597040R2)	<b>US\$ 153.98</b>	<a href="#">+Info</a>

469773	KEYSTONE -.120 (3mm) Soft EVA Round Clear 120mm - 25pk (#9611240R2)	US\$ 31.80	<a href="#">+Info</a>
469774	KEYSTONE .060 (1.5mm) Soft EVA Round Clear 120mm - 25pk (#9597090R2)	US\$ 26.32	<a href="#">+Info</a>
469775	KEYSTONE .040 (1mm) Soft EVA Round Clear 120mm - 25pk (#9596980R2)	US\$ 26.32	<a href="#">+Info</a>
469881	PRO-FORM - Soft EVA Round Clear .040 (1mm) 125mm - 300/Pkg (#9597040R1)	US\$ 153.98	<a href="#">+Info</a>

### Vacuum Forming Thermoplastic Materials

## Surgical Splint Material



### Splint Material - BUFFALO

SPLINT MATERIAL - BUFFALO High quality 5" x 5" sheet thermoplastic materials for creating a wide range of vacuum formed appliances and trays.

Pure grade resin sheets are latex-free and uniform in size, cut and thickness for improved performance. Resist bubbling, crazing and fracturing to promote superior adaptations every time. Clear soft Ethyl Vinyl Acetate (EVA) sheets designed to bond together for build-up and lamination of custom mouthguards



Features Rigid Clear tray and Splint Material

Buffalo Clear Material offers the best model adaptation and excellent rigidity. The result will be a clear, highly polished surface with no color dyes or bluing agents. The Extra Wide size stays secure in all frames.

Code	Name	Price	+
424500	.020 Tempo Spl Mat Clear 25pk 5x5	US\$ 20.99	<a href="#">+Info</a>
424979	.060" Rigid clear tr	US\$ 445.96	<a href="#">+Info</a>
424538	.080" Rigid clear tr	US\$ 23.18	<a href="#">+Info</a>
425003	.100" Rigid clear tr	US\$ 33.86	<a href="#">+Info</a>
424566	.040" Rigid clear tr	US\$ 14.70	<a href="#">+Info</a>
424744	.100" Rigid clear tr	US\$ 575.52	<a href="#">+Info</a>
424750	.040" Rigid clear tr	US\$ 283.80	<a href="#">+Info</a>
424858	.060" Rigid clear tr	US\$ 23.45	<a href="#">+Info</a>
424925	.080" Rigid clear tr	US\$ 510.77	<a href="#">+Info</a>
424930	.020" Temporary Spli	US\$ 188.69	<a href="#">+Info</a>

### SPLINT MATERIAL - KEYSTONE

#### SPLINT MATERIAL - KEYSTONE

This stiff clear material is ideal for hard splints, thin splints, surgical trays, orthodontic retainers and bleaching trays. The material sticks to acrylics for various dental appliances. The versatile splint material will surely be used in any lab or office.

#### Features:

- \* Stiff, Clear Material
- \* Acrylic Sticks To Material
- \* 5 x 5inch Sheets
- \* Packages of 25, 50 or 625
- \* For Hard splints, Surgical Trays , Orthodontic Retainers and More.





\* Vacuum Forming Material

\* Available in Five Thickness \* .020in \* .040in \* .060in \* .080in \* .10

### 25pk (Square)

Code	Name	Price	+
103924	PRO-FORM .020in Splint Material (.50mm) 5x5 (127x127mm) - 25pk #9615000	US\$ 10.97	+Info
122334	KEYSTONE Splint Material - .030in - ( 0.8mm ) (.75mm) Splint 5x5 (127x127mm) - ...	US\$ 23.25	+Info
103928	Splint Material - .040in (1mm) 5x5 (127x127mm) - 25pk #9615060	US\$ 16.28	+Info
103925	Splint Material - .060in ( 1.5mm ) Splint 5x5 (127x127mm) - 25pk #9614960	US\$ 20.78	+Info
103926	Splint Material - .080in ( 2mm ) Splint 5x5 (127x127mm) - 25pk #9614980	US\$ 24.53	+Info
103927	Splint Material - .100in Splint 5x5 (127x127mm) - 25pk #9615040	US\$ 23.44	+Info

### 50pk (Square)

Code	Name	Price	+
101184	PRO-FORM .020in Splint Material 5x5in 50pk Stiff CLEAR material is ideal - #...	US\$ 17.90	+Info
700105	PRO-FORM .030in Retainer Material. Stiff, CLEAR. 5x5in Sheets - 50pk - Acryli...	US\$ 27.22	+Info
101156	PRO-FORM .040in Splint Material, 5x5in 50pk Stiff CLEAR material is ideal - ...	US\$ 23.25	+Info
961478	PRO-FORM .060in Splint Material, 5x5in 50pk Stiff CLEAR material is ideal - ...	US\$ 37.65	+Info
961479	PRO-FORM .080in Splint Material, 5x5in 50pk Stiff CLEAR material is ideal - ...	US\$ 47.25	+Info
961482	Splint Material - .100in ( 2.5mm ) - 50pk - 5x5 sheets - ( K#9614820 ) ( 5X5in )	US\$ 36.55	+Info

### 625pk (Square)

Code	Name	Price	+
101157	PRO-FORM .020in Splint Material, 5x5 625pk Stiff CLEAR material is ideal - #...	US\$ 241.50	+Info
700107	PRO-FORM .030in (0.75mm) Splint Material, 5x5in (127 x 127mm) 625pk Stiff - ...	US\$ 273.63	+Info
961488	Splint Material - .040in ( 1mm ) - 5in x5in sheets - 600pk - ( K#9614880 ) ( 5X5...	US\$ 271.80	+Info
700037	KEYSTONE Splint Material - .060in ( 1.5mm ) - 600pk - 5in x5in - sheets - 9614840	US\$ 409.63	+Info
700041	Splint Material - .080in ( 2mm ) - 5in x5in sheets - 600pk - ( K#9614850 ) ( #15...	US\$ 370.25	+Info
700047	PRO-FORM .100in Splint Material, 5x5in 625pk Stiff CLEAR material is ideal - ...	US\$ 486.01	+Info

### Round 25pk

#### Round 120mm

Code	Name	Price	+
416706	PRO-FORM .020in (.50mm) Splint Material, Round 120mm, 25pk Stiff CLEAR materi...	US\$ 11.25	+Info
416708	PRO-FORM .030in (.75mm) Splint Material, Round 120mm, 25pk Stiff CLEAR materi...	US\$ 23.25	+Info

#### Round 125mm

Code	Name	Price	+
416705	PRO-FORM .020in (.50mm) Splint Material, Round 125mm, 25pk Stiff CLEAR materi...	US\$ 11.25	+Info
416707	KEYSTONE .030in (.75mm) Splint Material, Round 125mm, 25pk Stiff CLEAR - #961...	US\$ 23.25	+Info

### Surgical Splint Material - ASTRON

Surgical Splint Material - ASTRON

Code	Name	Price	+
115322	CLEARSplint Flexible Kit	US\$ 463.00	<a href="#">+Info</a>
115323	CLEARSplint Regular Kit	US\$ 463.00	<a href="#">+Info</a>

### Surgical Splint Material - DENTSPLY

Surgical Splint Material - DENTSPLY

Code	Name	Price	+
421103	ESSIX ACE 125 MM CIRCLE.040 100PK	US\$ 279.98	<a href="#">+Info</a>
421104	ESSIX ACE 125MM CIRCLE .035 100PK	US\$ 279.98	<a href="#">+Info</a>
421105	ESSIX ACE 125MM CIRCLE .030 100PK	US\$ 279.98	<a href="#">+Info</a>

### Surgical Splint Material - ERKODENT

Surgical Splint Material - ERKODENT

Code	Name	Price	+
420997	ERKODUR Splint Mat 0.2in 50pk	US\$ 470.39	<a href="#">+Info</a>

### Surgical Splint Material - JSP

Surgical Splint Material - JSP

Code	Name	Price	+
461867	Jsp - Splint Material - .080in - 25	US\$ 27.30	<a href="#">+Info</a>

### Surgical Splint Material - KEYSTONE

Surgical Splint Material - KEYSTONE

Code	Name	Price	+
470262	PRO-Splint Material .040in 25pk	US\$ 14.62	<a href="#">+Info</a>
470263	PRO-Splint Material .060in 25pk	US\$ 18.83	<a href="#">+Info</a>
470264	PRO-Splint Material .080in 25pk	US\$ 24.77	<a href="#">+Info</a>
470265	PRO-Splint Material .040in 625pk	US\$ 307.23	<a href="#">+Info</a>
470266	PRO-Splint Material .060in 625pk	US\$ 373.95	<a href="#">+Info</a>
470267	PRO-Splint Material .080in 625pk	US\$ 504.99	<a href="#">+Info</a>

### Surgical Splint Material - PRO-FORM

Surgical Splint Material - PRO-FORM

Code	Name	Price	+
122287-50	Splint Material -.080in 50pk	US\$ 32.55	<a href="#">+Info</a>

### Surgical Splint Material - PULPDENT

Surgical Splint Material - PULPDENT

Code	Name	Price	+
039-SPL-F	SplinMat S.Steel Grid 4mm	US\$ 101.17	<a href="#">+Info</a>
039-SPL-P	PerfoMat S.Steel Strip 3.5mm	US\$ 101.17	<a href="#">+Info</a>

### Vacuum Forming Thermoplastic Materials

## Thermoplastic Accesories



### MOUTHGUARD STRAPS - KEYSTONE

#### MOUTHGUARD STRAPS - KEYSTONE

Keystone has developed an easy-to-attach helmet strap used in the production of custom-made mouthguards. These straps have a safety release feature that allows the mouthguard to break away under 12 pounds of pressure.

- \* Available in black and clear.
- \* Sold in packages of 12 and 25.
- \* Mouthguard Straps compatible with all Mouthguards



#### 12PK

Code	Name	Price	+
34-9597400	PRO-FORM CLEAR Mouthguard Straps to be used with PRO-FORM mouthguards, pack - ...	US\$ 13.65	<a href="#">+Info</a>
34-9597340	PRO-FORM BLACK Mouthguard Straps to be used with PRO-FORM mouthguards, pack - ...	US\$ 13.65	<a href="#">+Info</a>

#### 25PK

Code	Name	Price	+
34-9597420	PRO-FORM - Mouthgards Straps CLEAR 25pk - Holder to retains mouthgard attacvhed ...	US\$ 21.83	<a href="#">+Info</a>
34-9597360	PRO-FORM - Mouthgards Straps BLACK 25pk - Holder to retains mouthgard attacvhed ...	US\$ 21.83	<a href="#">+Info</a>

### Thermoplastic Accesories - BUFFALO

Thermoplastic Accesories - BUFFALO

Code	Name	Price	+
424502	Hot Plate, Regular,	US\$ 14.83	<a href="#">+Info</a>

### Thermoplastic Accesories - DIGITECH

Thermoplastic Accesories - DIGITECH

Code	Name	Price	+
------	------	-------	---

123166	SHARK - Wheel Saw 16mm 5pk	US\$ 4.00	<a href="#">+Info</a>
123167	SHARK - Wheel Saw 19mm 5pk	US\$ 4.00	<a href="#">+Info</a>
123168	SHARK - Wheel Saw 22mm 5pk	US\$ 4.00	<a href="#">+Info</a>

### PUTTY-BLOCK - ADS

#### PUTTY-BLOCK - ADS

All new Putty-Block is ideal for boxing and beading impressions prior to pouring of any die material. When pressed into place, Putty-Block conforms tightly to the impression, closes open ends, and dams off any flow of liquids. Facial and posterior walls of impressions can be built up to permit easy, spill-free pouring. Releases cleanly from hardened die materials and impressions and is reusable. When molded into place, Putty-Block retains its shape, will not slump or run. Ready to use. Non-hazardous.



Code	Name	Price	+
118312	ADS Putty - Block 120 g Jar - #P800-1 ( Model Blocker )	US\$ 24.75	<a href="#">+Info</a>
118313	ADS Putty - Block 340 g Jar - #P800-2 ADS Putty - Block 120 g Jar - #P800-1 ( Mo...	US\$ 51.08	<a href="#">+Info</a>

### Thermoplastic Accesories - KEYSTONE

#### Thermoplastic Accesories - KEYSTONE

Code	Name	Price	+
034-1410010	838 Tray Handles- Alum. 100bx	US\$ 13.28	<a href="#">+Info</a>
123014	KEYSTONE - Wheel Saw	US\$ 20.85	<a href="#">+Info</a>
420354	PROCURE Block Out 1.2cc syr + tips	US\$ 44.78	<a href="#">+Info</a>

### Thermoplastic Accesories - PRO-FORM

#### Thermoplastic Accesories - PRO-FORM

Code	Name	Price	+
034-1009100	Centri-Fuse Laminate Kit	US\$ 151.30	<a href="#">+Info</a>
034-1009101	Centri-Fuse Refill Sheets 20pk	US\$ 35.05	<a href="#">+Info</a>
034-9575615	Heavy Duty Strap Box Blck 12pk	US\$ 10.50	<a href="#">+Info</a>
034-9575635	Heavy Duty Strap Box Gray 12pk	US\$ 10.50	<a href="#">+Info</a>

### Thermoplastic Accesories - SHARK

#### Thermoplastic Accesories - SHARK

Code	Name	Price	+
420850	SHARK - Galaxy cutting wheel 10mm	US\$ 36.46	<a href="#">+Info</a>

## Catalog Index / Indice

Base Strenghteners .....	1
Bleaching / Whitening Dual laminates .....	1
Coping Material .....	1
Crown & Bridge Material .....	2
Cutters and Knifes .....	3
Dual laminates 4 Bruxism .....	3
Mouthguards .....	5
Retainer Material .....	15
Soft Eva .....	16
Surgical Splint Material .....	18
Thermoplastic Accesories .....	21