

## Vacuum formers

### Vacuum forming Units



### ECCO VACUUM FORMING - ECCO

#### VACUUM FORMING - ECCO

- \* Vacuum Forming machine 220volts or 110volts.
- \* For all thermoplastic laminates.
- \* Square and round sheets suitable.
- \* Uniform heating distribution - 2 Year Warranty thermoformer vacuum forming.



**Shipping information:** 25?21?30cm - 6Kg - 15lb

Code	Name	Price	+
102348	BESQUAL Meta Dental Vacuum Forming Machine 110V. Heavy duty spring and latch mec...	US\$ 184.54	<a href="#">+Info</a>
102510	BESQUAL VAC - ESQUAL Meta Dental Vacuum Forming Machine 220V. Heavy duty spring ...	US\$ 203.68	<a href="#">+Info</a>

### PRO-FORM VACUUM FORMER - KEYSTONE

#### PRO-FORM VACUUM FORMER - KEYSTONE

The Pro-Form is an economical answer that can still excellently fabricate most thermoplastic appliances. The all-aluminum construction eliminates rusting and the heating element is symmetrically engineered for uniform heat distribution. The hexagonal post keeps the frame from rotating or sliding.

#### Specifications:

- \* Uses 5"x5" resin sheets
- \* Available in 110 Volt and 220 Volt
- \* Available in Single and Dual Chambered
- \* Symmetrical heating system
- \* Lightweight, rust-free aluminum frame
- \* Made in the U.S.A.
- \* 5-year warranty on parts and labor



#### PRO-FORM I

Code	Name	Price	+
101208	PRO-FORM Vacuum Former, 110 Volt, Single ChamberRED, Aluminum Construction, 5 - #...	US\$ 350.66	<a href="#">+Info</a>
110746	PRO-FORM Vacuum Former, 220 Volt for international plugs/outlets. #9594480	US\$ 350.66	<a href="#">+Info</a>

#### PRO-FORM II

Code	Name	Price	+
101209	KEYSTONE PRO - FORM Dual ChamberRED Vacuum Former - 110vol - made in USA - 5 yea...	US\$ 379.07	<a href="#">+Info</a>

110747	KEYSTONE PRO - FORM Dual ChamberRED Vacuum Former - 220V - made in USA - 5 year p...	US\$ 416.11	<a href="#">+Info</a>
--------	--	-------------	-----------------------

## PROLINE AIR VACUUM 2000 VENTURI - KEYSTONE

### PROLINE AIR VACUUM 2000 VENTURI - KEYSTONE

Don't let the size of this machine misguide you. The powerful pressure unit provides superior adaptation for dental plastics using laboratory or operatory internal air compressors. Because the unit has no motor, it is compact and very quiet. Also, the vacuum pressure from an office compressor is usually greater than the pressure generated from smaller vacuum units. The greater the pressure, the better definition and adherence the dental appliances will have.

#### Specifications:

- \* Can be used to layer appliances
- \* Very economical when compared to its European counterparts
- \* Uses 5"x5" resin sheets or round resin sheets
- \* 110 Volt
- \* Symmetrical heating system
- \* Lightweight, rust-free aluminum frame
- \* Made in the U.S.A.
- \* 5-year warranty on parts and labor

**Us Dental Depot Is An Authorized Keystone Dealer In US**



Code	Name	Price	+
034-7000470	KEYSTONE Proline Air Vac 2000 - 110vol - vacuum former - made in usa - 5 year w...	US\$ 545.94	<a href="#">+Info</a>
034-7000471	KEYSTONE PRoline Air Vacuum - 220volts vacuum former - made in usa - 5 year warr...	US\$ 537.54	<a href="#">+Info</a>

## THE MACHINE - KEYSTONE

### THE MACHINE - KEYSTONE

The Machine has been a staple in Keystone Industries' long-lasting line of highly recommended vacuum former machines.

The price, features and reliability make it the first choice in dental vacuum formers.

The machine has more power than any other dental vacuum former, which means the production of more accurate and definite dental thermoplastic appliance.

The aluminum construction is long-lasting, eliminating the possibility of rust. The heating element is symmetrically engineered for even heat distribution throughout the thermoplastic material.

This machine is used to process bleaching trays, sports mouthguards, implant stints, temporary splints, custom impression trays, denture bases, copings, temporary crowns and bridges and orthodontic appliances.



### THE MACHINE III

Code	Name	Price	+
034-1560010	KEYSTONE Machine III Vacuum Formers 110volts, - #7000330	US\$ 350.85	<a href="#">+Info</a>
034-7000331	KEYSTONE The Machine - III - 220 Volt vacuum former machine - vacuum forming -	US\$ 350.85	<a href="#">+Info</a>

### THE MACHINE IV

Code	Name	Price	+
34-7000338	KEYSTONE Machine IV Vacuum Former. Base of machine is made from super strong - #...	US\$ 395.56	<a href="#">+Info</a>
034-7000339	KEYSTONE The Machine IV - 220 Volt vacuum forming machine - (7000339)	US\$ 384.04	<a href="#">+Info</a>

### VACUUM FORMER UNITS - BIOART

Code	Name	Price	+
124059	PLASTVAC Vacuum Forming 110v	US\$ 374.83	<a href="#">+Info</a>
470731	PLASTVAC P7 Plus Vac Dig 110v	US\$ 401.35	<a href="#">+Info</a>
470732	PLASTVAC P7 Plus Vac Dig 220v	US\$ 380.24	<a href="#">+Info</a>
473333	PLASTVAC P7 Vac 110v OPEN BOX	US\$ 374.83	<a href="#">+Info</a>
473831	NEW PlastPress Vacuum	US\$ 841.27	<a href="#">+Info</a>

### Vacuum forming Units - BIOSTAR

Vacuum forming Units - BIOSTAR

Code	Name	Price	+
416696	MiniSTAR S W/ Scan Technology 220v	US\$ 2,749.99	<a href="#">+Info</a>
466302	MiniSTAR S W/ Scan Technology 110v	US\$ 3,449.94	<a href="#">+Info</a>

### Vacuum forming Units - BUFFALO

Vacuum forming Units - BUFFALO

Code	Name	Price	+
008-80165	Tray-Vac Vacuum Forming Machin	US\$ 407.60	<a href="#">+Info</a>
008-80186	Vac Form Heavy Duty 110V	US\$ 481.59	<a href="#">+Info</a>
008-80186-1	Vac Form Heavy Duty 220V	US\$ 512.33	<a href="#">+Info</a>
008-80190	Sta-Vac Premium Electric Vac	US\$ 629.85	<a href="#">+Info</a>
424738	Sta-Vac Mini-lab - Precision, High-	US\$ 763.10	<a href="#">+Info</a>
425216	Tray-Vac Vacuum Former, 220V AC, 7	US\$ 492.03	<a href="#">+Info</a>
425217	Sta-Vac II Mini-lab - Precision, Hi	US\$ 838.50	<a href="#">+Info</a>
425218	AccuVac Accu-Vac Vacuum Former, 220	US\$ 622.97	<a href="#">+Info</a>
425220	Sta-Vac II Mini-lab - Precision, Hi	US\$ 828.10	<a href="#">+Info</a>

### Vacuum forming Units - ERKODENT

Vacuum forming Units - ERKODENT

Code	Name	Price	+
420995	Erkoform 3D Motion 220V	US\$ 3,232.85	<a href="#">+Info</a>

420996	ERKODENT Occluform-3	US\$ 866.00	<a href="#">+Info</a>
--------	----------------------	-------------	-----------------------

## Vacuum forming Units - JSP

Vacuum forming Units - JSP

Code	Name	Price	+
461866	Jsp - Vacuum Forming Machine 220 Vo	US\$ 180.00	<a href="#">+Info</a>

## Vacuum forming Units - KEYSTONE

Vacuum forming Units - KEYSTONE

Code	Name	Price	+
100872	Pressure Dome 5x5 plastics	US\$ 132.80	<a href="#">+Info</a>

## Catalog Index / Indice

Vacuum forming Units ..... 1