

Vacuum formers

Vacuum forming Units

VACUUM FORMER UNITS - BI OART



PlastPress

PlastPress

Code	Name	Price	+
473831	NEW PlastPress Vacuum 110v	US\$ 940.24	+Info
473831-OB	NEW PlastPress Vacuum	US\$ 841.27	+Info

Plastvac P7

Plastvac P7

Plastvac P7

Plastvac P7

Code	Name	Price	+
124059	PLASTVAC Vacuum Forming 110v	US\$ 323.53	+Info
473333	PLASTVAC P7 Vac 110v OPEN BOX	US\$ 263.20	+Info

Plastvac P7 Plus

Plastvac P7 Plus

Code	Name	Price	+
470731	PLASTVAC P7 Plus Vac Dig 110v	US\$ 401.35	+Info
470732	PLASTVAC P7 Plus Vac Dig 220v	US\$ 384.04	+Info

Vacuum forming Units - KEYSTONE

Vacuum forming Units - KEYSTONE



THE MACHINE - KEYSTONE

THE MACHINE - KEYSTONE

The Machine has been a staple in Keystone Industries' long-lasting line of highly recommended vacuum former machines. The price, features and reliability make it the first choice in dental vacuum formers.

The machine has more power than any other dental vacuum former, which means the production of more accurate and definite dental thermoplastic appliance.

The aluminum construction is long-lasting, eliminating the possibility of rust. The heating element is symmetrically engineered for even heat distribution throughout the thermoplastic material.

This machine is used to process bleaching trays, sports mouthguards, implant stints, temporary splints, custom impression trays, denture bases, copings, temporary crowns and bridges and orthodontic appliances.



THE MACHINE III

Code	Name	Price	+
034-1560010	KEYSTONE Machine III Vacuum Formers 110volts, - #7000330	US\$ 384.06	+Info
034-7000331	KEYSTONE The Machine - III - 220 Volt vacuum former machine - vacuum forming -	US\$ 384.06	+Info

PRO-FORM VACUUM FORMER - KEYSTONE

PRO-FORM VACUUM FORMER - KEYSTONE

The Pro-Form is an economical answer that can still excellently fabricate most thermoplastic appliances. The all-aluminum construction eliminates rusting and the heating element is symmetrically engineered for uniform heat distribution. The hexagonal post keeps the frame from rotating or sliding.

Specifications:

- * Uses 5"x5" resin sheets
- * Available in 110 Volt and 220 Volt
- * Available in Single and Dual Chambered
- * Symmetrical heating system
- * Lightweight, rust-free aluminum frame
- * Made in the U.S.A.
- * 5-year warranty on parts and labor



PRO-FORM I

Code	Name	Price	+
101208	PRO-FORM Vacuum Former, 110 Volt, Single Chambered, Aluminum Construction, 5 - #...	US\$ 350.66	+Info
110746	PRO-FORM Vacuum Former, 220 Volt for international plugs/outlets. #9594480	US\$ 371.50	+Info

PROLINE AIR VACUUM 2000 VENTURI - KEYSTONE

PROLINE AIR VACUUM 2000 VENTURI - KEYSTONE

Don't let the size of this machine misguide you. The powerful pressure unit provides superior adaptation for dental plastics using laboratory or operatory internal air compressors. Because the unit has no motor, it is compact and very quiet. Also, the vacuum pressure from an office compressor is usually greater than the pressure generated from smaller vacuum units. The greater the pressure, the better definition and adherence the dental appliances will have.

Specifications:

- * Can be used to layer appliances
- * Very economical when compared to its European counterparts
- * Uses 5"x5" resin sheets or round resin sheets
- * 110 Volt
- * Symmetrical heating system
- * Lightweight, rust-free aluminum frame
- * Made in the U.S.A.
- * 5-year warranty on parts and labor



Us Dental Depot Is An Authorized Keystone Dealer In US

Code	Name	Price	+
------	------	-------	---

034-7000470	KEYSTONE Proline Air Vac 2000 - 110vol - vacuum former - made in usa - 5 year w...	US\$ 545.94	+Info
034-7000471	KEYSTONE PRoline Air Vacuum - 220volts vacuum former - made in usa - 5 year warr...	US\$ 537.54	+Info

Other Models

Other Models

Code	Name	Price	+
100872	Pressure Dome 5x5 plastics	US\$ 147.22	+Info

Vacuum forming Units - BUFFALO

Vacuum forming Units - BUFFALO



Econo-Vac

Econo-Vac

Code	Name	Price	+
008-80186	Vac Form Heavy Duty 110V	US\$ 431.61	+Info
008-80186-1	Vac Form Heavy Duty 220V	US\$ 520.01	+Info

Sta-Vac

Sta-Vac

Code	Name	Price	+
008-80190	Sta-Vac Premium Electric Vac	US\$ 671.23	+Info
424738	Sta-Vac Mini-lab - Precision, High-	US\$ 833.00	+Info

Tray-Vac

Tray-Vac

Code	Name	Price	+
008-80165	Tray-Vac Vacuum Forming Machin	US\$ 442.20	+Info
425216	Tray-Vac Vacuum Former, 220V AC, 7	US\$ 499.62	+Info

AccuVac

AccuVac

Code	Name	Price	+
425218	AccuVac Accu-Vac Vacuum Former, 220	US\$ 634.92	+Info

Sta-Vac II

Sta-Vac II

Code	Name	Price	+
425217	Sta-Vac II Mini-lab - Precision, Hi	US\$ 903.00	+Info
425220	Sta-Vac II Mini-lab - Precision, Hi	US\$ 903.00	+Info

ECCO VACUUM FORMING - BESQUAL

VACUUM FORMING - BESQUAL

- * Vacuum Forming machine 220volts or 110volts.
- * For all thermoplastic laminates.
- * Square and round sheets suitable.
- * Uniform heating distribution - 2 Year Warranty thermoformer vacuum forming.



Shipping information: 25?21?30cm - 6Kg - 15lb

Code	Name	Price	+
102348	BESQUAL Meta Dental Vacuum Forming Machine 110V. Heavy duty spring and latch mec...	US\$ 193.50	+Info
102510	BESQUAL VAC - ESQUAL Meta Dental Vacuum Forming Machine 220V. Heavy duty spring ...	US\$ 193.50	+Info

Vacuum forming Units - BIOSTAR

Vacuum forming Units - BIOSTAR

Code	Name	Price	+
416696	MiniSTAR S W/ Scan Technology 220v	US\$ 2,749.99	+Info
466302	MiniSTAR S W/ Scan Technology 110v	US\$ 3,449.94	+Info

Vacuum forming Units - ERKODENT

Vacuum forming Units - ERKODENT



Erkoform-3D

Erkoform-3D

Code	Name	Price	+
------	------	-------	---

420995	Erkoform 3D Motion 220V	US\$ 3,232.85	+Info
<p>Occluform occluform</p>			
Code	Name	Price	+
420996	ERKODENT Occluform-3	US\$ 866.00	+Info

Catalog Index / Indice

Vacuum forming Units 1