

Vacuum Mixers



Vacuum Mixers Units

TWISTER MIXING UNITS - RENFERT

TWISTER MIXING UNITS - RENFERT

Function & Performance

Twister and Twister Venturi are compact mixers with easy-to-use standard functions and powerful motor technology.

Advantages

- * Optimized vacuum mixing technology produces homogeneous mixing results.
- * Virtually maintenance-free components.
- * Automatic bowl coupling enables convenient single-hand operation.

Features

The units generally all have the same functions. The only difference is in the technology for generating the vacuum. While the Twister operates using a motor pump, the Twister venturi produces the vacuum with the aid of compressed air.

Direct manual control

All mixing parameters are entered easily and quickly.

Impressive motor power

The powerful motor guarantees homogeneous and bubble-free mixing with dental plasters, investments and silicones.

Vacuum level can be regulated

The vacuum level, which can be regulated between 70% (with Twister venturi 80%) and 100%, prevents boiling bubbles with certain plasters and duplicating silicones.

Universal use

The Twister units can be used for all the usual impression and model materials in dental technology.

#UsddRenfert



TWISTER MIXER - (vacuum with internal Pump)

Code	Name	Price	+
116580	RENFERT - TWISTER MIXER -Vacuum mixer 110-220volts (bi-Voltage) #1826-1000 #1...	US\$ 1,594.18	+Info

TWISTER VENTURI (vacuum with compressor)

Code	Name	Price	+
116579	RENFERT - TWISTER VENTURI Vacuum Mixer 110/220volts-Vacuum mixer-15 x 12 x 11in...	US\$ 1,335.77	+Info

TWISTER EVOLUTION VENTURI - (Programable + Compressor Vacuum)

Code	Name	Price	+
120699	RENFERT - TWISTER EVOLUTION VENTURI (DISCONTINUED) vacuum mixer 110/220v...	US\$ 930.08	+Info

TWISTER EVOLUTION MIXER - (Progamable + Internal Pump)

Code	Name	Price	+
116581	RENFERT - TWISTER EVOLUTION Mixer-vacuum mixer 110-220volts (bi-Voltage) #182...	US\$ 1,692.21	+Info

Vacuum Mixers Units - MESTRA

Vacuum Mixers Units - MESTRA



Code	Name	Price	+
419936	IRIS II Vacuum Mixer	US\$ 980.10	+Info
419937	IRIS II Evolution Inverter Vacuum	US\$ 1,105.65	+Info
419938	IRIS II Advanced Vacuum Mixer	US\$ 1,179.90	+Info

VPM2 VACUUM MIXING UNIT - WHIP MIX

VPM2 VACUUM MIXING UNIT - WHIP MIX

Features:

The VPM2 is designed for vacuum mixing all types of gypsums, plasters and investments. You can enter up to 30 separate programs using several parameters to custom mix a material. Utilizing multiple stages, the VPM2 can be programmed for time, speed and paddle direction.

Description:

- * Pre-mix stage blends powder and liquid using alternating paddle direction, eliminating the need for hand mixing
- * Twin chamber, oil-less vacuum pump offers a powerful vacuum quickly, providing a bubble-free mix
- * Vortex inducing paddle design quickly and efficiently incorporates mix
- * Inlet filter prevents material from entering vacuum pump
- * Display programmable in 5 languages (English, German, French, Italian and Spanish)
- * Compact, space saving design



VACUUM MIX WITH BRACKET

Code	Name	Price	+
045-29339	WHIP MIX VPM2 Complete Vacuum Mixing Unit - 110v With Bracket- 30 separate progr...	US\$ 2,972.65	+Info

VACUUM MIX WITH STAND

Code	Name	Price	+
045-29340	WHIP MIX VPM2 Programmable Vacuum Mixing Unit - 110v With Stand - 30 separate pr...	US\$ 3,169.31	+Info
045-29336	WHIP MIX VPM2 Complete Vacuum Mixing Unit - 220volts With Stand - 30 separate pr...	US\$ 3,169.31	+Info

**VENTURI VACUUM MIXERS - WORL
LDENT**

**VENTURI VACUUM MIXERS - WORL
DENT**

The VENTURI SYSEM Slim Line Spin Cycle vacuum mixer is a must have in any dental Lab. The VENTURI System features a compact smart design with hands free operation and a built in digital timer with auto tone for that extra convenience.

SPIN CYCLE VACUUM MIXER SLIM LINE EDITION



Code	Name	Price	+
100498	WORLDDENT SPIN CYCLE SL - VENTURI - Stainless Vacuum Mixer Slim Line (Mezcladora ...	US\$ 1,304.27	+Info
101666	WORLDDENT SPIN CYCLE GRAND - VENTURI - Stainless Vacuum Mixer with vibrator base ...	US\$ 1,382.81	+Info

ELITE MIXERS - ZHERMACK

ELITE MIXERS - ZHERMACK

Advantages:

- * Programming Fast, simple and intuitive thanks to the LCD display and the control lever
- * 10 customizable programs, 6 of these pre-configured for vacuum pump
- * Zhermack casts "oil free" without maintenance 2 year warranty

Features:

- * Electronic control by microprocessor with alphanumeric
- * LCD display control system multifunction control (MFC) 10 programs with the following configurable parameters: mixing time
- * Mixing speed of 100 to 450 rpm function pre-mix pre-vacuum connection interval for changing the direction of rotation of the fluted spreader or offset spatula in electronic vacuum control and double filter mounting on wall or on a bench by means of special support (available as accessory) glasses of mixture of 200 and 1000 cc (accessories)
- * Magnetic Sensor detection of the vessel System "safety" for manual removal of the vessel in case of power failure



Code	Name	Price	+
118507	ZHERMACK EliteMix (115 volt) - Z#C305122	US\$ 1,876.49	+Info

HATHO POLISHING KIT

Thermoplastics like *never* before



Packing Diamond D[®] Denture Acrylic with Tecnoflask



Keystone Industries: Laboratory

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



US DENTAL DEPOT

+1 954 874 6325
www.usdentaldepot.com



Vacuum Mixer units - AMANN GI RRBACH



SMARTMIX X2 - AMANN GIRRACH

AMANN GIRRACH - Smartmix X2 Vacuum
Mixer 110V - Vacuum Mixers Units
115700 - Mfg # 115700

Features:

Perfect mixture results thanks to patented
mixing blade geometry and amazingly simple
operating concept

The universal and futureproof
vacuum-mixing machine with the amazingly
simple operating concept.

Advantages!

- Multifunctional turn/ push button for simple
and quick programming and programme
selectio
- Easy operation - error-free working for
everyone
- Reproducible mixing ratio at the push of
a button, accurate to within a gram
- Integrated scales automatically measure
the dosage or for manual weighing
- Enables clean and dust free operation



When processing homogenous plaster, investment material or paste silicon and impression material, the Smartmix vacuum-mixing machine is an indispensable and productive piece of equipment.

Time, material and plaster savings of 20 to 25%. Finely sprinkled plaster using the new patented plaster dosage technology.

VACUUM MIXER unit

Code	Name	Price	+
115700	AMANN GIRRBACH - Smartmix X2 Vacuum Mixer 110V - Vacuum Mixers Units # 115700 - ...	US\$ 2,285.82	+Info

SMARTBOX X2 - AMANN GIRRBACH

smartbox(R) X2 The most homogenous plaster dispensing method yet Weighing and mixing plaster and water by hand is too imprecise and time-consuming by the high standards of the dental laboratory. Aside from ensuring a consistent high degree of dosage precision, the Smartbox is easy to operate and pays for itself within a few months.

Using the new dosage technology, the plaster is sprinkled very finely into the mixing beaker, which noticeably increases the homogeneity of the mixture.



Code	Name	Price	+
116170-2	Smartbox X2 Dispenser Mixer 115v	US\$ 3,645.97	+Info
116170	Smartbox X2 Dispenser Mixer 230v	US\$ 3,645.97	+Info

Vacuum Mixers

Vacuum mixers Accessories



Vacuum Mixers Accessories - WHIP-MI

X

WHIP MIX - Mixing Bowl 875ml - Gypsum
Vacuum Mixers Accessories
045-08834 - Mfg # 08834

Features:

The WHIP MIX VAC-U-MIXER allows thorough spatulating of stones, plasters, die materials, investments and alginates under vacuum to produce smooth, dense, detailed casts or impressions.

Includes:

- * 1X #6600 Mixing Bowl 875ml (800g Capacity)
- * 1x #6607 Paddle Assembly
- * 1x #6375 Drive Nut

Every VAC-U-MIXER size can be driven by both the COMBINATION UNIT or the VACUUM POWER MIXER



PLUS.

The VAC-U-MIXER is available in five sizes to satisfy the needs of any office or laboratory. These mixers provide a range of capacities from 25grams to 1000 grams. Sturdy construction ensures trouble-free operation; bowls are rigid clear plastic or stainless steel (#7600, 1200 cc). All models have stainless steel paddles and shafts.

Easy Operation

Simply insert drive nut in the drive chuck of your vacuum-mixing unit and the VAC-U-MIXER will do the rest. To avoid contamination, separate VAC-U-MIXERS should be used for mixing phosphate investments, alginates and gypsum materials.

VPM2 MIXER ACCESSORIES

Other Models

Other Models

Code	Name	Price	+
045-29343	VPM2 Vacuum Mixer Stand	US\$ 274.91	+Info

VPM2 MIXING BOWLS

Code	Name	Price	+
045-29348	VAC-U-MIXER #7685 Mixing Bowl 1200ml - (1000g capacity) - #29348	US\$ 324.47	+Info
045-29344	WHIP MIX - VAC-U-MIXER #4485 Mixing Bowl 300ml - (150g capacity) - #29344	US\$ 211.73	+Info
045-29345	VAC-U-MIXER #6585 Mixing Bowl 500ml - (350g capacity) - #29345	US\$ 290.04	+Info
045-29346	WHIP MIX - VAC-U-MIXER #6685 Mixing Bowl 875ml - (650g capacity) - #29346	US\$ 290.04	+Info

VAC-U-MIXER ACCESSORIES

VAC-U-MIXERS - WHIP MIX

The Whip Mix VAC-U-MIXER allows thorough spatulating of stones, plasters, die materials, investments and alginates under vacuum to produce smooth, dense, detailed casts or impressions. Every VAC-U-MIXER size can be driven by both the COMBINATION UNIT or the VACUUM POWER MIXER PLUS.

Easy Operation

Simply insert drive nut in the drive chuck of your vacuum-mixing unit and the VAC-U-MIXER will do the rest. To avoid contamination, separate VAC-U-MIXERS should be used for mixing phosphate investments, alginates and gypsum materials. The #4450, 300 mL size with a Teflon-coated paddle for easy cleanup is



recommended for alginates.

The VAC-U-MIXER is available in five sizes to satisfy the needs of any office or laboratory. These mixers provide a range of capacities from 25 grams to 1000 grams. Sturdy construction ensures trouble-free operation; bowls are rigid clear plastic or stainless steel (#7600, 1200 cc). All models have stainless steel paddles and shafts.

Provide a description of your present equipment and materials to be mixed when you order so that proper drive nuts will be supplied.

VAC-U-MIXER Complete Bowls

Mixing bowl complete with Paddle Assembly and Drive Nut

Code	Name	Price	+
045-08729	VAC-U-MIXER #4450 Mixing Bowl 300ml - (150g capacity) - with #4507 Paddle Assemb...	US\$ 243.21	+Info
045-08710	VAC-U-MIXER #4450 Mixing Bowl 300ml - (150g capacity) - 3 scoops of Alginate wit...	US\$ 247.16	+Info
115310	VAC-U-MIXER #4450 Mixing Bowl 300ml - (150g capacity) - with #4455 Paddle Assemb...	US\$ 243.21	+Info
045-08664	VAC-U-MIXER #6400 Mixing Bowl 200ml - 2 Scoop of Alginate with #6406 Paddle Asse...	US\$ 243.21	+Info
045-08672	VAC-U-MIXER #6400 Mixing Bowl 200ml with #4007 Paddle Assembly and #4109 Drive N...	US\$ 243.21	+Info
045-15652	VAC-U-MIXER #6615 MIXING BOWL 500ml - WITH #6506 PADDLE ASSEMBLY and #6393 Dri...	US\$ 274.14	+Info
045-08818	VAC-U-MIXER #6500 MIXING BOWL 500ml WITH #6508 FLAT PADDLE ASSEMBLY and #6375...	US\$ 243.21	+Info
115311	VAC-U-MIXER #6500 MIXING BOWL 500ml WITH #6506 PADDLE ASSEMBLY and #6375 Driv...	US\$ 243.21	+Info
045-15660	VAC-U-MIXER #6615 Mixing Bowl 875ml - (800g Capacity) with #6606 Paddle Assemb...	US\$ 274.14	+Info
045-08834	VAC-U-MIXER #6600 Mixing Bowl 875ml - (800g capacity) - with #6606 Paddle Assemb...	US\$ 243.21	+Info
115309	VAC-U-MIXER #7600 Mixing Bowl 1200ml - (1000g capacity) With (#7606 Paddle Asse...	US\$ 283.50	+Info

FLEXIBLE VAC-U-MIXER

Some items are accessories for VPM2

Code	Name	Price	+
120174	VAC-U-MIXER WHIP MIX #6750 Flexible Vac U Mixer w/6756 & 6865 875ml/650g - #15667	US\$ 475.81	+Info
120108	VAC-U-MIXER WHIP MIX #6790 Flexible Vac U Mixer Vpm2 875ml/650g - #29355	US\$ 294.98	+Info
120172	VAC-U-MIXER WHIP MIX #6720 Flexible Vac U Mixer w/6726 & 6865 5 Inch - #15665	US\$ 298.40	+Info
120176	VAC-U-MIXER WHIP MIX #6711 Spatula Plastic Round - #15680	US\$ 12.67	+Info

MIXING BOWL ONLY

Code	Name	Price	+
120203	WHIP MIX - VAC-U-MIXER - #6602 Plastic Bowl only 875ml - #09032	US\$ 97.09	+Info
045-09067	VAC-U-MIXER #7602 Stainless Steel Mixing Bowl 1200ml Size#7602 Stainless Steel B...	US\$ 105.92	+Info
045-8974	VAC-U-MIXER #6502 Plastic Bowl 500ml Capacity #08974 - Vac-U-Mixer For #6500- 50...	US\$ 79.94	+Info
045-07609	WHIP MIX VAC-U-MIXER #4452 Plastic Bowl 300ml for #4450 - 3inch wide x 3 inches ...	US\$ 72.24	+Info

PADDLE ASSEMBLY

Each Paddle comes with a different Drive Nut

Code	Name	Price	+
120201	VAC-U-MIXER WHIP MIX #7606 Paddle w/6865 Dn Rh - #09075	US\$ 183.46	+Info
120196	WHIP MIX #6606 Paddle w/6393 Dn Rh - #15644	US\$ 165.49	+Info
120198	VAC-U-MIXER - #6606 Paddle w/6375 Dn Rh - #09040	US\$ 165.49	+Info
045-15601	VAC-U-MIXER #6508 Paddle Assembly (Rh Thread) With Drive Nut #6393 For Use With ...	US\$ 148.32	+Info
045-8982	VAC-U-MIXER #6506 PADDLE ASSEMBLY FOR GYPSUMS With #6393 Drive Nut (Rh Thread) ...	US\$ 148.32	+Info
045-9008	VAC-U-MIXER #6508 PADDLE ASSEMBLY FLAT FOR ALGINATES With #6375 Drive Nut (Rh T...	US\$ 148.32	+Info
120193	VAC-U-MIXER WHIP MIX #6506 Paddle w/6393 Dn Rh - #15601	US\$ 148.32	+Info
120194	VAC-U-MIXER WHIP MIX #6508 Paddle w/6375 Dn Rh - #09008	US\$ 148.32	+Info
045-07617	VAC-U-MIX #4007 Paddle Assembly 200ml Size With Drive Nut LH #4109 #07617	US\$ 125.97	+Info
045-08907	VAC-U-MIXER #6406 Paddle Assembly for 200ml Size With Drive Nut RH #6375 #08907	US\$ 125.97	+Info
120162	VAC-U-MIXER WHIP MIX #6406 Paddle w/6375 Dn Rh - #08907	US\$ 125.97	+Info

OTHER VACUUM MIXER ACCESSORIES

Code	Name	Price	+
045-29342	WHIP MIX VPM2 Vacuum Mixer Wall Bracket #29342	US\$ 88.55	+Info
120130	WHIP MIX #6106 Vibrator Arm WHITE (new 12/2011) - #06548	US\$ 80.33	+Info
120158	WHIP MIX #5709 Transparent Tubing 5 Ft Piece - #09211	US\$ 73.56	+Info
120153	WHIP MIX #5701 Metal Trap Cap - #09172	US\$ 53.27	+Info
120233	WHIP MIX #6014-a Motor Switch Single Pole - #06246	US\$ 51.67	+Info
120152	WHIP MIX #5703 Gauze F/plastic Trap pk 6 - #09180	US\$ 14.36	+Info
120157	WHIP MIX #5704 Spring F/plastic Trap pk 1 - #09199	US\$ 5.48	+Info

VAC-U-CARE MAINTANCE KIT

Code	Name	Price	+
120169	WHIP MIX #4005 O Ring #3 pk 6 - #07382	US\$ 31.14	+Info
120134	WHIP MIX Vv Lubricant for all vacuum mixer units #56 4oz - #06475	US\$ 23.13	+Info
120170	WHIP MIX #4079 Tube Of Lubriplate pk 2 - #07420	US\$ 17.70	+Info
120148	WHIP MIX #6012-15 Felt Filter F/oiler Set 3 - #06408	US\$ 18.07	+Info
120122	WHIP MIX Vac U Care Maintenance Kit - #28398	US\$ 148.86	+Info

TWISTER VENTURI ACCESSORIES - RENFERT

TWISTER VENTURI ACCESSORIES - RENFERT

These accessories will make your dental units more efficient and improve de quality and function of your machine

#UsddRenfert



Code	Name	Price	+
113794	RENFERT - Stand Only for Twister Venturi - Stand Only Converts the Twister to a...	US\$ 332.89	+Info
120695	RENFERT - Mixing Spatula For Twister Bowls - #18210200	US\$ 21.28	+Info

IRIS Accessories | MESTRA

WHIP MIX - Mixing Bowl 875ml - Gypsum
Vacuum Mixers Accessories
045-08834 - Mfg # 08834

Features:

The WHIP MIX VAC-U-MIXER allows thorough spatulating of stones, plasters, die materials, investments and alginates under vacuum to produce smooth, dense, detailed casts or impressions.

Includes:

- * 1X #6600 Mixing Bowl 875ml (800g Capacity)
- * 1x #6607 Paddle Assembly
- * 1x #6375 Drive Nut



Every VAC-U-MIXER size can be driven by both the COMBINATION UNIT or the VACUUM POWER MIXER PLUS.

The VAC-U-MIXER is available in five sizes to satisfy the needs of any office or laboratory. These mixers provide a range of capacities from 25grsams to 1000 grams. Sturdy construction ensures trouble-free operation; bowls are rigid clear plastic or stainless steel (#7600, 1200 cc). All models have stainless steel paddles and shafts.

Easy Operation

Simply insert drive nut in the drive chuck of your vacuum-mixing unit and the VAC-U-MIXER will do the rest. To avoid contamination, separate VAC-U-MIXERS should be used for mixing phosphate investments, alginates and gypsum materials.

Code	Name	Price	+
419943	MESTRA IRIS Mixing Bowl for Vacuum 250cc #080523	US\$ 201.00	+Info
419944	MESTRA IRIS Mixing Bowl for Vacuum 550cc #080524	US\$ 226.50	+Info
419945	MESTRA IRIS Mixing Bowl #3 for Vacuum 550cc polycarbonate #080524P	US\$ 153.93	+Info
419939	MESTRA Suport for IRIS II Vacuum Mixer #080519	US\$ 228.00	+Info
419946	MESTRA IRIS Mixing Bowl #4 for Vacuum 750cc #080525	US\$ 252.00	+Info
419947	MESTRA IRIS Mixing Bowl for Vacuum 1000cc #080526	US\$ 252.00	+Info
419948	MESTRA IRIS Mixing Bowl for Vacuum 1250cc #080527	US\$ 252.00	+Info

Vacuum mixers Accesories - AMANN GIRRBACH

Vacuum mixers Accesories - AMANN GIRRBACH



FILTERS

Code	Name	Price	+
115701	AMANN GIRRBACH Smartmix X2 Spare Filter 5pcs #AG 115701	US\$ 69.84	+Info

SMARTMIX MIXING BOWLS - AMANN GIRRBACH

AMANN GIRRBACH Mixing Bowl with Blades

The Amann GIRRBAch Mixing bowls also know as Mixing Beakers are to be used with The Smartmix X2. Each mixing bowl come complete with blades.

The Amann GIRRBAch Vacuum Mixing bowls are also available with or without the blades.



Features:

Available in 5 Different Capacities
Mixing bowl with Blades

#UsddAmann

AMANN GIRRBACH is a trademark of Amann GIRRBAch AG

Code	Name	Price	+
115686	AMANN GIRRBACH Mixing Bowl with Blade 100ml Capacity #AG 115620 - Vacuum mixer 1...	US\$ 286.65	+Info
115687	AMANN GIRRBACH Mixing Bowl with Blade 250ml Capacity #AG 115630 - Vacuum mixer 2...	US\$ 257.98	+Info
115688	AMANN GIRRBACH Mixing Bowl with Blade 500ml Capacity #AG 115640 - Vacuum mixer 5...	US\$ 302.72	+Info
115689	AMANN GIRRBACH Beaker with Mixing Blade - 750ml Capacity #AG 115650 - Vacuum mix...	US\$ 320.14	+Info
115690	AMANN GIRRBACH Mixing Bowl with Blade 1000ml Capacity #AG 115660 - Vacuum mixer ...	US\$ 282.41	+Info
115694	AMANN GIRRBACH Mixing Bowl Only 1000ml Capacity #AG 115661 - Vacuum Mixer 1000ml...	US\$ 225.04	+Info
115693	AMANN GIRRBACH Mixing Bowl 750ml Capacity #AG 115651 - Vacuum Mixer 750ml Beaker	US\$ 203.59	+Info
115692	AMANN GIRRBACH Mixing Bowl 500ml Capacity #AG 115641 - Vacuum Mixer 500ml Beaker	US\$ 167.43	+Info
117201	AMANN GIRRBACH Rising Socket for Beakers - 2piece Set #AG 117201	US\$ 86.99	+Info

FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



**US
DENTAL
DEPOT**

+1 954 874 6325 
www.usdentaldepot.com 

Vacuum mixers Accesories - BEGO

Vacuum mixers Accesories - BEGO



Mixing Bowls MOTOVA



Code	Name	Price	+
10-16271	BEGO Mixing Bowl 250ml for Motova 1pk - Investment material and stone #16271 #16271	US\$ 435.08	+Info
10-16274	BEGO Mixing Bowl 1200ml for Motova 1pk - Investment material and stone #16274 #1...	US\$ 635.64	+Info
10-16275	BEGO Mixing Bowl 550ml for Motova 1pk - Silicone material #16275 #16275	US\$ 529.52	+Info
418965	BEGO 650ML REPLACEMENT BOWL- MOTOVA MP/MPV - #11097	US\$ 63.36	+Info
418981	BEGO LID FOR 650ML BOWL- MOTOVA MP/MPV - #12158	US\$ 92.33	+Info
419164	BEGO VACUUM BOWL 650ML FOR MOTOVA MP AND SLA - #14792	US\$ 551.23	+Info
419165	BEGO COVER w/ STIRRER FOR 650ML VACUUM BOWL- MOTOVA MP/MPV - #14796	US\$ 768.24	+Info
419674	BEGO O-RING - EASYMIX BOWLS- DUOSTAR NEW- EASYMIX- MOTOVA SLA - #35253 (SPECIAL...	US\$ 6.00	+Info

WORLDENT MIXING BOWLS

Worldent Mixing Bowls come in various selections...

Sizes:

- 200ml
- 500ml
- 900ml
- 1000ml

Set Up:

- Bowl Only
- Bowl and Spatula
- Spatula Only



Code	Name	Price	+
------	------	-------	---

106805	SPIN-CYCLE- Spin Cycle - Mixing Bowl Only - 500ml - replacement part #85B	US\$ 63.40	+Info
106806	WORLIDENT - SPIN-CYCLE- Spin Cycle - Mixing Bowl Only - 900ml - replacement part...	US\$ 62.77	+Info
106809	SPIN-CYCLE- Spin Cycle - Spatula Only - 900ml - replacement part - For Mixing Bo...	US\$ 156.75	+Info
116771	WORLIDENT Spin Cycles - Bowl & Spatula 200 MI (8200)	US\$ 166.78	+Info
116772	Spin Cycles - Bowl & Spatula 500 MI (8500)	US\$ 184.52	+Info
116773	WORLIDENT Spin Cycles - Bowl & Spatula 900 MI (8900)	US\$ 181.96	+Info
116789	Spin Cycles - O-Rings 200M1 Bowl (82Or)	US\$ 10.30	+Info
116790	Spin Cycles - O-Rings 500M1 Bowl (85Or)	US\$ 10.51	+Info
116791	Spin Cycles - O-Rings 900M1 Bowl (89Or)	US\$ 10.51	+Info
116777	Spin Cycles - Spatula Only 200 MI (82Sp)	US\$ 137.97	+Info
116778	Spin Cycles - Spatula Only 500 MI (85Sp)	US\$ 138.93	+Info

Vacuum Mixers

Dispenser Units



AQUASPENSE LIQUID DISPENSE R - WHIP MIX

AQUASPENSE LIQUID DISPENSER - WHIP MIX

The AquaSpense ensures proper water powder ratio's by automatically calculating the powder weight and then accurately dispensing the correct amount of water. Automating the dispensing of water to powder takes the human error out of the process creating more predictable working times, expansion, set times and hardness. The result is less re-work and increased profits.

Programmable for up to 15 various gypsums, the user will have the flexibility needed to standardize their lab processes for optimal product performance and results.

The AquaSpense SL unit dispenses both water and special liquid for your investing needs. This innovative unit takes the guess work out of measuring the powder and the water/liquid ratios by weighing the bowl with the powder on a calibrated scale and calculating the proper liquid for the volume. By dispensing the exact amount of water/liquid for the investment, inconsistencies and resulting remakes are virtually eliminated.

Predictable and accurate results save time and money, the AquaSpense and AquaSpense SL are sure to become indispensable in your lab.

Description:

- * Produces accurate casting sizes with less finishing time
- * Eliminates remakes due to human error and



variability in the measuring process

* Calculates the powder weight and accurately dispenses water and/or special liquid, resulting in a consistent mix every time.

* Programmable memory for up to 15 investments

* For the Model Department: Dispense liquids such as Gypsum Hardener or other enhancing liquids in place of the Special Liquids for Investments

AQUASPENSE LIQUID DISPENSER 115/230V

Code	Name	Price	+
045-09560	WHIP MIX AquaSpense Vacuum Mixer 110volts - 15 various gypsums Automatic powder ...	US\$ 2,788.85	+Info
045-09550	WHIP MIX AquaSpense SL Vacuum Mixer 110volts - Water and Special Liquid Dispensi...	US\$ 3,204.35	+Info

ACCESSORIES

Code	Name	Price	+
045-09552	WHIP MIX AquaSpense #5615A Bottle Cap Assembly - Coll Liquid #09552	US\$ 93.52	+Info
045-09556	WHIP MIX AquaSpense #5605A Bottle Cap Assembly - Water #09556	US\$ 20.56	+Info

Vacuum Mixers

Manual Vacuum Mixers



Manual Vacuum Mixers - KEYSTONE

Manual Vacuum Mixers - KEYSTONE

Code	Name	Price	+
034-1625000	Vac Mix Bowl Sm Comp 3-1/2in	US\$ 100.23	+Info
034-1625100	Vac Mix Bowl Med Comp 4-3/4in	US\$ 114.91	+Info
034-1625200	Vac Mix Bowl Lge Comp 5-1/2in	US\$ 123.06	+Info

Catalog Index / Indice

Vacuum Mixers Units	1
Vacuum mixers Accessories	6
Dispenser Units	13
Manual Vacuum Mixers	14