

Vibrators

Large Vibrators

VIBRATORS UNITS - BESQUAL

VIBRATORS UNITS - BESQUAL

- * The Vibrators are designed to promote easy pouring of plaster while preventing bubbles while pouring.
- * Ergonomic design
- * Long wave or short wave options
- * Easily detachable rubber top allows for convenient cleaning of vibration table



MEDIUM VIBRATOR

Code	Name	Price	+
101915	BESQUAL MED-Vibrator Rectangle 6 x 8 x 4.5in (medium) variable Speed - Heavy D...	US\$ 123.75	+Info

BIG VIBRATOR

Code	Name	Price	+
101914	BESQUAL BIG Vibrator 270x 220mm 3 speeds - (110vol ts) 2 frequencies - 300wats...	US\$ 250.96	+Info
114843	BESQUAL BIG Vibrator 270x 220mm 3 speeds - (220volts) 2 frequencies - 300wats ...	US\$ 253.50	+Info

BENCH VIBRATORS # 200 - BUFFALO

BENCH VIBRATORS # 200 - BUFFALO

This product is one of the largest, most powerful, variable speed bench vibrator for multiple flasks and dense castings.

Specifications:

- * Large, Powerful. Variable speed bench vibrator for multiple flasks and dense castings.
- * Large rubber-covered steel work platform improves convenience and spill control.
- * Improved electronics and fully variable speed control promote powerful 350 W, 60 Hz vibration.
- * Heavy beige-enameled cast aluminum housing provides stability and eases clean-up tasks.
- * Dimensions: 10 1/2" W x 7 1/2" D x 4 1/2" H
- * Weight: 19 lbs.



Code	Name	Price	+
008-84500	BUFFALO Vibrator #200 - 110v - Extra-Heavy Duty Vibrator, Variable - #84500	US\$ 318.28	+Info
008-84510	BUFFALO Vibrator #200 - 220v - Extra-Heavy Duty Vibrator, Variable - #84510	US\$ 383.46	+Info

EXPORT Vibrators - DENTAL FARM



Code	Name	Price	+
474519	MINIEXPORT	US\$ 735.02	+Info
474520	VIBRATOR EXPORT 2	US\$ 981.88	+Info
474521	MAXIEXPORT	US\$ 1,907.33	+Info

HEAVY DUTY VIBRATOR - WHIP MIX

HEAVY DUTY VIBRATOR - WHIP MIX

The Heavy Duty Vibrators is designed to promote easy pouring of plaster while preventing bubbles while pouring.

Description:

- * Strong vibration, with large vibrator and wide vibration table, enables simultaneous work on many rings
- * Vibration reduces bubbles
- * Ergonomic design.
- * Long wave or short wave options Easily detachable rubber top allows for convenient cleaning of vibration table



Code	Name	Price	+
115683	WHIP MIX Heavy Duty Vibrator 110volts 50 Hz #10675 - Dimensions - 11in Wide x 4...	US\$ 425.98	+Info

115684	WHIP MIX Heavy Duty Vibrator 220vol 50/60 Hz #10680 - Dimensions - 11in Wide x 4...	US\$ 393.36	+Info
--------	---	-------------	-----------------------

Large Vibrators - ESSENCE

Large Vibrators - ESSENCE

Code	Name	Price	+
124487	VIBRAMAXX Big Vibrator 110v	US\$ 160.00	+Info

Large Vibrators - SHOFU

Large Vibrators - SHOFU

Code	Name	Price	+
116917	Ceramosonic S	US\$ 678.83	+Info

Large Vibrators - WORLDDENT

Large Vibrators - WORLDDENT

Code	Name	Price	+
101759	Heavy Duty Vibrator Large 110v	US\$ 451.87	+Info
114844	Heavy Duty Vibrator Large 220v	US\$ 390.88	+Info

VIBRAX - RENFERT

VIBRAX - RENFERT

Function & Performance

The wide-ranging vibration spectrum optimally processes any material. The functionally sophisticated handling enables comfortable use even under difficult conditions.

Advantages

- Bubble-free flow behavior using two wave ranges with 4 levels of intensity

- Reduced transfer of the vibrations to the workbench as the housing is oscillation decoupled.

- Long service life using maintenance-free vibration magnet.

Details

- Uncomplicated adaptation of the intensity on the large lever switch using only one finger.

- Extremely quiet, ensuring a comfortable working atmosphere.

- Time-saving due to quickly cleaned elements .

- Non-tip and extremely stable even when high pressure is applied to the edges.

- No springing of the impression tray thanks to the soft tray rest.



US DENTAL DEPOT is an Authorized Renfert Dealer in US.

VIBRAX 110/220V

Code	Name	Price	+
23-18301000	RENFERT - Vibrax 110 Volts Vibrax- #1830-1000 #18301000-	US\$ 609.46	+Info
23-18300000	RENFERT - Vibrax 220 Volts-Large Vibrator- #1830-0000 #18300000 -for thin fluid...	US\$ 754.32	+Info

ACCESSORIES

Code	Name	Price	+
23-18300001	RENFERT - Vibration Ball for Vibrax- #1830-0001 #18300001	US\$ 138.12	+Info

Vibrators

Small Vibrators

GENERAL PURPOSE VIBRATOR

Code	Name	Price	+
424442	Buffalo Vibrator #1A	US\$ 123.07	+Info
424521	#2 Heavy Duty Vibrat	US\$ 244.46	+Info
424773	Buffalo Vibrator #1A 220v	US\$ 134.19	+Info
425257	#1B small vibrator w	US\$ 174.37	+Info
425258	#1B small vibrator w	US\$ 187.41	+Info
425259	#2 Heavy Duty Vibrat	US\$ 267.70	+Info

GENERAL PURPOSE VIBRATOR - WHIP MIX

GENERAL PURPOSE VIBRATOR - WHIP MIX

The General Purpose is designed to promote easy pouring of plaster while preventing bubbles while pouring.

Description:

- * Vibration reduces bubbles.
- * Ergonomic design.
- * Easily detachable rubber top allows for convenient cleaning of vibration table.



Code	Name	Price	+
120258	WHIP MIX Gen Purpose Vibrator 110v - #10650	US\$ 149.74	+Info

Small Vibrators - BESQUAL

Small Vibrators - BESQUAL

Code	Name	Price	+
101930	Vibrator Round Top 110V	US\$ 72.00	+Info

Small Vibrators - MESTRA

Small Vibrators - MESTRA

Code	Name	Price	+
419316	TEIDE CIMA Electrovibrator	US\$ 445.50	+Info

Small Vibrators - RAY FOSTER

Small Vibrators - RAY FOSTER

Code	Name	Price	+
464545	Vibrator Round Top 110V	US\$ 408.24	+Info
470920	DV34 Variable Speed Vibrator	US\$ 125.44	+Info

Catalog Index / Indice

Large Vibrators	1
Small Vibrators	4