

Waxers Electrical

Accessories for electric waxers



REPLACEMENT WAX TIPS - WHIP MIX

REPLACEMENT WAX TIPS - WHIP MIX

Whip Mix offer 9 different replacement tips for its new line of digital waxing pencils.



Code	Name	Price	+
120313	WHIP MIX Wax Carving Tips 6/pk - #10168	US\$ 191.88	+Info
120314	WHIP MIX Wax Carving Tips 2/pk - #10167	US\$ 58.23	+Info
120317	WHIP MIX Wax Tip #1 Offset Explorer - #10171	US\$ 29.28	+Info
120319	WHIP MIX Wax Tip #2 Pkt1 - #10172	US\$ 29.28	+Info
120312	WHIP MIX Wax Tip #3 Pkt2 - #10173	US\$ 29.28	+Info
120321	WHIP MIX Wax Tip #4 Discoid - #10174	US\$ 29.28	+Info
120322	WHIP MIX Wax Tip #5 Mini Beavertail - #10175	US\$ 29.28	+Info
120323	WHIP MIX Wax Tip #6 Beavertail - #10176	US\$ 29.28	+Info
120316	WHIP MIX Wax Tip #7 Large Beavertail - #10177	US\$ 28.99	+Info
120324	WHIP MIX Wax Tip #8 Pkt7 - #10178	US\$ 29.28	+Info
120320	WHIP MIX Wax Tip #9 Half Hollenback - #10179	US\$ 29.28	+Info

HANDPIECE FOR NT PRO AND EASY WAXER - YETI

HANDPIECE FOR NT PRO AND EASY WAXER YETI

Modeling tip to be pinned
 Heating pin
 Ergonomically shaped wooden handle
 Magnetic rings
 Plastic holder
 Cable/lead
 The handpiece is extremely light and well insulated against heat, remaining cool even after prolonged use
 The handpiece has a highly flexible connection cable/lead



Code	Name	Price	+
34-1860157	YETI Handpiece Pen for all waxers serise NT and pro with unipolar connector #411...	US\$ 123.20	+Info
34-1860158	YETI Handpiece Pen for pro waxer with angle 5 - pin connector #411-1000 Pen Only...	US\$ 123.20	+Info

MODELLING TIPS FOR IQ WAXER R NT LIGHT - YETI

MODELLING TIPS FOR IQ WAXER NT LIGHT YETI

Thanks to micro heating elements and the operation of modern laser technology, the heating was placed at the point of the modelling tip. Therefore, the hottest point of maximal heat of the instrument lies at the point of the modelling tip, so that:

The modelling tips heat up and cool down in fractions of a second.

There is no flow back of the wax drop.



Code	Name	Price	+
417050	YETI Accu Waxer Modelling tip for IQ Waxer NT light PKT 1 fine 10 mm #420-0100	US\$ 36.19	+Info
417051	YETI Accu Waxer Modelling tip for IQ Waxer NT light PKT. 2 thick 12 mm #420-0200	US\$ 33.32	+Info
417052	YETI Accu Waxer Modelling tip for IQ Waxer NT light Spear shaped #420-0400	US\$ 57.44	+Info

MODELLING TIPS FOR WAXER N T AND PROWAXER - YETI

MODELLING TIPS FOR WAXER NT AND PROWAXER YETI

We offer an assortment of six different tips for all our Waxers series NT, pro and light. The tips can be individually modified by grinding. The tips can be quickly changed thanks to the simple pin system. Their interior diameter is 2 mm.



Other Models

Other Models

Code	Name	Price	+
420172	Knife Shape P.K. t 4 Tip	US\$ 23.67	+Info
420173	Corkscrew thin P.K. t 5 Tip	US\$ 40.92	+Info
420174	Corkscrew thick P.K. t 6 Tip	US\$ 40.92	+Info

PK Thomas

PK Thomas

Code	Name	Price	+
115041	Yeti P.K. -1 Thomas Small Tip	US\$ 40.92	+Info
115042	Yeti P.K. -2 Thomas Tip	US\$ 40.81	+Info
115043	Hollenback Horizontal P.K. t 3 Tip	US\$ 40.92	+Info

WAXLECTRIC HANDLE - RENFERT

HANDPIECE REPLACEMENT PARTS

Feature

- * Maximum working comfort provided by heat-insulated, natural cork handle.
- * Heat-resistant cable protected against twisting.



US DENTAL DEPOT is an Authorized RENFERT Dealer in US.

Code	Name	Price	+
2151-1500	RENFERT - Kit for Waxlectric Holders- #2151-1500 #21511500	US\$ 67.04	+Info
2154-0001	RENFERT - Heating Handle (Handpiece only) YELLOW-For Electric Waxer- Phone P...	US\$ 106.18	+Info
2154-0002	RENFERT - Heating Handle (Handpiece only) red-For Electric Waxer- #2154-0002...	US\$ 102.56	+Info

WAXLECTRIC TIPS - RENFERT

Tip small - Intelligent tips

Function & Performance

The integrated heating element, customization option, easy and quick changing are only some of the special features of the sculpting tips.

Advantages

- ? 12 different versions.
- ? Very rapid heating and reheating rates thanks to the extreme thermal conductivity of the special alloy.
- ? Easy cleaning ensures long service life.



US DENTAL DEPOT is an Authorized RENFERT Dealer in US.

Waxlectric

Waxlectric

Code	Name	Price	+
2155-0101	Waxlectric Tip Probe Small	US\$ 56.87	+Info

Waxlectric Tips

Waxlectric Tips

Code	Name	Price	+
2155-0102	New Wax Electric Tip Probe Med	US\$ 56.87	+Info
2155-0103	Waxlectric Tip Probe Large	US\$ 56.87	+Info
2155-0104	Waxlectric Tip Blade Narrow	US\$ 56.87	+Info
2155-0105	Waxlectric Tip Blade Wide	US\$ 56.87	+Info
2155-0106	Waxlectric Tip Chamfered	US\$ 51.02	+Info

2155-0107	Waxlectric Tip Hollenback Hor	US\$ 56.87	+Info
2155-0108	Waxlectric Tip Hollenback Vert	US\$ 56.87	+Info
2155-0109	Waxlectric Tip Probe Mini Ang	US\$ 50.05	+Info
2155-0110	Waxlectric Tip Lance Blade	US\$ 56.87	+Info
2155-0111	Waxlectric Tip Beaver Tail	US\$ 56.87	+Info
2155-0112	Waxlectric Tip Large Spoon	US\$ 88.64	+Info

Accessories for electric waxers - BESQUAL

Accessories for electric waxers - BESQUAL

Code	Name	Price	+
471462	Carver Tip #3 for S500	US\$ 11.76	+Info
471463	Carver Tip #1 for S500	US\$ 12.00	+Info
471464	Carver Tip #2 for S500	US\$ 12.00	+Info
471465	Carver Tip #4 for S500	US\$ 12.00	+Info
471466	Carver Tip #5 for S500	US\$ 12.00	+Info
471472	Carver Tip #7 for S500	US\$ 14.90	+Info
471479	Carver Tip #6 for S500	US\$ 12.00	+Info

Accessories for electric waxers - DENTSPLY

Accessories for electric waxers - DENTSPLY

Code	Name	Price	+
418322	Waxpencil- Tip #7 (3Pk)	US\$ 47.52	+Info
418323	Waxpencil- Tip #6 (3Pk)	US\$ 47.52	+Info
418324	Waxpencil- Tip #8 (3Pk)	US\$ 47.52	+Info
418325	Waxpencil- Tip #2 (3Pk)	US\$ 47.52	+Info
418326	Waxpencil- Tip #3 (3Pk)	US\$ 47.52	+Info

Accessories for electric waxers - RENFERT

Accessories for electric waxers - RENFERT

Code	Name	Price	+
417585	Waxlectric tip Probe Small 1pk	US\$ 55.92	+Info
417586	Waxlectric tip Probe Large 1pk	US\$ 56.00	+Info
417587	Waxlectric tip Blade narrow 1pk	US\$ 56.00	+Info
417588	Waxlectric tip Blade Chamfer 1pk	US\$ 40.82	+Info
417589	Waxlectric tip Hollenback hor 1pk	US\$ 50.86	+Info
417590	Waxlectric tip Hollenback ve 1pk	US\$ 47.55	+Info
417591	Waxlectric Tip mini Probe 1pk	US\$ 55.92	+Info
417592	Waxlectric Tip blade Lance 1pk	US\$ 50.91	+Info
417593	Waxlectric Tip Blade beaver 1pk	US\$ 50.86	+Info
417640	O-ring for Waxlectric Tip 5pk	US\$ 6.38	+Info

Accessories for electric waxers - YETI

Accessories for electric waxers - YETI



Code	Name	Price	
417042	Clamp for handpieces	US\$ 9.05	+Info
417045	Wooden holder for Wxr NT / Light /	US\$ 7.44	+Info

Accessories for electric waxers - AMANN GIRRBACH

Accessories for electric waxers - AMANN GIRRBACH



AMANN GIRRBACH

Other Models

Other Models

Code	Name	Price	
467081	Circuit board for Smartwax Duo	US\$ 366.32	+Info
467082	Handpiece with cord	US\$ 147.01	+Info
467083	Probe small incl. handle	US\$ 66.15	+Info
467084	Probe large incl. handle	US\$ 66.15	+Info
467089	Waxjet for Smartwax	US\$ 270.76	+Info

Smartwax

smartwax

Code	Name	Price	
467087	Knife incl. handle	US\$ 66.15	+Info
467088	Spoon incl. handle	US\$ 66.15	+Info

Accessories for electric waxers - BESQUAL

Accessories for electric waxers - BESQUAL

Code	Name	Price	
471474	Carver Tip #8 for S500	US\$ 12.00	+Info

ACCESSORIES Smartwax Duo US

Accessories for use with the Amann Girschbach smartwax duo.



Code	Name	Price	
467085	Needle incl. handle	US\$ 66.15	+Info
467086	Beaver tail incl. handle	US\$ 66.15	+Info

Waxers Electrical
Electrical Waxer Units



EASY WAXER - YETI

EASY WAXER YETI

The development of the easy Waxer is based on the experiences of recent years and the goal to design an electric wax up unit every lab could afford and handle easily. With an analogue technology, a short heating time and a temperature range up to 250° it has a lot to offer. The highlight of the easy Waxer is definitely the LED indicator for an error diagnosis. In case the unit will show an error in some years dental technicians can find out the source of error easily.



Available Simple/ Duo

Code	Name	Price	+
417049	YETI EASY Waxer DUO 110vol #415-0115	US\$ 290.76	+Info

PRO WAXER - YETI

PRO WAXER YETI

Pro WAXER & pro Waxer duo
Handpieces and tips are not included in the warranty!
The cost-efficient alternative with analog technology.
High performance adjustable from 50 ° C - 225 ° C / 122 ° F - 437 ° F
Coarse display to indicate the setpoint temperature
The handpiece is ultralight, insensitive to temperature, remains pleasantly cool
The connection cable to the handpiece is highly flexible and absolutely insensitive to temperature



Available Simple/ Duo

Code	Name	Price	+
417046	YETI PRO Waxer 110vol #412-0115	US\$ 307.99	+Info
1860001	YETI PRO-Waxer-DUO - 110vol - #414-0115 (Keystone #1860001)	US\$ 331.25	+Info

WAXLECTRIC II - RENFERT

WAXLECTRIC II - RENFERT

Waxlectric. Waxing up in best time.

Renfert has all the answers if you want to sculpt wax stress-free and conveniently. The Renfert heat conductivity concept concentrates the temperature exactly where it is required. The sculpting tip forms the required heat center.



Temperature setting:

* Individual matching to any of your waxes thanks to a temperature setting of 50-200°C (122-392°F).

Sculpting tips:

- * Internal heating of the sculpting tip ensures direct transfer of heat to your wax.
- * Special alloy tips with optimum heat conductivity.

Handpiece:

- * Maximum working comfort provided by heat-insulated, natural cork handle.
- * Heat-resistant cable protected against twisting.

Housing:

- * 20° inclination allows view of the unit and optimum operation.
- * Hermetically sealed, easy to clean membrane board.
- * Increased work safety due to low-voltage technology. Power adaptor is outside the workstation.

#UsddRenfert

Code	Name	Price	+
2157-1000	RENFERT - Waxlectric II Programmable-110 Volts-Electric Wax Spatula of 2 pens- ...	US\$ 432.96	+Info
2157-0000	RENFERT - Waxlectric II Programmable-220 Volts-#2157-0000 #21570000 - (BI VOLT)	US\$ 418.03	+Info

WAXLECTRIC I - RENFERT

Renfert Waxlectric I - Waxing up in best time.

Renfert has all the answers if you want to sculpt wax stress-free and conveniently. The Renfert heat conductivity concept concentrates the temperature exactly where it is required. The sculpting tip forms the required heat center.

Temperature setting:

Individual matching to any of your waxes thanks to a temperature setting of 50-200°C (122-392°F).

Sculpting tips:

- Internal heating of the sculpting tip ensures direct transfer of heat to your wax.
- Special alloy tips with optimum heat conductivity.

Handpiece:

- Maximum working comfort provided by heat-insulated, natural cork handle.
- Heat-resistant cable protected against twisting.

Housing:

- 20° inclination allows view of the unit and optimum operation.
- Hermetically sealed, easy to clean membrane board.
- Increased work safety due to low-voltage technology. Power adaptor is outside the workstation.



#UsddRenfert

Code	Name	Price	+
2156-1000	RENFERT - Waxlectric I Programmable 110 Volts-Electric Wax Spatula of 1 pens-#2...	US\$ 323.52	+Info
2156-0000	RENFERT - Waxlectric I Programmable 220 Volts- #2156-0000 #21560000 Programmabl...	US\$ 309.53	+Info

WAXLECTRIC LIGHT I & II - RENFE RT

WAXLECTRIC LIGHT I & II - RENFERT

Function & Performance

These electric wax knives are the perfect introduction to efficient wax technology.

Advantages

- * Easy temperature control via a potentiometer.
- * Separately adjustable temperature with Waxlectric light II.
- * Compact, ergonomic design.

Features

Unlike the Waxlectric light I, the Waxlectric light II can be operated with two handpieces simultaneously.

Details

- * No contraction in the pattern, as the wax is neither overheated nor burned.
- * The technician can concentrate more on the pattern and less on the tools and instruments.



#UsddRenfert

Code	Name	Price	+
114204	RENFERT - Waxlectric Light I 110 Volts-1 channel-Included waxing-up tips- #2150...	US\$ 236.08	+Info
114206	RENFERT - Waxlectric Light II 110 Volts-2 channel-Included waxing-up tips- #215...	US\$ 237.74	+Info

Electrical Waxer Units - AMANN GIRRBACH

Electrical Waxer Units - AMANN GIRRBACH

Code	Name	Price	+

467080	Smartwax Duo US	US\$ 586.84	+Info
--------	-----------------	-------------	-----------------------

DIGITAL DUAL WAX CARVER UNITS - BESQUAL

DIGITAL DUAL WAX CARVER UNITS - BESQUAL

The Dental instrument is used to carve wax electrically and mold its form, featuring its availability for a variety of carving tips (6pcs), easy temperature control and precise waxing work.



Description:

- * Excellent heat conductivity, lasting service life, and easy use.
- * Less temperature deviation by the constant supply of the DC voltage.
- * Digitally visible temperature with choice of celsius or Fahrenheit.
- * Available to memorize the previous working temperature even after turning off power.

Code	Name	Price	+
114036	BESQUAL Digital Dual Wax Spatula Carver Unit 110vol with 2 pens and 6 tips #S50...	US\$ 238.50	+Info
114039	BESQUAL Digital Dual Wax Spatula Carver Unit 220v with 2 pens and 6 tips #S500N...	US\$ 188.41	+Info

DIGITAL WAX CARVING PENCIL - WHIP MIX

DIGITAL WAX CARVING PENCIL - WHIP MIX

The Digital Wax Carving Pencil is designed to carve and mold a variety of waxes used in the dental industry. The variety of carving tips and easy temperature control allow for precise waxing work.



Description:

- * Easy to use.
- * Excellent heat conductivity.
- * Long service life.
- * Easy-to-use, visible control buttons.
- * Accurate temperature maintenance with reset function.
- * Built-in digital CPU circuit accurately maintains desired temperature.
- * Digitally displays temperature in either Celsius or Fahrenheit.
- * Previous working temperature is kept in memory, even if the unit has been

turned off.

- * Usable with AC100-240V without additional adjustment.
- * Silicon handle reduces hand fatigue.
- * Translucent compartment in rear of unit for easy storage of carving tips.
- * Comes with 2 waxing tips (Offset Explorer and PKT 2)

DIGITAL WAX CARVING PENCIL

Code	Name	Price	+
120310	WHIP MIX Digital Wax Carving Pencil 110v - #10157	US\$ 314.81	+Info

ACCESSORIES

Code	Name	Price	+
120315	WHIP MIX Replacement Touch Carving Pencil W / Cord - #10165	US\$ 158.06	+Info
120318	WHIP MIX Replacement Wax Carving Pencil W / Cord - #10163	US\$ 145.29	+Info

DIGITAL WAX CARVING PENCILS-DUAL - WHIP MIX

DIGITAL WAX CARVING PENCILS-DUAL - WHIP MIX

The Digital Wax Carving Pencils-Dual is designed to carve and mold a variety of waxes used in the dental industry. The variety of carving tips for each pencil and easy temperature control allow for precise waxing work.

Description:

- * Easy to use.
- * Excellent heat conductivity.
- * Long service life.
- * Easy-to-use, visible control buttons.
- * Two-line circuit, designed for simultaneous work.
- * Minimal temperature deviation.
- * Accurate temperature maintenance with reset function.
- * Built-in digital CPU circuit accurately maintains desired temperature.
- * Digitally displays temperature in either Celsius or Fahrenheit.
- * Previous working temperature is kept in memory, even if the unit has been turned off.
- * Usable with AC100-240V without additional adjustment.
- * Silicon handle reduces hand fatigue.
- * Convenient magnetic cradle on surface or external cradle.
- * Translucent compartment in rear of unit for easy storage of carving tips.
- * Comes with 6 waxing tips (Offset Explorer, PKT 1, PKT 2, Discoid, Mini-Beaver Tail, Beaver Tail)



Code	Name	Price	+
120308	WHIP MIX Digital Wax Carving Pencil Dual 110v - #10150	US\$ 550.55	+Info

DIGITAL WAX POT - WHIP MIX

DIGITAL WAX POT - WHIP MIX

The Wax Pot is designed to melt a variety of wax used in the dental industry. The unit can heat wax quickly and maintain its liquid state. The setting temperature can be reduced easily to achieve optimum workability.

Description:

- * Capable of heating four different waxes simultaneously and maintaining their liquid state.
- * Built-in thermostat to accurately maintain desired temperature.
- * Excellent insulation and heat conductivity.
- * Easy to operate.
- * Heats wax quickly.
- * Built-in digital CPU circuit resets and maintains precise desired temperature.
- * Current and preset temperature readings are visible.
- * LED indicates automatic heater power and turns off when unit has reached the reset temperature.
- * Precise centigrade temperature sensor



Code	Name	Price	+
120301	WHIP MIX Wax Pot 220vol - #10258	US\$ 488.91	+Info

Electrical Waxer Units - BUFFALO

Electrical Waxer Units - BUFFALO

Code	Name	Price	+
425323	Dental Waxing Unit Aluminum 220v	US\$ 114.59	+Info
425324	Dental Waxing Unit Aluminum 110v	US\$ 118.04	+Info

Electrical Waxer Units - DENTSPLY

Electrical Waxer Units - DENTSPLY

Code	Name	Price	+
418320	Wax Pencil Pro- 120V	US\$ 817.76	+Info
418321	Wax Pencil Pro- 220V	US\$ 852.80	+Info

Electrical Waxer Units - KEYSTONE

Electrical Waxer Units - KEYSTONE

Code	Name	Price	+
446927	Yeti Waxer NT Digital Display Shows	US\$ 710.91	+Info

ELECTRONIC WAXING SYSTEM - VANIMAN

ELECTRONIC WAXING SYSTEM - VANIMAN

Rapid heat up, Rapid cool down. Numerous Easy-on/Easy-off golden tips available. Small PKT tip comes with each unit. Digital microprocessor electronics. Temperature sensing for accurate control. Automatic selectable temperatures. Sleek modern ergonomic design. Ergonomic spatula/handpiece. Lightweight, detachable cord. Operates in ?F or ?C.



ELECTRONIC WAXING SYSTEM

Code	Name	Price	+
122653	VANIMAN Vector Electronic Waxer - Includes 1-small Pkt Tip - Electronic Waxing S...	US\$ 423.18	+Info

ACCESSORIES

Other Models

Other Models

Code	Name	Price	+
122654	Needle Point Electr Wxn Sys Ea	US\$ 20.81	+Info
122656	Lg PKT. Elctr Wxn System Ea	US\$ 32.96	+Info
122657	Hollenback Elctr Wxn Sys Ea	US\$ 34.57	+Info
122658	Beavertail Elctr Wxn Sys Ea	US\$ 32.96	+Info
122659	Sm #7 Electronic Wxn Sys Ea	US\$ 32.96	+Info
122660	Dent. Spn Tp Elctr Wxn Sys Ea	US\$ 38.38	+Info
122662	Extra Spat Assy Elctr Wxn Sys	US\$ 66.61	+Info
122663	Extra Spat Assy w/ Tp Cmbnt Ea	US\$ 66.61	+Info
122870	Cushion Grips Vector Waxer 3pk	US\$ 5.71	+Info
122871	Repl Pwr Sply f/ Vector Wxr Ea	US\$ 54.32	+Info
122872	Spat Assy Hnp w/ Cord Vctr Wxr	US\$ 36.25	+Info
122873	Spatula Assy Hnp Vector Wxr Ea	US\$ 66.61	+Info
122874	Spat Assy Hnp w/ Tp Vctr Wxr Ea	US\$ 66.61	+Info
122875	Spat Repl Cord Vector Waxer Ea	US\$ 26.27	+Info

PK Thomas

PK Thomas

Code	Name	Price	+
122661	Tip Assmt Kt -Elctr Wxn Sys Ea	US\$ 141.76	+Info

Catalog Index / Indice

Accessories for electric waxers	1
Electrical Waxer Units	6