

Abrasive Wheels unMonuted






Grinding wheels unmounted



DURA-GREEN UNMOUNTED GRINDING WHEELS - SHOFU

ABRASIVES UNMOUNTED FINISHING WHEELS - SHOFU

Use Dura-Green Wheels, Chipless Wheels (silicon carbide) for: porcelains and precious alloys
 Use Cut-Off Disks (resin-bonded) and SprueCut Wheels (silicon carbide) for: non-precious alloys and model castings

SHAPE CATEGORY	No. 7	No. 7a	No. 8	Cut-Off Disk	SpruCut No. 8
SHAPE					
ISO #	190	190	190	250	195
THICKNESS	1.7mm	2.4mm	3mm	-	3mm
DURA-GREEN					
(12 pk) PN	0150	0151	0152	-	-
CHIPLESS					
(12 pk) PN	0455	0456	0457	-	-
UNM. WHEELS					
(12 pk) PN	-	-	-	0441 (100 pk)	0160

US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.

FINISHING DISCS

Code	Name	Price	+
111521	SHOFU DURA-GREEN UnMounted Stones Wheels #8 12pk - #0152	US\$ 19.08	+Info
111519	SHOFU DURA-GREEN UnMounted Stones Wheels #7 12pk - #0150	US\$ 19.08	+Info
111520	SHOFU DURA-GREEN UnMounted Stones Wheels #7A 12pk - #0151	US\$ 19.08	+Info
111522	SHOFU Spru Cut UnMounted Stones Wheels 12pk - #0160	US\$ 18.43	+Info
111705	SHOFU Cutoff discs 100pk UnMounted Stones Finshing wheels - #0441	US\$ 26.73	+Info
111710	SHOFU UnMounted Stones Chipless Wheel #7 12pk - #0455	US\$ 18.43	+Info
111711	SHOFU UnMounted Stones Chipless Wheel #8 12pk - #0457	US\$ 18.43	+Info

Grinding wheels unmounted - DEDECO

Grinding wheels unmounted - DEDECO

Code	Name	Price	+
420199	Koolies Wheels White #1 (50)	US\$ 36.01	+Info

UNMOUNTED STONE WHEELS - BESQUAL

UNMOUNTED STONE WHEELS - BESQUAL

Cut Off Disk, Pin Cutters & Roughing Disks



Code	Name	Price	+
116091	BESQUAL #3 PINK Unmounted Stone Wheel P303 - 100pk -19mm x 4mm Piedra Verde Rued...	US\$ 15.98	+Info
116092	BESQUAL #2 PINK Unmounted Stone Wheel P302 - 100pk -19mm x 2mm Piedra Verde Rued...	US\$ 15.98	+Info

501-0001	BESQUAL #2 GREEN Unmounted Stone Wheel G302 - 100pk -19mm x 3mm Piedra Verde Rue...	US\$ 15.98	+Info
501-0002	BESQUAL #3 GREEN Unmounted Stone Wheel G303 - 100pk -19mm x 4mm Piedra Verde Rue...	US\$ 15.98	+Info
501-0003	BESQUAL #2 WHITE Unmounted Stone Wheel W302 - 100pk -19mm x 3mm Piedra Blanca Ru...	US\$ 15.98	+Info
501-0004	BESQUAL #3 WHITE Unmounted Stone Wheel W303 - 100pk -19mm x 4mm Piedra Blanca Ru...	US\$ 15.98	+Info

UNMOUNTED STONE WHEELS - KEYSTONE

UNMOUNTED STONE WHEELS - KEYSTONE

These Red Engine Wheels are made of aluminum oxide grit that performs a superior job on chrome or gold. The wheels come in square edge and knife edge in 7/8" (22.2 mm) diameter or 5/8" (15.8 mm) diameters.



Description:

- * Made of aluminum oxide
- * Available in square edge and knife edge
- * Available in 7/8" (22.2 mm) and 5/8" (15.8 mm) diameters
- * Available in 12/Pk

Code	Name	Price	+
034-1900380	KEYSTONE RED Engine Stone Wheels - 0.88in x 0.12in #1900380 -	US\$ 12.25	+Info
034-1900420	KEYSTONE R - 128 RED Engine Stone Wheels 0.62in - Knife Edge - 12pk #1900420 -	US\$ 12.25	+Info
034-1900430	KEYSTONE R - 129 RED Engine Stone Wheels - 0.88in - Knife Edge - 12pk #1900430 -	US\$ 12.25	+Info

Abrasive Wheels unMonuted

Heatless Wheels

Heatless Grinding Wheels - BESQUAL

#1
#1

Code	Name	Price	+
471485	Heatless Wheel HG1 25x4.8mm	US\$ 17.52	+Info
474238	Heatless Wheel HG1 25x4.8mm	US\$ 17.52	+Info

#2
#2

Code	Name	Price	+
471486	Heatless Wheel HG2 25x3.1mm	US\$ 17.52	+Info
474239	Heatless Wheel HG2 25x3.1mm	US\$ 15.42	+Info

#3 #3

Code	Name	Price	+
471487	Heatless Wheel HG3 25x2.4mm	US\$ 17.52	+Info
474240	Heatless Wheel HG3 25x2.4mm	US\$ 17.52	+Info

#4 #4

Code	Name	Price	+
471488	Heatless Wheel HG4 22x4.8mm	US\$ 17.52	+Info
474241	Heatless Wheel HG4 22x4.8mm	US\$ 17.52	+Info

#5 #5

Code	Name	Price	+
471489	Heatless Wheel HG5 22x3.1mm	US\$ 17.52	+Info
474242	Heatless Wheel HG5 22x3.1mm	US\$ 17.52	+Info

#6 #6

Code	Name	Price	+
471490	Heatless Wheel HG6 2x2.4mm	US\$ 17.52	+Info
474237	Heatless Wheel HG1 25x2.4mm	US\$ 16.29	+Info

Heatless Wheels - JSP

Heatless Wheels - JSP

Code	Name	Price	+
461590	Jsp - Heat Free (Heatless)	US\$ 18.11	+Info

MIZZY HEATLESS WHEEL REGULAR - KEYSTONE

KEYSTONE - MIZZY HEATLESS WHEELS REGULAR GRIT

Features:

Heatless original wheels (grey) are made of a silicone carbide abrasive to make a unique, non-loading binder for rapid reduction through the final shaping of non-precious, semi-precious, precious alloys, porcelain, acrylic, denture teeth and natural dentition, without generating excessive heat. These wheels come in 15 different sizes and in regular grit. They also withstand a maximum of 20,000 rpm.

The white heatless wheels are non-contaminating, aluminum oxide wheels that can be used on either porcelain or precious and semi-precious metal. They cut fast, run true, wear evenly, and won't ever clog. They also come in 15 different sizes.

- * Made of silicone carbide
- * White heatless wheels made of aluminum oxide
- * Available in pack for 50 or 5 pieces



50pk

#1
#1

Code	Name	Price	+
400121	Mizzy Heatless Regular Grit #1 50pk	US\$ 23.70	+Info

#10

#10

Code	Name	Price	+
400130	Mizzy Heatless Regular Grit #10 50pk	US\$ 23.70	+Info

#11

#11

Code	Name	Price	+
400131	Mizzy Heatless Regular Grit #11 50pk	US\$ 23.70	+Info

#12

#12

Code	Name	Price	+
400132	Mizzy Heatless Regular Grit #12 50pk	US\$ 23.70	+Info

#13

#13

Code	Name	Price	+
400133	Mizzy Heatless Regular Grit #13 50pk	US\$ 23.70	+Info

#14

#14

Code	Name	Price	+
400134	Mizzy Heatless Regular Grit #14 50pk	US\$ 23.70	+Info

#15

#15

Code	Name	Price	+
400135	Mizzy Heatless Reg Grit #15 50pk	US\$ 23.70	+Info

#2

#2

Code	Name	Price	+
400122	Mizzy Heatless Regular Grit #2 50pk	US\$ 23.70	+Info

#3

#3

Code	Name	Price	+
400123	Mizzy Heatless Regular Grit #3 50pk	US\$ 23.70	+Info

#4

#4

Code	Name	Price	+
400124	Mizzy Heatless Regular Grit #4 50pk	US\$ 23.70	+Info

#5
#5

Code	Name	Price	+
400125	Mizzy Heatless Regular Grit #5 50pk	US\$ 23.23	+Info

#6
#6

Code	Name	Price	+
400126	Mizzy Heatless Regular Grit #6 50pk	US\$ 23.70	+Info

#7
#7

Code	Name	Price	+
400127	Mizzy Heatless Regular Grit #7 50pk	US\$ 23.70	+Info

#8
#8

Code	Name	Price	+
400128	Mizzy Heatless Regular Grit #8 50pk	US\$ 22.51	+Info

#9
#9

Code	Name	Price	+
400129	Mizzy Heatless Regular Grit #9 50pk	US\$ 23.70	+Info

5pk

#1
#1

Code	Name	Price	+
400150	Mizzy Heatless Regular Grit #1 5pk	US\$ 4.22	+Info

#10
#10

Code	Name	Price	+
400159	Mizzy Heatless Regular Grit #10 5pk	US\$ 4.22	+Info

#2
#2

Code	Name	Price	+
400151	Mizzy Heatless Regular Grit #2 5pk	US\$ 4.22	+Info

#3
#3

Code	Name	Price	+
------	------	-------	---

400152	Mizzy Heatless Regular Grit #3 5pk	US\$ 4.22	+Info
#4			
#4			
Code	Name	Price	+
400153	Mizzy Heatless Regular Grit #4 5pk	US\$ 4.22	+Info
#5			
#5			
Code	Name	Price	+
400154	Mizzy Heatless Regular Grit #5 5pk	US\$ 3.84	+Info
#6			
#6			
Code	Name	Price	+
400155	Mizzy Heatless Regular Grit #6 5pk	US\$ 4.22	+Info
#7			
#7			
Code	Name	Price	+
400156	Mizzy Heatless Regular Grit #7 5pk	US\$ 4.22	+Info
#8			
#8			
Code	Name	Price	+
400157	Mizzy Heatless Regular Grit #8 5pk	US\$ 3.29	+Info
#9			
#9			
Code	Name	Price	+
400158	Mizzy Heatless Regular Grit #9 5pk	US\$ 4.22	+Info

MIZZY HEATLESS WHEELS WHITE - KEYSTONE

* The White Mizzy Heatless Wheels are non-contaminating, aluminum oxide wheels that can be used on either porcelain or precious and semi-precious metal.

* The White Mizzy Heatless Wheels cut fast, run true, wear evenly, and won't ever clog.

* They also come in 15 different sizes.

* Available in pack of 50 or 5 pieces



50pk

#10
#10

Code	Name	Price	+
034-6151000	#10 Wht Htless 0.6inx4.8mm 50p	US\$ 23.70	+Info

#11
#11

Code	Name	Price	+
034-6151100	#11 Wht Htless 0.6inx3.1mm 50p	US\$ 23.70	+Info

#12
#12

Code	Name	Price	+
034-6151200	#12 Wht Htless .62inx2.4mm 50p	US\$ 23.70	+Info

#13
#13

Code	Name	Price	+
034-6151300	#13 Wht Htless 0.5inx4.8mm 50p	US\$ 23.70	+Info

#14
#14

Code	Name	Price	+
034-6151400	#14 Wht Htless 0.5inx3.1mm 50p	US\$ 23.70	+Info

#15
#15

Code	Name	Price	+
034-6151500	#15 Wht Htless 0.5inx2.4mm 50p	US\$ 25.88	+Info

#2
#2

Code	Name	Price	+
034-6150200	#2 Wht Htless 1inx3.1mm 50p	US\$ 23.70	+Info

#3
#3

Code	Name	Price	+
034-6150300	#3 Wht Htless 1inx2.4mm 50p	US\$ 23.70	+Info

#4
#4

Code	Name	Price	+
034-6150400	#4 Wht Htless 0.88inx4.8mm 50p	US\$ 23.70	+Info

#5
#5

Code	Name	Price	+
034-6150500	#5 Wht Htless 0.88inx3.1mm 50p	US\$ 23.70	+Info

#6

#6

Code	Name	Price	+
034-6150600	#6 Wht Htless 0.88inx2.4mm 50p	US\$ 23.70	+Info

#7

#7

Code	Name	Price	+
034-6150700	#7 Wht Htless 0.75inx4.8mm 50p	US\$ 23.70	+Info

#8

#8

Code	Name	Price	+
034-6150800	#8 Wht Htless 0.75inx3.1mm 50p	US\$ 22.51	+Info

#9

#9

Code	Name	Price	+
034-6150900	#9 Wht Htless 0.75inx2.4mm 50p	US\$ 23.70	+Info

5pk

#10

#10

Code	Name	Price	+
400144	Mizzy Heatles White #10 5pk	US\$ 4.22	+Info

#11

#11

Code	Name	Price	+
400145	Mizzy Heatles White #11 5pk	US\$ 4.22	+Info

#12

#12

Code	Name	Price	+
400146	Mizzy Heatles White #12 5pk	US\$ 4.22	+Info

#13

#13

Code	Name	Price	+
400147	Mizzy Heatles White #13 5pk	US\$ 4.22	+Info

#14
#14

Code	Name	Price	+
400148	Mizzy Heatles White #14 5pk	US\$ 4.22	+Info

#15
#15

Code	Name	Price	+
400149	Mizzy Heatles White #15 5pk	US\$ 4.63	+Info

#2
#2

Code	Name	Price	+
400136	Mizzy Heatles White #2 5pk	US\$ 4.22	+Info

#3
#3

Code	Name	Price	+
400137	Mizzy Heatles White #3 5pk	US\$ 4.22	+Info

#4
#4

Code	Name	Price	+
400138	Mizzy Heatles White #4 5pk	US\$ 4.22	+Info

#5
#5

Code	Name	Price	+
400139	Mizzy Heatles White #5 5pk	US\$ 4.22	+Info

#6
#6

Code	Name	Price	+
400140	Mizzy Heatles White #6 5pk	US\$ 4.22	+Info

#7
#7

Code	Name	Price	+
400141	Mizzy Heatles White #7 5pk	US\$ 4.22	+Info

#8
#8

Code	Name	Price	+
400142	Mizzy Heatles White #8 5pk	US\$ 3.25	+Info

#9
#9

Code	Name	Price	+
400143	Mizzy Heatles White #9 5pk	US\$ 4.28	+Info

RED HEATLESS WHEELS - KEYSTONE

RED HEATLESS WHEELS - KEYSTONE

The Red Heatless Wheels have a ceramic-coated abrasive grain that ensures a smooth, free, cutting wheel at the coolest temperatures. This is due to the unusual micro-chipping action of the crystals as they wear.

The wheels are made of non-contaminating pure alpha aluminum oxide hexagonal crystals with a tough ceramic coating. Great for use on all non-precious and semi-precious metals, porcelains, cercons and lava ? there?s no need to change back to a silicon carbide wheel. The wheels are long lasting and exhibits excellent uniform wear and available in the most popular Mizzy sizes

* Made of a pure alpha aluminum oxide hexagonal crystals

* Available in the most popular Mizzy sizes: 1, 2, 3, 5, 8



#1

#1

Code	Name	Price	+
101128	#1 red Htless 4.8 x 25.4 mm 50p	US\$ 23.78	+Info

#3

#3

Code	Name	Price	+
101129	#3 red Htless 2.4x25.4 50p	US\$ 23.70	+Info

#5

#5

Code	Name	Price	+
6131605	#5 red Htless 3.1x22 50p	US\$ 23.70	+Info

#8

#8

Code	Name	Price	+
034-6131608	#8 red Htless 3.1 x 19 mm 50p	US\$ 25.62	+Info

Abrasive Wheels unMonuted

Knock Down Wheels



Knock Down Wheels - BUFFALO

Knock Down Wheels - BUFFALO

Code	Name	Price	+
424670	#21 Abrasive Lathe Wheel	US\$ 12.07	+Info
424692	Blue #461 Abrasive L	US\$ 11.09	+Info

Knock Down Wheels - KEYSTONE

Knock Down Wheels - KEYSTONE



Code	Name	Price	+
034-1900080	White Knock D. Wheels 3x0.37in	US\$ 9.69	+Info
034-1900090	White Knock D. Wheels 4x0.5in	US\$ 11.64	+Info
034-1900060	Wheels Grinding Black 3x 0.5in	US\$ 10.85	+Info

Catalog Index / Indice

Grinding wheels unmounted	1
Heatless Wheels	2
Knock Down Wheels	10