

Wires Clasps & Arch wires

Ball Clasps

Ball Clasps - KEYSTONE

Ball Clasps - KEYSTONE

Code	Name	Price	+
446026	Ball Clasp SS Small .028 dia 10pk	US\$ 4.30	+Info

CLASPS - KEYSTONE

CLASPS - KEYSTONE

* Anchorage of removable appliances, frequently applied on incisors and small posterior teeth.

* The ball clasp can be bent in different manners under the condition that the ball is always inserted gingivally in the interdental area.

* Its application is normally with short clinical crowns or with teeth that don't have retentive undercuts (for example lower incisors).

* It should lie in the interdental area, gingival of the contact point. For this, it is necessary to slightly carve the papila.



Our Stainless Steel Ball Clasps come in two different sizes:

* Small - .028" dia. (.7mm)

* Large - .032" dia. (.8mm)

Comes in either: 10/Pkg or 100/Pkg.

10PK

Code	Name	Price	+
034-1270004-10	KEYSTONE 0.024in / 0.6 mm Ball Clasps 10pk - (MFG #1270004-10)	US\$ 3.60	+Info
034-1270005	KEYSTONE 0.028in / 0.7mm Stainless steel Ball Clasps 10pk - (MFG #1270005-10)	US\$ 4.08	+Info
034-1270010	KEYSTONE 0.032in / 0.8mm Ball Clasps 10pk - (MFG #1270010-10)	US\$ 4.08	+Info
034-1270021-10	KEYSTONE 0.036in / 0.9mm Ball Clasps 10pk - (MFG #1270021-10)	US\$ 6.37	+Info
034-1270022-10	KEYSTONE 0.040in / 1.0 mm Ball Clasps 10pk -	US\$ 8.51	+Info

100PK

Code	Name	Price	+
034-1270004	KEYSTONE 0.024in / 0.6 mm Ball Clasps 100pk #127004 -	US\$ 33.72	+Info
034-1270006	KEYSTONE Stainless Steel Ball Clasps Small .028in (.7mm) dia, Package of 100 -...	US\$ 38.32	+Info
034-1270020	KEYSTONE Stainless Steel Ball Clasps Large .032in (.8mm) dia, Package of 100 -...	US\$ 38.32	+Info
034-1270021	KEYSTONE Stainless Steel Ball Clasps .036in (.9mm) dia, Package of 100 - #1270021	US\$ 38.32	+Info
034-1270022	KEYSTONE Stainless Steel Ball Clasps .040in (1.0mm) dia, Package of 100 - #127...	US\$ 38.32	+Info

Wires Clasps & Arch wires

Clasp Wires



CLASP WIRE - KEYSTONE

STAINLESS STEEL CLASP WIRE - KEYSTONE

These Stainless Steel Clasp Wires are available in half-round and round configurations.



Features:

- * Stainless Steel Clasp Wire - Gauge Half Round
- * Stainless Steel Clasp Wire - Gauge Round

HALF-ROUND CLASP WIRE - KEYSTONE

Code	Name	Price	+
034-1900720	KEYSTONE Stainless Steel Orthodontic Clasp Wire, 12ga .08in2mm , Half Round, 1 ...	US\$ 6.72	+Info
034-1900730	KEYSTONE Stainless Steel Orthodontic Clasp Wire, 14ga .064in , Half Round, 1 ft...	US\$ 6.72	+Info
034-1900750	KEYSTONE Stainless Steel Orthodontic Clasp Wire, 16ga .051in , Half Round, 1 ft...	US\$ 7.48	+Info

ROUND CLASP WIRE - KEYSTONE

Code	Name	Price	+
034-1900760	KEYSTONE 17ga Orthodontic Clasp Wire - Stainless Round 10pk - 1.10mm - 110 300m...	US\$ 7.48	+Info
034-1900770	KEYSTONE Stainless Steel Orthodontic Clasp Wire, 18ga, Round, Straight 10 x 1ft...	US\$ 6.72	+Info
034-1900780	KEYSTONE Stainless Steel Orthodontic Clasp Wire, 19ga, Round, 10 x 1 ft Section...	US\$ 8.19	+Info
034-1900790	KEYSTONE 20ga Orthodontic Clasp Wire - Stainless Round 10pk x - 0.80mm - 080 30...	US\$ 7.48	+Info

Clasp Wires - JSP

Clasp Wires - JSP

Code	Name	Price	+
437650	JSP Clasp Wire - 1/2 Round 16 Ga. (US\$ 6.00	+Info
437656	JSP Clasp Wire - Round 17 Ga. (1.1m	US\$ 6.00	+Info
437664	JSP Clasp Wire - Round 20 Ga. (.81m	US\$ 6.00	+Info
437668	JSP Clasp Wire - 1/2 Round 12 Ga. (US\$ 9.55	+Info

Clasp Wires - KEYSTONE

Clasp Wires - KEYSTONE

Code	Name	Price	+
446489	Keystone Stainless Steel Clasp Wire	US\$ 6.72	+Info

Clasp Wires - TICONIUM

Clasp Wires - TICONIUM

Code	Name	Price	+
420989	18ga Clasp Wire Round 5pk	US\$ 25.00	+Info

Wires Clasps & Arch wires

Retainer Wire



RETAINER WIRE - ECCO

RETAINER WIRE - ECCO

- * Diamond drawn medium-hard temper for retainers and functional appliances
- * Easy to bend all loops and will not fracture
- * Fifty or ten 14 lights per tube
- * Also 50gm and 500gm spools



10pk

Code	Name	Price	+
106977	ECCO - 30gauge Ligature Retainer Wire .25mm x14 inch log - 10pk (Alambre de ort...	US\$ 3.50	+Info
106979	ECCO - 28ga Ligature Retainer Wire .30mm x14 inch log - 10pk (Alambre de ortod...	US\$ 4.75	+Info
106981	ECCO - 27gauge Ligature Retainer Wire .36mm x14 inch log - 10pk (Alambre de ort...	US\$ 6.00	+Info
106983	ECCO - 25gauge Ligature Retainer Wire .45mm x14 inch log - 10pk (Alambre de ort...	US\$ 6.00	+Info

50pk

Code	Name	Price	+
106970	ECCO - 30ga Ligature Retainer Wire .25mm size x14 inch log - 50pk (Alambre de ...	US\$ 16.00	+Info
106972	ECCO - 28gauge Ligature Retainer Wire .30mm size x14 inch log - 50pk (Alambre d...	US\$ 20.00	+Info
106974	ECCO - 27ga Ligature Retainer Wire .36mm size x14 inch log - 50pk (Alambre de ...	US\$ 19.80	+Info
106976	ECCO - 25ga Ligature Retainer Wire .45mm size x14 inch log - 50pk (Alambre de ...	US\$ 21.60	+Info

Wires Clasps & Arch wires

ROACH CLASPS




ROACH CLASPS - KEYSTONE

STAINLESS STEEL WIRE CLASPS

This Stainless Steel Wire Clasp is bicuspid medium. Items come in straight, left, and right.

Available in package of 10

Size: Medium



Code	Name	Price	+
034-1271210	KEYSTONE Stainless Steel Wire Roach Clasp, Bicuspid Medium, Right, pack of 10 -...	US\$ 22.00	+Info
034-1271160	KEYSTONE Stainless Steel Wire Roach Clasp, Anterior Small, Straight, pack of 10...	US\$ 12.38	+Info
034-1271190	KEYSTONE Stainless Steel Wire Roach Clasp, Bicuspid Medium, Straight, pack of 1...	US\$ 22.00	+Info
446911	KEYSTONE Stainless Steel Wire Roach Clasp, Molar Large, Right, pack of 100 - #1...	US\$ 139.71	+Info
446912	KEYSTONE Stainless Steel Wire Roach Clasp, Bicuspid Medium, Straight, 100 pack ...	US\$ 181.17	+Info
446913	KEYSTONE Stainless Steel Wire Roach Clasp, Molar Large, Left, pack of 100 - #12...	US\$ 139.71	+Info
446914	KEYSTONE Stainless Steel Wire Roach Clasp, Molar Large, Straight, pack of 100 -...	US\$ 139.71	+Info
034-1271200	KEYSTONE Stainless Steel Wire Roach Clasp, Bicuspid Medium, Left, pack of 10 - ...	US\$ 22.00	+Info

Wires Clasps & Arch wires

Spooled Ligat & Retainer Wire



Spooled Ligat & Retainer Wire - ORTHOTECHNOLOGY

Spooled Ligat & Retainer Wire - ORTHOTECHNOLOGY

Code	Name	Price	+
421034	Spool Twist Wire 0175 30 per spool	US\$ 44.52	+Info

SPOOLED LIGATURE WIRE - BESQUAL

SPOOLED LIGATURE WIRE - BESQUAL

High quality stainless steel ligature wire is super smooth for easy tying. Soft stainless steel with a bright, diamond die high quality finish. Spools for the office that prefers to form their own ligatures.

Available in 0.8in; 0.9in and 1.0in

Weight: 50gr



Code	Name	Price	+
123011	BESQUAL Spooled Ligature Wire .008in Diameter - 20GA 0.8mm (0.032in) x 11m - (W...	US\$ 12.51	+Info
123012	BESQUAL Spooled Ligature Wire .008in Diameter - 19GA 0.9mm (0.036in) x 9m - (We...	US\$ 12.51	+Info
123013	BESQUAL Retainer Wire - 18GA. 1 spool, 1.0 mm (0.040 in) x 7 m, stainless - #55...	US\$ 12.51	+Info

SPOOLED LIGATURE WIRE - ECCO

SPOOLED LIGATURE WIRE - ECCO

Bright Finish with Soft Temper for Easy Use

Weight: 1pound/455grams



Code	Name	Price	+
------	------	-------	---

113556	ECCO - 32ga Spooled Ligature Wire .008in Diameter T304/304L #6 - (0.20mm Diame...	US\$ 35.04	+Info
113557	ECCO - 31gauge Spooled Ligature Wire .009in Diameter T304/304L # 7 - (0.23mm Di...	US\$ 47.00	+Info
113558	ECCO - 30gauge Spooled Ligature Wire .010in Diameter T304/304L # 8- (0.25mm Dia...	US\$ 47.00	+Info
113559	ECCO - 29gauge Spooled Ligature Wire .011in Diameter - T304/304L # 9 (0.28 mm D...	US\$ 47.00	+Info
113560	ECCO - 28gauge Spooled Ligature Wire .012in Diameter T304/304L # 10 - (0.31mm D...	US\$ 47.00	+Info
113561	ECCO - 26gauge Spooled Ligature Wire .014in Diameter T304/304L # 11 - (0.36mm D...	US\$ 47.00	+Info

Wires Clasps & Arch wires

Spring Tensile Wire



SPRING TENSILE ROUND WIRE - ECCO

SPRING TENSILE ROUND WIRE - ECCO



Code	Name	Price	+
113562	JSP .012in 28ga Round Orthodontic Wire, 14in Long Spring Tensile 10pk - #OD11...	US\$ 7.00	+Info
113563	JSP .014in 27ga Round Orthodontic Wire, 14in Long Spring Tensile 10pk - #OD114	US\$ 7.00	+Info
113564	JSP .016in 26ga Round Orthodontic Wire, 14in Long Spring Tensile 10pk # OD116	US\$ 7.00	+Info
113565	JSP .018in 25ga Round Orthodontic Wire, 14in Long Spring Tensile 10pk - #OD118...	US\$ 6.09	+Info
113566	JSP .020in Round Orthodontic Wire, 14in Long Spring Tensile 10/Pack - #OD120 #...	US\$ 7.00	+Info
113567	JSP .022in Round Orthodontic Wire, 14in Long Spring Tensile 10/Pack - #OD122 #...	US\$ 7.00	+Info
113568	JSP .024in Round Orthodontic Wire, 14in Long Spring Tensile 10/Pack - #OD124 #...	US\$ 7.00	+Info
113569	ECCO - Spring Tensile Round Wire .025Diam T304 SPR 14 - 10pk - Round Orthodontic...	US\$ 7.00	+Info

Catalog Index / Indice

Ball Clasps	1
Clasp Wires	2
Retainer Wire	3
ROACH CLASPS	3
Spoiled Ligat & Retainer Wire	4
Spring Tensile Wire	5