


CAD/CAM abrasives

CAD/CAM Burs

AMANN GIRRBACH DIAMOND CAD/CAM MILLING BURS - ATTRI TOR

AMANN GIRRBACH DIAMOND CAD/CAM
MILLING BURS - ATTRITOR

Available in 0.4mm - 1mm - 1.4mm - 1.8mm



Code	Name	Price	+
124124	ATTRITOR - Amann Girrbach Roto Diamond Cad Cam Bur 0.4mm for disilicate and feld...	US\$ 25.29	+Info
124125	ATTRITOR - Amann Girrbach Roto Diamond Cad Cam Bur 1mm for disilicate and feldsp...	US\$ 25.29	+Info
124126	ATTRITOR - Amann Girrbach Roto Diamond Cad Cam Bur 1.4mm for disilicate and feld...	US\$ 34.18	+Info
124127	ATTRITOR - Amann Girrbach Roto Diamond Cad Cam Bur 1.8mm for disilicate and feld...	US\$ 25.29	+Info

AMANN GIRRBACH DLC CAD/CAM MILLING BURS - ATTRITOR

AMANN GIRRBACH DLC CAD/CAM MILLING BURS - ATTRITOR


Features

- * Ultra-precisi
- * Hard, wear-resistant material
- * Last longer
- * Milling Bur 2.5mm - 1.0mm - 0.6mm

Can be applied to milling the Zirconia disc, PMMA disc, Acrylic disc for CAD/CAM

#UsddAmann

AMANN GIRRBACH is a trademark of Amann Girrbach AG




Code	Name	Price	+
115577	ATTRITOR - AMANN GIRRBACH Cad/Cam DLC Coated Milling Bur 0.6mm - each - for Cera...	US\$ 18.00	+Info
416754	ATTRITOR - AMANN GIRRBACH Cad/Cam Coated Milling Bur 0.3mm - each - for Ceramill...	US\$ 17.82	+Info
115571	ATTRITOR- AMANN GIRRBACH Cad/Cam DLC Coated Milling Bur 1mm - each - for Ceramil...	US\$ 18.00	+Info
115565	ATTRITOR- AMANN GIRRBACH Cad/Cam DLC Coated Milling Bur 2.5mm - each - for Ceram...	US\$ 18.00	+Info

CEREC COMPATIBLE COMPACT MILLING BURS - ATTRITOR

CEREC COMPACT COMPATIBLE MILLING BURS- ATTRITOR

The Cerec Inlab Sirona Milling Instruments Are The Machine Parts, Known As Burs, That Do The Actual Milling Of The Composite Blocks In Their Respective Unit.

For Cerec Inlab Compact Cad/Cam Milling Diamond



Code	Name	Price	+
113474	CEREC compatble Compact Cylinder Pointed Bur - 6pk - #5855734 SGH for Cerec Inla...	US\$ 135.73	+Info
113475	CEREC compatble Compact Step Bur #12mm - 6 pk -5855734 for Cerec Inlab Compact C...	US\$ 135.73	+Info

CEREC MCXL COMPATIBLE MILLING BURS - ATTRITOR


CEREC MCXL COMPATIBLE MILING BURS - ATTRITOR

Details:

- * Full Compatible Sinter Diamond Burs For Cerec Sirona
- * Available For Ceramic And Zirconia

CEREC is a registered trademark of Sirona Dental Systems.

#UsddCerec



Cylinder Burs

Code	Name	Price	+
112250	CEREC MCXL Cylinder Pointed Bur #CMP 12S for MCXL - 1pk - #6240159 Cad/Cam Milli...	US\$ 22.63	+Info
116130	CEREC MCXL compatble Cylinder Pointed Step Bur #20 - 1pk - for MCXL #6259589 Cad...	US\$ 22.63	+Info

Step Burs

Code	Name	Price	+
116222	* CEREC compatble MCXL Step Bur #12 - 1pk - for MCXL # 626025 Cad/Cam Milling...	US\$ 19.55	+Info
112251	CEREC MCXL compatble Step Bur #12S - 1pk - for MCXL #6240167 Cad/Cam Milling Dia...	US\$ 22.63	+Info
116131	CEREC compatble MCXL Step Bur #20 - 1pk - for MCXL #6259597 Cad/Cam Milling Diam...	US\$ 22.63	+Info


CEREC X5 COMPATIBLE MILLING BURS - ATTRIT

OR

CEREC X5 COMPATIBLE MILLING BURS - ATTRITOR

Dental milling bur compatible with Sirona CEREC inLab MC X5 for milling zirconia.

Special coating to ensure long and accurate milling.



Code	Name	Price	+
123456	CEREC - MC X5 Bur 0.5mm BLUE Ring Composite # SIC05MM Compatible# 6478171	US\$ 42.00	+Info
123455	CEREC - MC X5 Bur 1.0mm BLUE Ring Composite # SIC10MM Compatible# 6478163	US\$ 42.00	+Info
123454	CEREC - MC X5 Bur 2.5mm BLUE Ring Composite # SIC25MM Compatible# 6478155	US\$ 43.50	+Info
123453	CEREC X5 Bur 0.5mm RED Ring PMMA - WAX # SIP05MM Compatible# 6478114	US\$ 42.00	+Info
123452	CEREC X5 Bur 1.0mm RED PMMA-WAX # SIP10MM Compatible# 6478106	US\$ 51.17	+Info
123451	CEREC - MC X5 Bur 2.5mm RED Ring PMMA- WAX # SIP25MM Compatible# 6478098	US\$ 55.01	+Info
123450	CEREC - MC X5 Bur 0.5mm YELLOW Ring ZrO2 (Zirconio) # SI05MM Compatible# 6478056	US\$ 42.00	+Info
123449	CEREC - MC X5 Bur 1.0mm YELLOW Ring ZrO2 (Zirconio) # SI10MM Compatible# 6478049	US\$ 42.00	+Info
123448	CEREC - MC X5 Bur 2.5mm YELLOW Ring ZrO2 (Zirconio) # SI25MM Compatible# 6478031	US\$ 42.00	+Info

DIAMOND COATED (DC) CAD/CAM MILLING BURS - ATTRITOR

DIAMOND COATED (DC) CAD/CAM MILLING BURS - ATTRITOR

Features:

Compatible with AMANN GIRRBACH - VHF - WIELAND - ROLAND and ZIRKONZAHN



AMANN GIRRBACH

Code	Name	Price	+
124070	ATTRITOR - AMANN GIRRBACH Cad/Cam Bur DC round 0.6mm head - 3mm Shank - 47 mm le...	US\$ 32.00	+Info
124071	ATTRITOR - AMANN GIRRBACH Cad/Cam Bur DC round 1mm head - 3mm Shank - 47 mm leng...	US\$ 32.00	+Info
124072	ATTRITOR - AMANN GIRRBACH Cad/Cam Bur DC round 2.5mm head - 3mm Shank - 47 mm le...	US\$ 32.00	+Info

ROLAND

Code	Name	Price	+
124076	ATTRITOR - ROLAND DC Cad/Cam Bur DIAMOND COATED round 0.6mm head - 4mm Shank - ...	US\$ 32.00	+Info
124077	ATTRITOR - ROLAND DC Cad/Cam Bur DIAMOND COATED round 1mm head - 4mm Shank - 50...	US\$ 32.00	+Info
124078	ATTRITOR - ROLAND DC Cad/Cam Bur DIAMOND COATED round 2mm head - 4mm Shank - 50...	US\$ 32.00	+Info

VHF

Code	Name	Price	+
124063	ATTRITOR - VHF & YETI Cad/Cam Bur DC DIAMOND COATED 0.6mm head - 3mm Shank - 35 ...	US\$ 49.21	+Info
124064	ATTRITOR - VHF & YETI Cad/Cam Bur DC DIAMOND COATED 1mm head - 3mm Shank - 35 mm...	US\$ 49.21	+Info
124065	ATTRITOR - VHF & YETI Cad/Cam Bur DC DIAMOND COATED 2mm head - 3mm Shank - 35 mm...	US\$ 49.21	+Info

WIELAND

Code	Name	Price	+
124067	ATTRITOR - WIELAND Cad/Cam Bur DIAMOND COATED 0.7mm head - 3mm Shank - 35 mm len...	US\$ 61.85	+Info
124068	ATTRITOR - WIELAND Cad/Cam Bur DIAMOND COATED 1mm head - 3mm Shank - 35 mm lenh...	US\$ 61.85	+Info
124069	ATTRITOR - WIELAND Cad/Cam Bur DIAMOND COATED 2-5mm head - 3mm Shank - 35 mm len...	US\$ 61.85	+Info

ZIRKONZAHN

Code	Name	Price	+
124073	ATTRITOR - ZIRKONZAHN Cad/Cam Bur DIAMOND COATED 0.6mm head - 3mm Shank - 47 mm ...	US\$ 19.20	+Info
124074	ATTRITOR - ZIRKONZAHN Cad/Cam Bur DIAMOND COATED 1mm head - 3mm Shank - 47 mm le...	US\$ 32.00	+Info
124075	ATTRITOR - ZIRKONZAHN Cad/Cam Bur DIAMOND COATED 2mm head - 3mm Shank - 47 mm le...	US\$ 32.00	+Info

ROLAND COMP 52D CAD/CAM MILLING BURS - ATTRITOR

ROLAND COMP 52D CAD/CAM MILLING BURS - ATTRITOR

Dental Milling Burs for its advanced D50 - D51 zirconia milling machine.
Each tool is specially engineered for precise performance, maximum durability and longer life.

Features:

- *Compatible with Roland DWX-52D, DWX-52DC, DWX-52DCi
- * Available in 4 Different Sizes
- * Carbide Milling Burs
- * Compatible with roland existing Carbide Tools

Rolands milling burs are available in four different sizes allowing dental professionals to choose an ideal tool for every milling strategy. Because they are fully compatible with Rolands existing carbide tools, users can switch to the new tools without any changes in the CAM software, making the upgrade a true "plug-n-play" solution.



SET

Code	Name	Price	+
------	------	-------	---

420987	ATTRITOR - ROLAND DLC Cad/Cam Bur KIT - 13pcs - 50 mm length - 2 blades #MBD-R1016	US\$ 188.50	+Info
--------	--	-------------	-----------------------

D50 - D51

Code	Name	Price	+
115578	ATTRITOR - ROLAND Cad/Cam DLC Milling Burs 0.8mm ea. for DWX-52D-D51-D50-52DCI	US\$ 16.27	+Info
115583	ATTRITOR - ROLAND Cad/Cam DLC Milling Burs 0.6mm ea. for DWX-52D-D51-D50-52DCI	US\$ 18.00	+Info
115572	ATTRITOR - ROLAND Cad/Cam DLC Milling Burs 1mm ea. - MBR1016 - for	US\$ 18.00	+Info
115566	ATTRITOR - ROLAND Cad/Cam DLC Milling Burs 2mm ea. - MBR2016 - for DWX-52D-D51-D...	US\$ 18.00	+Info

ROLAND COMP DIAMOND CAD/CAM MILLING BURS - ATTRITOR

ROLAND COMP DIAMOND CAD/CAM MILLING BURS - ATTRITOR

Features:

A standard milling tool for finishing and adding detail to dental prosthetics, is designed for use on the Roland DWX-4W and DWX-42W wet mill. It is approved for grinding glass ceramics and composite resins.

* Roland Compatible

Available in 1mm - 2.5mm - 0.6mm



Code	Name	Price	+
124122	ATTRITOR - Roland Disilicate Diamond Cad Cam Bur 0.6mm for disilicate and feldsp...	US\$ 34.18	+Info
124121	ATTRITOR - Roland Disilicate Diamond Cad Cam Bur 1mm for disilicate and feldspar...	US\$ 34.18	+Info
124123	ATTRITOR - Roland Disilicate Diamond Cad Cam Bur 2.5mm for disilicate and feldsp...	US\$ 34.18	+Info

VHF & YETI DLC CAD/CAM MILLING BURS - ATTRITOR

VHF & YETI DLC CAD/CAM COATED MILLING BUR - ATTRITOR

Features:

- * Ultra-precision.
- * Hard, wear-resistant material.
- * Last longer.
- * Can be applied to milling the Zirconia disc, PMMA disc, Acrylic disc for CAD/CAM

Quantity: 1 pc

Sizes: 0.6 mm - 1mm -2mm

Made in Germany



Set

Code	Name	Price	+
416999	VHF Cad/Cam DLC Coated Milling Bur - 40mm long 15 pcs set = 7 x 2mm - 5 1mm & 3 ...	US\$ 679.67	+Info

35mm 4 Axes

Code	Name	Price	+
115582	VHF & YETI Cad/Cam DLC Coated Milling Bur 0.6mm - 35mm long Each - Made in Germa...	US\$ 41.88	+Info
115575	VHF & YETI Cad/Cam DLC Coated Milling Bur 1mm - 35mm long Each - Made in Germany...	US\$ 41.88	+Info
115569	VHF & YETI Cad/Cam DLC Coated Milling Bur 2mm - 35mm long Each - Made in Germany...	US\$ 41.88	+Info

40mm 5 Axes

Code	Name	Price	+
124307	VHF & YETI Cad/Cam DLC Coated Milling Bur 0.6mm - 40mm long Each - Made in Germa...	US\$ 41.04	+Info
124306	VHF & YETI Cad/Cam DLC Coated Milling Bur 1mm - 40mm long Each - Made in Germany...	US\$ 40.62	+Info
124305	VHF & YETI Cad/Cam DLC Coated Milling Bur 2.5mm - 40mm long Each - Made in Germa...	US\$ 34.76	+Info

WIELAND DLC CAD/CAM MILLING BURS - ATTRIT

OR

ATTRITOR - WIELAND DLC CAD/CAM MILLING BURS


Length 35mm & 40mm

WIELAND MINI Cad Cam Coated Bur. The Wieland Mini Cad / Cam Burs are sold individually and for use with the Ceramill Motion - Made in Germany

Features:

- * Super Accuracy.
- * Hard Ware-Resistant Material.
- * Longer Life.
- * Can Be Applied To Milling Zirconia Disc, Wax Block Disc, Acrylic Disc CAD / CAM Use.

Quantity: 1 Pc
Made In Germany



35MM LENGH

Code	Name	Price	+
115579	WIELAND MINI Cad/Cam Coated Milling Bur 0.6mm x 35mm lengh - each - MBV060 - Mad...	US\$ 11.04	+Info
115573	WIELAND MINI Cad/Cam Coated Milling Bur 1mm x 35mm lengh - each - for Ceramill M...	US\$ 15.36	+Info
115567	WIELAND MINI Cad/Cam Coated Milling Bur 2.5mm x 35mm lengh - each - for Ceramill...	US\$ 15.36	+Info

40MM LENGH

Code	Name	Price	+
420161	WIELAND MINI Cad/Cam Coated Milling Bur 0.6mm x 40mm lengh - each - Made in Germany	US\$ 16.00	+Info
420162	WIELAND MINI Cad/Cam Coated Milling Bur 1mm x 40mm lengh - each - for Ceramill M...	US\$ 16.00	+Info
420163	WIELAND MINI Cad/Cam Coated Milling Bur 2.5mm x 40mm lengh - each - for Ceramill...	US\$ 16.00	+Info

WIELAND DLC CAD/CAM MILLING BURS - ATTRIT OR

WIELAND DLC CAD/CAM MILLING BURS - ATTRITOR


Length 40mm

WIELAND Cad Cam Coated Bur.

Features:

- * Milling Bur 2.0mm - 1.0mm - 0.6mm
- * Super Accuracy.
- * Hard ware-resistant material.
- * Longer life.
- * Can be applied to milling zirconia disc, wax block disc, acrylic disc CAD / CAM use.
- * Quantity: 1 pc

Made In Germany



Code	Name	Price	+
115580	WIELAND Cad/Cam Coated Milling Bur 0.6mm x 40mm length - each -for Wieland & - Ma...	US\$ 11.04	+Info
115574	WIELAND Cad/Cam Coated Milling Bur 1mm x 40mm length - each - for Wieland & - Mad...	US\$ 15.84	+Info
115568	WIELAND Cad/Cam Coated Milling Bur 2.0mm x 40mm length - each - for Wieland & - M...	US\$ 16.00	+Info

YETI & VHF GLASS CERAMIC MI LLING BURS 4 & 5 Axes - ATTRITO R



Code	Name	Price	+
417087	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 2.6mm FLAT END for di...	US\$ 22.97	+Info
417088	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 2.4mm BALL END for di...	US\$ 38.28	+Info
417081	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 1mm BALL END for disi...	US\$ 38.28	+Info
417086	ATTRITOR - VHF & YETI Disilicate Sinter Diamond CadCam Bur 0.6mm BALL END for di...	US\$ 22.97	+Info

ZIRKONZAHN DLC M5 CAD/CAM MILLING BURS - ATTRITOR

ZIRKONZAHN DLC M5 CAD/CAM MILLING BURS - ATTRITOR

M5 - TEC Milling - Bur 0.6mm

This M5 Zirkonzahn carbide burs are round end milling burs, and are designed to mill pre-sintered Zirconia with CAD/CAM Milling Unit M5 (Precise Milling)

Features:

- * Round End
- * Carbide Milling Bur for M5
- * Available in 3 Different Sizes 0.6mm - 1mm - 2mm

ZIRKONZAHN is a trademark of ZIRKONZAHN GMBH



Code	Name	Price	+
115581	ZirkonZahn M5 Cad/Cam DLC Milling Bur 0.5mm - each - Made in Germany #951-310506	US\$ 37.39	+Info
115576	ZirkonZahn M5 Cad/Cam DLC Milling Bur 1mm - each - Made in Germany #951-311012	US\$ 37.39	+Info
115570	ZirkonZahn M5 Cad/Cam DLC Milling Bur 2mm - each - Made in Germany #951-312018	US\$ 37.39	+Info

AMANN GIRRBACH DLC CAD/CAM MILLING BURS - AMANN GIRRBACH

DLC CAD/CAM MILLING BURS - AMANN GIRRBACH

Ceramill Roto Motion Burs with Magnet Holder

Features:

- * Ultra-precision.
- * Hard, wear-resistant material.
- * Last longer.
- * Can be applied to milling the Zirconia disc, PMMA disc, Acrylic disc for CAD/CAM

Available in Size

- * 0.9mm
- * 1.0mm
- * 1.2mm
- * 2.0mm
- * 2.5mm

#UsddAmann

AMANN GIRRBACH is a trademark of Amann Girrbach AG



Code	Name	Price	+
114391	AMANN GIRRBACH Ceramill Roto Motion 1.0mm bur with magnetic holder 1pk #760669	US\$ 96.22	+Info
114394	AMANN GIRRBACH Ceramill Roto Motion 0.6mm bur with magnetic holder 1pk #760606	US\$ 100.05	+Info
114395	AMANN GIRRBACH Ceramill Roto Motion 0.3mm bur with magnetic holder 1pk #760607	US\$ 100.05	+Info

114392	AMANN GIRRBACH Ceramill Roto Motion 2.5mm bur with magnetic holder 1pk #760668	US\$ 96.22	+Info
468252	AMANN GIRRBACH - Roto RFID 1,0 SC - # 760681	US\$ 101.87	+Info
468268	AMANN GIRRBACH - Roto RFID 2,0 Model - # 760698	US\$ 59.96	+Info
468253	AMANN GIRRBACH - Roto RFID 0,6 SC - # 760682	US\$ 101.87	+Info
468254	AMANN GIRRBACH - Roto RFID 0,3 SC - # 760683	US\$ 106.03	+Info
468255	AMANN GIRRBACH - Roto RFID 2,5 DMB DC - # 760684	US\$ 198.59	+Info
468256	AMANN GIRRBACH - Roto RFID 1,0 DMB DC - # 760685	US\$ 198.59	+Info
468257	AMANN GIRRBACH - Roto RFID 0,6 DMB DC - # 760686	US\$ 198.59	+Info
468258	AMANN GIRRBACH - Roto RFID 0,3 DMB DC - # 760687	US\$ 180.93	+Info
468259	AMANN GIRRBACH - Roto RFID 2,5 CoCr - # 760688	US\$ 68.59	+Info
468260	AMANN GIRRBACH - Roto RFID 1,5 CoCr - # 760689	US\$ 68.59	+Info
468261	AMANN GIRRBACH - Roto RFID 1,0 CoCr - # 760690	US\$ 68.59	+Info
468262	AMANN GIRRBACH - Roto RFID 2,5 PMMA Denture - # 760692	US\$ 184.03	+Info
468263	AMANN GIRRBACH - Roto RFID 1,5 PMMA Denture - # 760693	US\$ 184.03	+Info
468264	AMANN GIRRBACH - Roto RFID 1,2 Drill - # 760694	US\$ 52.45	+Info
468265	AMANN GIRRBACH - Roto RFID 1,0 Wax Denture - # 760695	US\$ 68.59	+Info
468266	AMANN GIRRBACH - Roto RFID 3,0 Wax Denture - # 760696	US\$ 68.59	+Info
468267	AMANN GIRRBACH - Roto RFID 2,5 Model - # 760697	US\$ 59.96	+Info
468239	AMANN GIRRBACH Roto RFID 2,5 Sintron - #760664	US\$ 78.16	+Info
467774	AMANN GIRRBACH - Calibration pin Roto RFID - # 181205	US\$ 129.95	+Info
468118	AMANN GIRRBACH - Roto RFID 1.2 ST - # 751001	US\$ 101.71	+Info
468119	AMANN GIRRBACH - Roto RFID 0.6 ST - # 751002	US\$ 101.71	+Info
468120	AMANN GIRRBACH - Roto RFID 1.5 T-Shape - # 751003	US\$ 114.56	+Info
468235	AMANN GIRRBACH - Roto RFID 2,5 ZI - # 760660	US\$ 84.19	+Info
468236	AMANN GIRRBACH - Roto RFID 1,0 ZI - # 760661	US\$ 84.19	+Info
468237	AMANN GIRRBACH - Roto RFID 0,6 ZI - # 760662	US\$ 84.19	+Info
468238	AMANN GIRRBACH - Roto RFID 0,3 ZI - # 760663	US\$ 84.19	+Info
468251	AMANN GIRRBACH - Roto RFID 2,5 SC - # 760680	US\$ 101.87	+Info
468240	AMANN GIRRBACH Roto RFID 1,0 Sintron - #760665	US\$ 78.16	+Info
468241	AMANN GIRRBACH Roto RFID 0,6 Sintron - #760666	US\$ 78.16	+Info
468242	AMANN GIRRBACH Roto RFID 0,3 Sintron - #760667	US\$ 78.16	+Info
468243	AMANN GIRRBACH - Roto RFID 0,6 PMMA/Wax - # 760670	US\$ 86.27	+Info
468244	AMANN GIRRBACH - Roto RFID 0,3 PMMA/Wax - # 760671	US\$ 86.27	+Info
468245	AMANN GIRRBACH - Roto RFID 2,0 Ti - # 760672	US\$ 68.59	+Info
468246	AMANN GIRRBACH - Roto RFID 1,0 Ti - # 760673	US\$ 68.59	+Info

Cad Cam Burs - UP3D

DC Coated Zirconia burs

Code	Name	Price	+
474321	UP3D - P52/53 Cad Bur 2mm	US\$ 35.64	+Info
474322	UP3D - P52/53 Cad Bur 1mm	US\$ 35.64	+Info
474323	UP3D - P52/53 Cad Bur 0.6mm	US\$ 35.64	+Info
474324	UP3D - P52/53 Cad Bur 0.3mm	US\$ 35.64	+Info

DLC PMMA Cad Cam Burs

Code	Name	Price	+
------	------	-------	---

474325	UP3D - PMMA Bur 2.0mm	US\$ 23.76	+Info
474326	UP3D - PMMA Bur 1.0mm	US\$ 23.76	+Info
474327	UP3D - PMMA Bur 0.6mm	US\$ 23.76	+Info

CAD/CAM Burs - AMANN GIRRBACH

CAD/CAM Burs - AMANN GIRRBACH

Code	Name	Price	+
464820	Roto Diamond Burs 1.8mm	US\$ 64.90	+Info
469827	Roto RFID Calibration Pin	US\$ 144.94	+Info

CAD/CAM Burs - BESQUAL

CAD/CAM Burs - BESQUAL

Code	Name	Price	+
417130	Zircon Reader for 2mm di	US\$ 19.20	+Info

CAD/CAM Burs - DENTSPLY

CAD/CAM Burs - DENTSPLY

Code	Name	Price	+
473682	ORIGINAL Step Bur #12S 6pk	US\$ 219.99	+Info
473683	ORIGINAL Cylinder Ptd Bur #12S 6pk	US\$ 260.00	+Info

CAD/CAM Burs - UP3D

CAD/CAM Burs - UP3D

Code	Name	Price	+
474339	UP3D - P52/53 Cad Bur 0.3mm	US\$ 35.64	+Info
474342	UP3D - P52/53 Cad Bur 0.6mm	US\$ 35.64	+Info
474343	UP3D - P52/53 Cad Bur 1.0mm	US\$ 35.64	+Info
474344	UP3D - P52/53 Cad Bur 2.0mm	US\$ 35.64	+Info

CERAMILL MOTION ORIGINAL - AMANN GIRRBACH

CERAMILL MOTION ORIGINAL - AMANN GIRRBACH

High Precision, long life sensing and milling tools for cost effective machining of ceramill blanks.

Available sizes:

- * 1.2mm
- * 0.2mm
- * 0.3mm and additional special patterns.

Material removal and dust removal are optimally controlled due to the tool geometry tailored to the material.

Extensive assortment of milling cutters for the management of a wide range of



indications. Custom selection for Zirconia Abutments selected from Amann Girrbach Burs and Readers
Sold by unit or as a set

#UsddAmann

AMANN GIRRBACH is a trademark of Amann Girrbach AG

Code	Name	Price	+
760751	AMANN GIRRBACH Zirconia Bur - F0.6 - #760751 (for item T0.6) Used for Implants...	US\$ 79.26	+Info
760709	AMANN GIRRBACH Zirconia Bur - F0.9mm #10 - #760709 (for item T0.9mm) Used for ...	US\$ 94.34	+Info
760711	AMANN GIRRBACH Zirconia Bur - F1.2 - #760711 (for item T1.2) Used for Implants...	US\$ 121.99	+Info
760721	AMANN GIRRBACH Zirconia Bur - F2 - #760721 (for item T2) Used for Implants & A...	US\$ 105.50	+Info
760720	AMANN GIRRBACH Zirconia Bur - UF2mm #28 - #760720 (for item UT2mm) Used for Im...	US\$ 84.38	+Info
760730	AMANN GIRRBACH Zirconia Bur UF3mm - #760730(for item UT3mm) Conical tools with...	US\$ 84.38	+Info
760701	AMANN GIRRBACH Zirconia Bur - F4 - #760701 (for item T4) Used for Implants & A...	US\$ 121.99	+Info
760740	AMANN GIRRBACH Zirconia Bur - CF2mm #14 - #760740 (for item CT2mm) Used for Im...	US\$ 84.38	+Info
760761	AMANN GIRRBACH Zirconia Bur - SF2mm - #760761(for item ST2mm) Whool without fa...	US\$ 64.13	+Info
760771	AMANN GIRRBACH Zirconia Reader SF1.2mm - #760771(for item ST1.2mm) Shank-type ...	US\$ 84.38	+Info
760781	AMANN GIRRBACH Zirconia Bur SF0.6mm - #760781 (for item ST0.6mm) Thin shank-ty...	US\$ 79.26	+Info
760791	AMANN GIRRBACH Zirconia Bur KF2.5mm #760791- (for item KT2.5mm) Ball tool for ...	US\$ 79.26	+Info
760892	AMANN GIRRBACH Zirconia Reader - KT2mm #17 (KT1.8mm) - #760892 (for item KF2m...	US\$ 94.34	+Info
760881	AMANN GIRRBACH Zirconia Reader ST0.6mm - #760881 (for item SF0.6mm) Thin shank...	US\$ 84.78	+Info
760891	AMANN GIRRBACH Zirconia Reader KT2.5mm - #760891(for item KF2.5mm) Ball tool f...	US\$ 84.78	+Info
760871	AMANN GIRRBACH Zirconia Reader ST1.2mm - #760871 (for item SF1.2mm) Shank-type...	US\$ 80.71	+Info
760861	AMANN GIRRBACH Zirconia Reader ST2mm - #760861 (for item SF2mm) Whool without ...	US\$ 80.71	+Info
760851	AMANN GIRRBACH Zirconia Reader - T0.6 - #760851 (for item F0.6) Used for Impla...	US\$ 84.78	+Info
760840	AMANN GIRRBACH Zirconia Reader - CT2mm #13 - #760840 (for item CF2mm) Used for...	US\$ 80.71	+Info
760831	AMANN GIRRBACH Zirconia Reader - T3 - #760831 (for item F3) Used for Implants ...	US\$ 105.50	+Info
760830	AMANN GIRRBACH Zirconia Reader UT3mm - #760830 (for item UF3mm) Conical tools ...	US\$ 80.71	+Info
760821	AMANN GIRRBACH Zirconia Reader - T2 - #760821 (for item F2) Used for Implants ...	US\$ 71.70	+Info
760820	AMANN GIRRBACH Zirconia Reader - UT2mm #27 - #760820 (for item UF2mm) Used for...	US\$ 80.71	+Info
760811	AMANN GIRRBACH Zirconia Reader - T1.2 - #760811 (for item F1.2) Used for Impla...	US\$ 105.50	+Info
760809	AMANN GIRRBACH Zirconia Reader - T0.9mm #9 - #760809 (for item F0.9mm) Used fo...	US\$ 94.34	+Info
760801	AMANN GIRRBACH Zirconia Reader - T4 - #760801 (for item F4) Used for Implants ...	US\$ 72.92	+Info
760792	AMANN GIRRBACH Zirconia Bur - KF1.8mm (KF2mm) #18 - #760792 (for item KT1.8mm...	US\$ 105.50	+Info
420319	AMANN GIRRBACH #F3 Ceramill Roto Bur #760731	US\$ 105.50	+Info

MANUAL ZIRCONIA BURS - ECCO

MANUAL ZIRCONIA BURS - ECCO

- * Compatible with Zircozhan
- * Comatible with Amann Girrbach Ceramill Manual
- * Compatiblewith Ecco-Mill System
- * 03 mm shank
- * 47 mm Lengh
- * Made in Germany

If you change bur brands we strongly recomend to buy the



correspondent Reader the first time.

KIT

Code	Name	Price	+
102214	ECCO - Zircon Drill bur/read Kit - 10 pieces - 5 readers / 5 burs. 3mm shank for...	US\$ 275.00	+Info

ECCO

Code	Name	Price	+
102352	ECCO - Zircon Bur 1mm diam x 8 mm length - each - 3mm shank Total length 47mm (...	US\$ 35.64	+Info
102351	ECCO - Zircon Bur 1mm diam x 15 mm length - each - 3mm shank - Total length 47mm...	US\$ 35.64	+Info
102350	ECCO - Zircon Bur 2mm diam - each - 3mm shank - Total length 47mm (compatible w...	US\$ 35.64	+Info
103404	ECCO - Zircon Bur 3mm diam - each - 3mm shank Total length 47mm (compatible wit...	US\$ 35.64	+Info
102349	ECCO - Zircon Bur 4mm diam - each - 3mm shank Total length 47mm (compatible wit...	US\$ 35.64	+Info

READERS

Code	Name	Price	+
102355	ECCO - Zircon Reader for 1mm diam x 12 mm length (0.8mm diam) each - 3mm shank...	US\$ 19.20	+Info
103405	ECCO - Zircon Reader for 3mm diam -(2.4mm diam) each - 3mm shank Total length ...	US\$ 19.20	+Info
102353	ECCO - Zircon Reader for 4mm diam -(3.2mm diam) each - 3mm shank Total length ...	US\$ 19.20	+Info

HANSONG CAD CAM BURS



Amann Girrbach Burs

Code	Name	Price	+
436910	Hansong Milling Burs Amann Girrbach	US\$ 27.89	+Info
436913	Hansong Milling Burs Amann Girrbach	US\$ 28.87	+Info
436926	Hansong Milling Burs Amann Girrbach	US\$ 31.89	+Info

Roland Burs

Code	Name	Price	+
436911	Hansong Milling Burs Roland Type-2	US\$ 42.74	+Info
436916	Hansong Milling Burs Roland Type- 2	US\$ 37.88	+Info
436919	Hansong Milling Burs Roland Type- 2	US\$ 37.88	+Info
436920	Hansong Milling Burs Roland Type- 2	US\$ 37.88	+Info
436922	Hansong Milling Burs Roland Type- 2	US\$ 37.88	+Info
436931	Hansong Milling Burs Roland Type-2O	US\$ 41.58	+Info
436939	Hansong Milling Burs Roland Type- 2	US\$ 37.88	+Info
436941	Hansong Milling Burs Roland Type-2	US\$ 41.58	+Info

Sirona Type Burs

Code	Name	Price	+
436927	Hansong Milling Burs Sirona Type- 2	US\$ 37.88	+Info
436930	Hansong Milling Burs Sirona Type- 2	US\$ 42.74	+Info
436934	Hansong Milling Burs Sirona Type- 2	US\$ 40.43	+Info
436936	Hansong Milling Burs Sirona Type- 2	US\$ 42.74	+Info
436937	Hansong Milling Burs Sirona Type- 2	US\$ 42.74	+Info
436943	Hansong Milling Burs Sirona Type- 2	US\$ 38.11	+Info

VHF Burs

Code	Name	Price	+
436912	Hansong Milling Burs VHF Type-2 VTR	US\$ 37.88	+Info
436914	Hansong Milling Burs VHF Type-2 VTR	US\$ 37.88	+Info
436923	Hansong Milling Burs VHF Type-2 VDR	US\$ 38.11	+Info
436924	Hansong Milling Burs VHF Type-2 VDR	US\$ 38.11	+Info
436929	Hansong Milling Burs VHF Type-2 VDR	US\$ 40.43	+Info
436932	Hansong Milling Burs VHF Type-2 VDR	US\$ 42.74	+Info
436933	Hansong Milling Burs VHF Type-2 VDR	US\$ 34.65	+Info
436935	Hansong Milling Burs VHF Type-2 VDR	US\$ 41.58	+Info
436938	Hansong Milling Burs VHF Type-2 VDR	US\$ 38.11	+Info
436944	Hansong Milling Burs VHF Type-2 VDR	US\$ 42.74	+Info

Wieland Type Burs

Code	Name	Price	+
436915	Hansong Milling Burs Wieland Type-	US\$ 41.58	+Info
436918	Hansong Milling Burs Wieland Type-	US\$ 42.74	+Info
436928	Hansong Milling Burs Wieland Type-	US\$ 38.11	+Info
436940	Hansong Milling Burs Wieland Type-	US\$ 41.58	+Info
436942	Hansong Milling Burs Wieland Type-	US\$ 45.05	+Info

ZirkonZahn Burs

Code	Name	Price	+
436917	Hansong Milling Burs Zirkonzahn Typ	US\$ 45.05	+Info
436921	Hansong Milling Burs Zirkonzahn Typ	US\$ 50.82	+Info
436925	Hansong Milling Burs Zirkonzahn Typ	US\$ 45.05	+Info

ORIGINAL ROLAND CAD/CAM MILLING BURS



CARBIDE

BALL END

Code	Name	Price	+
420215	ROLAND 0.3mm Carbide Ball End Mill - High Quality ZCB-15D-US	US\$ 103.39	+Info
420214	ROLAND 0.6mm Carbide Ball End Mill ZCB-30D-US	US\$ 103.39	+Info
420213	ROLAND 1mm Carbide Ball End Mill - High Quality ZCB-50D-US	US\$ 103.39	+Info
420212	ROLAND 2mm Carbide Ball End Mill - High Quality ZCB-100D-US	US\$ 103.39	+Info

FLAT END

Code	Name	Price	+
420300	ROLAND 0.6mm Carbide Flat End Mill #PT06-10 Extremo Chato	US\$ 36.97	+Info
420301	ROLAND 1mm Carbide Flat End Mill #PT1- 20 Extremo Chato	US\$ 38.51	+Info
420302	ROLAND 2mm Carbide Flat End Mill #PT2- 20 Extremo Chato	US\$ 38.51	+Info
420306	ROLAND 3mm Flat End Carbide Milling Bur #ZRB-100D Extremo Chato	US\$ 115.78	+Info

TORIC

Code	Name	Price	+
420334	ROLAND Rapid Razor-Sharp 3mm Carbide Toric End Bur SKU#RSB-150D-US Made in Germany	US\$ 121.71	+Info

DIAMOND COATED

Code	Name	Price	+
420218	ROLAND 0.6mm DIAMOND COATED DC Ball End Mill #ZDB-30D-US	US\$ 144.46	+Info
420217	ROLAND 1mm DIAMOND COATED DC Ball End Mill ZDB-50D-US	US\$ 144.46	+Info
420216	ROLAND 2mm DIAMOND COATED DC Ball End Mill ZDB-100D-US	US\$ 144.46	+Info

DIAMOND GRINDING BURS

For Glass Ceramics

LONG DIAMOND GRINDING BURS

Code	Name	Price	+
420316	ROLAND 0.5mm Long Diamond Grinding Bur #ZGB2- 25D ZBG2- 25D Wet Milling Cad Burs...	US\$ 86.71	+Info

420315	ROLAND 1mm Long Diamond Grinding Bur ZGB2-50D ZBG2-50D Wet Milling Cad Burs for	US\$ 86.71	+Info
420314	ROLAND 2.5mm Long Diamond Grinding Bur #ZGB2-125D ZBG2-125D Wet Milling Cad Burs...	US\$ 86.71	+Info

REGULAR DIAMOND GRINDING BURS

Code	Name	Price	+
420312	ROLAND 0.5mm Regular Diamond Grinding Bur #ZGB- 25D Wet Milling Cad Burs for Dis...	US\$ 103.39	+Info
420313	ROLAND 1mm Regular Diamond Grinding Bur #ZGB-50D Wet Milling Cad Burs for Disili...	US\$ 103.39	+Info
420311	ROLAND 2.5mm Regular Diamond Grinding Bur #ZGB-125D Wet Milling Cad Burs for Dis...	US\$ 103.39	+Info

SHORT DIAMOND GRINDING BURS

Code	Name	Price	+
420318	ROLAND 2.5mm Short Diamond Grinding Bur #ZGB2-125DS ZBG2-125DS Wet Milling Cad B...	US\$ 86.71	+Info

LLOLIPOP

Code	Name	Price	+
420333	ROLAND 2mm Carbide Lollipop End Mill SKU#LOLLIPOP Made in Germany	US\$ 34.70	+Info

SPECIAL COATED

Code	Name	Price	+
420221	ROLAND 0.6mm Special Coated Ball End Mill-Ccomposites Hybrid #ZDB-30D-HY	US\$ 169.96	+Info
420220	ROLAND 1mm Special Coated Ball End Mill-Ccomposites Hybrid #ZDB-50D-HY	US\$ 169.96	+Info
420219	ROLAND 2mm Special Coated Ball End Mill-Ccomposites Hybrid #ZDB-100D-HY	US\$ 169.96	+Info

ROTO DIAMOND BURS - AMANN GIRRBACH

ROTO DIAMOND BURS - AMANN GIRRBACH

Available en 3 sizes 1.0mm - 1.4mm - 1.8mm

Made in Germany



Code	Name	Price	+
420952	AMANN GIRRBACH Roto Diamond Burs (Wet milling disilicate) Diameter 1.0mm #76...	US\$ 60.26	+Info
420951	AMANN GIRRBACH Roto Diamond Burs (Wet milling disilicate) Diameter 1.4mm #76...	US\$ 60.26	+Info
420953	AMANN GIRRBACH Roto Diamond Burs (Wet milling disilicate) Diameter 1.4mm #76...	US\$ 61.34	+Info
420950	AMANN GIRRBACH - Roto Diamond Burs (Wet milling disilicate) Diameter 1.8mm	US\$ 39.10	+Info

YETI & VHF GLASS CERAMIC MILLING BURS 4 & 5 Axes - YETI

YETI - Milling Cutter Glass Ceramic CAD/CAM Burs

for 4 & 5 Axes machines



Code	Name	Price	+
417062	YETI YETI digi MILL Milling Cutter M200-R4-32 (K/4) # 481-1121	US\$ 57.38	+Info
417064	YETI DIGI-Milling Cutters - Glass Ceramics (LiSi) CYLINDER BALL END (Round) - ...	US\$ 59.70	+Info
417063	YETI DIGI-Milling Cutter Glass Ceramic (LiSi) CONICAL - Radius 0.60 mm Lenght...	US\$ 59.70	+Info
417091	YETI DIGI-Milling Cutter Glass Ceramic (LiSi) CYLINDER BALL END (Flat) - Toru...	US\$ 59.70	+Info
417119	YETI DIGI-Milling Cutters - Glass Ceramics (LiSi) CYLINDER FLAT END - 1.20 mm...	US\$ 59.70	+Info
417120	YETI DIGI-Milling Cutters - Glass Ceramics (LiSi) CYLINDER FLAT END - 2.60 mm...	US\$ 59.70	+Info
417121	YETI DIGI-Milling Cutter Glass Ceramic SET . 9 Pcs (wet diamond Sintered for Di...	US\$ 430.11	+Info

YETI DLC CAD/CAM MILLING BURS - YETI

YETI - DLC CAD/CAM MILLING BURS YETI



Code	Name	Price	+
417122	YETI DIGI-Milling Cutter Zirconia PMMA Wax 0.60 mm Lenght: 40 mm #481-0506 U06...	US\$ 67.03	+Info
417123	YETI DIGI-Milling Cutter Zirconia PMMA Wax 2.0 mm Lenght: 40 mm #481-0420 Z200...	US\$ 69.24	+Info
417124	YETI DIGI-Milling Cutter Zirconia PMMA Wax 1.0 mm Lenght: 40 mm #481-0410 Z100...	US\$ 69.24	+Info


Z-GRINDERS ZIRCONIA GRINDERS - ATTRITOR

Z-GRINDERS ZIRCONIA GRINDERS - ATTRITOR

Proud of this range of Sinter Diamond stones specially designed for working on all types of **Ceramics & Zirconia materials**.

With our unique combination of **Aluminum Oxide**, **Glass**, **Diamond Particle** and a range of ingredients that create this spectacular range of abrasive instruments we have created a product that is second to none in the market.

Our unique product far exceeds the competitors in cutting ability and life span.



Z-GRINDERS

Diamond impregnated polishers for ceramics and zirconia materials. Superior in performance and much more economical due to its durability.

Code	Name	Price	+
112947	ATTRITOR - Z-Grinder for Zirconia Wheel 15mm x 3.5mm - medium Grit - 1pk- Sin...	US\$ 20.00	+Info
112228	ATTRITOR - Z-Grinder for Zirconia Wheel 22mm x 2mm - medium Grit - 1pk- Sinte...	US\$ 20.00	+Info
100754	ATTRITOR - Z-Grinder for Zirconia Wheel 22mm x 4.5mm- medium Grit - 1pk- Sint...	US\$ 20.00	+Info
112948	ATTRITOR - Z-Grinder for Zirconia Taper 3.5mm x 11mm - medium Grit - 1pk- S...	US\$ 12.50	+Info
100756	ATTRITOR - Z-Grinder for Zirconia Large Barrel 6mm x 13mm - medium Grit - 1pk...	US\$ 15.00	+Info
100755	ATTRITOR - Z-Grinder for Zirconia Medium Cylinder bARREL 5mm x13mm - medium Gr...	US\$ 12.50	+Info

Z-ZIRCON ZIRCONIA GRINDERS - ATTRITOR

Z-GRINDERS ZIRCONIA GRINDERS - ATTRITOR

Proud of this range of Sinter Diamond stones specially designed for working on all types of **Ceramics & Zirconia materials**.

With our unique combination of **Aluminum Oxide**, **Glass**, **Diamond Particle** and a range of ingredients that create this spectacular range of abrasive instruments we have created a product that is second to none in the market.

Our unique product far exceeds the competitors in cutting ability and life span.



Z-ZIRCON

Code	Name	Price	+
115225	ATTRITOR - Zirconia Grinder Z-Zircon Medium INverted Cone - medium Grit - 1pk-...	US\$ 22.50	+Info
114356	ATTRITOR - Zirconia Grinder Z-Zircon discs 22mm x 0.7mm- medium Grit - 1pk- Si...	US\$ 19.80	+Info
114177	ATTRITOR - Zirconia Grinder Z-Zircon Wheel 22mm x 4.5mm- medium Grit - 1pk- Si...	US\$ 20.00	+Info
114176	ATTRITOR - Zirconia Grinder Z-Zircon Wheel 22mm x 2mm - medium Grit - 1pk- Sin...	US\$ 19.60	+Info
114175	ATTRITOR - Zirconia Grinder Z-Zircon Wheel 15mm x 3.5mm - medium Grit - 1pk- S...	US\$ 30.00	+Info

ABRASIVE DISCS - MPF

ABRASIVE DISCS - MPF

Features:

Highly durable discs for ceramic and alloy restoration, without the risk of chipping or cracking the surface.

Quantity: 100pk

Available in: Small - Medium - Large



Code	Name	Price	+
123935	MPF BRUSH Aluminum Oxide Abrasive Disc small (22 x 0-2-0-25 x 1-8) 100 pcs 123...	US\$ 39.86	+Info
123933	MPF BRUSH Aluminum Oxide Abrasive Disc medium (25 x 0-6 x 1-8) 100 pcs 123-0001...	US\$ 32.85	+Info
123934	MPF BRUSH - Aluminum Oxide Abrasive Disc large (38 x 0-6 x 1-8) 100 pcs 123-00...	US\$ 46.50	+Info

BESQUAL ZIR-STONES ABRASIV ES

Premium diamond abrasive grinding stones for zirconia and all ceramic materials.



ZIR-STONES

Code	Name	Price	+
469856	Zircon Diamond Stone KIT 7pk	US\$ 57.33	+Info

469857	Zircon Diam Stone Concave #9 5pk	US\$ 42.63	+Info
469858	Zircon Diam Stone Knife-Edge #10	US\$ 42.63	+Info
469859	Zircon Diam Stone Wheel #11 5pk	US\$ 42.63	+Info
469871	Zircon Diam Stone Taper #20 5pk	US\$ 42.63	+Info
469872	Zircon Diam Stone Barrel #39 5pk	US\$ 42.63	+Info
469873	Zircon Diam Stone Inv Cone #42 5pk	US\$ 42.63	+Info
469874	Zircon Diam Stone Bullet #13	US\$ 42.63	+Info

COREMASTER - SHOFU

COREMASTER - SHOFU

Features:

CoreMaster is designed for adjusting marginal areas and the thickness of aluminum oxide and zirconia copings.

- * Diamond impregnated polishers eliminate microfractures
- * Two grit sizes and three shapes
- * Contouring and finishing
- * No water cooling needed
- * Use CoreMaster for: zirconia frameworks, alumina frameworks and porcelains



US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.

Standard Grit



Code	Name	Price	+
106767	SHOFU HP Coremaster KN7 HP - 2pk - #0139	US\$ 46.18	+Info
106772	SHOFU HP Coremaster WH6 HP - 2pk - #0140C	US\$ 38.44	+Info
106773	SHOFU HP Coremaster Standard Unmounted Wheel-1pk - #0141	US\$ 59.71	+Info

WHITE PEARLS ZIRCONIA GRINDERS - MPF Brush Co

WHITE PEARLS ZIRCONIA GRINDERS - MPF Brush Co



Features:

Zirconia has become very popular in the dental industry in the past few years and its growing day by day.

Due to zirconias hardness and sensativity to heat, there was a great need for a rotary instrument that could cut and abrade this material with a minimum of heat generation as possible.

We have created a rotary instrument with a very high content of fine diamond particle which has

been bonded with special binders that hold all the elements together, creating a superb Zirconia cutting material with exceptional life span, 2-3 times more life than any product on the market today. With our high content of diamond particle our product can cut and abrade with very little pressure needed from the handpiece, thus much less heat generation to the zirconia substructure.

Utilizing stainless iron mandrels for extra strength we guarantee a very high end looking and functioning product. We have laser printed onto our shanks for specific marking & branding for this exclusive product line.

#UsddMpf

Code	Name	Price	+
123210	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Inverted Small Cone 6mm x 8mm - 1p...	US\$ 43.38	+Info
123209	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Inverted Cone 12mm x 6mm - 1pk- #...	US\$ 39.04	+Info
123208	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Cone 3.5mm x 11mm - 1pk- #140-0007...	US\$ 43.38	+Info
123207	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Large Barrel 6.5mm x 13mm - 1pk- ...	US\$ 43.38	+Info
123205	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Thin discs 22mm x 0.7mm - 1pk- #1...	US\$ 54.85	+Info
123206	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Small Barrel 4.8mm x 13mm - 1pk- ...	US\$ 43.38	+Info
123203	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Small discs 15mm x 3.5mm - 1pk- #1...	US\$ 54.85	+Info
123204	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel discs 22mm x 2mm - 1pk- #140-0003 ...	US\$ 54.85	+Info
123202	MPF BRUSH -WHITE PEARL Zirconia Grinder Wheel Large discs 22mm x 4.5mm - 1pk- #1...	US\$ 63.79	+Info

ZIRCOCERA ZIRCONIA GRINDER

S - MPF Brush Co

ZIRCOCERA ZIRCONIA GRINDERS - MPF Brush Co



Features:

The GSS Series are a diamond impregnated rotary instrument of the highest quality, they boast a very high content of industrial diamond, bound together with a special ceramic binder. The GSS Series rotary instruments have extraordinary abrasive ability, they can cut through Zirconia with a minimal amount of pressure, and this alleviates the chances of micro fractures within the Zirconia substructure.

GSS Series have the ability to cut all types of ceramic material with great ease, and due to the fine grain.

Industrial diamond and special binder, the ceramic surface is left with a smooth, clean finish.

Due to the composition of the diamond to binder ratio, the GSS Series boasts a very high service life, ergonomically it is by far the best solution for the ceramist today!

#UsddMpf

Code	Name	Price	+
123201	MPF BRUSH -ZIRCO CERA GSS Zirconia Grinder Wheel Inverted Small Cone 6mm x 8mm -...	US\$ 37.64	+Info
123200	MPF BRUSH -ZIRCO CERA GSS Zirconia Grinder Wheel Inverted Cone 12mm x 6mm - 1pk-...	US\$ 37.64	+Info

123199	MPF BRUSH - ZIRCO CERA GSS Zirconia Grinder Wheel Cone 3.5mm x 11mm - 1pk- #121-0...	US\$ 37.64	+Info
123198	MPF BRUSH - ZIRCO CERA GSS Zirconia Grinder Wheel Large Barrel 6.5mm x 13mm - 1pk...	US\$ 37.64	+Info
123197	MPF BRUSH ZIRCO CERA GSS Zirconia Grinder Wheel Small Barrel 4.8mm x 13mm - 1pk-...	US\$ 36.89	+Info
123196	MPF BRUSH - ZIRCO CERA GSS Zirconia Grinder Wheel Thin discs 22mm x 0.7mm - 1pk- ...	US\$ 47.84	+Info
123194	MPF BRUSH - ZIRCO CERA GSS Zirconia Grinder Wheel Small discs 15mm x 3.5mm - 1pk-...	US\$ 47.84	+Info
123195	MPF BRUSH - ZIRCO CERA GSS Zirconia Grinder Wheel discs 22mm x 2mm - 1pk- #121-0...	US\$ 47.84	+Info
123192	MPF BRUSH - ZIRCO CERA GSS Zirconia Grinder Wheel Large discs 22mm x 4.5mm - 1pk-...	US\$ 55.49	+Info

Zirconia & Emax Grinders - ATTRITOR

Zirconia & Emax Grinders - ATTRITOR

Code	Name	Price	+
416991	ZIRCONIA GRINDER KIT 8pcs	US\$ 164.90	+Info

Zirconia Contouring - BESQUAL

BESQUAL - Zirconia Countouring KIT

This contouring Kit is ideal use for pre-sintered (unsintered) Zirconia. Switching burs while doing contouring work is such a hassle, this contouring 2-in-1 abrasive/polisher will solve just that. From quick reduction to final smooth finishing, you can accomplish great result in no time. This kit includes three other shapes for those hard to reach ea



Code	Name	Price	+
469842	Zirconia Contouring Kit	US\$ 36.38	+Info
469843	Z-KC51 Pre-Sint Md Cn 5.6x13mm	US\$ 7.65	+Info
469844	Zirconia Medium Disc 15x0.3mm	US\$ 7.65	+Info
469845	Dual Grit Fine/Medium Egg 7x19mm	US\$ 13.20	+Info
469846	Dual Grit Fine/Medium Egg 10x23mm	US\$ 12.80	+Info
469847	Medium Flame 7x13mm	US\$ 7.65	+Info

ZMAX ZIRCONIA GRINDERS - DIAMOND IMPREGNATED DISCS - MPF Brush Co

ZMAX ZIRCONIA GRINDERS - DIAMOND IMPREGNATED DISCS - MPF Brush Co



Features:

Zirconia has become very popular in the dental industry in the past few years and its growing day by day.

Due to zirconias hardness and sensativity to heat, there was a great need for a rotary instrument that could cut and abrade this material with a minimum

of heat generation as possible.

We have created a rotary instrument with a very high content of fine diamond particle which has been bonded with special binders that hold all the elements together, creating a superb Zirconia cutting material with exceptional life span, 2-3 times more life than any product on the market today. With our high content of diamond particle our product can cut and abrade with very little pressure needed from the handpiece, thus much less heat generation to the zirconia substructure.

Utilizing stainless iron mandrels for extra strength we guarantee a very high end looking and functioning product. We have laser printed onto our shanks for specific marking & branding for this exclusive product line.

The ZMax come in 9 different shapes

#UsddMpf

Code	Name	Price	+
123186	MPF BRUSH -ZMAX Zirconia Grinder Wheel Thin 22mm x 0.7mm - 1pk- # 120-0004 1200004	US\$ 55.49	+Info
123185	MPF BRUSH -ZMAX Zirconia Grinder Wheel 22mm x 2mm - 1pk- #120-0003 1200003	US\$ 55.49	+Info
123184	MPF BRUSH -ZMAX Zirconia Grinder Wheel Small 15mm x 3.5mm - 1pk- #120-0002 1200002	US\$ 55.49	+Info
123183	MPF BRUSH -ZMAX Zirconia Grinder Wheel Large 22mm x 4.5mm - 1pk- # 120-0001 120...	US\$ 45.45	+Info
123190	MPF BRUSH -ZMAX Zirconia Grinder Wheel Inverted Cone 12mm x 6mm - 1pk- #120-000...	US\$ 47.84	+Info
123188	MPF BRUSH -ZMAX Zirconia Grinder Wheel Small Barrel 6.5mm x 13mm - 1pk- #120-00...	US\$ 32.53	+Info
123191	MPF BRUSH -ZMAX Zirconia Grinder Wheel Inverted Small Cone 6mm x 8mm - 1pk- #12...	US\$ 47.84	+Info
123187	MPF BRUSH -ZMAX Zirconia Grinder Wheel Small Barrel 4.8mm x 13mm - 1pk- #120-00...	US\$ 47.84	+Info
123189	MPF BRUSH -ZMAX Zirconia Grinder Wheel Cone 3.5mm x 11mm - 1pk- #120-0007 1200007	US\$ 47.84	+Info

CAD/CAM abrasives

Zirconia & Emax Polishers



RAINBOW CERAMIC POLISHER - MPF Brush Co

RAINBOW CERAMIC POLISHER - MPF Brush Co

Features:

Premium diamond abrasive wheels for perfectly sealed ceramic or zirconia restoration.

For minimal reduction on zirconia marginal regions and all types of ceramic materials.

Available in color red, blue, gray and yellow



Code	Name	Price	+
123929	MPF BRUSH -Rainbow- Ceramic Polisher- Coarse- RED 122-0001 - #MPF BRUSH co	US\$ 32.60	+Info
123930	MPF BRUSH Rainbow- Ceramic Polisher- Medium- BLUE 122-0002 - #MPF BRUSH co	US\$ 46.57	+Info
123931	MPF BRUSH Rainbow- Ceramic Polisher- Grey- Fine 122-0003 - #MPF BRUSH co	US\$ 32.60	+Info
123932	MPF BRUSH Rainbow- Ceramic Polisher- YELLOW- High Shine 122-0004 - #MPF BRUSH co	US\$ 44.24	+Info

Zirconia & Emax Polishers - AMANN GIRRBACH

Zirconia & Emax Polishers - AMANN GIRRBACH

Code	Name	Price	+
468122	Ceramill Zolid FX Starter Kit	US\$ 3,221.73	+Info
468865	Zolid green-state finishing Kit	US\$ 162.34	+Info
468873	Zolid preparation dent kit	US\$ 85.97	+Info
468880	Zolid sinter-state polishing kit	US\$ 191.46	+Info
468893	Zolid polishing dent kit update	US\$ 104.34	+Info

Zirconia & Emax Polishers - ATTRITOR


Zirconia & Emax Polishers - ATTRITOR

Code	Name	Price	+
416558	ZIRCONIA POLISHER KIT 6pcs	US\$ 125.00	+Info
470735	NANO Ceramic POLISHER KIT 5pcs	US\$ 17.50	+Info
470762	Zirconia Luster Polishers 6pcs	US\$ 97.50	+Info
470835	Zirc / Disil Polish Adjust 5 pcs	US\$ 32.34	+Info
470836	Zirc / Disil Polish Adjust 5 pcs	US\$ 30.36	+Info
470837	Zirc / Disil Polish Adjust 5 pcs	US\$ 32.34	+Info
470838	Zirc / Disil Polish Adjust 5 pcs	US\$ 31.68	+Info
470839	Zirc / Disil Polish Adjust 5 pcs	US\$ 31.68	+Info

ZIRCONIA POLISHERS POINT - ATTRITOR

ZIRCONIA POLISHERS POINT - ATTRITOR

- * For fast shaping
- * Scratch removal
- * Leaves satin finish on all porcelains and final stage gloss finish on cercon
- * Alumina and Zirconia copins



Code	Name	Price	+
100933	ATTRITOR - Dura Point Extra Fine - BLUE-Silicone Diamond point - for fast shapin...	US\$ 18.00	+Info
100932	ATTRITOR - Dura Point Fine PINK -Silicone Diamond point - for fast shaping -sca...	US\$ 14.40	+Info
100934	ATTRITOR - Dura Point Coarse - GREEN Silicone Diamond point - for fast shaping -...	US\$ 19.38	+Info

HATHO POLISHING KIT

Thermoplastics
like *never* before



Packing Diamond D®
Denture Acrylic with
Tecnoflask




Keystone Industries: Laboratory



ZIRCONIA POLISHERS WHEELS - ATTRITOR

ZIRCONIA POLISHERS WHEELS - ATTRITOR


- * Diamond Impregnated Silicone Discs Two tone (22 x3mm)
- * Diamond Impregnated Rubber Discs Mono Tone
- * Diamond Impregnated Silicone Wheels



HIGH SHINE SILICON WHEELS

These Mono tone color SILICONE WHEELS (Square edge) have a high content of refined diamond particle, which is bonded together with odorless Japanese Silicone material.

- * EXTRA COARSE - BROWN
- * COARSE - GRAY
- * MEDIUM - BLUE
- * FINE - PINK
- * EXTRA FINE - WHITE



Code	Name	Price	+
100930	ATTRITOR - Dura Wheel - SC Super Coarse BROWN - Silicone Diamond Wheel - for fa...	US\$ 15.36	+Info
100723	ATTRITOR - Dura Wheel - Coarse- GRAY -Silicone Diamond Wheel - For high shining ...	US\$ 15.36	+Info
100722	ATTRITOR - Dura Wheel Medium - BLUE - Silicone Diamond Wheel - for slow shaping ...	US\$ 15.36	+Info
100763	ATTRITOR - Dura Wheel - Fine PINK -Silicone Diamond Wheel - for shining -scratch...	US\$ 15.36	+Info
115563	ATTRITOR - Dura Wheel - Super Fine WHITE -Silicone Diamond Wheel - for shining -...	US\$ 11.83	+Info

HIGH SHINE SILICONE DISCS

These two tone color SILICONE DISCS have a very high content of refined diamond particle, which is bonded together with odorless Japanese Rubber material.



This product line is extremely popular with technicians around the world; the dimensions are 27mm x 0.7mm disc shape.

These wheels can be used to abrade and polish ZIRCONIA AND ALL TYPES OF CERAMICS We advise that the RPM are maintained to our specifications so as to obtain the optimum result and allows for a longer life span of the product.

Code	Name	Price	+
112956	ATTRITOR - Zirconia High Shine YELLOW Super Fine - 1pk- 27mm x 0.7mm- Bicolor Di...	US\$ 15.36	+Info
112957	ATTRITOR - Zirconia High Shine Grey Fine - 1pk- 27mm x 0.7mm- Bicolor Diamond Im...	US\$ 17.28	+Info
112958	ATTRITOR - Zirconia High Shine BLUE Medium - 1pk-27mm x 0.7mm- Bicolor Diamond I...	US\$ 15.36	+Info
112959	ATTRITOR - Zirconia High Shine RED Coarse - 1pk-27mm x 0.7mm- Bicolor Diamond Im...	US\$ 15.21	+Info

POLISHING RUBBER DISCS

These Mono tone color rubber discs have a very high content of refined diamond particle, which is bonded together with odorless Japanese Rubber material.

This product line is extremely popular with technicians around the world; the dimensions are 27mm x 0.7mm disc shape. These wheels can be used to abrade and polish zirconia and all types of dental ceramic.

We advise that the RPM are maintained to our specifications so as to obtain the optimum result and allows for a longer life span of the product.

There abrasive grade range from Super coarse, coarse, medium, fine and superfine.



Code	Name	Price	+
112953	ATTRITOR - Zirconia polisher Light GRAY Super Fine - 1pk-27mm x 0.7mm- Diamond I...	US\$ 17.28	+Info
112952	ATTRITOR - Zirconia polisher RED Fine - 1pk- 27mm x 0.7mm- Diamond Impregnated R...	US\$ 15.36	+Info
112951	ATTRITOR - Zirconia polisher BLUE Medium - 1pk- 27mm x 0.7mm- Diamond Impregnate...	US\$ 15.36	+Info
112950	ATTRITOR - Zirconia polisher GREEN Coarse - 1pk- 27mm x 0.7mm- Diamond Impregnat...	US\$ 15.36	+Info
112949	ATTRITOR - Zirconia polisher GRAY SC Super Coarse - 1pk- 27mm x 0.7mm -Diamond ...	US\$ 17.28	+Info

Catalog Index / Indice

CAD/CAM Burs	1
Zirconia & Emax Grinders	16
Zirconia & Emax Polishers	22