

Grafting Material

Resorbable Bone Grafting Mater



BIOVIVA Allograft Mineralized Bone - 3D Dent

Code	Name	Price	+
473996	BIOVIVA BIOVIVA Allograft Mineraliz	US\$ 108.80	+Info
473997	BIOVIVA BIOVIVA Allograft Mineraliz	US\$ 200.00	+Info
473998	BIOVIVA BIOVIVA Allograft Mineraliz	US\$ 54.45	+Info
473999	BIOVIVA BIOVIVA Allograft Mineraliz	US\$ 91.89	+Info

ESSENTIALS IN IMPLANTOLOGY- OUSIA

ESSENTIALS IN IMPLANTOLOGY- OUSIA

OUSIA - Essentials in Implantology OUSIA - we are excited to introduce the latest products to the Dental Implantology market. This product offering combines both the latest advances in science and technology allowing for optimal results.



BIOXEN

An Osteoconductive Matrix Ousia bioXEN is a naturally sourced bovine bone mineral. The unique manufacturing process of bioXEN results in a bone grafting particulate with the crystalline structure of carbonate apatite, which is a structure similar to intact human bone. In addition, the manufacturing process for bioXEN allows for a combination of a porous structure for vascular access with a roughened surface that is favorable to cell adhesion and subsequent bone matrix deposition.

Available: 1cc- 2.5cc- 4cc- 5cc/ Cancellous Particulate- Xenograft Part

CALCIUM

Medical Crade Calcium Sulfate Kit

Available: 1g- 5g

GRAF

The best of cortical and cancellous features for the ideal mixture
Remodels completely
Provides scaffold for cell ingrowth and acts as a guide for rapid bone regeneration

Available: 1cc- 2.5cc- 5cc

Bio Cal

Code	Name	Price	+
123164	OUSIA - BioCal .5g medical grade calcium sulfate Kit - #OUSBC1	US\$ 53.30	+Info
123165	OUSIA - BioCal 1g medical grade calcium sulfate Kit - #OUSBC2	US\$ 63.42	+Info

Bio Xen

Code	Name	Price	+
123155	OUSIA - BioXen .5cc - Mineralized Xenograft .25-1mm Cancellous Particulate - #OU...	US\$ 97.70	+Info
123159	OUSIA - BioXen 1cc - Mineralized Xenograft 1-2mm Particulate - #OUSM5	US\$ 145.20	+Info
123156	OUSIA - BioXen 1cc - Mineralized Xenograft .25-1mm Cancellous Particulate - #OUSM2	US\$ 161.33	+Info
123160	OUSIA - BioXen 2.5cc - Mineralized Xenograft 1-2mm Particulate - #OUSM6	US\$ 251.66	+Info
123157	OUSIA - BioXen 2.5cc - Mineralized Xenograft .25-1mm Cancellous Particulate - #O...	US\$ 297.60	+Info
123158	OUSIA - BioXen 4cc - Mineralized Xenograft .25-1mm Particulate - #OUSM4	US\$ 388.73	+Info

Bio-Graf Allograft Particulate

Code	Name	Price	+
123161	OUSIA - BioGraf .5cc - Mineralized Natural Blend Allograft .25-1mm - #OUSNB1	US\$ 70.30	+Info
123162	OUSIA - BioGraf 1cc - Mineralized Natural Blend Allograft .25-1mm - #OUSNB2	US\$ 107.56	+Info
123163	OUSIA - BioGraf 2.5cc - Mineralized Natural Blend Allograft .25-1mm - #OUSNB3	US\$ 243.10	+Info

Bio-Graft Cortical Allograft Chips

Code	Name	Price	+
468944	Mineralized Cancellous allograft- .	US\$ 124.06	+Info
468945	Mineralized Cancellous allograft- .	US\$ 221.53	+Info
468946	Mineralized Cortical allograft- .25	US\$ 76.19	+Info
468947	Mineralized Cortical allograft- .25	US\$ 124.06	+Info
468948	Mineralized Cortical allograft- .25	US\$ 232.62	+Info
468949	Demineralized Cortical allograft- 1	US\$ 49.60	+Info
468950	Demineralized Cortical allograft- 1	US\$ 76.19	+Info
468951	Demineralized Cortical allograft- 1	US\$ 124.06	+Info
468952	Demineralized Cortical allograft- 1	US\$ 324.34	+Info

Grafting Material

Resorbable Membrane



BIO ADAP BARRIER MEMBRANE- OUSIA


BIO ADAP BARRIER MEMBRANE- OUSIA

For Successful Bone Regeneration

The use of barrier membranes for guided bone regeneration, has significantly impacted the dental industry in the last 20 years. Findings early on revealed that by excluding fibrous connective tissue, a favorable environment for osteoblastic activity was created. Since that time, research has shown that the protection provided by a barrier membrane results in a highly significant increase in bone regeneration. (1)

With the growing demand for dental implant therapy, the increased number of patients needing guided bone regeneration demands that a clinician have an option for a predictable, high quality, barrier membrane for guided bone and tissue regeneration. Now you do.

Available Size: 15 mm x 20 mm- 20 mm x 30 mm- 30 mm x 40 mm



Code	Name	Price	+
123152	OUSIA - BioAdapt barrier membrane- 15x20 - 1 per box - #OUSB1	US\$ 164.58	+Info
123153	OUSIA - BioAdapt barrier membrane- 20x30 - 1 per box - #OUSB2	US\$ 197.38	+Info
123154	OUSIA - BioAdapt barrier membrane- 30x40 - 1 per box - #OUSB3	US\$ 272.11	+Info

BIO COLLAGEM- OUSIA

BIO COLLAGEM- OUSIA

Fully resorbable
 Ideal for applications in minor oral wounds
 Provides an excellent matrix to assist in wound healing
 Protects the oral wound bed
 Controls minor bleeding due to the intrinsic hemostatic properties of collagen fibers



Available Size: PLUG 1cm x 2cm- KOTE 2cm x 4cm- STRIP 2.5cm x 7.5cm

Code	Name	Price	+
123146	OUSIA - BioPlug collagen plug- resorbable 10 per box - #OUSPGS	US\$ 151.65	+Info
123148	OUSIA - BioKote collagen mat- resorbable 10 per box - #OUSKTS	US\$ 171.55	+Info
123147	OUSIA - BioStrip collagen tape- resorbable 10 per box - #OUSTP	US\$ 221.95	+Info

INTEGRA COLLAGEN PRODUCTS - MILTEX

INTEGRA COLLAGEN PRODUCTS - MILTEX

HeliMend collagen membrane is an absorbable collagen membranes made from purified bovine tendon that provides wound stabilization and creates space during guided tissue regeneration procedures where a longer absorption rate is desired.

* Absorbed in 4-8 weeks

Available: 15 X 20mm- 20 X 30mm- 30 X 40mm

Absorbable wound dressing made of collagen obtained from bovine deep flexor (Achilles) tendon. These soft, white, pliable and nonfriable sponge plugs are used to help control bleeding, develop blood clots and protect the wound site in order for the healing process to begin.

* Reabsorbed by the body in 10 to 14 days after placement

* Packaged in individual clam shells with an easy-to-peel backing

* Sterile

Available: 1 X 2cm- 2.5 X 7.5cm- 2 X 4cm



Absorbable Collagen Wound Dressings

MILTEX
MILTEX

HeliTape

HeliTape

Code	Name	Price	+
401969	Collagen Dress 2.5x7.5cm 10pk	US\$ 221.97	+Info

Other Models

Other Models

Code	Name	Price	+
401967	Collagen Dress 2x4cm 10pk	US\$ 157.95	+Info
401968	Collagen Dress 1x2cm Ster 10pk	US\$ 161.14	+Info

MILTEX-INTEGRA

MILTEX-INTEGRA

HeliCote

HeliCote

Code	Name	Price	+
445081	HeliCote Collagen Wound Dressing, 0	US\$ 171.22	+Info

HeliTape

HeliTape

Code	Name	Price	+
445069	HeliTape Collagen Wound Dressing, 1	US\$ 239.46	+Info

Other Models

Other Models

Code	Name	Price	+
445065	HeliPlug 3/8" x 3/4" Absorbable Col	US\$ 158.52	+Info

Advanced Collagen Membrane

MILTEX

MILTEX

Code	Name	Price	+
401970	Collagen Membrane 15x20mm	US\$ 132.62	+Info
401971	Collagen Membrane 20x30mm	US\$ 154.76	+Info
401972	Collagen Membrane 30x40mm	US\$ 226.79	+Info
401973	Absor Collagen Membr 15x20mm	US\$ 113.30	+Info
401974	Absor Collagen Membr 20x30mm	US\$ 139.26	+Info
401975	Absor Collagen Membr 30x40mm	US\$ 220.37	+Info

MILTEX-INTEGRA

MILTEX-INTEGRA

Code	Name	Price	+
445066	HeliMend Collagen Membrane 15 x 20	US\$ 123.26	+Info
445067	HeliMend Collagen Membrane 20 x 30	US\$ 145.87	+Info
445068	HeliMend Advanced Collagen Membrane	US\$ 243.09	+Info

445072	HeliMend Advanced Collagen Membrane	US\$ 138.91	+Info
445073	HeliMend Collagen Membrane 30 x 40	US\$ 236.15	+Info
445075	HeliMend Advanced Collagen Membrane	US\$ 166.67	+Info

Resorbable Membrane - OUSIA

Resorbable Membrane - OUSIA

Code	Name	Price	+
447598	BioMaintain Collagen Bovine Membran	US\$ 142.63	+Info
447599	BioMaintain Collagen Bovine Membran	US\$ 213.95	+Info
447600	BioMaintain Collagen Bovine Membran	US\$ 116.04	+Info

Catalog Index / Indice

Resorbable Bone Grafting Mater	1
Resorbable Membrane	2