

Cosmetic Accessories

Accessories

CURING SHIELDS - PLASDENT

HAND MIRRORS - PLASDENT



CURING HAND SHIELDS

Code	Name	Price	+
001-1860	PLASDENT Curing Shield- (4in W x 10in L)- #1860-Each	US\$ 16.99	+Info
001-1860TH	PLASDENT Tooth Curing Shield-(5.5 W x 11inL)-#1860TH-Each	US\$ 15.08	+Info

ULTRAVIOLET COVERS - PLASDENT

ULTRAVIOLET COVERS - PLASDENT



Code	Name	Price	+
001-1404-S	PLASDENT Small Ultraviolet Covers-#1404-S-Each	US\$ 10.06	+Info
001-1404-M	PLASDENT Medium Ultraviolet Covers-#1404-M-Each	US\$ 13.87	+Info

Cosmetic Accessories

Applicator Needles / Dispensing Tips

APPLICATOR NEEDLES - PLASDENT

APPLICATOR NEEDLES - PLASDENT



BENT LUER-LOCK TIPS

Code	Name	Price	+
001-BNT-100	PLASDENT Bent Luer-Lock Tips-#BLT-100-Color: BLACK-(100 Pcs/Bag)-Pre-Bent	US\$ 6.78	+Info

BENT NEEDLES TIPS

Code	Name	Price	+
------	------	-------	---

001-BNT-18	PLASDENT Bent Needle Tips-#BNT-18-18Ga-Color: PINK-(100 Pcs/Bag)-Pre-Bent	US\$ 5.23	+Info
001-BNT-20B	PLASDENT Bent Needle Tips-#BNT-20B-20ga-Color: BLACK-(100 Pcs/Bag)-Pre-Bent	US\$ 5.23	+Info
001-BNT-20Y	PLASDENT Bent Needle Tips-#BNT-20Y-20ga-Color: YELLOW-(100 Pcs/Bag)-Pre-Bent	US\$ 5.23	+Info
001-BNT-22	PLASDENT Bent Needle Tips-#BNT-22-22Ga-Color: GRAY (100 Pcs/Bag)-Pre-Bent	US\$ 5.23	+Info
001-BNT-25	PLASDENT Bent Needle Tips-#BNT-25-25ga-Color: BLUE (100 Pcs/Bag)-Pre-Bent	US\$ 5.23	+Info

Application Tips - MARK3



Liner

Infusor

Infusor

Code	Name	Price	+
425446	MARK3 Infusor Tips Bent 100pk, Pre	US\$ 39.83	+Info

Pre-Bent

Pre-Bent

Code	Name	Price	+
425448	Pre-Bent Needle Appl Tip 19ga 100pk	US\$ 6.29	+Info
425515	Pre-Bent Needle Appl Tip 18ga 100pk	US\$ 6.29	+Info
425530	Pre-Bent Needle Appl Tip 20ga 100pk	US\$ 6.29	+Info
425598	Pre-Bent Needle Appl Tip 22ga 100pk	US\$ 6.29	+Info
425707	Pre-Bent Needle Appl Tip 25ga 100pk	US\$ 6.29	+Info

Applicator Needles / Dispensing Tips - 3M ESPE

Applicator Needles / Dispensing Tips - 3M ESPE

Code	Name	Price	+
427695	Filtek Flow Dispensing Tips, Packag	US\$ 9.09	+Info

Applicator Needles / Dispensing Tips - BD

Applicator Needles / Dispensing Tips - BD

Code	Name	Price	+
425890	BD Vacutainer 25 G x 3/4" needle x	US\$ 128.12	+Info

Applicator Needles / Dispensing Tips - GC AMERICA

Applicator Needles / Dispensing Tips - GC AMERICA

Code	Name	Price	+
422356	G-aenial Dispensing Tips III - Meta	US\$ 49.96	+Info

Applicator Needles / Dispensing Tips - IVOCLAR VIVADENT

Applicator Needles / Dispensing Tips - IVOCLAR VIVADENT

Code	Name	Price	+
442649	Helioseal F Black applicator tips,	US\$ 0.00	+Info

Applicator Needles / Dispensing Tips - PAC-DENT

Applicator Needles / Dispensing Tips - PAC-DENT

Code	Name	Price	+
448859	Pac-Dent Blue 25 gauge Etchant Tips	US\$ 9.77	+Info

Applicator Needles / Dispensing Tips - PULPDENT

Applicator Needles / Dispensing Tips - PULPDENT



Multi-Cal

Multi-Cal

Code	Name	Price	+
039-DCN	Pre-bent Tips D.Blue 22ga 20pk	US\$ 19.62	+Info

Other Models

Other Models

Code	Name	Price	+
039-DCL_A	Pre-Bent Tips D.Blue 22ga 100pk	US\$ 29.77	+Info
039-ETN	Prebent Tips 25ga L.Blue 20pk	US\$ 9.40	+Info
039-ETN-50	Prebent Tips 25ga L.Blue 50pk	US\$ 31.06	+Info
457795	Pulpdent Applicator Tips 22gau	US\$ 16.69	+Info

457826	Pre-Bent applicator tips 22g 20pk	US\$ 24.66	+Info
457841	Seal-Rite Applicator Tips 15/Pk. Bl	US\$ 21.89	+Info

Applicator Needles / Dispensing Tips - TAUB

Applicator Needles / Dispensing Tips - TAUB

Code	Name	Price	+
122397	Auto Mix Tips For 7gm Syr 25pk	US\$ 39.74	+Info
122398	Fusion Zr Ven Cem Disp Tips 25	US\$ 12.72	+Info
122399	Fusion Zr Ven Cem Disp Tips 10	US\$ 41.87	+Info

Applicator Needles / Dispensing Tips - ULTRADENT

Applicator Needles / Dispensing Tips - ULTRADENT

Code	Name	Price	+
463908	Micro Tip Black 22ga 20/Pk.	US\$ 13.37	+Info

PRE BENT APPLICATOR NEEDLES - ECCO

PRE BENT APPLICATOR NEEDLES - ECCO

- * Applicator needles
- * 100pk
- * Gauges Available: 15, 18, 20, 22, 25
- * Ideal for applying chemicals and resins from 1.2 & 5cc syringes



Code	Name	Price	+
557047100	PreBent Applicator Needle Tips -ga 25ga - 100pk - YELLOW - Ideal for applying ch...	US\$ 3.48	+Info

Cosmetic Accessories

Bonding & MicroBrushes Applicators



Bonding & MicroBrushes Applicators - 3M ESPE

Bonding & MicroBrushes Applicators - 3M ESPE

Code	Name	Price	+
427191	3M ESPE Bendable Applicator Brushes	US\$ 41.64	+Info
427270	3M Mini-Sponge Applicator 200pk	US\$ 42.58	+Info
427328	3M ESPE Applicator Brush Set: 3 Ass	US\$ 61.10	+Info

Bonding & MicroBrushes Applicators - COLTENE WHALADENT

Bonding & MicroBrushes Applicators - COLTENE WHALADENT

Code	Name	Price	+
426900	One Coat self-etch bond microbrushe	US\$ 19.43	+Info

433597	One Coat self-etch bond One Coat Se	US\$ 0.00	+Info
--------	-------------------------------------	-----------	-----------------------

Bonding & MicroBrushes Applicators - ECCO

Bonding & MicroBrushes Applicators - ECCO

Bendable Brush

Bendable Brush

Code	Name	Price	+
112434	Applicator Brush Pink 100pk	US\$ 4.24	+Info

Handles

Handles

Code	Name	Price	+
557612	Composite Brush Handle 5pk	US\$ 0.93	+Info

Round Cylinder

Round Cylinder

Code	Name	Price	+
115950	Micro brush Round Cylinder 3mm	US\$ 3.75	+Info

Bonding & MicroBrushes Applicators - HOUSE BRAND

Bonding & MicroBrushes Applicators - HOUSE BRAND

Code	Name	Price	+
470704	Bonding Brush Tips 100Pk	US\$ 3.38	+Info

Bonding & MicroBrushes Applicators - KEYSTONE

Bonding & MicroBrushes Applicators - KEYSTONE

Bonding Brushes

Bonding Brushes

Code	Name	Price	+
445800	PIP Mizzy Pressure Indicator Paste	US\$ 11.72	+Info
445893	Bosworth Applicator Brushes, Packag	US\$ 17.33	+Info

Micro-Applicators

Micro-Applicators

Code	Name	Price	+
446224	Small Micro-applic fine tip-400pk	US\$ 15.66	+Info
446225	XSmall Micro-applic fine tip-400pk	US\$ 15.66	+Info
446305	Medium Micro-applicat PINK 400pk	US\$ 15.66	+Info

Bonding & MicroBrushes Applicators - MARK3

Bonding & MicroBrushes Applicators - MARK3

Code	Name	Price	+
425787	MARK3 Cargus Regular Brush Tips, Pa	US\$ 8.30	+Info

Bonding & MicroBrushes Applicators - MICROBRUSH

Bonding & MicroBrushes Applicators - MICROBRUSH

MicroBrush Applicators refills

MicroBrush Applicators refills

Fine 1.5mm

Fine 1.5mm

Code	Name	Price	+
400756	Microbrush App Fine Asstd 400pk	US\$ 36.38	+Info
400757	Microbrush Appl Fine Yellow 400pk	US\$ 36.38	+Info
401031	Appli Refill Fine Pink 400pk	US\$ 56.03	+Info
401032	Appl Refill Fine Pink/Yell 400pk	US\$ 56.03	+Info
401033	Appl Refill Fine Yellow 400pk	US\$ 56.03	+Info
401048	Microbrush Fine Pink 400pk tube	US\$ 56.03	+Info

Regular 2.0mm

Regular 2.0mm

Code	Name	Price	+
400758	Microbrush Appli Reg Asstd 400pk	US\$ 52.04	+Info
400759	Microbrush Appli Reg Green 400pk	US\$ 36.38	+Info
400760	Microbrush Appli Reg Purple 400pk	US\$ 36.38	+Info
401034	Appl Refill Regular Asstd 400pk	US\$ 56.03	+Info
401035	Appl Refill Regular Blue 400pk	US\$ 56.03	+Info
401036	Appl Refill Regular Green 400pk	US\$ 56.03	+Info
401037	Appl Refill Regular Peach 400pk	US\$ 56.03	+Info
401038	Appl Refill Regular Purple 400pk	US\$ 56.03	+Info
401052	Microbrush Reg Blue 400pk	US\$ 56.03	+Info
401054	Microbrush Reg Peach 400pk	US\$ 56.03	+Info

Super Fine 1.0mm

Super Fine 1.0mm

Code	Name	Price	+
400761	Microbrush Appli SFine Wte 400pk	US\$ 36.38	+Info
401039	Refill Supfine White 400pk	US\$ 56.03	+Info
401584	Microbr XThin Black 100pk	US\$ 28.23	+Info

MicroBrush Dispenser

MicroBrush Dispenser

Fine 1.5mm

Fine 1.5mm

Code	Name	Price	+
401027	Dispenser KIT Fine Pink/Yell 400pk	US\$ 76.08	+Info
401028	Dispenser Kit Fine Yellow 50pk	US\$ 15.54	+Info

Regular 2.0mm

Regular 2.0mm

Code	Name	Price	+
401029	Dispenser Kit Reg Asstd 400pk	US\$ 76.08	+Info

Super Fine 1.0mm

Super Fine 1.0mm

Code	Name	Price	+
401030	Dispenser Kit Sup fine White 400pk	US\$ 76.08	+Info
401583	Microbr Disp KIT XThin Black 100pk	US\$ 39.80	+Info

MicroBrush Plus

MicroBrush Plus

Fine 1.5mm

Fine 1.5mm

Code	Name	Price	+
400762	Microbr Plus Appli Fine Asstd 400pk	US\$ 36.38	+Info

Regular 2.0mm

Regular 2.0mm

Code	Name	Price	+
400763	Microbr Plus Appl Reg Asstd 400pk	US\$ 36.38	+Info

Super Fine 1.0mm

Super Fine 1.0mm

Code	Name	Price	+
400764	Microbr Plus Appli SFine White 400pk	US\$ 36.38	+Info

UltraBrush Dispenser

UltraBrush Dispenser

Fine 1.5mm

Fine 1.5mm

Code	Name	Price	+
401579	Ultrabrush Applic Fine Bl/Or 200pk	US\$ 45.98	+Info
401580	Ultrabrush Fine Blue w/disp 100pk	US\$ 37.51	+Info

Regular 2.0mm

Regular 2.0mm

Code	Name	Price	+
------	------	-------	---

401581	Ultrabrush Appl Reg Yel/Gr 200pk	US\$ 45.98	+Info
401582	Ultrabrush Reg Green w/disp 100pk	US\$ 37.51	+Info

Bonding & MicroBrushes Applicators - PAC-DENT

Bonding & MicroBrushes Applicators - PAC-DENT

Code	Name	Price	+
448908	Pac-Dent Disposable all-purpose Mic	US\$ 28.30	+Info
448956	Pac-Dent Disposable all-purpose Mic	US\$ 21.20	+Info
448981	Pac-Dent Disposable all-purpose Mic	US\$ 21.20	+Info
449005	Pac-Dent Disposable all-purpose Mic	US\$ 21.20	+Info

Bonding & MicroBrushes Applicator s - PLASDENT

Bonding & MicroBrushes Applicators - PLASDENT



Bonding Brushes

Bonding Brushes

Code	Name	Price	+
001-8404BRH	Bonding Brush Tips 100Pk Ass	US\$ 5.19	+Info
001-8404KIT	Brsh & Micro Tps Mx Kt 50tp5hn	US\$ 5.50	+Info

Brush Applicators

Brush Applicators

Code	Name	Price	+
423542	Plasdent Disposable Brush Applicato	US\$ 15.28	+Info
423729	Plasdent Disposable Brush Applicato	US\$ 15.28	+Info
424092	Plasdent Disposable Brush Applicato	US\$ 15.28	+Info
424166	Plasdent Disposable Brush Applicato	US\$ 15.28	+Info

Handles

Handles

Code	Name	Price	+
001-8404HND	Univ. Brush Handle - WHITE - Ea	US\$ 0.80	+Info
001-8404HND-B	Univ. Brush Handle BLUE - Ea	US\$ 0.71	+Info
001-8404HND-G	Univ. Brush Handle GREEN - Ea	US\$ 0.48	+Info
001-8404HND-R	Univ. Brush Handle red - Ea	US\$ 0.48	+Info
001-8404HND-Y	Univ. Brush Handle YELLOW - Ea	US\$ 0.48	+Info

MAXBRUSH

MAXBRUSH

Code	Name	Price	+
001-8101ST-A	Maxbrush App 4In 100pk - Ass	US\$ 5.25	+Info
001-8101ST-R	Maxbrush App 4In 100pk red	US\$ 5.25	+Info

MaxMicro

MaxMicro

Extended

Extended

Code	Name	Price	+
001-600-EX-11	Maxext Micro App - Ext - Blck 100p	US\$ 5.42	+Info

Fine 1.5mm

Fine 1.5mm

Code	Name	Price	+
001-600-F	MaxMicro App 1.5Mm Fn Yll 100p	US\$ 4.81	+Info

Regular 2.0mm

Regular 2.0mm

Code	Name	Price	+
001-600-R	MaxMicro App - 2.0Mm Rgl Pu 100p	US\$ 4.63	+Info
001-600-R-12P	MaxMicro App - 2.0Mm Rgl Pch 100	US\$ 4.63	+Info
001-600-R-2	MaxMicro App - 2.0Mm - Blue - 100pk	US\$ 4.81	+Info
001-600-R-4	MaxMicro App - 2.0Mm - Green - 100pk	US\$ 4.63	+Info
600-R-4	MaxMicro App - 2.0Mm Rgl Pu 400p	US\$ 17.34	+Info

Super Fine 1.0mm

Super Fine 1.0mm

Code	Name	Price	+
001-600SF-1	MaxMicro App - 1.0Mm SF Wh 100p	US\$ 4.81	+Info

Ultra Fine 0.5mm

Ultra Fine 0.5mm

Code	Name	Price	+
001-600-UF-11	MaxMicro App - 0.5Mm UF Blck 100	US\$ 4.81	+Info

Micro App Bendable Tip

Micro App Bendable Tip

Code	Name	Price	+
423828	Plasdent Micro Applicators - Blue 1	US\$ 7.91	+Info
424121	Plasdent Micro Applicators - Black	US\$ 7.91	+Info
424122	Plasdent Micro Applicators - red 14	US\$ 7.91	+Info
424163	Plasdent Micro Applicators - Purple	US\$ 7.91	+Info
424203	Plasdent Micro Applicators - Yellow	US\$ 7.91	+Info

424216	Plasdent Micro Applicators - Green	US\$ 7.91	+Info
Micro Applicators Micro Applicators			
Fine 1.5mm Fine 1.5mm			
Code	Name	Price	+
423939	Fine Tip Micro Applicat PINK 100pk	US\$ 3.92	+Info
Super Fine 1.0mm Super Fine 1.0mm			
Code	Name	Price	+
423567	MaxMicro Tips 1.0mm SF Wte 100pk	US\$ 4.81	+Info
MICRO-TIPS MICRO-TIPS			
Code	Name	Price	+
001-8404TIP	Micro Tips - 100 pk	US\$ 4.28	+Info

Bonding & MicroBrushes Applicators - TPC

Bonding & MicroBrushes Applicators - TPC



Micro Applicators

Micro Applicators

Code	Name	Price	+
464083	TPC Advanced Technology Disposable	US\$ 14.00	+Info
464085	TPC Advanced Technology Disposable	US\$ 14.00	+Info
464090	Micro Applicators Fine 4x100pk	US\$ 14.00	+Info
464097	TPC Advanced Technology Disposable	US\$ 14.00	+Info

Other Models

Other Models

Code	Name	Price	+
464120	TPC Advanced Technology Disposable	US\$ 0.00	+Info

Cosmetic Accessories

Cheek Retractor

Cheek Retractor - HAGER

Cheek Retractor - HAGER



Mirahold

Mirahold

Code	Name	Price	+
434981	Open-Ex Silicone mouth props 12pk	US\$ 84.52	+Info

Mirasuto HVE Tips

Mirasuto HVE Tips

Code	Name	Price	+
434884	Hager Worldwide Aspirator Tips with	US\$ 68.33	+Info

Cheek Retractor - OSUNG

Cheek Retractor - OSUNG

Code	Name	Price	+
449955	Columbia Cheek & Lip Retractor	US\$ 18.60	+Info
449981	Osung Lip Retractor, CRL, Used for	US\$ 55.80	+Info
450010	Osung Plastic lip widener, Small, u	US\$ 8.80	+Info
450046	Osung Plastic lip widener, Large, u	US\$ 8.80	+Info
450275	Osung Lip and cheek retractor with	US\$ 10.85	+Info

Cheek Retractor - PLASDENT

Cheek Retractor - PLASDENT

Code	Name	Price	+
423490	EXTND Double Ended Cheek Retractors	US\$ 7.19	+Info
423498	EXTND Cheek Retractor - Hand-held A	US\$ 6.66	+Info
423506	EXTND Cheek Retractor - Hand-free A	US\$ 7.19	+Info
423523	EXTND Cheek Retractor Hand-held 2pk	US\$ 6.66	+Info
423546	EXTND Cheek Retractor - Hand-free A	US\$ 6.54	+Info
423629	EXTND Double Ended Cheek Retractors	US\$ 7.19	+Info
423639	EXTND Cheek Retractor - Hand-free C	US\$ 7.19	+Info
423669	EXTND Disposable Intraoral Cheek Re	US\$ 7.39	+Info
423761	EXTND Intraoral Cheek Retractors -	US\$ 7.45	+Info
423781	EXTND Cheek Retractor - Hand-held A	US\$ 7.19	+Info
423807	EXTND Wire Type Cheek Retractor, CI	US\$ 7.72	+Info
423813	EXTND Intraoral Cheek Retractors -	US\$ 7.45	+Info
423814	EXTND Cheek Retractor - Hand-held C	US\$ 7.19	+Info
423822	EXTND Cheek Retractor - Hand-free C	US\$ 7.19	+Info

423874	EXTND Intraoral Cheek Retractors -	US\$ 7.45	+Info
423877	EXTND Disposable Intraoral Cheek Re	US\$ 6.80	+Info

Cheek Retractor - SHARK

Cheek Retractor - SHARK

Code	Name	Price	+
917003	Cheek Retr Medi w/Tong Ret 10pk	US\$ 3.52	+Info

CHEEK RETRACTOR - SHARK

CHEEK RETRACTOR - SHARK

Features:

- * Easy to insert
- * One size fits all
- * Convenient tongue guard
- * No Sharp Edges
- * Transparent plastic
- * Oval shape retractor
- * Allows comfortable tongue placement



Code	Name	Price	+
917001	SHARK - 2x Cheek Retractor - Adult Size - Disposable -10pk 5x2pcs - #748002 A	US\$ 4.00	+Info
917002	SHARK - 2x Cheek Retractor - Child Size - Disposable - 10pk 5x2pcs - #748002 C	US\$ 4.00	+Info

Cosmetic Accessories

Composite Gun Dispensers

Composite Gun Dispensers - 3M ESPE

Composite Gun Dispensers - 3M ESPE

Code	Name	Price	+
427182	3M ESPE Restorative Caps Dispenser.	US\$ 0.00	+Info
427831	Aplicap Applier	US\$ 85.24	+Info
427892	Aplicap Activator	US\$ 76.69	+Info
428071	Filtek Restoratives Dispenser Gun	US\$ 62.34	+Info

Composite Gun Dispensers - CENTRIX

Composite Gun Dispensers - CENTRIX

Code	Name	Price	+
109-130600	Centrix Pistola Plastica Plast	US\$ 93.86	+Info

Composite Gun Dispensers - MARK3

Composite Gun Dispensers - MARK3

Code	Name	Price	+
425461	MARK3 Composite Gun - Compules Disp	US\$ 25.58	+Info

Composite Gun Dispensers - PAC-DENT

Composite Gun Dispensers - PAC-DENT

Code	Name	Price	+
448852	Pac-Dent Plastic composite dispensi	US\$ 28.23	+Info
448912	Pac-Dent Composite Dispensing Gun w	US\$ 65.51	+Info

Composite Gun Dispensers - PLASDENT

Composite Gun Dispensers - PLASDENT

ACUPUSH

ACUPUSH

Code	Name	Price	+
113737	Composite Gun Light Gray	US\$ 11.03	+Info
423626	AcuPush Unidose Dispens Gun Lt Gray	US\$ 12.02	+Info
423633	Dispenser Gun Unidose Black	US\$ 11.61	+Info
423638	AcuPush Unidose Dispens Gun Blue Oc	US\$ 11.61	+Info
423770	AcuPush Unidose Dispens Gun Orange	US\$ 11.61	+Info
423776	Dispenser Gun Unidose Green	US\$ 11.61	+Info
423852	Dispenser Gun Unidose Purple Gala	US\$ 11.49	+Info

GUN-CAPS

GUN-CAPS

Code	Name	Price	+
423623	Capsule Applier Universal	US\$ 34.60	+Info

Cosmetic Accessories

Composite Heaters

Composite Heaters - COLTENE

Composite Heaters - COLTENE

Code	Name	Price	+
430790	Ultrafil 3D Introductory Kit with H	US\$ 0.00	+Info

Composite Heaters - PEGASUS

Composite Heaters - PEGASUS

Code	Name	Price	+
------	------	-------	---

466398	Anesthetic & Composite warmer	US\$ 133.13	+Info
--------	-------------------------------	-------------	-----------------------

Composite Heaters - PREMIER

Composite Heaters - PREMIER

Code	Name	Price	+
452983	Cartridge Warmer Complete Unit	US\$ 79.89	+Info

Cosmetic Accessories

Lip Expanders



Lip Expanders - IVOCLAR VIVADENT

Cheek Retractor - IVOCLAR VIVADENT



Code	Name	Price	+
442202	OptraGate 3D Retractors - Adult Reg	US\$ 159.45	+Info
442208	OptraGate 3D OptraGate Retractors -	US\$ 159.30	+Info
442225	OptraGate 3D OptraGate Retractors -	US\$ 159.30	+Info
442304	OptraGate 3D OptraGate Retractors -	US\$ 164.59	+Info
442447	OptraGate 3D OptraGate Retractors -	US\$ 112.50	+Info
442503	OptraGate 3D OptraGate Retractors -	US\$ 112.50	+Info

Lip Expanders - EASY-SMILE

Lip Expanders - EASY-SMILE

Code	Name	Price	+
466266	EASYSMILE Mouth Opener	US\$ 7.72	+Info

LIP EXPANDERS - ECCO

LIP EXPANDERS - ECCO

- * To hold lips and cheeks away during etching and bonding procedures
- * Autoclavable plastic Soft/Har
- * Comfortable for patient
- * 1pk
- * Autoclave cycle 121C x 10 minutes



ANTERIOR

- * Lip expander
- * One size fits all
- * Used for bonding
- * 1pk

Code	Name	Price	+
822350	Anterior Lip Expander - BLUE - each One size fits all - One Size fits all - Ster...	US\$ 4.00	+Info
822366	Anterior Lip Expander - CLEAR - each One size fits all - One Size fits all - Ste...	US\$ 4.00	+Info

EXPANDER CLEAR

Code	Name	Price	+
822666	Autoclavable Lip Expander CLEAR Mini size (Child - Pedo) - Each - #X-130 CRO2 ...	US\$ 3.88	+Info
822144	Autoclavable Lip Expander CLEAR Medium Size each - #JG13#M	US\$ 4.00	+Info
822444	Autoclavable Lip Expander CLEAR Large (Adult Standard size) CLEAR - Each #X-13...	US\$ 4.00	+Info

Lip Expanders - IVOCLAR VIVADENT

Lip Expanders - IVOCLAR VIVADENT

Code	Name	Price	+
442522	OptraGate 3D OptraGate Retractors -	US\$ 0.00	+Info
442533	OptraGate 3D OptraGate Retractors -	US\$ 0.00	+Info

Lip Expanders - OSUNG

Lip Expanders - OSUNG



Code	Name	Price	+
449965	Osung Plastic lip widener, Medium,	US\$ 10.00	+Info
450005	Osung Plastic lip widener, Small, u	US\$ 10.00	+Info
450006	Osung Plastic lip widener, Large, u	US\$ 10.00	+Info

Lip Expanders - PREMIER

Cheek Retractor - PREMIER



Code	Name	Price	+
------	------	-------	---

453046	Comfortview Small Size - Lip and Ch	US\$ 71.36	+Info
453529	Comfortview Cushion 10pk	US\$ 27.82	+Info

Cosmetic Accessories

Mixing Wells



BONDWELL-II DISPOSABLE MIXING WELLS - PLASDENT

BONDWELL-II DISPOSABLE MIXING WELLS - PLASDENT



2 WELL

Code	Name	Price	+
001-901MW-2	Bndwll-II Dspb Mxng White-2w 200pk	US\$ 7.90	+Info
423478	Bondwell 2-Well Disposab Wte 200pk	US\$ 14.55	+Info
423504	Bondwell 2-Well Disposable Mixing W	US\$ 34.50	+Info

4 WELL

Code	Name	Price	+
001-901MW-4	Bndwll - II Dspb Mxng Wh -	US\$ 8.96	+Info
423668	Bondwell 4-Well Disposable Mixing W	US\$ 16.75	+Info
423670	Bondwell 4-Well Disposable Mixing W	US\$ 36.54	+Info
423974	Bondwell 4-Well Disposable Mixing W	US\$ 16.31	+Info

Mixing Wells - 3D-DENTAL

Mixing Wells - 3D-DENTAL

Code	Name	Price	+
474033	3D Dental Disposable Mixing Well 2-	US\$ 15.20	+Info
474034	3D Dental Disposable Mixing Well 4-	US\$ 15.75	+Info

Mixing Wells - 3M ESPE

Mixing Wells - 3M ESPE

Code	Name	Price	+
427235	3M Mixing Wells Disposable 96pk	US\$ 30.65	+Info
427346	Clinpro Sealant Mixing Well 4pk	US\$ 44.72	+Info

Mixing Wells - IVOCAR VIVADENT

Mixing Wells - IVOCAR VIVADENT

Code	Name	Price	+
------	------	-------	---

442316	Vivapad Plastic Mixing Well w/Cover	US\$ 23.32	+Info
--------	-------------------------------------	------------	-----------------------

Mixing Wells - MARK3

Mixing Wells - MARK3

Code	Name	Price	+
425570	MARK3 Cargus 2-Well Disposable Mixi	US\$ 18.23	+Info

Mixing Wells - PLASDENT

Mixing Wells - PLASDENT

Code	Name	Price	+
423832	AcuWell Universal Material Well - B	US\$ 5.27	+Info
423959	AcuWell Univer Mater Well Lt Gray	US\$ 5.42	+Info
424002	AcuWell Universal Material Well - B	US\$ 5.47	+Info
424003	AcuWell Universal Material Well - P	US\$ 5.47	+Info
424004	AcuWell Universal Material Well - O	US\$ 5.47	+Info
424299	AcuWell Univers Material Well Green	US\$ 5.47	+Info

Mixing Wells - PULPDENT

Mixing Wells - PULPDENT

Code	Name	Price	+
039-MW-2	Disp Mixing Wells 2-Well 480pk	US\$ 95.88	+Info
039-MW-4	Disp Mixing Well 4-Well 420pk	US\$ 95.88	+Info
457701	Pulpdent 4-Well Disposable Mixing W	US\$ 75.24	+Info
457708	Pulpdent 2-Well Disposable Mixing W	US\$ 78.22	+Info
457787	Pac-Dent 2-Well Disposable Mixing W	US\$ 33.74	+Info

Mixing Wells - SHARK

Mixing Wells - SHARK

Code	Name	Price	+
115838	Mixing Well 4 Cavity	US\$ 8.73	+Info

Cosmetic Accessories

Mouth Props

Mouth Props - DENTSPLY

Mouth Props - DENTSPLY

Code	Name	Price	+
469189	EZ-Prop Mouse Props 100/Bx. Disposa	US\$ 0.00	+Info

Mouth Props - PLASDENT

Mouth Props - PLASDENT

Extend

Extend

Code	Name	Price	+
423562	Disposab Mouth Props Small Wte 48pk	US\$ 8.57	+Info
423641	Disposable Mouth Props MED 48pk	US\$ 8.04	+Info
423647	Disposab Mouth Props Large Wte 48pk	US\$ 8.04	+Info

Mouth Props

Mouth Props

Rubber

Rubber

Code	Name	Price	+
423527	EXTND Rubber Mouth Props - Medium (	US\$ 7.31	+Info
423563	EXTND Rubber Mouth Props - Small (P	US\$ 7.31	+Info
423933	EXTND Rubber Mouth Props Large 2pk	US\$ 7.31	+Info

Silicone

Silicone

Code	Name	Price	+
423480	EXTND Silicone Mouth Props - Medium	US\$ 8.71	+Info
423501	EXTND Silicone Mouth Props - Large	US\$ 8.71	+Info
423508	EXTND Silicone Mouth Props - Small	US\$ 8.71	+Info
423645	EXTND Silicone Mouth Props - Large	US\$ 8.71	+Info
423679	EXTND Silicone Mouth Props - Medium	US\$ 8.71	+Info
423848	EXTND Silicone Mouth Props - Small	US\$ 8.71	+Info

Mouth Props - SHARK

Mouth Props - SHARK

Code	Name	Price	+
113698	Mouth Props Medium S 3pk	US\$ 3.96	+Info

Cosmetic Accessories

Preformed Angle Formers

Preformed Angle Formers - 3M ESPE

Preformed Angle Formers - 3M ESPE

Code	Name	Price	+
427192	3M ESPE #1 Upper Left Central Pedia	US\$ 32.50	+Info

427239	3M ESPE #1 Upper Right Lateral Pedi	US\$ 32.50	+Info
427236	3M ESPE #1 Upper Right Central Pedi	US\$ 32.50	+Info
427229	3M ESPE #1 Upper Left Lateral Pedia	US\$ 32.50	+Info
427222	3M ESPE #4 Upper Left Lateral Pedia	US\$ 32.50	+Info
427208	3M ESPE #4 Upper Left Central Pedia	US\$ 32.50	+Info
427204	3M ESPE #4 Upper Right Lateral Pedi	US\$ 32.50	+Info
427203	3M ESPE #2 Upper Right Lateral Pedi	US\$ 32.50	+Info
427197	3M ESPE #4 Upper Right Central Pedi	US\$ 32.50	+Info
427144	3M ESPE #3 Upper Right Lateral Pedi	US\$ 32.50	+Info
427189	3M ESPE #3 Upper Left Central Pedia	US\$ 32.50	+Info
427188	3M ESPE #2 Upper Left Lateral Pedia	US\$ 32.50	+Info
427175	3M ESPE #2 Upper Left Central Pedia	US\$ 32.50	+Info
427174	3M ESPE #2 Upper Right Central Pedi	US\$ 32.50	+Info
427171	3M ESPE Intro Kit - Pediatric Strip	US\$ 761.18	+Info
427148	3M ESPE #3 Upper Right Central Pedi	US\$ 32.50	+Info
427146	3M ESPE #3 Upper Left Lateral Pedia	US\$ 32.50	+Info

COFORM - DIRECTA

COFORM - DIRECTA

CoForm, matrix system, is a unique set of pre-formed transparent matrices, made of celluloid plastic. There are specifically designed to deal with composite Class IV restorations as well as various tooth fractures.

The method and advantages of celluloid plastic matrices.

The cavity is filled in layers in the usual manner. The only difference is that when the last layer is placed, instead of shaping and curing, the matrix is pressed onto the uncured composite giving the restoration the anatomy of a healthy tooth. The composite is then cured through the clear matrix.



ASSORTED KIT

Code	Name	Price	+
121243	DIRECTA - CoForm 128-Piece Matrix Kit: 2 each of 4 sizes for 8 Mesial Corners an...	US\$ 221.00	+Info
121242	DIRECTA - CoForm 64-Piece Matrix Kit: Sizes 11, 12, 21, 22, Mesial & Distal - #6...	US\$ 112.11	+Info

DISTAL

Code	Name	Price	+
121245	DIRECTA - CoForm Matrix Refill #11 Distal 8/Pk. Innovative aid for forming resto...	US\$ 18.05	+Info
121247	DIRECTA - CoForm Matrix Refill #21 Distal 8/Pk. Innovative aid for forming resto...	US\$ 18.05	+Info

Preformed Angle Formers - DIRECTA

Preformed Angle Formers - DIRECTA

Code	Name	Price	+
------	------	-------	---

417386	Full Form Strip off crown 60pc kit	US\$ 106.01	+Info
--------	------------------------------------	-------------	-----------------------

Preformed Angle Formers - HENRY SCHEIN

Preformed Angle Formers - HENRY SCHEIN

Code	Name	Price	+
420275	Strip Crown Asstd Anterior 60pk	US\$ 58.81	+Info
420276	Strip Crown Asstd Posterior 60pk	US\$ 51.60	+Info

Cosmetic Accessories

Veneers Carriers



Sticky Tips

Sticky Tips



MicroBrush

MicroBrush

Code	Name	Price	+
401577	Micro Tip Appl Blue 64pk	US\$ 43.66	+Info

PLASDENT

PLASDENT

Code	Name	Price	+
001-8303-1	Sticky Applicator White - 60 pk	US\$ 15.87	+Info

PULPDENT

PULPDENT

Code	Name	Price	+
039-PIE	Pic-N-Stic Stiky Applicator 60pk	US\$ 29.48	+Info

Veneers Carriers - MPF BRUSH

Veneers Carriers - MPF BRUSH

Code	Name	Price	+
------	------	-------	---

123937

Crown & Veneer Holder KIT

US\$ 35.14

[+Info](#)

VENEER ME - SMILE LINE

SMILE LINE - VENEER ME

VENEER ME: MAXIMUM SECURITY FOR YOUR VENEERS

Thanks to veneerME, you'll be able to carry and ship and also manipulate your hyper-delicate veneers in the highest conditions of security !

The sticky polyurethane gel-pad

The base of veneerME box is covered with a sticky pad out of polyurethane gel. Veneers are placed on this pad on their buccal surface, and adhere naturally. All the preparation process of the internal part of the veneer (acid etching, silanization) can be done in this position without any fear of a veneer escaping or falling on the floor. The pad's material is super resistant and can be re-used several times (it is easily replaceable in case of deterioration).

The security membrane-frame

Once the veneers are secured on the pad, the membrane-frame is clipped on the base of veneerME box for extra fixation. The xx-thin membrane, made out of nylon, delicately adapts to the shape of the veneers. This done, you just need to close the lid and your master work is ready for delivery.

The membrane-frame can be replaced -if needed- after having been used several times.

Huge advantage : the membrane is permeable ! In other words the clinician can proceed with rinsing and drying after having replaced the membrane-frame on the base. This way the clinician can work in the best conditions and in total security.



Code	Name	Price	+
420194	SMILELINE - VeneerME set of 3 complete boxes 3 spare PU gel pads and 3 spare mem...	US\$ 284.00	+Info
420195	SMILELINE - VeneerME Mini-set of 1 complete box 1 spare PU gel pad and 1 spare m...	US\$ 118.50	+Info
473362	SMILE LINE - Tinto Seal for Transportation (10 pcs) - # 1300-SEAL	US\$ 16.17	+Info
420196	SMILELINE -VeneerME spare PU pads (10 pcs) #2610	US\$ 73.00	+Info
420197	SMILELINE - VeneerME spare membrane-frames (10 pcs) #2620	US\$ 124.96	+Info

Cosmetic Accessories

Wedges



ACUWEDGES PLASTIC ASSORTED - PLASDENT

ACUWEDGES PLASTIC ASSORTED - PLASDENT



Code	Name	Price	+
001-WG-A	PLASDENT Acuwedges Plastic Wedges- #WG-A- ASSORTED 4 Sizes (100 Pcs/Box)	US\$ 6.46	+Info
001-WG-16	PLASDENT Acuwedges Plastic Wedges-#WG-16-16Mm-Color: PURPLE (100 Pcs/Box)	US\$ 6.03	+Info
001-WG-14	PLASDENT Acuwedges Plastic Wedges-#WG-14-14Mm-Color: YELLOW (100 Pcs/Box)	US\$ 6.46	+Info
001-WG-12	PLASDENT Acuwedges Plastic Wedges-#WG-12-12Mm-Color: GREEN (100 Pcs/Box)	US\$ 6.43	+Info
001-WG-10	PLASDENT Acuwedges Plastic Wedges- #WG-10- 10mm- Color: BLUE (100 Pcs/Box)	US\$ 6.46	+Info

BARMANS ANATOMICAL WEDGES - DIRECTA

BARMAN'S ANATOMICAL WEDGES - DIRECTA

The barman wedges achieve safer procedures thanks to its anatomical shape. They are available in two sizes. small (natural) and large (green). The sides of the wedge are shaped according to the interdental space



Code	Name	Price	+
121221	DIRECTA - Directa Barman's Wedges Anatomical, Mixed 200/Pk. Progressively thic...	US\$ 32.15	+Info

FENDER PRIME - DIRECTA

FENDER PRIME - DIRECTA

FenderPrime is an interproximal plastic wedge with an inbuilt steel plate. The "boat-shaped" tip allows the gingiva to be compressed without damage and reduces the risk for bleeding. FenderPrime Green (long) can be shaped to the contour of the tooth.

FenderPrime is available in two sizes, long and short. FenderPrime green (the long wedge) is bendable. FenderPrime yellow has a shorter, more rigid plate for very tight interproximal spaces.



KITS

Code	Name	Price	+
121235	DIRECTA - Fenderfamily Complete Assorted Kit - Kit - 6028-10 - #602810	US\$ 435.63	+Info

121237	DIRECTA - FenderMate Prime Assortment Kit 36/Pk. Includes 18 each of short and l...	US\$ 80.65	+Info
--------	---	------------	-----------------------

VALUE PACK

Code	Name	Price	+
121232	DIRECTA - FenderMate Prime Short, NEON YELLOW 100/Value Pack. Unique combination...	US\$ 142.38	+Info
121233	DIRECTA - FenderMate Prime Long, GREEN 100/Value Pack. Unique combination - #602809	US\$ 172.96	+Info

WEDGES ONLY

Code	Name	Price	+
121238	DIRECTA - FenderMate Prime Short, YELLOW 18/Pk. Unique combination of a preparat...	US\$ 48.43	+Info

FENDER WEDGE - DIRECTA

FENDER WEDGE - DIRECTA

FenderWedge is a powerful combination of a wedge and a protective stainless steel plate thin enough to accommodate most narrow spaces easily. When inserted, FenderWedge pre-separates the teeth by a few tenths of a millimeter, which aids in the final building of the contact point.



- * Protects adjacent teeth during Class II preparation
- * Facilitates easy placement of matrix and establishment of contact points
- * As easy to insert as a wedge
- * Jams in and stays in position

KIT

Code	Name	Price	+
121210	DIRECTA - FenderWedge Assortment Kit 144pk 36 wedges of each size. Designed to ...	US\$ 163.33	+Info

SMALL

Code	Name	Price	+
121211	DIRECTA - FenderWedge Small 1.4mm, ORANGE 36pk Designed to protect adjacent tee...	US\$ 57.73	+Info

FENDER-PREP - DIRECTA

FENDER-PREP - DIRECTA

FenderPrep is a plastic wedge with an inbuilt steel plate specially designed for larger preparations such as crowns, veneers and inlays. FenderPrep improves the quality of preparations and saves time thanks to increased vision and safety.

FenderPrep has a thin, "boat-shaped" tip which enables it to be inserted into tight proximal spaces at the same time as it compresses the gingiva without damage, avoiding bleeding and delay of impression taking.

The steel plate is perforated thus allowing it to be bent before insertion to lean against the adjacent tooth, providing better access and vision.



Code	Name	Price	+
121236	DIRECTA - FenderWedge Prep FenderWedge Prep, GRAY - Value Pack - 100 wedges. is ...	US\$ 236.04	+Info

FENDERMATE - DIRECTA

FENDERMATE - DIRECTA

FenderMate the fastest sectional matrix in the world. Establish tight marginal seal and create contact point in seconds

A matrix that does not perfectly adapt to the cavity margin under the contact point may cause overhang and a control examination with a probe or floss may not detect this. Over a period of time occlusal pressure causes fracturing of unbonded excess material, wich creates a trap for food impaction and plaque retention. Sectional matrix systems consisting of a matrix, wedge and ring, may create a risk of leakage due to lighter pressure of the wedge against the matrix when a retention ring is applied to separate the teeth.



FENDERMATE KIT

Code	Name	Price	+
121215	DIRECTA - FenderMate Assortment Kit: One hand operated organizer box with 72 pie...	US\$ 273.48	+Info

FENDERMATE VALUE PACK

Code	Name	Price	+
121234	DIRECTA - FenderMate Temp FenderMate-Temp, Lilac Value Pack - 100 Wedges.Stainle...	US\$ 223.44	+Info
121228	DIRECTA - FenderMate Left Regular, Dark GREEN 100/Value Pack. Combines a pre-cur...	US\$ 349.42	+Info
121229	DIRECTA - FenderMate Left Narrow, Light GREEN 100/Value Pack. Combines a pre-cur...	US\$ 349.42	+Info
121230	DIRECTA - FenderMate Right Regular, Dark BLUE 100/Value Pack. Combines a pre-cur...	US\$ 349.42	+Info
121231	DIRECTA - FenderMate Right Narrow, LIGHT BLUE 100/Value Pack. Combines a pre-cur...	US\$ 349.42	+Info

Wedges - BOSWORTH

Wedges - BOSWORTH

COLORED

COLORED

Code	Name	Price	+
445849	Wonder Wedges 12mm Small Blue, 100/	US\$ 16.33	+Info
446105	Wonder Wedges Assorted, Colored, 10	US\$ 18.15	+Info
446163	Wonder Wedges 13mm Medium Green	US\$ 19.97	+Info
446205	Wonder Wedges 15mm Large Yellow, 10	US\$ 18.15	+Info
446604	Wonder Wedges 17mm XL Bleached	US\$ 19.97	+Info

NATURAL

NATURAL

Code	Name	Price	+
123385	WONDER WEDGES Box Assorted Nat	US\$ 18.15	+Info
445864	Wonder Wedges 12mm Small, Natural.	US\$ 16.33	+Info
445962	Wonder Wedges 15mm Large, Natural.	US\$ 18.15	+Info
446057	Wonder Wedges Assorted, Natural, 10	US\$ 18.15	+Info
446248	Wonder Wedges 17mm Extra-Large, Nat	US\$ 18.15	+Info
446252	Wonder Wedges 13mm MD Natural	US\$ 18.15	+Info

Wedges - DIGITECH

Wedges - DIGITECH

FenderWedge

FenderWedge

Code	Name	Price	+
474822	FenderWedge S-1.9mm-LT GREEN 50pk	US\$ 14.93	+Info
474823	FenderWedge M-1.9mm-LT ORANGE 50pk	US\$ 14.16	+Info
474824	FenderWedge L-2.8mm-LT PURPLE 50pk	US\$ 15.23	+Info
474825	Fdwgd Sm Md Lg - Or 150pcs	US\$ 38.76	+Info

Silicone Wedge

Silicone Wedge

Code	Name	Price	+
474826	Silicone Wedge 20pcs	US\$ 9.19	+Info

WEDGES - MICROBRUSH

WEDGES - MICROBRUSH

Features:

The ContactPro? Wedges are anatomically shaped to provide optimum adaptation along the tooth cervical margin. Made of Maple Wood from sustainably managed forests, these wedges come in four sizes and colors to work as a system with the ContactPro? Matrix Ring.



Code	Name	Price	+
401555	MICROBRUSH ContactPro Wedge- Medium- GREEN- 300pk #291748 -MIC 291748	US\$ 40.25	+Info
401554	MICROBRUSH ContactPro Wedge- Large- YELLOW- 300pk #291749 -MIC 291749	US\$ 40.25	+Info

Wedges - PREMIER

Wedges - PREMIER

Code	Name	Price	+
------	------	-------	---

400810	Sycamore Wedges Asstd 400pk	US\$ 50.49	+Info
400811	Sycamore Wdgs UT Orange 400pk	US\$ 50.49	+Info

Wedges - PULPDENT

Wedges - PULPDENT

Code	Name	Price	+
039-FIN	Wedges Fintrec Transp 250pk	US\$ 65.25	+Info

Wedges - TOR

Wedges - TOR

Code	Name	Price	+
400287	Wooden Wedges 6 types 400pk	US\$ 13.75	+Info

WIZARD WEDGES - WATERPIK

WIZARD WEDGES - WATERPIK

Features

Waterpik Brand Dental Items Supplies, Reorder Code 061215-000

These proven wedges ensure better access and stable holding of matrix systems during your critical restorative procedures. Specially crafted from birch wood for greater tooth separation and faster, better tooth position recovery upon removal. Available in four sizes.

- * Birchwood
- * Triangular Shape
- * Assorted



Code	Name	Price	+
400854	WaterPik - Wizard Wedges ASSORTED 500pk #061215	US\$ 51.10	+Info
400857	WaterPik - Wizard Wedges Small 500pk - #061211	US\$ 53.29	+Info
400855	WaterPik - Wizard Wedges Large 500pk - #061212	US\$ 47.25	+Info
400858	WaterPik - Wizard Wedges Wide 500pk - WaterPik - #061214	US\$ 53.29	+Info
400856	WaterPik - Wizard Wedges Slim Jim 500pk - #061213	US\$ 51.10	+Info

Catalog Index / Indice

Accessories	1
Applicator Needles / Dispensing Tips	1
Bonding & MicroBrushes Applicators	4
Cheek Retractor	10
Composite Gun Dispensers	12
Composite Heaters	13
Lip Expanders	14
Mixing Wells	16
Mouth Props	17
Preformed Angle Formers	18
Veneers Carriers	20
Wedges	22