

Electro Surgery Units

Electrosurgery Units



ELECTRO- SURGERY UNIT- BONART

ELECTRO- SURGERY UNIT- BONART

The Bonart Electro Surgery Unit comes with an easy to user User friendly control panel. Safely designed electronics. 10 Level power setting. 7 piece Various standard electrodes.

Technical Specifications:

- * Maximum power output: 50W (400 ohm load)
- * Power supply 110v: 110V±5% 50/60Hz 92VA
- * Power supply 220v: 220V±5% 50/60Hz 92VA
- * Work frequency: 1.5~1.7MHz±5%

Includes:

- * Cutting unit comes with 7 piece variety electrodes set
- * Foot Switch
- * Indifference Plate
- * 1 Handpiece
- * Electrode set (item number #)
- * AC Power Cord
- * 2yr warranty on the PC Board
- * 1 yr on the accessories
- * Users Manual
- * Cut - Cut Coagulation
- * Including Large Loop- Small Loop
- * Heavy Wire
- * Fine Wire
- * Heavy Ball- Diamond & Oval Shaped

REPLACEMENT SET TIPS : # 112212

Weight & Dimentions:

- * Dimenion: 23(L) x 22.5 (W) x 8.5(H) cm
- * Weight: 4 Kg
- * Handpiece Cable: 190cm
- * ind. Plate Cable: 190cm
- * FDA approved

Available: 110 Volts- 220 Volts



ELECTRON

ELECTRON

Code	Name	Price	+
101228	Electro-Surge Unit ART-E1 110v	US\$ 846.45	+Info

[Other Models](#)

Other Models

Code	Name	Price	+
118773	Electro-Surge Unit ART-E1 220v	US\$ 777.60	+Info

Catalog Index / Indice

Electrosurgery Units 1