

Etching Materials

Phosphoric Acid



ETCHANT GEL 37% - PRO-LINE

[ETCHANT GEL 37% - PRO-LINE](#)

Features:

- * Gel Prime-Dent Etchant Gel 37% phosphoric is a medium viscosity etchant gel used for etching enamel and dentin prior to bonding procedures.
- * Prime-Dent Etchant Gel is thixotropic, stays where placed and will not run.
- * Does not dry out and has a long shelf life when properly stored.
- * Washes off clean.



Code	Name	Price	+
008-0360	PRO-LINE - Etching Gel Phosphoric Acid 37% Syringe 12gr + application tips (10...	US\$ 6.60	+Info
402369	Pro-line Etching Gel Phosphoric Acid 37% Syringe 60grs (5x12grs) + 25 applica...	US\$ 18.15	+Info

ETCH-RITE ETCHING GEL - PULPDENT

[ETCH-RITE ETCHING GEL - PULPDENT](#)

Features:

38% Phosphoric Acid

Etch-Rite is a soft, thixotropic gel with handling characteristics most preferred by clinicians. It dispenses through small gauge needles, stays where placed, washes off with ease, and provides the optimal etch pattern on dentin and enamel surfaces to ensure mechanical retention of bonding agents, restorative resins and resin cements. The blue gel is easily visible and the material is available in a wide variety of packaging options.



The world standard for more than two decades and millions of applications

- * Washes off easily without leaving any residue
- * Stays where placed. Does not run onto exposed tissue
- * 15-second etch
- * Dispenses through small gauge tips for precise placement
- * Packaged in syringes for convenient application

There are many packaging options, including 1.2, 3, 6, 9.4, 12 and 25 mL syringes. Economical bulk packs are also available.

Code	Name	Price	+
039-ET-24	PULPDENT Etch-Rite Etching Gel 24/Pk. 38% Phosphoric Acid, Bulk Kit: 24pk Syring...	US\$ 92.57	+Info
039-ET-50	PULPDENT Etch-Rite Jumbo Kit - 38% Phosphoric Acid Etching Gel: 2 - 25 ml Syring...	US\$ 67.00	+Info
039-ET-50R	PULPDENT Etch-Rite Jumbo Syringe Refill - 38% Phosphoric Acid Etching Gel: 2pk x...	US\$ 55.64	+Info


039-ETCH	PULPDENT Etch-Rite etching gel - 38% Phosphoric Acid, Kit: 4 - 1.2 ml Syringes a...	US\$ 20.62	+Info
039-ET-TWIN	PULPDENT Etch-Rite etching gel - 38% Phosphoric Acid, Twin Pack: 2 - 3 ml Syring...	US\$ 27.42	+Info
039-ET-12	PULPDENT Etch-Rite etching gel - 38% Phosphoric Acid, 12 ml Syringe Only - #ET-12	US\$ 19.49	+Info
039-ET-6	PULPDENT Etch-Rite etching gel - 38% Phosphoric Acid, Single 6 ml Syringe - #ET-6	US\$ 14.17	+Info
039-SY-ET3	PULPDENT Etch-Rite Dental Etching Gel - 38% Phosphoric Acid 3ml empty Refill Syr...	US\$ 11.31	+Info

ETCHING MATERIALS 37% - IVOCLAR VIVADENT



Code	Name	Price	+
442579	Total Etch 37% Phosphoric Acid Etch	US\$ 31.37	+Info

ETCHING MATERIALS GEL 37% - MARK3



Code	Name	Price	+
425522	Etchant Gel 37% Phosphoric 1pk 50ml	US\$ 22.88	+Info
425537	Etchant Gel 37% Bulk 25pk	US\$ 44.48	+Info
425704	Etchant Gel 37% Posph 4pk 1.2ml	US\$ 9.38	+Info

PHOSPHORIC ACID (85%) REFILL - ADS

PHOSPHORIC ACID (85%) REFILL - ADS

High purity, Premium grade. Generally used to etch enamel.



Code	Name	Price	+
118276	ADS Phosphoric Acid 0.85 4oz. - #P883-7	US\$ 16.97	+Info
118277	ADS Phosphoric Acid 0.85 pt 16oz. - #P883-5	US\$ 52.23	+Info

Phosphoric Acid - 3M ESPE

Phosphoric Acid - 3M ESPE

Code	Name	Price	+
427640	Scotchbond Etching Gel Refill - 37%	US\$ 79.43	+Info
427718	Scotchbond Export Package - Blue et	US\$ 110.46	+Info
427785	Scotchbond Universal Etchant, 2 - 3	US\$ 58.50	+Info
427913	Ketac Conditioner Liquid - 25% Phos	US\$ 22.04	+Info
428041	Scotchbond Etching Gel Dispensing T	US\$ 19.49	+Info
428172	Concise Etching Liquid 37% Phosphor	US\$ 76.93	+Info

Phosphoric Acid - COLTENE WHALADENT

Phosphoric Acid - COLTENE WHALADENT

Code	Name	Price	+
426559	Etchant Gel S, 35% (blue) Phosphori	US\$ 40.24	+Info
432052	One Coat self-etch bond Introducator	US\$ 265.06	+Info

Phosphoric Acid - DENMAT

Phosphoric Acid - DENMAT

Code	Name	Price	+
428513	Tenure UniBond Complete Kit with GI	US\$ 0.00	+Info

Phosphoric Acid - DENTSPLY

Phosphoric Acid - DENTSPLY

Code	Name	Price	+
469565	Caulk Tooth Conditioner Gel- 34% Ph	US\$ 69.86	+Info

Phosphoric Acid - IVOCLAR VIVADENT

Phosphoric Acid - IVOCLAR VIVADENT

Code	Name	Price	+
442577	Enamel Prep Blue Etching Gel, 37% P	US\$ 0.00	+Info

Phosphoric Acid - KEYSTONE

Phosphoric Acid - KEYSTONE

Code	Name	Price	+
445989	Acid-Etch 35% Blue Etching Gel, Sta	US\$ 50.85	+Info
446212	Keystone Etching Gel for Metals, 15	US\$ 26.18	+Info

Phosphoric Acid - PREMIER

Phosphoric Acid - PREMIER



Code	Name	Price	
45327	Spot On Acid Dentin/Enamel Etch Gel	US\$ 38.91	+Info
453892	Premier Etch 37% Phosphoric Acid Ge	US\$ 41.26	+Info
454129	Spot On Acid Dentin/Enamel Etch Gel	US\$ 35.16	+Info

Phosphoric Acid - PULPDENT

Phosphoric Acid - PULPDENT

Code	Name	Price	
039-ET-144	Etch-Rite Clinic Kit 1.2ml 144pk	US\$ 0.00	+Info

Phosphoric Acid - ULTRADENT

Phosphoric Acid - ULTRADENT

Code	Name	Price	
463886	Ultra-Etch IndiSpense Syringe Refil	US\$ 107.87	+Info
463903	Ultra-Etch Kit: 4 x 1.2 ml Syringes	US\$ 51.43	+Info

Etching Materials

Hydrofluoric Acid for Porcelain



HydroFluoric Acid for Porcelain - ULTRADENT

Fluorhydric Acid for Porcelain - ULTRADENT



Code	Name	Price	
463901	Porce Etch Acid 9% 2pk	US\$ 58.10	+Info
463909	Porcelain Etch Kit Buffered 9.6%	US\$ 91.78	+Info

HydroFluoric Super Stripper - ADS

ADS - Super Stripper Refill
Fluorhydric Acid for Porcelain
#118403 - Mfg # A864-52



Features:

- * Powerful stripping strength
- * Shipped only in pints
- * Cleans and devests castings and strips porcelain like no other Hydrofluoric Acid substitute available

Quantity: 4 pints - No neutralizer

Code	Name	Price	+
118402	ADS Super Stripper refills 2 - 16oz. btl - #A 864-60	US\$ 37.83	+Info
118403	ADS Super Stripper refills 4 - 16oz. btl - #A 864-52	US\$ 54.97	+Info
118404	ADS Super Stripper (hydrofluoric acid replacement formula) 8oz Stripper (Acid...	US\$ 42.39	+Info
118405	ADS Super Stripper with neutralizer kit - #S 222	US\$ 25.16	+Info

HydroFluoric Acid for Porcelain - KEYSTONE

Fluorhydric Acid for Porcelain - KEYSTONE



Code	Name	Price	+
123871	Ceramic Etch Gel 15ml	US\$ 26.18	+Info

446176	Strip-it Acid Substitute, 4 Quarts	US\$ 42.23	+Info
--------	------------------------------------	------------	-----------------------

HydroFluoric Acid for Porcelain - PRO-LINE

Fluorhydic Acid for Porcelain - PRO-LINE



Code	Name	Price	+
008-050	PRO-LINE - Porc Etch Gel Syr 3gr	US\$ 15.20	+Info

HydroFluoric Acid for Porcelain - PULPDENT

Fluorhydic Acid for Porcelain - PULPDENT



Code	Name	Price	+
039-PEG	Porcelain Etch Gel Kit 4x1.2gr syr	US\$ 57.56	+Info
101498	Porc Etch Gel 4x1.2ml Syr Only	US\$ 57.56	+Info
101497	Porc Prep Etch Kit 4x2.1ml syr	US\$ 78.57	+Info

Hydrofluoric Acid for Porcelain - ULTRADENT

Hydrofluoric Acid for Porcelain - ULTRADENT

Code	Name	Price	+
470493	Porcelain Etch Acid 9% 4pk	US\$ 104.26	+Info

PORCELAIN WASH - ADS

ADS - Porcelain Wash - Gallon
Fluorhydric acid for Porcelain
#118305 - Mfg #P996-8

Features:

- * Microscopically cleans and decontaminates all soiled porcelain restorations
- * Chemical scrubbing action dissolves organic debris trapped in the pores of dental porcelain
- * Restorations that have been in the mouth can now be refired with no fear of blow-ups, discoloration or out-gassing
- * High purity and potency
- * Excellent for removal of finishing and grinding impurities that, even if not visible to the naked eye, would contaminate ceramic work pieces
- * Eliminates costly remakes due to exploded, cracked or discolored restorations
- * Will not affect porcelain shades



Quantity: 1 gallon

Code	Name	Price	+
118304	ADS Porcelain Wash qt - #P996-6	US\$ 22.26	+Info
118305	ADS Porcelain Wash gal - #P996-8	US\$ 60.97	+Info

Etching Materials

Hydrochloric Acid for Porcelain

HATHO
POLISHING KIT

Thermoplastics
like *never* before



Packing Diamond D®
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory



HYDROCHLORIC ACID - ADS

ADS -HYDROCHLORIC ACID

Features

- High purity, Premium grade. Generally used as a cleaner and pickle.
- All ADS, Inc. Acid Kits are concentrated acids in premeasured bottles.
- All acids are high purity, premium grade.
- Shipped by standard methods. No hazard fees.



CONCENTRATED

Code	Name	Price	+
118087	ADS Hydrochloric Acid concentrate refill 6 - 1 oz. btls - #A200-2	US\$ 27.77	+Info
118090	ADS Hydrofluoric Acid concentrate refill 6 - 1 oz. btls - #A210-2 49%	US\$ 27.77	+Info

KITS

Code	Name	Price	+
118086	ADS Hydrochloric Acid concentrate kit - #A200-1	US\$ 47.50	+Info
118089	ADS Hydrofluoric Acid concentrate kit - #A210-1 6 x 1oz acide bottles - glasse...	US\$ 48.81	+Info

PLUS

Code	Name	Price	+
118088	ADS Hydrochloric Acid PLUS qt - #H886-6	US\$ 32.16	+Info

Etching Materials

Acid Neutralizer



Acid Neutralizer - ECCO

Acid Neutralizer - ECCO

Code	Name	Price	+
567320	Acid Neutralizer 4x1.2cc Syrin	US\$ 14.40	+Info

Etching Materials

Eugenol



Eugenol - KEYSTONE

Eugenol - KEYSTONE

Code	Name	Price	+
446910	Keystone Eugenol, USP Grade, 1 oz.	US\$ 8.74	+Info

EUGENOL - ADS

EUGENOL - ADS

Used in various restorative, endodontic and prosthodontics procedures as an analgesic and antiseptic.



Code	Name	Price	+
117954	ADS Eugenol 4oz. - #S883-5-4	US\$ 21.62	+Info
117955	ADS Eugenol 8oz. - #S883-5-5	US\$ 29.23	+Info

EUGENOL USP - PULPDENT

EUGENOL USP - PULPDENT

Features:

Used in restorative and prosthodontic applications as a antiseptic and analgesic.
Highly purified eugenol for a multitude of dental applications.
Available in 2, 4 and 16 oz bottles.



Code	Name	Price	+
039-E-02	PULPDENT Pulpdent Eugenol, 2 oz. Bottle - #E-02	US\$ 37.31	+Info
039-E-04	PULPDENT Pulpdent Eugenol, 4 oz. Bottle - #E-04	US\$ 38.70	+Info
039-E-16	PULPDENT Eugenol USP 16Oz (472mL) #E-16	US\$ 0.00	+Info

Catalog Index / Indice

Phosphoric Acid -----	1
Hydrofluoric Acid for Porcelain -----	4
Hydrochloric Acid for Porcelain -----	7
Acid Neutralizer -----	8
Eugenol -----	8