

Evacuation & Suction Systems

Portable Suction Evacuation systems



Portable Suction Evacuation systems - ASEPTIC



**ASEPTICO - Vacuum Saliva Ejector - Portable Units
Vacuum Systems
115333 - Mfg # AA-20**

According to the Royal College of Dental Surgeons, the Portable Dental Saliva Ejector is one of the most-used auxiliary powered backups for suction. Powered by a rechargeable battery pack, wall outlet, or 12V DC vehicle source, it can be used as a battery backup for sedation dentistry. This portable dental suction unit features adjustable 80-550 mmHg vacuum range and an 800cc waste canister with float shut-off, and it packs into its included carrying case when it is no longer needed.

Features

- * Powered via rechargeable battery pack, wall outlet or 12 VDC vehicle source
- * Recharge battery via 110 to 220 V wall outlet or 12 VDC car lighter power adapter (included)
- * Adjustable 80 - 550 mm Hg vacuum range
- * Includes padded carrying case
- * 800 cc waste canister with float shut-off

Dimensions: 8.3" x 8" x 8.5"

Weight: 6.6 lbs (3 kg)

Content per box:

- * Suction Unit and Tubing (AA-20)
- * Soft Case



Code	Name	Price	+
115333	Vacuum Saliva Ejector Portable	US\$ 742.22	+Info
417370	Aseptico High Volume Suction	US\$ 669.41	+Info
417371	Aseptico Low Volume Suction	US\$ 665.98	+Info

**Portable Suction Evacuation system
s - PEGASUS**

Portable Suction Evacuation systems - PEGASUS



Code	Name	Price	+
118715	Vacuum Saliva Ejector PORT-S2	US\$ 302.90	+Info

Portable Suction Evacuation systems - TPC

Portable Suction Evacuation systems - TPC



Code	Name	Price	+
117368	Mobile-Vac Portable Evacuation Unit	US\$ 1,187.50	+Info

Vacuum Systems - TECHWEST

Vacuum Systems - TECHWEST



ECO-VAC

ECO-VAC

1-4 Chairs

1-4 Chairs

Code	Name	Price	+
VPD2S2	ECO-VAC 2 Chairs 1 x1 HP 220v	US\$ 8,683.99	+Info
VPD4D2	ECO-VAC 4 Chairs 2 x1 HP 220v	US\$ 15,599.48	+Info

5-12 Chairs

5-12 Chairs

Code	Name	Price	+
VPD10D2	ECO-VAC 12 Chairs 2x2 HP 220v	US\$ 18,321.66	+Info
VPD5S2	ECO-VAC 7 Chairs 1x2 HP 220v	US\$ 9,676.97	+Info

GOLDENVAC

GOLDENVAC

No Rec/No Ssep

No Rec/No Ssep

Code	Name	Price	+
VPL3SS	GV 3 Cha Pump 1x1.5HP 110/220v	US\$ 2,477.48	+Info
VPL4S2	GV 4 Chairs Pump 1x2 HP 220v	US\$ 2,953.86	+Info
VPL6D2	GV 6 Chairs Pump 2x1.5 HP 220v	US\$ 6,448.48	+Info

No Rec/Yes Sep

No Rec/Yes Sep

Code	Name	Price	+
VPL12T2	GV 12 Chairs Pump 3x2 HP 220v	US\$ 9,806.29	+Info
VPL4D2	GV 4 Chairs Pump 2x1 HP 220v	US\$ 4,927.65	+Info
VPL8D2	GV 8 Chairs Pump 2x2 HP 220v	US\$ 6,566.40	+Info
VPL9T2	GV 9 Chairs Pump 3x1.5 HP 220v	US\$ 9,075.10	+Info

Yes Rec/No Sep

Yes Rec/No Sep

Code	Name	Price	+
474335	GV 4 Chairs Pump 2x1 HP 220v	US\$ 3,315.68	+Info
VPL3SSR	GV 3 Cha Pump 1x1.5HP 110/220v	US\$ 2,517.21	+Info

Vacuum Systems - SOLMETEX

Vacuum Systems - SOLMETEX



Code	Name	Price	+
475015	Solmetex DryVac Tankless Vacuum	US\$ 6,443.55	+Info

Accessories for DryVac

Code	Name	Price	+
475016	SOLMETEX DryVac Tankless MAINTENANCE KIT # DRYVC-NXT-MKT	US\$ 213.00	+Info

Catalog Index / Indice

Portable Suction Evacuation systems	1
Vacuum Systems	2