

## Military Dental Products



### Dental units

#### Dental units - ASEPTICO

Dental units - ASEPTICO

Code	Name	Price	+
111241	Dental Cart Duovolt	US\$ 26,437.50	<a href="#">+Info</a>

#### PORTABLE FIELD SINK - ASEPTICO

PORTABLE FIELD SINK - ASEPTICO

- \* Hands-free washing using included foot controls
- \* Integrated electric water heater
- \* Stainless steel side tray mounts to either side
- \* Use standard military water containers for water supply and waste disposal
- \* Water and waste lines designed to protect against cross contamination
- \* Designed for easy maintenance
- \* Dual-voltage, auto-switching: 120/220V - 50/60 Hz



#### PORTABLE FIELD SINK

Code	Name	Price	+
111234	ASEPTICO Command Air Portable Field Sink 6520-01-456-5380	US\$ 6,656.25	<a href="#">+Info</a>
111235	ASEPTICO Portable Field Scrub Sink Large	US\$ 8,484.38	<a href="#">+Info</a>

#### ACCESSORIES

Code	Name	Price	+
111229	ASEPTICO Extended Spout For Adu-40Cf Assy #AA-04EXS	US\$ 451.50	<a href="#">+Info</a>
111230	ASEPTICO Soap Dispenser For Adu-40Cf #AA-04SD	US\$ 514.50	<a href="#">+Info</a>

## Catalog Index / Indice

Dental units ..... 1