

Pins and Posts (Dentine)

ParaPulpal Screw Pin / Posts ALL BRANDS



PARAPULPAL DENTINE SCREW POSTS COLTENE



ACCESSORIES

ACCESSORIES

Code	Name	Price	+
426940	TMS B-59-3 Pin hand wrench used wit	US\$ 23.69	+Info
427062	TMS B-59-5 Pin hand wrench used wit	US\$ 29.84	+Info

BONDENT Titanium Gold Plated Pins

Screw Post Gold Plated - COLTENE WHALADENT



BONDENT Drills

Code	Name	Price	+
431589	COLTENE Bondent Drills (Depth Limiting), package of 6 drills - # DB-91	US\$ 131.46	+Info

Other Models

Other Models

Code	Name	Price	+
430792	Bondent .425 mm x 2.5 mm titanium a	US\$ 197.68	+Info
431596	Bondent .425 mm x 2.5 mm titanium a	US\$ 568.97	+Info
436184	Bondent .425 mm x 2.0 mm titanium a	US\$ 197.68	+Info

Drills Kodex K-97 ParaPulpal (Common for all Parapulpal Pins Coltene)

A parapulpal pin in dentistry refers to a small metal pin that is placed in the dentin, near the pulp, to enhance the retention and support of a dental restoration. These pins can be used in situations where tooth structure is compromised and there's not enough retention for restorations like fillings or crowns.

Here's a brief overview of how they're used:

- 1.Preparation: After removing decay and preparing the tooth for the restoration, the dentist drills a small hole in the dentin where the pin will be placed.
- 2.Insertion: The pin is then threaded or tapped into the hole. Once in place, it provides additional support for the dental material that will be used to restore the tooth.
- 3.Restoration: Dental material (like amalgam, composite, etc.) is then placed over and around the pin, effectively locking it in place and using it as



an anchor.

It's essential for the pin to be placed in a manner that avoids the pulp (the nerve and blood vessels in the center of the tooth) to prevent any possible injury. Proper placement ensures that the restoration is both functional and comfortable for the patient.

However, with advancements in dental materials and techniques, the use of parapulpal pins has become less common. Modern adhesive dentistry techniques often allow for sufficient retention of restorations without the need for such pins. Still, they may be used in certain situations where added support is required.

Code	Name	Price	+
426720	Kodex K-97 Gold .027" x .197" / .67	US\$ 89.04	+Info
426729	Kodex K-93 Silver .021"x.157" 6pk	US\$ 78.92	+Info
426762	Kodex K-91 red .017" X .059" / .425	US\$ 128.63	+Info
426824	Kodex K-96 Gold .027" x .079" / .67	US\$ 140.32	+Info
426971	Kodex K-89 Pink .0135" x .051" / .3	US\$ 147.03	+Info
433445	Kodex K-90 red .017" x .059" / .425	US\$ 147.90	+Info
436463	Kodex K-98 Green .028" x .197" / .7	US\$ 149.43	+Info
426493	Kodex K-92 Silver .021" x .079" / .	US\$ 103.95	+Info

MAX Pin Titanium Pins

Titanium Post - COLTENE



Drills Max Pin

Code	Name	Price	+
426632	COLTENE Max M-92 PURPLE Drills .021in x .095in / .525 mm x 2.4 mm (depth limit...	US\$ 217.46	+Info
426865	COLTENE Max M-91 BLUE Drills .017in x .075in / .425 mm x 1.9 mm (depth limitin...	US\$ 227.08	+Info

Other Models

Other Models

Code	Name	Price	+
426586	Max M-23 Bulk Pack .021" / .525 mm	US\$ 124.79	+Info
426611	Max M-21 Complete Kit .021" / .525	US\$ 154.63	+Info

426668	Max M-13 Bulk Pack .017" / .425 mm	US\$ 280.05	+Info
426745	Max M-11 Complete Kit .017" / .425	US\$ 269.35	+Info
436179	Bondent .425 mm x 2.0 mm titanium a	US\$ 579.94	+Info

TMS DENTINE RETENTION PINS

TMS pins are available as alpha-beta titanium alloy or gold plated stainless steel and four different sizes. The system is completed with the colour coded KODEX-drills to ensure an optimal preparation of the canal. Sizes:

Minuta Pink: 0.350 mm / .0135"
 Minikin Red: 0.425 mm / .017"
 Minim Silver: 0.525 mm / .021"
 Regular Gold: 0.675 mm / .027"



TMS (original)



Minikin Pins

Code	Name	Price	+
430776	TMS Self Shearing TMS T-109 Minikin	US\$ 487.71	+Info
431578	TMS Self Shearing TMS T-108 Minikin	US\$ 292.30	+Info
431580	TMS Self Shearing TMS T-107 Minikin	US\$ 147.13	+Info

Minim Pins

Code	Name	Price	+
426506	TMS T-04 Minim .021" / .525 mm Gold	US\$ 91.64	+Info
426510	TMS T-05 Minim .021" / .525 mm Gold	US\$ 152.48	+Info
426675	TMS T-10 Minim .021" / .525 mm Gold	US\$ 133.73	+Info
426717	TMS T-12 Minim.021" / .525 mm Gold-	US\$ 522.62	+Info
430779	TMS T-11 Minim .021" / .525 mm Gold	US\$ 216.28	+Info
430799	TMS T-06 Minim .021" / .525 mm Gold	US\$ 274.46	+Info

Minuta Pins

Code	Name	Price	+
433714	Coltene/Whaledent TMS Minuta Hand W	US\$ 15.15	+Info
432149	TMS Self Shearing TMS T-111 Minuta	US\$ 146.85	+Info
436183	TMS Self Shearing TMS T-112 Minuta	US\$ 292.30	+Info

Tms Regular Pins

Code	Name	Price	+
426590	TMS T-02 Regular .027" / .675 mm Go	US\$ 152.48	+Info
426672	TMS T-01 Regular .027" / .675 mm Go	US\$ 99.70	+Info

426834	TMS T-07 Regular .027" / .675 mm Go	US\$ 145.70	+Info
431562	TMS T-08 Regular .027" / .675 mm Go	US\$ 216.28	+Info
431607	TMS T-03 Regular .027" / .675 mm Go	US\$ 274.46	+Info
431657	TMS Self Shearing TMS T-13 Regular	US\$ 147.11	+Info
432164	TMS Self Shearing TMS T-14 Regular	US\$ 292.30	+Info
436181	TMS Self Shearing TMS T-15 Regular,	US\$ 407.78	+Info

TMS Bending Tool



Code	Name	Price	+
432118	TMS B-61-M MINIM Bending Tool, sing	US\$ 40.24	+Info
432199	TMS B-61-R REGULAR Bending Tool, si	US\$ 40.34	+Info

TMS Link Plus



Minim Titanium Pins

Bulk

Code	Name	Price	+
426503	COLTENE TMS Link Plus L-722 Minim Silver, Bulk .021in / .525 mm Gold-Plated sta...	US\$ 164.22	+Info
426589	COLTENE TMS Link Plus L-732 Minim 2/1 Silver, Bulk .021in / .525 mm Gold-Plated...	US\$ 227.85	+Info
426609	COLTENE TMS Link Plus L-842 Minim 2/1 Silver, Bulk .021in / .525 mm Titanium Al...	US\$ 227.85	+Info
426661	COLTENE TMS Link Plus L-822 Minim Silver, Bulk .021in / .525 mm Titanium Alloy ...	US\$ 164.22	+Info

Other Models

Other Models

Code	Name	Price	+
426453	TMS Link Plus L-721 Minim Silver .0	US\$ 93.24	+Info
426593	TMS Link Plus L-731 Minim 2/1 Silve	US\$ 124.05	+Info

426665	TMS Link Plus L-821 Minim Silver .0	US\$ 93.24	+Info
426676	TMS Link Plus L-841 Minim 2/1 Silve	US\$ 124.05	+Info
436182	TMS Link Plus Minim .021" / .525 mm	US\$ 332.15	+Info

Regular SS Gold Plated Pins

Code	Name	Price	+
426520	TMS Link Plus L-751 Regular 2/1 Gol	US\$ 98.87	+Info
426634	TMS Link Plus L-741 Regular Gold .0	US\$ 93.15	+Info
431565	TMS Link Plus L-752 Regular 2/1 Gol	US\$ 227.85	+Info
431591	TMS Link Plus L-742 Regular Gold- B	US\$ 164.22	+Info

TMS Link Series



Minikin Pins

Code	Name	Price	+
426504	TMS Link L-512 Minikin red, Bulk .0	US\$ 199.17	+Info
426598	TMS Link L-511 Minikin red .017" /	US\$ 87.12	+Info
426646	TMS Link L-811 Minikin red .017" /	US\$ 90.17	+Info
426647	TMS Link L-812 Minikin red, Bulk .0	US\$ 215.69	+Info

Minim Titanium Pins

Code	Name	Price	+
426499	TMS Link L-532 Minim 2/1 Silver, Bu	US\$ 226.18	+Info
426517	TMS Link L-522 Minim Silver, Bulk .	US\$ 172.59	+Info
426629	TMS Link L-531 Minim 2/1 Silver .02	US\$ 126.18	+Info
426879	TMS Link L-521 Minim Silver .021" /	US\$ 86.89	+Info

Other Models

Other Models

Code	Name	Price	+
431608	TMS Link L-501 Minuta Pink .0135" /	US\$ 132.10	+Info
436180	TMS Link L-502 Minuta Pink, Bulk .0	US\$ 280.60	+Info

Regular SS Gold Plated Pins

Code	Name	Price	+
426645	TMS Link L-551 Regular 2/1 Gold .02	US\$ 123.70	+Info
431586	TMS Link L-541 Regular Gold .027" /	US\$ 86.79	+Info
432227	TMS Link L-542 Regular Gold, Bulk .	US\$ 172.59	+Info
436177	TMS Link L-552 Regular 2/1 Gold, Bu	US\$ 226.18	+Info

PARAPULPAL PLATED PINS ECCO



SCREW POST GOLD PLATED SCREW POST GOLD PLATED - ECCO

Golden cylindrical gold plated screw posts with head for Cross Key.
Golden Composite cylindrical gold plated screw posts with special composite head slot for increased mechanical retention.



Code	Name	Price	+
107562	ECCO - Dental Gold Plated Screw Posts Assortment Kit 120pk (10x12 Posts- 2 Keys)...	US\$ 30.94	+Info
107563	ECCO - Dental Gold Plated Screw Posts Assortment Kit 240pk (20x12 Posts- 2 Keys)...	US\$ 35.64	+Info

SCREW POSTS STAINLESS STEEL SCREW POSTS STAINLEES STEEL - ECCO



Code	Name	Price	+
107548	ECCO - Dental Stainless Steel Screw Posts Assortment Kit 120pk (10x12 Posts- 2 K...	US\$ 17.42	+Info
107549	ECCO - Dental Stainless Steel Screw Posts Assortment Kit 240pk (20x12 Posts- 2 K...	US\$ 36.96	+Info

TITANIUM SCREW POSTS Titanium Screw Posts HOUSE BRAND

Code	Name	Price	+
118716	TITANIUM Scw Pst As.120pk	US\$ 35.20	+Info

ParaPulpal Screw Pin / Posts ALL BRANDS - COLTENE

ParaPulpal Screw Pin / Posts ALL BRANDS - COLTENE

Code	Name	Price	+
426711	Bondent complete kit of Titanium Al	US\$ 327.53	+Info
436178	TMS Self Shearing TMS T-113 Minuta	US\$ 0.00	+Info

Screw Post Gold Plated - HAGER

Screw Post Gold Plated - HAGER

Pins and Posts (Dentine)

ParaPulpal Screw Pin / Posts ALL BRANDS

Code	Name	Price	+
434973	Mirapost Classic Gold Plated Screw	US\$ 74.91	+Info
434974	Mirapost Classic Gold Plated Screw	US\$ 74.91	+Info

Catalog Index / Indice

ParaPulpal Screw Pin / Posts ALL BRANDS 1